






















FastGene™ RNA Premium Kit

培養細胞および組織等からのトータルRNAの精製とゲノムDNA除去
 FG-81006 (6preps), FG-81050 (50preps), FG-81250 (250preps)

| ステップ | スタンダードプロトコル | ラージインプットプロトコル |
|----------------------------|--|---|
| サンプルの準備と量の確認 | < 5×10 ⁶ 培養細胞 < 10 mg 組織 | < 1×10 ⁷ 培養細胞 < 20 mg 組織 |
| 細胞の溶解とホモジナイズ |  350 μL バッファー RL ^{*1} 添加後十分にホモジナイズ S |  600 μL バッファー RL ^{*1} 添加後十分にホモジナイズ S |
| ライセートの清澄化 |   FastGene™ RNA filter column にライセート添加 ≥ 10,000 x g (室温: 20 ~ 25°C) 1 min FastGene™ RNA filter column 廃棄後 ろ液を回収 |  600 μL バッファー RL ^{*1} 添加後十分にホモジナイズ S |
| カラム結合条件の調整 |  350 μL 70% エタノール ピペティングで混合 |  600 μL 70% エタノール ピペティングで混合 |
| カラム結合 |   FastGene™ RNA binding column に 最大 700 μL までのサンプル溶液を添加 ≥ 10,000 x g (室温: 20 ~ 25°C) 1 min ろ液廃棄後 カラムを元のコレクションチューブ (2.0mL) に戻す | サンプル溶液が なくなるまで繰り返す |
| メンブレン洗浄 1 (タンパク除去) |  600 μL バッファー RW1 ≥ 10,000 x g (室温: 20 ~ 25°C) 30 s カラムを新しいコレクションチューブ (2.0mL) に移す | |
| メンブレン洗浄 2 (塩類の除去) |  700 μL バッファー RW2 ^{*1} ≥ 10,000 x g (室温: 20 ~ 25°C) 30 s カラムを新しいコレクションチューブ (2.0mL) に移す | |
| メンブレン乾燥 |  フルスPEEDで遠心 (室温: 20 ~ 25°C) 1min カラムを新しいコレクションチューブ (1.5mL) に移す | |
| 溶出 |  50 μL バッファー RE (注: メンブレンの中央に添加) ≥ 10,000 x g (室温: 20 ~ 25°C) 1 min FastGene™ RNA binding column 廃棄後 溶出液を回収 S | |
| DNase I 反応条件の調整 |  5 μL 10x DNase I reaction buffer | |
| DNase I 反応 (DNAの分解) |  1 μL DNase I ^{*1} ピペティングで混合 (ピペットを 50 μL にセット) インキュベート (室温: 20 ~ 25°C) 10 min | |
| カラム結合条件の 再調整 |  250 μL バッファー RBD ^{*1} ピペティングで混合 | |
| カラム結合 |   FastGene™ RNA mini-elute column にサンプル溶液を添加 ≥ 10,000 x g (室温: 20 ~ 25°C) 1 min ろ液廃棄後 カラムを元のコレクションチューブ (2.0mL) に戻す | |
| メンブレン洗浄 3 (塩類と分解DNAの除去) |  700 μL バッファー RW2 ^{*1} ≥ 10,000 x g (室温: 20 ~ 25°C) 30 s カラムを新しいコレクションチューブ (2.0mL) に移す | |
| メンブレン乾燥 |  フルスPEEDで遠心 (室温: 20 ~ 25°C) 1 min カラムを新しいコレクションチューブ (1.5mL) に移す | |
| 溶出 |  適量 ^{*2} の バッファー RE (注: メンブレンの中央に添加) ≥ 10,000 x g (室温: 20 ~ 25°C) 1 min FastGene™ RNA mini-elute column 廃棄後 溶出液を回収 | |
























※1: これらの試薬は事前調整が必要です。 ※2: スタンダードプロトコル 20 μL (10 ~ 50 μL)、ラージインプットプロトコル 50 μL (20 ~ 50 μL)

S Safety Stopping Point. この操作後、-70°C以下の保存も可能です。


FastGene™ RNA Premium Kit

For purification of total RNA (gDNA removed) from animal cells/tissues

FG-81006 (6preps) , FG-81050 (50preps) , FG-81250 (250preps)

| Step | Standard protocol | Large input protocol |
|--|--|---|
| Sample quantity | < 5×10 ⁶ cultured animal cells < 10 mg animal tissues | < 1×10 ⁷ cultured animal cells < 20 mg animal tissues |
| Resuspension/homogenization by cell lysis |  350 μL buffer RL ^{※1} Vortex vigorously  |  600 μL buffer RL ^{※1} Vortex vigorously  |
| Filtration of cellular debris |   Transfer lysate into a FastGene™ RNA filter column Centrifuge at ≥ 10,000 x g (RT : 20 ~ 25°C) 1 min Discard filter column, Harvest flow-through | |
| Optimize RNA binding conditions |  Add 350 μL 70 % ethanol Mix thoroughly by pipetting |  Add 600 μL 70 % ethanol Mix thoroughly by pipetting |
| RNA binding |   Load mix (up to 700 μL) onto FastGene™ RNA binding column Centrifuge at ≥ 10,000 x g (RT : 20 ~ 25°C) 1 min Discard flow-through, Re-insert binding column in collection tube (2.0 mL) | Repeat the procedure for larger volume. |
| Protein elimination |  Add 600 μL of buffer RW1 Centrifuge at ≥ 10,000 x g (RT : 20 ~ 25°C) 30 s Transfer binding column to new collection tube (2.0 mL) | |
| Desalination |  Add 700 μL of buffer RW2 ^{※1} Centrifuge at ≥ 10,000 x g (RT : 20 ~ 25°C) 30 s Transfer binding column to new collection tube (2.0 mL) | |
| Removal of RW2 |  Centrifuge at full speed (RT : 20 ~ 25°C) 1 min Transfer binding column to new collection tube (1.5 mL) | |
| Elution of RNA |  Add 50 μL of buffer RE to membrane center Centrifuge at ≥ 10,000 x g (RT : 20 ~ 25°C) 1 min Discard binding column, Harvest eluted solution  | |
| Optimize DNase I conditions |  Add 5 μL 10x DNase I reaction buffer | |
| DNA digestion |  Add 1 μL of DNase I ^{※1} to the mixture Mix thoroughly by pipetting Incubate for 10 min (at room temp. 20-25°C) | |
| RNA rebinding optimization |  Add 250 μL of buffer RBD ^{※1} to the mixture Mix thoroughly by pipetting | |
| RNA binding |   Transfer mixture into FastGene™ RNA mini-elute column Centrifuge at ≥ 10,000 x g (RT : 20 ~ 25°C) 1 min Discard flow-through, Re-insert binding column in collection tube (2.0 mL) | |
| Desalination/ Elimination of digested DNA |  Add 700 μL buffer RW2 ^{※1} Centrifuge at ≥ 10,000 x g (RT : 20 ~ 25°C) 30 s Transfer mini-elute column in new collection tube (2 mL) | |
| Removal of RW2 |  Centrifuge at full speed (RT : 20 ~ 25°C) 1 min Transfer mini-elute column in new collection tube (1.5 mL) | |
| Elution of RNA |  Add appropriate volume ^{※2} of buffer RE to the membrane center Centrifuge at ≥ 10,000 x g (RT : 20 ~ 25°C) 1 min Discard mini-elute column, Harvest eluted solution | |

※1 : need preparation before use. ※2 : Standard protocol 20 μL (10 ~ 50 μL). Large input protocol 50 μL (20 ~ 50 μL)

 Safety Stopping Point. Sample after homogenize step or 1st elution step can be stored at -70°C.