


KAPA HyperCap Workflow v3.4 Guide to Success



Using KAPA HyperPrep, KAPA HyperPlus or KAPA EvoPlus Kits

Sample prep workflow using superior probe technology with optimized Library Preparation and Target Enrichment.



Get to chopping
Fragmentation of high quality gDNA has been optimized to 180 – 220 bp average fragment length using 22 min at 37°C. Look for further optimization instructions in the KAPA HyperPlus Kit TDS.

It's not a typo!
Ensure that you are adding the correct volume of KAPA Frag Buffer (5 µL) and KAPA Frag Enzyme (10 µL) to each reaction.

Component	Volume
Double-stranded DNA (with KAPA Frag Conditioning Solution, if needed)	35 µL
KAPA Frag Buffer (10X)	5 µL
KAPA Frag Enzyme	10 µL
Total volume	50 µL

Got EDTA?

- If using the KAPA HyperPlus Kits (proceed to step 2B), the enzymatic fragmentation reaction is sensitive to EDTA.
- If your DNA contains EDTA, please refer to the KAPA HyperPlus Kit Technical Data Sheet (TDS) on how to easily remove it with KAPA HyperPure Beads, or mitigate it (in case you know its concentration) by using the Conditioning Solution provided in the kit.
- If using the KAPA EvoPlus Kits, no need to worry – the enzymatic fragmentation reaction is **insensitive to EDTA (up to 2 mM)**, as well as numerous other inhibitors.

The KAPA FragTail ReadyMix is viscous and requires special attention during pipetting.
Fragmentation at 37°C for 18 minutes followed by A-tailing at 55°C for 30 minutes.

Which End Repair and A-Tailing Enzyme Mix to use?

The KAPA HyperCap Workflow v3.4 has been validated with the updated HyperPlus ERAT Enzyme Mix (orange cap) for improved performance.

Do not use the HyperPlus ERAT Enzyme Mix (orange cap) if you are using Covaris sonicated or mechanically sheared gDNA.

Incubation Conditions
Remember the additional first incubation step when using KAPA HyperPrep Kits. 30 min at 20°C first followed by 30 mins at 65°C.

How many cycles of amplification in total?

KAPA HyperPrep Kit	8 cycles
KAPA HyperPlus Kit	6 cycles
KAPA EvoPlus Kit	7 cycles

So, did we make a good library?

Determine the concentration (using Qubit), the size distribution and quality (using the Agilent Bioanalyzer) of the amplified library. Is the library good enough to proceed to target capture?

We can help...
That's what we're here for

Visit us at sequencing.roche.com/support

HyperDesign Tool

Intuitive online custom design tool with Roche's renowned probe-design algorithm in a simple, easy-to-use interactive web interface.

www.hyperdesign.com

8 TARGET CAPTURE HYBRIDIZATION-BASED

How much Hybridization Component H to add?

Hybridize your KAPA HyperCap Probes to your single/multiplex DNA sample library, where the amount of Hybridization Component H depends on the capture target size.

<40 Mbp Capture Target Size	12 µL
≥40 Mbp Capture Target Size	9 µL

9 WASH & RECOVERY

55°C for 5min

Prepare your wash buffers to create working solutions before you start the day.

Prepare your Capture Beads based on the capture target size.

How many total cycles of amplification are required?

Capture target*	Total cycles
<100 kb	18
100 kb – 2 Mb	16
>2 Mb – 40 Mb	10
>40 Mb	8

Capture Target Size

<40 Mbp	50 µL per rxn
≥40 Mbp	100 µL per rxn

Great! Need more KAPA Sample Prep Reagents?

Contact your local sales representative to order more.

* Check the "Coverage summary" file (HyperDesign) for the "Total capture space" value.

Get the latest Operators Manuals (Method Sheets, Users Guides and Instructions for Use) at navifyportal.roche.com



