Application

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KAPA Library Preparation Kit
KAPA Real-Time Library Amplification Kit
KAPA Library Quantification Kit

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[LIMprep Workflow]

<Start material>	10ng fragmented mouse genome DNA (n=9)
<KAPA library preparation kit>	End-repair, dA-tail, adaptor ligation
<KAPA real-time library amplification kit>	Determination of PCR cycle
<KAPA library amplification kit>	PCR enrichment
<Resulted Library DNA>	Quantified by pico-green Average: 73.8±4.8 ng, CV: 0.066
<Quantified by qPCR>	1ng Library DNA / gene locus (n=9) Detected 8 gene locus (A-H)
<KAPA library quantification kit>	Final QC
<MiSeq Run>	SE50bp


* LIMprep: Ligation-based Illumina Multiplex library PREPreparation method for low amount DNA using TruSeq adaptor and KAPA library preparation kit

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[LIMprep Workflow]

<Start material>	10ng fragmented mouse genome DNA (n=9)
<KAPA library preparation kit>	End-repair, dA-tail, adaptor ligation
<KAPA real-time library amplification kit>	Determination of PCR cycle ①
<KAPA library amplification kit>	PCR enrichment
<Resulted Library DNA>	Quantified by pico-green Average: 73.8±4.8 ng, CV: 0.066 ②, ③
<Quantified by qPCR>	1ng Library DNA / gene locus (n=9) Detected 8 gene locus (A-H)
<KAPA library quantification kit>	Final QC
<MiSeq Run>	SE50bp

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

<KAPA library amplification kit>

PCR enrichment

Quantified by pico-green Average: 73.8±4.8 ng, CV: 0.066

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Michael A Quail et al., Optimal enzymes for amplifying sequencing libraries, Nature Methods 9, 10-11 (2012)

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[LIMprep Workflow]

<Start material>	10ng fragmented mouse genome DNA (n=9)	
<KAPA library preparation kit>	End-repair, dA-tail, adaptor ligation	
<KAPA real-time library amplification kit>	Determination of PCR cycle	
<KAPA library amplification kit>	PCR enrichment	
<Resulted Library DNA>	Quantified by pico-green Average: 73.8±4.8 ng, CV: 0.066	
<Quantified by qPCR>	1ng Library DNA / gene locus (n=9) Detected 8 gene locus (A-H)	4
<KAPA library quantification kit>	Final QC	
<MiSeq Run>	SE50bp	5

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