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New Design

DRAFT

① Start design — ② Configure design — ③ Add targets — ④ Review targets — ⑤ Submit design — ⑥ Done

- ① Start design
- ② Configure design
- ③ Add targets
- ④ Review targets
- ⑤ Submit design
- ⑥ Done

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HyperDesign Roche Sequencing Solutions Education and Support Log in Sign up Roche

## Welcome to the new HyperDesign custom design tool.

1.

New probe technology combined with our renowned expertise in probe design

More intuitive user interface from start to finish

Simple and streamlined custom design process

Already have a NimbleDesign account? All your data is here! Simply log in using your NimbleDesign credentials. [Log in](#)

### Sign up now

First name \* Last name \*

Country \*

Email \*

Choose a Password \*

I agree to the [Terms of Service](#) and [Privacy Policy](#).

2. [Sign up](#)

[Learn more about the tool](#) [KAPA library prep products](#)

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Online Technical Support

## Your recent designs

Custom Design Name	Product	Created	Status
hello	KAPA HyperChoice	28 Jan 2020	Draft
HyperDesign	KAPA HyperChoice	28 Jan 2020	Ready for review
test123	KAPA HyperChoice	28 Jan 2020	Draft
forMarkos	KAPA HyperChoice	27 Jan 2020	Approved
sadsad	KAPA HyperChoice	27 Jan 2020	Ready for review

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## Messages

Clear all messages

Design HyperDesign has a completed probe selection. The job is ready for you to review.

Design forMarkos has a completed probe selection. The job is ready for you to review.

HyperDesign, a new design tool is coming soon, to complement our new target enrichment technology!

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New Design

DRAFT

1 Start design — 2 Configure design — 3 Add targets — 4 Review targets — 5 Submit design — 6 Done

## Product

KAPA Target Enrichment Probes

## Select organism(s)

Human ? Other

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## Name your design

Custom Design Name

A unique Internal Reference Number (IRN) will also be automatically assigned at time of submission or approval. Custom Design Name and IRN are printed on the tube label. 0 / 14

Description

Eg. design targeting transcripts from a cancer gene panel, version 3

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## SUMMARY

Product KAPA Target Enrichment Probes

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The screenshot shows the 'Configure design' step in the HyperDesign workflow. The progress bar indicates the current step is '2 Configure design', with previous steps '1 Start design' and '3 Add targets' completed, and '4 Review targets', '5 Submit design', and '6 Done' remaining. The 'Advanced options' section contains five checkboxes with descriptions in Japanese and English. The 'Next' button is highlighted in blue. A 'SUMMARY' panel on the right shows design details: Custom Design Name, Product (KAPA Target Enrichment Probes), and Genomic build (Human (hg38/GRCh38)).

**HyperDesign** Home Your designs **+ New design** Help

**DRAFT**

1 Start design — **2 Configure design** — 3 Add targets — 4 Review targets — 5 Submit design — 6 Done

**Advanced options**

If you select any of the options below, completing your design will require collaboration with a Roche Expert Designer.

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Nothing above applies to your design?  
Simply click Next

**Back** **Next**

**SUMMARY**

Custom Design Name

Product **KAPA Target Enrichment Probes**

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The screenshot shows the HyperDesign web interface. At the top, there is a navigation bar with 'HyperDesign', 'Home', 'Your designs', and a '+ New design' button. On the right, there are links for 'Help', a user profile icon, and the Roche logo. The main content area has a green header with a document icon and a 'DRAFT' label. Below this is a progress bar with six steps: 'Start design', 'Configure design', 'Add targets' (the current step), 'Review targets', 'Submit design', and 'Done'. The 'Add targets' step is highlighted with a green circle. Below the progress bar, there is a text prompt: 'Use any or all of the following input types to define your targets.' Three options are presented in light blue boxes: 'Manual input' (with a pencil icon), 'Add file' (with a folder icon and 'Upload file size limit 150 MB'), and 'Add link' (with a link icon and 'Provide link to external file location (eg. URL, FTP server)'). At the bottom left is a 'Back' button and at the bottom right is a 'Next' button. On the right side of the interface, there is a 'SUMMARY' panel with a right arrow. It contains the following information: 'Custom Design Name' (blurred), 'Product' (KAPA Target Enrichment Probes), and 'Genomic build' (Human (hg38/GRCh38)).

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Type or paste gene IDs below or search by gene name using the panel on the right. Your input to the target regions will be mapped automatically.

**Your input**

Type or paste Gene IDs, for example:  
SMA14  
NM\_005359  
4089

Or enter genomic coordinates, for example  
chrN (tab) start (tab) stop (tab) comment

You must enter a value

**Gene search**

mlh

MLH1 HGNC symbol(s) Add gene

MLH1 WikiGene Name(s) Add gene

MLH3 HGNC symbol(s) Add gene

MLH3 WikiGene Name(s) Add gene

Format Select format

**Target options**

Protein coding transcript BioType regions only ?

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The screenshot shows the 'Add targets' step in the Roche Expert Designer workflow. A progress bar at the top indicates the current step is '3 Add targets', with previous steps 'Start design', 'Configure design', and 'Review targets' completed, and subsequent steps 'Submit design' and 'Done' pending. Below the progress bar, there are three input options: 'Manual input', 'File upload', and 'Add link'. The 'Add link' option is selected, and a dialog box is open. The dialog box has the title 'Add link - are you sure?' and contains the following text: 'Remember that your design request will be processed by a Roche Expert Designer since the HyperDesign automated tool is not capable of inputting data from an external link. For quicker probe selection, consider adding targets directly, using the file upload or manual input methods.' Below this text is a text input field with the placeholder text 'Provide link to file location (URL, FTP server) or describe how you will deliver your targets later on.' At the bottom right of the dialog box are 'Cancel' and 'Add' buttons. A black arrow points from the 'Add link' option in the background to the 'Add' button in the dialog box. In the bottom right corner of the overall image, there is orange text: 「Add」ボタンを押して次の画面に進みます.

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## ④ Review targets

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The screenshot displays a web application interface for reviewing targets. At the top, a green header bar contains a document icon, a 'DRAFT' label, and a navigation menu with six steps: 'Start design', 'Configure design', 'Add targets', 'Review targets' (highlighted with a green circle), 'Submit design', and 'Done'. Below the header, there is a 'Download view as BED' button and a 'View Consolidated regions' toggle switch. The main content area features a table with the following columns: a selection checkbox, a right-pointing arrow, 'Chr', 'Start', 'Stop', 'Gene', 'HGNC ID', 'Transcript', 'Transcript BioType', and 'Tags'. The table contains five rows of data, each with a selection checkbox, an arrow, and various numerical values for coordinates and identifiers. At the bottom of the table, there is a 'Back' button on the left, a 'Next' button on the right, and a summary of region counts: 'Initial regions count: 6120 (1,170,937 bp)' and 'Final regions count after consolidation: 1567 (478,394 bp)'. The interface also includes a 'Items per page' dropdown set to 5 and a pagination indicator '1 - 5 of 6120'.

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Choose probe stringency

Before submitting your design for selection, please define **maximum close matches** (?) in genome, depending how stringent or relaxed do you expect your probes to be.

Define maximum close matches in genome (?)

Stringent probes      Relaxed probes

From 1 to 4 - Best for unique regions      From 17 to 20 - Best for repetitive regions

Specify overhang in genome (?)

0 bp Overhang      125 bp Overhang

From 0 to 30 - optimal

Back      Submit your regions

**SUMMARY**

Custom Design Name: 20200507Desian

Product: KAPA Target Enrichment Probes

Genomic build: Human (hg19/GRCh37)

Initial intake: 1 files

Initial regions count: 6120 (1,170,937 bp)

Final regions count after consolidation: 1567 (478,394 bp)

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The screenshot shows a web interface for submitting a design. At the top, a green header contains the document name '20200507Design' and a 'DRAFT' status. Below the header is a progress bar with six steps: 'Start design', 'Configure design', 'Add targets', 'Review targets', '5 Submit design', and '6 Done'. The 'Submit design' step is currently active. On the left, there are sections for 'Choose probe stringency' and 'Define maximum close matches in genome', each with a slider. A 'Back' button is at the bottom left. On the right, a 'SUMMARY' panel lists design details. A white confirmation dialog box is centered on the screen, with a close button (X) in the top right corner. The dialog contains the following text:

Are you sure you want to submit the design for selection?

- Probe selection may take up to 24 hours and you will get an email notification when it's ready.
- After probe selection is complete, you will be able to review the estimated coverage over your target regions.
- Manufacturing will not begin until you approve the coverage and place an order.

Confirmation will be sent to the email provided in your profile settings:  
[hisano@genetics-n.co.jp](mailto:hisano@genetics-n.co.jp)

If it's not correct, please update your profile settings, before you submit your design.

At the bottom of the dialog are two buttons: 'Cancel' and 'Submit for selection'. An arrow points from the 'Submit for selection' button in the dialog to the 'Submit your regions' button on the main interface.

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The screenshot shows the Roche Expert Designer interface. At the top, there is a green header with a document icon, a "SUBMITTED" button, and a window icon. The main content is divided into two panels. The left panel, titled "SUBMITTED DESIGN", lists design parameters: Custom Design Name (blurred), Product (KAPA HyperChoice), Submitted (blurred), Genomic build (Human (hg19/GRCh37)), Stringency (Maximum close matches - 20), Probe overhang (30 bp), Initial intake (1 files), Initial regions count (6120 (1,170,937 bp)), and Final regions count after consolidation (1567 (478,394 bp)). The right panel, titled "SELECTION RESULTS", features a green DNA double helix icon with arrows. Below the icon, it states: "Your design has been successfully submitted for probe selection. You will be notified via Email when selection results are available here for your review." It also notes: "This design is no longer editable" and "To restart the design flow with the same input regions and submit it for selection again, please use the  Clone design option." At the bottom of the right panel, there are two paragraphs of Japanese text explaining the design process and notification timeline.

**SUBMITTED DESIGN**

Custom Design Name [blurred]

Product **KAPA HyperChoice**

Submitted [blurred]

Genomic build **Human (hg19/GRCh37)**

Stringency **Maximum close matches - 20**

Probe overhang **30 bp**

**Initial intake** 1 files

**Initial regions count** 6120 (1,170,937 bp)

**Final regions count after consolidation** 1567 (478,394 bp)

**SELECTION RESULTS**

Your design has been successfully submitted for probe selection.  
You will be notified via Email when selection results are available here for your review.

This design is no longer editable

To restart the design flow with the same input regions and submit it for selection again, please use the  Clone design option.

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**Initial intake** 1 files

**Initial regions count** 6120 (1,170,937 bp)

**Final regions count after consolidation** 1567 (478,394 bp)

### SELECTION RESULTS

7 May 2020, 02:12 PM

Selection for **20200507Design** design is completed and **ready for your review**. See the [Design Review and Approval Guidelines](#) for further information about the design review and approval process.

**Full report** **Probe selection results**

**Selection results summary**

Final sequence total	478,394 bp
Estimated coverage	452,782 bp
% target based covered	94.6%

**Estimated % coverage**  
**94.6%**

**Approve design** **Clone design**

**Ready to order?** Approve this design and work with your Local Roche Representative to get pricing and place an order.  
**Looking to make changes?** Clone your design and make adjustments to your design targets or probe stringency.

Ask your Local Roche Representative about the following **material numbers** and **pack sizes** corresponding to your design:

KAPA HyperChoice MAX 0.5Mb T1, 12 rxn (9052143001)  
KAPA HyperChoice MAX 0.5Mb T1, 24 rxn (9052151001)  
KAPA HyperChoice MAX 0.5Mb T1, 48 rxn (9052160001)

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HyperDesign

Dear [redacted]

You have approved your design [redacted]

Reach out to your Local Roche Representative [redacted] to choose pack sizes and place an order. Manufacturing will begin once your order is placed.

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