
Roche Applied Science

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The image shows a blue tray filled with small vials, likely containing PCR probes, and a single larger vial in the foreground. To the right, a computer screen displays the Roche Applied Science website interface for the Universal ProbeLibrary. The screen shows the following information:

Roche Applied Science

Universal ProbeLibrary for Human

ProbeFinder has designed the optimal real-time PCR assay for:

[NM_0011312](#) Homo sapiens actin, beta (ACTB), mRNA.

Assay details:

Use Universal ProbeLibrary probe: #11, cat.no. 04685105001
(Formerly Exon ProbeLibrary probe: Human#11)

Primer	Length	Position	Tm	%GC	Sequence
Left Primer	18	821 - 838	59	50	attggcaatgagcggttc
Right Primer	19	878 - 896	60	58	ggatgccacaggactccat

Probe (78 nt)

gagcgggtccgctgcctgaggcactcttcagccttccttcctggcatg
ggcatcc

Download peak meet Download report Order probes or set

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
Home Products Special Interest Sites Support & Resources News

Home > Special Interest Sites > Genomic Systems > Real-Time PCR Systems > **Universal ProbeLibrary System**

Universal ProbeLibrary

- ▼ Real-Time PCR Systems
 - ▶ LightCycler® Carousel-Based System
 - ▶ LightCycler® 480 System
 - ▼ **Universal ProbeLibrary System**
 - ▶ System Description
 - ▶ Technology
 - ▶ **Assay Design Center**
 - ▶ User Statements and Application
 - ▶ Assay List
 - ▶ Performance
 - ▶ Product List
 - ▶ Support
 - ▶ Literature and References
 - ▶ Multimedia Presentations
 - ▶ Product Information and Pack Inserts

Gene Expression Quantification with Real-Time PCR - Simple and Fast

- ◆ Design real-time qPCR assays online in seconds.
- ◆ Rely on just **165 prevalidated probes** for over five million qPCR assays for a large variety of organisms.
- ◆ Reduce the cost of gene expression analysis by performing **multiplex qPCR assays** with Universal ProbeLibrary Reference Gene Assays.
- ◆ Perform real-time PCR on any real-time PCR instrument with standard protocols.

✉ [Sign up](#) for E-Mail Services to stay updated on "Real-Time PCR Systems"

Ⓜ [View](#) Universal ProbeLibrary Multi Media Presentation

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The screenshot shows the 'Universal ProbeLibrary Assay Design Center' page from Roche Applied Science. On the left is a navigation menu with 'Assay Design Center' selected. The main content area includes a description of the ProbeFinder software and instructions to select an organism. A dropdown menu is open, listing various organisms such as Homo sapiens (Human), Mus musculus (Mouse), Rattus norvegicus (Rat), Anopheles gambiae (Mosquito), Caenorhabditis elegans (C. elegans), Danio rerio (Zebrafish), Drosophila melanogaster (Fruitfly), Pan troglodytes (Primates), Arabidopsis thaliana, Oryza sativa (Rice), Zea mays (Maize), Saccharomyces cerevisiae (Baker's yeast), and 'Other organism (paste sequence only)'.

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The screenshot shows the 'Universal ProbeLibrary for Human' target specification form. It asks the user to 'Specify your target(s):' and provides two options: 'By sequence ID, gene name or keyword' and 'By sequence'. The 'By sequence ID, gene name or keyword' option includes an example: 'e.g. ENST00000331789, NM_001101 or X00351 or beta-actin'. The 'By sequence' option includes an example: 'e.g. >part of X00351 Human mRNA for beta-actin' followed by a long DNA sequence: 'CACGGCATCGTCACCAACTGGGACGACATGGAGAAAAATCTGGCACCACACCTTCTACAAT GAGCTGCGTGTGGCTCCCGAGGAGCACCCCGTGCTGCTGACCGAGGCCCCCTGAACCCC AAGGCCAACCCGAGAGATGACCCAGATCATGTTTGAGACCTTCAACACCCCAGCCATG TACGTTGCTATCCAGGCTGTGCTATCCCTGTACGCCCTCTGGCCGTACCACTGGCATCGTG ATGGACTCCGGTGACGGGGTCACCCACACTGTGCCCATCTACGAGGGGTATGCCCTCCCC'.

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Automatically select an intron spanning assay. Design multiplex PCR with reference gene.

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Universal Probelibrary for Human Roche Applied

Multiple sequences were identified based on your submission.
Please choose the sequence(s) you would like to continue with. You can select up to 10 sequences.
Displaying number 1 to 30 of the 399 genes matched by your query.

Name	Length	Description
<input type="checkbox"/> ENST00000117192.2	1840	TBBE_HUMAN Uniprot@SWISSPROT Tubulin beta-6 chain [Source: Uniprot@SWISSPROT; Acc: G9BUJ5]
<input type="checkbox"/> ENST00000134295.3	2064	TP052_HUMAN Uniprot@SWISSPROT Tubulin polyglutamylase complex subunit 2 (P042). [Source: Uniprot@SWISSPROT; Acc: G66CL5]
<input type="checkbox"/> ENST00000183058.1	1935	TP052_HUMAN Uniprot@SWISSPROT Tubulin polyglutamylase complex subunit 2 (P042). [Source: Uniprot@SWISSPROT; Acc: G66CL5]
<input type="checkbox"/> ENST00000171025.4	2123	G550V5_HUMAN Uniprot@SPTREMBL tubulin, beta 8 [Source: RefSeq_peptide; Acc: NP_817124]
<input checked="" type="checkbox"/> ENST00000109812.2	1436	NP_817124.1 RefSeq_peptide tubulin, beta 8 [Source: RefSeq_peptide; Acc: NP_817124]
<input type="checkbox"/> ENST00000128874.2	1440	NP_817124.1 RefSeq_peptide tubulin, beta 8 [Source: RefSeq_peptide; Acc: NP_817124]
<input type="checkbox"/> ENST00000132108.5	820	NP_817124.1 RefSeq_peptide tubulin, beta 8 [Source: RefSeq_peptide; Acc: NP_817124]
<input type="checkbox"/> ENST00000108418.1	1798	NP_879078.1 RefSeq_peptide tubulin, alpha-like 3 [Source: RefSeq_peptide; Acc: NP_879078]
<input type="checkbox"/> ENST00000115811.3	1381	NP_879078.1 RefSeq_peptide tubulin, alpha-like 3 [Source: RefSeq_peptide; Acc: NP_879078]
<input type="checkbox"/> ENST00000175939.1	2629	NP_001008409.1 RefSeq_peptide tubulin bromine igase-like family, member 9 [Source: RefSeq_peptide; Acc: NP_001008409]
<input type="checkbox"/> ENST00000175934.1	3261	NP_001008409.1 RefSeq_peptide tubulin bromine igase-like family, member 9 [Source: RefSeq_peptide; Acc: NP_001008409]
<input type="checkbox"/> ENST00000175934.1	2967	NP_001008409.1 RefSeq_peptide tubulin bromine igase-like family, member 9 [Source: RefSeq_peptide; Acc: NP_001008409]
<input type="checkbox"/> ENST00000175932.1	3038	NP_001008409.1 RefSeq_peptide tubulin bromine igase-like family, member 9 [Source: RefSeq_peptide; Acc: NP_001008409]

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The screenshot shows a software interface with two checked options: "Automatically select an intron spanning assay." and "Design multiplex PCR with reference gene." Below these is a dropdown menu for "Preferred Reference Gene:" with "Any" selected. The dropdown list includes "Any", "ACTB", "HPRT", "PBGD", and "PGK1". A "Design" button is visible to the left of the dropdown. At the bottom left, the text "ProbeFinder Version: 2.40" is displayed.

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Real-Time PCR Systems Roche Applied Science

Figure legend

Score (mouse over for details)	
All criteria met	
All criteria met	Scoring details for probe no. 9
All criteria met	Amplicon length (76 nt):
All criteria met	Insilico PCR:
All criteria met	

Sequence The light-purple bar illustrates the input sequence and the numbers refer to the nucleotide positions in the input sequence.

Amplicons The length and positions of the PCR amplicons are depicted as bars along the sequence. Enlarged views also show the positions of probes and primers in the amplicons. The numbers refer to the ProbeLibrary probe matching the individual amplicons.

Introns The wedges indicate the exon-exon splice junctions with the intron sizes depicted in the enlarged view. Mouse over an intron to see information about the position.

SNP The magenta color illustrates that SNP (Single Nucleotide Polymorphism) is found at the marked position. Run the mouse over magenta box for more details on this SNP. Mouse over a SNP to see extra information.

Overview The example shows multiple PCR amplicons spanning a 110 nucleotide intron and their positions relative to the input sequence. Mouse over each amplicon to see scoring information.


Scoring The example shows the mouse over the scoring column. This provides additional information about the scoring.

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Probe no. 77 scoring details	
Amplicon length (71 nt):	Green
Intron spanning:	Orange
In silico PCR:	Green
SNP:	Green
Overall:	Orange

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Universal ProbeLibrary for Human *Roche Applied Science*

Please choose the sequence(s) you would like to continue with.

	Name	Length	Description
❑	NM_002751.4	2222	Homo sapiens mitogen-activated protein kinase 11 (MAPK11), transcript variant 1, mRNA.
❑	NM_138993.1	1530	Homo sapiens mitogen-activated protein kinase 11 (MAPK11), transcript variant 2, mRNA.

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ProbeFinder has designed the optimal real-time PCR assays for the 2 genes :

[NM_002751.5](#) Homo sapiens mitogen-activated protein kinase 11 (MAPK11), transcript variant 1, mRNA.
[NM_138993.2](#) Homo sapiens mitogen-activated protein kinase 11 (MAPK11), transcript variant 2, mRNA.

Assay details:

Use Universal ProbeLibrary probes: #1(Human#01) and #2(Human#02)

Input	Probe ID	Primer	Length	Position	Tm	%GC	Sequence
NM_002751.5	#1 cat.no. 04684974001	Left	21	647 - 667	59	38	tggatgcattacaaccaaca
		Right	21	731 - 751	59	52	ctggtcaatgtagtgcgttcc
		Amplicon (105 nt) tggatgcattacaaccaacagtggatctctggtccgtgggctgcatcatggctgagctg ctccagggaagccctctcccggaagcgactacattgaccag					
NM_138993.2	#2 cat.no. 04684982001	Left	20	1241 - 1260	59	55	cccacactcctatgggtcaca
		Right	20	1313 - 1332	60	55	tccgaggttactggatcagc
		Amplicon (92 nt) cccacactcctatgggtcagactctggcctaggaccctcgcttcaggagaatctac acgcatgatggagctgatccagtaacctcgga					

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ProbeFinder has designed real-time PCR assays **differentiating** between the 2 input sequences :

[NM_002751.5](#) Homo sapiens mitogen-activated protein kinase 11 (MAPK11), transcript variant 1, mRNA.
[NM_138993.2](#) Homo sapiens mitogen-activated protein kinase 11 (MAPK11), transcript variant 2, mRNA.

Use Universal ProbeLibrary probe(s): **#27(Human#27)** and **#2(Human#02)**

Input	Probe ID	Primer	Length	Position	Tm	%GC	Sequence
NM_002751.5	#27 cat.no. 04687582001	Left	19	1988 - 2006	60	53	aggggctttctccttgat
		Right	20	2093 - 2112	59	55	gtccgagtccaagtccacat
Amplicon (125 nt) aggggctttctccttgatgtgggaccacagcagcttctgaggctgggggtgggtggg tgggtggtttggccttgaggacgctagggcaggcagcacacctggatgtggacttggact cggac							
NM_138993.2	#2 cat.no. 04684982001	Left	20	1241 - 1260	59	55	cccacactcctatggtcaca
		Right	20	1313 - 1332	60	55	tccgaggttactggatcagc
Amplicon (92 nt) cccacactcctatggtcacagacttctggcctaggaccctcgcttcaggagaatctac acgcatgatggagctgatccagttaacctcgga							

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ProbeFinder has designed real-time PCR assays that are **common** for the 2 genes :

[NM_002751.5](#) Homo sapiens mitogen-activated protein kinase 11 (MAPK11), transcript variant 1, mRNA.
[NM_138993.2](#) Homo sapiens mitogen-activated protein kinase 11 (MAPK11), transcript variant 2, mRNA.

Use Universal ProbeLibrary probe: #1(cat.no.04684974001)
(Formerly Exiqon ProbeLibrary Probe: Human#01)

Input	Probe ID	Primer	Length	Position	Tm	% GC	Sequence
NM_002751.5	#1 cat.no. 04684974001	Left	21	647 - 667	59	38	tggatgcattacaaccaaca
		Right	21	731 - 751	59	52	ctggtcaatgtagtcgcttcc
Amplicon (105 nt) tggatgcattacaaccaaacagtggatctctggtccgctggctgcatcatggctgagctg ctccagggcaaggccctcttcccgggaagcgactacattgaccag							

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(例: ENST00000158302, NM_001101, AB062273 など)

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Zhihua Xie, Susan Westmoreland, Mary E Bahn, Guo-Lin Chen, Hong Yang, Eric Vallender, Wei-Dong Yao, Bertha K Madras, and Gregory M. Miller (2007). Rhesus monkey trace amine-associated receptor 1 signaling: enhancement by monoamine transporters and attenuation by the D2 autoreceptor *in vitro*. *J Pharmacol Exp Ther.* 321(1):116-27.

Christoph Leucht and Laure Bally-Cuif (2007). The Universal ProbeLibrary – A versatile Tool for Quantitative Expression Analysis in the Zebrafish. *Biochemica.2*, 16-16

Ina Horst and Christoph Peterhänsel (2007). Quantification of Zea mays mRNAs by Real-Time PCR Using the Universal ProbeLibrary. *Biochemica 1*, 8-10

Li-Hui Wang, Seok-Hyung Kim, Jung Hyun Lee, Yoon-La Choi, Young Chul Kim, Tae Sung Park, Yun-Chul Hong, Chun-Fu Wu and Young Kee Shin (2007). Inactivation of SMAD4 Tumor Suppressor Gene During Gastric Carcinoma Progression. *Clinical Cancer Research 13*: 102-110

Nien-Yi Lin, Chung-Tien Lin, Yu-Ling Chen, Ching-Jin Chang (2007). Regulation of tristetraprolin during differentiation of 3T3-L1 preadipocytes. *FEBS Journal 274 (3)*, 867-878

Maldonado-Saldivia J, van den Bergen J, Krouskos M, Gilchrist M, Lee C, Li R, Sinclair AH, Surani MA, Western PS (2007). Dppa2 and Dppa4 are closely linked SAP motif genes restricted to pluripotent cells and the germ line. *Stem Cells 25*: 19-28

Aguilar-Martínez JA, Poza-Carrión C, Cubas P. Arabidopsis BRANCHED1 acts as an integrator of branching signals within axillary buds. *Plant Cell. (2007) ;19(2):458-72.*

Kimberly Johung, Edward C. Goodwin, and Daniel DiMaio (2007). Human Papillomavirus E7 Repression in Cervical Carcinoma Cells Initiates a Transcriptional Cascade Driven by the Retinoblastoma Family, Resulting in Senescence. *Journal of Virology 81, No. 5*, 2102-2116

Seyed H. Hosseini, James J. Kohler, Chad P. Haase, Nina Tioleco, Tami Stuart, Erin Keebaugh, Tomika Ludaway, Rodney Russ, Elgin Green, Robert Long, Liya Wang, Staffan Eriksson and William Lewis (2007). Targeted Transgenic Over-expression of Mitochondrial Thymidine Kinase (TK2) Alters Mitochondrial DNA (mtDNA) and Mitochondrial Polypeptide Abundance. *American Journal of Pathology 170:865-874.*

Rong-Mei Wu, Marion Wood, Anthony Thrush, Eric F. Walton, and Erika Varkonyi-Gasic (2007). Real-time PCR quantification of plant miRNAs using Universal ProbeLibrary technology. *Biochemica 2*, 12-15

Esther Carrera, Tara Holman, Anne Medhurst, Wendy Peer, Heike Schmutz, Steven Footitt, Frederica L. Theodoulou, and Michael J. Holdsworth (2007). Gene Expression Profiling Reveals Defined Functions of the ATP-Binding Cassette Transporter COMATOSE Late in Phase II of Germination. *Plant Physiol.* 143(4): 1669–1679.

M. Julieta Del Prete, Rolando Vernal, Helmut Dolznig, Ernst W. Müllner, and Jose A. Garcia-Sanz (2007). Isolation of polysome-bound mRNA from solid tissues amenable for RT-PCR and profiling experiments. *RNA* 13: 414–421.

Sahin O, Löbke C, Korf U, Appelhans H, Sülthmann H, Poustka A, Wiemann S, Artl D (2007). Combinatorial RNAi for quantitative protein network analysis. *Proc Natl Acad Sci USA* 104(16):6579–84.

Kuo, W.P., et. al. (2006). A sequence-oriented comparison of gene expression measurements across different hybridization-based technologies. *Nature Biotechnology* 24 (7): 832–840.

Ulrike Bichelmeier, Thorsten Schmidt, Jeannette Hübener, Jana Boy, Lukas Rüttiger, Karina Häbig, Sven Poths, Michael Bonin, Marlies Knipper, Werner J. Schmidt, Johannes Wilbertz, Hartwig Wolburg, Franco Laccone, and Olaf Riess. Nuclear Localization of Ataxin-3 Is Required for the Manifestation of Symptoms in SCA3: *In Vivo* Evidence *The Journal of Neuroscience* (2007) 27(28), 7418 –7428

Nabanita Nag, Katherine Peterson, Keith Wyatt, Sonja Hess, Sugata Ray, Jack Favor, Debora Bogani, Mary Lyon, and Graeme Wistow. Endogenous retroviral insertion in Cryge in the mouse No3 cataract mutant *Genomics* (2007) 89, 512–520.

Liping Qiao, Paul S. MacLean, Hanning You, Jerome Schaack, and Jianhua Shao
Knocking Down Liver CCAAT/Enhancer-Binding Protein alpha by Adenovirus-Transduced Silent Interfering Ribonucleic Acid Improves Hepatic Gluconeogenesis and Lipid Homeostasis in db/db Mice. *Endocrinology* (2006) 147(6), 3060–3069.

Annemiek Beverdam and Peter Koopman (2006). Expression profiling of purified mouse gonadal somatic cells during the critical time window of sex determination reveals novel candidate genes for human sexual dysgenesis syndromes. *Human Molecular Genetics* 15(3):417–431.

Jalila Chagraoui, Sherry L. Niessen, Julie Lessard, Simon Girard, Philippe Coulombe, Martin Sauvageau, Sylvain Meloche and Guy Sauvageau (2007). E4F1: a novel candidate factor for mediating BMI1 function in primitive hematopoietic cells. *Genes & Dev.* 20: 2110–2120.

Amy K. Jassen, Hong Yang, Gregory M. Miller, Elizabeth Calder, and Bertha K. Madras (2006). Receptor Regulation of Gene Expression of Axon Guidance Molecules: Implications for Adaptation. *Mol Pharmacol.* 70(1):71–7.

Lind K, Ståhlberg A, Zoric N, Kubista M.(2006). Combining sequence-specific probes and DNA binding dyes in real-time PCR for specific nucleic acid quantification and melting curve analysis. *Biotechniques* 40(3): 315-9

Guccione, E., Martinato, F., Finocchiaro, G., Luzi, L., Tizzoni, L., Dall' Olio, V., Zardo, G., Bernard, L., Amati, B. (2006). Myc-binding-site recognition in the human genome is determined by chromatin context. *Nature Cell Biology* 8(7): 764-70

Ann M. O'Hara, Padraig O'Regan, Áine Fanning, Caitlin O'Mahony, John MacSharry, Anne Lyons, John Bienenstock, Liam O'Mahony, and Fergus Shanahan
Functional modulation of human intestinal epithelial cell responses by *Bifidobacterium infantis* and *Lactobacillus salivarius*.
Immunology (2006) 118(2), 202-215

Rein, R.S., Mauritz, R.P. and Geyer M. (2006). Universal ProbelLibrary: a new concept for streamlining gene expression analysis with qPCR. *Nature Methods Application Notes*

Onno Bakker (2006). Comparison of qPCR Assays Using the LightCycler® 2.0 Instrument and either SYBR Green I Intercalation or the Universal ProbelLibrary as Detection Format. *Biochemica* 3, 8-10.

Granado, M., Priego, T., Martin, A.I., Villanua, M.A., Lopez-Calderon, A.(2005). Ghrelin receptor agonist GHRP-2 prevents arthritis-induced increase in E3 ubiquitin-ligating enzymes MuRF1 and MAFbx gene expression in skeletal muscle. *Am J Physiol Endocrinol Metab.* 289(6).

Taranger, C.K., Agate Noer, Sørensen, A.L., Håkelién, AM., Boquest, A.C., Collas, P.(2005). Induction of Dedifferentiation, Genome-wide Transcriptional Programming and Epigenetic Reprogramming by Extracts of Carcinoma and Embryonic Stem Cells. *Molecular Biology of the Cell* 16: 5719 - 5735.

Viola Gesellchen, David Kuttenkeuler, Michael Steckel, Nadège Pelte & Michael Boutros (2005). An RNA interference screen identifies Inhibitor of Apoptosis Protein 2 as a regulator of innate immune signalling in *Drosophila*. *EMBO Rep.* 6(10):979-84.

Kruhøffer, M, Jensen, J.L., Laiho, P., Dyrskjot, L., Salovaara, R., Arango, D., Birkenkamp-Demtroder, K., Sorensen, F.B., Christensen, L.L., Buhl, L., Mecklin, J.P., Jarvinen, H., Thykjaer, T., Wikman, F.P., Bech-Knudsen, F., Juhola, M., Nupponen, N.N., Laurberg, S., Andersen, C.L., Aaltonen, L.A., Ørntoft, T.F. (2005). Gene expression signatures for colorectal cancer microsatellite status and HNPCC. *Br J Cancer.* 20;92(12):2240-8.

Mouritzen, P., Noerholm, M., Nielsen, P.S., Jacobsen, N., Lomholt, C., Pfundheller H.M. and Tolstrup N.(2005). ProbelLibrary: A new method for faster design and execution of quantitative real-time PCR. *Nature Methods* 2: 313-17

Magnusson, N.E., Larsen, A., Rungby, J., Kruhøffer, M., Ørntoft, T.F., Stoltenberg, M.(2005). Gene expression changes induced by bismuth in a macrophage cell line. *Cell Tissue Res* 321: 195210.

Mauritz, R.P., Mouritzen, P., Pfundheller, H.M., Tolstrup, N., and Lomholt, C. (2005). Universal ProbeLibrary Set: One Transcriptome – One Kit. *Biochemica* 2, 22-24.

Effi Rees, Bianca Horeis, Nadja Gernold, Oliver Heil, and Bernhard Korn (2005). Gene Silencing Using esiRNA – Efficient, Robust and Not Influenced by Positional Effects. *Biochemica*.3, 26–29.

Michael Steckel and Michael Boutros (2005). Rapid Development of Real-Time RT-PCR Assays Using Universal ProbeLibrary: Applications for Dissecting Signaling Pathways by RNA Interference. *Biochemica* .3, 17–19.

Mouritzen, P., Nielsen, P.S, Jacobsen, N., Noerholm, M., Lomholt, C., Pfundheller, H.M., Ramsing, N.B., Kauppinen, S., and Tolstrup, N. (2004). The ProbeLibrary™ – Expression profiling 99% of all human genes using only 90 dual-labeled real-time PCR Probes. *Biotechniques* (2004) 37:492–495.

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Disclaimer of License

A license to perform the patented 5' Nuclease Process for research is obtained by the purchase of (i) both Licensed Probe and Authorized 5' Nuclease Core Kit, (ii) a Licensed 5' Nuclease Kit, or (iii) license rights from Applied Biosystems.

This product contains Licensed Probe. Its purchase price includes a limited, non-transferable immunity from suit under U.S.

Patents Nos. 6,214,979 and 5,804,375 (claims 1-12 only) and corresponding patent claims outside the United States. The

purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using

only this amount of product for the purchaser's own internal research. The right to use this product in the 5' Nuclease

Process under the applicable claims of US Patents Nos. 5,210,015 and 5,487,972, and corresponding patent claims outside the

United States, can be obtained through purchase of an Authorized 5' Nuclease Core Kit. Except under separate license

rights available from Applied Biosystems, no right under any other patent claim, or to perform commercial services of any kind,

including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, or to

sublicense, repackage with other products, or resell in any form, is conveyed expressly, by implication, or by estoppel. This product

is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on

purchasing licenses may be obtained from the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City,

California 94404, USA.

The technology used for the LightCycler® System is licensed from Idaho Technology Inc., Salt Lake City, UT, USA.

ProbeLibrary is covered by US and other patent applications owned by Exiqon A/S. Locked Nucleic Acids (LNA) are covered by

U.S. Patents No US 6,794,499, US 6,670,461, US 6,268,490 & US 6,770,748 and other patents and patent applications owned

by Exiqon A/S and Prof. Takeshi Imanishi. The quencher used in the probes is covered by patent applications owned by Exiqon

A/S.

Trademarks

LIGHTCYCLER, LC, FASTSTART, TAQMAN, MAGNA PURE, and HIGH PURE are trademarks of Roche.

PROBELIBRARY and LNA are registered trademarks of Exiqon A/S, Vedbaek, Denmark.

FAM, TAMRA, and ROX are trademarks of Applied Biosystems Corporation, Norwalk, CT, USA.

Other brand or product names are trademarks of their respective holders.

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