MagNA Pure DNA Tissue Lysis Buffer

REF 06 640 702 001

 **Version 06**
 Content version: August 2021
 Store at +15 to +25°C

1. Intended Use

The MagNA Pure 96 and 24 Systems are automated nucleic acid purification systems consisting of the instrument, software, control unit (only for MagNA Pure 96 System), consumables, and reagents. The MagNA Pure 96 and 24 Systems are intended for use by professional users for the purification of nucleic acids from biological samples for *in vitro* diagnostic purposes.

The MagNA Pure DNA Tissue Lysis Buffer is for use with the MagNA Pure 96 and 24 Systems.

2. Explanation of the Reagent

The MagNA Pure DNA Tissue Lysis Buffer is designed for:

- Lysis of many different types of tissues.
- Stabilization of DNA in tissue lysates.
- Purification of genomic DNA from tissue sample lysates using the MagNA Pure 96 System and MagNA Pure 24 System.

3. Reagent Principle/ Summary

The MagNA Pure DNA Tissue Lysis Buffer is used to preprocess tissue samples for isolation of genomic DNA using the MagNA Pure 96 Instrument and MagNA Pure 24 Instrument.

4. Reagents - Working Solutions

| Vial/Cap | Label | Contents/Function |
|---------------|------------------------------------|--|
| 1 (colorless) | MagNA Pure DNA Tissue Lysis Buffer | <ul style="list-style-type: none"> ▪ one bottle with 100 ml ▪ ready-to-use reagent for lysis of tissue ▪ < 6 M Urea, < 1 M Na-Citrat, < 100 mM NaCl, < 500 mM K₂HPO₄/KH₂PO₄ |

5. Precautions and Warnings

5.1 Precautions

 The MagNA Pure DNA Tissue Lysis Buffer contains no hazardous substances in concentrations that need to be declared. Nevertheless, always wear gloves and follow standard safety precautions to minimize contact when handling.

5.2 Handling Instructions

- Do not use the reagent after its expiration date has passed.
- Efficient disruption and homogenization of the sample material is essential for isolation of genomic DNA from tissues. Incomplete tissue disruption will result in significantly reduced DNA yields. On the other hand, too extensive disruption and homogenization will lead to shearing of high-molecular genomic DNA. Therefore, the tissue lysis procedure shall be validated. In addition, to minimize the risk of carryover contamination which may result in false positive results, follow the guidelines listed below:

- Perform sample preparation, PCR/RT-PCR setup and PCR/RT-PCR in physically-separated locations.
 - Discard pipette tips in sealed containers to prevent airborne contamination.
- Nuclease contaminated reagents and reaction vessels will degrade template nucleic acid. Please follow these guidelines to minimize the risk of contamination:
- Avoid touching surfaces or materials that could cause nuclease carryover.
 - Clean and decontaminate work areas and instruments, including pipettes, with commercially available decontamination reagents.
 - Use only new nuclease-free aerosol-blocking pipette tips and microcentrifuge tubes.
 - Use a work area specifically designed for RNA work, and if possible use reaction vessels and pipettors dedicated only for work with template RNA.

5.3 Laboratory Procedures

- All human sourced material and all resulting waste should be considered potentially infectious. Thoroughly clean and disinfect all work surfaces with disinfectants recommended by the local authorities.
- As the sensitivity and titer of potential pathogens in the sample material can vary, the operator must optimize pathogen inactivation and follow the appropriate measures according to local safety regulations.
- Do not eat, drink or smoke in the laboratory work area.
- Do not pipette by mouth.
- Wear protective disposable gloves, laboratory coats, and eye protection when handling specimens and reagents.
- Avoid microbial and nuclease contamination of reagents when removing aliquots from reagent bottles. Use sterile disposable pipette tips.
- Wash hands thoroughly after handling specimens and reagent.

5.4 Waste Handling

- For US: Material Safety Data Sheets (MSDS) are available on www.usdiagnostics.roche.com, or upon request from the local Roche office.
- For all other countries: Material Safety Data Sheets (MSDS) are available online at www.e-labdoc.roche.com, or upon request from the local Roche office.
- Dispose unused reagents and waste in accordance with country, federal, state, and local regulations.

6. Storage and Stability (Reagents)

- Store the MagNA Pure DNA Tissue Lysis Buffer at +15°C to +25°C.
- When properly handled, the MagNA Pure DNA Tissue Lysis Buffer is stable until the expiration date printed on the label.

7. Materials

7.1 Materials Provided

See Reagents - Working Solutions

7.2 Materials and Devices Required but not Provided

- Standard laboratory equipment
- Vortex mixer
- MagNA Pure 96 DNA and Viral NA Large Volume Kit (Cat. No. 06 374 891 001)
- MagNA Pure 96 DNA and Viral NA Small Volume Kit (Cat. No. 06 543 588 001)
- MagNA Pure 96 Instrument (Cat. No. 06 541 089 001)
- MagNA Pure 96 System Fluid (Internal) (Cat. No. 06 430 112 001)
- MagNA Pure 96 System Fluid (External) (Cat. No. 06 640 729 001)
- MagNA Pure 24 Total NA Isolation Kit (Cat. No. 07 658 036 001)
- MagNA Pure 24 Instrument (Cat. No. 07 290 519 001)
- Optional, Proteinase K, PCR grade, Activity (+37°C) ≥ 0.6 U/μl (e.g., Cat. No. 03 115 828 001)
- Optional, MagNA Lyser Instrument (Cat. No. 03 358 968 001 as of SN 40467540, Cat. No. 03 358 976 001 as of SN 40405218) or similar device for homogenization of fresh-frozen tissue
- Optional, MagNA Lyser Green Beads (Cat. No. 03 358 941 001)

8. Assay Procedures

8.1 General Remarks

It is the user's responsibility to validate system performance for any procedures used in their laboratory.

8.2 Purification Protocol

For a detailed description regarding the assay procedure, please refer to the Instructions for Use/Method Sheets for **MagNA Pure 96 DNA and Viral NA Small Volume Kit**, **MagNA Pure 96 DNA and Viral NA Large Volume Kit** and **MagNA Pure 24 Total NA Isolation Kit**.

9. Quality Control

 Always run appropriate controls.

To control the entire process, starting from sample preparation to analysis, perform the following controls:

- Positive control**, using a sample material positive for target.
- Negative control**, using PBS in place of the sample.
- Extraction control**, using a sample material negative for target.
- Internal control (IC)**, by adding a defined amount of a control template to all samples to be purified.

For applications that could produce false negative results, such as the detection of pathogens, the use of an appropriate internal control (IC) is mandatory. The IC is added during nucleic acid isolation, preferably using the automated IC function of the MagNA Pure 96 System and MagNA Pure 24 System. The IC can also be added manually to the sample. In this case, the IC must be stable in the sample material, and a nuclease-sensitive IC, such as unprotected RNA, should not be used for this purpose.



10. Limitations and Interferences

- Reliable results are dependent on appropriate specimen collection, transport, storage and processing procedures.
- Use of this product should be limited to personnel trained in techniques of nucleic acid purification and isolation and PCR.
- False negative results may occur if a specimen is improperly collected, transported or handled. False negative results may also occur if insufficient amount of template is present in the specimen.
- Any IVD application using the sample preparation procedure in conjunction with any downstream IVD nucleic acid testing should be evaluated with regard to the individual IVD parameter.
- To minimize the risk of a negative impact on the results, adequate controls should be used.
- Storage conditions (temperature, time) for lysates shall be validated with regard to the individual IVD parameter.

11. Supplementary Information

11.1 Symbols

In this Instruction Manual, the following symbols are used to highlight important information:

| Symbol | Description |
|---|--|
|  | Important Note |
|  | Information Note |
|  | For <i>in vitro</i> diagnostic use. |
|  | The reagent complies with the requirements of the IVDR Regulation (EU) 2017/746. |
|  | Catalogue Number |
|  | Global Trade Item Number |
|  | Unique Device Identifier |
|  | Batch Code |
|  | Use-by date |
|  | Date of Manufacture |
|  | Content of Kit |
|  | Temperature Limit |
|  | Consult instructions for use |
|  | Distributed by |
|  | Manufacturer |
|  | Authorized representative in the European Community |
|  | Importer |

11.2 Changes to previous version

- Change of the name of the product.
- Update to comply with the requirements of the IVDR Regulation (EU) 2017/746.
- Add reference to MagNA Pure 24 System.

12. Trademarks

MAGNA PURE is a trademark of Roche.

All other product names and trademarks are the property of their respective owners.

13. Regulatory Disclaimer

For *in vitro* diagnostic use.

Roche Molecular Systems, Inc.
1080 US Highway 202 South
Branchburg, NJ 08876 USA
Manufactured in Germany

Roche Diagnostics GmbH
Sandhofer Strasse 116,
68305 Mannheim, Germany
+49 621 759 0

Distributed in USA by Roche Diagnostics, Indianapolis, IN, USA
US Customer Technical Support: 1-800-526-1247

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Das Reagenz MagNA Pure DNA Tissue Lysis Buffer dient zur:

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