# MagNA Pure External Lysis Buffer

REF 06 374 913 001

**Version 08**  
Content version: August 2021  
Store at +15 to +25°C

## 1. Intended Use

The MagNA Pure 96 and 24 Systems are automated nucleic acid purification systems consisting of the instrument, software, control unit (only for MagNA Pure 96 System), consumables, and reagents. The MagNA Pure 96 and 24 Systems are intended for use by professional users for the purification of nucleic acids from biological samples for *in vitro* diagnostic purposes.

The MagNA Pure External Lysis Buffer is for use with the MagNA Pure 96 and 24 Systems.

## 2. Explanation of the Reagent

The MagNA Pure External Lysis Buffer is designed for:

- Lysis of different/various sample materials, e.g.,
  - Whole blood
  - Serum
  - Plasma, EDTA/citrate-preserved
- Stabilization of nucleic acids with lysates.
- Nucleic acid purification using the MagNA Pure 96 System and the MagNA Pure 24 System.

## 3. Reagent Principle/ Summary

The MagNA Pure 96 External Lysis Buffer is used to lyse cells and inactivate enzymes, such as nucleases.

## 4. Reagents - Working Solutions

Vial/ Cap	Label	Contents/Function
green cap	MagNA Pure External Lysis Buffer	<ul style="list-style-type: none"> <li>▪ one bottle with 100 ml Blue Reagent</li> <li>▪ ready-to-use reagent for cell lysis and binding of total nucleic acids</li> <li>▪ &lt; 6 M GuSCN, &lt; 30 % Triton X-100, &lt; 100 mM Tris-HCl, 0.01 % Bromophenol blue</li> </ul>

## 5. Precautions and Warnings

### 5.1 Precautions

- ⚠ The MagNA Pure External Lysis Buffer contains guanidine isothiocyanate which is an irritant. Always wear gloves and follow standard safety precautions to minimize contact when handling.
- ⚠ Do not allow this buffer to touch the skin, eyes, or mucous membranes. If contact does occur, wash the affected area immediately with large amounts of water. If the reagents are spilled, dilute the spill with water before wiping it up.
- ⚠ Do not allow reagents containing guanidine thiocyanate to contact sodium hypochlorite (bleach) solution or acids. These mixtures produce a highly toxic gas.

- ⚠ Do not use the buffer if it contains precipitates. If a precipitate is visible, place the bottle at +37°C and mix from time to time until the precipitate is completely dissolved. Do not warm the buffer longer at +37°C than is actually needed for complete dissolution of the precipitate. Ensure that the MagNA Pure External Lysis Buffer is equilibrated to +15 to +25°C before use.

- ⚠ Handle all samples as if potentially infectious. Follow standard laboratory safety procedures for hazardous materials.

### 5.2 Handling Instructions

- Wear disposable gloves and change them frequently.
  - Do not use the reagent after its expiration date has passed.
- In addition, to minimize the risk of carryover contamination which may result in false positive results, follow the guidelines listed below:
- Perform sample preparation, PCR/RT-PCR setup and PCR/RT-PCR in separate locations.
  - Discard pipette tips in sealed containers to prevent airborne contamination.

Nuclease contaminated reagents and reaction vessels will degrade template nucleic acid. Please follow these guidelines to minimize the risk of contamination:

- Avoid touching surfaces or materials that could cause nuclease carryover.
- Clean and decontaminate work areas and instruments, including pipettes, with commercially available decontamination reagents.
- Use only new nuclease-free aerosol-blocking pipette tips and microcentrifuge tubes.
- Use a work area specifically designed for RNA work, and if possible use reaction vessels and pipettors dedicated only for work with template RNA.

### 5.3 Laboratory Procedures

- All human sourced material and all resulting waste should be considered potentially infectious. Thoroughly clean and disinfect all work surfaces with disinfectants, as recommended by the local authorities.
- As the sensitivity and titer of potential pathogens in the sample material can vary, the operator must optimize pathogen inactivation and follow the appropriate measures according to local safety regulations.
- Do not eat, drink or smoke in the laboratory work area.
- Do not pipette by mouth.
- Wear protective disposable gloves, laboratory coats, and eye protection when handling specimens and reagent.
- Avoid microbial and nuclease contamination of reagents when removing aliquots from reagent bottles. Use sterile disposable pipette tips.
- Wash hands thoroughly after handling specimens and reagents.

## 5.4 Waste Handling

- For US: Material Safety Data Sheets (MSDS) are available on [www.usdiagnostics.roche.com](http://www.usdiagnostics.roche.com), or upon request from the local Roche office.
- For all other countries: Material Safety Data Sheets (MSDS) are available online at [www.e-labdoc.roche.com](http://www.e-labdoc.roche.com), or upon request from the local Roche office.
- Dispose unused reagents and waste in accordance with country, federal, state, and local regulations.

## 5.5 For customers in the European Economic Area

Contains SVHC: octyl/nonylphenol ethoxylates. For use as part of an IVD method and under controlled conditions only – acc. to Art. 56.3 and 3.23 REACH Regulation.

## 6. Storage and Stability (Reagents)

- The MagNA Pure External Lysis Buffer is shipped at ambient temperature.
- Store the MagNA Pure External Lysis Buffer at +15°C to +25°C.
- When properly handled, the MagNA Pure External Lysis Buffer is stable until the expiration date printed on the label.

## 7. Materials

### 7.1 Materials Provided

see Reagents - Working Solution

### 7.2 Materials and Devices Required but not Provided

- Standard laboratory equipment
- Vortex mixer
- MagNA Pure 96 DNA and Viral NA Large Volume Kit (Cat. No. 06 374 891 001)
- MagNA Pure 96 DNA and Viral NA Small Volume Kit (Cat. No. 06 543 588 001)
- MagNA Pure 96 Instrument (Cat. No. 06 541 089 001)
- MagNA Pure 96 System Fluid (Internal) (Cat. No. 06 430 112 001)
- MagNA Pure 96 System Fluid (External) (Cat. No. 06 640 729 001)
- MagNA Pure 24 Total NA Isolation Kit (Cat. No. 07 658 036 001)
- MagNA Pure 24 Instrument (Cat. No. 07 290 519 001)

## 8. Assay Procedures

### 8.1 General Remarks

It is the user's responsibility to validate system performance for any procedures used in their laboratory.

### 8.2 Purification Protocol

For a detailed description regarding the assay procedure, please refer to the Instructions for Use/Method Sheets for **MagNA Pure 96 DNA and Viral NA Small Volume Kit**, **MagNA Pure 96 DNA and Viral NA Large Volume Kit** and **MagNA Pure 24 Total NA Isolation Kit**.

## 9. Quality Control

 Always run appropriate controls.

To control the entire process, starting from sample preparation to analysis, perform the following controls:

- Positive control**, using a sample material positive for target
- Negative control**, using PBS in place of the sample
- Extraction control**, using a sample material negative for target
- Internal control (IC)**, by adding a defined amount of a control template to all samples to be purified

For applications that could produce false negative results, such as the detection of pathogens, the use of an appropriate internal control (IC) is mandatory. The IC is added during nucleic acid isolation, preferably using the automated IC function of the MagNA Pure 96 System and the MagNA Pure 24 System. The IC can also be added manually to the sample. In this case, the IC must be stable in the sample material, and a nuclease-sensitive IC, such as unprotected RNA, should not be used for this purpose.

## 10. Limitations and Interferences

- Reliable results are dependent on appropriate specimen collection, transport, storage and processing procedures.
- Use of this product should be limited to personnel trained in techniques of nucleic acid purification and isolation and PCR.
- False negative results may occur if a specimen is improperly collected, transported or handled. False negative results may also occur if insufficient amount of template is present in the specimen.
- Any IVD application using the sample preparation procedure in conjunction with any downstream IVD nucleic acid testing should be evaluated with regard to the individual IVD parameter.
- To minimize the risk of a negative impact on the results, adequate controls should be used.
- Storage conditions (temperature, time) for lysates shall be validated with regard to the individual IVD parameter.

## 11. Supplementary Information

### 11.1 Symbols

In this Instruction Manual, the following symbols are used to highlight important information:

Symbol	Description
	Important Note
	Information Note
	For <i>in vitro</i> diagnostic use.
	The reagent complies with the requirements of the IVDR Regulation (EU) 2017/746.
	Catalogue Number
	Global Trade Item Number
	Unique Device Identifier
	Batch Code
	Use-by date
	Date of Manufacture
	Content of Kit
	Temperature Limit
	Consult instructions for use
	Distributed by
	Manufacturer
	Authorized representative in the European Community
	Importer

## 11.2 Changes to previous version

- Change of the name of the product.
- Update to comply with the requirements of the IVDR Regulation (EU) 2017/746.
- Add reference to MagNA Pure 24 System.

## 12. References

- 1 Sambrook J., Fritsch E.F., & Maniatis T. (1989). Molecular Cloning: A Laboratory Manual. Cold Spring Harbor Laboratory Press, Cold Spring Harbor, NY.
- 2 Farrell R.E. (1993). RNA Methodologies: A Laboratory Guide for Isolation and Characterization, Academic Press, San Diego

## 13. Trademarks

MAGNA PURE is a trademark of Roche.

All other product names and trademarks are the property of their respective owners.

## 14. Regulatory Disclaimer

For *in vitro* diagnostic use.

Roche Molecular Systems, Inc.  
1080 US Highway 202 South  
Branchburg, NJ 08876 USA  
Manufactured in Germany

Roche Diagnostics GmbH  
Sandhofer Strasse 116,  
68305 Mannheim, Germany  
+49 621 759 0

Distributed in USA by Roche Diagnostics, Indianapolis, IN, USA  
US Customer Technical Support: 1-800-526-1247

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