MagNA Pure Bacterial Lysis Buffer

REF 06 374 921 001

Version 06
Content version: August 2021
↓ Store at +15 to +25°C

1. Intended Use

The MagNA Pure 96 and 24 Systems are automated nucleic acid purification systems consisting of the instrument, software, control unit (only for MagNA Pure 96 System), consumables, and reagents. The MagNA Pure 96 and 24 Systems are intended for use by professional users for the purification of nucleic acids from biological samples for *in vitro* diagnostic purposes.

The MagNA Pure Bacterial Lysis Buffer is for use with the MagNA Pure 96 and 24 Systems.

2. Explanation of the Reagent

The MagNA Pure Bacterial Lysis Buffer is designed for:

- Lysis of different/various sample materials, *e.g.*,
 - Urine
 - Bronchoalveolar lavage (BAL)
 - Sputum*
 - Cerebrospinal fluid (CSF)*
 - Swabs
 - Stool
 - Bacterial cultures
 - Blood
 - Plasma
 - Serum
- Stabilization of nucleic acids with lysates.
- Nucleic acid purification using the MagNA Pure 96 System and MagNA Pure 24 System.

⊗ * Sample materials not validated with MagNA Pure 24 System.

3. Reagent Principle/ Summary

The MagNA Pure Bacterial Lysis Buffer is used to lyse bacteria.

4. Reagents - Working Solutions

Vial/ Cap	Label	Contents/Function
1 (colorless)	MagNA Pure Bacterial Lysis Buffer	<ul style="list-style-type: none"> ▪ one bottle, 20 ml ▪ ready-to-use reagent for lysis of bacteria ▪ < 5 % Laurylsarcosine, ▪ < 200 mM Tris-HCl, ▪ < 1 M Na-citrate

5. Precautions and Warnings

5.1 Precautions

⚠ The MagNA Pure Bacterial Lysis Buffer contains no hazardous substances in concentrations to be declared. Nevertheless, always wear gloves and follow standard safety precautions to minimize contact when handling.

⚠ Do not let this buffer touch your skin, eyes or mucous membranes. If contact does occur wash the affected area immediately with large amounts of water. If you spill the reagent, dilute the spill with water before wiping it up.

⚠ Handle all samples as if potentially infectious. Follow standard laboratory safety procedures for hazardous materials. Since the infectivity and titer of potential pathogens in the sample material can vary, always use effective pathogen inactivation methods and take the appropriate measures prescribed by local safety regulations. Work in a dedicated safety cabinet until you are sure the pathogens are inactivated. Note that even a boiling step cannot guarantee 100% inactivation of all kinds of pathogens.

5.2 Handling Instructions

- Wear disposable gloves and change them frequently.
- Do not use the reagent after its expiration date has passed.

In addition, to minimize the risk of carryover contamination which may result in false positive results, follow the guidelines listed below:

- Perform sample preparation, PCR/RT-PCR setup and PCR/RT-PCR in separate locations.
 - Discard pipette tips in sealed containers to prevent airborne contamination.
- Nuclease contaminated reagents and reaction vessels will degrade template nucleic acids. Please follow these guidelines to minimize the risk of contamination:
- Avoid touching surfaces or materials that could cause nuclease carryover.
 - Clean and decontaminate work areas and instruments, including pipettes, with commercially available decontamination reagents.
 - Use only new nuclease-free aerosol-blocking pipette tips and microcentrifuge tubes.
 - Use a work area specifically designed for RNA work, and if possible use reaction vessels and pipettors dedicated only for work with template RNA.

5.3 Laboratory Procedures

- All human sourced material, and all resulting waste, should be considered potentially infectious. Thoroughly clean and disinfect all work surfaces with disinfectants recommended by the local authorities.
- As the sensitivity and titer of potential pathogens in the sample material can vary, the operator must optimize pathogen inactivation, and follow the appropriate measures according to local safety regulations.
- Do not eat, drink or smoke in the laboratory work area.
- Do not pipette by mouth.
- Wear protective disposable gloves, laboratory coats, and eye protection when handling specimens and reagents.
- Avoid microbial and nuclease contamination of reagents when removing aliquots from reagent bottles. Use sterile disposable pipette tips.
- Wash hands thoroughly after handling specimens and reagent.

5.4 Waste Handling

- For US: Material Safety Data Sheets (MSDS) are available on www.usdiagnostics.roche.com, or upon request from the local Roche office.
- For all other countries: Material Safety Data Sheets (MSDS) are available online at www.e-labdoc.roche.com, or upon request from the local Roche office.
- Dispose unused reagents and waste in accordance with country, federal, state, and local regulations.

6. Storage and Stability (Reagents)

- Store the MagNA Pure Bacterial Lysis Buffer at +15°C to +25°C.
- When properly handled, the MagNA Pure Bacterial Lysis Buffer is stable until the expiration date printed on the label.

7. Materials

7.1 Materials Provided

see Reagents - Working Solution

7.2 Materials and Devices Required but not Provided

- Standard laboratory equipment
- Vortex mixer
- Thermomixer
- Centrifuge
- MagNA Pure 96 DNA and Viral NA Large Volume Kit (Cat. No. 06 374 891 001)
- MagNA Pure 96 DNA and Viral NA Small Volume Kit (Cat. No. 06 543 588 001)
- MagNA Pure 96 Instrument (Cat. No. 06 541 089 001)
- MagNA Pure 96 System Fluid (Internal) (Cat. No. 06 430 112 001)
- MagNA Pure 96 System Fluid (External) (Cat. No. 06 640 729 001)
- MagNA Pure 24 Total NA Isolation Kit (Cat. No. 07 658 036 001)
- MagNA Pure 24 Instrument (Cat. No. 07 290 519 001)
- Proteinase K, PCR grade, Activity (+37°C) ≥ 0.6 U/μl (e.g., Cat. No. 03 115 828 001)

8. Assay Procedures

8.1 General Remarks

- It is the user's responsibility to validate system performance for any procedures used in their laboratory.
- Due to the great variety of sample materials, no single universally applicable procedure is possible.
- Different sample materials may need adapted procedures, such as additional lysis steps to increase the DNA recovery from some bacterial species. These procedures shall be validated with regard to the individual IVD parameter.
- The preparation of a semi-liquid sample (BAL, sputum, stool, etc.) for nucleic acid isolation will depend on the type of sample material, sample viscosity, and particle type and content.

8.2 Purification Protocol

For a detailed description regarding the assay procedure, please refer to the Instructions for Use/Method Sheets for **MagNA Pure 96 DNA and Viral NA Small Volume Kit**, **MagNA Pure 96 DNA and Viral NA Large Volume Kit** and **MagNA Pure 24 Total NA Isolation Kit**.

9. Quality Control

⚠ Always run appropriate controls.

To control the entire process, starting from sample preparation to analysis, perform the following controls:

- Positive control**, using a sample material positive for target.
- Negative control**, using PBS in place of the sample.
- Extraction control**, using a sample material negative for target.
- Internal control (IC)**, by adding a defined amount of a control template to all samples to be purified.

⚠ For applications that could produce false negative results, such as the detection of pathogens, the use of an appropriate internal control (IC) is mandatory. The IC is added during nucleic acid isolation, preferably using the automated IC function of the MagNA Pure 96 System and MagNA Pure 24 System. The IC can also be added manually to the sample. In this case, the IC must be stable in the sample material, and a nuclease-sensitive IC, such as unprotected RNA, should not be used for this purpose.

10. Limitations and Interferences

- Reliable results are dependent on appropriate specimen collection, transport, storage and processing procedures.
- Use of this product should be limited to personnel trained in techniques of nucleic acid purification and isolation, and PCR.
- False negative results may occur if a specimen is improperly collected, transported or handled. False negative results may also occur if insufficient amount of template is present in the specimen.
- Any IVD application using the sample preparation procedure in conjunction with any downstream IVD nucleic acid testing should be evaluated with regard to the individual IVD parameter.
- To minimize the risk of a negative impact on the results, adequate controls should be used.
- Storage conditions (temperature, time) for lysates shall be validated with regard to the individual IVD parameter.

11. Supplementary Information

11.1 Symbols

In this Instruction Manual, the following symbols are used to highlight important information:

Symbol	Description
⚠	Important Note
ℹ	Information Note
IVD	For <i>in vitro</i> diagnostic use.
CE	The reagent complies with the requirements of the IVDR Regulation (EU) 2017/746.
REF	Catalogue Number
GTIN	Global Trade Item Number
UDI	Unique Device Identifier
LOT	Batch Code
🕒	Use-by date
🏭	Date of Manufacture
CONTENT	Content of Kit
🌡	Temperature Limit
📖	Consult instructions for use
D	Distributed by
🏭	Manufacturer
EC REP	Authorized representative in the European Community
🌐	Importer

11.2 Changes to previous version

- Change of the name of the product.
- Update to comply with the requirements of the IVDR Regulation (EU) 2017/746.
- Add reference to MagNA Pure 24 System.

12. References

- 1 Sambrook J., Fritsch E.F., & Maniatis T. (1989). Molecular Cloning: A Laboratory Manual. Cold Spring Harbor Laboratory Press, Cold Spring Harbor, NY.
- 2 Farrell R.E. (1993). RNA Methodologies: A Laboratory Guide for Isolation and Characterization, Academic Press, San Diego

13. Trademarks

MAGNA PURE is a trademark of Roche.

All other product names and trademarks are the property of their respective owners.

14. Regulatory Disclaimer

For *in vitro* diagnostic use.

Roche Molecular Systems, Inc.
1080 US Highway 202 South
Branchburg, NJ 08876 USA
Manufactured in Germany

Roche Diagnostics GmbH
Sandhofer Strasse 116,
68305 Mannheim, Germany
+49 621 759 0

Distributed in USA by Roche Diagnostics, Indianapolis, IN, USA
US Customer Technical Support: 1-800-526-1247

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- © * Provmaterialen är inte validerade med MagNA Pure 24 System.