Twist DX Probe & Primerの選定方法

2018.05.07

TwistAmp® Amp Probe作成

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Flip Reverse Primer: ☐

About In-Silico PCR

In-Silico PCR searches a sequence database with a pair of PCR primers, using an indexing strategy for fast performance. See an example [video](#) on our YouTube channel.

Configuration Options

Genome and Assembly - The sequence database to search.
Target - If available, choose to query transcribed sequences.
Forward Primer - Must be at least 15 bases in length.
Reverse Primer - On the opposite strand from the forward primer. Minimum length of 15 bases.
Max Product Size - Maximum size of amplified region.
Min Perfect Match - Number of bases that match exactly on 3' end of primers. Minimum match size is 15.
Min Good Match - Number of bases on 3' end of primers where at least 2 out of 3 bases match.
Flip Reverse Primer - Invert the sequence order of the reverse primer and complement it.

Output

When successful, the search returns a sequence output file in fasta format containing all sequence in the database that lie between and include the primer pair. The fasta header describes the region in the database and the primers. The fasta body is capitalized in areas where the primer sequence matches the database sequence and in lower-case elsewhere. Here is an example from human:

```
>chr22:31000551+31001000 TAACAGATTGATGATGCATGAAATGGG CCCATGAGTGGCTCCTAAAGCAGCTGC
TtACAGATTGATGATGCATGAAATGGGgggtggccagggggtgggggtga
gactgcagagaaagcagggctggttcataacaagctttgtgctcccaa
tatgacagctgaagttttccaggggctgatggtgagccagtgagggtaa
tacacagaacatcctagagaaacccctcattccttaaagattaaaaataa
gacttgctgctgtgaagggttggtattcctatttgagaaattctgtta
tccagaatggcttaccccaaatgctgaaagtgtgtacgtaattctcaa
agcaagctcctcctcagacagagaacaccagcggtcacaggaagcaaa
aatttgcttcacttttaaggtgaatccagaacccagatgtcagagctcc
aagcactttgctcagctccacGCAGCTGCTTTAGGAGCCACTCATGaG
```

The + between the coordinates in the fasta header indicates this is on the positive strand.

Author

In-Silico PCR was written by [Jim Kent](#). Interactive use on this web server is free to all. Sources and executables to run batch jobs on your own server are available free for academic, personal, and non-profit purposes. Non-exclusive commercial licenses are also available. Contact Jim for details.

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
Primer Melting Temperatures

Forward: 81.1 C cctcaccgagcgcggctacagcttcaccac
Reverse: 78.8 C ttgctcgaagtcaggcgacgtagcacag
The temperature calculations are done assuming 50 mM salt and 50 nM annealing oligo concentration. The code to calculate the melting temp comes from [Primer3](#).

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
UCSC In-Silico PCR

Genome: Assembly: Target:

Max Product Size: Min Perfect Match: Min Good Match:

Forward Primer: Reverse Primer:

Flip Reverse Primer: ☐

UCSC In-Silico PCR

>[chr7:5568024-5568138](#) 115bp CCTCACCGAGCGCGGCTACAGCTTCACCAC AGCCGTGGCCATCTCTTGCTCGAAGTCCAG
CCTCACCGAGCGCGGCTACAGCTTCACCACcaccggccgagcgggaaatcg
tgcgtagacattaaggagaagctgtgctacgtcgccCTGGACTTCGAGCAA
GAGATGGCCACGGCT

Primer Melting Temperatures

Forward: 81.1 C cctcaccgagcgcggctacagcttcaccac
Reverse: 78.9 C agccgtggccatctcttgctcgaagtccag
The temperature calculations are done assuming 50 mM salt and 50 nM annealing oligo concentration. The code to calculate the melting temp comes from [Primer3](#).

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