## Software Release notes: Pippin version 6.20

For previous software release notes, contact us at [support@sagescience.com](mailto:support@sagescience.com)

### Should you upgrade to v. 6.20?

PacBio / BluePippin users or users who are collecting large DNA should upgrade their instruments. This software version turns off the LED detectors when they are not required. This eliminates concerns over possible DNA nicking with prolonged exposure to blue light. It also includes a new high-pass protocol which collects fragments over 40kb.

Pippin Prep users may want to forgo an upgrade unless there is a concern about exposure to green light.

### Note ethidium bromide users:

Ethidium bromide cassettes will have LEDs turned off in the non-marker lanes. This will turn off the visualization of the DNA in the sample lanes, which will be different than expected. If users wish to monitor the DNA in sample lanes, users will need to manually turn on the LEDs in the protocol editor screen and save or re-save the protocol.

### What's new from v.6.13 to v.6.20?

#### Changes in 6.20:

The LED detectors turn will now turn off after detecting calibration markers. When internal standard protocol are used, the LEDs will turn off after the markers have been detected. With external standard protocols, the sample lane LEDs will be off by default.

#### Pippin Prep™

##### New Cassette Definitions or Improvements

None.

#### BluePippin™

##### New Cassette Definitions or Improvements

#### 40kb High-Pass Definition (“0.75%DF Marker U1 High Pass 40”):

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