

Sest Gene® Restriction Enzyme Hpa I



Cat.# FG-Hpal 5

Size 500 units Conc. 4 units/µl

Store at -20°C

Supplied with: 10X FastGene® Buffer IV (FG-REB4) 10X FastGene® FastCut Buffer (FG-REBHF) 6X DNA Loading Buffer Sterile water

Recognition site



For Research Use Only. Not for use in diagnostic procedures.

Heat Inactivation

No.

Methylation sensitivity

dam methylation: Not sensitive *dcm* methylation: Not sensitive CpG methylation: Conditionally sensitive

Prolonged incubation

A minimum amount of enzyme required to digest 1 μg substrate DNA for 16 hr; 0.25 U.

Relative activity in FastGene® Buffers

FastGene®	Buffer I:	0%
FastGene®	Buffer II:	50%
FastGene®	Buffer III:	25%
FastGene®	Buffer IV:	100%
FastGene®	FastCut Buffer:	100%

Note

Cleavage of mammalian genomic DNA is blocked by CpG methylation that partially overlaps its recognition sequence. Reaction condition with excess enzyme, excess glycerol (>5%) or high pH (>8.0) may result in star activity. It is relatively unstable during incubation. Thus, use a higher level of enzyme in order to achieve complete cleavage in 1 hr. It is sensitive to impurities in DNA and its activity can be reduced by desiccation in a frost free freezer. Source: Haemophilus parainfluenzae

Reaction conditions

1X FastGene[®] Buffer IV 37℃ 1X FastGene[®] FastCut Buffer, 37℃

FastGene[®] FastCut Buffer

 $\mathsf{FastGene}^{\circledast}$ restriction enzyme can cut substrate DNA in 5-15 min with $\mathsf{FastGene}^{\circledast}$ $\mathsf{FastCut}$ Buffer.

1X FastGene® Buffer IV

20 mM Tris-acetate (pH 7.9 at 25°C) 50 mM potassium acetate 10 mM magnesium acetate 100 μg/ml BSA

Unit definition

One unit is defined as the amount of enzyme required for complete digestion of 1 μg bacteriophage λ at 37°C for 1 hr in 50 μl reaction mixtures.

Quality control

- Unit definition assay
- Overdigestion assay
- Endonuclease assay
- Extreme pure assay

Dilution buffer

FastGene® Diluent A

Standard reaction condition

- Normal	protocol
----------	----------

Component	Final Conc.	Volume
Substrate DNA	1 µg	Χ μΙ
10X FastGene [®] Buffer IV	1 X	5 µl
Hpa I	4 unit	1 µl
Sterile water		up to 50 µl
→ Incubate at 37°C for 1 hr		

Fast protocol

Component	Final Conc.	Volume
Substrate DNA	1 µg	Xμl
10X FastGene [®] FastCut Buffer	1 X	5 µl
Hpa I	4 unit	1 µl
Sterile water		up to 50 µl
Incubate at 27°C for 1E mir		

 \rightarrow Incubate at 37°C for 15 min

% We recommend 5-10 units of enzyme per μg DNA and 10-20 units for genomic DNA in a 1 h digest.