

# Restriction Enzyme Fok I



Cat.# FG-Fokl Size 1,000 units Conc. 5 units/µl

Store at -20℃

**Supplied with:** 10X FastGene® Buffer IV (FG-REB4)

10X FastGene® FastCut Buffer (FG-REBHF) 6X DNA Loading Buffer

Sterile water

# Recognition site

For Research Use Only. Not for use in diagnostic procedures.

**ISO**9001

## Dilution buffer:

FastGene® Diluent A

# **Heat Inactivation**

Fok I can be inactivated at 65°C for 20 min.

# Methylation sensitivity

dam methylation: Not sensitive

dcm methylation: Conditionally sensitive

CpG methylation: Conditionally sensitive

# Prolonged incubation

A minimum amount of enzyme required to digest 1  $\mu g$  substrate DNA for 16 hr; 1 U.

# Relative activity in FastGene® Buffers

 FastGene® Buffer I:
 100%

 FastGene® Buffer II:
 100%

 FastGene® Buffer III:
 10%

 FastGene® Buffer III:
 100%

 FastGene® FastCut Buffer:
 100%

#### Note

Cleavage is inhibited by *dcm* methylation and CpG methylation partially overlapping its cleavage site.

Overdigestions of >5 units of Fok I per µg of DNA and incubation times >2 hours are not recommended.

Source: Flavobacterium okeanokoites

# Reaction conditions

1X FastGene® Buffer IV 37°C 1X FastGene® FastCut Buffer, 37°C

# FastGene® FastCut Buffer

FastGene® restriction enzyme can cut substrate DNA in 5-15 with FastGene® FastCut Buffer.

# 1X FastGene® Buffer IV

20 mM Tris-acetate (pH 7.9 at 25°C) 50 mM potassium acetate 10 mM magnesium acetate 100 μα/ml BSA

## Unit definition

One unit is defined as the amount of enzyme required for complete digestion of 1  $\mu$ g bacteriophage  $\lambda$  at 37°C for 1 hr in 50  $\mu$ l reaction mixtures.

# Quality control

- Unit definition assay
- Overdigestion assay
- Endonuclease assay
- Extreme pure assay

## Standard reaction condition

- Normal protocol

Component	Final Conc.	Volume
Substrate DNA	1 μg	Χ μΙ
10X FastGene® Buffer IV	1 X	5 μΙ
Fok I	5 unit	1 μΙ
Sterile water		up to 50 μl

- → Incubate at 37°C for 1 hr
- Fast protocol

Component	Final Conc.	Volume
Substrate DNA	1 μg	Χ μΙ
10X FastGene® FastCut Buffer	1 X	5 μΙ
Fok I	5 unit	1 μΙ
Sterile water		up to 50 μl
In a class at 27%C for 15 main	_	

→ Incubate at 37°C for 15 min

 $\ensuremath{\mathbb{X}}$  We recommend 5-10 units of enzyme per  $\mu g$  DNA and 10-20 units for genomic DNA in a 1 h digest.