

# Restriction Enzyme Afl II



Cat.# FG-AfIII

Size 2,000 units

Conc. 4 units/µl

Store at -20°C

Supplied with: 10X FastGene® Buffer IV (FG-REB4)

10X FastGene® FastCut Buffer (FG-REBHF)

6X DNA Loading Buffer

Sterile water

#### Recognition site

For Research Use Only. Not for use in diagnostic procedures.

**ISO**9001

## Dilution buffer:

FastGene® Diluent A

#### **Heat Inactivation**

Afl II can be inactivated at 65°C for 20 min.

## Methylation sensitivity

dam methylation: Not sensitive dcm methylation: Not sensitive CpG methylation: Not sensitive

#### Prolonged incubation

A minimum amount of enzyme required to digest 1  $\mu g$  substrate DNA for 16 hr; 0.13 U.

## Relative activity in FastGene® Buffers

 FastGene® Buffer I:
 75%

 FastGene® Buffer II:
 100%

 FastGene® Buffer III:
 75%

 FastGene® Buffer IV:
 100%

 FastGene® FastCut Buffer:
 100%

#### Note

It is not sensitive to *dam, dcm* or CpG methylation. Less than 50% of cleaved fragments are ligated in 20  $\mu$ l reaction mixtures containing 100-500 units of T4 DNA ligase.

Source: Anabaena flos-aquae

## Reaction conditions

1X FastGene® IV 37℃

1X FastGene® FastCut Buffer, 37°C

### FastGene® FastCut Buffer

FastGene® restriction enzyme can cut substrate DNA in 5-15 with FastGene® FastCut Buffer.

#### 1X FastGene® Buffer IV

20 mM Tris-acetate (pH 7.9 at 25°C) 50 mM potassium acetate 10 mM magnesium acetate 100 μg/ml BSA

#### Unit definition

One unit is defined as the amount of enzyme required for complete digestion of 1  $\mu$ g bacteriophage  $\lambda$  (Hind III digestion) at 37°C for 1 hr in 50  $\mu$ l reaction mixtures.

## Quality control

- Unit definition assay
- Overdigestion assay
- Endonuclease assay
   Extreme pure assay

#### Standard reaction condition

- Normal protocol

- Normai protocoi		
Component	Final Conc.	Volume
Substrate DNA	1 μg	Χ μΙ
10X FastGene® Buffer IV	1 X	5 μΙ
Afl II	4 unit	1 μΙ
Sterile water	up to 50 μl	

→ Incubate at 37°C for 1 hr

- Fast protocol

Component	Final Conc.	Volume
Substrate DNA	1 μg	ΧμΙ
10X FastGene® FastCut Buffer	1 X	5 µl
Afl II	4 unit	1 μΙ
Sterile water		up to 50 μl

→ Incubate at 37°C for 15 min