

Restriction Enzyme Hinc II



Cat.# FG-HincII Size 1,000 units Conc. 10 units/µl

Store at -20℃

Supplied with: 10X FastGene® Buffer IV (FG-REB4)

10X FastGene® FastCut Buffer (FG-REBHF) 6X DNA Loading Buffer

Sterile water

Recognition site



For Research Use Only. Not for use in diagnostic procedures.

ISO9001

Dilution buffer

FastGene® Diluent B

Heat Inactivation

Hinc II can be inactivated at 65°C for 20 min.

Methylation sensitivity

dam methylation: Not sensitive dcm methylation: Not sensitive CpG methylation: Conditionally sensitive

Prolonged incubation

A minimum amount of enzyme required to digest 1 μ g substrate DNA for 16 hr; 0.13 U.

Relative activity in FastGene® Buffers

 FastGene® Buffer I:
 75%

 FastGene® Buffer II:
 50%

 FastGene® Buffer III:
 50%

 FastGene® Buffer IV:
 100%

 FastGene® FastCut Buffer:
 100%

Note

Cleavage of mammalian genomic DNA is blocked by CpG methylation partially overlapping its recognition sequence. It is sensitive to impurities present in quick DNA preparations.

Source: Haemophilus influenzae Rc

Reaction conditions

1X FastGene® Buffer IV 37°C 1X FastGene® FastCut Buffer, 37°C

FastGene® FastCut Buffer

FastGene® restriction enzyme can cut substrate DNA in 5-15 min with FastGene® FastCut Buffer.

1X FastGene® Buffer IV

20 mM Tris-acetate (pH 7.9 at 25°C) 50 mM potassium acetate 10 mM magnesium acetate 100 μg/ml BSA

Unit definition

One unit is defined as the amount of enzyme required for complete digestion of 1 μ g bacteriophage λ at 37°C for 1 hr in 50 μ l reaction mixtures.

Quality control

- Unit definition assay
- Overdigestion assay
- Endonuclease assay - Extreme pure assay

Standard reaction condition

- Normal protocol

Component	Final Conc.	Volume
Substrate DNA	1 μg	Xμl
10X FastGene® Buffer IV	1 X	5 µl
Hinc II	10 unit	1 µl
Sterile water		up to 50 μl

→ Incubate at 37°C for 1 hr

- Fast protocol

Component	Final Conc.	Volume
Substrate DNA	1 μg	Χ μΙ
10X FastGene® FastCut Buffer	1 X	5 μΙ
Hinc II	10 unit	1 μΙ
Sterile water		up to 50 μl
→ Incubate at 37°C for 15 mir	1	

* We recommend 5-10 units of enzyme per μg DNA and 10-

20 units for genomic DNA in a 1 h digest.