


I N S T R U C T I O N S



Velveteen Squares

Catalog No. H37848-0001 (White) & H37848-0002 (Black)



These 100% cotton, sterile Velveteen Replica Plating Squares are especially suited for use in conjunction with the Scienceware® Replica Plating Tool #F37848-0000.

The Velveteen Squares are supplied sterilized and ready for use. These may be re-used a number of times after careful washing and sterilization. To maintain sterility do not open the remaining bags until they are ready to be used. Seal tightly when not in use to reduce airborne contamination.

Bel-art Products

PEQUANNOCK, NJ 07440-1992 USA
TEL: 1-800-4BEL-ART • FAX: 973 694-7199 • www.belart.com

937848016


Scienceware® is a reg. TM of Bel-Art Products
Bel-Art Products assumes no obligation or liability for any advice furnished by it, or for results obtained with respect to these products. All such advice is given and accepted at the buyer's risk.

11/06

OTHER FINE SCIENCEWARE® PRODUCTS FOR YOUR LABORATORY

Replicate! Replicate! Replicate!
**The Bel-Blotter™ 96 Wells
Replicating Tool**

Catalog No. F37876-0002



- Efficient, Cost Effective, Precise
- Capillary Action Insures That Each Well Will Pick Up 10ml Per Tip
- 96 "Open Ended" Pipette Tips Mate Perfectly With Any 96 Well Plate
- Molded of Autoclaveable Polypropylene



Bel-art Products

TEL: 1-800-4BEL-ART • FAX: 973 694-7199 • www.belart.com