

Smart sample storage solutions

Microtube Storage Boxes

Designed for easy storage of 1.5/2.0 mL microtubes and cryovials using the latest plastics technology

Our polypropylene microtube storage boxes feature a clear alphanumeric well identifiers (except black versions), uniform wells and a sleek wave design finish that also indicates the boxes capacity on the outside.



81 well boxes have a keyed friction-fit lid, with a frosted labeling patch and will fit into standard 2" freezer racking. Suitable for temperatures down to -80 ° C.



50- and 100-well boxes feature a hinged, snap clasp for extra-secure sample storage.



Precision molding ensures all wells are uniform and can accommodate tubes



Choose opaque black for light-sensitive samples

A autoclavable



Heathrow Scientific, LLC
620 Lakeview Parkway, Vernon Hills IL 60061 USA
Tel.: +1-847-816-5070 | Fax: +1-847-816-5072
www.heathrowscientific.com | E-mail: info@heatsci.com

High transparency,
yet strong coloration

Microtube Storage Boxes

6 bright colors facilitate color coding



50-well boxes



81-well boxes



100-well boxes

81-Well Microtube Storage Boxes

- Measure 130 x 130 x 47 mm
- Friction-fit lid; pack of 5

part no.	tubes	wells	color
HS120034	1.5/2.0 mL	81	Natural
HS120035	1.5/2.0 mL	81	Blue
HS120036	1.5/2.0 mL	81	Green
HS120037	1.5/2.0 mL	81	Pink
HS120038	1.5/2.0 mL	81	Yellow
HS120039	1.5/2.0 mL	81	Orange
HS120040	1.5/2.0 mL	81	*
HS120041	1.5/2.0 mL	81	Black

50-Well Microtube Storage Boxes

- Measure 142 x 89 x 56 mm
- Snap clasp lid; pack of 5

part no.	tubes	wells	color
HS120032	1.5/2.0 mL	50	Natural
HS120033	1.5/2.0 mL	50	*

100-Well Microtube Storage Boxes

- Measure 141 x 151 x 57 mm
- Snap clasp lid; pack of 5

part no.	tubes	wells	color
HS120042	1.5/2.0 mL	100	Natural
HS120043	1.5/2.0 mL	100	Black
HS120044	1.5/2.0 mL	100	*

*Assorted colors: blue, green, purple, yellow, and orange.

Ask for these and other innovative Heathrow Scientific products:

True North® Cryogenic Vials

Innovative design saves freezer space!

Made from superior puncture- and rupture-resistance polypropylene with a silicon lid seal, our vials include a printed writing area and graduations. For mechanical freezers and vapor phase LN₂ (-196°); include a traceable sterility certificate.



**RNase/DNase -
Pyrogen-FREE!**

Microtubes

High-quality tubes for critical applications

Made from 100% virgin polypropylene and are certified sterile, RNase/DNase and pyrogen free. Autoclavable and can withstand high centrifugal forces.



A autoclavable

D disposable