

4titude 一次元バーコードオプション



1. 標準バーコード（プレート貼付 or バーコード単品）

・4titude社製フルスカートおよびセミスカート全製品対応可能

仕様

- ・標準のバーコード形式は*****（例：12345678）で、ほぼすべてのバーコードリーダーで読み取り可能なCode128を採用。
- ・ラベル右側にヒューマンリーダブルコード（人が目視で読める形式）を印字。常に確実なサンプルの認識が可能。
- ・摩擦に強く、温度耐性は-20℃～+100℃、DMSOなどの溶剤にも耐性のあるラベル素材。

ご注文方法

プレート1枚あたりの貼付枚数により、型番末尾に /SBC、/DBC、/TBC、/QBC いずれかを加えてご注文ください。バーコード追加料金が自動的に付加されます。

※プレート枚数が500枚未満の場合：バーコードセットアップ料が必要

Cat.No.	単位	価格
4ti-####/SBC	シングルバーコード印字 (1箇所) (右図A)	プレート1枚分
4ti-####/DBC	ダブルバーコード印字 (2箇所) (右図A⑧)	プレート1枚分
4ti-####/TBC	トリプルバーコード印字 (3箇所) (右図A⑧⑨)	プレート1枚分
4ti-####/QBC	クワッドラプルバーコード印字 (4箇所) (右図A⑧⑨⑩)	プレート1枚分

Cat.No.	単位	価格
4ti-007	バーコードセットアップ料	プレート500枚未満 お問い合わせください

バーコード単品

バーコードシール単品（ロール）の販売も行っております。

750枚までは一定価格、751枚以上は追加料金でご購入いただけます。

Cat.No.	単位	価格
4ti-004	バーコード単品（標準）	バーコードラベル750枚まで/ロール お問い合わせください

2. 特注バーコード（プレート貼付 or バーコード単品）

・標準仕様でカバーされていない形式（例：任意の番号、貼付位置など）をご希望の場合に特注対応を承ります。

ご注文方法

お好みのプレートの型番+バーコード追加用型番にてご注文ください。専用オーダーフォームがございます。

※プレート枚数が500枚未満の場合：バーコードセットアップ料が必要

Cat.No.	価格	Cat.No.	単位	価格
4ti-####	お問い合わせください	4ti-006	シングルバーコード印字 (1箇所)	プレート1枚分
		4ti-005	ダブルバーコード印字 (2箇所)	プレート1枚分
		4ti-008	トリプルバーコード印字 (3箇所)	プレート1枚分
		4ti-0011	クワッドラプルバーコード印字 (4箇所)	プレート1枚分

Cat.No.	単位	価格
4ti-007	バーコードセットアップ料	プレート500枚未満 お問い合わせください

バーコード単品

バーコードシール単品（ロール）のみでの販売も行っております。



750枚までは一定価格、751枚以上は追加料金でご購入いただけます。

Cat.No.	単位	価格
4ti-003	バーコード単品（特注）	バーコードラベル750枚まで/ロール お問い合わせください

お見積依頼およびご注文方法

専用オーダーフォームにご記入の上で注文ください。

	標準バーコードプレート	特注バーコードプレート
プレート型番	4ti-####/SBC、/DBC、/TBC、/QBC、1プレートあたりの貼付枚数により末尾のコードを変更してください。	4ti-####
バーコード追加用型番	不要です。上記型番でご注文いただく事により、バーコード料金は自動付加されます。	4ti-006、005、008 or 0011 1プレートあたりの貼付枚数により異なる
バーコードセットアップ料型番	4ti-007（プレート500枚未満）	4ti-007（プレート500枚未満）
専用オーダーフォーム	不要	必要



サイズ	最小	最大
A ラベル長	-	46mm
B コード長	-	25mm
C クワイエットゾーン	3mm	-
D ラベル幅	5mm	6mm
E コード幅	5mm	6mm

1. Linear Barcode Options


1.1 Standard Barcoded Plates and Barcodes, Supplied on a Roll

4titude® standard barcoded plates are available off-the-shelf, single to quadruple barcoded for all our fully-skirted and semi-skirted plates. The standard barcode format is ##### (e.g. 12345678), printed in Code 128 which is compatible with most barcode readers. In addition, human readable format is printed to the right of the label, ensuring samples can always be identified. The labels are highly scratch resistant and can withstand temperatures from -20°C to 100°C and solvents, such as DMSO.

The single barcode is positioned opposite the numbers (south side of the plate), the double barcode is additionally positioned opposite the letters (east side) and the triple barcode is additionally positioned north side of the plate.

To order either version, please add the following codes to the ordering code of the respective plate: /SBC (for single bar code), /DBC (for double bar code), /TBC (for triple bar code) and /QBC (for quadruple bar code). The price will include your current plate price plus an additional surcharge to cover the cost of the barcoding. Please refer to section 1.3 to select the corresponding product codes in our price list or [contact 4titude](#).

4titude® can also supply the same barcode labels, as detailed above, loose supplied on a roll. Our standard charge (4ti-004) applies for a quantity of up to 750 labels, above 750 labels, additional charges apply. Please [contact 4titude](#) for details.



Position of Standard Barcode Labels



	Dimension	Min.	Max.
A	Label length	-	46 mm
B	Code length	-	25 mm
C	Quiet zone	3 mm	-
D	Label width	5 mm	6 mm
E	Code width	5 mm	6 mm

Label Specification

1.2 Custom Barcoded Plates

Should you require other barcode options not covered by our standard barcoding service, e.g. specified barcodes or label positions, our custom service can apply them in your preferred configuration or format, with a starting sequence of your specification.

Pricing will depend on your custom barcoding requirements, please [contact 4titude](#). For any custom order of less than 500 plates/10 boxes, an additional set-up charge (4ti-007) applies.

4titude® can also supply custom linear barcode labels, as detailed above, loose supplied on a roll. Our standard charge (4ti-003) applies for a quantity of up to 750 labels, above 750 labels, additional charges apply. Please [contact 4titude](#) for details.

Please complete this **Custom Linear Barcode Request Form** (see section 2) and submit it with every order of custom barcoded plates or custom barcode labels.

1.3 Ordering Information

Ordering Information		
Code	Description	Quantity
4ti-006	Single Barcode applied to a plate	Surcharge for 1 plate
4ti-005	Double Barcode applied to a plate	Surcharge for 1 plate
4ti-008	Triple Barcode applied to a plate	Surcharge for 1 plate
4ti-0011	Quadruple Barcode applied to a plate	Surcharge for 1 plate
4ti-007	Barcoding set-up charge	For orders below 500 plates
4ti-004	Standard linear barcode labels, supplied on a roll	Charge for a qty of <750 barcode labels, for a qty of >750 labels, additional charges apply
4ti-003	Custom linear barcode labels, supplied on a roll	Charge for a qty of <750 barcode labels, for a qty of >750 labels, additional charges apply

Information required	Standard Barcoded Plates	Custom Barcoded Plates
Product code <i>plate</i>	4ti-####/SBC, /DBC, /TBC, /QBC, depending on the number of barcodes per plate	4ti-#### without addition
Product code <i>barcode surcharge</i>	Not required, surcharge will be added automatically	4ti-006, 005, 008 or 0011, depending on the number of barcodes per plate
Product code <i>barcoding set-up charge</i>	4ti-007, applies to orders below 500 plates	4ti-007, applies to orders below 500 plates
<i>Custom Linear Barcode Request Form</i>	Not required	Required , please refer to section 2 of this document

Please note

Barcode pricing is available on request. Please use this **Custom Linear Barcode Request Form** to submit your price request, please send the completed form to orders@4ti.co.uk, your account manager or your local distributor.

Incomplete or incorrect details may result in product delay. Please [contact 4titude](#) if you require assistance in completing this form. Please refer to the next page for the 4titude® **Custom Linear Barcode Request Form** (section 2).

2. Custom Linear Barcode Request Form

Use this form for an order of **Custom Linear Barcoded Plates** or **Custom Linear Barcode Labels** in conjunction with a standard purchase order or for a **Custom Linear Barcode Price Request**. 4titude® custom barcoded plates and barcode labels supplied on a roll will be manufactured to your specific requirements. Barcodes will be printed in Code 128, for other options, please [contact 4titude®](#). Use "Additional Information" (section 2.4.2) if specific details are required. Incomplete or incorrect details may result in product delay.

2.1 Kind of Request/Contact Details

I am Making a price request Submitting this form with an order, PO number _____

End user
Name
Company/Institute
Department
Street
City
State/Province/County
Country
Tel
Email

Distributor
Distributor
Contact Name
Country
Tel
Email
4titude® Account Number

2.2 Product Requirements *(Note: applies to Custom Linear Barcoded Plates, only)*

Part number	Description of item to be barcoded	Quantity per case	Number of cases required

2.3 Barcode Print Options

2.3.1 Sequence information

4titude® can produce linear barcodes using any combination of numbers or letters up to a maximum of 10 digits (e.g. 12345678, ABC1234567, 76543ACB12).

I wish to use the **4titude® standard code**, and accept the next available sequence of the format ##### (e.g. 12345678).

I wish to use a **custom code**. Please specify sequence **START** code:

Please specify sequence **END** code:

Note: We recommend to include only necessary information in the barcode and leave out fillers. Unnecessary information will increase the complexity of the code which in turn may reduce its readability.

2.3.2 Number of barcodes per plate *(Note: applies to Custom Linear Barcoded Plates, only)*

I require **one barcode per plate**.

I require **multiple barcodes per plate**, all barcodes on the plate to have **the same** sequence.

barcodes on the plate to **differ in** sequence.

Please [contact 4titude®](#) to discuss options, additional charges apply for the application of varying sequences per plate.

2.3.3 Position human readable code

With the standard barcode format, the human readable code is printed to the right in reference to the barcode. Alternatively, it can be printed to the left. Please specify the position of the human readable code in reference to the barcode.

Human readable code **right** (standard)



Human readable code **left**



2.4 Barcode Positioning *(Note: applies to Custom Linear Barcoded Plates, only)*

The 4titude® standard barcode label is applied opposite the numeric grid reference (Opposite Numeric), positioned centrally in horizontal and vertical direction.



2.4.1 Position on the plate

I wish to use the 4titude® standard position (Opposite Numeric).

I wish to use a **custom position on the plate**.


Please check **one box** below to indicate the position required. For Double/Triple Bar Codes please check **two/three boxes**.

N (Numeric)

A (Alpha)

OA (Opposite Alpha)

ON (Opposite Numeric)



2.4.2 Additional information/Special requirements

Please specify

2.5 Confirmation of Requirements

2.5.1 Price request

Please send the completed form to orders@4ti.co.uk, your account manager or your local distributor to receive your estimate. Please save this form for your future order.

2.5.2 Submitting this form with an order

I agree that all the information presented in this form is accurate and complete, I herewith authorise production.

Full Name _____

Date _____

Please send the completed form **together with your purchase order** to orders@4ti.co.uk, your account manager or your local distributor.