

Nanobind[®] HT 1 mL blood kit RT

Product PN: 102-763-100

Description


The Nanobind HT 1 mL blood kit RT contains reagents for components of the Nanobind HT 1 mL blood kit (102-762-800) that are to be stored at room temperature. This kit is used in conjunction with Nanobind HT 1 mL blood kit 4C (102-763-200) to complete automated HMW DNA extraction from 1 mL blood (96 reactions).

Kit configuration & layout

Configuration

#	Component	Part number	Qty	Color	Volume
1	Proteinase K	102-782-400	1	clear	12 mL
2	Nanobind disks	102-782-900	1	clear	96 each
3	Buffer EB	102-782-800	2	clear	30 mL
4	Buffer BL3	102-782-500	1	clear	100 mL
5	Buffer CW1	102-782-600	1	clear	100 mL
6	Buffer CW2	102-782-700	2	clear	100 mL

Layout



Quantities

96 reactions

Storage and handling requirements

Product should be stored ambient between 15 to 30 °C.

Refer to protocols for additional recommendations for reagent handling.

Shelf life

Refer to expiration date on kit label.

Intended use

For research use only. Not for use in diagnostic procedures.

Support

For additional product and technical information, such as safety data sheets (SDS) or procedures, visit our [document website](#).

Legal

Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable license terms at pacb.com/license. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, and SBB are trademarks of PacBio. All other trademarks are the sole property of their respective owners.