

Cryovial Rack

(#61050)

1. Description

SPL Cryovial Racks are designed for handling and storage of Cryovials.

2. Features

- ❖ Numeric labeling
- ❖ 3 colors: Blue, Green, Red
- ❖ Autoclavable


3. Product Range

Cat. No.	Material	Description	Packaging
61050	PP	5 x 10 (50 Holes)	10

4. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

5. Product Photographs



SPL-TDS-CRYRA


For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
 © SPL Life Sciences Co., Ltd., Version 3, November 20th, 2018

6. Specifications & Technical Drawing

- Specifications

Dimensions		Material (Resin)	Packaging
Length:	193.00 mm	Polypropylene (PP)	• Pieces/Case: 10
Width:	85.93 mm		
Height:	27.55 mm		

- Technical Drawing



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr

To place an order, contact your local distributor or

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr

SPL-TDS-CRYRA

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd., Version 3, November 20th, 2018