



Henke-Sass, Wolf GmbH
Keltenstrasse 1 · 78532 Tuttlingen · Germany

Henke-Sass, Wolf GmbH
Keltenstrasse 1
78532 Tuttlingen · Germany
Telefon: +49 7462 9466-137
Telefax: +49 7462 9466-5137
heiko.schwoerer@henkesasswolf.de
www.henkesasswolf.de

To whom it may concern

Tuttlingen, 05.09.2017

Inner diameter / Outer diameter - disposable syringes

Item / Artikel:	Inner diameter (mm):	Outer diameter (mm):
HSW NORM-JECT® 1 ML, LDS, Luer Slip	4,7	6,39
HSW NORM-JECT® 2 ML / 3 ML, Luer Slip and Luer Lock	9,7	10,8
HSW NORM-JECT® 5 ML (6 ML), Luer Slip and Luer Lock	12,5	13,7
HSW NORM-JECT® 10 ML (12 ML), Luer Slip and Luer Lock	15,9	17,3
HSW NORM-JECT® 20 ML (24 ML), Luer Slip and Luer Lock	20,1	21,6
HSW NORM-JECT® 30 ML, Luer Slip and Luer Lock	22,6	24,1
HSW NORM-JECT® 50 ML (60 ML), Luer Slip and Luer Lock	29,0	30,9

HSW SOFT-JECT® 1 ML, Luer Slip	4,7	6,39
HSW SOFT-JECT® 1 ML, LDS, Luer Slip and Luer Lock	4,7	6,4
HSW SOFT-JECT® 2 ML / 3 ML, Luer Slip and Luer Lock	9,7	10,8
HSW SOFT-JECT® 5 ML, Luer Slip and Luer Lock	12,5	13,7
HSW SOFT-JECT® 10 ML (12 ML), Luer Slip and Luer Lock	15,9	17,3
HSW SOFT-JECT® 20 ML, Luer Slip and Luer Lock	20,1	21,6
HSW SOFT-JECT® 30 ML, Luer Slip and Luer Lock	22,0	24,0
HSW SOFT-JECT® 50 ML (60 ML), Luer Slip	27,75	30,5
HSW SOFT-JECT® 50 ML (60 ML), Luer Lock (BD)	26,6	29,1
HSW SOFT-JECT® 50 ML (60 ML), Catheter	27,8	30,45

NON-CONTROLLED COPY

Sitz der Gesellschaft: Tuttlingen
Amtsgericht Stuttgart · HRB 450056
Geschäftsführung:
Kathrin McKenna, Nina Stackmann,
Peter Decker, Armin Lekitsch
Steuer-Nr. 21101/33306
USt.-ID-Nr. DE142927433

Deutsche Bank Tuttlingen
BLZ 653 700 75
Kto. 211494000
IBAN: DE89 6537 0075 02114940 00
BIC: DEUTDESS653

Kreissparkasse Tuttlingen
BLZ 643 500 70
Kto. 8509297
IBAN: DE57 6435 0070 0008 5092 97
BIC: SOLADES1TUT

BW-Bank Stuttgart
BLZ 600 501 01
Kto. 4879048
IBAN: DE06 6005 0101 0004 8790 48
BIC: SOLADEST

Commerzbank Tuttlingen
BLZ 694 440 007
Kto. 240260000
IBAN: DE15 6944 0007 0240 2600 00
BIC: COBADEFFXXX