

## MD-MINI Blocks Information

Cat. No.	Block Picture	Well Shape	Description	Well diameter	Well depth	Block dimension
MD-MINI-B01			For 0.2ml tube (PCR Strip tube), 32 wells	6.35mm	17mm	L71 x W47 x H32mm
MD-MINI-B02			For 1.5ml tube, 12 wells	10.88mm	30mm	L71 x W47 x H32mm
MD-MINI-B03			For 15ml centrifuge tube, 6 wells	17.3 mm	70 mm	L71 x W47 x H75mm
MD-MINI-B04			For 50ml tube, 2 wells	29 mm	70 mm	L71 x W47 x H75mm
MD-MINI-B05			For 0.5ml tube, 12 wells	7.9 mm	25 mm	L71 x W47 x H32mm
MD-MINI-B06			For 1.5 or 2.0ml tube, 12 wells	11 mm	30 mm	L71 x W47 x H32mm
MD-MINI-B07			For 1.5 or 2.0ml tube, 12 wells	11 mm	30 mm	L71 x W47 x H32mm