





	96 Well μltra Plate	96 Well Fast Plate	Semi- Skirt96 Raised Rim Plate	Skirted96 Plate	96-Well MB Plate	96-Well 480 Plate	384-Well NX Plate	384-Well 480 Plate
代表型番	21970	38800	35800	23080	28110	38870	39620	38820
<b>PCR 装置</b>								
<b>Eppendorf</b>								
MasterCycler	●	●	●	●	●			
MasterCycler ep	●	●		●	●			
MasterCycler M384				●			●	
PTC-100	●	●	●	●	●		●	
PTC-200	●	●	●		●		●	
PCT- 225 Tetrad	●	●	●	●	●		●	
Dynad/Disciple	●	●	●	●	●		●	
RoboCycler	●	●						
TouchGene Gradient (TC512)	●	●	●	●	●		●	
G enius (TC412)	●	●	●	●				
TouchGene X	●	●	●	●			●	
<b>Flexi</b>								
Gene	●		●	●	●		●	
Genius	●		●	●	●			
<b>ThermoHybaid</b>								
PCR Sprint			●	●				
PCR Express	●	●	●	●	●		●	
Multi Block System	●	●	●	●	●		●	
Touchdown	●	●	●	●	●		●	
Omnigene	●	●	●	●	●		●	
Omn-E	●	●	●	●	●			
pxe	●	●	●	●	●		●	
px2	●	●	●	●	●		●	
<b>TaKaRa</b>								
TP240				●				
TP3000	●			●	●			
<b>シーケンサー</b>								
<b>Applied Biosystems (Life Technologies)</b>								
310	●		●					
3100	●		●	●	●		●	
3130	●			●	●		●	
3700	●			●	●		●	
3730/3730x	●			●	●		●	
<b>Amersham</b>								
Megabace 500				●	●			
Megabace 1000					●			
MegaBace 4000							●	
<b>MJ Research</b>								
BaseStation				●				
<b>Transgenomics</b>								
Wave System				●				