

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価（税抜）	2024年4月1日より 新定価（税抜）	備考
10035	Λ°トリエ イソシ 35x10mm,9.4cm2 ST 20x25	¥15,100	¥16,400	
10040	1000ul マルチフィットチップ° BL 1000/bag	¥4,300	¥4,600	
10060	Λ°トリエ イソシ 60x15mm,21.50cm2 ST20x25	¥16,900	¥18,300	
10080	1000ul マルチフィットチップ° BL 100X10R	¥9,200	¥9,900	
10090	Λ°トリエ イソシ 90x15mm,57.50cm2 ST10x50	¥17,900	¥19,400	
10091	Λ°トリエ イソシ 90x20mm,58cm2 ST 10x20	¥11,400	¥12,400	
10095	Λ°トリエ イソシ 90x15mm 57.50cm2 ST10x20	¥9,000	¥9,800	
10100	Λ°トリエ イソシ 100x15mm,72.30cm2ST10x50	¥25,900	¥28,000	
10130	1000ul マルチフィットチップ° BL S100X10R	¥10,400	¥11,100	
10150	Λ°トリエ イソシ 150x20mm,148cm2 ST 10x12	¥18,300	¥19,800	
10190	1000ul マルチフィットチップ° NT 1000/bag	¥4,500	¥4,800	
10230	1000ul マルチフィットチップ° NT S100X10R	¥10,700	¥11,500	
10245	スクリプト イソシ 243x243x19.30 ST 5x4袋	¥16,700	¥18,100	
10310	10ul ワタツチチップ° S96X10R	¥16,200	¥16,700	
10320	10ul ワタツチ フィルターチップ° S96X10R	¥24,900	¥25,600	
10340	200ul ワタツチチップ° S96X10R	¥16,200	¥16,700	
10350	200ul ワタツチ フィルターチップ° S96X10R	¥24,900	¥25,600	
10370	1000ul ワタツチチップ° S96X8R	¥15,500	¥16,700	
10380	1000ul ワタツチ フィルターチップ° S96X8R	¥24,900	¥25,600	
10470	200ul マルチフィットチップ° NT 1000/bag	¥3,300	¥3,600	
10520	200ul マルチフィットチップ° NT S96X10R	¥8,200	¥8,800	
10590	200ul マルチフィットチップ° YL 1000/bag	¥3,300	¥3,600	
10660	200ul マルチフィットチップ° YL S96X10R	¥8,200	¥8,800	
11035	セルカチヤーテ° イソシ 35mm ST 10x50袋 ExG	¥16,200	¥17,500	
11060	セルカチヤーテ° イソシ 60mm ST 10x50袋 ExG	¥19,700	¥21,400	
11090	セルカチヤーテ° イソシ 90mm ST 10x50袋	¥21,400	¥23,200	
11125	スクリプト イソシ 126.4x126.4x20mm ST10x24	¥32,800	¥35,500	
11150	0.65ml マイクロチューブ° NT 1000/bag	¥6,500	¥7,000	
11300	0.65ml LBマイクロチューブ° NT 500/bag	¥5,200	¥5,500	
11510	1.7ml マイクロチューブ° NT 500/bag	¥3,000	¥3,300	
12000	2ml マイクロチューブ° NT 400/bag	¥3,700	¥4,000	
12160	2.0ml LBマイクロチューブ° NT 200/bag	¥2,900	¥3,100	
13760	200ul マルチフレックスチップ° 平0.4mm 200X4R	¥43,100	¥46,300	
13770	200ul マルチフレックスチップ° 平0.4mm ST200X4R	¥57,400	¥61,800	
13790	200ul マルチフレックスチップ° 丸 200X4R	¥21,500	¥23,200	
13810	200ul マルチフレックスチップ° 丸 ST200X4R	¥24,800	¥26,700	
13980	10ul Eチップ° NT 1000/bag	¥4,100	¥4,400	
14200	1000ul マルチフィルターチップ° NT S100x10R	¥19,500	¥21,000	
14210	20ul マルチフィルターチップ° NT S96x10R	¥18,400	¥19,800	
14220	200ul マルチフィルターチップ° NT S96x10R	¥18,300	¥19,600	
15020	10ul マルチフィルター Eチップ° NT S96x10R	¥18,400	¥19,800	
15100	10ul ミニフレックスチップ° 丸 200X4R	¥23,000	¥24,700	
15110	10ul ミニフレックスチップ° 丸 ST200X4R	¥24,400	¥26,300	
15170	48well 0.5-0.65ml用 ウェルラック 6R/box	¥16,400	¥17,700	
15190	48well 1.5-2.0ml用 ウェルラック 6R/box	¥16,400	¥17,700	
15370	200ul ソフトフィットチップ° YL S96X10R	¥8,200	¥8,800	
15670	200ul ソフトフィットチップ° NT S96X10R	¥8,200	¥8,800	
15720	200ul ソフトフィットチップ° NT 1000/bag	¥3,300	¥3,600	
15730	200ul ソフトフィットチップ° YL 1000/bag	¥3,300	¥3,600	
16140	2.0ml ドルフィンチューブ° NT 1000/bag	¥10,200	¥10,900	
16249	200ml/バレル付ステンレス製フィルター- (->5L)1	¥127,800	¥148,000	
16307	ガラス製フィルター- 250ml 1	¥36,000	¥40,000	
16309	フィルター- (ガラス) 16309	¥72,000	¥80,000	
16460	1000ul用 空ラック 10R/box	¥4,600	¥4,900	
16510	ポリカーボネート製フィルター- 16510 1	¥18,500	¥18,900	
16757	イーボートMD8 本体/フィルター- /充電器	¥845,000	¥921,100	
17360	10ul ミニフレックスチップ° 平0.2mm ST200X4R	¥70,400	¥75,700	
17655	アルミ製フィルター- 1	¥92,500	¥94,400	
17803	エアリア° ラー・アクセリア°-/BACTair用アダプター	¥26,000	¥26,500	
20009	カバ°スリッパ° 9mm,48ウェル用 ST 100x2袋	¥30,600	¥33,100	
20012	カバ°スリッパ° 12mm,24ウェル用 ST 100x2袋	¥36,300	¥39,300	
20035	セルカチヤーテ° イソシ 35mm ST 10x50袋 ExG	¥18,500	¥18,900	
20060	セルカチヤーテ° イソシ 60mm ST 10x50袋 ExG	¥22,000	¥22,900	
20100	セルカチヤーテ° イソシ 90mm ST 10x20袋	¥13,000	¥14,100	
20150	セルカチヤーテ° イソシ 150x25mm ST 10x12袋	¥25,000	¥26,000	
20151	セルカチヤーテ° イソシ 150x20mm ST 10x12袋	¥21,600	¥23,400	
21008	8チャンネルリサ°-バ°- PS 7mL ST 1x100袋	¥17,500	¥18,900	
21012	12チャンネルリサ°-バ°- PS 4mL ST 1x100袋	¥18,100	¥19,500	
21750	1000XT用 空ヒンジ°ラック 8R/box	¥5,400	¥5,800	
21970	96well ウェルラック° PCR° レート NT 25/box	¥17,600	¥18,900	
21980	96well ウェルラック° PCR° レート RD 25/box	¥17,600	¥18,900	
21990	96well ウェルラック° PCR° レート BL 25/box	¥17,600	¥18,900	
22000	96well ウェルラック° PCR° レート GR 25/box	¥17,600	¥18,900	
22010	96well ウェルラック° PCR° レート YL 25/box	¥17,600	¥18,900	
22050	1チャンネルリサ°-バ°- PS 50mL ST 1x50袋	¥7,500	¥8,100	
22170	96well ウェルラック° PCR° レート PL 25/box	¥17,600	¥18,900	
22210	10ul Eチップ° NT S96X10R	¥10,600	¥11,400	
23080	96well スカ°テット° PCR° レート NT 25/box	¥17,600	¥18,900	
23570	10ul ウェルラック° NT S 96X10R	¥11,600	¥12,500	
23580	10ul ウェルラック° NT 1000/bag	¥4,200	¥4,500	
26035	SPL3Dウェルラック° イソシ 35mm 10枚x2袋	¥53,200	¥57,600	
26060	SPL3Dウェルラック° イソシ 60mm 10枚x2袋	¥59,100	¥63,900	
26410	10ul ロング°リーチ NT ST ヒンジ° 96x10R	¥13,000	¥13,900	
26430	200ul ソフトフィットチップ° NT ST ヒンジ° 96x10R	¥8,800	¥9,500	
26500	リアルタイムPCR用圧着シール 100/box	¥55,700	¥60,000	

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
26520	96wellPCRプレート用シロケット 10/box	¥12,800	¥13,700	
27730	200ulチップステーションリフトチップ YL ST96X50	¥40,300	¥43,400	
27770	200ulチップステーションリフトチップ NT ST96X50	¥40,300	¥43,400	
27790	10ulチップステーションリフトチップ NT ST96X50	¥48,000	¥51,600	
28090	300ul NXチップ用空トラック 10/box	¥5,400	¥5,800	
28200	10ulマルチフィルターチップ NT S96x10R	¥18,900	¥20,400	
28310	200ul/10ulリフトチップ用空トラック10R/box	¥5,600	¥6,000	
28360	10ulチップ用空トラック 10R/box	¥5,600	¥6,000	
28480	200ulマルチウェルチップ丸 1000/bag	¥7,500	¥8,100	
30004	カルチャープレート4well TC処理 ST 4枚x30	¥19,700	¥21,300	
30006	カルチャープレート6well TC処理 ST 50枚	¥12,100	¥12,700	
30012	カルチャープレート12well TC処理 ST 50枚	¥12,100	¥12,700	
30024	カルチャープレート24well TC処理 ST 50枚	¥12,600	¥13,200	
30048	カルチャープレート48well TC処理 ST 50枚	¥12,600	¥13,400	
30096	カルチャープレート96well TC処理 ST 50枚	¥12,100	¥12,800	
30196	カルチャープレート96well白平底 ST 50枚	¥13,200	¥14,300	
30260	200ul用空トラック 10R/box	¥2,900	¥3,200	
30296	カルチャープレート96well黒平底 ST 50枚	¥13,200	¥14,300	
30340	Sharp 10ulフィルター NT ST 96x10R	¥19,100	¥20,600	
30350	96well FastPCRプレート WH 25/box	¥20,400	¥21,900	
30396	カルチャープレート96well白平底 50	¥14,400	¥15,600	
30430	200ulクイックトラックリフトチップ NT 96X50	¥33,500	¥36,000	
30440	200ulクイックトラックリフトチップ YL 96X50	¥35,200	¥37,800	
30450	200ulクイックトラックマルチチップ NT 96X50	¥33,500	¥36,000	
30460	200ulクイックトラックマルチチップ YL 96X50	¥33,500	¥36,000	
30470	10ulクイックトラックチップ NT 96X50	¥34,500	¥37,100	
30550	200NXフィルターチップ NT S96X10R	¥18,400	¥19,800	
30570	Sharp 10ul NT ST 96x10R	¥13,000	¥13,900	
30800	300ul XTチップ NT 1000/bag	¥4,200	¥4,500	
30930	300ul XTチップ NT S96X10R	¥9,800	¥10,500	
31024	カルチャープレート24well 50枚	¥13,100	¥14,200	
31096	カルチャープレート96well 10枚x10袋	¥24,900	¥27,000	
31610	1000ul XTチップ NT 768/bag	¥5,100	¥5,500	
31770	300ul NXチップ NT 1000/bag	¥4,100	¥4,500	
31790	300ul NXチップ NT S96X10R	¥9,400	¥10,100	
32006	カルチャープレート6well ST 50枚	¥12,500	¥13,500	
32012	カルチャープレート12well ST 50枚	¥12,500	¥13,500	
32024	カルチャープレート24well 平底 ST 50枚	¥13,100	¥14,200	
32048	カルチャープレート48well ST 1x50袋	¥13,700	¥14,800	
32096	カルチャープレート96well 平底 ST 50枚	¥13,500	¥14,700	
33196	フラックプレート96well ST 1x20袋	¥203,100	¥219,800	
33212	フラックプレート12well TC処理 ST 1x10袋	¥84,700	¥91,700	
33224	フラックプレート24well TC処理 ST 1x10袋	¥84,700	¥91,700	
33296	フラックプレート96well TC処理 ST 1x20袋	¥148,900	¥161,200	
33312	フラックプレート12well TC処理 ST 1x10袋	¥32,600	¥35,300	
33324	フラックプレート24well TC処理 ST 1x10袋	¥32,600	¥35,300	
33384	384well HTプレート, 10x4袋	¥25,200	¥25,500	
33396	フラックプレート96well TC処理 ST 5x5袋	¥43,200	¥46,700	
33696	初回プレート96well TC処理 ST 5枚x5袋	¥39,500	¥42,700	
34000	1000XTフィルターチップ NT S96X8R	¥19,700	¥21,200	
34096	カルチャープレート96well 丸底 ST 1x50袋	¥13,800	¥14,900	
34596	ティッププレート96well/0.8mL V底 5x5袋	¥26,800	¥29,000	
34696	ティッププレート96well/2mL U底ST 5x5袋	¥31,000	¥33,600	
34760	1000ul XTチップ NT S96X8R	¥10,700	¥11,500	
34770	1000XTクイックトラック NT 96X40	¥35,500	¥38,100	
35224	Insert™ リング 24well ST 12x4袋	¥37,000	¥40,000	
35524	Insert™ スタテリング 24well ST 12x4袋	¥37,000	¥40,000	
35810	96well セミスケレットPCRプレート GR 25/box	¥20,800	¥22,300	
35820	96well セミスケレットPCRプレート RD 25/box	¥20,800	¥22,300	
35830	96well セミスケレットPCRプレート YL 25/box	¥20,800	¥22,300	
35840	96well セミスケレットPCRプレート BL 25/box	¥21,000	¥22,600	
35850	96well セミスケレットPCRプレート PL 25/box	¥20,800	¥22,300	
36060	100ulマルチフィルターチップ NT S96x10R	¥18,100	¥19,500	
36590	リアルタイムPCR用粘着シール 100/box	¥44,000	¥47,300	
37384	カルチャープレート384wellTCヨリ ST 10x4袋	¥61,800	¥66,900	
37390	96well セミスケレットPCRプレート WH 25/box	¥21,000	¥22,600	
37650	10ul マイクロリフトチップ NT 1000/bag	¥4,700	¥5,000	
37660	10ul マイクロリフトチップ NT S96X10R	¥12,300	¥13,300	
37820	10ul クイックトラックリフトチップ NT 96X50	¥35,500	¥38,200	
37870	96well カムトラックPCRプレート WH 25/box	¥18,500	¥19,900	
37990	8連チップ(フラットチップ) NT 250ストリップ	¥12,200	¥13,100	
38000	10ul フィルターリフトチップ NT S96x10R	¥19,600	¥21,100	
38780	8連チップ(トームチップ) NT 250ストリップ	¥11,500	¥12,400	
38800	96well FastPCRプレート NT 25/box	¥14,500	¥15,600	
39024	カルチャープレート24well ST 5枚	¥21,300	¥23,000	
39096	カルチャープレート96well ST 5枚	¥21,300	¥23,000	
39224	カルチャープレート24well ST 1x5袋	¥21,300	¥23,000	
39620	384well NX PCRプレート NT 50/box	¥41,800	¥45,000	
40005	テストチューブ 5ml, 12x75mm, EOG滅菌 25x20	¥13,200	¥15,600	
40205	テストチューブ 5ml, 12x75mm, 100本x10袋	¥12,000	¥14,900	
43012	SPLクイックイアル 1.8mL ST 50x10袋 709	¥19,600	¥21,200	
43021	SPLクイックイアル 1.2mL ST 50x10袋 709	¥25,200	¥27,200	
43022	SPLクイックイアル 1.8mL ST 50x10袋 709	¥24,200	¥25,900	
43023	SPLクイックイアル 5.0mL ST 50x10袋 709	¥38,300	¥41,400	
43032	チップインサート100x5袋/5色	¥6,400	¥6,900	
43111	SPLクイックイアル 1.2mL ST 50x10袋 709	¥25,200	¥27,200	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
43112	SPLクワイバ [®] イル 1.8mL ST 50x10袋 アウター	¥25,900	¥28,100	
43113	SPLクワイバ [®] イル 5.0mL ST 50x10袋 アウター	¥38,300	¥41,400	
44022	SPLクワイバ [®] イル 1.8mL ST 50x10袋 インナー	¥31,600	¥34,200	
44112	SPLクワイバ [®] イル 1.8mL ST 50x10袋 アウター	¥31,700	¥34,300	
45002	スクリューキャップ [®] チューブ 2.0mL 500本x4袋	¥84,600	¥91,600	
50015	遠心チューブ 15mL ST 50本x10袋	¥16,800	¥18,200	
50050	遠心チューブ 50mL ST 25本x20袋	¥21,800	¥23,600	
50115	遠心チューブ 15mL ST 25本x20パック	¥19,200	¥20,800	
50150	遠心チューブ 50mL ST 25本x12パック	¥16,300	¥17,600	
50215	スクリューキャップ [®] チューブ 15ml, エコカ, ST 50x10袋	¥17,900	¥18,800	
50250	遠心チューブ 50mL スカート ST 25本x20袋	¥32,400	¥35,000	
50351	バ [®] イルアウター, 50ml, ST 3x30袋	¥18,300	¥19,800	
50650	スクリューキャップ [®] チューブ 50ml, スカート, ST 25x20袋	¥29,100	¥31,500	
51005	スクリューキャップ [®] チューブ エコカ 5mL 100本x2袋	¥17,100	¥18,500	
51105	スクリューキャップ [®] チューブ エコカ5mLST100本x2袋	¥20,500	¥22,200	
52005	5ml スクリューキャップ [®] アダプター, アセラル, 10本	¥7,100	¥7,700	
52050	遠心チューブ 50ml用ラック 25本立て 20パック	¥6,000	¥6,500	
52202	2ウェル遠心チューブ用ラック 20個	¥3,900	¥4,300	
53005	スクリューキャップ [®] チューブ 自立 5mL 100本x2袋	¥19,700	¥21,300	
53105	スクリューキャップ [®] チューブ 自立5mL ST100本x2袋	¥23,500	¥25,500	
54015	遠心チューブ 15mL ST 50本x10袋	¥38,600	¥41,800	
54115	遠心チューブ 15mL ST 50本x4袋	¥18,100	¥19,500	
55085	ローホ [®] トル 850cm2 TC処理 ST 2本x20袋	¥47,500	¥51,400	
56000	スクリューキャップ [®] チューブ 1000ml, ST 12本x2袋	¥18,200	¥19,700	
56005	スクリューキャップ [®] チューブ エコカ 5mL 100本x2袋	¥16,500	¥17,800	
56105	スクリューキャップ [®] チューブ エコカ5mLST100本x2袋	¥19,400	¥20,900	
56125	スクリューキャップ [®] チューブ 125ml, ST 24x4袋	¥34,600	¥37,400	
56250	スクリューキャップ [®] チューブ 250ml, ST 30x2袋	¥25,900	¥28,100	
56500	スクリューキャップ [®] チューブ 500ml, ST 24x2袋	¥24,300	¥26,300	
61050	クワイバ [®] イルラック (50ウェル) 10個	¥6,900	¥7,400	
62015	マイクロ遠心チューブ 1.5mL 200本x25袋	¥32,800	¥35,400	
62017	マイクロ遠心チューブ 1.7mL 500本x10袋	¥32,800	¥35,400	
70012	セルカルチャーフラスコ 25cm2 TC処理 5x40袋	¥36,600	¥39,600	
70025	セルカルチャーフラスコ 25cm2 TC処理 5x40袋	¥22,300	¥24,100	
70075	セルカルチャーフラスコ 75cm2 TC処理 ST 5x20袋	¥21,600	¥23,400	
70112	セルカルチャーフラスコ 25cm2 TC処理 5x40袋	¥59,800	¥64,700	
70125	セルカルチャーフラスコキャップ [®] 付 25cm2 ST 5x40袋	¥20,900	¥22,600	
70175	セルカルチャーフラスコキャップ [®] 付 75cm2 ST 5x20袋	¥20,600	¥22,300	
70212	セルカルチャーフラスコ 25cm2 5x40袋	¥59,800	¥64,700	
70225	セルカルチャーフラスコキャップ [®] 付 25cm2 ST 5x40袋	¥24,000	¥26,000	
70275	セルカルチャーフラスコキャップ [®] 付 75cm2 ST 5x20袋	¥22,600	¥24,400	
70312	セルカルチャーフラスコ 25cm2 5x40袋	¥63,200	¥68,400	
70325	セルカルチャーフラスコキャップ [®] 付 25cm2 ST 5x40袋	¥25,100	¥27,200	
70375	セルカルチャーフラスコキャップ [®] 付 75cm2 ST 5x20袋	¥23,700	¥25,600	
70500	10ul リド [®] リド [®] NT 96x40	¥31,900	¥34,300	
70510	200ul リド [®] リド [®] NT 96x40	¥29,900	¥32,200	
70520	200ul リド [®] リド [®] YL 96x40	¥29,900	¥32,200	
70530	300NX リド [®] リド [®] NT 96x40	¥38,600	¥41,600	
70540	1000XT リド [®] リド [®] NT 96x32	¥38,400	¥41,300	
70600	10ul リド [®] リド [®] LB 96x40	¥34,500	¥37,100	
70610	200ul リド [®] リド [®] LB 96x40	¥32,500	¥35,000	
70620	300NX リド [®] リド [®] LB 96x40	¥41,000	¥44,100	
70700	10ul リド [®] リド [®] スターキットNT 96x10	¥11,400	¥12,300	
70710	200ul リド [®] リド [®] スターキットNT 96x10	¥10,600	¥11,400	
70720	200ul リド [®] リド [®] スターキットYL 96x10	¥10,600	¥11,400	
70730	300ul リド [®] リド [®] スターキットNT 96x10	¥11,700	¥12,600	
70740	1000ul リド [®] リド [®] スターキットNT 96x8	¥11,500	¥12,400	
71000	10ul 200ul 300ul 空ラック 24R/box	¥13,000	¥14,000	
71010	1000XT用空ラック 16R/box	¥11,700	¥12,600	
71175	セルカルチャーフラスコ 175cm2 TC処理 ST 5x8袋	¥17,000	¥18,400	
72175	セルカルチャーフラスコキャップ [®] 付 175cm2 ST 5x8袋	¥15,900	¥17,200	
73250	セルカルチャーフラスコ 250ml ST 1x8袋	¥9,200	¥9,900	
73500	セルカルチャーフラスコ 500ml ST 1x8袋	¥11,500	¥12,500	
74020	リザ [®] -バ [®] -セル [®] -SBS BK 10/box	¥12,800	¥13,700	
74030	200ml 96wellリザ [®] -バ [®] -深 WH ST 50/box	¥31,900	¥34,300	
74250	セルカルチャーフラスコ 250mlフィルターキャップ [®] ST1x8袋	¥10,100	¥10,900	
74310	96well セミスケレット [®] PCR [®] レート NT25/box	¥21,000	¥22,600	
74320	96well セミスケレット [®] PCR [®] レート GR25/box	¥21,000	¥22,600	
74330	96well セミスケレット [®] PCR [®] レート RD25/box	¥21,000	¥22,600	
74340	96well セミスケレット [®] PCR [®] レート YL25/box	¥21,000	¥22,600	
74350	96well セミスケレット [®] PCR [®] レート BL25/box	¥21,000	¥22,600	
74360	96well セミスケレット [®] PCR [®] レート PL25/box	¥21,000	¥22,600	
74370	96well セミスケレット [®] PCR [®] レートNT25x4/box	¥83,500	¥89,800	
74380	96well セミスケレット [®] PCR [®] レート BK25/box	¥23,700	¥25,500	
74390	96well セミスケレット [®] PCR [®] レート WH25/box	¥23,700	¥25,500	
74500	セルカルチャーフラスコ 500mlフィルターキャップ [®] ST1x8袋	¥12,700	¥13,700	
75250	セルカルチャーフラスコ 250ml 1x8袋	¥8,000	¥8,700	
76250	セルカルチャーフラスコ 250ml フィルターキャップ [®] 1x8袋	¥8,800	¥9,500	
80025	クワイバ [®] ユス 5x5(25穴) 1x20袋	¥11,700	¥12,600	
80081	クワイバ [®] ユス 9x9(81穴) 1x10袋	¥11,700	¥12,600	
80281	クワイバ [®] ユス 5ml 9x9(81穴) 1x6袋	¥10,300	¥11,100	
90001	ループ [®] &コントロール 1uL, 195mm ST 10x50袋	¥7,500	¥8,200	
90010	ループ [®] &コントロール 10uL, 195mm ST 10x50袋	¥7,500	¥8,200	
90020	セルスルーバ [®] - ST 1x100袋	¥10,500	¥11,300	
90021	セルスルーバ [®] - ST 1x100袋	¥10,500	¥11,300	
90030	セルスルーバ [®] - ST 1x100袋	¥11,800	¥12,800	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
90031	セルスルー- ST 1x100袋	¥11,800	¥12,800	
90032	セルスルー- ST 1x100袋	¥26,400	¥28,600	
90040	セルスルー- ST 1x100袋	¥8,800	¥9,500	
90050	スリット- ST 10x50袋	¥10,500	¥11,300	
91001	細胞シカヒ [®] 1mL ST 200本x4箱	¥30,900	¥33,400	
91002	細胞シカヒ [®] 2mL ST 150本x4箱	¥25,300	¥27,400	
91005	細胞シカヒ [®] 5mL ST 100本x4箱	¥21,900	¥23,700	
91010	細胞シカヒ [®] 10mL ST 100本x4箱	¥23,700	¥25,700	
91025	細胞シカヒ [®] 25mL ST 50本x4箱	¥18,000	¥19,500	
91050	細胞シカヒ [®] 50mL ST 40本x4箱	¥54,100	¥58,600	
91110	カッティングトッパ [®] 10mL ST 100本x4箱	¥29,400	¥31,900	
93001	細胞シカヒ [®] 1mL ST 50本x20袋	¥34,900	¥37,800	
93002	細胞シカヒ [®] 2mL ST 50本x20袋	¥38,400	¥41,600	
93005	細胞シカヒ [®] 5mL ST 50本x10袋	¥24,800	¥26,800	
93010	細胞シカヒ [®] 10mL ST 50本x10袋	¥26,800	¥29,100	
93025	細胞シカヒ [®] 25mL ST 25本x8袋	¥16,300	¥17,600	
93040	セルスルー- 40μm ST 1x50袋	¥13,600	¥14,700	
93050	細胞シカヒ [®] 50mL ST 25本x8袋	¥61,500	¥66,600	
93070	セルスルー- 70μm ST 1x50袋	¥13,000	¥14,100	
93100	セルスルー- 100μm ST 1x50袋	¥14,200	¥15,400	
94020	マルチセルスルー-20μm ST 1x50袋	¥35,500	¥38,500	
94030	マルチセルスルー-30μm ST 1x50袋	¥35,500	¥38,500	
94040	マルチセルスルー-40μm ST 1x50袋	¥34,600	¥37,400	
94070	マルチセルスルー-70μm ST 1x50袋	¥34,100	¥36,900	
94100	マルチセルスルー-100μm ST 1x50袋	¥35,900	¥38,800	
95001	細胞シカヒ [®] 1mL ST 200本x4袋	¥29,400	¥31,900	
95002	細胞シカヒ [®] 2mL ST 150本x4袋	¥24,100	¥26,100	
95005	細胞シカヒ [®] 5mL ST 100本x4袋	¥20,800	¥22,600	
95010	細胞シカヒ [®] 10mL ST 100本x4袋	¥22,600	¥24,500	
95025	細胞シカヒ [®] 25mL ST 50本x4袋	¥17,200	¥18,600	
95050	細胞シカヒ [®] 50mL ST 40本x4袋	¥51,800	¥56,100	
96000	SPL粘着プレート 143x79mm 100枚	¥19,100	¥20,600	
97003	透析チャンバ - 3mL MWCO 3.5K 1x10袋	¥32,000	¥34,600	
97103	透析チャンバ - 15mL MWCO 3.5K 1x8袋	¥41,800	¥45,300	
101043	ヒト [®] 2.9 mm 1枚	¥9,300	¥9,600	
104026	16-Base B 300個	¥69,900	¥71,900	
104143	6-Tube B 1536個	¥89,500	¥92,000	
190022	Loading Solution 500μL	¥1,800	¥2,000	
190039	PippinElectrophoresisBuffer 60ml 1	¥9,200	¥9,800	
201008	リチウムイオン [®] 40mg @	¥7,800	¥8,000	
201207	インキ [®] 25ml	¥7,800	¥8,000	
201924	スクラパー 24ウェル用 先端幅0.5mm ST 20袋	¥8,800	¥9,500	
201925	スクラパー 24ウェル用 先端幅1mm ST 20袋	¥8,800	¥9,500	
210001	ルートドライ [®]	¥73,500	¥78,400	
310061	アソシエート 100	¥10,700	¥11,600	
310062	アソシエート 100	¥29,100	¥31,500	
310074	インキ [®] 30x2袋	¥3,700	¥4,000	
310100	プレートカッター [®] 100x40mm, ST 5x40	¥24,000	¥26,000	
310103	プレートカッター [®] 100x50mm, 5x40	¥38,900	¥40,000	
558250	遠心チューブ 250mL ST 6本x9袋	¥23,800	¥25,800	
562000	バ [®] 2L, ST 4本x20袋	¥156,500	¥169,300	
565000	バ [®] 5L, ST 1本x6袋	¥98,200	¥106,300	
651463	1.2mL Cryoアンプル 144/box	¥37,600	¥40,100	
651483	1.2mL Cryoアンプル(Pre-scored)144/box	¥45,600	¥48,700	
651486	2.0mL Cryoアンプル(Pre-scored)144/box	¥52,200	¥55,700	
651490	1.0mL Cryoアンプル用BOX 15/box	¥68,800	¥73,400	
651492	2.0mL Cryoアンプル用BOX 15/box	¥52,800	¥56,300	
651502	1.0mL 凍乾アンプル(Pre-scored)144/box	¥59,800	¥63,800	
651905	5.0mL Vacule [®] 200/box	¥103,400	¥110,300	
651907	10mL Vacule [®] 200/box	¥98,000	¥104,600	
651954	4mL Vacule [®] 200/box	¥73,800	¥78,800	
700225	セルスルー- 225cm ² TC処理 ST 5x5袋	¥23,400	¥25,400	
701225	セルスルー- 225cm ² ST 5x5袋	¥22,300	¥24,100	
720040	BluePippin/PippinPrep用 リンクセット	¥9,200	¥9,800	
901998.1	MMA30用秤量皿(ガラス製) 200	¥7,500	¥8,200	2024/3/14 追加
906940	MA45用秤量皿(ガラス製) 80枚	¥9,200	¥10,200	2024/3/14 追加
906941	ガラス製-製フィルタ [®] 軟質 200枚入	¥23,200	¥25,800	2024/3/14 追加
6980115	ガラス製ベース 30ml	¥10,800	¥11,600	
6980116	ガラス製ベース	¥9,200	¥9,800	
6980178	フィルタ [®] -用シリコン [®] O-リング 42x3mm	¥1,200	¥1,300	
6980228	上フタ(PC製フィルタ [®] -用)	¥2,500	¥2,600	
6980230	上部フタ 250ml(PC製フィルタ [®] -用)	¥3,700	¥3,800	
6980233	フィルタ [®] -用シリコン [®] O-リング	¥3,500	¥3,600	
6980235	シリコン [®] O-リング 14x2mm(PC製フィルタ [®] -用)	¥2,500	¥3,000	
6980236	シリコン [®] 製キャップ(PC製フィルタ [®] -用)10	¥3,500	¥3,600	
6980569	シリコン [®] 製キャップ(PC製フィルタ [®] -用)10	¥3,000	¥3,200	
6980570	シリコン [®] 製キャップ(PC製フィルタ [®] -用)10	¥2,100	¥2,200	
6980737	フィルタ [®] -用PTFE加工アング [®] ドレインφ47	¥8,000	¥8,400	
6981090	シリコン [®] 製キャップ(PC製フィルタ [®] -用)	¥4,700	¥4,800	
6988094	シリコン [®] 製キャップ 1	¥10,600	¥11,200	
3023982001	O-RING	¥1,600	¥2,000	
3053563001	TENSION GAUGE FOR DRIVE BELTS	¥330,000	¥466,000	
4468325001	LID HOUSING UPPER	¥214,000	¥272,000	
4468333001	SWITCH THERMO	¥4,000	¥5,000	
4468350001	SENSOR MODULE	¥55,000	¥73,000	
4468368001	PCB MAIN	¥74,000	¥97,000	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
4468384001	HOUSING	¥155,000	¥197,000	
4468414001	CABLE SENSOR	¥900	¥2,000	
4468422001	FAN INTERN	¥2,300	¥3,000	
4468481001	SPACER SHORT	¥5,800	¥6,000	
4468490001	CHAMBER MEASURE	¥340,000	¥350,000	
4468511001	DAMPER LID	¥31,000	¥32,000	
4468520001	PHOTOMETER LC1.3	¥2,100,000	¥2,110,000	
4468546001	SCREW SET LC 1.3	¥500	¥1,000	
4468554001	FLANGE	¥116,000	¥120,000	
4468589001	LID LOCK	¥340,000	¥434,000	
4468597001	CHAMBER LID	¥320,000	¥330,000	
4468619001	FAN WHEEL	¥18,000	¥19,000	
4468627001	STEPMOTOR ROTOR	¥130,000	¥140,000	
4468635001	WINDOW OPTIC	¥20,000	¥21,000	
4685954001	LOADING UNIT	¥1,460,000	¥1,910,000	
4685989001	READER BARCODE INT.	¥130,000	¥177,000	
4685997001	POWER SUPPLY	¥439,000	¥546,000	
4686012001	PCB STATUS BOARD	¥20,000	¥21,000	
4686039001	PCB BLOCK TC CONTROL	¥190,000	¥276,000	
4686047001	PCB DETECTION CONTROL	¥80,000	¥106,000	
4686055001	PCB DATAREDUCTION	¥215,000	¥291,000	
4686063001	PANEL LEFT SIDE	¥90,000	¥132,000	
4686071001	PANEL RIGHT SIDE	¥88,000	¥131,000	
4686098001	BLOCK CYCLER DOOR W. FAN	¥297,000	¥413,000	
4686128001	LC 480 Ventilation Dust Filter 4入	¥2,000	¥3,000	
4686144001	MOTOR MWP LOADING	¥106,000	¥151,000	
4686152001	MOTOR COVER LIFT	¥170,000	¥210,000	
4686179001	CABLE FAN	¥1,600	¥2,000	
4686187001	CABLE MWP DETECTOR	¥1,600	¥2,000	
4686195001	LIGHT BARRIER L350	¥5,600	¥7,000	
4686217001	SWITCH MICRO	¥3,000	¥4,000	
4686225001	SOLENOID	¥36,000	¥49,000	
4686233001	DETECTOR MWP	¥38,000	¥52,000	
4686268001	FUSE 10A T 250V 5X20 SET OF 10	¥200	¥1,000	
4720610001	SYSTEMCOMPONENTS LC 2.0	¥280,000	¥290,000	
4793897001	FAN 24VDC 60X60 ASSY	¥7,000	¥9,000	
4824121001	O-RING CO-RE, SET OF 8	¥15,000	¥16,000	
4871731001	RETAINER MULTI WELL PLATE	¥120,000	¥162,000	
4899555001	PROTECTIVE CAPS FOR OPTICS	¥2,000	¥3,000	
4900375001	CABLE MAINS CONNECTION USA	¥5,000	¥7,000	
4914899001	DISK DRIVE HD SATA DC7600SFF	¥44,000	¥45,000	
4914902001	DISK DRIVE FLOPPY DC7600CMT	¥60,000	¥76,000	
4914929001	DISK DRIVE FLOPPY DC7600SFF	¥32,000	¥40,000	
4914937001	DVD DRIVE +/-RW DC7600	¥87,000	¥109,000	
4914945001	KEYBOARD PS/2 US INT.	¥8,000	¥10,000	
5065976001	LID ROTOR MAGNA LYSER	¥264,000	¥321,000	
5091268001	CABEL GROUNDING	¥2,000	¥3,000	
5128234001	ROPE TRACTION LOADING	¥21,000	¥31,000	
5128242001	FAN 24V ASSY	¥28,000	¥29,000	
5206243001	KIT MAINTENANCE MAGNA LYSER	¥58,000	¥87,000	
5214947001	CERTIFICATE IQ LC480	¥7,600	¥8,000	
5215013001	CERTIFICATE OQ LC480	¥8,300	¥9,000	
5360021001	TOOL BELT TENSION MEASURING V2	¥445,000	¥568,000	
5382785001	CABLE RS232 INTERNAL	¥27,000	¥34,000	
5382807001	KIT FOR CHAMBER LID	¥105,000	¥110,000	
5449596001	PCB R-HEAD	¥94,000	¥95,000	
5449766001	PCB INTERLOCK	¥75,000	¥76,000	
5517133001	LAN PCIE FH NIC RP5700/DC7800CMT	¥20,000	¥25,000	
5523800001	TRAY WASTE INCLUDING FUNNEL	¥66,000	¥87,000	
5571740001	TOOL SCREW DRIVER BIT 1/4in 10ヶ付	¥2,300	¥3,000	
5587131001	CABLE SET LOADING LC480	¥34,000	¥42,000	
5652766001	TOOL BRUSH 12X134MM SET OF 10	¥7,200	¥10,000	
5667976001	PUMP-MEMBRAN 24VDC/18W	¥240,000	¥300,000	
5667984001	PUMP-MEMBRAN 24VDC/7.8W	¥120,000	¥146,000	
5667992001	BLOCK CONNECTING FLUID-CHANNEL	¥20,000	¥29,000	
5668018001	SOLENOID SHUTTER 600MM MOLEX 1.25	¥16,000	¥22,000	
5668026001	RACK BALANCE (CUBE) FSR ONLY	¥90,000	¥119,000	
5675693001	FUSE 2.5AT H 250V 5X20 SET OF 10	¥200	¥1,000	
5694973001	SWITCH MICRO	¥7,000	¥8,000	
5741424001	FOAM BLOCK SHIPPING TC LC480	¥6,300	¥7,000	
5772940001	COVER INSTRUMENT LC480 II	¥150,000	¥224,000	
6309828001	KIT RETROFIT EXTERNAL FLUID	¥608,000	¥804,000	
6407498001	TOOL CALIBRATION WEIGHT 10G OWS	¥73,000	¥95,000	
6407501001	TOOL WEIGHING PAN 5ML OWS	¥81,000	¥110,000	
6407528001	POWER ADAPTER EL-BOX MT OWS	¥110,000	¥120,000	
6436170001	KIT OF ALL CLAMPS	¥60,000	¥72,000	
6436188001	SWITCH MAGNETIC	¥7,000	¥9,000	
6436196001	DAMPER LOAD FLAP	¥37,000	¥38,000	
6436218001	QUICK-CONNECTOR MALE PLCD D8	¥13,000	¥17,000	
6449964001	CASE WEIGHING SYSTEM ONBOARD	¥231,000	¥290,000	
6548300001	COVER HEATER ASSY	¥200,000	¥285,000	
6548318001	DETECTION ASSY	¥1,350,000	¥1,460,000	
6548326001	MOTOR STEP ASSY	¥108,000	¥140,000	
6548334001	CYCLER ASSY	¥1,170,000	¥1,080,000	
6548342001	POWER SUPPLY AC/DC	¥68,000	¥85,000	
6548385001	DISPLAY ASSY	¥82,000	¥102,000	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
6548458001	LIGHT BARRIER	¥8,000	¥11,000	
6591043001	PCB RT CPU 168 SERVICE V2	¥316,000	¥411,000	
6593437001	COVER DRAWER L/R	¥9,000	¥12,000	
6621619001	FUSE 8A T 5X20 250VAC/150VDC UL	¥450	¥1,000	
6621643001	FUSE 3.15A T 5X20 250VAC/300VDC UL	¥1,800	¥2,000	
6652697001	FLAP MAINTENANCE	¥289,000	¥290,000	
6652719001	HOOD SERVICE	¥339,000	¥340,000	
6652751001	UPGRADEKIT LC480 GSCAN SW1.5.1	¥11,000	¥12,000	
6656447001	DAMPING SLIDERS R	¥5,600	¥6,000	
6656463001	COVER TOP	¥38,000	¥51,000	
6656480001	COVER RIGHT	¥33,000	¥42,000	
6713033001	COVER FILTER PAD	¥6,700	¥7,000	
6727867001	DETECTION UNIT COMPLETE IV	¥2,580,000	¥2,690,000	
6727891001	KIT UPGRADE LED	¥19,000	¥23,000	
6778909001	TOOL TEMPERATURE QUALIFICATION 96	¥1,300,000	¥1,640,000	
6821235001	TRANSPORT LOCKING DEVICE	¥1,600	¥2,000	
6889964001	ADHESIVE TAPE TESA 4287	¥200	¥1,000	
6889972001	ADHESIVE TAPE KAPTON	¥200	¥1,000	
6889999001	COVER CYCLER	¥42,000	¥55,000	
6920683001	TRANSPORT LOCK - FAILURE STATE	¥1,300	¥2,000	
6920691001	COMPACT FLASH CARD SW1.1	¥118,000	¥120,000	
7030576001	CHASSIS UPGRADE KIT	¥315,000	¥409,000	
7105444001	SW OSW7 Recovery DVD RP5800 V2.1.1	¥5,200	¥6,000	
7317727001	COLLAR EXTERNAL WASTE	¥35,000	¥45,000	
7317832001	SHIMS LPS SET OF 4	¥8,500	¥11,000	
7320582001	SW LC2.0/1.x V5.0.0 (LSR only) CD	¥5,000	¥7,000	
7379773001	TIP RETAINING PLATE	¥65,000	¥84,000	
7391145001	Drive Blu-ray USB	¥59,000	¥75,000	
7396171001	SPACER (SET OF 4 PCS)	¥5,800	¥6,000	
7657544001	CABLE FFC POWER CYCLER CONTROL	¥15,000	¥16,000	
7657579001	CABLE FFC DISPLAY	¥9,600	¥10,000	
7657587001	CABLE LCD 4-WIRE	¥12,000	¥16,000	
7699000001	CIRCUIT BREAKER RESETTABLE 4A	メーカー価格	¥3,000	
7743556001	BLOCK CYCLER 96 S LC480	¥900,000	¥1,540,000	
7743564001	BLOCK CYCLER 384 S LC480	¥960,000	¥820,000	
7802455001	POWER SUPPLY 500W 24V	¥184,000	¥221,000	
7802463001	PROCESSING STATION	¥600,000	¥610,000	
7802471001	COOLING STATION	¥523,000	¥530,000	
7928475001	GUIDE RAIL T L=530	¥101,000	¥131,000	
7928505001	GUIDE RAIL U L=530	¥103,000	¥131,000	
7928530001	Z-AXIS MODULE	¥490,000	¥469,000	
7928556001	MOTOR DC INCL. PULLY	¥101,000	¥110,000	
7928564001	CAM-DISCS SET OF 2	¥5,200	¥6,000	
7928572001	DAMPER COVER FRONT	¥171,000	¥218,000	
7928599001	TRAY DECK	¥12,000	¥15,000	
7930879001	PCB INTERLOCK	¥148,000	¥190,000	
7930887001	PCB 3 DCM2TEMP	¥150,000	¥201,000	
7930895001	MOTOR DCX35L	¥101,000	¥129,000	
7931018001	COVER FRONT	¥330,000	¥400,000	
7931026001	TOOL PIPETTER RELEASER	¥110,000	¥120,000	
7931034001	SOCKET UV TUBE LIGHT	¥27,000	¥33,000	
7931069001	SCANNER CCD LATERAL	¥117,000	¥120,000	
7931077001	HANDHELD BARCODE READER	¥91,000	¥121,000	
7931085001	TOUCHPANEL MODULE COMPL.	¥420,000	¥576,000	
7943300001	BARCODE READER SERVICE	¥97,000	¥118,000	
7972369001	CABLE COOLING UNIT SET OF 2	¥25,000	¥48,000	
7988001001	LIGHT SOURCE 2.1 ASSY	¥208,000	¥210,000	
7995580001	SWITCH POWER	¥7,800	¥8,000	
8017018001	CABLE FFC COMMUNICATION	¥9,800	¥10,000	
8020523001	PCB MAINBOARD	¥310,000	¥320,000	
8020531001	KIT UPDATE WASTE UNIT INTERNAL	¥241,000	¥323,000	
8021422001	TRAY LOWER STAGE INTERNAL	¥130,000	¥171,000	
8021473001	AXIS-Z R-HEAD CHANNEL 2 AND 4	¥472,000	¥480,000	
8021481001	MOTOR DC ASSY WITH PULLY Z15 S2M	¥117,000	¥120,000	
8022283001	FILTER METAL 1/4 INCH THREADED	¥100,000	¥128,000	
8022542001	MOTOR 24V ASSY WITH PULLY Z20/T2.5	¥203,000	¥210,000	
8022569001	AXIS-Z R-HEAD CHANNEL 1 AND 3	¥460,000	¥470,000	
8022917001	FLAP LOADING	¥241,000	¥291,000	
8022933001	SENSOR LPS COMPL.	¥670,000	¥976,000	
8022968001	FILTER COOLING UNIT DRAWER MODULE	¥73,000	¥74,000	
8022976001	PCB ROBOTICS CONTROL 2.1V	¥250,000	¥260,000	
8022984001	PUMP MICRO ANNULAR GEAR MZR-4622	¥680,000	¥690,000	
8022992001	HEAD PROCESS 96 CHANNELS V2	¥4,090,000	¥7,130,000	
8023000001	PROCESS STATION COMPL.	¥1,570,000	¥1,760,000	
8023018001	COOLING STATION COMPL	¥920,000	¥1,240,000	
8023026001	POWER BOX COMPL	¥410,000	¥420,000	
8023034001	SWITCH CAPACITIVEPROXIMITY BALLUFF	¥97,000	¥148,000	
8023069001	MOTOR ASSY DROP CATCHER	¥406,000	¥410,000	
8033994001	TRAY PROCESSING STATION	¥23,000	¥24,000	
8034010001	MOTOR DC, 2W, 21V	¥114,000	¥120,000	
8034028001	MOTOR 4,5W, 24V INCL.GEAR	¥201,000	¥210,000	
8034036001	CABLE HANDHELD BARCODE READER USB	¥6,300	¥9,000	
8035288001	TOOL MP24 MOUNTING KIT OWS	¥12,000	¥34,000	
8035318001	SET FUSE MP24	¥1,600	¥2,000	
8035334001	COVER SIDE, LEFT AND RIGHT	¥324,000	¥444,000	
8035342001	COVER TOP	¥133,000	¥164,000	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
8035385001	CABLE SET POWER SUPPLY	¥109,000	¥138,000	
8035393001	CABLE RIBBON PIPETTER, SET OF 2	¥6,900	¥7,000	
8042187001	TRAY TIP PARK PEN SET OF 3	¥3,400	¥4,000	
8053472001	INRUSH CURRENT LIMITERKNFPUMP MP96	¥25,000	¥31,000	
8099022001	mSATA 64GB	¥55,000	¥70,000	
8153922001	FLAT BAND CABLE SET OF 3	¥183,000	¥190,000	
8153949001	SIDE COVER TOP	¥69,000	¥86,000	
8153965001	SIDE COVER LEFT	¥94,000	¥116,000	
8153981001	SIDE COVER RIGHT	¥94,000	¥116,000	
8154015001	GUIDE RAIL (SET OF 2)	¥325,000	¥392,000	
8229279001	CCD-SCANNER DCR200i - RS232	¥430,000	¥569,000	
8234515001	KEYBOARD US TRACKPAD USB	¥158,000	¥160,000	
8274754001	KIT UPGRADE IPC-BAYTRAIL	¥32,000	¥39,000	
8274762001	IPC-BAYTRAIL LC96	¥234,000	¥302,000	
8274789001	PCB INFRASTRUCTURE LC96	¥87,000	¥122,000	
8316937001	PCBA INFRASTRUCTURE	¥99,000	¥126,000	
8417458001	USB FLASH DRIVE 32GB/HW ENC	¥72,000	¥94,000	
8493634001	FIXTURE LAN CONNECTOR	¥900	¥1,000	
8735913001	PCB BATTERY CASE, CR2032	¥40,000	¥49,000	
8913498001	FUSE 16A T 250V 5X20 H SET OF 10	¥1,100	¥2,000	
9050795001	LAMP UNIT V2	¥670,000	¥685,000	
9050825001	CONTROL UNIT LC 480 ENGAGE FLEX	¥794,000	¥958,000	
9240667001	IPC GENE-BT05 WITHOUT MSATA V2	¥233,000	¥309,000	
9252070001	CONTROLUNIT MP96 HP ENGAGE FLX PRO	¥700,000	¥710,000	
9257403001	THRUST ROD VERTICAL	¥63,500	¥64,000	
9257438001	SHAFT WITH LEVER AND FLAGS	¥94,000	¥137,000	
9257446001	TRANSPORT LOCK SCREW	¥4,200	¥6,000	
9263390001	PCB INDICATOR	¥104,000	¥110,000	
9383310001	TOOL CUBE ADAPTER FOR OWS	¥521,000	¥530,000	
9412883001	ALIGNMENT BAR TRANSFER HEAD	¥26,000	¥33,000	
9412891001	PIPETTER MP24 V2	¥834,000	¥1,020,000	
9412921001	STRIPPER SLEEVE SET OF 8	¥11,600	¥12,000	
9412930001	TOOL CONNECTING SLEEVE HOLDER	¥145,000	¥150,000	
9417249001	PCB Z-AXIS MODULE CONNECTOR	¥42,000	¥58,000	
9417265001	CABLE RIBBON LCD MONITOR	¥15,300	¥16,000	
9420517001	LIGHT SOURCE 2.2 ASSY	¥345,000	¥350,000	
9429123001	COMPRESSION MOUNT COVER HEATER	¥71,000	¥92,000	
9439811001	PCB LED CONTROL RELAX V2	¥27,000	¥34,000	
9586067001	TOOL TORQUE DRIVER 0.1-1.2NM	¥144,000	¥150,000	
9825363001	TOOL TORX BIT TX8	¥5,400	¥6,000	
9888098001	USB FLASH DRIVE 32GB/HW ENC BOOT	¥81,600	¥109,000	
9946764001	TOOL BALANCE OWS 2.0	¥3,863,000	¥4,940,000	
9946772001	CASE INCL. ACCESSORIES OWS 2.0	¥1,090,000	¥1,370,000	
12023989001	TOOTH BELT	¥14,000	¥18,000	
12023997001	WINDOW PHOTOMETER CPL	¥20,000	¥21,000	
26213006001	FUSE 5x20 6.3AT 250V ULR/IEC	¥200	¥1,000	
26350020001	FAN 24VDC 119x119x25	¥43,000	¥65,000	
28006687001	TUBE PE NA D4/6	¥540	¥1,000	
28006695001	TUBE PE NA D6/8	¥400	¥1,000	
28009643001	ADAPTER MALE G1/8 INCH	¥1,100	¥2,000	
28014574001	FILTER ASSEMBLY 120 MM	¥4,900	¥5,000	
28035792001	Filter for Fan 120mm	¥700	¥1,000	
28057168001	FUSE 5x20 1.6AT 250V ULR/IEC	¥200	¥1,000	
28062625001	FUSE 5x20 3.15AT 250V ULR/IEC	¥200	¥1,000	
28095493001	ADAPTER D4/6, MALE G1/8INCH ELB	¥3,600	¥4,000	
28103348001	COUPLING QUICK MALE 4/6MM	¥6,500	¥8,000	
28108285001	CONTAINER 10 L HD-PE YELLOW	¥10,000	¥14,000	
28117802001	MOTOR DC E/24V/90W	¥114,000	¥145,000	
28127328001	Lamp UVC 15W/24VDC	¥3,000	¥4,000	
28133484001	LUBRICANT POLYLUBE GLY 151	¥4,000	¥5,000	
28134847001	ADAPTER FEMALE G1/8 INCH WHITE	¥5,800	¥6,000	
28134863001	ADAPTER FEMALE G1/4 INCH WHITE	¥4,900	¥5,000	
28151989001	CABLE MAINS POWER (EU)	¥4,300	¥5,000	
28166935001	LUBRICANT GREASE HIGH VACUUM DC	¥1,600	¥3,000	
28167826001	MOTOR DC E/24V/90W ASSY MA P	¥130,000	¥165,000	
28174563001	TOOL ADJUSTING BOLT, pack of 2	¥9,000	¥14,000	
28175454001	FITTING TUBE LIGHT	¥37,000	¥38,000	
28181055001	ADAPTER MALE D PLUGGABLE D6	¥7,000	¥10,000	
083.MSP17.20X	Spin Tips 96本	¥3,900	¥4,000	
083.MWP01.20X	2.2 ml 96 deep well plate 100枚	¥77,700	¥79,900	
10010-111	防護服(タイツ)製上下	¥125,900	¥134,800	
10010-113	防護服(タイツ)製上下	¥125,900	¥134,800	
10010-115	防護服(タイツ)製上下	¥125,800	¥134,700	
10010-117	防護服(タイツ)製上下	¥138,200	¥148,000	
10010-119	防護服(タイツ)製上下	¥146,500	¥156,800	
10010-121	防護服(タイツ)製上下	¥141,500	¥151,500	
10010-123	防護服(タイツ)製上下	¥152,800	¥163,600	
10010-136	腕カバー- テレホンアロケティブアルルタイツ	¥86,600	¥92,700	
10010-238	防護服(カバー-オール)付き密封縫製	¥132,000	¥141,300	
10010-240	防護服(カバー-オール)付き密封縫製	¥132,000	¥141,300	
10010-242	防護服(カバー-オール)付き密封縫製	¥132,000	¥141,300	
10010-244	防護服(カバー-オール)付き密封縫製	¥139,900	¥149,800	
10010-246	防護服(カバー-オール)付き密封縫製	¥148,300	¥158,800	
10010-248	防護服(カバー-オール)付き密封縫製	¥157,000	¥168,100	
10010-252	防護服(カバー-オール)付き密封縫製	¥121,800	¥130,400	
10010-254	防護服(カバー-オール)付き密封縫製	¥121,800	¥130,400	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
10010-256	防護服(カバーオール/フード付き密封縫製)	¥120,300	¥128,700	
10010-258	防護服(カバーオール/フード付き密封縫製)	¥130,000	¥139,200	
10010-260	防護服(カバーオール/フード付き密封縫製)	¥137,800	¥147,500	
10010-262	防護服(カバーオール/フード付き密封縫製)	¥144,500	¥154,700	
10010-786	電極 PH ORION GELFILL	¥128,500	¥137,800	
10011-002	マット(保存およびPCR用)	¥120,600	¥127,200	
10011-006	マット Flat Mat 10	¥122,300	¥129,000	
10011-117	シーリングフィルム SealPlate 100	¥53,000	¥55,900	
10011-119	シーリングフィルム ST 100枚入	¥108,100	¥114,000	
10011-122	384-Well Sealing Mat	¥111,600	¥117,800	
10011-124	マット(保存およびPCR用シリコン製)	¥120,200	¥126,800	
10011-126	96ウェルシーリングマット ラウト	¥109,800	¥115,800	
10011-130	マット(保存およびPCR用シリコン製)	¥122,400	¥129,200	
10011-164	PCRプレート 24ウェル(PP製)	¥65,400	¥68,200	
10011-180	PCRプレート 32ウェル(PP製)	¥90,200	¥94,200	
10011-194	PCRプレート 384ウェル(PP製) 10入	¥45,300	¥47,300	
10011-228	PCRプレート 96ウェル Full Skirt	¥248,400	¥259,400	
10011-300	Axy アルチップ 1000ul 1000本x5入	¥185,100	¥196,300	
10011-326	マイクロリットルチップ 200ul ST 960	¥77,100	¥81,800	
10011-488	凍結保存チューブ 0.5ml ST	¥108,400	¥113,200	
10011-508	凍結保存チューブ 0.5ml(自立型)	¥91,900	¥95,900	
10011-510	凍結保存チューブ 0.5ml(自立型)ST	¥103,600	¥108,200	
10011-512	凍結保存チューブ 0.5ml(自立型)	¥100,600	¥105,100	
10011-520	凍結保存チューブ 0.5ml(自立型)	¥103,600	¥108,200	
10011-542	凍結保存チューブ	¥92,400	¥96,500	
10011-544	凍結保存チューブ	¥102,600	¥107,200	
10011-566	凍結保存チューブ 1.5ml(自立型)	¥103,000	¥107,600	
10011-598	凍結保存チューブ 2.0ml エコ グリ	¥97,500	¥101,800	
10011-600	凍結保存チューブ 2.0ml エコ グリ ST	¥102,400	¥107,000	
10011-620	凍結保存チューブ 2.0ml(自立型)	¥87,700	¥91,600	
10011-624	凍結保存チューブ 2.0ml(自立型)	¥95,800	¥100,000	
10011-628	凍結保存チューブ 2.0ml(自立型)	¥97,600	¥102,000	
10011-654	凍結保存チューブ(キャップレス)	¥166,400	¥173,700	
10011-656	凍結保存チューブ(キャップレス)自立 0.5ml	¥193,500	¥202,000	
10011-658	凍結保存チューブ(キャップレス)エコ 0.5ml	¥195,100	¥203,800	
10011-660	凍結保存チューブ(キャップレス)	¥159,900	¥166,900	
10011-662	凍結保存チューブ(キャップレス) 500	¥27,300	¥28,500	
10011-664	凍結保存チューブ(キャップレス)自立 1.5ml	¥183,100	¥191,300	
10011-668	凍結保存チューブ(キャップレス)	¥170,300	¥177,700	
10011-670	凍結保存チューブ(キャップレス)	¥168,700	¥176,100	
10011-672	凍結保存チューブ 自立 2.0ml 500	¥30,000	¥31,300	
10011-674	凍結保存チューブ(キャップレス)エコ 2.0ml	¥189,600	¥197,900	
10011-680	マイクロ遠心チューブ 0.6ml 1000個	¥45,400	¥47,400	
10011-682	マイクロ遠心チューブ 0.6ml ST 500本	¥41,000	¥42,800	
10011-696	マイクロ遠心チューブ 1.5ml(耐熱性) 500本	¥22,300	¥23,300	
10011-700	マイクロ遠心チューブ 1.5ml(耐熱性) 500本	¥22,900	¥23,900	
10011-702	マイクロ遠心チューブ 1.5ml ST 250本	¥25,300	¥26,400	
10011-706	マイクロ遠心チューブ 1.5ml(キャップレス) 500	¥24,800	¥25,900	
10011-718	マイクロ遠心チューブ 0.6ml/1.7ml	¥26,900	¥28,100	
10011-724	マイクロ遠心チューブ 1.7ml 250本	¥25,800	¥27,000	
10011-742	マイクロ遠心チューブ 2.0ml(耐熱性) 500本	¥33,900	¥35,400	
10011-744	マイクロ遠心チューブ 2.0ml 250本	¥25,500	¥26,700	
10011-748	マイクロ遠心チューブ 2.0ml(キャップレス) 500	¥25,200	¥26,300	
10011-760	PCR 0.2ml8連チューブ 詰合パック (Capタ)	¥109,700	¥114,500	
10011-764	PCR 0.2ml8連チューブ NT (Capタ)	¥86,300	¥90,100	
10011-766	PCR 0.2ml8連チューブ グリーン (Capタ)	¥110,300	¥115,100	
10011-774	PCR 0.2ml 8連チューブ 100- (Capタ)	¥107,600	¥112,400	
10011-776	0.2ml PCRチューブ トームアソート 1000	¥62,400	¥65,100	
10011-780	0.2ml PCRチューブ フラット グリ 1000個	¥60,800	¥63,600	
10011-782	PCR 8連チューブ ストリップ 詰合(Tubeタ)	¥24,000	¥25,100	
10011-786	PCR8連チューブ ストリップ グリ 125本	¥25,900	¥27,000	
10011-798	0.2mlPCRチューブ トームアソート 1000	¥63,700	¥66,500	
10011-802	0.2mlPCRチューブ トーム グリ 1000個	¥61,700	¥64,500	
10011-816	0.2mlPCRチューブ キャップレス グリ 1000	¥56,300	¥58,800	
10011-934	Axy 200ul円筒チップ 1000本x20入	¥621,500	¥658,600	
10011-940	96ウェル保存/培養用 96ウェルプレート-1.0mL	¥181,900	¥191,800	
10011-942	96ウェル保存/培養用 96ウェルプレート-1.6mL	¥169,000	¥178,300	
10022-150	Plastic Capillary Tubes 1000	¥126,400	¥138,700	
10024-258	サンプリングテナプレート	¥70,700	¥77,600	
10025-682	VWR Centrifuge Tubes ST 50ml 50	¥24,500	¥28,300	
10029-324	InsulatedEnvelope&GelWrapSet 32	¥136,600	¥153,400	
10029-390	Large NanoCool Shipping system 2入	¥232,700	¥287,200	
10032-802	Flowmi Cell Strainers 40um 50	¥71,300	¥76,300	
10035-976	Poxygrid Bag Holder Kit	¥35,400	¥37,900	
100491-546	Rodent Brain Matrix Mouse 30g	¥679,300	¥650,800	
100492-848	Acrylic RBM Adult Mouse Coronal	¥131,500	¥125,800	
100498-482	血球計算盤(Petroff hausser)	¥887,300	¥973,100	
100498-520	Petroff Hausser Cover Glass 1	¥51,500	¥56,500	
100500-224	Disposable Embedding Molds 288	¥185,300	¥225,600	
100500-232	Glass Coplin Staining Jar 4	¥103,300	¥113,300	
100500-422	Sample Vials w/White Caps ST 500	¥85,500	¥93,700	
100501-004	POLYSINE ADHESIOON SLIDE 72入	¥62,100	¥68,100	
100502-208	Round Cover Slips 12mm	¥129,100	¥141,600	
100503-222	Tri-Pour Beaker Caps 1000mL 250	¥26,600	¥29,200	
10053-766	Handle HDPE 3Ft. 1	¥36,600	¥39,800	
10053-772	Beaker Double Spout HDPE 500ml 1	¥27,100	¥29,500	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
10062-868	Cell Culture Flasks 200	¥112,400	¥118,000	
10062-892	VWR Cell Culture Plate 6well 100	¥82,800	¥86,900	
10062-894	VWR Cell Culture Plate 12well 100	¥96,400	¥101,200	
10062-896	VWR Cell Culture Plate 24well 100	¥110,100	¥115,500	
10062-898	VWR Cell Culture Plate 48well 100	¥127,600	¥133,900	
10062-900	VWR Cell Culture Plate 96well 100	¥132,900	¥139,500	
10081-038	VWR Plant Tissue Culture Container	¥25,200	¥27,100	
10123-792	VWRペトリプレート 1000枚/巻	¥4,300	¥4,700	
10124-110	Volumetric Pipets Class A 17mL 3入	¥74,400	¥80,800	
10124-430	Reclosable Biohazard Bags 1000	¥28,900	¥31,700	
10126-410	VWR Erlenmeyer Flasks 3L Vented	¥123,700	¥133,100	
10126-438	VWR Erlenmeyer Flasks 2L Vented	¥104,900	¥112,900	
10127-460	Lead-Free Autoclave Tape 1個	¥2,800	¥2,900	
10128-916	Polycarbonate Shaker Flask 24	¥148,200	¥159,500	
10128-920	PC製フラスコ 250ml 83x142.3mm 12入	¥83,300	¥98,600	
10128-922	Polycarbonate Shaker Flask 12個	¥83,300	¥89,600	
10128-924	PC製フラスコ 500ml 105x165.3 mm 12入	¥127,400	¥150,700	
10146-290	VWR Ice Buckets with Lids 2.5L 1	¥41,100	¥44,200	
10153-660	Kettle Brush Large	¥86,400	¥93,400	
10153-666	No-Wire Racks 60本立 1個入	¥21,800	¥24,000	
10190T	1000ulマルチチップ LB 1000/bag	¥5,300	¥5,700	
10199-874	Glass Coverslips 30mm 100枚入	¥40,200	¥43,600	
10204-924	CellStrainers for1000ul PipetTips	¥63,600	¥69,300	
102096-460	Accu-Punche 8.0mm 25	¥79,800	¥87,500	
102096-464	Accu-Punche 10.0mm 25	¥105,800	¥116,000	
102096-468	Accu-Punche 12.0mm 25	¥112,000	¥122,800	
10210-300	C-FLEX TUBING Roll of 50	¥173,200	¥185,400	
10210-326	C-FLEX TUBING 50Feet Roll	¥73,000	¥78,100	
102-208-200	ShortReadEliminatorKit XS 15	¥39,800	¥43,500	
102-208-300	Short Read Eliminator Kit 15回用	¥39,800	¥43,500	
102-208-400	ShortReadEliminatorKit XL 15	¥39,800	¥43,500	
102-301-900	Nanobind CBB Big DNA Kit 24回用	¥78,000	¥86,900	
102-302-000	NanobindPlantNucleiBigDNA Kit 24回	¥78,000	¥86,900	
102-302-100	Nanobind Tissue Big DNA Kit 24回	¥78,000	¥86,900	
10230T	1000ul マルチチップ LB S100X10R	¥11,700	¥12,500	
102-762-700	NanoBind HT CBB Kit 96回	¥254,000	¥278,100	
102-762-800	NanoBind HT 1mL Blood Kit 96回	¥508,000	¥556,200	
10279-140	集材(ストレッツ)	¥382,000	¥521,900	
10324-328	Disposable Loop 10uL	¥82,400	¥90,300	
10497-176	Zip Style Reclosable Transport Bag	¥105,500	¥118,500	
10520T	200ul マルチチップ LB S96X10R	¥9,200	¥9,900	
10600-036	Rack no wire cuvette 42穴 8個入	¥14,300	¥15,700	
10600-044	No-Wire Test Tube Half Racks 1	¥14,200	¥15,700	
10717-702	小動物用ケージ(コリアニマルケージ) 1	¥32,400	¥34,000	
10717-704	小動物用ケージリフト(コリアニマルケージ) 1	¥74,900	¥78,600	
10717-715	小動物用ケージ(コリアニマルケージ)	¥269,900	¥283,200	
10718-010	小動物用手術台 28.x11.4cm 1	¥144,400	¥151,500	
10718-012	小動物用手術台 41x21cm 1	¥143,000	¥150,100	
10718-014	Rabbit Surgery Board	¥262,800	¥275,800	
10718-015	Rabbit Surgery Board	¥299,900	¥314,800	
10718-028	小動物注射用固定器(円錐型) 1	¥293,700	¥308,300	
10718-046	Rodent Animal Restrainer	¥81,600	¥85,600	
10718-048	Rodent Animal Restrainer	¥89,800	¥94,200	
10718-052	Rodent Animal Restrainer	¥103,200	¥108,300	
10718-054	小動物用手術台	¥82,300	¥86,300	
10718-056	小動物用固定器具 14x3.8cm	¥89,400	¥93,800	
10718-064	給水ビュ&先管	¥715,800	¥1,264,000	
10718-076	Rodent Animal Restrainer	¥86,200	¥90,400	
10718-330	SIPPER TUBE ECONO 2.5 STRAIGHT	¥6,900	¥7,300	
10754-774	VWR Crystallizing Dishes 340mL台	¥89,400	¥106,900	
10769-152	Stopcock Apiezon Grease H 25g	¥140,700	¥172,800	
10769-154	Grease Apiezon L 25g 1入	¥92,300	¥113,700	
10769-156	Stopcock Apiezon Grease M 25g 1入	¥59,100	¥72,800	
10769-158	Apiezon Grease M 100g 1入	¥104,700	¥129,000	
10769-160	Apiezon Grease N 25g 1入	¥267,600	¥329,600	
10769-164	Apiezon Grease T 25g 1入	¥178,900	¥241,700	
10769-582	Disposable Samplers Pack of 20	¥158,200	¥178,100	
10769-584	Disposable Samplers Pack of 20	¥245,900	¥276,900	
10790-798	Space Saver Vacuum Desiccators	¥79,500	¥85,500	
10798-786	Vapor Line Indicator 100	¥19,100	¥20,500	
10799-210	Disposable Borosilicate Glass 1000	¥224,500	¥269,300	
10803-732	バケツ(樹脂製安全バケツ)2かごのフナシ 1	¥15,100	¥16,100	
10803-992	VWR SpecimenContainers 滅菌 120mL	¥77,700	¥149,100	
10805-018	塗布棒(木製アプリアクター)1000本x20	¥175,700	¥188,100	
10805-154	綿棒(滅菌済み)	¥98,500	¥104,400	
10805-165	Puritan Single Swab 100本	¥15,000	¥17,000	
10805-905	ワイルド(ケイタイEX-L/クリーン用) 119入	¥14,600	¥15,300	
10806-000	Cotton-Tipped Applicators 1000本	¥18,900	¥21,900	
10806-001	Puritan Swabs Cotton Head 152.4mm	¥13,200	¥14,200	
10806-005	Cotton-Tipped Applicators 1000	¥146,200	¥164,400	
10806-020	綿棒L°アスワブ(発泡材、広葉樹材使用)	¥213,700	¥228,800	
10806-021	綿棒L°アスワブ(発泡材、広葉樹材使用)	¥195,400	¥209,200	
10806-023	綿棒L°アスワブ(発泡材、広葉樹材使用)	¥283,400	¥303,300	
10806-024	綿棒L°アスワブ(発泡材、広葉樹材使用)	¥172,500	¥184,600	
10806-025	綿棒L°アスワブ(発泡材、広葉樹材使用)	¥166,600	¥178,300	
10806-026	綿棒L°アスワブ(発泡材、広葉樹材使用)	¥179,900	¥192,600	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
10806-027	綿棒 [®] 17577 [®] (発泡材、広葉樹材使用)	¥176,400	¥188,800	
10806-032	綿棒 [®] 17577 [®] (発泡材、広葉樹材使用)	¥280,000	¥299,700	
10806-286	Animal Ear Punch Thumb-Style 1	¥6,700	¥7,800	
10806-290	Animal Ear Punch, Plier-Style	¥13,000	¥14,800	
10806-294	Animal Ear Punch	¥11,200	¥12,900	
10810-322	インナーローブ [®] クリアBCR(ナイロン/MF製)	¥175,600	¥183,400	
10810-782	CleanWIPE Wipers PK25	¥16,500	¥17,300	
10810-942	MiraWIPE Wipers 150	¥179,400	¥195,700	
10815-702	カバーオールクリーニング GP(ホリワ [®] 0 [®] レ)	¥182,000	¥195,700	
10815-704	カバーオールクリーニング GP(ホリワ [®] 0 [®] レ)	¥138,400	¥153,900	
10815-732	LabCoat Non-Strile 初付 M 25枚入	¥151,800	¥162,500	
10815-734	LabCoat Non-Strile 初付 L 25枚入	¥151,800	¥162,500	
10815-968	カバーオールクリーニング GP(ホリワ [®] 0 [®] レ)	¥195,800	¥210,500	
10815-976	腕カバークリーニング GP(ホリワ [®] 0 [®] レ繊維)	¥84,000	¥84,300	
10816-246	ティッシュ [®] サンプル [®] 115 21cm ST 100	¥86,200	¥92,700	
10821-313	Disposable Samplers sterile 20	¥102,300	¥115,300	
10821-343	Disposable Sampler 100ml 滅菌 20入	¥153,400	¥172,700	
10824-045	透明 [®] 100 [®] (ビニール製/5 [®] 用)	¥233,300	¥273,900	
10824-067	透明 [®] 100 [®] (ビニール製/5 [®] 用)	¥244,000	¥284,800	
10824-089	透明 [®] 100 [®] (ビニール製/5 [®] 用)	¥403,700	¥416,100	
10824-125	透明 [®] 100 [®] (ビニール製/5 [®] 用)	¥204,700	¥235,900	
10825-048	透明 [®] 100 [®] (ビニール製/5 [®] 用)	¥259,500	¥302,200	
10826-007	100 [®] (PE製/ティッシュ [®] サンプル)	¥139,000	¥148,200	
10832-028	カバーオール (クラス100クリーニング対応) ST M	¥936,100	¥884,200	
10832-036	カバーオール (クラス100クリーニング対応) ST 3X	¥1,050,500	¥992,300	
10836-556	バケツ(樹脂製安全バケツ)用 2 [®] 2 [®] 1	¥2,700	¥3,000	
10837-262	クリーニング無塵紙 250	¥619,300	¥663,000	
10837-322	カバーオールクリーニング GP(ホリワ [®] 0 [®] レ)	¥150,400	¥161,000	
10837-326	カバーオールクリーニング GP(ホリワ [®] 0 [®] レ)	¥150,400	¥161,900	
10837-328	カバーオールクリーニング GP(ホリワ [®] 0 [®] レ)	¥182,000	¥195,700	
10837-332	カバーオールクリーニング GP(ホリワ [®] 0 [®] レ)	¥189,900	¥204,200	
10837-376	100 [®] コンフォート [®] (PP製/ティッシュ [®])	¥68,200	¥73,000	
10837-378	カバーオールクリーニング [®] 2 [®]	¥114,400	¥122,800	
10837-380	カバーオールクリーニング [®] 2 [®]	¥114,400	¥122,800	
10837-382	カバーオールクリーニング [®] 2 [®]	¥114,400	¥122,800	
10837-384	カバーオールクリーニング [®] 2 [®]	¥114,400	¥122,800	
10837-386	カバーオールクリーニング [®] 2 [®]	¥108,600	¥120,600	
10837-390	カバーオールクリーニング [®] 2 [®] (3層構造素材)	¥130,700	¥139,900	
10843-053	Bouffant Cap 500入	¥127,500	¥136,500	
10854-122	Vial Rack for Headspace Vials 1	¥22,500	¥24,900	
10854-124	Vial Rack for Headspace Vials 5	¥94,300	¥105,100	
10861-554	VWR Tissue culture plate 6well 100	¥75,200	¥78,900	
10861-558	VWR Tissue culture plate 24well 100	¥78,600	¥82,500	
11000-266	耐熱性 [®] 100 [®] (高温作業/羊毛裏地付)	¥76,000	¥81,400	
11000-270	耐熱性 [®] 100 [®] (高温作業用)	¥94,800	¥103,100	
11000-566	Glove Box Glove 15 [®] 1 [®] 1 [®] 1 [®]	¥264,200	¥285,600	
11000-568	Glove Box Glove 15 [®] 1 [®] 1 [®] 1 [®]	¥276,200	¥298,600	
11000-570	Glove Box Glove 15 [®] 1 [®] 1 [®] 1 [®]	¥180,600	¥193,300	
11000-572	Glove Box Glove 15 [®] 1 [®] 1 [®] 1 [®]	¥260,300	¥281,400	
11000-628	靴カバー [®] ユニバーサル [®] Universal 1	¥13,700	¥15,200	
11000-632	100 [®] (シリコン [®] /4H製 [®] 1 [®] 型) XLサイズ	¥61,000	¥67,500	
11000-634	100 [®] (シリコン [®] /4H製 [®] 1 [®] 型) Lサイズ	¥75,600	¥83,700	
11000-636	100 [®] (シリコン [®] /4H製 [®] 1 [®] 型) Mサイズ	¥63,400	¥70,200	
11000-638	100 [®] (シリコン [®] /4H製 [®] 1 [®] 型) Sサイズ	¥70,600	¥78,200	
11000-640	100 [®] (シリコン [®] /4H製 [®]) 1枚	¥19,200	¥21,300	
11000-642	SilverShield/4H NorthSafety サイズ [®] 11	¥47,400	¥51,300	
11000-644	SilverShield/4H NorthSafety サイズ [®] 10	¥57,000	¥30,700	
11000-646	SilverShield/4H NorthSafety サイズ [®] 9	¥66,000	¥61,700	
11000-648	SilverShield/4H NorthSafety サイズ [®] 8	¥56,400	¥51,700	
11000-650	SilverShield/4H NorthSafety サイズ [®] 7	¥57,200	¥52,300	
11000-652	SilverShield/4H NorthSafety サイズ [®] 6	¥49,400	¥53,500	
11010-084	ハントステルス [®] スチール [®] コイル [®] (1 [®] 1 [®] 1 [®] 付) 12.1cm	¥75,300	¥80,600	
11010-138	Spiral Coated Capsule Sinkers	¥35,200	¥37,600	
11010-142	Spiral No Coated Capsule Sinkers 1	¥23,100	¥24,800	
11010-270	Porous Dissolution Filter 35um 1000	¥579,000	¥619,800	
11010-382	Sinker 3 Prong w/Magnetic Legs 1	¥44,800	¥48,000	
11019-554	Kord Valmark Petri Dishes 100	¥107,600	¥115,800	
11019-556	Kord Valmark Petri Dishes ST 500	¥139,800	¥150,400	
11019-566	Kord Valmark Petri Dishes 500	¥101,200	¥108,900	
11105-050N	11105-050N 111 0.65um/50mm	¥17,600	¥19,800	
11106-047N	CA φ47mm 0.45um 100pc	¥17,600	¥19,800	
11106-050N	11106-050N 111 0.45um/50mm	¥17,600	¥19,800	
11106-90G	低吸着 11106-90G 111 0.45u/90	¥13,800	¥15,600	
11107-047ACN	CA φ47mm 0.2um 100pc	¥27,500	¥28,200	
11107-047N	CA φ47mm 0.2um 100pc	¥17,600	¥19,800	
11150SPL	11150SPL 150mm ST 10x12袋	¥22,600	¥24,400	
11200-000	11200-000 洗浄液 1	¥4,100	¥4,900	
11200-002	11200-002 洗浄液 3.78Lx4本	¥188,900	¥206,300	
11200-008	11200-008 クリーナー(ペットタイプ)	¥63,800	¥77,700	
11200-010	11200-010 クリーナー(ペットタイプ)	¥49,300	¥59,100	
11211-908	11211-908 テープ(クリーニング用/ホリワ [®] 0 [®] 1 [®] 製) 1	¥22,000	¥23,600	
11211-910	11211-910 シールテープ(水容器用/ホリワ [®] 0 [®] 1 [®] 製)	¥406,800	¥436,300	
11211-912	11211-912 テープ(クリーニング用/ホリワ [®] 0 [®] 1 [®] 製) 1	¥12,400	¥13,300	
11211-914	11211-914 テープ(クリーニング用/ビニール製)	¥12,500	¥13,400	
11211-916	11211-916 テープ(クリーニング用/ビニール製)	¥16,400	¥17,600	
11211-922	11211-922 シールテープ(水容器用/ホリワ [®] 0 [®] 1 [®] 製)	¥382,800	¥410,600	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
11211-926	シールテープ (水容器シール用/ホリフン製)	¥382,800	¥410,600	
11211-928	シールテープ (水容器シール用/ホリフン製)	¥382,800	¥410,600	
11211-930	シールテープ (水容器シール用/ホリフン製)1	¥7,200	¥7,700	
11211-932	シールテープ (水容器シール用/ホリフン製)	¥382,800	¥410,600	
11211-934	シールテープ (水容器シール用/ホリフン製)	¥380,900	¥408,600	
11211-936	テープ (クリールム用ビニール製)	¥11,900	¥12,800	
11211-938	テープ (クリールム用ビニール製)DarkBlue	¥11,900	¥14,000	
11211-940	テープ (クリールム用ビニール製)	¥11,900	¥12,800	
11211-942	テープ (クリールム用ビニール製)DarkGreen	¥11,700	¥12,600	
11211-944	テープ (クリールム用ビニール)LightGreen	¥11,900	¥12,800	
11211-946	テープ (クリールム用ビニール)透明	¥11,800	¥12,600	
11211-950	テープ (クリールム用ビニール製)	¥12,100	¥13,000	
11211-952	テープ (クリールム用ビニール製)	¥11,900	¥12,700	
11212-036	クリールムテープ ホリフン 幅5.1cmx55m 3	¥185,700	¥214,900	
11212-098	クリールム用テープ (ワイドテープ) 6本	¥85,000	¥93,700	
11213-021	滅菌用バスタブ (蒸気/ガス滅菌兼用) 200枚	¥17,700	¥18,600	
11213-023	滅菌用バスタブ (蒸気/ガス滅菌兼用) 200枚	¥21,300	¥22,300	
11213-025	滅菌用バスタブ (蒸気/ガス滅菌兼用) 200枚	¥27,600	¥29,000	
11213-027	滅菌用バスタブ (蒸気/ガス滅菌兼用) 200枚	¥53,600	¥58,800	
11213-029	滅菌用バスタブ (蒸気/ガス滅菌兼用) 200枚	¥35,600	¥37,400	
11213-031	蒸気/ガス滅菌用バスタブ 200枚	¥34,200	¥35,900	
11213-033	滅菌用バスタブ (蒸気/ガス滅菌兼用) 200枚	¥47,000	¥49,300	
11213-035	滅菌用バスタブ (蒸気/ガス滅菌兼用) 200枚	¥57,900	¥60,800	
11213-036	滅菌用バスタブ (蒸気/ガス滅菌兼用) 100枚	¥45,900	¥48,100	
11213-038	滅菌用バスタブ (蒸気/ガス滅菌兼用) 100枚	¥59,300	¥62,200	
11213-040	滅菌用バスタブ (蒸気/ガス滅菌兼用) 100枚	¥45,400	¥47,600	
11213-042	滅菌用バスタブ (蒸気/ガス滅菌兼用) 100枚	¥53,400	¥56,100	
11213-044	滅菌用バスタブ (蒸気/ガス滅菌兼用) 100枚	¥112,000	¥117,500	
11213-062	滅菌用バスタブ (蒸気/ガス滅菌兼用) 100	¥61,800	¥64,200	
11213-230	滅菌用バスタブ (蒸気/ガス滅菌兼用)	¥76,300	¥80,000	
11213-236	滅菌用バスタブ (蒸気/ガス滅菌兼用)	¥345,700	¥352,100	
11213-237	滅菌用バスタブ (蒸気/ガス滅菌兼用) 200枚	¥36,500	¥37,900	
11213-241	滅菌用バスタブ (蒸気/ガス滅菌兼用) 200枚	¥29,400	¥30,900	
11213-243	滅菌用バスタブ (蒸気/ガス滅菌兼用) 200枚	¥80,500	¥83,700	
11213-245	滅菌用バスタブ (蒸気/ガス滅菌兼用) 200枚	¥83,500	¥86,800	
11213-268	滅菌用バスタブ (ガス滅菌用)	¥1,115,800	¥1,223,800	
11213-800	Tubular Rollstock 41cm x 38m 1	¥961,400	¥1,080,700	
11214-301	ホリフンバスタブ 118.3ml PK100	¥46,200	¥52,000	
11214-338	ホリフンバスタブ 709.8ml PK50	¥40,300	¥45,300	
11214-367	ホリフンバスタブ 950ml 40入	¥47,000	¥52,800	
11214-469	ホリフンバスタブ 100枚	¥94,400	¥106,100	
11214-491	ホリフンバスタブ 480ml 75入	¥91,800	¥103,200	
11214-527	ホリフンバスタブ 950ml 40入	¥94,400	¥106,100	
11214-708	バイオハザード用廃棄容器 (オートクレーブ用)	¥157,900	¥165,800	
11214-711	バイオハザード用廃棄容器 (オートクレーブ用)	¥152,900	¥168,300	
11214-715	ビペットキラー	¥148,000	¥163,100	
11215-002	バイオバック 1	¥26,500	¥29,000	
11215-202	バイオバック (厚手タイプ)	¥16,700	¥17,800	
11215-235	バイオバック (厚手タイプ)	¥21,400	¥22,800	
11215-246	バイオバック (厚手タイプ)	¥25,200	¥26,900	
11215-257	バイオバック (厚手タイプ)	¥26,400	¥28,200	
11215-280	バイオバック (厚手タイプ)	¥28,100	¥29,900	
11215-302	洗浄瓶 (細口)250mL 6本入	¥20,900	¥22,400	
11215-348	バイオバック (厚手タイプ)	¥39,100	¥41,600	
11215-359	バイオバック (厚手タイプ)	¥45,600	¥49,800	
11215-392	バイオバック (厚手タイプ)	¥50,800	¥54,100	
11215-508	バイオハザードバック用長方形バックホルダー	¥21,800	¥23,200	
11215-550	SCIENCEWARE Flaskup Flask Holder 3	¥19,900	¥21,200	
11215-680	バイオハザードバック (PP製厚さ2ミル)	¥265,400	¥291,100	
11215-684	VWR Biohazard Specimen Bags 500枚	¥67,400	¥74,000	
11215-824	オートクレーブバック (廃棄処理用)	¥196,600	¥215,600	
11215-836	バイオハザードバック 14x33cm	¥175,400	¥187,800	
11215-847	バイオハザードバック用スタック 1個	¥22,900	¥24,500	
11215-853	オートクレーブバック (廃棄処理用)	¥15,600	¥16,700	
11215-855	オートクレーブバック (廃棄処理用)	¥24,200	¥25,900	
11215-857	バイオハザードバック (HDPE製)	¥15,900	¥17,000	
11215-862	バイオハザードバック	¥97,800	¥105,700	
11215-880	脱臭バット	¥119,000	¥127,800	
11215-904	バイオハザードバック (PP製厚さ2ミル)	¥146,300	¥160,400	
11215-983	脱臭バット	¥92,600	¥114,500	
11215-985	脱臭バット	¥92,600	¥99,700	
11215-987	脱臭バット	¥92,200	¥103,600	
11216-000	バイオバック	¥77,300	¥81,000	
11216-002	バイオバック	¥34,200	¥35,800	
11216-012	バイオバック	¥72,000	¥75,600	
11216-056	バイオバック	¥108,300	¥113,700	
11216-089	バイオバック	¥113,300	¥130,800	
11216-090	バイオバック	¥135,300	¥156,200	
11216-092	FilterBags for HomogenizerBlenders	¥279,600	¥293,500	
11216-280	ホリフンバック 15x23cm 500枚入	¥147,600	¥170,400	
11216-373	バイオバック	¥58,800	¥61,700	
11216-409	バイオバック	¥135,300	¥156,200	
11216-420	バイオバック	¥155,000	¥162,700	
11216-422	Whirl-Pak バイオバック 500	¥196,200	¥226,500	
11216-424	ホリフンバック 15x23cm 500枚入	¥142,600	¥163,600	
11216-500	バイオバック	¥222,500	¥233,600	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
11216-515	Whirl-Pak Filter Bags 1.63L	¥276,600	¥318,900	
11216-520	Whirl-Pak Filter Bags 710mL 250	¥187,200	¥196,500	
11216-759	チオ硫酸トリウム入りワールパック 100	¥41,900	¥44,000	
11216-772	ザップパック	¥77,100	¥80,900	
11216-774	ザップパック	¥78,200	¥82,100	
11216-776	ザップパック	¥110,100	¥124,800	
11216-777	ホリパック 19x30cm 500枚入	¥206,900	¥217,200	
11216-780	チオ硫酸トリウム入りワールパック	¥1,189,200	¥1,248,100	
11216-790	ザップパック (環境試料用)	¥125,100	¥131,300	
11216-796	ホリパック 13x19cm 500枚入	¥95,500	¥100,300	
11216-890	ホリパック	¥120,800	¥133,200	
11216-892	ホリパック	¥135,400	¥147,800	
11216-898	ホリパック	¥95,200	¥104,000	
11216-900	ホリパック	¥37,500	¥43,800	
11216-902	ホリパック	¥57,400	¥62,700	
11216-920	ザップパック (Stand-up)710mL 500枚	¥128,200	¥134,500	
11217-082	診断用検体/感染性物質輸送容器	¥129,700	¥147,100	
11217-086	コンテナ (Speci-Freez耐熱輸送容器) 1	¥28,700	¥31,400	
11217-088	ミニスローバル (感染性物質用)	¥7,300	¥8,000	
11217-100	ビトラン液体密閉試験パック 250枚	¥74,600	¥81,400	
11217-102	ビトラン液体密閉試験パック	¥302,600	¥330,400	
11217-104	ビトラン液体密閉試験パック	¥481,200	¥525,500	
11217-106	ビトラン液体密閉試験パック	¥494,300	¥539,900	
11217-108	ビトラン液体密閉試験パック	¥777,000	¥848,600	
11217-110	ビトラン液体密閉試験パック	¥817,500	¥892,900	
11217-114	Infecon3000交換用フィルタ 12入	¥79,100	¥86,400	
11217-120	ビトラン液体密閉試験パック 250枚	¥71,900	¥78,600	
11217-122	ビトラン液体密閉試験パック	¥324,200	¥354,000	
11217-124	ビトラン液体密閉試験パック	¥376,100	¥410,800	
11217-126	ビトラン液体密閉試験パック	¥519,800	¥624,400	
11217-128	ビトラン液体密閉試験パック	¥653,100	¥713,300	
11217-132	Infecon3000 (コンテナ) 用梱包材	¥82,500	¥90,100	
11217-134	Infecon3000 (コンテナ) 用梱包材	¥24,800	¥27,000	
11217-136	Boxes+labels for INF3000 12束入	¥46,900	¥51,200	
11217-138	Infecon3000 (コンテナ) 用容器 1L	¥137,800	¥150,500	
11217-190	Infecon Specimen Transport Bag 250	¥86,100	¥94,100	
11217-192	Infecon Specimen Transport Bag 250	¥128,400	¥140,300	
11217-194	Infecon Specimen Transport Bag 250	¥178,600	¥195,100	
11217-508	パック (シールド付) 100枚x10	¥64,600	¥70,400	
11217-510	パック (シールド付)	¥45,100	¥49,300	
11217-512	パック (シールド付)	¥45,200	¥49,400	
11217-514	パック (シールド付)	¥63,700	¥69,500	
11217-516	パック (シールド付) 100	¥6,000	¥6,500	
11217-524	吸収シート (コンテナ-Infecon3000用) 50ml	¥8,000	¥8,700	
11217-526	吸収シート (コンテナ-Infecon3000用) 100ml	¥10,900	¥11,900	
11217-532	Shipping Declaration Form 150入	¥12,400	¥13,500	
11217-542	スローバル (ドライブ用)	¥4,300	¥4,700	
11217-544	Speci-Freez返却用バル	¥4,500	¥4,900	
11217-548	ビトラン液体密閉試験パック	¥90,100	¥98,400	
11217-550	ビトラン液体密閉試験パック	¥200,500	¥219,000	
11217-552	ビトラン液体密閉試験パック	¥153,400	¥167,600	
11217-554	ビトラン液体密閉試験パック	¥332,700	¥363,700	
11217-556	インフェコン試験輸送パック (ホリパック製)	¥201,500	¥220,100	
11217-558	パック (防護用エアレーション付)	¥24,700	¥26,900	
11217-560	インフェコン試験輸送パック (ホリパック製)	¥87,300	¥94,400	
11217-570	コンテナ (テストテープ・診断検体輸送用)	¥62,600	¥68,300	
11217-620	バル (感染性物質用) 50枚入	¥6,400	¥7,000	
11217-630	Bitran Liquid-tight Formalin Bags	¥219,400	¥239,700	
11217-634	Bitran Liquid-tight Formalin Bags	¥179,400	¥196,000	
11217-636	Infecon Specimen Transport Bag 250	¥145,400	¥158,800	
11217-638	Infecon Specimen Transport Bag 250	¥129,400	¥141,300	
11217-660	コンテナ (検体・感染性物質輸送用)	¥31,400	¥34,400	
11217-664	Infecon5500用交換フィルタ 4入	¥58,300	¥63,600	
11301-025N	CN φ25mm 8um 100pc	¥12,200	¥14,000	
11301-047N	CN φ47mm 8um 100pc	¥17,800	¥20,400	
11302-047ACN	CN φ47mm 3um 100pc	¥27,500	¥28,600	
11302-047N	CN φ47mm 3um 100pc	¥17,800	¥20,400	
11303-047N	CN φ47mm 1.2um 100pc	¥17,800	¥20,400	
11304-047N	CN φ47mm 0.8um 100pc	¥17,800	¥20,400	
11304-050N	CN φ50mm 0.8um 100pc	¥17,800	¥20,400	
11306-013N	CN φ13mm 0.45um 100pc	¥11,800	¥13,800	
11306-025N	CN φ25mm 0.45um 100pc	¥12,200	¥14,000	
11306-047N	CN φ47mm 0.45um 100pc	¥17,800	¥20,400	
11306-100G	ニトリル-スチレン113 0.45um/100mm 25	¥14,600	¥16,800	
11306-25-ACN	CN スリットフィルタ φ25mm 0.45um 100	¥16,000	¥16,600	
11310-688	ポリ口瓶PBT CAP	¥56,100	¥60,700	
11310-690	ポリ口瓶PBT CAP GL32 10個入	¥30,700	¥33,200	
11310-694	ポリ口瓶シリコンキャップ	¥46,300	¥50,000	
11310-696	Silicone Septa GL32 10個入	¥15,600	¥16,300	
11311-162	TruCLEAN Mops (ステルスバル) 1本	¥115,400	¥123,600	
11311-648	Vial 40cc 24mm GREEN 144	¥137,100	¥151,900	
11311-656	ガラスバール用スリットバル 8OZ 24	¥32,300	¥34,800	
11311-658	ガラスバール用スリットバル 16OZ 24	¥41,100	¥45,500	
11376-412	効能補給剤	¥47,100	¥51,000	
11376-414	効能補給剤	¥46,400	¥50,200	
11376-416	効能補給剤	¥58,900	¥63,700	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
11376-418	効補剤	¥56,600	¥61,300	
11403-047ACN	CN W/B φ47mm 1.2um 100pc	¥19,000	¥19,800	
11404-80ALN	ニトロセルロースメンブレン 0.8um/80mm 100	¥57,000	¥59,000	
11406-047ACN	CN W/B φ47mm 0.45um 100pc	¥19,000	¥19,800	
11803-013N	PTFE φ13mm 1.2um 100pc	¥19,000	¥22,000	
11803-100G	PTFE φ100mm 1.2um 25pc	¥32,000	¥36,800	
11806-047N	PTFE φ47mm 0.45um 100pc	¥28,000	¥32,800	
11842-050N	PTFE φ50mm 5um 100pc	¥28,000	¥32,800	
12000-206	PCR 0.2ml8連フューネル 125	¥90,900	¥95,000	
12000-568	防護服(タイツ製/フード&フック付き)	¥143,100	¥153,200	
12000-570	防護服(タイツ製/フード&フック付き)	¥143,100	¥153,200	
12000-572	防護服(タイツ製/フード&フック付き)	¥153,000	¥163,700	
12000-574	防護服(タイツ製/フード&フック付き)	¥151,700	¥162,500	
12000-576	防護服(タイツ製/フード&フック付き)	¥160,800	¥172,200	
12000-578	防護服(タイツ製/フード&フック付き)	¥173,700	¥186,000	
12175-001	秤量皿(アルミ)	¥97,600	¥112,400	
12177-266	秤量皿(アルミ)	¥49,500	¥59,800	
12177-299	秤量皿(アルミ)	¥57,200	¥62,700	
12177-301	秤量皿(アルミ)60ml 50枚	¥12,100	¥14,400	
12177-305	秤量皿(アルミ)	¥70,400	¥85,000	
12243-062	Transparent Poly Tubing Rolls	¥243,100	¥267,900	
12243-068	Transparent Poly Tubing Rolls	¥311,500	¥343,300	
12243-126	Pink Antistatic Poly Tubing Rolls	¥248,000	¥273,300	
12243-134	Pink Antistatic Poly Tubing Rolls	¥577,700	¥636,600	
12243-246	Breather Bags	¥266,700	¥293,900	
12243-278	Breather Tubing Rolls 30.5cm x 30m	¥93,700	¥103,300	
12444-044	テラカ(ワット温度計)	¥175,500	¥190,700	
12444-120	テラカ(温度計)	¥25,000	¥29,400	
12567-330	天秤台	¥798,500	¥1,444,000	
12570-378	卓上除震台	¥747,700	¥1,386,000	
12574-064	ヘッドアップ顕微鏡	¥156,000	¥173,600	
12576-934	Disposable Polystyrene Forceps 100	¥138,700	¥78,000	
12577-060	秤量皿(アルミ/超小型)	¥41,000	¥44,300	
12577-062	秤量皿(アルミ/超小型)	¥21,800	¥23,900	
12577-064	秤量皿(アルミ/超小型)	¥45,700	¥50,100	
12577-066	秤量皿(アルミ/超小型)250枚	¥42,300	¥51,100	
12577-068	秤量皿(アルミ/超小型)	¥20,300	¥21,900	
12577-070	秤量皿(アルミ/超小型)	¥48,300	¥58,300	
12577-072	秤量皿(アルミ/超小型)0.12mL	¥50,600	¥55,400	
12577-081	VWR Aluminium Crinkle Dish(付)	¥50,800	¥58,400	
12578-103	Glass Fiber Pad 1	¥740,800	¥808,600	
12578-121	VWR計量紙 7.6x7.6cm 500枚入	¥17,800	¥19,200	
12578-165	VWR計量紙 10.2x10.2cm 500枚入	¥16,600	¥18,000	
12578-201	VWR計量紙 15.2x15.2cm 500枚入	¥16,800	¥18,200	
12602-47ALK	ゼラチン φ47mm 3um 50pc	¥39,000	¥43,500	
12602-47ALN	ゼラチン φ47mm 3um 100pc	¥66,500	¥74,200	
12602-50ALK	ゼラチン φ50mm 3um 50pc	¥39,000	¥43,500	
12602-80ALK	ゼラチンフィルター 50枚	¥43,500	¥48,500	
12620-874	Microtube&Microplate Holder Kit	¥64,600	¥69,700	
12620-878	Microplate Holder Set of 2	¥37,300	¥39,100	
12621-138	VWR Modular Heating Block 6穴	¥82,900	¥149,000	
12621-246	VWR Talon S-Connector(スルス)1個	¥19,900	¥20,900	
12777-200	Accessories for BRAND seripettor	¥115,400	¥124,800	
12777-838	トータル温度計(冷蔵/冷凍庫用)	¥55,700	¥61,400	
12985-068	VWR Taron Support Stands	¥36,700	¥38,500	
13004-047N	CN B/W φ47mm 0.8um 100pc	¥18,500	¥19,200	
13005-047ACN	CN B/W φ47mm 0.65um 100pc	¥19,000	¥19,800	
13006-025N	格子入メンブレンフィルター-130 25mm0.45um 10	¥13,000	¥13,500	
13006-047ACN	CN B/W φ47mm 0.45um 100pc	¥19,000	¥19,800	
13196-048	実験台保護シート(吸水性)	¥60,400	¥67,800	
13196-058	Glascrite Pen	¥24,100	¥25,600	
13196-070	Plastic Caliper with Dial	¥78,700	¥142,000	
13197-116	バレルドック 61x76cm 100枚入	¥96,600	¥103,400	
13197-118	バレルドック 100枚	¥115,600	¥123,800	
13197-128	SCIENCEWARE Micro Funnels	¥15,300	¥16,500	
13197-150	バレルドック	¥132,300	¥146,200	
13197-170	Centrifuge Tubes Rimless 45mL	¥25,500	¥27,400	
13197-176	Tissue Grinder-Conical 6個入	¥13,100	¥14,100	
13197-284	廃液保存ボトル(HDPE製)	¥149,900	¥160,400	
13197-366	Large Weighted Safety Shields 76.2	¥510,600	¥568,100	
13197-374	Quaverette Spatula	¥41,900	¥45,300	
13197-390	薬さし	¥25,400	¥26,900	
13197-392	薬さし	¥27,400	¥29,000	
13197-394	薬さし	¥294,200	¥311,900	
13197-398	薬さし	¥56,900	¥60,300	
13197-414	StainlessSteelLandle&Reservoir 1	¥335,200	¥359,500	
13197-506	ワイルドシステム(ウェットスク/水性薬液使用)	¥115,100	¥127,200	
13229E	ヒパルシドフィルター 5000MWCO CTA12入	¥15,800	¥16,400	
13239E	ヒパルシドフィルター 10000MWCO CTA12入	¥15,800	¥16,400	
13249E	ヒパルシドフィルター 20000MWCO CTA12入	¥15,800	¥16,400	
13270-650	ウォーター用温度プローブ/冷却コイル	¥142,300	¥259,000	
13270-655	ウォーター用温度プローブ/冷却コイル	¥153,700	¥280,000	
13270-660	ウォーター用温度プローブ/冷却コイル	¥181,100	¥329,000	
13270-792	シリコバス	¥811,800	¥1,477,000	
13270-794	シリコバス	¥820,700	¥1,493,000	
13270-868	ウォーター用温度プローブ/冷却コイル	¥100,000	¥182,000	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
13271-060	恒温槽 室温+5℃-70℃	¥852,800	¥1,551,000	
13271-500	クーラー(循環型/投込式)	¥1,690,600	¥3,075,000	
13271-502	クーラー(循環型/投込式)40℃-40℃	¥1,197,600	¥3,093,000	
13271-504	クーラー(循環型/投込式)40℃-40℃	¥1,657,700	¥3,015,000	
13272-034	シリカガラス	¥132,000	¥240,000	
13279E	ヒパルシ°ンフィルター 300000MWCO PES12入	¥15,800	¥16,400	
13303-730	ウォーターバス(振盪機能付)	¥94,800	¥173,000	
13304-100	ウォーターバス(振盪機能付)	¥81,900	¥149,000	
13304-132	ウォーターバス(振盪機能付)	¥40,600	¥74,000	
13304-165	ウォーターバス(振盪機能付)	¥81,300	¥148,000	
13304-223	ウォーターバス(振盪機能付)	¥65,100	¥119,000	
13304-256	ウォーターバス(振盪機能付)	¥61,900	¥113,000	
13304-289	ウォーターバス(振盪機能付)	¥102,300	¥186,000	
13306-181	ウォーターバス(振盪機能付)	¥218,600	¥398,000	
13306-202	ウォーターバス(振盪機能付)	¥118,800	¥217,000	
13306-203	ウォーターバス(振盪機能付)	¥167,500	¥305,000	
13306-309	ウォーターバス(振盪機能付)	¥421,400	¥767,000	
13307-196	循環恒温槽(糞便性大腸菌検査用)	¥125,300	¥229,000	
13307-425	循環恒温槽(糞便性大腸菌検査用)	¥278,000	¥507,000	
13330-187	循環恒温槽(デジタル式過熱循環型)	¥328,700	¥598,000	
13400-047Q	GF φ47mm 500pc	¥24,800	¥28,800	
13440-047Q	GFnb φ47mm 500pc	¥52,000	¥57,200	
13450-196	耐熱性ガラス(Temp-Tec) サイズ9	¥45,300	¥48,500	
13450-310	シリカガラスコートケール付ヒューリアクワイア	¥162,300	¥176,300	
13452-318	DuPontフイルコート(2°ケット) XL 30枚入	¥83,000	¥88,800	
13453-578	防護服(タイツ製カール/ゴム入り)	¥159,800	¥171,100	
13485-118	ウォーターバス(リング付き) 40℃-100℃	¥1,070,900	¥1,948,000	
13485-119	ウォーターバス(リング付き) 40℃-100℃	¥1,606,800	¥2,923,000	
13500-390	カンパルバッグ(Stand-up)118mL 500枚	¥70,700	¥74,200	
13500-392	カンパルバッグ(Stand-up)532mL 500枚	¥106,900	¥120,200	
13500-422	7700ルーフマスク(バックル付)M	¥91,100	¥99,700	
13500-498	ThermoSafe VIP Insulated Shipper	¥242,600	¥260,000	
13500-582	ThermoSafe VIP Insulated Shipper	¥165,700	¥177,500	
13500-857	保護メガネ Uvextreme 10	¥86,400	¥93,100	
13500-922	溶存酸素計 500	¥148,400	¥273,000	
13501-912	VWR Tray Pad Set	¥66,800	¥71,100	
13501-946	ガラス片 100mL 12入	¥75,500	¥81,600	
13501-966	溶存酸素計(マルチター)	¥827,800	¥1,521,000	
13501-968	溶存酸素計(マルチター)	¥685,100	¥1,258,000	
13501-970	溶存酸素計(マルチター)	¥443,600	¥660,000	
13501-972	溶存酸素計(マルチター)	¥741,700	¥1,362,000	
13501-974	溶存酸素計(マルチター)	¥264,600	¥486,000	
13501-976	溶存酸素計(マルチター)	¥200,000	¥368,000	
13501-980	溶存酸素計(マルチター)	¥19,300	¥35,000	
13501-986	溶存酸素計(マルチター)	¥19,700	¥37,000	
13502-131	安全ヘルメット	¥267,500	¥296,300	
13502-139	安全ヘルメット	¥184,300	¥54,700	
13502-141	安全ヘルメット	¥184,300	¥54,700	
13502-182	ヘルメットライナー	¥81,100	¥45,600	
13502-193	安全ヘルメット	¥134,200	¥148,700	
13502-200	安全ヘルメット	¥89,800	¥48,400	
13502-274	フェイスシールド	¥10,900	¥12,100	
13502-275	フェイスシールド	¥28,600	¥31,800	
13502-278	フェイスシールドアーム	¥6,200	¥6,900	
13502-281	フェイスシールドアーム	¥200,100	¥222,300	
13503-188	UNIPLATE Microplates 10mL 25枚	¥204,000	¥242,900	
13503-204	Microplate 48well PP 25枚	¥211,300	¥231,200	
13890-024	ビーカー(ホリワ®)1L	¥59,200	¥69,400	
13890-046	ビーカー(ホリワ®)1L	¥70,200	¥82,300	
13890-068	ビーカー(ホリワ®)1L	¥97,200	¥125,400	
13890-080	ビーカー(ホリワ®)1L	¥107,000	¥125,400	
13890-104	ビーカー(ホリワ®)1L	¥156,400	¥183,500	
13890-126	ビーカー(ホリワ®)1L	¥140,600	¥164,900	
13890-148	ビーカー(ホリワ®)1L	¥86,900	¥101,900	
13890-160	ビーカー(ホリワ®)1L	¥113,200	¥132,800	
13890-181	ビーカー(ホリワ®)1L 4000ml 4個入	¥113,400	¥132,900	
13912-604	YREX ビーカー 1L 6本	¥123,600	¥133,600	
13914-177	ビーカー(PTFE)	¥72,100	¥77,600	
13914-199	ビーカー(PTFE)	¥87,400	¥101,300	
13914-213	ビーカー(PTFE)	¥97,600	¥109,200	
13914-235	ビーカー(PTFE)	¥139,500	¥156,000	
13914-257	ビーカー(PTFE)	¥137,800	¥154,100	
13914-279	ビーカー(PTFE)	¥262,700	¥293,900	
13914-290	ビーカー(PTFE)	¥283,300	¥317,000	
13915-985	ビーカー(ホリワ®)	¥45,100	¥54,400	
13915-996	ビーカー(ホリワ®)	¥54,800	¥65,900	
13916-004	ビーカー(ホリワ®)	¥69,200	¥83,100	
13916-015	ビーカー(ホリワ®)	¥56,400	¥61,800	
13916-037	ビーカー(ホリワ®)	¥36,200	¥39,700	
13916-060	ビーカー(ホリワ®)	¥33,500	¥40,400	
13916-707	ビーカー(ハンドル付き)	¥12,600	¥13,900	
13916-730	ビーカー(ハンドル付き)	¥12,400	¥13,600	
13916-762	ビーカー(ハンドル付き)	¥12,400	¥13,600	
13916-795	ビーカー(ハンドル付き)	¥23,300	¥25,700	
13916-820	ビーカー(ハンドル付き)	¥36,600	¥40,300	
13916-830	36本掛けチューブホルダー	¥33,600	¥61,000	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
13916-832	10本(15ml)+16本(5ml)掛けチューブ 緑	¥33,400	¥61,000	
13916-834	6本掛けチューブ 緑	¥33,600	¥61,000	
13917-580	ビーカー(テロ) 50mL	¥49,800	¥55,800	
13917-582	ビーカー(テロ) 125mL	¥91,200	¥102,100	
13917-584	ビーカー(テロ) 250mL	¥139,800	¥156,400	
13917-586	ビーカー(テロ) 500mL	¥147,700	¥165,200	
13920-012	KIMAX Berzelius TM -カ- 100mL 12本	¥67,100	¥79,000	
13975-025	ビーカー(ステン)ス	¥28,500	¥33,400	
13975-081	ビーカー(ステン)ス	¥25,200	¥31,600	
13977-010	ビーカー(ステン)ス	¥37,400	¥44,800	
13977-042	ビーカー(ステン)ス	¥39,800	¥46,800	
13977-064	ビーカー(ステン)ス	¥42,000	¥49,800	
13977-097	ビーカー(ステン)ス	¥52,600	¥61,800	
13977-122	ビーカー(ステン)ス	¥78,200	¥91,800	
13977-144	ビーカー(ステン)ス	¥79,100	¥94,700	
13977-166	ビーカー(ステン)ス	¥86,500	¥101,700	
13978-067	ビーカー(ハンドル付き ステン)ス	¥89,100	¥104,600	
13978-089	ビーカー(ハンドル付き ステン)ス	¥92,500	¥108,600	
13980T	10ul Eチップ LB 1000/bag	¥4,600	¥4,900	
13983-023	ビーカー ハンドルつき(ステン)ス	¥58,100	¥68,300	
13983-067	ビーカー ハンドルつき(ステン)ス	¥59,700	¥70,100	
13983-089	ビーカー ハンドルつき(ステン)ス	¥104,400	¥122,700	
14002-650	ワイパー(ワイパー-ル30Iコマ付)	¥79,700	¥83,700	
14002-957	TruCLEAN Mop(Swivel Mop Frame)	¥128,200	¥134,600	
14003-078	温度計 Type-K Datalogging	¥340,500	¥369,900	
14003-580	シーカー	¥230,500	¥427,000	
14004-140	クリーナー CIP300 5kg	¥203,600	¥218,500	
14007-066	VistaLabリサ TM -ル TM -(7kg付)ST 100ml 50枚	¥84,900	¥91,600	
14100-404	発泡スロ-ルコンテナ(ホ-リスポンフォーム・蓋付)12	¥72,000	¥78,900	
14100-410	発泡スロ-ルコンテナ(ホ-リスポンフォーム・蓋付)	¥125,100	¥130,000	
14100-412	発泡スロ-ルコンテナ(ホ-リスポンフォーム・蓋付)	¥21,600	¥22,400	
14100-414	発泡スロ-ルコンテナ(ホ-リスポンフォーム・蓋付)	¥21,400	¥22,200	
14100-428	輸送コンテナ(厚み4.4cm)6個入	¥82,900	¥86,100	
14100-430	発泡スロ-ルコンテナ(ホ-リスポンフォーム・蓋付)	¥136,700	¥142,000	
14100-432	発泡スロ-ルコンテナ/ホ-リスポンフォーム/7kg付 4	¥78,200	¥81,200	
14100-462	Insulated Shipper case of 3	¥177,900	¥190,600	
14100-464	ThermoSafe Insulated Shippers CS 2	¥263,800	¥282,700	
14100-466	ThermoSafe Insulated Shippers 1	¥77,000	¥82,500	
14100-468	ThermoSafe Insulated Shippers	¥119,500	¥128,000	
14100-478	PolarPack保冷剤 24個入	¥10,400	¥11,200	
14183-500	カスケット	¥77,000	¥83,600	
14183-521	カスケット	¥108,100	¥117,400	
14183-565	カスケット	¥121,600	¥167,100	
14200T	1000ul マルチフィルターチップ LB S100x10R	¥19,600	¥21,200	
14202-676	Dounce Tissue Grinders 2mL	¥127,800	¥139,500	
14203-100	Bottle Specific Gravity 25ml	¥109,000	¥120,100	
14203-342	Bottle Specific Gravity 100ml	¥133,000	¥146,600	
14204-636	Bottle Specific Gravity 50ml	¥128,300	¥141,400	
14210-316	Separatory Funnels Squibb 125mL 1	¥83,800	¥91,500	
14210-386	Bottle Specific Gravity 5mL	¥100,800	¥124,000	
14210-660	Bottle Specific Gravity 10ml	¥103,400	¥114,000	
14210T	20ul マルチフィルターチップ LB S96x10R	¥18,400	¥19,800	
14215-224	シーカー	¥40,700	¥74,000	
14215-226	シーカー	¥65,000	¥115,000	
14215-228	シーカー	¥83,700	¥551,000	
14215-230	シーカー	¥98,400	¥177,000	
14215-232	シーカー	¥109,000	¥196,000	
14215-236	シーカー	¥75,800	¥136,000	
14215-238	シーカー	¥85,600	¥154,000	
14215-240	シーカー	¥247,700	¥448,000	
14215-242	シーカー	¥243,200	¥440,000	
14215-244	シーカー	¥231,100	¥420,000	
14215-760	KIMTECH PURE G3初作トリル L 100	¥31,900	¥34,100	
14215-762	KIMTECH PURE G3初作トリル M 100	¥31,900	¥34,100	
14215-764	KIMTECH PURE G3初作トリル S 100	¥31,900	¥34,100	
14215-766	KIMTECH PURE G3初作トリル XL 100	¥31,900	¥34,100	
14215-768	KIMTECH PURE G3初作トリル XS 100	¥31,900	¥34,100	
14216-006	保護カ TM genesi (UV保護用)	¥79,000	¥82,900	
14216-236	褐色ガラス瓶 100ml	¥77,500	¥84,500	
14218-194	真空 TM (オレリ)	¥1,011,700	¥1,840,000	
14220-000	バ TM (PP製厚さ1.5mm) 100	¥16,400	¥18,000	
14220-002	バ TM (PP製厚さ1.5mm) 100	¥30,400	¥33,300	
14220-004	バ TM (PP製厚さ1.5mm) 100	¥74,700	¥81,900	
14220-006	バ TM (PP製厚さ1.5mm) 100	¥118,700	¥130,100	
14220-008	バ TM (PP製厚さ1.5mm)	¥288,900	¥316,900	
14220-010	バ TM (PP製厚さ1.5mm)	¥282,900	¥310,300	
14220-012	バ TM (PP製厚さ1.5mm)	¥334,000	¥366,300	
14220-014	バ TM (PP製厚さ1.5mm)	¥428,900	¥470,400	
14220-016	バ TM (PP製厚さ1.5mm) 100	¥10,100	¥11,100	
14220-018	バ TM (PP製厚さ1.5mm) 100	¥12,400	¥13,600	
14220-020	バ TM (PP製厚さ1.5mm)	¥271,100	¥297,300	
14220-022	バ TM (PP製厚さ1.5mm)	¥313,400	¥343,700	
14220-024	バ TM (PP製厚さ1.5mm)	¥356,700	¥391,200	
14220-026	バ TM (PP製厚さ2mm)	¥94,900	¥104,100	
14220-028	バ TM (PP製厚さ2mm)	¥70,600	¥77,500	
14220-030	バ TM (PP製厚さ2mm)	¥143,200	¥168,800	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
14220-032	バ イザ ト バ ック (PP製厚さ2ミル)	¥272,300	¥298,600	
14220-034	バ イザ ト バ ック (PP製厚さ2ミル)	¥158,500	¥173,800	
14220-036	オートクレープ バ ック (廃棄処理用)20.3x30.5	¥159,600	¥175,100	
14220-038	オートクレープ バ ック (廃棄処理用)35.6x48.3	¥122,300	¥134,100	
14220-040	オートクレープ バ ック (廃棄処理用)48.3x58.4	¥166,400	¥182,500	
14220-042	オートクレープ バ ック (廃棄処理用)63.5x88.9	¥294,100	¥322,600	
14220-044	オートクレープ バ ック (廃棄処理用)	¥247,300	¥271,200	
14220-046	バ イザ ト バ ック (PP製厚さ2ミル)	¥153,400	¥168,300	
14220-048	バ イザ ト Bag (PP/厚2ミル)35.6x48.3cm	¥74,000	¥81,200	
14220-050	バ イザ ト バ ック (PP製厚さ2ミル)	¥153,800	¥168,700	
14220-052	バ イザ ト バ ック (PP製厚さ2ミル)	¥238,400	¥261,500	
14220-054	バ イザ ト バ ック (PP製厚さ2ミル)	¥269,500	¥295,500	
14220-066	バ イザ ト バ ック (PP製厚さ2ミル)	¥116,900	¥128,200	
14220-068	バ イザ ト バ ック (PP製厚さ2ミル)	¥148,900	¥163,300	
14220-070	バ イザ ト バ ック (PP製厚さ2ミル)	¥150,700	¥165,300	
14220-082	バ イザ ト バ ック (PP製厚さ3ミル)	¥124,200	¥136,300	
14220-084	バ イザ ト バ ック (PP製厚さ3ミル)	¥273,300	¥299,800	
14220-086	バ イザ ト バ ック (PP製厚さ3ミル)	¥310,100	¥340,100	
14220-088	バ イザ ト バ ック (PP製厚さ3ミル)	¥127,100	¥139,500	
14220-090	バ イザ ト バ ック (PP製厚さ2ミル)	¥100,400	¥110,100	
14220-092	バ イザ ト バ ック (PP製厚さ2ミル)	¥229,400	¥251,600	
14220-094	バ イザ ト バ ック (PP製厚さ2ミル)	¥271,500	¥297,800	
14220-096	バ イザ ト バ ック (PP製厚さ2ミル)	¥271,500	¥297,800	
14220-098	バ イザ ト バ ック 96.5x121.9cm	¥209,600	¥229,900	
14220-100	バ イザ ト バ ック (PP製厚さ2ミル)	¥81,700	¥89,600	
14220-102	バ イザ ト バ ック (PP製厚さ2ミル)	¥211,900	¥232,500	
14220-104	バ イザ ト バ ック (PP製厚さ2ミル)	¥243,800	¥267,400	
14220-106	バ イザ ト バ ック (PP製/焼却処理用)	¥141,500	¥155,200	
14220-108	バ イザ ト バ ック (PP製/焼却処理用)	¥138,700	¥152,100	
14220-110	バ イザ ト バ ック (PP製/焼却処理用)	¥142,100	¥155,800	
14220-112	バ イザ ト バ ック (PP製/焼却処理用)	¥220,200	¥241,600	
14220-114	バ イザ ト バ ック (PP製/焼却処理用)	¥201,100	¥220,600	
14220-729	VistaLab Ovation Tip 1400ul 1000本	¥48,300	¥51,700	
14220-735	Tip Ovation 35ul バ ル 1000本入	¥39,400	¥42,100	
14220-737	Tip Ovation 200ul バ ル 1000本入	¥38,800	¥41,500	
14220T	200ul マルチフィルチップ LB S96x10R	¥18,700	¥20,100	
14221-792	バ ック(再密閉可能)	¥21,500	¥23,100	
14221-806	Weighing Funnels 20個	¥72,400	¥77,900	
14221-808	Weighing Funnels 50個	¥76,500	¥82,300	
14221-810	Weighing Funnels 50個	¥57,800	¥62,200	
14221-826	Weighing Funnels 6個	¥73,600	¥79,100	
14221-972	ホルテックス エ-用ア ットホ-ムハット	¥37,400	¥69,000	
14222-300	CURITY GauzeSpongesST 8-Ply 3000枚	¥203,300	¥217,700	
14224-534	Immerision General-Purpos Probe	¥49,300	¥53,600	
14226-876	防火キャビネット(シングルドア)45.4L	¥834,300	¥1,518,000	
14226-878	防火キャビネット(シングルドア)45.4L	¥975,000	¥1,774,000	
14226-880	防火キャビネット(シングルドア)60.6L	¥968,700	¥1,762,000	
14226-882	防火キャビネット(シングルドア)60.6L	¥1,116,600	¥2,031,000	
14227-676	Exhaust Filter(真空ホ-ウ用) 1個	¥285,500	¥529,000	
14228-062	ユニラック(UniRack)5色x各2個	¥36,000	¥38,800	
14228-076	Lab System w/racks 3000本	¥4,216,200	¥9,960,000	
14230-628	真空ホ-ウ (オイルフリー)	¥1,293,800	¥2,353,000	
14231-358	Tenbroeck Tissue Grinders 1mL 2本	¥141,600	¥154,400	
14231-364	GrindingChambar & SerratedPestle2ml	¥187,800	¥204,800	
14232-324	TITAN2 Syringe Filters 30mm/PTFE 5	¥886,700	¥963,300	
14232-384	SunFlo Dissolution Filter 1000個	¥738,400	¥798,200	
14232-798	ガラスバ イル 5x32mm クラ 1000本	¥340,700	¥368,300	
14233-756	ワイル-(82003-820)専用テ イス ッカー	¥20,000	¥21,600	
14301-110-K	空のア-レート 50	¥845,000	¥921,100	
14320-110ACD	トリアイ寒天培地 (φ110mm、10個入)	¥8,500	¥9,500	
14321-110ACD	ワ-寒天培地 (φ110mm、10個入) ④	¥8,500	¥9,500	
14555-300	フェア-(産業用)	¥408,200	¥437,100	
14555-328	フェア-(産業用加-ムめつき仕上げ)	¥412,100	¥444,800	
14670-001	トランスファ-ヒ-ア-ット(ス-イト) 500本	¥31,000	¥34,400	
14670-003	トランスファ-ヒ-ア-ット(ス-イト) 500	¥44,700	¥49,500	
14670-006	トランスファ-ヒ-ア-ット(ス-イト) 500	¥40,900	¥45,300	
14670-103	トランスファ-ヒ-ア-ット(ス-イト) 500	¥31,000	¥34,300	
14670-114	トランスファ-ヒ-ア-ット(ス-イト) 500	¥39,900	¥44,200	
14670-147	トランスファ-ヒ-ア-ット(ス-イト)	¥209,200	¥231,800	
14670-149	トランスファ-ヒ-ア-ット(ス-イト)	¥157,200	¥174,200	
14670-169	トランスファ-ヒ-ア-ット(ス-イト/細口タイプ) 500	¥37,600	¥41,700	
14670-170	トランスファ-ヒ-ア-ット(ス-イト/細口タイプ)	¥244,300	¥270,600	
14670-180	トランスファ-ヒ-ア-ット(ス-イト/細口タイプ)	¥61,200	¥67,800	
14670-191	トランスファ-ヒ-ア-ット(ス-イト/細口) 500本	¥43,700	¥48,500	
14670-196	トランスファ-ヒ-ア-ット(ス-イト/細口タイプ)	¥218,100	¥242,000	
14670-200	トランスファ-ヒ-ア-ット(ス-イト) 500入	¥28,000	¥31,000	
14670-201	トランスファ-ヒ-ア-ット(ス-イト)	¥205,300	¥227,200	
14670-205	トランスファ-ヒ-ア-ット(ス-イト) 500本	¥28,900	¥32,000	
14670-250	トランスファ-ヒ-ア-ット(ス-イト) 500	¥52,800	¥58,600	
14670-282	トランスファ-ヒ-ア-ット(ス-イト) ST 500入	¥64,600	¥71,500	
14670-307	トランスファ-ヒ-ア-ット(ス-イト) 500入	¥25,100	¥27,900	
14670-345	トランスファ-ヒ-ア-ット(ス-イト/細口タイプ) 500	¥26,700	¥29,500	
14670-364	トランスファ-ヒ-ア-ット(ス-イト) 500	¥44,400	¥49,200	
14670-366	トランスファ-ヒ-ア-ット(ス-イト)	¥203,400	¥225,200	
14670-398	トランスファ-ヒ-ア-ット(ス-イト)400本入	¥79,100	¥87,700	
14670-399	トランスファ-ヒ-ア-ット(ス-イト)	¥208,600	¥230,800	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
14670-400	トランスアピロト(スプレッド) 100本	¥72,600	¥80,400	
14670-401	トランスアピロト(スプレッド)22.9cm 400本	¥43,900	¥48,700	
14670-405	トランスアピロト(スプレッド) 100本	¥63,500	¥70,300	
14670-417	トランスアピロト(スプレッド)	¥133,700	¥163,000	
14670-418	トランスアピロト(スプレッド)	¥295,400	¥327,300	
14670-419	トランスアピロト(スプレッド)	¥190,800	¥211,400	
14670-420	トランスアピロト(スプレッド)	¥218,200	¥241,500	
14670-430	トランスアピロト(スプレッド/ボトル付)	¥179,100	¥198,200	
14670-439	トランスアピロト(スプレッド/ボトル付)	¥148,200	¥164,200	
14670-454	トランスアピロト(スプレッド)	¥136,600	¥151,200	
14670-455	トランスアピロト(スプレッド)	¥233,700	¥259,100	
14670-468	トランスアピロト(スプレッド)	¥186,200	¥206,100	
14672-200	ピスツルビロト(ガラス製) 14.6cm	¥69,800	¥75,800	
14672-380	ピスツルビロト(ガラス製)22.9 1000	¥46,300	¥50,300	
14672-400	ピスツルビロト(綿栓付 ガラス製)	¥68,600	¥74,500	
14672-410	ピスツルビロト(綿栓付 ガラス製)	¥75,100	¥81,600	
14672-412	ピスツルビロト(綿栓付 ガラス製)	¥75,400	¥81,900	
14672-608	ピスツルビロト(ガラス製)	¥86,900	¥103,900	
14673-010	ピスツルビロト(オートタイプ ガラス製)2ml	¥48,000	¥52,200	
14673-043	ピスツルビロト(オートタイプ ガラス製)2ml	¥58,800	¥63,800	
14705-003	キビトリチュフ 70ul 200入	¥9,700	¥10,600	
14705-008	キビトリチュフ 70ul 200入	¥7,100	¥7,700	
14900-532	コンパクトリップ XS 100枚	¥10,300	¥10,800	
14905-000	温度記録計(温度幅2種 選択可能)	¥467,300	¥922,000	
14905-014	USBケーブル付付録付録付録	¥73,000	¥79,300	
14905-016	USBケーブル付付録付録付録	¥157,100	¥170,700	
14905-018	温湿度計用ラジウムカード	¥22,100	¥24,700	
14905-020	温湿度計用ラジウムカードキット	¥41,100	¥45,900	
15001-372	スリッパ	¥56,400	¥73,400	
15001-374	スリッパ - カート	¥23,200	¥25,300	
15001-376	スリッパ - カート	¥24,600	¥26,700	
15001-378	スリッパ - カート	¥27,600	¥29,900	
15001-380	スリッパ - カート	¥29,900	¥32,600	
15001-382	スリッパ - カート	¥35,700	¥38,800	
15001-384	スリッパ - カート	¥37,700	¥40,900	
15001-664	安全キビネット(2ドア)	¥1,180,600	¥2,147,000	
15001-668	安全キビネット(2ドア)	¥1,711,800	¥3,114,000	
15020T	10ul 2mlフィルター Eタイプ LB S96x10R	¥18,800	¥20,300	
15145-508	SWAB FOAMTIP #120 PK50	¥45,600	¥48,800	
15146-027	Tissue Freezing Media(70ul) 4本入	¥87,800	¥96,300	
15147-802	ミカトム替刃	¥133,900	¥146,900	
15148-031	Tissue Freezing Media	¥185,600	¥203,600	
15159-409	Paraplast Embedding Media 1kg x 8	¥103,500	¥112,900	
15159-464	Paraplast Plus 1kgx8	¥107,100	¥117,400	
15159-486	Paraplast X-tra 1kg x 8PK	¥123,100	¥135,000	
15160-099	パイロカプテイスル Truncated 8x8mm	¥222,600	¥181,700	
15160-157	パイロカプテイスル Truncated 12x12mm	¥252,800	¥286,500	
15160-215	パイロカプテイスル Square 22x22mm	¥220,500	¥255,500	
15160-270	パイロカプテイスル 288入	¥221,500	¥243,300	
15160-339	パイロカプテイスル	¥223,700	¥246,800	
15170-048	血球計算盤	¥391,900	¥434,000	
15170-079	血球計算盤(位相差顕微鏡用)	¥243,100	¥269,200	
15170-081	血球計算盤	¥235,900	¥261,200	
15170-089	血球計算盤	¥235,600	¥260,900	
15170-090	Nageotte計算盤	¥398,500	¥441,200	
15170-168	血球計算盤(Complete Chamber Set)	¥285,500	¥316,100	
15170-170	血球計算盤(サイズ 75x32mm)	¥208,100	¥230,400	
15170-172	血球計算盤 0.1mm	¥218,000	¥241,400	
15170-173	血球計算盤	¥180,200	¥199,500	
15170-208	血球計算盤	¥177,900	¥196,900	
15170-230	血球計算盤	¥254,200	¥281,500	
15170-263	血球計算盤(位相差顕微鏡用)	¥198,800	¥220,100	
15170-274	血球計算盤	¥167,600	¥185,600	
15401-402	マイクロピロト(ガラス製)	¥158,900	¥169,200	
15401-413	マイクロピロト(ガラス製)	¥179,700	¥191,400	
15401-560	毛細管(マイクロ入り/目盛付)	¥13,700	¥14,600	
15404-355	WintrobeTube Non-Plugged 115mm 200	¥76,900	¥83,900	
15406-025N	タイプ 154滅菌ろ過用ポリエチレン膜 100	¥12,800	¥14,800	
15406-047N	タイプ 154滅菌ろ過用ポリエチレン膜 100	¥17,800	¥20,800	
15406-142G	フィルター膜, PES, φ142mm, 0.45um, 25	¥34,200	¥37,600	
15407-47-MIN	PES製マイクロフィルター 47mm 100	¥17,800	¥20,800	
15425-024	Indicator Strips 250入	¥13,100	¥14,400	
15431-360	小動物用ケージ(エコノミケージ) 1	¥55,500	¥58,300	
15551-016	VWR Antistatic StaticAway Brush	¥24,400	¥27,400	
15552-194	SPILFYTER 汚染防止用スプレー 500L	¥57,300	¥61,300	
15555-702	ふるい振とう器	¥340,100	¥382,000	
15555-710	ふるい振とう器	¥335,100	¥376,400	
15555-714	ふるい振とう器	¥323,500	¥363,400	
15555-716	ふるい振とう器	¥306,800	¥346,100	
15555-718	ふるい振とう器	¥335,000	¥376,400	
15555-720	ふるい振とう器	¥301,400	¥340,000	
15555-722	ふるい振とう器	¥329,800	¥370,500	
15555-724	ふるい振とう器	¥303,700	¥342,700	
15555-726	ふるい振とう器	¥335,100	¥376,400	
15555-728	ふるい振とう器	¥305,200	¥344,400	
15555-730	ふるい振とう器	¥335,100	¥376,400	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
15555-732	ふるい振とう器	¥355,100	¥400,700	
15555-734	ふるい振とう器	¥359,000	¥403,300	
15555-740	ふるい振とう器	¥373,700	¥421,700	
15555-742	ふるい振とう器	¥403,700	¥453,500	
15555-746	ふるい振とう器	¥1,539,000	¥2,922,000	
15555-748	ふるい振とう器	¥311,300	¥351,200	
15555-750	ふるい振とう器	¥335,000	¥376,400	
15555-752	ふるい振とう器	¥303,800	¥342,700	
15555-754	ふるい振とう器	¥335,100	¥376,500	
15555-756	ふるい振とう器	¥295,400	¥333,300	
15555-758	ふるい振とう器	¥340,400	¥382,400	
15555-760	ふるい振とう器	¥303,700	¥342,700	
15555-764	ふるい振とう器	¥303,700	¥342,600	
15555-766	ふるい振とう器	¥328,500	¥369,100	
15555-768	ふるい振とう器	¥303,700	¥342,700	
15555-772	ふるい振とう器	¥326,300	¥368,100	
15555-774	ふるい振とう器	¥357,500	¥401,600	
15555-776	ふるい振とう器	¥327,500	¥369,600	
15555-782	ふるい振とう器	¥376,100	¥422,500	
15555-786	ふるい振とう器	¥403,100	¥452,800	
15555-790	ふるい振とう器	¥388,700	¥436,600	
15555-792	ふるい振とう器	¥670,700	¥1,318,000	
15566-024	XRF Thin-Film Sample Support	¥76,700	¥82,500	
15566-030	XRF Thin-Film Sample Support	¥105,000	¥113,000	
15566-036	XRF Thin-Film Sample Support	¥54,600	¥59,300	
15566-050	XRF Thin-Film Sample Support 500	¥67,000	¥72,800	
15566-052	XRF Thin-Film Sample Support	¥67,900	¥73,700	
15566-058	XRF Thin-Film Sample Support	¥140,000	¥152,100	
15670T	200ul リトフットチップ LB S96X10R	¥9,600	¥10,300	
15704-112	試料コップ	¥211,800	¥232,300	
15704-122	試料コップ	¥345,100	¥378,500	
15704-126	ホジナゲ	¥186,400	¥204,400	
15713-507	発泡スロ-ルコンテナ(ホ-リスレンフォーム・蓋付)	¥16,100	¥16,700	
15713-531	発泡スロ-ルコンテナ(ホ-リスレンフォーム・蓋付)	¥33,200	¥34,500	
15713-543	Thermosafe輸送コップ(カート付)3個入	¥97,700	¥101,500	
15713-551	発泡スロ-ルコンテナ(ホ-リスレンフォーム・蓋付) 3	¥87,200	¥90,600	
15713-607	発泡スロ-ルコンテナ(ホ-リスレンフォーム・蓋付)	¥80,600	¥88,400	
15713-608	発泡スロ-ルコンテナ(ホ-リスレンフォーム・蓋付)	¥120,700	¥132,700	
15713-610	発泡スロ-ルコンテナ(ホ-リスレンフォーム・蓋付)	¥85,700	¥89,300	
15714-465	Thermosafe輸送コップ(カート付)8個入	¥135,200	¥140,500	
15714-500	発泡スロ-ルコンテナ(ホ-リスレンフォーム・蓋付)	¥32,700	¥34,000	
15714-522	発泡スロ-ルコンテナ(ホ-リスレンフォーム・蓋付)	¥147,200	¥166,000	
15714-704	発泡スロ-ルコンテナ ThermoSafe 2個入	¥63,300	¥65,800	
15715-103	U-TEK 保冷剤	¥32,400	¥34,800	
15715-105	U-TEK 保冷剤	¥35,500	¥38,100	
15715-127	U-TEK 保冷剤 12袋入	¥22,600	¥24,200	
15715-130	U-TEK 保冷剤	¥24,000	¥25,700	
15715-132	U-tek保冷剤(-23℃) 6個入	¥76,700	¥82,200	
15715-137	U-tek保冷剤(-23℃) 36個入	¥42,400	¥45,500	
15715-152	U-tek保冷剤(-23℃) 16個入	¥24,500	¥26,300	
15715-156	U-tek保冷剤(-23℃) 8個入	¥20,500	¥21,900	
15720T	200ul リトフットチップ LB 1000/bag	¥4,000	¥4,300	
16001-170	トランスファーピ-ット(ホ-ット)	¥94,600	¥99,800	
16001-172	トランスファーピ-ット(ホ-ット) 500	¥12,200	¥14,200	
16001-174	トランスファーピ-ット(ホ-ット) 500	¥11,900	¥13,800	
16001-176	トランスファーピ-ット(ホ-ット)	¥87,400	¥101,500	
16001-178	トランスファーピ-ット(ホ-ット) 500	¥14,600	¥16,900	
16001-180	トランスファーピ-ット(ホ-ット) 500	¥13,000	¥13,700	
16001-182	トランスファーピ-ット(ホ-ット) 500	¥14,700	¥15,500	
16001-186	トランスファーピ-ット(ホ-ット) 500	¥16,100	¥18,300	
16001-188	トランスファーピ-ット(ホ-ット)7.0mL 500	¥13,100	¥15,200	
16001-190	トランスファーピ-ット(ホ-ット) 400	¥9,200	¥10,300	
16001-192	トランスファーピ-ット(ホ-ット)	¥103,000	¥108,100	
16001-194	トランスファーピ-ット(ホ-ット) 500	¥14,700	¥16,600	
16001-196	トランスファーピ-ット(ホ-ット) 100本入	¥11,000	¥11,700	
16001-588	MPC Series Plastic Couplings	¥550,400	¥595,000	
16002-560	HistoGel Cooling Block 1個	¥156,100	¥168,800	
16002-920	15L Complete Reactor System J	¥15,908,600	¥28,937,000	
16003-352	温度モニター DIG Type T	¥695,900	¥1,266,000	
16003-446	Extension Coad Stan-mini Yell	¥40,200	¥43,400	
16003-472	KWIK MOUNT KLAMP UNIVERSAL 1個	¥47,900	¥52,100	
16004-096	VWR VistaVision加-ゲ-ス No.1 100	¥7,500	¥8,200	
16004-300	VWR VistaVision加-ゲ-ス No.1 18mm	¥13,100	¥14,400	
16004-304	VWR VistaVision加-ゲ-ス No.2	¥55,800	¥61,000	
16004-314	VWR VistaVision加-ゲ-ス No.1 100	¥7,200	¥8,000	
16004-322	VWR VistaVision加-ゲ-ス No.1.5	¥8,200	¥9,000	
16004-344	VWR VistaVision加-ゲ-ス No.1.5	¥7,900	¥8,600	
16004-406	VWR VistaVision HistoBondスライト 72	¥54,900	¥60,200	
16004-422	VWR VistaVisionスライト 72	¥19,400	¥21,300	
16004-430	VWR VistaVisionスライト 72	¥17,400	¥20,700	
16007-000	VWR 血球計算盤(加-ゲ-ス)20x26mm 10	¥26,700	¥29,300	
16007-050	ホ-リ-トリック広口コップ 20ml 1本	¥35,800	¥38,900	
16053-962	広口コップ	¥116,200	¥126,200	
16054-224	Sample Vials 40mL 144本	¥409,000	¥452,900	
16058-757	ホ-リ-トリックコップ(ホ-リスレン)	¥56,000	¥58,800	
16058-780	ホ-リ-トリックコップ(ホ-リスレン)	¥74,100	¥77,700	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
16058-815	ネジ口サンワルボトル(ホリソ製)	¥87,800	¥92,100	
16058-848	ネジ口サンワルボトル(ホリソ製)	¥61,200	¥64,200	
16058-870	ネジ口サンワルボトル(ホリソ製)	¥46,400	¥50,000	
16059-590	250ml広口ボトル(ホリソ製) 250個	¥211,900	¥226,600	
16059-592	500ml広口ボトル(ホリソ製) 125個	¥141,100	¥151,000	
16059-594	1000ml広口ボトル(ホリソ製) 50個	¥101,000	¥108,100	
16087-042	ガラスバール	¥117,400	¥130,000	
16087-044	ガラスバール	¥144,700	¥160,200	
16087-048	ガラスバール	¥142,700	¥158,000	
16087-052	ガラスバール	¥148,400	¥164,400	
16087-056	ガラスバール	¥184,800	¥204,700	
16087-058	ガラスバール	¥179,100	¥198,300	
16087-060	ガラスバール	¥221,600	¥245,400	
16087-062	ガラスバール	¥88,800	¥98,300	
16087-064	ガラスバール	¥93,000	¥103,000	
16087-066	ガラスバール	¥80,400	¥89,000	
16087-068	ガラスバール	¥95,000	¥105,200	
16087-070	ガラスバール	¥107,700	¥119,300	
16087-072	ガラスバール	¥92,700	¥102,600	
16087-256	ガラスバール用キャップ Unlined 100	¥21,000	¥22,500	
16087-260	ガラスバール用キャップ Unlined 100	¥23,200	¥25,000	
16120-753	角型ボトル PP 1000ml 6本	¥67,700	¥76,400	
16121-183	ホリソ製角型ボトル 250mL 6個	¥70,000	¥75,600	
16125-720	広口ホリソ	¥84,900	¥89,100	
16125-752	広口ホリソ	¥84,100	¥88,200	
16125-785	広口ホリソ	¥100,400	¥105,400	
16125-810	広口ホリソ	¥135,700	¥156,700	
16125-843	広口ホリソ	¥81,600	¥94,200	
16125-876	広口ホリソ	¥58,500	¥61,400	
16125-901	広口ホリソ	¥55,500	¥58,300	
16125-903	ホリソ(キャップ付)広口丸瓶	¥59,600	¥62,600	
16125-967	広口ホリソ	¥34,900	¥36,600	
16125-970	液体リテナー 箱付き	¥155,300	¥163,800	
16125-984	ホリソ・角型(広口・高密度ホリソ製)	¥406,400	¥437,300	
16125-985	ホリソ・角型(広口・高密度ホリソ製)	¥32,100	¥34,500	
16125-986	ホリソ・角型(広口・高密度ホリソ製)	¥211,200	¥227,200	
16125-987	ホリソ・角型(広口・高密度ホリソ製)	¥41,400	¥44,600	
16125-988	ホリソ・角型(広口・高密度ホリソ製)	¥246,700	¥265,400	
16125-989	ホリソ・角型(広口・高密度ホリソ製)	¥39,600	¥42,600	
16125-990	ホリソ・角型(広口・高密度ホリソ製)	¥179,700	¥193,300	
16125-991	ホリソ・角型(広口・高密度ホリソ製)	¥55,000	¥59,100	
16125-992	ホリソ・角型(広口・高密度ホリソ製)	¥245,300	¥263,900	
16125-993	ホリソ・角型(広口・高密度ホリソ製)	¥32,600	¥35,100	
16125-996	ホリソ・角型(広口・高密度ホリソ製)	¥516,800	¥553,200	
16125-997	ホリソ・角型(広口・高密度ホリソ製)	¥32,800	¥35,300	
16126-176	広口ホリソ	¥35,100	¥38,300	
16126-198	広口ホリソ	¥42,300	¥46,200	
16126-212	広口ホリソ	¥50,700	¥55,300	
16126-234	広口ホリソ 250ml	¥56,400	¥60,400	
16126-278	広口ホリソ	¥40,600	¥44,400	
16126-314	広口ホリソ	¥69,600	¥74,800	
16126-358	広口ホリソ	¥46,100	¥51,800	
16126-446	リテナー 5ml 6	¥24,000	¥25,700	
16126-490	バアルラック 1	¥23,800	¥26,100	
16126-492	バアルラック	¥162,900	¥177,500	
16126-494	バアルラック	¥133,500	¥147,400	
16126-668	スライド用ホックス(杓付/ピン)	¥24,800	¥26,500	
16129-061	Nalgeneボトル 1L 24個入	¥159,500	¥180,200	
16149-151	I-LOC CLSRS Black 100個入	¥53,100	¥57,900	
16150-821	広口ガラスボトル	¥132,100	¥143,400	
16150-824	広口ガラスボトル	¥90,100	¥97,900	
16150-827	広口ガラスボトル	¥131,500	¥142,900	
16150-921	広口ガラスボトル	¥134,400	¥146,000	
16150-924	広口ガラスボトル	¥98,900	¥107,400	
16150-927	広口ガラスボトル	¥132,600	¥144,100	
16150-954	広口ガラスボトル	¥138,200	¥150,100	
16151-027	ガラス瓶(細口・ホリソ型・キャップ付)	¥87,100	¥93,200	
16151-049	ガラス瓶(細口・ホリソ型・キャップ付)	¥52,000	¥56,400	
16151-060	ガラス瓶(細口・ホリソ型・キャップ付)	¥54,600	¥59,300	
16151-082	ガラス瓶(細口・ホリソ型・キャップ付)	¥74,700	¥81,200	
16151-107	ガラス瓶(細口・ホリソ型・キャップ付)	¥51,700	¥56,200	
16151-129	ガラス瓶(細口・ホリソ型・キャップ付)	¥72,300	¥78,600	
16151-220	ガラス瓶(細口・ホリソ型・キャップ付)	¥108,200	¥117,500	
16151-242	ガラス瓶(細口・ホリソ型・キャップ付)	¥59,500	¥64,600	
16151-286	ガラス瓶(細口・ホリソ型・キャップ付)	¥60,000	¥65,200	
16151-300	ガラス瓶(細口・ホリソ型・キャップ付)	¥74,100	¥80,500	
16151-322	ガラス瓶(細口・ホリソ型・キャップ付)	¥48,900	¥58,400	
16151-344	ガラス瓶(細口・ホリソ型・キャップ付)	¥68,300	¥74,200	
16151-650	ガラス瓶(細口・ホリソ型・キャップ付)	¥97,100	¥104,000	
16151-672	ガラス瓶(細口・ホリソ型・キャップ付)	¥50,700	¥55,100	
16151-719	ガラス瓶(細口・ホリソ型・キャップ付)	¥51,000	¥55,400	
16151-730	ガラス瓶(細口・ホリソ型・キャップ付)	¥61,200	¥66,500	
16151-752	ガラス瓶(細口・ホリソ型・キャップ付)	¥43,000	¥46,700	
16151-774	ガラス瓶(細口・ホリソ型・キャップ付)	¥61,200	¥66,500	
16151-865	ガラス瓶(細口・ホリソ型・キャップ付)	¥99,500	¥106,500	
16151-901	ガラス瓶(細口・ホリソ型・キャップ付)	¥56,000	¥60,900	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
16151-923	ガラス瓶(細口・ホスト型・キャップ付)	¥69,000	¥75,000	
16151-945	ガラス瓶(細口・ホスト型・キャップ付)	¥54,900	¥58,700	
16151-967	ガラス瓶(細口・ホスト型・キャップ付)	¥80,200	¥87,100	
16152-020	ガラス瓶 褐色(細口・ホスト型・キャップ付)	¥51,800	¥56,200	
16152-041	ガラス瓶 褐色(細口・ホスト型・キャップ付)	¥59,300	¥64,400	
16152-063	ガラス瓶 褐色(細口・ホスト型・キャップ付)	¥64,500	¥70,100	
16152-085	ガラス瓶 褐色(細口・ホスト型・キャップ付)	¥54,200	¥58,900	
16152-100	ガラス瓶 褐色(細口・ホスト型・キャップ付)	¥75,500	¥82,100	
16152-405	ガラス瓶 褐色(細口・ホスト型・キャップ付)	¥51,300	¥55,700	
16152-427	ガラス瓶 褐色(細口・ホスト型・キャップ付)	¥57,100	¥62,000	
16152-449	ガラス瓶 褐色(細口・ホスト型・キャップ付)	¥62,800	¥65,900	
16152-460	ガラス瓶 褐色(細口・ホスト型・キャップ付)	¥47,000	¥51,000	
16152-482	ガラス瓶 褐色(細口・ホスト型・キャップ付)	¥65,800	¥71,500	
16152-642	ガラス瓶 褐色(細口・ホスト型・キャップ付)	¥51,100	¥55,500	
16152-664	ガラス瓶 褐色(細口・ホスト型・キャップ付)	¥59,600	¥64,800	
16152-686	ガラス瓶 褐色(細口・ホスト型・キャップ付)	¥67,000	¥70,400	
16152-700	ガラス瓶 褐色(細口・ホスト型・キャップ付)	¥61,500	¥66,800	
16152-722	ガラス瓶 褐色(細口・ホスト型・キャップ付)	¥78,000	¥84,700	
16153-102	ガラス瓶 褐色(細口・ホスト型・キャップ付)	¥59,900	¥65,100	
16153-135	ガラス瓶 褐色(細口・ホスト型・キャップ付)	¥62,000	¥66,400	
16153-168	ガラス瓶 褐色(細口・ホスト型・キャップ付)	¥76,000	¥79,700	
16153-190	ガラス瓶 褐色(細口・ホスト型・キャップ付)	¥53,900	¥64,400	
16153-226	ガラス瓶 褐色(細口・ホスト型・キャップ付)	¥67,900	¥73,700	
16159-128	角型テフロンボトル 60ml 96本入	¥185,800	¥210,800	
16161-337	キャップ(テフロン付 耐薬品性)	¥201,600	¥219,000	
16161-360	キャップ(テフロン付 耐薬品性)	¥197,300	¥214,300	
16161-428	キャップ(テフロン付 耐薬品性)	¥299,200	¥314,100	
16161-519	キャップ(テフロン付 耐薬品性)	¥31,100	¥33,800	
16161-541	キャップ(テフロン付 耐薬品性)	¥82,000	¥89,100	
16161-933	Kimble BTL Serum Glass Clear 20ml	¥227,900	¥248,500	
16161-935	BOTTLE SERUM CLR 50ML 288本入	¥361,200	¥393,900	
16161-937	BOTTLE SERUM CLR 100ML 144本入	¥268,900	¥293,200	
16161-938	BOTTLE SERUM CLEAR 125ML 144本	¥297,200	¥324,100	
16171-284	ガラスバケツ	¥57,400	¥60,200	
16171-285	ガラスバケツ	¥54,200	¥56,900	
16171-288	ガラスバケツ	¥102,600	¥111,900	
16171-293	ガラスバケツ	¥76,300	¥80,100	
16171-651	ストッパ(ストレートラック) 1000個入	¥949,300	¥1,035,200	
16171-885	ビーズ(栓) 外径13mm	¥422,300	¥460,500	
16179-988	広口茶瓶	¥35,400	¥38,400	
16180-003	広口茶瓶	¥58,800	¥63,900	
16180-025	広口茶瓶	¥70,100	¥76,200	
16180-047	広口茶瓶	¥92,100	¥100,000	
16180-058	広口茶瓶 250mL 84	¥179,500	¥195,000	
16180-207	広口茶瓶 30mL(キャップ28-400)24本入	¥99,000	¥107,600	
16180-229	広口茶瓶 60mL(キャップ33-400)24本入	¥86,500	¥94,000	
16180-240	広口茶瓶	¥98,000	¥106,500	
16180-262	広口茶瓶	¥131,500	¥142,900	
16180-284	広口茶瓶	¥72,300	¥78,500	
16180-309	広口茶瓶	¥102,200	¥111,000	
16180-681	広口茶瓶	¥36,700	¥39,900	
16180-706	広口茶瓶	¥48,400	¥52,600	
16180-728	広口茶瓶	¥47,500	¥51,600	
16180-852	広口茶瓶	¥36,700	¥39,800	
16180-874	広口茶瓶	¥49,500	¥53,800	
16180-903	広口ガラスボトル	¥215,000	¥225,600	
16180-904	広口ガラスボトル	¥106,500	¥115,700	
16180-906	広口ガラスボトル	¥122,700	¥133,300	
16188-027	広口角ガラスビン	¥97,800	¥106,200	
16188-049	広口角ガラスビン	¥114,700	¥122,800	
16188-082	広口角ガラスビン	¥71,700	¥77,900	
16188-107	広口角ガラスビン	¥96,000	¥102,800	
16188-129	広口角ガラスビン	¥121,000	¥130,100	
16188-286	広口角ビン	¥361,500	¥392,700	
16188-300	広口角ビン	¥221,700	¥237,300	
16188-504	広口角ガラスビン	¥103,600	¥112,600	
16188-526	広口角ガラスビン	¥121,300	¥131,700	
16188-548	広口角ガラスビン	¥107,900	¥117,200	
16188-560	広口角ガラスビン	¥86,300	¥93,800	
16188-581	広口角ガラスビン	¥99,700	¥108,300	
16188-708	広口角ガラスビン	¥99,200	¥106,200	
16188-741	広口角ガラスビン	¥43,700	¥47,400	
16188-854	広口角ガラスビン	¥65,100	¥69,700	
16188-876	広口角ガラスビン 480 24入	¥125,400	¥134,900	
16188-899	広口角ガラスビン	¥112,300	¥117,900	
16189-008	広口角ガラスビン	¥119,700	¥130,000	
16189-030	広口角ガラスビン	¥156,200	¥169,600	
16189-096	広口角ガラスビン	¥94,300	¥102,400	
16189-121	広口角ガラスビン	¥169,200	¥183,800	
16189-154	広口角ガラスビン	¥174,500	¥186,900	
16189-188	広口角ガラスビン	¥149,800	¥157,200	
16193-810	広口ガラス瓶	¥78,800	¥85,600	
16193-843	広口ガラス瓶	¥59,100	¥64,200	
16193-876	広口ガラス瓶	¥48,800	¥53,000	
16193-901	広口ガラス瓶	¥47,700	¥51,800	
16193-919	広口ガラスボトル	¥120,400	¥130,800	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
16194-008	広口ガラス瓶	¥91,600	¥99,500	
16194-020	広口ガラス瓶	¥94,800	¥103,000	
16194-041	広口ガラス瓶	¥58,100	¥63,700	
16194-052	広口ガラス瓶	¥263,700	¥289,200	
16194-100	保存ビン(ガラス製)	¥92,900	¥101,000	
16194-256	広口ガラス瓶	¥120,700	¥131,100	
16194-278	広口ガラス瓶	¥142,100	¥171,500	
16194-290	広口ガラス瓶	¥97,000	¥106,400	
16194-314	広口ガラス瓶 240ml(58-400) 24本入	¥119,400	¥129,700	
16194-609	広口ガラス瓶	¥74,200	¥80,600	
16194-620	広口ガラス瓶	¥82,200	¥89,300	
16194-722	広口ガラス瓶	¥66,900	¥72,700	
16194-744	広口ガラス瓶	¥83,800	¥91,000	
16194-857	広口ガラス瓶	¥99,600	¥108,200	
16194-879	広口ガラス瓶	¥108,800	¥118,100	
16194-909	広口ガラスボトル	¥186,800	¥203,000	
16195-022	保存ビン(ガラス製)	¥117,900	¥140,900	
16195-088	保存ビン(ガラス製)	¥98,500	¥107,000	
16195-102	保存ビン(ガラス製)	¥69,500	¥75,500	
16195-306	保存ビン(ガラス製)	¥170,200	¥184,900	
16195-328	保存ビン(ガラス製)	¥100,600	¥109,300	
16195-340	保存ビン(ガラス製)	¥106,800	¥116,000	
16195-361	保存ビン(ガラス製)	¥149,000	¥161,800	
16195-383	保存ビン(ガラス製)	¥128,900	¥140,000	
16195-408	保存ビン(ガラス製)	¥137,200	¥149,100	
16195-667	保存ビン(ガラス製)	¥73,800	¥80,200	
16195-689	保存ビン(ガラス製)	¥65,900	¥71,600	
16195-703	保存ビン(ガラス製)	¥52,100	¥56,600	
16195-725	保存ビン(ガラス製)	¥101,000	¥109,700	
16195-747	保存ビン(ガラス製)	¥48,700	¥52,900	
16195-769	保存ビン(ガラス製)	¥62,000	¥67,400	
16195-882	保存ビン(ガラス製)	¥316,800	¥344,200	
16195-884	保存ビン(ガラス製)	¥486,800	¥528,800	
16196-025	広口ボトル(ホリリテレン)	¥31,700	¥35,800	
16196-047	広口ボトル(ホリリテレン)	¥46,200	¥49,700	
16196-069	広口ボトル(ホリリテレン)	¥28,000	¥30,100	
16196-080	広口ボトル(ホリリテレン)	¥19,900	¥21,400	
16196-207	広口ガラス瓶	¥98,400	¥107,900	
16196-229	広口ガラス瓶	¥129,600	¥140,800	
16196-240	広口ガラス瓶 960mL	¥82,400	¥89,600	
16196-262	広口ガラス瓶	¥58,300	¥63,300	
16196-284	広口ガラス瓶	¥55,700	¥60,500	
16198-910	ハイルキャップ	¥64,800	¥70,700	
16198-911	ハイルキャップ	¥72,700	¥79,300	
16198-912	ハイルキャップ	¥91,500	¥99,800	
16198-913	ハイルキャップ	¥125,500	¥136,800	
16198-915	遠心チューブ (ガラス製加工ホリリテレン)	¥144,000	¥157,000	
16198-918	ハイルキャップ	¥131,300	¥143,200	
16198-938	ハイルキャップ	¥251,200	¥273,900	
16198-960	ハイルキャップ	¥116,200	¥126,700	
16198-993	ハイルキャップ	¥107,500	¥117,200	
16199-001	キャップ(ヒニールラテ付)	¥273,600	¥297,200	
16199-023	キャップ(ヒニールラテ付)	¥191,300	¥207,800	
16199-045	キャップ(ヒニールラテ付)	¥194,400	¥211,200	
16199-056	キャップ(ヒニールラテ付)	¥234,700	¥254,900	
16199-089	キャップ(ヒニールラテ付)	¥232,000	¥252,000	
16199-103	キャップ(ヒニールラテ付)	¥250,900	¥272,600	
16199-125	キャップ(ヒニールラテ付)	¥267,200	¥290,200	
16199-136	キャップ(ヒニールラテ付)	¥302,800	¥329,000	
16199-169	広口角ビン	¥11,000	¥12,000	
16199-180	キャップ(ヒニールラテ付)	¥426,000	¥462,800	
16199-205	キャップ(ヒニールラテ付)	¥484,700	¥526,600	
16199-249	キャップ(ヒニールラテ付)	¥514,800	¥559,200	
16199-260	キャップ(ヒニールラテ付)	¥295,100	¥320,500	
16199-282	キャップ(ヒニールラテ付)	¥330,100	¥358,600	
16199-307	キャップ(ヒニールラテ付)	¥372,900	¥405,100	
16199-329	キャップ(ヒニールラテ付)	¥223,500	¥242,800	
16200-166	キャップ(ホリリテレンラテ付)	¥10,400	¥11,300	
16200-188	キャップ(ホリリテレンラテ付)	¥14,100	¥15,300	
16200-202	キャップ(ホリリテレンラテ付)	¥15,400	¥18,400	
16200-224	キャップ(ホリリテレンラテ付)	¥16,700	¥18,100	
16200-280	キャップ(ホリリテレンラテ付)	¥24,300	¥26,400	
16200-304	キャップ(ホリリテレンラテ付)	¥28,200	¥30,600	
16200-326	キャップ(ホリリテレンラテ付)	¥51,200	¥55,600	
16200-444	PTFE ティッシュ (43-400) 500	¥748,600	¥785,700	
16200-448	PTFE ティッシュ (48-400) 500	¥894,800	¥939,200	
16214-998	キャップ(シリコンラテ付)	¥7,200	¥7,800	
16215-017	キャップ(シリコンラテ付)	¥8,400	¥8,800	
16215-040	キャップ(シリコンラテ付)	¥10,100	¥11,000	
16215-050	キャップ(シリコンラテ付)	¥10,300	¥11,200	
16215-061	キャップ(シリコンラテ付)	¥10,600	¥11,500	
16215-083	キャップ(シリコンラテ付)	¥10,700	¥11,700	
16215-094	キャップ(シリコンラテ付)	¥16,000	¥19,100	
16215-119	キャップ(シリコンラテ付)	¥20,300	¥22,000	
16215-130	キャップ(シリコンラテ付)	¥22,900	¥24,800	
16216-982	キャップ(ホリリテレンラテ付)	¥7,500	¥8,100	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
16217-002	キャップ (ホ°リシ-ルライ-付)	¥7,800	¥8,500	
16217-023	キャップ (ホ°リシ-ルライ-付)	¥351,900	¥382,300	
16217-024	キャップ (ホ°リシ-ルライ-付)	¥8,000	¥8,700	
16217-046	キャップ (ホ°リシ-ルライ-付)	¥8,700	¥9,500	
16217-068	キャップ (ホ°リシ-ルライ-付)	¥9,200	¥10,000	
16217-090	キャップ (ホ°リシ-ルライ-付)	¥9,100	¥9,900	
16217-095	キャップ (ホ°リシ-ルライ-付)	¥352,900	¥383,400	
16217-104	キャップ (ホ°リシ-ルライ-付)	¥9,200	¥10,000	
16217-126	キャップ (ホ°リシ-ルライ-付)	¥12,200	¥13,300	
16217-148	キャップ (ホ°リシ-ルライ-付)	¥20,000	¥21,700	
16218-124	ガラスバ-イアル	¥107,100	¥118,600	
16218-126	ガラスバ-イアル	¥88,400	¥97,900	
16218-128	ガラスバ-イアル	¥98,600	¥109,100	
16218-130	ガラスバ-イアル	¥126,500	¥140,100	
16218-132	ガラスバ-イアル	¥124,900	¥138,300	
16218-134	ガラスバ-イアル	¥162,000	¥179,300	
16218-136	ガラスバ-イアル	¥147,700	¥163,600	
16218-186	Glass Jugs Without Cap 3785ml 4本	¥31,100	¥33,800	
16218-238	Glass Jugs without Cap 1893ml 6本	¥40,600	¥44,100	
16299-280	お°口瓶入りキャップ°/PTFE/45mm/10個入	¥76,900	¥83,800	
16300-062	LDPE Dropper Tips 8mm ピ°ツク	¥19,200	¥20,900	
16300-072	LDPE Dropper Tips 13mm ピ°ツク	¥25,300	¥27,600	
16300-093	LDPE Dropper Tips 20mm ピ°ツク	¥26,700	¥29,200	
16300-134	PP Dropper Tip Caps(20-410)100入	¥28,800	¥31,400	
16300-342	LDPE Dropping Bottels(20-410)30ml	¥132,000	¥144,700	
16344-101	ス°レ°ホ°トル Qorpak 500mL 1	¥10,800	¥11,600	
16354-024	ス°ト°ホ°トル/ホ°ト°ラッ° (ホ°リシ-ル製)15mL	¥24,100	¥25,800	
16354-046	ス°ト°ホ°トル/ホ°ト°ラッ° (ホ°リシ-ル製)30mL	¥25,400	¥27,300	
16354-068	ス°ト°ホ°トル/ホ°ト°ラッ° (ホ°リシ-ル製)60mL	¥24,000	¥25,700	
16358-236	スト°ル-	¥117,400	¥128,000	
16362-102	ス°ト°ホ°トル(褐色) 30ml 12本入	¥144,300	¥156,700	
16401-4707ACK	バ-イザ°ト100 モ°タ- 0.2um	¥29,000	¥30,300	
16455-000	細口ガラスピ°	¥232,200	¥254,600	
16455-241	コ°テナ- (塗料缶輸送用)	¥129,000	¥140,600	
16455-262	発泡ス°ト°コ°テナ- (PS°ホ°ム°蓋付)	¥242,300	¥254,300	
16465-001	コ°テナ- (危険性物質輸送用)	¥173,900	¥214,300	
16465-003	コ°テナ- (危険性物質輸送用)	¥155,500	¥208,200	
16465-005	コ°テナ- (危険性物質輸送用)	¥160,300	¥205,900	
16465-007	コ°テナ- (危険性物質輸送用)	¥169,200	¥218,500	
16465-009	コ°テナ- (危険性物質輸送用)	¥255,100	¥267,800	
16465-011	コ°テナ- (危険性物質輸送用)	¥268,000	¥281,200	
16465-013	コ°テナ- (危険性物質輸送用)	¥283,500	¥297,600	
16465-015	コ°テナ- (危険性物質輸送用)	¥275,900	¥291,200	
16465-043	コ°テナ- (危険性物質輸送用)250mL	¥189,000	¥223,700	
16465-047	コ°テナ- (危険性物質輸送用)	¥259,300	¥292,400	
16465-049	コ°テナ- (危険性物質輸送用)	¥168,100	¥193,300	
16465-051	コ°テナ- (危険性物質輸送用)	¥382,500	¥401,500	
16465-160	メ°ラ-シ°テム (郵送用梱包器材セット)	¥140,700	¥150,800	
16466-030	凍結保存チューブ° (SuperClear) 500本	¥29,800	¥31,300	
16466-032	凍結保存チューブ° (SuperClear) 500本	¥29,700	¥31,200	
16466-034	凍結保存チューブ° (SuperClear) 500本	¥30,000	¥31,500	
16466-036	凍結保存チューブ° (SuperClear) 500本	¥29,700	¥31,100	
16466-038	凍結保存チューブ° (SuperClear) 500本	¥30,000	¥31,400	
16466-040	凍結保存チューブ° (SuperClear) 500本	¥29,600	¥31,000	
16466-042	凍結保存チューブ° (SuperClear)2.0mL ST	¥104,000	¥109,200	
16466-044	凍結保存チューブ° (SuperClear)2.0ml	¥81,400	¥85,500	
16466-046	凍結保存チューブ° (SuperClear)1.5mL ST	¥103,100	¥108,200	
16466-048	凍結保存チューブ° (SuperClear)1.5ml	¥80,300	¥84,300	
16466-050	凍結保存チューブ° (SuperClear)0.5mL ST	¥103,000	¥108,100	
16466-052	凍結保存チューブ° (SuperClear)0.5ml	¥79,800	¥83,800	
16466-054	凍結保存チューブ° (SuperClear)0.5mL ST	¥103,900	¥109,100	
16466-056	凍結保存チューブ° (SuperClear)0.5ml	¥81,300	¥85,400	
16466-058	Micro Centrifuge Tube 500	¥103,500	¥108,600	
16466-060	凍結保存チューブ° (SuperClear)2.0ml	¥80,300	¥84,200	
16466-062	凍結保存チューブ° (SuperClear)1.5mL ST	¥104,000	¥109,200	
16466-064	凍結保存チューブ° (SuperClear)1.5ml	¥81,100	¥85,100	
16466-068	キャップ (ス°ル°-ク°リア凍結保存チューブ°用)	¥54,300	¥57,000	
16466-070	キャップ (ス°ル°-ク°リア凍結保存チューブ°用)	¥54,800	¥57,500	
16466-072	キャップ (ス°ル°-ク°リア凍結保存チューブ°用)	¥54,800	¥57,500	
16466-074	キャップ (ス°ル°-ク°リア凍結保存チューブ°用)	¥54,800	¥57,500	
16466-076	キャップ (ス°ル°-ク°リア凍結保存チューブ°用)	¥54,500	¥57,200	
16466-078	キャップ (ス°ル°-ク°リア凍結保存チューブ°用)	¥54,300	¥57,000	
16466-080	キャップ (ス°ル°-ク°リア凍結保存チューブ°用)	¥54,500	¥57,200	
16466-082	キャップ (ス°ル°-ク°リア凍結保存チューブ°用)	¥54,000	¥56,700	
16466-084	キャップ (ス°ル°-ク°リア凍結保存チューブ°用)	¥54,200	¥56,900	
16508B	PC製フィルタ°ホ°タ°- φ47/50mm	¥28,400	¥32,800	
16514E	PC製フィルタ°ホ°タ°- φ13mm 12	¥10,200	¥12,000	
16517E	PC製フィルタ°ホ°タ°- φ25mm	¥16,400	¥17,800	
16532GUK	ミ°ニ°サ°ル°ト°ハ°イ°ロ°- 0.2u 28mm LL型 γ線ST	¥20,000	¥21,800	
16532K	ミ°ニ°サ°ル°ト°ハ°イ°ロ°- 0.2u 28mm ル°ア°/ロ°ツク EOG	¥18,000	¥19,800	
16532Q	ミ°ニ°サ°ル°ト°ハ°イ°ロ°- 未滅菌 0.2u 28mm LL型	¥150,000	¥160,000	2024/3/14 追加
16533GUK	ミ°ニ°サ°ル°ト°ハ°イ°ロ°- 0.45u 28mm LL型 γ線ST	¥20,000	¥21,800	
16533K	ミ°ニ°サ°ル°ト°ハ°イ°ロ°- 0.45u28mm ル°ア°/ス°リッ° EOG	¥18,000	¥19,800	
16533Q	ミ°ニ°サ°ル°ト°ハ°イ°ロ°- 未滅菌 0.45u 28mm LS	¥150,000	¥160,000	2024/3/14 追加
16537K	ミ°ニ°サ°ル°ト°ハ°イ°ロ°- 0.45u 28mm ル°ア°/ロ°ツク S	¥18,000	¥19,800	
16537Q	ミ°ニ°サ°ル°ト°ハ°イ°ロ°- 未滅菌 0.45u 28mm LL	¥150,000	¥160,000	2024/3/14 追加

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
16541K	ミニバルブ0-0.2u 28mm ルア-スリップ S	¥18,000	¥19,800	
16541Q	ミニバルブ0-未滅菌 0.2u 28mm LS型	¥150,000	¥160,000	2024/3/14 追加
16550-101	Transpette Dispアルチトランスファーピペット	¥43,000	¥46,900	
16550-103	Transpette Dispアルチトランスファーピペット	¥181,300	¥209,300	
16550-104	96ウェルティープウェルプレート 2.5ml 24入	¥160,600	¥177,000	
16550-106	96ウェルティープウェルプレート 2.5ml ST	¥181,400	¥199,900	
16550-108	96ウェルティープウェルプレート用カバー(未滅菌)	¥43,300	¥46,400	
16550-110	Side Loader Collars(未滅菌)	¥62,300	¥68,700	
16553K	ミニバルブ0-0.1u 28mm ルア-ロック S	¥18,000	¥19,800	
16592GUK	ミニバルブ 0.8u 28mm LL型 γ線ST 50	¥18,200	¥20,000	
16592K	ミニバルブ 0.8u 28mm ルア-ロック EOG滅菌	¥16,200	¥18,000	
16592Q	ミニバルブ 0.8u 28mm ルア-ロック 500入	¥142,000	¥152,000	
16596HNC	ミニバルブ 0.2um 50個	¥21,000	¥23,100	
16628-029	比重ビン	¥229,400	¥250,200	
16628-041	比重ビン	¥207,900	¥226,700	
16628-062	比重ビン	¥215,100	¥234,600	
16628-067	比重ビン	¥124,200	¥135,400	
16649-900	洗浄瓶(広口)	¥18,100	¥21,300	
16649-901	洗浄瓶(広口)	¥23,500	¥27,700	
16649-905	洗浄瓶(広口)	¥18,100	¥19,400	
16649-910	洗浄瓶(広口)	¥18,100	¥19,400	
16649-911	洗浄瓶(広口)	¥18,100	¥21,300	
16649-915	洗浄瓶(広口)	¥21,800	¥23,300	
16649-917	洗浄瓶(広口)	¥18,100	¥19,400	
16649-918	洗浄瓶(広口)	¥24,000	¥25,700	
16649-925	洗浄瓶(広口)	¥18,000	¥19,300	
16649-935	洗浄瓶(広口)	¥18,000	¥19,400	
16649-940	洗浄瓶(広口)	¥17,900	¥19,100	
16649-943	洗浄瓶(広口)	¥17,800	¥19,100	
16649-945	洗浄瓶(広口)	¥18,100	¥19,400	
16649-950	洗浄瓶(広口)	¥18,100	¥19,400	
16649-955	洗浄瓶(広口)	¥18,100	¥19,400	
16649-966	洗浄瓶(広口)	¥18,100	¥19,400	
16649-980	薬品トレー	¥25,200	¥27,000	
16650-049	スプレーボトル 1000mL	¥34,200	¥36,700	
16650-082	洗浄瓶(広口)	¥101,100	¥108,700	
16650-093	洗浄瓶(広口)	¥101,100	¥108,700	
16650-107	洗浄瓶(広口)	¥101,100	¥108,100	
16650-162	洗浄瓶(広口)	¥80,500	¥86,200	
16650-173	洗浄瓶(広口)	¥80,500	¥86,600	
16650-184	洗浄瓶(広口)	¥80,500	¥86,200	
16650-195	洗浄瓶(広口)	¥78,900	¥84,500	
16650-242	洗浄瓶(広口)	¥62,100	¥66,500	
16650-253	洗浄瓶(広口)	¥70,400	¥75,300	
16650-264	洗浄瓶(広口)	¥62,100	¥73,100	
16650-275	洗浄瓶(広口) 6	¥23,400	¥25,100	
16650-322	洗浄瓶(広口)	¥42,700	¥45,700	
16650-333	洗浄瓶(広口)	¥42,700	¥45,700	
16650-344	洗浄瓶(広口)	¥42,700	¥45,700	
16650-355	洗浄瓶(広口)	¥42,700	¥45,700	
16650-537	洗浄瓶(広口) 500mL 6本入	¥21,700	¥23,200	
16650-550	広口角ボトル	¥162,100	¥178,700	
16650-552	広口角ボトル	¥77,200	¥85,100	
16650-554	広口角ボトル	¥151,600	¥167,100	
16650-556	広口角ボトル	¥129,900	¥143,100	
16650-558	広口角ボトル	¥112,700	¥124,200	
16650-560	広口角ボトル	¥88,200	¥97,200	
16650-562	広口角ボトル	¥200,700	¥221,200	
166MP-4	マイクロピペット イソエット 吸引ポンプ 1	¥205,000	¥209,100	
16746-01-COM	IPソケット MD8ソケット-111エット 1	¥845,000	¥921,100	
16746-01-EXC	接続ケーブル(30m) 1	¥93,000	¥101,400	
16746-01-SHB	ソケットアダプター (パイネット用) 1	¥900,000	¥981,000	
16746-01-SHT	ソケットアダプター (サニタリーフランジ用) 1	¥900,000	¥981,000	
16746SHBCOM	パイネットアダプター 1	¥1,750,000	¥1,907,500	2024/3/14 追加
16746SHTCOM	サニタリーフランジアダプター 1	¥2,080,000	¥2,267,200	2024/3/14 追加
17003-118	VALVE HIVAC 900 0-4MM 1個	¥65,800	¥71,400	
17003-120	VALVE HIVAC 900 0-8MM 1個	¥77,500	¥84,200	
17003-122	VALVE HIVAC 900 0-12MM 1個	¥90,200	¥98,000	
17003-144	Valve Hivac 1-ARM 0-8MM	¥83,800	¥91,000	
17003-150	0-8 Chem-Cap Plug with knob	¥49,900	¥54,300	
17006-592	Pie-Block System w/holder	¥482,900	¥878,000	
17018-022	フラスコ(テストチューブ用)	¥12,500	¥14,200	
17018-044	フラスコ(テストチューブ用)	¥13,500	¥15,400	
17027-023	フラスコ(テストチューブ用 先丸型)	¥7,600	¥8,600	
17027-045	フラスコ(テストチューブ用 先丸型)	¥8,100	¥9,200	
17027-067	フラスコ(テストチューブ用 先丸型)	¥8,900	¥10,200	
17027-089	フラスコ(テストチューブ用 先丸型)	¥10,000	¥11,400	
17028-026	フラスコ(テストチューブ用 先丸型)	¥12,700	¥14,400	
17028-048	フラスコ(テストチューブ用 先丸型)	¥13,400	¥15,300	
17028-060	フラスコ(テストチューブ用 先丸型)	¥19,300	¥22,000	
17030-006	フラスコ(チューブ用)	¥8,100	¥9,300	
17041-002	フラスコ(シリンダー-ボトル用)	¥10,800	¥12,100	
17050-003	フラスコ(乳脂計用)	¥11,600	¥13,200	
17065-000	フラスコ(ボトル用)	¥9,600	¥11,000	
17066-058	フラスコ(ボトル用)	¥12,900	¥14,600	
17066-080	フラスコ(ボトル用)	¥25,000	¥31,400	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
17066-089	ブラシ(ボトル用)	¥29,300	¥33,500	
17068-009	ブラシ(乳脂計用)	¥13,300	¥16,700	
17110-025	Buret and Pipet Brushes 5本	¥6,400	¥7,300	
17130-022	ブラシ(テストチューブ用)	¥8,000	¥9,200	
17130-044	ブラシ(テストチューブ用)	¥13,800	¥15,700	
17130-066	ブラシ(テストチューブ用)	¥13,500	¥15,400	
17140-004	ブラシ(ボトル用 先丸型)	¥16,600	¥7,300	
17145-020	ブラシ(ボトル洗浄用)	¥28,000	¥18,400	
17145-041	ブラシ(ボトル洗浄用)	¥25,800	¥22,400	
17145-063	ブラシ(ボトル洗浄用)	¥30,200	¥25,300	
17150-020	ブラシ(フラスコ・ボトル用)	¥13,800	¥15,600	
17150-041	ブラシ(フラスコ・ボトル用)	¥15,300	¥17,400	
17155-001	ブラシ(フラスコ・ボトル用)	¥22,800	¥11,200	
17160-001	ブラシ(ファン用)5本入	¥23,400	¥26,600	
17161-004	ブラシ(円錐形)	¥103,300	¥117,800	
17170-005	ブラシ(ビーカー用)	¥61,400	¥70,100	
17190-002	ブラシ(ビーカー用)	¥9,000	¥12,100	
17190-007	ブラシ(ビーカー用)	¥44,900	¥51,200	
17210-008	ブラシ(平型)	¥71,100	¥81,000	
17216-026	ブラシ(らくだ毛)	¥8,900	¥10,100	
17216-048	ブラシ(らくだ毛)	¥11,400	¥13,000	
17227-022	ワイヤクリーナーブラシ 直径0.3cm/全長15m	¥29,300	¥36,800	
17227-044	ワイヤクリーナーブラシ 直径0.5cm/全長15m	¥33,700	¥38,400	
17230-027	ブラシ(やぎ毛)	¥74,700	¥84,500	
17230-049	ブラシ(やぎ毛)	¥86,700	¥98,700	
17230-060	ブラシ(やぎ毛)	¥123,500	¥140,800	
17235-020	Paint/Dust Brushes	¥26,000	¥29,700	
17235-042	Paint/Dust Brushes ブラシ幅5.1cm	¥3,400	¥3,800	
17250-002	ブラシ(テーブル用)	¥9,200	¥10,500	
17263-004	ワイヤブラシ	¥75,400	¥85,900	
17302-117	バケツ(耐腐食性)	¥410,100	¥439,000	
17370T-X	200ul フィルターカート 3チップ LB S96X10R	¥20,400	¥22,000	
17370-X	200ul フィルターカート 3チップ NT S96X10R	¥19,700	¥21,200	
17528-80ACD	ディスク・サンプラー・チップ・フィルター タイプ 175	¥26,000	¥29,000	
17528-80BZD	ディスク・サンプラー・チップ・フィルター 10リ	¥27,000	¥30,100	
17558K	ミニサンプラーSRP 0.2u 15mm スピンドル 50	¥12,600	¥13,800	
17558Q	ミニサンプラーSRP 0.2u 15mm スピンドル 500	¥102,000	¥114,000	
17559K	ミニサンプラーSRP 0.45u 15mm スピンドル 50	¥12,600	¥13,800	
17559Q	ミニサンプラーSRP 0.45u 15mm スピンドル 500	¥102,000	¥114,000	
17573ACK	ミニサンプラーSRP15 0.2u 25mm 50	¥18,000	¥19,000	
17573K	ミニサンプラーSRP 0.2u 15mm ルアー/スリッパ 50	¥12,600	¥13,800	
17573Q	ミニサンプラーSRP 0.2u 15mm ルアー/スリッパ 500	¥102,000	¥114,000	
17574K	ミニサンプラーSRP 0.45u15mm ルアー/スリッパ 50	¥12,600	¥13,800	
17574Q	ミニサンプラーSRP 0.45u 15mm ルアー/スリッパ 500	¥102,000	¥114,000	
17575K	ミニサンプラーSRP 0.2u 25mm ルアー/スリッパ 50	¥26,600	¥29,260	
17575Q	ミニサンプラーSRP 0.2u 25mm ルアー/スリッパ 500	¥210,000	¥231,000	
17576K	ミニサンプラーSRP 0.45u25mm ルアー/スリッパ 50	¥26,600	¥28,800	
17576Q	ミニサンプラーSRP 0.45u 25mm ルアー/スリッパ 500	¥210,000	¥230,000	
1757HK	ミニサンプラー PES 0.2u 25mm LS型 50	¥25,400	¥26,800	
17593K	ミニサンプラー 1.2u 28mm ルアー/ロック S 50入	¥16,200	¥18,000	
17593Q	ミニサンプラー 1.2u 28mm ルアー/ロック 500入	¥142,000	¥152,000	
17594Q	ミニサンプラー 5.0u 28mm ルアー/ロック 500入	¥142,000	¥152,000	
17683-060	Support Stands SCIENCEWARE	¥73,200	¥78,300	
17683-258	VWR Talon Double Buret Clamp	¥36,200	¥39,600	
17761ACK	ミニサンプラーRC15-AC 0.2u15mm ルアー/ロック S	¥14,000	¥15,000	
17761K	ミニサンプラーRC 0.2u 15mm ルアー/スリッパ 50	¥10,200	¥11,600	
17761Q	ミニサンプラーRC 0.2u 15mm ルアー/スリッパ 500	¥86,000	¥98,800	
17762K	ミニサンプラーRC 0.45u 15mm ルアー/スリッパ 50	¥10,200	¥11,600	
17762Q	ミニサンプラーRC 0.45u 15mm ルアー/スリッパ 500	¥86,000	¥98,800	
17764ACK	ミニサンプラーRC25-AC 0.2u25mm ルアー/ロック S	¥18,400	¥19,800	
17764K	ミニサンプラーRC 0.2u 25mm ルアー/スリッパ 50	¥12,400	¥14,000	
17764Q	ミニサンプラーRC 0.2u 25mm ルアー/スリッパ 500	¥105,800	¥118,000	
17765ACK	ミニサンプラーRC25-AC 0.45u25mm ルアー/ロック S	¥18,400	¥19,800	
17765K	ミニサンプラーRC 0.45u 25mm ルアー/スリッパ 50	¥12,400	¥14,000	
17765Q	ミニサンプラーRC 0.45u 25mm ルアー/スリッパ 500	¥105,800	¥118,000	
1776BK	ミニサンプラーNY 0.2u 15mm ルアー/スリッパ 50	¥12,400	¥13,600	
1776CK	ミニサンプラーNY 0.45u 15mm ルアー/スリッパ 50	¥12,400	¥13,600	
1776DACK	ミニサンプラーPES 0.2u 15mm ルアー/スリッパ 50	¥17,000	¥18,600	
17788-102	白衣 サイズ 56 1	¥31,500	¥35,800	
17788-106	Reusable Lab Coat 89cm 34サイズ	¥25,600	¥25,800	
17820K	ミニサンプラーSRP 0.45u 4mm ルアー/スリッパ 50	¥8,000	¥9,000	
17820Q	ミニサンプラーSRP 0.45u 4mm ルアー/スリッパ 500	¥70,000	¥78,000	
17821K	ミニサンプラーRC 0.2u 4mm ルアー/スリッパ 50	¥7,000	¥8,000	
17821Q	ミニサンプラーRC 0.2u 4mm ルアー/スリッパ 500	¥58,800	¥68,000	
17822K	ミニサンプラーRC 0.45u 4mm ルアー/スリッパ 50	¥7,000	¥8,000	
17822Q	ミニサンプラーRC 0.45u 4mm ルアー/スリッパ 500	¥58,800	¥68,000	
17823K	ミニサンプラーラズ 0.2u 28mm ルアー/ロック S	¥32,200	¥35,420	
17823Q	ミニサンプラーラズ 0.2u 28mm ルアー/ロック 500	¥290,600	¥319,660	
17824K	ミニサンプラーGF 28mm ルアー/ロック 50入	¥19,000	¥20,900	
17824Q	ミニサンプラーGF 28mm ルアー/ロック 500入	¥182,000	¥200,200	
17825Q	ミニサンプラー 1.2u 28mm ルアー/ロック 500入	¥290,600	¥319,660	
17829K	ミニサンプラーラズ 0.45u 28mm ルアー/ロック S	¥32,200	¥35,420	
17829Q	ミニサンプラーラズ 0.45u 28mm ルアー/ロック 500	¥290,600	¥319,660	
17845ACK	ミニサンプラーNY 0.2u/25mm ルアー/スリッパ S 50	¥17,600	¥18,600	
17845Q	ミニサンプラーNY 0.2u/25mm ルアー/スリッパ 500	¥158,000	¥168,000	
17846ACK	ミニサンプラーNY 0.45u/25mm ルアー/スリッパ S 50	¥17,600	¥18,600	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
17846Q	ミニサルトNY 0.45u/25mm ルア-/スリッ° 500	¥158,000	¥168,000	
1784C-K	ミニサルトNY 0.45u/25mm ルア-/スリッ° 50	¥22,600	¥24,800	
17966-622	実験室用 ラック	¥68,900	¥74,500	
17966-624	実験室用 ラック	¥70,900	¥76,700	
17966-626	実験室用 ラック	¥83,500	¥90,300	
180-22E	サルトホスト 150mL 12	¥4,800	¥5,600	
180-23E	サルトホスト 250mL 12	¥6,600	¥7,800	
180-24E	サルトホスト 500mL 12	¥11,600	¥13,000	
180-25E	サルトホスト 1000mL 12	¥14,000	¥16,000	
18075D	サルトホ P20 ホー/エ-テルスルホ ルア-/ロツク S	¥18,600	¥19,200	
18089D	サルトホ P20 0.2um HB/HB 10	¥19,600	¥20,800	
180D03E	サルトホ RF/レシ-バ°-付 0.1um 250mL 12	¥15,800	¥18,600	
180D05E	サルトホ RF/レシ-バ°-付 0.1um 1L 12	¥32,000	¥37,400	
180E01-2	サルトホ RF50 0.22um 50mL 24	¥23,200	¥26,800	
180E02E	サルトホ RF/レシ-バ°-付0.22um 150ml 12	¥14,000	¥16,200	
180E03E	サルトホ RF/レシ-バ°-付0.22um 250ml 12	¥15,800	¥18,600	
180E04E	サルトホ RF/レシ-バ°-付0.22um 500ml 12	¥25,600	¥29,800	
180E05E	サルトホ RF/レシ-バ°-付 0.22um 1L 12	¥32,000	¥37,400	
180E12E	サルトホ BT/レシ-バ°-なし0.22u 150ml 12	¥12,800	¥14,800	
180E13E	サルトホ BT/レシ-バ°-なし0.22um 250ml 12	¥13,600	¥15,800	
180E14E	サルトホ BT/レシ-バ°-なし0.22um500ml 12	¥16,600	¥19,200	
180E15E	サルトホ BT/レシ-バ°-なし 0.22um 1L 12	¥25,600	¥29,800	
180F01-2	サルトホ RF50 0.45um 50mL 24	¥23,200	¥26,800	
180F02E	サルトホ RF/レシ-バ°-付0.45um 150ml 12	¥14,000	¥16,200	
180F03E	サルトホ RF/レシ-バ°-付0.45um 250ml 12	¥15,800	¥18,600	
180F04E	サルトホ RF/レシ-バ°-付0.45um 500ml 12	¥25,600	¥29,800	
180F05E	サルトホ RF/レシ-バ°-付0.45um 1000ml 12	¥32,000	¥37,400	
180F15E	サルトホ BT/レシ-バ°-なし0.45um 1000mL12	¥25,600	¥29,800	
18343-051	滅菌トチャンバ°	¥234,700	¥439,000	
18385-220	Laboratory Gas Lighter 1入	¥22,400	¥24,200	
18385-223	Replacement Flint Tips 40個入	¥3,600	¥4,100	
18406-047N	RC φ47mm 0.45um 100pc	¥22,000	¥25,400	
18407-047N	RC φ47mm 0.2um 100pc	¥22,000	¥25,400	
18700-115	収納トレ-	¥325,800	¥352,300	
18700-148	収納トレ-	¥304,700	¥329,400	
18700-150	収納トレ-	¥326,200	¥352,600	
18889-002	PROSAT ヲバ°-(ホ°リ°ル°レ°) 30	¥9,900	¥10,500	
18900-012	メソ°レン°フィルタ° 3um 未滅菌 10枚入	¥60,300	¥71,700	
18900-014	メソ°レン°フィルタ° 0.2um 未滅菌 10枚入	¥70,000	¥75,600	
18900-282	メソ°レン°フィルタ° 3um 未滅菌 1枚	¥10,800	¥11,700	
18900-292	Wall Mounting Unit	¥37,000	¥40,000	
18900-296	シリコ°ア°タ°フ°タ°+PTFE°フ°イ°ツカ°ル°バ°	¥21,200	¥22,900	
18901-114	BRAND DischargeTube 1個	¥94,400	¥108,900	
18901-220	BRAND seripettor Dispensers	¥129,200	¥139,700	
19001-003	VWR シャ°フ° スリ°チ°タ°ス° M レ°ト° 3.8L	¥79,200	¥83,100	
19420-164	BOMB GENERAL PURPOSE 125ML 1	¥1,641,200	¥1,836,000	
19420-166	PTFE Cup with Cover 23mL	¥159,900	¥175,200	
19472-088	Cup Cover PTFE 23ml	¥160,100	¥172,700	
19757-082	カ°ス°ホ°ン°ハ° 用カ°ト	¥236,000	¥259,700	
19766-014	ユ°テ°リ°テ°イ°カ°ト	¥259,300	¥280,300	
19766-016	ユ°テ°リ°テ°イ°カ°ト	¥340,800	¥368,400	
19766-018	多目的カ°ト	¥390,600	¥422,300	
19766-022	多目的カ°ト	¥424,600	¥459,000	
19766-026	多目的カ°ト	¥412,300	¥445,700	
19766-032	多目的カ°ト	¥391,500	¥423,200	
19766-040	多目的カ°ト	¥547,000	¥591,400	
19766-068	実験室用 ラック	¥26,000	¥28,100	
19766-070	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥46,700	¥50,500	
19766-072	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥143,100	¥154,700	
19766-074	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥148,400	¥160,400	
19766-076	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥49,100	¥53,100	
19766-080	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥48,500	¥52,400	
19766-082	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥76,300	¥82,500	
19766-084	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥159,100	¥172,000	
19766-086	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥101,700	¥109,900	
19766-088	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥200,200	¥216,400	
19766-090	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥100,800	¥109,000	
19766-092	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥258,100	¥279,000	
19766-094	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥79,400	¥85,900	
19766-098	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥180,600	¥195,200	
19766-106	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥49,100	¥53,100	
19766-114	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥54,000	¥58,300	
19766-118	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥62,500	¥67,600	
19766-122	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥254,100	¥274,700	
19766-134	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥71,800	¥77,600	
19766-136	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥72,700	¥78,600	
19766-140	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥274,100	¥296,300	
19766-142	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥135,700	¥146,700	
19766-158	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥61,000	¥66,000	
19766-162	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥53,600	¥58,000	
19766-166	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥71,400	¥77,100	
19766-168	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥211,200	¥228,300	
19766-176	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥65,000	¥70,300	
19766-182	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥71,400	¥77,200	
19766-192	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥83,400	¥90,200	
19766-196	実験器具棚 ス°バ°-レ°ク°タ°リ°ス°	¥143,400	¥155,000	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
19766-200	実験器具棚 ス-IP-レクタリ-ス	¥100,000	¥108,100	
19766-204	ネットカート 61x152.4	¥362,400	¥391,800	
19766-208	実験器具棚 ス-IP-レクタリ-ス	¥109,500	¥118,300	
19766-222	実験室用 ラック	¥18,000	¥19,400	
19766-226	実験室用 ラック	¥64,000	¥69,100	
19766-236	キャスター(ストップ付)	¥30,600	¥33,100	
19766-238	キャスター(ストップ付)	¥25,000	¥27,000	
19766-240	キャスター(ストップ付)	¥33,200	¥35,900	
19766-242	MetroMAX Stem Casters InterMetro	¥26,100	¥28,300	
19766-244	キャスター(ストップ付)	¥55,000	¥59,500	
19766-246	キャスター(ストップ付)	¥49,700	¥53,700	
19766-248	MetroMAX Stem Casters InterMetro	¥43,000	¥46,500	
19766-250	キャスター(ストップ付)	¥37,000	¥40,000	
19766-252	MetroMAX Stem Casters InterMetro	¥47,800	¥51,600	
19766-254	MetroMAX Stem Casters InterMetro	¥25,000	¥27,000	
19766-256	実験室用 ラック	¥20,000	¥21,600	
19766-260	実験室用 ラック	¥48,900	¥52,800	
19766-264	実験室用 ラック	¥20,500	¥22,200	
19766-274	実験室用 ラック	¥34,400	¥37,200	
19766-282	実験室用 ラック	¥28,000	¥30,300	
19766-328	実験室用 ラック	¥21,500	¥23,200	
19766-334	実験器具棚 ス-IP-レクタリ-ス	¥66,300	¥71,700	
19766-364	実験器具棚 ス-IP-レクタリ-ス	¥80,900	¥87,500	
19766-386	実験器具棚 ス-IP-レクタリ-ス	¥372,300	¥402,400	
19766-400	実験器具棚 ス-IP-レクタリ-ス	¥81,000	¥87,600	
19766-414	実験器具棚 ス-IP-レクタリ-ス	¥136,300	¥147,400	
19766-416	実験器具棚 ス-IP-レクタリ-ス	¥97,900	¥105,900	
19766-428	実験器具棚 ス-IP-レクタリ-ス	¥118,200	¥127,800	
19766-474	キャスター(ストップ付)	¥22,700	¥24,500	
19766-476	MetroMAX Stem Casters InterMetro	¥40,500	¥43,800	
19766-478	MetroMAX Stem Casters InterMetro	¥34,800	¥37,700	
19766-480	ネットカート	¥125,400	¥135,500	
19766-482	ネットカート	¥118,400	¥128,000	
19766-486	MetroMAX Stem Casters InterMetro	¥45,100	¥48,700	
19766-492	キャスター(ストップ付)	¥31,600	¥34,200	
19766-494	MetroMAX Stem Casters InterMetro	¥31,600	¥34,200	
19766-504	実験器具棚 ス-IP-レクタリ-ス	¥129,700	¥140,200	
19766-506	実験器具棚 ス-IP-レクタリ-ス	¥160,300	¥173,300	
19766-508	実験器具棚 ス-IP-レクタリ-ス	¥182,300	¥197,100	
19766-510	実験器具棚 ス-IP-レクタリ-ス	¥269,700	¥291,600	
19766-512	実験器具棚 ス-IP-レクタリ-ス	¥188,700	¥204,000	
19766-516	実験器具棚 ス-IP-レクタリ-ス	¥324,600	¥350,900	
19766-518	実験器具棚 ス-IP-レクタリ-ス	¥210,300	¥227,300	
19766-520	実験器具棚 ス-IP-レクタリ-ス	¥184,100	¥199,000	
19766-522	実験器具棚 ス-IP-レクタリ-ス	¥210,600	¥227,700	
19766-524	実験器具棚 ス-IP-レクタリ-ス	¥211,200	¥228,300	
19766-526	実験器具棚 ス-IP-レクタリ-ス	¥280,800	¥303,500	
19766-528	実験器具棚 ス-IP-レクタリ-ス	¥283,000	¥305,900	
19766-530	実験器具棚 ス-IP-レクタリ-ス	¥332,000	¥358,900	
19766-532	実験器具棚 ス-IP-レクタリ-ス	¥192,600	¥208,200	
19766-534	実験器具棚 ス-IP-レクタリ-ス	¥195,800	¥211,700	
19766-536	ネットカート 61x91.4	¥302,900	¥327,500	
19766-538	実験器具棚 ス-IP-レクタリ-ス	¥271,300	¥293,300	
19766-540	ネットカート 61x121.9	¥422,100	¥456,300	
19766-542	実験器具棚 ス-IP-レクタリ-ス	¥288,600	¥312,000	
19766-544	実験器具棚 ス-IP-レクタリ-ス	¥527,100	¥959,000	
19766-620	実験器具棚 ス-IP-レクタリ-ス	¥84,100	¥90,900	
19766-622	多目的カート	¥389,000	¥420,500	
19766-624	多目的カート	¥425,200	¥459,700	
19766-626	多目的カート	¥459,200	¥496,400	
19766-628	多目的カート	¥422,100	¥456,300	
19766-630	多目的カート	¥457,200	¥494,200	
19766-632	多目的カート	¥503,500	¥544,300	
19766-634	多目的カート	¥507,100	¥548,200	
19766-636	多目的カート	¥474,400	¥512,900	
19766-638	多目的カート	¥598,800	¥647,300	
19766-940	実験室用 ラック	¥189,200	¥204,500	
19766-942	実験室用 ラック	¥187,800	¥203,100	
19766-944	実験室用 ラック	¥187,800	¥203,100	
19766-946	実験室用 ラック	¥222,800	¥240,900	
19766-948	実験室用 ラック	¥222,800	¥240,900	
19766-950	実験室用 ラック	¥281,400	¥304,200	
19766-952	実験室用 ラック	¥282,500	¥305,400	
19766-954	実験室用 ラック	¥336,800	¥364,200	
19766-956	実験室用 ラック	¥199,800	¥216,000	
19766-958	実験室用 ラック	¥203,900	¥220,500	
19766-960	実験室用 ラック	¥233,200	¥252,100	
19766-962	実験室用 ラック	¥278,600	¥301,200	
19766-964	実験室用 ラック	¥275,900	¥298,300	
19766-966	実験室用 ラック	¥297,000	¥321,100	
19766-968	実験室用 ラック	¥352,800	¥381,400	
19766-970	実験室用 ラック	¥425,800	¥460,300	
19806-001	Adhesive Permabond 910 3g	¥6,300	¥6,800	
19806-007	Adhesive Permabond 910FS 1oz	¥24,800	¥26,500	
19872-004	SPORTROL Spore Strips 100入	¥241,900	¥265,400	
19872-038	インジケータラベル Inprinted 5000枚/ロール	¥123,900	¥135,900	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
1ZAD0025	サリ-フワン`変換アダプター	¥845,000	¥921,100	
1ZAD0026	47mm GMF専用アダプター	¥845,000	¥921,100	
1ZAS-0007	シリコ製チューブ (1m/肉厚2mm) 1	¥2,200	¥2,420	
1ZPX-D0002	培地カバー	¥5,900	¥6,490	
20001-396	軽作業グローブ (革製) M 1双	¥10,800	¥11,500	
20068-542	ニードル(液体混合用)固定バーなし	¥211,000	¥221,400	
20068-546	ニードル(液体混合用)固定バー付	¥316,000	¥331,600	
20068-548	ニードル(液体混合用)固定バー付	¥309,200	¥324,500	
20068-550	ニードル(液体混合用)固定バー付	¥333,700	¥350,200	
20068-552	ニードル(液体混合用)固定バー付	¥333,700	¥350,200	
20068-554	ニードル(液体混合用)固定バー付	¥329,300	¥345,600	
20068-556	ニードル(液体混合用)固定バーなし	¥252,000	¥264,500	
20068-558	ニードル(液体混合用)固定バーなし	¥252,000	¥264,500	
20068-606	給餌針/挿管針	¥523,500	¥549,400	
20068-608	給餌針/挿管針	¥501,500	¥526,400	
20068-610	給餌針/挿管針	¥541,400	¥568,200	
20068-612	給餌針/挿管針	¥543,700	¥570,600	
20068-614	給餌針/挿管針	¥545,800	¥572,900	
20068-616	給餌針/挿管針	¥520,300	¥546,100	
20068-620	給餌針/挿管針	¥522,600	¥548,500	
20068-622	給餌針/挿管針	¥519,700	¥545,400	
20068-624	給餌針/挿管針	¥522,600	¥548,500	
20068-628	給餌針/挿管針	¥527,200	¥553,300	
20068-630	給餌針/挿管針	¥522,200	¥548,000	
20068-636	給餌針/挿管針	¥654,200	¥686,700	
20068-638	給餌針/挿管針	¥582,700	¥611,500	
20068-640	給餌針/挿管針	¥609,700	¥639,900	
20068-642	給餌針/挿管針	¥612,700	¥643,000	
20068-644	給餌針/挿管針(小動物用/curved)	¥580,600	¥609,400	
20068-646	給餌針/挿管針	¥496,800	¥521,400	
20068-648	給餌針/挿管針	¥496,800	¥521,400	
20068-650	給餌針/挿管針	¥576,500	¥605,100	
20068-652	給餌針/挿管針	¥533,700	¥560,100	
20068-666	給餌針(滅菌済/デイスル-ザブル)	¥735,100	¥1,298,000	
20068-668	給餌針(滅菌済/デイスル-ザブル)	¥315,900	¥331,600	
20068-670	給餌針(滅菌済/デイスル-ザブル)	¥311,400	¥326,800	
20068-684	ニードル 16G 5.1cm 12本入	¥200,500	¥224,300	
20068-686	ニードル 16G 10.2cm 12本入	¥200,300	¥224,100	
20068-688	ニードル 14G 10.2cm 12本入	¥313,800	¥351,100	
20068-690	Laboratory Pipetting Needle 14G	¥224,100	¥250,700	
20069-038	ニードル	¥525,000	¥587,400	
201256TAN	溶解バフター, ボトル, 20ml	¥8,000	¥8,300	
20170-002	8連/12連PCRチューブ&キャップ	¥79,400	¥80,900	
20170-010	PCRシカケルチューブ&ト-ムキャップ	¥34,100	¥34,700	
20170-012	PCRシカケルチューブ&ワットキャップ 1000本	¥34,600	¥35,200	
20170-022	マイクロ遠心チューブ 1.5ml 500本	¥13,400	¥15,500	
20170-024	マイクロ遠心チューブ/耐遠心/1.5ml/500x10	¥100,200	¥105,000	
20170-026	マイクロ遠心チューブ/耐遠心/1.5ml/500x10	¥100,200	¥105,000	
20170-028	マイクロ遠心チューブ/耐遠心/1.5ml/500x10	¥99,700	¥104,500	
20170-030	マイクロ遠心チューブ/耐遠心/1.5ml/500x10	¥100,200	¥105,000	
20170-032	マイクロ遠心チューブ/耐遠心/1.5ml/500x10	¥100,200	¥105,000	
20170-034	マイクロ遠心チューブ/耐遠心/1.5ml/500x10	¥99,700	¥104,500	
20170-036	マイクロ遠心チューブ 1.5ml 500本入	¥13,000	¥13,600	
20170-038	マイクロ遠心チューブ 1.5ml 500本	¥13,000	¥13,600	
20170-084	マイクロ遠心チューブ 2.0ml 500本	¥15,700	¥16,500	
20170-086	マイクロ遠心チューブ 2.0ml 500本	¥15,700	¥18,100	
20170-088	マイクロ遠心チューブ 2.0ml 500本	¥15,700	¥16,500	
20170-090	マイクロ遠心チューブ 2.0ml 500本	¥15,700	¥16,500	
20170-092	マイクロ遠心チューブ 2.0ml 500本	¥15,700	¥16,500	
20170-094	マイクロ遠心チューブ 2.0ml 500本	¥15,700	¥16,500	
20170-096	マイクロ遠心チューブ 2.0ml 500本	¥15,700	¥16,400	
20170-098	マイクロ遠心チューブ 2.0ml 500本	¥17,600	¥18,500	
20170-170	マイクロ遠心チューブ 2.0mL NT 500本	¥16,500	¥17,300	
20170-688	マイクロチューブ(凍結保存チューブ用)	¥100,000	¥111,200	
20901-540	TUBE MICROCNTRFG NAT 2ML CS500	¥47,700	¥51,000	
20901-541	TUBE CENT SAFLOC CS500 AMB 2ML 1	¥59,300	¥63,900	
20901-905	Microcentrifuge Floating Rack 4	¥28,200	¥30,200	
20907-938	Rack Grip No-Wire 60穴 1個	¥19,200	¥20,900	
21001-000	PEEK Tubing チューブ 1	¥41,600	¥45,000	
21008-000	遠心チューブ 15ml 100本x10箱	¥82,000	¥86,100	
21008-022	遠心チューブ 15ml キャップ 500個x2袋	¥29,900	¥31,400	
21008-103	遠心チューブ 15ml(ス-ル-クア PP製)	¥92,300	¥96,800	
21008-177	ス-ル-クア遠心チューブ 50ml 500入	¥70,300	¥73,800	
21008-178	遠心チューブ 50ml(ス-ル-クア PP製)ST	¥114,400	¥120,000	
21008-189	遠心チューブ 15ml(ホ-リスル製) 1000本	¥141,900	¥149,000	
21008-197	遠心チューブ 15ml(ホ-リスル製)	¥95,200	¥99,900	
21008-202	遠心チューブ 15ml(ホ-リスル製) 50x10パック	¥91,900	¥96,400	
21008-212	遠心チューブ 15ml(ホ-リスル製) 50x10袋	¥89,700	¥103,500	
21008-214	ス-ル-クア遠心チューブ 10ml 50本x10	¥72,200	¥75,800	
21008-216	ス-ル-クア遠心チューブ 15ml S 50本x10	¥71,100	¥74,600	
21008-240	ス-ル-クア遠心チューブ 50ml 50本x10袋	¥110,200	¥115,600	
21008-242	ス-ル-クア遠心チューブ 50ml ST 50本x10袋	¥112,500	¥118,000	
21008-485	チューブ用ラック 50mLチューブ用140x76x89mm	¥405,400	¥477,400	
21008-486	チューブ用ラック 50mLチューブ用140x140x89mm	¥387,200	¥413,800	
21008-487	チューブ用ラック 50mLチューブ用210x210x89mm	¥561,600	¥600,000	
21008-918	Falcon遠心管 PP ST 15mL	¥193,800	¥186,800	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
21008-959	TUBE MICROCNTRFG NAT 1.5ML 500	¥43,700	¥46,900	
21009-342	Nalgene Oak Ridge遠心チューブ 50mL 10	¥93,600	¥109,800	
21010-402	Nalgen Square Media Bottles 500mL	¥119,800	¥137,300	
21020-607	遠心チューブ (ガラス製加工済)	¥83,900	¥91,500	
21020-640	遠心チューブ (ガラス製加工済) 125	¥85,800	¥93,600	
21020-684	遠心チューブ (ガラス製加工済)	¥90,900	¥99,100	
21020-695	遠心チューブ (ガラス製加工済)	¥78,900	¥86,000	
21020-720	遠心チューブ (ガラス製加工済) 5ml	¥67,300	¥73,300	
21020-742	遠心チューブ (ガラス製加工済) 10mL	¥73,200	¥79,900	
21020-764	遠心チューブ (ガラス製加工済) 15ml	¥78,500	¥85,600	
21021-422	丸底遠心チューブ (Plain Top) 30ml	¥127,900	¥139,500	
21021-424	Plain-Top Tubes 15mL 6本入	¥76,300	¥83,200	
21060-000	Graduated Hopkins Vaccine CentTube	¥423,500	¥461,800	
21150-404	Multi-Cap Tool Black	¥36,700	¥40,800	
21432-739	加圧チューブ (TLC展開用チューブ) 7分付	¥268,900	¥292,100	
21432-761	加圧チューブ (滑り止め7分付チューブ)	¥261,900	¥284,500	
21432-794	加圧チューブ (ガラス製の)	¥56,800	¥61,700	
21432-819	加圧チューブ (滑り止めの)	¥114,600	¥124,500	
21475-524	定規 (蛍光タイプ) 25cm	¥15,500	¥17,100	
21475-529	UVトランスイルミネーター	¥43,800	¥81,000	
21503-123	KIMAX Chromatographic Columns 1	¥123,800	¥133,100	
21511-004	加圧チューブ (HPLC用フィルター)	¥23,100	¥25,000	
21511-006	加圧チューブ (HPLC用フィルター)	¥28,600	¥30,900	
21511-012	Tubing Cutter Plastic	¥27,800	¥30,100	
21511-380	PEEK Tubing 1分付	¥29,500	¥31,900	
21511-448	Inline Solvent Filter 0.5µm	¥91,600	¥99,000	
21511-494	PEEK Tubing Sleeve 1	¥4,400	¥4,700	
21511-544	チューブカッター用替刃 (A327)	¥6,900	¥7,400	
21511-890	Sleeve Micrtight	¥2,400	¥2,600	
21571-000	Bunsen Round Jaw Extension Clamps	¥24,500	¥25,800	
21571-258	Bunsen Round Jaw Extension Clamps	¥34,300	¥36,100	
21571-305	Bunsen Round Jaw Extension Clamps	¥42,500	¥44,700	
21572-501	VWR Taron Regular Clamp Holder	¥11,400	¥12,500	
21572-658	Talon All-Position Clamp Holder	¥37,700	¥41,100	
21572-805	VWR Talon II Extension Clamps 1	¥24,600	¥27,000	
21573-265	VWR Talon Single Buret Clamp	¥25,700	¥28,100	
21573-708	VWR Talon ThreeProng Swivel Clamp1	¥26,300	¥28,700	
21716-102	ホースコック	¥12,500	¥13,500	
21733-024	ステンレス製クリップ サイズ6 10個入	¥64,200	¥69,700	
21733-032	ステンレス製クリップ サイズ10 10個入	¥64,500	¥70,100	
21733-040	ステンレス製クリップ サイズ16 10個入	¥65,900	¥71,600	
21800-000	バット (ナイロ製)	¥22,600	¥24,500	
21800-014	ライトメーター (Lux/Fc両用露出計)	¥188,200	¥350,000	
21815-443	Safety Plunger Cans	¥55,900	¥59,900	
21834-051	粉末洗剤 (非イオン系/洗浄機用) 1.8kg	¥48,500	¥54,600	
21834-075	粉末洗剤 (非イオン系/洗浄機用) 11.3kg	¥171,400	¥183,500	
21834-084	粉末洗剤 (非イオン系/洗浄機用) 22.7kg	¥324,200	¥347,100	
21834-109	錠剤型洗剤 (低発泡性)	¥146,100	¥164,500	
21834-120	Cleaner Luminox LW-FM 1kg 4本	¥274,400	¥293,700	
21899-716	洗浄剤 (Cage-Klenz 100) 5kg	¥153,200	¥168,800	
21902-985	ワipe (コット製) 100枚入	¥33,000	¥35,300	
21903-005	ワipe (ゲイラーEX-L/クリーン用)	¥131,900	¥138,400	
21903-010	ワipe (スコット製)	¥207,000	¥217,300	
21903-012	ワipe (システム/ウェット/水性薬液使用)	¥89,500	¥106,700	
21903-014	ワipe (システム/ウェット/水性薬液使用)	¥163,300	¥171,400	
21903-020	ワipe (ワipe-ルX60/エシ加工済)	¥117,100	¥122,900	
21903-021	ワipe (ゲイラーEX-L/クリーン用)	¥226,100	¥237,300	
21903-022	ワipe (ワipe-ルX60/エシ加工済)	¥100,600	¥110,000	
21903-025	ワipe (プレジジョン/実験室用) 280枚	¥6,700	¥7,000	
21903-026	ワipe (プレジジョン/実験室用)	¥165,300	¥189,800	
21903-040	ワipe (ワipe-ルX60/エシ加工済)	¥49,500	¥51,900	
21905-005	ワipe (ワipe-ル30/エシ加工済) 100枚	¥17,800	¥19,400	
21905-011	ワipe EX-L 12x12cm 196枚	¥10,600	¥11,500	
21905-026	ワipe EX-L/半導体向け対応 280枚	¥5,200	¥5,400	
21905-049	ワipe (ワipe EX-L/半導体向け対応)	¥13,400	¥16,000	
21908-000	ワipe (ワipe-ル40)	¥199,100	¥236,400	
21908-006	ワipe (ワipe-ル30/エシ加工済)	¥165,000	¥173,100	
21908-008	ワipe (ワipe-ル30/エシ加工済)	¥87,900	¥94,100	
21908-023	ワipe (ワipe-ル10)	¥212,400	¥223,000	
21908-205	KIMTECH PURE ワipe	¥124,900	¥131,100	
21908-673	ワipe 無塵紙	¥611,700	¥654,800	
21908-674	ワipe 無塵紙	¥548,400	¥587,100	
21908-763	Haemo-Sol Regular Cleaners	¥24,800	¥27,700	
21908-765	Haemo-Sol Non-sudsing Cleaners	¥152,200	¥185,000	
21909-140	Disposable Biopsy Punches 4.0mm 50	¥113,800	¥119,400	
21909-142	Disposable Biopsy Punches 5.0mm 50	¥118,900	¥88,700	
21909-292	Metzenbaum Scissors 14cm	¥71,500	¥75,100	
21909-404	ワipe ストレート ステンレス製 11.4 cm	¥51,500	¥54,000	
21909-666	Disposable Scalpels ST 10	¥20,600	¥21,600	
21910-009	VWR Premoistened Clean-Wipes 100枚	¥11,900	¥13,000	
21910-105	ワipe (コット製)	¥48,000	¥53,200	
21910-107	ワipe (コット製) 10.1x10.1cm	¥33,500	¥37,400	
21910-111	ワipe (ワipe/ウェット/トワイ)	¥14,300	¥15,800	
21910-118	Cheesecloth Wipes BBA Nonwovens	¥1,091,000	¥1,207,200	
21910-121	Cheesecloth Wipes BBA Nonwovens	¥1,137,400	¥1,808,600	
21912-042	ワipe 用ワipe (MF&E加工済/低発塵)	¥312,800	¥334,800	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
21912-046	クリーニング用ワイパ®-(ホ°リステル製/洗浄済)	¥32,200	¥34,500	
21912-158	Lensx 90 Lens Tissues 1000入	¥33,700	¥36,100	
21912-262	卓上除震台	¥735,000	¥1,363,000	
21913-125	クリーニング用ワイパ®- 500	¥39,700	¥42,500	
21913-147	クリーニング用ワイパ®- 500	¥59,600	¥63,800	
21913-169	クリーニング用ワイパ®- 500	¥129,800	¥139,000	
21913-214	クリーニング用ワイパ®-(ホ°リステル製/洗浄済)	¥29,100	¥31,100	
21913-216	クリーニング用ワイパ®-(ホ°リステル製/洗浄済)	¥49,900	¥53,400	
21913-812	かまろワイパ® 300 30.5x30.5cm 10枚	¥35,200	¥37,800	
21914-000	Ultra-Seal 3000 Wipers 100枚入	¥115,400	¥123,500	
21914-002	クリーニング用ワイパ®-(MF/洗浄済) 100枚	¥197,100	¥211,000	
21914-758	クリーニング用ワイパ®-(MF&セルロ-ス製) 300	¥20,000	¥21,400	
21914-772	クリーニング用ワイパ®-(MF製/滅菌済)クラス10	¥271,600	¥290,700	
21914-777	クリーニング用ワイパ®-(MF&セルロ-ス/滅菌済)	¥238,900	¥255,700	
21921-844	Hazmart Boot Covers, Large 1双	¥13,700	¥15,700	
21923-220	靴カバー/ショート	¥102,600	¥108,600	
22001SPL	1リットル用リザーバ®-PS 100mL ST 1x50袋	¥10,000	¥10,900	
22179-004	シリコナイズ®G-チューブ®(遠心チューブ)1.5ml	¥34,600	¥38,500	
22210T	10ul Eチップ LB S96X10R	¥11,600	¥12,500	
22221-075	テープ®(クリーニング用ビニ-ル製)DarkGreen	¥11,700	¥12,500	
22221-077	テープ®(クリーニング用ビニ-ル製)LightGreen	¥11,700	¥12,500	
22221-091	テープ®(クリーニング用ビニ-ル製)透明	¥11,700	¥12,500	
22228-000	作業用フロアマット・静電気制御品	¥123,900	¥132,600	
22228-002	作業用フロアマット・静電気制御品	¥174,400	¥186,700	
22228-004	作業用フロアマット・静電気制御品	¥261,500	¥280,000	
22230-000	マット(作業用フロアマット)	¥108,800	¥116,500	
22230-002	マット(作業用フロアマット)	¥148,600	¥159,100	
22230-004	マット(作業用フロアマット)	¥230,600	¥246,900	
22231-010	マット(作業用フロアマット)	¥102,900	¥110,100	
22231-012	マット(作業用フロアマット)	¥226,200	¥242,200	
22231-014	マット(作業用フロアマット)	¥159,100	¥170,300	
22231-016	マット(作業用フロアマット)	¥212,100	¥227,000	
22231-018	マット(作業用フロアマット)	¥329,400	¥352,700	
22231-020	マット(作業用フロアマット)	¥292,600	¥313,200	
22231-022	マット(作業用フロアマット)	¥397,400	¥425,500	
22231-024	マット(作業用フロアマット)	¥150,900	¥161,600	
22231-026	マット(作業用フロアマット)	¥176,900	¥189,300	
22231-027	マット(作業用フロアマット)	¥292,600	¥313,200	
22231-028	マット(作業用フロアマット)	¥424,100	¥454,100	
22231-030	マット(作業用フロアマット)	¥265,100	¥283,800	
22231-032	マット(作業用フロアマット)	¥424,100	¥454,100	
22231-034	マット(作業用フロアマット)	¥353,400	¥378,400	
22231-036	マット(作業用フロアマット)	¥1,490,300	¥1,595,500	
22231-038	マット(作業用フロアマット)	¥993,600	¥1,063,700	
22231-050	作業用フロアマット・静電気制御品	¥178,900	¥191,500	
22231-052	作業用フロアマット・静電気制御品	¥238,500	¥255,400	
22231-106	マット(作業用フロアマット)	¥662,400	¥709,100	
22231-108	マット(作業用フロアマット)	¥1,987,100	¥2,127,300	
22231-110	マット(作業用フロアマット)	¥496,800	¥531,800	
22231-112	マット(作業用フロアマット)	¥636,200	¥681,100	
22231-114	マット(作業用フロアマット)	¥993,600	¥1,063,700	
22231-116	マット(作業用フロアマット)	¥2,980,700	¥3,191,000	
22231-120	マット(作業用フロアマット)	¥529,900	¥567,300	
22231-122	マット(作業用フロアマット)	¥662,400	¥709,100	
22231-124	マット(作業用フロアマット)	¥794,900	¥850,900	
22231-126	マット(作業用フロアマット)	¥1,324,700	¥1,418,200	
22231-128	マット(作業用フロアマット)	¥1,987,100	¥2,127,300	
22231-130	マット(作業用フロアマット)	¥3,974,200	¥4,254,600	
22233-176	実験器具棚 ス-パ-エレクタリ-ス®	¥62,100	¥67,100	
22233-178	実験器具棚 ス-パ-エレクタリ-ス®	¥73,400	¥79,300	
22233-494	実験器具保管棚 シェルフ	¥18,900	¥20,400	
22233-496	実験器具保管棚 シェルフ	¥25,300	¥27,300	
22233-498	実験器具保管棚 シェルフ	¥26,700	¥28,800	
22233-504	実験器具保管棚 シェルフ	¥36,900	¥39,900	
22233-640	MetroMAX Stem Casters InterMetro	¥57,600	¥62,200	
22233-642	MetroMAX Stem Casters InterMetro	¥63,600	¥68,700	
22233-644	MetroMAX Stem Casters InterMetro	¥55,400	¥59,800	
22233-690	多目的カート	¥1,000,700	¥1,820,000	
22233-692	多目的カート	¥1,253,500	¥2,280,000	
22233-694	多目的カート	¥1,448,100	¥2,634,000	
22233-696	ネットカート	¥2,556,500	¥2,763,700	
22233-698	ネットカート	¥3,083,200	¥3,333,000	
22233-700	ネットカート	¥2,739,100	¥2,961,100	
22233-716	実験室用 ラック	¥48,800	¥52,800	
22233-718	実験室用 ラック	¥56,300	¥60,900	
22233-720	実験室用 ラック	¥56,300	¥60,900	
22233-722	実験室用 ラック	¥67,800	¥73,300	
22233-724	実験室用 ラック	¥68,700	¥74,300	
22233-726	実験室用 ラック	¥79,900	¥86,400	
22233-728	実験室用 ラック	¥78,900	¥85,300	
22233-730	実験室用 ラック	¥92,600	¥100,100	
22233-732	実験室用 ラック	¥68,900	¥74,500	
22233-734	実験室用 ラック	¥70,800	¥76,600	
22233-736	実験室用 ラック	¥70,800	¥76,500	
22233-738	実験室用 ラック	¥81,800	¥88,400	
22233-740	実験室用 ラック	¥82,500	¥89,200	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
22233-742	実験室用 ラック	¥95,200	¥102,900	
22233-744	実験室用 ラック	¥96,400	¥104,300	
22233-746	実験室用 ラック	¥113,800	¥123,000	
22234-042	マイク遠心チューブ 0.6ml(低吸着) 500個	¥41,100	¥42,900	
22234-044	マイク遠心チューブ 1.5ml(低吸着) 250個	¥26,200	¥27,400	
22234-046	マイク遠心チューブ 1.7ml(低吸着) 250個	¥25,600	¥26,700	
22234-048	マイク遠心チューブ 2.0ml(低吸着) 250個	¥26,100	¥27,200	
22877-252	加二用シリンダ ティスク 6x8	¥51,100	¥53,300	
22877-254	加二用シリンダ ティスク 8x8	¥51,100	¥53,300	
22939-956	TruCLEAN Replacement Bucket レッド	¥58,900	¥63,100	
22939-958	TruCLEAN Replacement Bucket 緑	¥58,600	¥62,800	
22939-960	TruCLEAN Replacement Bucket ブルー	¥58,600	¥62,700	
22939-964	TruCLEAN Stainless Steel Busket	¥154,300	¥165,200	
22939-966	TruCLEAN Supply Pail 8L ブルー	¥36,400	¥38,900	
22939-968	TruCLEAN Supply Pail 7L 緑	¥15,600	¥16,600	
22939-970	TruCLEAN Supply Pail 6L レッド	¥36,700	¥39,300	
22940-008	TruCLEAN(Double-Bucket System)	¥1,512,600	¥1,619,300	
22940-014	TruCLEAN(Triple-Bucket System)	¥1,678,300	¥1,796,700	
22940-020	TruCLEAN Yarn Mop Head 12本	¥186,200	¥199,300	
22940-023	TruCLEAN Mop(Sponge) 25本	¥390,600	¥418,100	
22940-024	ブラシ(ナイロン)	¥15,600	¥16,700	
22940-032	ドライワイパー	¥28,400	¥30,400	
22940-034	ドライワイパー	¥34,700	¥37,200	
22940-036	ドライワイパー 1個	¥40,100	¥42,900	
22940-046	スカーフ	¥6,000	¥6,400	
22940-058	シヨール(ホリフビ製)	¥54,600	¥58,400	
22940-062	シヨール(ホリフビ製)	¥28,300	¥30,300	
22940-064	シヨール(ホリフビ製)	¥28,300	¥30,300	
22940-070	ブラシ(ナイロン)	¥6,800	¥7,300	
22940-072	ツルハカ	¥34,200	¥36,600	
22940-074	ツルハカ	¥6,600	¥7,100	
22940-076	ブラシ(ナイロン)	¥11,500	¥12,300	
22940-077	ブラシ(ナイロン)	¥11,700	¥12,600	
22940-088	ブラシ(ほうき)	¥17,900	¥19,200	
22940-090	ブラシ(ほうき)	¥17,700	¥19,000	
22940-092	ブラシ(ほうき)	¥17,900	¥19,200	
22940-098	Telescopic Extensionホリフビ	¥38,600	¥41,300	
22940-114	ちりとり(ホリフビ製) ブルー	¥16,200	¥17,300	
22940-120	ちりとり(ホリフビ製) レッド	¥15,900	¥17,000	
22940-154	ブラシ(ほうき)	¥17,900	¥19,200	
22940-158	ブラシ(ほうき)	¥17,900	¥19,200	
22940-182	ブラシ(ナイロン)	¥15,600	¥16,700	
22940-186	ブラシ(ナイロン)	¥15,600	¥16,700	
22940-191	TruCLEAN Cover for Sponge Mop Head 25	¥178,600	¥191,200	
22940-304	TruCLEAN Mop(Hydro-Sorb) 25本	¥415,200	¥444,500	
22940-308	TruCLEAN Hydro-Sorb Sponges 25個	¥131,600	¥140,900	
22940-326	TruCLEAN Mop(Fixed Mop Frame)1本	¥114,700	¥122,800	
23006-047N	PTE φ47mm 0.4um 100pc	¥27,200	¥29,000	
23201-012	導電率計用電極 YSI 3400シリーズ	¥412,700	¥758,000	
23201-045	導電率計用電極 YSI 3400シリーズ	¥341,200	¥627,000	
23201-078	導電率計用電極 YSI 3400シリーズ	¥341,200	¥627,000	
23201-103	導電率計用電極 YSI 3400シリーズ	¥233,900	¥430,000	
23201-205	導電率計用電極 YSI 3400シリーズ	¥348,100	¥640,000	
23270-004	コンスタメーター	¥476,700	¥907,000	
23420-027	コルク栓 100個入	¥19,300	¥28,800	
23420-038	コルク栓	¥54,000	¥80,200	
23420-050	コルク栓	¥65,200	¥96,900	
23420-071	コルク栓	¥54,000	¥80,300	
23420-082	コルク栓	¥17,700	¥19,900	
23420-093	コルク栓	¥69,700	¥78,000	
23420-107	コルク栓	¥24,100	¥28,400	
23420-118	コルク栓	¥91,300	¥107,600	
23420-129	コルク栓	¥19,600	¥23,100	
23420-130	コルク栓	¥64,400	¥76,000	
23420-140	コルク栓	¥18,400	¥21,700	
23420-151	コルク栓	¥63,800	¥75,200	
23420-162	コルク栓	¥23,700	¥28,600	
23420-173	コルク栓	¥71,800	¥86,800	
23420-184	コルク栓	¥35,800	¥44,300	
23420-195	コルク栓	¥75,100	¥93,000	
23420-209	コルク栓	¥24,000	¥29,300	
23420-210	コルク栓	¥80,200	¥97,800	
23420-220	コルク栓	¥34,400	¥39,000	
23420-231	コルク栓	¥146,500	¥165,500	
23420-242	コルク栓	¥41,000	¥50,800	
23420-264	コルク栓	¥59,500	¥70,000	
23420-286	コルク栓	¥65,500	¥73,900	
23420-300	コルク栓	¥89,100	¥105,200	
23420-322	コルク栓	¥75,400	¥85,100	
23420-366	コルク栓	¥136,500	¥178,100	
23420-399	コルク栓	¥129,300	¥155,000	
23420-413	コルク栓	¥203,500	¥226,000	
23420-435	コルク栓	¥218,800	¥246,500	
23420-479	コルク栓	¥257,700	¥319,300	
23420-504	コルク栓	¥58,700	¥65,600	
23420-515	コルク栓	¥291,800	¥326,100	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
23420-526	コルク栓	¥78,700	¥88,000	
23420-559	コルク栓	¥495,000	¥553,300	
23420-708	Cork Sheet 101x305mm	¥35,900	¥39,900	
23457-060	コルクロー	¥125,700	¥135,200	
23570T	10ul ウルトラGチップ LB S 96X10R	¥12,700	¥13,600	
23580T	10ul ウルトラGチップ LB 1000/bag	¥4,700	¥5,000	
23609-102	携帯用計数カウンタ	¥27,200	¥30,000	
23609-172	VWR Digital Surface Thermometer	¥47,300	¥52,200	
23609-182	スプレートル 240ml PP 3	¥19,200	¥20,700	
23609-196	VWR 99M/59S Timer	¥20,800	¥22,900	
23609-230	デジタル温度計(2チャンネル)	¥153,100	¥168,200	
23649-006	Howard Mold Counting Chamber	¥256,700	¥284,200	
23649-030	血球計算盤	¥38,100	¥42,100	
23649-032	血球計算盤	¥38,900	¥43,000	
23984-000	灰吸法ト	¥91,000	¥95,500	
23984-066	灰吸法ト	¥155,600	¥163,300	
24065-000	Walter Crucible Holder 6個	¥54,800	¥58,900	
24301-704	Quick Serve Tap(液体リフト用)	¥4,400	¥4,800	
24631-446	加圧チップ 96孔 16x125mm 1000	¥256,300	¥279,500	
24780-369	双シリンダ(ホリオン製)2000ml 1本	¥42,400	¥45,400	
24780-371	双シリンダ(ホリオン製)1000ml 1本	¥21,200	¥22,700	
24980-005	デシケータ(ステンレス製)	¥746,500	¥1,358,000	
24980-107	デシケータ(ステンレス製)	¥565,000	¥998,000	
24981-008	デシケータ	¥1,800,300	¥3,275,000	
24982-010	デシケータ(バキューム用、ガラス扉強化材)	¥249,600	¥196,600	
24982-012	デシケータ(バキューム用、ガラス扉強化材)	¥161,700	¥124,000	
24983-331	デシケータ(積み重ね型)	¥307,800	¥565,000	
24987-048	デシケータ 250MM	¥225,900	¥446,000	
24988-164	デシケータ(小型/真空式)	¥97,000	¥104,400	
24988-233	デシケータ(小型/真空式)	¥62,500	¥67,200	
24988-255	デシケータ(小型/真空式)	¥81,000	¥87,200	
24988-277	デシケータ(小型/真空式)	¥106,100	¥114,200	
25000-134	VWR CryoPro紙製チップス 81本立	¥12,600	¥13,600	
25006-047N	PA φ47mm 0.45um 100pc	¥22,800	¥27,000	
25007-047N	PA φ47mm 0.2um 100pc	¥22,800	¥27,000	
25020-020	デシケータ	¥95,900	¥105,700	
25020-041	デシケータ	¥116,100	¥128,000	
25337-059	蒸発皿(テフロン樹脂)	¥37,100	¥41,500	
25337-070	蒸発皿(テフロン樹脂)	¥60,000	¥69,500	
25384-090	VWR® 96ウェル イソスタック Semi-Stackable ST 50	¥63,000	¥66,100	
25384-092	96ウェル イソスタック ST 60x15mm 500	¥68,000	¥71,300	
25384-093	96ウェル イソスタック 96ウェルプレートST60x15mm 500	¥85,100	¥98,300	
25384-144	試料コンテナ(滅菌済 100個包装)	¥21,600	¥25,200	
25384-146	試料コンテナ(滅菌済 10個包装)	¥22,600	¥26,300	
25384-148	試料コンテナ(未滅菌 100個包装)	¥57,600	¥67,100	
25384-150	ピペタ(ホリオン製) 50ml 100個入	¥28,600	¥30,800	
25384-152	ピペタ(ホリオン製) 100ml 100個入	¥31,500	¥36,600	
25384-154	ピペタ(ホリオン製) 250ml 100個入	¥61,200	¥65,100	
25384-156	ピペタ(ホリオン製) 400ml 100個入	¥56,900	¥65,400	
25384-158	ピペタ(ホリオン製) 800ml 100個入	¥68,500	¥75,200	
25384-160	ピペタ(ホリオン製) 1000ml 100個入	¥106,500	¥116,800	
25384-302	96ウェル イソスタック Slippable 100x15 500個入	¥84,000	¥88,100	
25384-304	Petri Dishes Segmented ST 20	¥99,600	¥104,600	
25384-308	VWR® 96ウェル イソスタック ST X-Plate 20枚x25袋	¥105,000	¥110,300	
25384-326	VWR® 96ウェル イソスタック Slippable ST 100入	¥61,500	¥64,500	
25384-342	96ウェル イソスタック ST 100x15mm 500	¥83,600	¥87,700	
25395-251	Chemware Chemfluor PFA Petri Dish	¥67,600	¥75,600	
25395-380	複製用ツール Replica-Plating Tool	¥60,700	¥65,600	
25395-382	複製用ツール(綿拭き) 36枚入	¥63,000	¥67,500	
25405-225	96ウェル イソスタック用チップ	¥27,900	¥29,900	
25405-228	96ウェル イソスタック用チップ	¥74,100	¥79,400	
25415-061	アルミニウムイソスタック(蓋付き)	¥161,000	¥173,200	
25415-083	アルミニウムイソスタック(蓋付き)	¥168,000	¥180,700	
25420-028	アルミニウムイソスタック(蓋付き)	¥87,500	¥94,200	
25420-040	アルミニウムイソスタック(蓋付き)	¥111,100	¥119,500	
25420-061	アルミニウムイソスタック(蓋付き)	¥135,400	¥145,600	
25420-072	アルミニウムイソスタック(蓋付き)	¥135,700	¥145,900	
25420-083	アルミニウムイソスタック(蓋付き)	¥138,400	¥148,900	
25433-008	秤量皿(アルミ)	¥110,300	¥132,800	
25433-010	秤量皿(アルミ)	¥74,200	¥86,800	
25433-014	秤量皿(アルミ/丸型) 70mL	¥10,300	¥11,900	
25433-016	秤量皿(アルミ/丸型) 20mL	¥8,500	¥10,000	
25433-018	秤量皿(アルミ/丸型) 500mL	¥12,900	¥14,100	
25433-020	秤量皿(アルミ/丸型) 75mL	¥8,400	¥10,200	
25433-022	秤量皿(アルミ/丸型) 200mL	¥9,000	¥10,600	
25433-030	VWR秤量皿 140mL 100枚	¥25,500	¥30,400	
25433-052	秤量皿(アルミ) 20mL 100入	¥10,400	¥11,400	
25433-085	秤量皿(アルミ)	¥104,500	¥125,900	
25433-089	秤量皿(アルミ)	¥113,600	¥132,200	
25460-907	染色用スライド Case of 4	¥60,800	¥67,600	
25461-010	Staining Dish 20-Slide 3/Case	¥86,900	¥93,500	
25461-016	Staining Dish 16/30/40-Slide	¥196,000	¥210,600	
25485-031	連続分注チップ	¥25,400	¥27,400	
25500-800	マイクロチップ 7673A 0.5ul	¥135,900	¥146,900	
25500-802	マイクロチップ 7673A 1.0ul	¥92,700	¥100,200	
25500-806	マイクロチップ 7673A 1.0ul	¥92,300	¥99,700	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
25500-808	マイクロシリンジ 7673A 0.5ul	¥135,900	¥146,900	
25500-810	マイクロシリンジ 7673A 10ul	¥314,200	¥339,700	
25500-812	マイクロシリンジ 7673A 10ul	¥54,600	¥59,100	
25607-109	ピンセットストレート 200mm 1	¥6,300	¥6,600	
25607-120	ピンセット 200mm 1	¥14,400	¥15,100	
25607-186	ピンセット(Micro-Adson) 120mm 1本	¥5,000	¥5,200	
25607-302	止血鉗子(Hemostatic)12.5cm 1	¥10,000	¥10,500	
25607-380	止血鉗子/Rochester-Ochsner14cm 1	¥14,700	¥15,500	
25607-470	止血鉗子(Rochester-Pean)16cm 1	¥17,000	¥17,800	
25607-517	止血鉗子(Rochester-Pean)16cm 1	¥9,700	¥10,200	
25607-608	ピンセット(組織用) 145mm ステンレス製 1	¥3,800	¥4,000	
25607-925	スカルパノーム #7 140mm ステンレス製 1本	¥4,600	¥4,800	
25607-947	スカルパノーム #3 ステンレス製 102mm 1	¥3,400	¥3,600	
25608-065	刃用替刃 #10 2.2cm ST 個包装 100	¥40,200	¥44,300	
25608-087	刃用替刃 #11 20mm ST 個包装 100	¥35,100	¥38,600	
25608-203	ハサミ ストレート 115mm 1	¥7,300	¥7,700	
25608-225	ハサミ Curved 115mm 1	¥4,500	¥4,700	
25608-280	ハサミ ストレート ステンレス製 約130mm 1	¥12,200	¥12,800	
25608-316	ハサミ ストレート ステンレス製 153mm 1	¥13,200	¥13,800	
25608-360	ハサミ ストレート ステンレス製 127mm 1	¥13,400	¥14,100	
25608-509	ハサミ 反刃 ステンレス製 130mm 1	¥5,000	¥5,200	
25608-575	ハサミ(スプリングハルト式) ステンレス 98mm 1	¥20,100	¥21,100	
25608-586	ハサミ(スプリングハルト式) ステンレス 114mm 1	¥59,300	¥62,200	
25608-868	Tissue-Tek Slide Holders 24place 1	¥30,100	¥32,600	
25608-874	スライドメイカー	¥93,200	¥100,900	
25608-916	Tissue-TekCryomoldMoldsAdapter 100	¥22,300	¥24,200	
25608-922	Tissue-TekCryomoldMoldsAdapter 100	¥22,500	¥24,400	
25608-944	マイクロトームクリーニングシステム	¥138,100	¥149,700	
25608-946	マイクロトームクリーニングシステム	¥76,900	¥83,100	
25609-006	スカルパノーム 2inch 1	¥72,000	¥79,700	
25609-008	スカルパノーム #61 ST 20枚	¥86,700	¥96,000	
25609-010	スカルパノーム #62 ST 20枚	¥76,800	¥85,000	
25610-001	解剖用具セット(ステンレス製) 6品セット	¥12,400	¥13,100	
25640-002	解剖用具セット(ステンレス製) 9品セット	¥9,900	¥10,400	
25672-100	ピンセット(計量用)	¥3,500	¥3,800	
25718-008	ピンセット(解剖用) 150mm 5本入	¥25,900	¥27,100	
25718-012	ピンセット(ストレートタイプ 解剖用) 117mm	¥11,200	¥12,000	
25719-044	ピンセット(標本用) 114mm 1	¥3,000	¥3,100	
25719-066	ピンセット(標本用) 140mm 1	¥5,400	¥5,600	
25729-081	ピンセット 120mm	¥6,000	¥6,300	
25729-387	ピンセット(標本用) 305mm 1	¥12,900	¥13,500	
25729-445	ピンセット(標本用) 152mm 1	¥4,800	¥5,000	
25729-503	ピンセット(標本用) 203mm 1	¥6,700	¥7,100	
25729-569	ピンセット(標本用) 254mm 1	¥9,700	¥10,200	
25729-627	ピンセット(標本用 用型) 258mm 1	¥7,500	¥7,900	
25778-000	ニードル(解剖用) 144mm ステンレス製 12本	¥2,200	¥2,400	
25872-008	ハサミ ストレート ステンレス製 115mm 1	¥2,100	¥2,300	
25874-003	ハサミ FineBent ステンレス製 115mm 1	¥3,600	¥3,800	
25874-105	ハサミ ストレート ステンレス製 140mm 1	¥8,800	¥9,200	
25876-009	ハサミ ストレート ステンレス製 約115mm 1	¥1,800	¥1,900	
25877-103	ハサミ ストレート ステンレス製 140mm 1	¥10,300	¥10,800	
25878-026	ハサミ ストレート ステンレス製 1	¥1,600	¥1,700	
25878-048	ハサミ ストレート ステンレス製 140mm 1	¥6,500	¥6,900	
25878-106	ハサミ ストレート ステンレス製 165mm 1	¥7,700	¥8,100	
25881-009	止血鉗子(解剖用)140mm ステンレス製 1	¥4,500	¥4,700	
25940-006	シカーフローブ 150mm 3本入	¥4,900	¥5,200	
25942-001	シカーフローブ 125mm 3本	¥4,800	¥5,000	
26316-694	フィルター(フィルターホルダー)	¥36,300	¥39,500	
26316-698	フィルター(フィルターホルダー)	¥29,700	¥32,300	
26316-702	フィルター(フィルターホルダー)	¥7,100	¥7,700	
26316-710	フィルター(フィルターホルダー)	¥46,800	¥51,100	
26316-712	フィルター(フィルターホルダー)	¥81,500	¥94,400	
26316-726	フィルター(フィルターホルダー)	¥148,100	¥161,500	
26316-740	フィルター(フィルターホルダー)	¥110,600	¥120,600	
26396-506	ガラスビーズ	¥14,700	¥16,400	
26396-508	ガラスビーズ	¥11,200	¥12,500	
26396-521	ガラスビーズ	¥10,600	¥11,800	
26396-563	ガラスビーズ 4mm 1本セット	¥11,300	¥12,600	
26396-596	ガラスビーズ	¥10,500	¥11,600	
26396-621	ガラスビーズ	¥9,900	¥11,100	
26396-630	ガラスホール	¥87,800	¥98,000	
26396-633	ガラスホール	¥86,900	¥96,900	
26396-636	ガラスホール	¥80,700	¥90,000	
26396-639	ガラスホール	¥78,700	¥87,800	
26397-103	PTFE Boiling Stones	¥115,700	¥129,500	
26397-409	多孔質沸騰石(226g容器入り)	¥12,200	¥13,300	
26405-110	Element For HG301A	¥30,200	¥32,300	
26410T	10ul ロングリチ NT LB ST ヒンジ 96x10R	¥13,800	¥14,800	
26430T	200ul ソフトフィットチップ NT LB ST ヒンジ 96x10R	¥10,300	¥11,100	
26674-451	Flask 24/40 1000mL 1個	¥405,400	¥311,500	
26675-230	フラスコ(ファストフリーズ)600ml	¥129,100	¥108,800	
26675-319	Fast-Freeze Flasks 900mL 1個	¥263,200	¥217,500	
26678-008	真空ホールド(オレガリ二段型)	¥1,181,800	¥2,150,000	
26678-012	真空ホールド(オレガリ二段型)	¥1,308,600	¥2,679,000	
26678-020	真空ホールド(オレガリ二段型)	¥2,138,500	¥4,378,000	
26678-056	真空ホールド(オレガリ二段型)	¥470,300	¥1,221,000	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
27148-022	マット(粘着タイプ) 24x36 60枚重ねx4入	¥225,000	¥240,900	
27373-406	Owl Gel Staining Box	¥50,600	¥57,100	
27770T	200ul チップステーション リーチアップ LB ST96X50	¥43,900	¥47,200	
27790T	10ul チップステーション リーチアップ LB ST96X50	¥51,500	¥55,500	
28130-025	ガラスヤすり	¥77,400	¥82,900	
28130-047	ガラスヤすり	¥84,100	¥91,300	
28130-069	ガラスヤすり	¥82,300	¥88,100	
28130-081	ガラスヤすり	¥124,700	¥133,500	
28200T	10ul 汎用フィルター Gチップ LB S96x10R	¥19,600	¥21,100	
28297-082	フィルター(定性濾紙)	¥18,100	¥19,600	
28297-289	ファイルーグラス	¥24,600	¥28,000	
28297-500	ファイルーグラス	¥33,600	¥36,500	
28297-504	ファイルーグラス	¥21,600	¥23,400	
28297-506	ファイルーグラス	¥32,400	¥35,100	
28297-508	ファイルーグラス	¥30,600	¥33,200	
28297-510	ファイルーグラス	¥39,700	¥43,100	
28297-826	フィルター(定性濾紙)	¥9,500	¥10,300	
28297-846	フィルター(定性濾紙)	¥24,900	¥27,000	
28297-868	フィルター(定性濾紙) 5.5cm 50枚入	¥2,600	¥2,800	
28297-930	フィルター(定量濾紙)	¥11,400	¥12,400	
28297-974	ファイルーグラス	¥34,200	¥38,900	
28297-976	ファイルーグラス	¥12,100	¥13,100	
28297-978	ファイルーグラス	¥9,600	¥10,400	
28297-982	ファイルーグラス	¥12,200	¥13,200	
28297-984	ファイルーグラス	¥21,500	¥23,300	
28297-986	フロテイングプレート 13x20cm	¥75,400	¥81,700	
28297-988	フロテイングプレート 16x18cm	¥61,800	¥67,000	
28297-990	フロテイングプレート 14x20cm	¥45,500	¥49,300	
28297-992	フロテイングプレート 14x16cm	¥26,900	¥29,100	
28297-994	フロテイングプレート 11x18cm	¥32,600	¥35,300	
28297-996	フロテイングプレート 11x14cm	¥30,400	¥32,900	
28297-998	フロテイングプレート 11x13cm 厚さ2.6mm	¥21,800	¥23,700	
28298-010	フロテイングプレート 10x14cm	¥31,000	¥33,700	
28298-014	フロテイングプレート 7x10cm	¥22,200	¥24,000	
28298-016	フロテイングプレート 7x8cm	¥29,200	¥31,700	
28298-018	フロテイングプレート	¥24,500	¥26,500	
28298-020	フロテイングプレート - 46x57cm 100枚	¥104,900	¥113,700	
28298-022	フロテイングプレート - 20x20cm 100枚入	¥42,900	¥46,500	
28298-026	フロテイングプレート - 14x16cm 100枚入	¥18,600	¥20,100	
28298-028	フロテイングプレート -	¥23,300	¥25,300	
28298-030	フロテイングプレート - 7x10cm 100枚入	¥14,700	¥15,900	
28298-032	フロテイングプレート 20x20cm	¥92,600	¥100,400	
28306-051	フィルター(定量濾紙)	¥5,700	¥6,200	
28306-095	フィルター(定量濾紙)	¥9,300	¥10,100	
28306-110	フィルター(定量濾紙)	¥10,900	¥11,900	
28306-131	フィルター(定量濾紙)	¥13,300	¥14,400	
28306-153	フィルター(定量濾紙)	¥16,700	¥18,100	
28306-175	フィルター(定量濾紙)	¥19,400	¥21,000	
28306-346	フィルター(定量濾紙)	¥7,700	¥8,300	
28306-368	フィルター(定量濾紙)	¥10,000	¥10,800	
28306-380	フィルター(定量濾紙)	¥10,800	¥11,700	
28306-404	フィルター(定量濾紙)	¥13,700	¥14,800	
28307-134	フィルター(定量濾紙)	¥12,400	¥13,500	
28309-967	フィルター(定性濾紙)	¥2,700	¥3,000	
28309-989	フィルター(定性濾紙)	¥2,900	¥3,100	
28310-015	フィルター(定性濾紙)	¥3,100	¥3,400	
28310-026	フィルター(定性濾紙)	¥3,500	¥3,800	
28310-048	フィルター(定性濾紙) 9.0cm	¥3,000	¥3,200	
28310-081	フィルター(定性濾紙)	¥3,300	¥3,600	
28310-106	フィルター(定性濾紙)	¥4,300	¥4,700	
28310-128	フィルター(定性濾紙)	¥5,400	¥6,200	
28310-140	フィルター(定性濾紙)	¥9,300	¥10,000	
28310-161	フィルター(定性濾紙)	¥7,900	¥9,400	
28310-172	フィルター(定性濾紙)	¥14,700	¥16,000	
28310-208	フィルター(定性濾紙)	¥16,100	¥19,200	
28311-029	フィルター(定性濾紙)	¥144,200	¥172,000	
28313-046	フィルター(定性濾紙)	¥2,600	¥2,800	
28313-057	フィルター(定性濾紙)	¥3,300	¥3,500	
28313-068	フィルター(定性濾紙)	¥3,800	¥4,100	
28313-080	フィルター(定性濾紙)	¥4,700	¥5,100	
28313-104	フィルター(定性濾紙)	¥5,400	¥5,800	
28313-126	フィルター(定性濾紙)	¥6,900	¥7,500	
28313-148	フィルター(定性濾紙)	¥7,000	¥7,600	
28313-160	フィルター(定性濾紙)直径25cm 50枚入	¥11,200	¥12,100	
28313-181	フィルター(定性濾紙)	¥15,100	¥16,400	
28313-240	フィルター(定性濾紙)	¥34,400	¥41,000	
28317-472	UNIPLATE Microplates 2mL 25枚	¥162,000	¥177,300	
28320-019	フィルター(定性濾紙)	¥2,900	¥3,500	
28320-020	フィルター(定性濾紙)	¥3,100	¥3,300	
28320-041	フィルター(定性濾紙)	¥3,200	¥3,500	
28320-085	フィルター(定性濾紙)	¥3,900	¥4,200	
28320-100	フィルター(定性濾紙) 100入	¥5,100	¥5,800	
28320-121	フィルター(定性濾紙)	¥6,200	¥6,700	
28320-143	フィルター(定性濾紙)	¥9,000	¥10,700	
28320-187	フィルター(定性濾紙)	¥10,700	¥11,600	
28320-201	フィルター(定性濾紙)	¥13,900	¥15,100	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
28320-223	フィルター(定性濾紙)	¥20,200	¥24,100	
28320-406	抽出用円筒ろ紙	¥83,200	¥91,700	
28320-500	抽出用円筒ろ紙	¥83,100	¥91,600	
28320-506	抽出用円筒ろ紙	¥95,000	¥104,600	
28320-510	抽出用円筒ろ紙	¥122,500	¥135,000	
28321-077	フィルター(定性濾紙)	¥4,400	¥4,800	
28321-113	フィルター(定性濾紙)	¥5,800	¥6,200	
28321-135	フィルター(定性濾紙)	¥8,600	¥9,400	
28331-026	フィルター(折りたたみ定性濾紙)	¥24,400	¥26,500	
28331-048	フィルター(折りたたみ定性濾紙)	¥24,600	¥26,600	
28331-060	フィルター(折りたたみ定性濾紙)	¥30,100	¥32,600	
28331-081	フィルター(折りたたみ定性濾紙)	¥30,700	¥33,200	
28331-106	フィルター(折りたたみ定性濾紙)	¥33,300	¥36,200	
28331-128	フィルター(折りたたみ定性濾紙)	¥57,700	¥62,500	
28333-021	フィルター(折りたたみ定性濾紙)	¥22,400	¥24,300	
28333-043	フィルター(折りたたみ定性濾紙)	¥24,400	¥26,400	
28333-065	フィルター(折りたたみ定性濾紙) 18.5cm	¥26,300	¥30,000	
28333-087	フィルター(折りたたみ定性濾紙)	¥29,400	¥31,900	
28333-101	フィルター(折りたたみ定性濾紙)	¥32,200	¥37,300	
28333-123	フィルター(折りたたみ定性濾紙)	¥35,700	¥38,700	
28333-125	ファイバークラス	¥55,000	¥59,600	
28333-127	ファイバークラス	¥40,300	¥46,200	
28333-129	ファイバークラス 100入	¥12,600	¥13,700	
28333-133	ファイバークラス	¥39,500	¥45,500	
28333-135	ファイバークラス	¥20,800	¥22,600	
28333-137	ファイバークラス	¥12,600	¥13,700	
28333-139	ファイバークラス 100入	¥10,600	¥11,500	
28334-172	プロテイングバット 14x16 50	¥33,900	¥37,400	
28334-174	プロテイングバット 16x18cm	¥39,300	¥43,400	
28342-010	FilterPaper 9cm BLK PK100	¥26,000	¥28,600	
28456-043	Grade3 Filter Paper 径7cm 100枚入	¥14,200	¥15,600	
28496-174	フィルター(ガラス繊維濾紙)	¥47,300	¥52,100	
28496-176	フィルター(ガラス繊維濾紙)	¥64,200	¥70,700	
28496-178	フィルター(ガラス繊維濾紙)	¥117,800	¥129,800	
28496-180	フィルター(ガラス繊維濾紙)	¥133,000	¥146,600	
28496-182	フィルター(ガラス繊維濾紙)	¥191,300	¥210,800	
28496-184	フィルター(ガラス繊維濾紙)	¥347,500	¥382,900	
28511-007	顕微鏡用アセタリ(ワイル)	¥52,400	¥60,700	
28641-001	水流アシレーター	¥76,600	¥137,000	
28641-010	保持リング	¥20,100	¥19,100	
28641-012	保持リング	¥27,000	¥26,300	
28641-014	保持リング	¥37,000	¥35,000	
29442-396	96well Microtitration Plate 360uL	¥144,800	¥152,000	
29630-000	KIMAX Volumetric Flasks 110mL	¥111,800	¥121,900	
29700-004	ワスリ用重し(鉛製) 200g	¥16,000	¥17,600	
29700-026	ワスリ用重し(鉛製) 500g	¥17,200	¥19,000	
29700-048	ワスリ用重し(鉛製) 700g	¥20,400	¥22,400	
29700-060	ワスリ用重し(鉛製) 900g	¥21,900	¥24,100	
29700-066	ワスリ用重し(リットル型) 1	¥29,300	¥31,700	
29876-001	Flow Monitor Water	¥26,300	¥33,000	
30002-014	Pipets 1ml GRAD 500本入	¥63,900	¥71,500	
30002-116	L字型 Plastic Spreaders 1/Peel pk	¥111,700	¥124,900	
30002-118	L字型 Plastic Spreaders 5/Peel pk	¥171,900	¥192,200	
30002-120	L字型 Plastic Spreaders 10/Peel pk	¥62,900	¥70,400	
30180-030	VWR ホスト丸底 PC製 500mL 12	¥71,500	¥77,000	
30231-080	漏斗(ろ過用)	¥100,900	¥109,000	
30249-530	Urbanti High-Speed Filter Funnel 6	¥20,200	¥22,100	
30252-951	粉末漏斗(PP製)	¥28,200	¥36,300	
30252-953	粉末漏斗(PP製)	¥32,600	¥38,200	
30252-955	粉末漏斗(PP製)	¥20,300	¥23,800	
30252-957	粉末漏斗(PP製)	¥17,100	¥22,000	
30256-004	粉末漏斗(PP製 高耐久性)	¥49,900	¥59,300	
30256-005	粉末漏斗(PP製 高耐久性)	¥65,700	¥70,300	
30256-006	粉末漏斗(PP製 高耐久性)	¥40,000	¥42,800	
30256-007	粉末漏斗(PP製 高耐久性)	¥44,600	¥51,000	
30286-026	漏斗(ステンレス製)	¥41,600	¥48,900	
30286-048	漏斗(ステンレス製)	¥45,200	¥53,200	
30286-060	漏斗(ステンレス製) 660mL	¥47,300	¥55,600	
30286-081	漏斗(ステンレス製)	¥51,600	¥60,600	
30286-503	漏斗(PE製 工業用)	¥72,700	¥78,200	
30287-029	秤量皿漏斗(ガラス製)	¥59,400	¥73,300	
30287-040	秤量皿漏斗(ガラス製)	¥99,000	¥108,600	
30306-407	Buchner Funnels, Two Piece, 315mL	¥18,500	¥19,800	
30308-090	漏斗(フナー型 ろ過用 陶器製)	¥298,600	¥321,800	
30308-092	漏斗(フナー型 ろ過用 陶器製)	¥737,900	¥790,000	
30308-094	漏斗(フナー型 ろ過用 陶器製)	¥1,248,700	¥1,376,100	
30309-000	漏斗(フナー型 ろ過用 陶器製)	¥94,600	¥101,300	
30309-002	漏斗(フナー型 ろ過用 陶器製)	¥174,500	¥186,800	
30309-004	漏斗(フナー型 ろ過用 陶器製)	¥223,500	¥242,400	
30360T	Sharp100ulマイクロフィルター-NT LB ステンレス 96x10	¥19,100	¥20,600	
30370T	20ul マイクロフィルター- NT LB ST ステンレス 96x10R	¥19,100	¥20,600	
30380T	100ul V マイクロフィルター- NT LB ST ステンレス 96x10R	¥19,100	¥20,600	
30390T	200ulVマイクロフィルター- NT LB ST ステンレス 96x10	¥19,100	¥20,600	
30430T	200ul クイックラック リフトチップ LB 96X50	¥38,600	¥41,500	
30450T	200ul クイックラック リフトチップ LB 96X50	¥37,300	¥40,100	
30470T	10ul クイックラック Gチップ LB 96X50	¥38,300	¥41,200	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
30510T-X	300XT フィルターチップ° LB S96X10R	¥19,100	¥19,600	
30510-X	300XT フィルターチップ° NT S96X10R	¥18,400	¥18,900	
30534-905	分液ロトトラック	¥219,500	¥235,000	
30550T	200NX フィルターチップ° LB S96X10R	¥19,100	¥20,600	
30587-000	ロト台	¥199,600	¥218,200	
30617-002	Clear Protective Tape Rolls 3入	¥41,200	¥44,500	
30617-012	Write-On 3ハ°リンクテープ°	¥19,000	¥20,400	
30617-020	Write-On 3ハ°リンクテープ° 3本	¥18,600	¥20,000	
30617-024	Write-On 3ハ°リンクテープ° 6本	¥47,800	¥51,400	
30617-028	Write-On 3ハ°リンクテープ° 6本	¥48,300	¥51,900	
30617-030	Write-On 3ハ°リンクテープ°	¥48,300	¥51,900	
30617-036	Write-On 3ハ°リンクテープ°	¥42,700	¥46,000	
30617-048	Write-On 3ハ°リンクテープ°	¥40,100	¥43,100	
30617-608	Clavies Autoclavable Bags 100枚	¥55,100	¥58,700	
30617-644	Poxygrid 50mL Centrifuge Tube Rack	¥50,600	¥54,200	
30617-648	ガラスディスク°サドル°ボックス	¥46,600	¥52,600	
30620-678	FLASK BRUSH NYLON 1000ML 1本	¥4,900	¥20,500	
30621-204	Clear Straight Sided Jars 24	¥57,300	¥62,300	
30623-042	Vented Screw Caps 未滅菌 Chemglass	¥4,700	¥5,600	
30624-022	STOPPER BUTYL RUBBER PK100	¥54,000	¥64,600	
30624-048	液体洗剤 7X 1ガ°ロ×4入	¥126,600	¥165,400	
30800T	300ul XTチップ° LB 1000/bag	¥5,000	¥5,400	
30930T	300ul XTチップ° LB S96X10R	¥10,700	¥11,500	
31759-021	U-Tube気圧計	¥318,100	¥345,600	
31759-043	U-Tube気圧計	¥347,700	¥377,700	
31759-065	U-Tube気圧計	¥382,200	¥415,200	
31759-167	U-Tube気圧計	¥681,900	¥1,246,000	
31759-189	U-Tube気圧計	¥130,100	¥141,400	
31759-203	U-Tube気圧計	¥130,100	¥141,400	
31759-225	U-Tube気圧計	¥151,900	¥165,000	
31759-280	U-Tube気圧計	¥225,800	¥245,300	
31759-327	U-Tube気圧計	¥354,700	¥385,300	
31770T	300ul NXチップ° LB 1000/bag	¥5,000	¥5,400	
31790T	300ul NXチップ° LB S96X10R	¥10,400	¥11,200	
32800-012	デシケーター	¥396,300	¥735,000	
32800-042	デシケーター	¥288,800	¥536,000	
32800-044	デシケーター	¥480,000	¥692,000	
32800-046	デシケーター	¥424,000	¥786,000	
32829-024	シリコンチューブ°	¥81,400	¥85,400	
32829-026	シリコンチューブ°	¥86,700	¥91,000	
32829-028	シリコンチューブ°	¥96,000	¥100,700	
32829-030	シリコンチューブ°	¥111,000	¥116,500	
32829-032	シリコンチューブ°	¥124,000	¥130,100	
32829-034	シリコンチューブ°	¥162,300	¥170,400	
32829-036	シリコンチューブ°	¥134,300	¥141,000	
32829-038	シリコンチューブ°	¥138,300	¥145,100	
32829-040	シリコンチューブ°	¥159,800	¥167,800	
32829-044	シリコンチューブ°	¥410,900	¥431,300	
32829-062	シリコンチューブ°	¥109,800	¥115,200	
32829-104	シリコンチューブ°	¥189,200	¥198,600	
32829-182	Platinum-Cured Silicone Tubing 1	¥51,500	¥54,000	
32860-010	交換ホイル(ガラス管カッター用) 10入	¥164,100	¥215,800	
32861-011	ガラスカッター	¥146,000	¥158,600	
32861-012	ガラスカッター	¥86,000	¥93,400	
32885-780	耐熱/耐寒グロ°フ° ORP ホ°リフ(クリ°ン°ル°ム)	¥352,500	¥370,000	
32885-805	耐熱性グロ°フ° ハ°イ°ル°グ°ラ°ト° 33cm 1双	¥29,200	¥32,800	
32887-904	耐薬品製手袋(7°ポ°ン°製)	¥28,400	¥30,400	
32887-906	耐薬品製手袋(7°ポ°ン°製)	¥29,700	¥31,800	
32887-908	耐薬品製手袋(7°ポ°ン°製)	¥29,700	¥31,800	
32887-910	耐薬品製手袋(7°ポ°ン°製)	¥31,600	¥33,800	
32887-912	耐薬品製手袋(7°ポ°ン°製)	¥30,200	¥32,300	
32887-914	耐薬品製手袋(7°ポ°ン°製)	¥23,000	¥24,700	
32887-916	耐薬品製手袋(7°ポ°ン°製)	¥24,000	¥25,700	
32887-918	耐薬品製手袋(7°ポ°ン°製)	¥22,900	¥24,500	
32887-920	耐薬品製手袋(7°ポ°ン°製)	¥24,900	¥26,700	
32887-921	耐薬品製手袋(7°ポ°ン°製)	¥28,600	¥31,600	
32887-922	耐薬品製手袋(7°ポ°ン°製)	¥22,300	¥23,900	
32887-924	耐薬品製手袋(7°ポ°ン°製)	¥31,800	¥35,300	
32887-926	耐薬品製手袋(7°ポ°ン°製)	¥25,400	¥24,900	
32887-928	耐薬品製手袋(7°ポ°ン°製)	¥33,500	¥37,100	
32887-929	耐薬品製手袋(7°ポ°ン°製)	¥22,000	¥24,000	
32887-930	耐薬品製手袋(7°ポ°ン°製)	¥33,500	¥37,100	
32887-931	耐薬品製手袋(7°ポ°ン°製)	¥22,000	¥24,000	
32887-932	耐薬品製手袋(7°ポ°ン°製)	¥33,500	¥37,100	
32887-935	耐薬品製手袋(7°ポ°ン°製)	¥28,100	¥31,100	
32887-936	耐薬品製手袋(7°ポ°ン°製)	¥22,000	¥21,500	
32887-938	耐薬品製手袋(7°ポ°ン°製)	¥44,300	¥47,900	
32887-940	耐薬品製手袋(7°ポ°ン°製)	¥45,000	¥48,700	
32887-942	耐薬品製手袋(7°ポ°ン°製)	¥46,200	¥49,900	
32887-944	耐薬品製手袋(7°ポ°ン°製)	¥43,200	¥46,700	
32887-946	耐薬品製手袋(7°ポ°ン°製)	¥38,600	¥41,800	
32887-948	耐薬品製手袋(7°ポ°ン°製)	¥42,300	¥46,900	
32887-949	耐薬品製手袋(7°ポ°ン°製)	¥39,900	¥44,200	
32887-952	耐薬品製手袋(7°ポ°ン°製)	¥45,300	¥44,400	
32887-954	耐薬品製手袋(7°ポ°ン°製)	¥47,300	¥51,700	
32887-956	耐薬品製手袋(7°ポ°ン°製)	¥62,100	¥68,800	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
32887-958	耐薬品製手袋(フッ素樹脂製)	¥61,200	¥67,800	
32887-962	耐薬品製手袋(フッ素樹脂製)	¥38,900	¥38,100	
32887-970	耐薬品製手袋(フッ素樹脂製)	¥41,200	¥45,000	
32887-972	耐薬品製手袋(フッ素樹脂製)	¥54,000	¥59,800	
32888-020	CHEMSOFT CEニトリルグローブ XS 100枚	¥27,900	¥29,800	
32888-024	CHEMSOFT CEニトリルグローブ M 100枚	¥27,100	¥29,000	
32888-202	AK NaturalLatexCleanroom Gloves 10	¥79,900	¥86,500	
32888-204	AK NaturalLatexCleanroom Gloves 10	¥47,600	¥51,600	
32888-206	AK Natural Latex Cleanroom Gloves	¥291,300	¥315,400	
32888-208	AK Natural Latex Cleanroom Gloves	¥437,700	¥487,200	
32888-210	AK NaturalLatexCleanroom Gloves 10	¥56,800	¥61,500	
32888-244	耐薬品性グローブ ノースーフイ(ニトリル製)	¥291,800	¥130,000	
32888-246	耐薬品性グローブ ノースーフイ(ニトリル製)	¥283,300	¥313,700	
32888-248	耐薬品性グローブ ノースーフイ(ニトリル製)	¥157,100	¥78,100	
32888-272	耐薬品性グローブ ノースーフイ(ニトリル製)	¥190,700	¥94,300	
32888-274	耐薬品性グローブ ノースーフイ(ニトリル製)	¥188,100	¥184,200	
32888-278	耐薬品性グローブ ノースーフイ(ニトリル製)	¥190,700	¥94,300	
32888-282	耐薬品性グローブ ノースーフイ(ニトリル製)	¥205,100	¥200,900	
32888-284	耐薬品性グローブ ノースーフイ(ニトリル製)	¥390,000	¥432,000	
32888-286	耐薬品性グローブ ノースーフイ(ニトリル製)	¥205,100	¥200,900	
32888-306	耐薬品性グローブ ノースーフイ(テック製)	¥45,500	¥50,400	
32888-350	耐薬品性グローブ ホック用グローブ SZ9.75	¥168,800	¥203,300	
32888-355	Glove Box Gloves(ハイブリッド)	¥162,200	¥177,300	
32889-013	Glove Mid Temp. 12インチ Woman	¥426,800	¥447,900	
32889-289	品質検査グローブ メジック (木綿レス)	¥37,700	¥40,300	
32889-400	軽作業用グローブ メジック (コットンレス)	¥15,000	¥16,000	
32889-989	耐薬品性グローブ スーパー XLサイズ	¥5,600	¥6,000	
32890-052	耐薬品性グローブ サイズ 7-7.5 12双	¥68,800	¥73,700	
32890-054	耐薬品性グローブ サイズ 8-8.5 12双	¥68,800	¥73,700	
32890-056	耐薬品性グローブ サイズ 9-9.5 12双	¥68,800	¥73,700	
32890-058	耐薬品性グローブ サイズ 10 12双	¥68,800	¥73,700	
32890-060	耐薬品性グローブ サイズ 11 12双	¥68,800	¥73,700	
32890-118	SOL-VEXニトリルグローブ SZ8-8.5 12双	¥104,400	¥111,700	
32890-120	SOL-VEXニトリルグローブ SZ9-9.5 12双	¥104,400	¥111,700	
32890-122	SOL-VEXニトリルグローブ SZ10 12双	¥104,400	¥111,700	
32890-124	SOL-VEXニトリルグローブ SZ11 12双	¥104,400	¥111,700	
32890-341	耐薬品性グローブ スーパー サイズ	¥5,600	¥6,000	
32890-343	耐薬品性グローブ スーパー Mサイズ	¥5,600	¥6,000	
32890-345	耐薬品性グローブ スーパー Lサイズ	¥5,500	¥5,900	
32890-594	耐薬品性グローブ ケミプロ サイズ 8 12双	¥31,300	¥33,600	
32890-596	耐薬品性グローブ ケミプロ サイズ 9 12双	¥31,400	¥33,700	
32891-470	BCRグローブ (ハーフフィンガータイプ) 10双	¥17,200	¥17,900	
32891-840	ニトリルグローブ ハンズマニファクチャリング 50	¥16,600	¥17,500	
32891-841	ニトリルグローブ ハンズマニファクチャリング 50	¥16,600	¥17,500	
32891-842	ニトリルグローブ ハンズマニファクチャリング 50	¥16,600	¥17,500	
32891-844	ニトリルグローブ ハンズマニファクチャリング 50	¥16,600	¥17,500	
32892-027	耐薬品性グローブ ニトリル (ニトリル製)	¥190,600	¥204,000	
32892-028	耐薬品性グローブ ニトリル (ニトリル製)	¥190,600	¥204,000	
32892-029	耐薬品性グローブ ニトリル (ニトリル製)	¥190,600	¥204,000	
32892-031	耐薬品性グローブ ニトリル (ニトリル製)	¥190,600	¥204,000	
32892-032	耐薬品性グローブ ニトリル (ニトリル製)	¥209,700	¥224,500	
32892-033	耐薬品性グローブ ニトリル (ニトリル製)	¥209,700	¥224,500	
32892-034	耐薬品性グローブ ニトリル (ニトリル製)	¥209,700	¥224,500	
32892-035	耐薬品性グローブ ニトリル (ニトリル製)	¥209,700	¥224,500	
32892-036	耐薬品性グローブ ニトリル (ニトリル製)	¥209,700	¥224,500	
32892-043	耐薬品性グローブ ニトリル (ニトリル製)	¥209,700	¥224,500	
32892-044	耐薬品性グローブ ニトリル (ニトリル製)	¥209,700	¥224,500	
32892-045	耐薬品性グローブ ニトリル (ニトリル製)	¥209,700	¥224,500	
32892-046	耐薬品性グローブ ニトリル (ニトリル製)	¥209,700	¥224,500	
32892-047	耐薬品性グローブ ニトリル (ニトリル製)	¥209,700	¥224,500	
32892-048	耐薬品性グローブ (ハイブリッド製/裏地無し)	¥193,100	¥209,900	
32892-049	耐薬品性グローブ (ハイブリッド製/裏地無し)	¥193,100	¥209,900	
32892-050	耐薬品性グローブ (フッ素樹脂製)	¥59,500	¥76,000	
32892-051	耐薬品性グローブ (フッ素樹脂製)	¥59,500	¥76,000	
32892-052	耐薬品性グローブ (フッ素樹脂製)	¥59,500	¥76,000	
32892-053	耐薬品性グローブ クロムブラック(耐薬品製)	¥44,500	¥47,700	
32892-054	耐薬品性グローブ クロムブラック(耐薬品製)	¥44,500	¥47,700	
32892-055	耐薬品性グローブ クロムブラック(耐薬品製)	¥44,500	¥47,700	
32892-056	耐薬品性グローブ クロムブラック(耐薬品製)	¥44,500	¥47,700	
32892-057	耐薬品性グローブ ケムスター 耐薬品/天然ゴム	¥30,000	¥32,200	
32892-058	耐薬品性グローブ ケムスター 耐薬品/天然ゴム	¥30,000	¥32,200	
32892-059	耐薬品性グローブ ケムスター 耐薬品/天然ゴム	¥30,000	¥32,200	
32892-060	耐薬品性グローブ ケムスター 耐薬品/天然ゴム	¥30,000	¥32,200	
32892-076	耐薬品性グローブ ハンズマニ (テック製)	¥25,400	¥27,200	
32892-077	耐薬品性グローブ ハンズマニ (テック製)	¥25,400	¥27,200	
32892-078	耐薬品性グローブ ハンズマニ (テック製)	¥25,400	¥27,200	
32916-100	耐酸性グローブ トリチック サイズ 6 12	¥27,800	¥30,000	
32916-102	耐酸性グローブ トリチック サイズ 7 12	¥28,200	¥30,500	
32916-104	耐酸性グローブ トリチック サイズ 8 12	¥27,800	¥30,000	
32916-108	耐酸性グローブ トリチック サイズ 10 12	¥27,800	¥30,000	
32916-229	耐酸性グローブ トリチック 12	¥32,700	¥38,700	
32916-234	耐酸性グローブ トリチック 12	¥32,700	¥35,200	
32916-235	耐酸性グローブ トリチック 12	¥32,700	¥38,700	
32916-237	耐酸性グローブ トリチック 12	¥32,700	¥38,700	
32916-240	耐酸性グローブ トリチック サイズ 6 12双	¥33,800	¥36,500	
32916-241	耐酸性グローブ トリチック 12	¥32,700	¥35,200	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
32916-250	耐酸性グローブトリック 12入	¥31,700	¥34,200	
32916-271	耐酸性グローブトリック サイズ8 12双	¥31,700	¥37,700	
32916-293	耐酸性グローブトリック サイズ9 12双	¥31,700	¥37,700	
32916-318	耐酸性グローブトリック サイズ10 12双	¥31,700	¥34,200	
32916-330	耐酸性グローブトリック サイズ11 12双	¥31,700	¥34,200	
32916-350	GLOVE TRIPOLY 521 サイズ8	¥55,500	¥60,300	
32916-804	Derma Free ビニールグローブ M 100	¥7,700	¥8,300	
32917-203	耐薬品グローブワートン 衬あり/天然ゴム	¥21,300	¥22,800	
32917-206	耐薬品グローブワートン 衬あり/天然ゴム	¥21,300	¥22,800	
32917-209	耐薬品グローブワートン 衬あり/天然ゴム	¥21,300	¥22,800	
32917-212	耐薬品グローブワートン 衬あり/天然ゴム	¥21,300	¥22,800	
32917-220	耐薬品性グローブテクニク	¥36,500	¥39,100	
32917-223	耐薬品性グローブテクニク	¥36,500	¥39,100	
32917-224	耐薬品性グローブテクニク	¥36,500	¥39,100	
32917-229	耐薬品性グローブテクニク	¥42,900	¥45,900	
32917-230	耐薬品性グローブテクニク	¥42,900	¥45,900	
32917-231	耐薬品性グローブテクニク	¥42,900	¥45,900	
32917-232	耐薬品性グローブテクニク	¥42,900	¥45,900	
32917-233	耐薬品性グローブテクニク	¥43,000	¥46,000	
32917-240	耐薬品性グローブテクニク	¥65,400	¥70,000	
32917-241	耐薬品性グローブテクニク	¥65,400	¥70,000	
32917-242	耐薬品性グローブテクニク	¥65,400	¥70,000	
32917-707	StanSolv Flock-Lined NitrileGloves	¥125,500	¥74,800	
32917-718	StanSolv Flock-Lined NitrileGloves	¥125,500	¥134,400	
32917-729	StanSolv Flock-Lined NitrileGloves	¥125,500	¥134,400	
32917-730	StanSolv Flock-Lined NitrileGloves	¥125,500	¥134,400	
32917-740	StanSolv Flock-Lined NitrileGloves	¥131,100	¥140,300	
32917-810	StanSolv Unlined Nitrile Gloves	¥23,100	¥24,700	
32917-820	StanSolv Unlined Nitrile Gloves	¥23,100	¥24,700	
32917-831	StanSolv Unlined Nitrile Gloves	¥24,300	¥26,000	
32917-842	StanSolv Unlined Nitrile Gloves	¥23,100	¥24,700	
32917-853	StanSolv Unlined Nitrile Gloves	¥25,700	¥27,500	
32917-855	耐薬品性グローブスタンソルブ(ニトリル製)	¥195,500	¥208,400	
32917-856	耐薬品性グローブスタンソルブ(ニトリル製)	¥195,500	¥208,400	
32917-857	耐薬品性グローブスタンソルブ(ニトリル製)	¥195,500	¥208,400	
32917-858	耐薬品性グローブスタンソルブ(ニトリル製)	¥195,500	¥208,400	
32917-859	耐薬品性グローブスタンソルブ(ニトリル製)	¥195,500	¥208,400	
32917-860	耐薬品性グローブスタンソルブ(ニトリル製)	¥195,500	¥208,400	
32920-012	Temp-Tec Insulated Neoprene Gloves	¥23,000	¥24,700	
32920-017	耐熱性グローブ(Temp-Tec)サイズ10	¥23,400	¥25,000	
32920-020	グローブ スタンソルブ(衬あり/ニトリル/ニット裏地)	¥61,600	¥66,200	
32920-042	グローブ スタンソルブ(衬あり/ニトリル/ニット裏地)	¥61,600	¥66,200	
32920-064	グローブ スタンソルブ(衬あり/ニトリル/ニット裏地)	¥61,600	¥66,200	
32920-086	グローブ スタンソルブ(衬あり/ニトリル/ニット裏地)	¥61,600	¥66,200	
32920-270	グローブ スタンソルブ(ニトリル/ニット裏地付)	¥439,900	¥470,900	
32920-272	グローブ スタンソルブ(ニトリル/ニット裏地付)	¥439,300	¥470,200	
32920-274	グローブ スタンソルブ(ニトリル/ニット裏地付)	¥439,300	¥470,200	
32920-276	グローブ スタンソルブ(ニトリル/ニット裏地付)	¥439,300	¥470,200	
32928-088	StanSolv Unlined Nitrile Gloves	¥25,500	¥27,300	
32928-102	StanSolv Unlined Nitrile Gloves	¥25,500	¥12,200	
32928-124	StanSolv Unlined Nitrile Gloves	¥27,100	¥29,000	
32928-146	StanSolv Unlined Nitrile Gloves	¥25,500	¥27,300	
32929-604	StanSolv Unlined Nitrile Gloves	¥70,500	¥75,400	
32929-808	StanSolv Unlined Nitrile Gloves	¥70,500	¥75,400	
32929-820	StanSolv Unlined Nitrile Gloves	¥67,200	¥40,900	
32929-900	Blue-Grip Natural Rubber Gloves	¥23,900	¥25,600	
32929-901	Blue-Grip Natural Rubber Gloves	¥23,900	¥25,600	
32929-902	Blue-Grip Natural Rubber Gloves	¥23,900	¥25,600	
32929-903	Blue-Grip Natural Rubber Gloves	¥23,900	¥25,600	
32930-054	グローブ アディティブ(ニトリル/天然ゴム)	¥116,600	¥124,900	
32930-056	グローブ アディティブ(ニトリル/天然ゴム)	¥116,600	¥124,900	
32930-058	グローブ アディティブ(ニトリル/天然ゴム) 12	¥22,100	¥23,700	
32930-060	グローブ アディティブ(ニトリル/天然ゴム)	¥116,600	¥124,900	
32930-110	StanSolv Flock-Lined NitrileGloves	¥259,300	¥157,600	
32930-111	StanSolv Flock-Lined NitrileGloves	¥259,500	¥277,800	
32930-112	StanSolv Flock-Lined NitrileGloves	¥260,200	¥278,600	
32930-788	GLOVE TRIPOLY 514+ サイズ6	¥49,500	¥53,000	
32930-791	GLOVE TRIPOLY 514+ サイズ7	¥49,500	¥53,000	
32930-792	GLOVE TRIPOLY 514+ サイズ8	¥49,500	¥53,000	
32931-490	グローブ スタンソルブ(衬あり/ニトリル/ニット裏地)	¥88,200	¥48,200	
32931-515	グローブ スタンソルブ(衬あり/ニトリル/ニット裏地)	¥88,200	¥94,900	
32931-537	グローブ スタンソルブ(衬あり/ニトリル/ニット裏地)	¥88,200	¥94,900	
32932-270	品質検査用グローブ(収縮性ナイロン製)	¥26,800	¥28,700	
32932-272	品質検査用グローブ(収縮性ナイロン製)	¥25,500	¥27,300	
32935-392	インナーグローブ ウェルスマットCut-Te	¥85,700	¥89,700	
32935-394	インナーグローブ ウェルスマットCut-Te	¥99,400	¥103,700	
32935-396	インナーグローブ ウェルスマットCut-Te	¥113,500	¥118,600	
32935-398	インナーグローブ ウェルスマットCut-Te	¥117,700	¥122,900	
32935-448	Cowhide Drivers Gloves サイズ	¥9,200	¥9,600	
32935-449	Cowhide Drivers Gloves	¥9,200	¥9,600	
32935-450	Cowhide Drivers Gloves	¥9,200	¥9,600	
32935-460	品質検査用グローブ(木綿レース製)	¥7,200	¥7,600	
32935-461	品質検査用グローブ(木綿レース製)	¥7,000	¥7,400	
32935-481	インナーグローブ ウェルスマットスパータ	¥134,100	¥143,500	
32935-483	インナーグローブ ウェルスマットスパータ	¥134,100	¥143,500	
32935-485	インナーグローブ ウェルスマットスパータ	¥134,100	¥143,500	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
32935-487	インナーグローブ ウェルスメント スパータ	¥134,100	¥143,500	
32935-509	インナーグローブ ウェルスメント スパータ	¥125,300	¥134,200	
32935-581	インナーグローブ ウェルスメント(ナイロン製) S	¥29,900	¥32,000	
32935-583	インナーグローブ ウェルスメント(ナイロン製) M	¥29,900	¥32,000	
32935-585	インナーグローブ ウェルスメント(ナイロン製) L	¥29,900	¥32,000	
32935-587	インナーグローブ ウェルスメント(ナイロン製) XL	¥29,900	¥32,000	
32940-073	Glove Box Glove 15ミル アクリル製 1双	¥293,800	¥317,600	
33000-159	消防用ヘルメット	¥153,000	¥172,100	
33000-506	保護メガネ astro OTG3001(UV保護用) 1	¥13,900	¥15,000	
33000-508	保護メガネ astro OTG 3001(UV保護用)	¥109,600	¥117,900	
33000-511	保護メガネ astro OTG3001(UV保護用) 1	¥15,100	¥16,300	
33000-518	保護メガネ astro OTG 3001(UV保護用)	¥103,100	¥110,900	
33000-538	保護メガネ astro OTG 3001(UV保護用)	¥183,600	¥197,700	
33001-580	耳栓 E-A-R1	¥53,500	¥57,300	
33002-036	保護メガネ astrospec3000(UV保護用) 1	¥11,000	¥11,800	
33002-049	保護メガネ astrospec 3000(UV保護用)	¥83,700	¥89,900	
33002-067	保護メガネ(UV保護用)	¥22,400	¥24,600	
33002-078	保護メガネ Blak-Ray1(UV保護用)	¥22,400	¥24,000	
33002-210	保護メガネ astrospec 3000(UV保護用)	¥9,700	¥10,400	
33002-332	保護メガネ astro OTG 3001(UV保護用)	¥15,700	¥16,900	
33012-000	保護メガネ Visorgog グリーン	¥6,400	¥7,000	
33012-005	保護メガネ Visorgog グリーン	¥6,400	¥7,000	
33012-006	保護メガネ Visorgog ブルー	¥6,400	¥7,000	
33012-007	保護メガネ Visorgog ホワイト	¥7,000	¥7,600	
33012-008	保護メガネ Visorgog レッド	¥6,400	¥7,000	
33012-010	保護メガネ用 Replacement Lens	¥5,900	¥6,500	
33500-238	Adapter Thermocouple 24/40JT	¥53,100	¥57,700	
33501-430	Stopper Polyethylene 14/20 100個	¥10,500	¥11,400	
33501-432	Stopper Polyethylene 19/22 100個	¥12,100	¥13,200	
33501-434	Stopper Polyethylene 24/25 100個	¥13,500	¥14,700	
33501-438	Stopper Polyethylene 29/42 100個	¥14,600	¥15,800	
33501-674	Tru-Stir Stirrer Shaft Coupling	¥212,800	¥230,000	
33501-844	FUNNEL SYSTEM HDPE 41N 2L BTL	¥252,700	¥278,500	
33501-892	Funnels FLASH Chromatography 60ml	¥30,400	¥44,000	
33501-894	Funnels FLASH Chromatography 120ml	¥25,900	¥37,500	
33501-896	Funnels FLASH Chromatography 240ml	¥25,900	¥37,500	
33501-898	Funnels FLASH Chromatography 960ml	¥26,200	¥37,900	
33501-900	Funnels FLASH Chromatography 960ml	¥26,100	¥37,800	
33503-124	StainlessSteel用バケツ 19L	¥272,100	¥330,600	
33674-000	硬度計	¥807,100	¥1,567,000	
33676-005	硬度計	¥804,100	¥1,470,000	
33729-065	Heating Tape(幅12.7mm/長さ2.4m)	¥166,300	¥184,300	
33742-206	カット用カッター	¥820,600	¥895,800	
33750-265	Spillproof Mantle Emx Series	¥634,300	¥1,188,000	
33994-118	StirringElementsノドール 3枚羽 PTFE	¥589,900	¥1,133,000	
33994-130	StirringElementsノドール 8/1.5-6/35	¥175,700	¥264,000	
33994-132	StirringElementsノドール 8/2.4-10/55	¥235,500	¥357,000	
33994-136	StirringElementsノドール 8/7/55/800	¥196,500	¥375,000	
34104-737	Electrode Cable 15FT	¥179,400	¥192,700	
34105-859	VWR pH Electrode Arm	¥95,900	¥106,800	
34559-903	比重計	¥19,900	¥23,800	
34559-925	比重計	¥19,900	¥21,600	
34559-947	比重計	¥19,800	¥21,500	
34559-969	比重計	¥20,000	¥21,700	
34559-980	比重計	¥20,000	¥21,700	
34560-006	比重計	¥20,500	¥22,200	
34565-000	比重計	¥20,000	¥21,700	
34570-000	比重計	¥12,700	¥13,800	
34575-003	比重計	¥13,200	¥14,400	
34580-003	比重計	¥12,900	¥14,000	
34585-007	比重計	¥12,800	¥13,900	
34610-002	比重計	¥14,100	¥15,600	
34615-006	比重計	¥14,200	¥15,600	
34620-006	比重計	¥14,300	¥15,700	
34620-162	比重計	¥13,100	¥14,500	
34622-001	比重計	¥14,000	¥15,400	
34622-158	比重計	¥14,000	¥15,400	
34622-205	比重計	¥14,000	¥15,400	
34622-250	比重計	¥13,900	¥15,100	
34623-106	比重計	¥13,800	¥15,000	
34623-150	比重計	¥13,900	¥15,300	
34623-208	比重計	¥13,900	¥15,100	
34623-252	比重計	¥13,700	¥15,100	
34623-300	比重計	¥13,700	¥15,100	
34627-027	比重計	¥17,900	¥19,400	
34627-050	比重計	¥17,900	¥19,700	
34627-082	比重計	¥16,900	¥18,300	
34627-118	比重計	¥16,900	¥18,600	
34627-130	比重計	¥16,900	¥18,600	
34627-173	比重計	¥16,900	¥18,600	
34627-174	比重計	¥16,900	¥18,600	
34627-231	比重計	¥16,900	¥18,600	
34627-253	比重計	¥16,900	¥18,600	
34627-297	比重計	¥16,900	¥18,600	
34627-311	比重計	¥16,900	¥18,600	
34627-355	比重計	¥16,900	¥18,600	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
34627-388	比重計	¥16,800	¥18,500	
34627-413	比重計	¥16,800	¥18,500	
34627-446	比重計	¥16,700	¥18,400	
34627-479	比重計	¥16,900	¥18,600	
34627-504	比重計	¥17,900	¥19,400	
34627-537	比重計	¥17,500	¥19,000	
34627-560	比重計	¥17,000	¥18,700	
34627-628	比重計	¥17,900	¥19,700	
34627-650	比重計	¥17,500	¥19,300	
34630-000	比重計	¥165,900	¥180,300	
34640-207	比重計	¥13,800	¥15,000	
34640-262	比重計	¥13,900	¥15,300	
34640-284	比重計	¥14,600	¥16,200	
34640-309	比重計	¥13,900	¥15,000	
34640-320	比重計	¥13,600	¥14,800	
34640-342	比重計	¥14,300	¥15,800	
34670-106	比重計	¥13,800	¥15,200	
34670-128	比重計	¥14,500	¥16,000	
34670-140	比重計	¥14,100	¥15,500	
34670-161	比重計	¥13,900	¥15,000	
34780-002	比重計	¥23,000	¥24,400	
34785-006	比重計	¥17,300	¥18,300	
34787-089	比重計	¥17,200	¥19,000	
34787-114	比重計	¥17,000	¥18,700	
34787-147	比重計	¥17,400	¥19,100	
34787-170	比重計	¥16,900	¥18,600	
34787-205	比重計	¥14,200	¥15,700	
34787-238	比重計	¥15,200	¥16,700	
34787-260	比重計	¥15,100	¥16,700	
34787-329	比重計	¥17,400	¥19,200	
34808-021	比重計	¥15,800	¥17,200	
34808-101	比重計	¥17,700	¥19,200	
34808-156	比重計	¥18,700	¥20,400	
34808-214	比重計	¥18,700	¥20,400	
34808-269	比重計	¥17,600	¥19,100	
34808-327	比重計	¥17,600	¥20,300	
34808-382	比重計	¥18,100	¥19,600	
34808-418	比重計	¥18,200	¥21,900	
34808-462	比重計	¥17,400	¥18,900	
35100-100	Sample Vials 144	¥112,900	¥125,100	
35100-104	Vial Borosilicate 3.75ml 144	¥157,000	¥173,800	
35123-324	比重計	¥19,000	¥20,600	
35123-326	比重計	¥149,900	¥181,000	
35123-346	比重計	¥20,500	¥22,300	
35123-348	比重計	¥149,900	¥181,000	
35123-368	比重計	¥19,000	¥20,600	
35123-370	比重計	¥146,700	¥177,100	
35123-380	比重計	¥18,700	¥20,300	
35123-382	比重計	¥145,500	¥154,200	
35123-404	比重計	¥18,800	¥20,400	
35123-406	比重計	¥149,900	¥181,000	
35123-426	比重計	¥18,800	¥20,400	
35123-428	比重計	¥177,900	¥205,400	
35123-450	比重計	¥173,600	¥188,600	
35123-460	比重計	¥18,300	¥19,900	
35123-462	比重計	¥151,600	¥183,000	
35125-649	比重計	¥162,500	¥196,200	
35125-671	比重計	¥159,900	¥193,000	
35125-682	比重計	¥162,100	¥195,600	
35125-707	比重計	¥190,500	¥206,900	
35125-729	比重計	¥161,200	¥194,600	
35125-751	比重計	¥148,700	¥179,500	
35125-762	比重計	¥159,100	¥192,000	
35125-784	比重計	¥159,100	¥192,000	
35125-809	比重計	¥159,100	¥192,000	
35400-100	ウォーター入(振盪機能付)	¥37,200	¥68,000	
35790-084	Rocker II Platform Mixer Boekel	¥35,300	¥64,000	
35915-116	リソクター 96ウェルプレート用	¥267,000	¥486,000	
36060T	100ul マルチフィルターチップ LB S96x10R	¥18,400	¥19,800	
36312-004	保存ビン(ステンレス製) 1.2L	¥74,900	¥89,600	
36312-005	保存ビン(ステンレス製) 1.8L	¥81,300	¥95,500	
36312-006	保存ビン(ステンレス製) 2.8L	¥76,000	¥95,300	
36312-007	保存ビン(ステンレス製) 4.5L	¥85,300	¥102,100	
36318-852	保存ビン(ポリプロピレン製)	¥79,400	¥85,400	
36318-854	保存ビン(ポリプロピレン製)	¥95,100	¥102,300	
36318-856	保存ビン(ポリプロピレン製)960ml 24本入	¥189,000	¥203,300	
36318-858	保存ビン(ポリプロピレン製)	¥105,000	¥113,000	
36318-860	保存ビン(ポリプロピレン製)	¥132,500	¥142,500	
36318-862	保存ビン(ポリスチレン製)	¥220,400	¥237,200	
36318-864	保存ビン(ポリスチレン製)	¥127,500	¥137,200	
36318-866	保存ビン(ポリスチレン製)	¥141,700	¥152,500	
36318-868	保存ビン(ポリスチレン製)	¥93,400	¥100,500	
36318-870	保存ビン(ポリスチレン製)	¥99,500	¥107,100	
36319-268	広口茶瓶 30mL	¥70,600	¥76,700	
36319-343	ポリ瓶(キャップ付)広口丸瓶	¥56,300	¥59,100	
36319-728	広口褐色ガラス瓶(キャップ付)30ml	¥39,200	¥42,600	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
36319-729	広口褐色ガラス瓶(キャップ付)60ml	¥60,300	¥65,500	
36319-730	広口褐色ガラス瓶(キャップ付)120ml	¥71,100	¥77,200	
36319-731	広口褐色ガラス瓶(キャップ付)250ml	¥97,000	¥105,400	
36319-732	広口褐色ガラス瓶(キャップ付)500ml	¥54,600	¥59,300	
36319-733	広口褐色ガラス瓶(キャップ付)950ml	¥74,100	¥80,500	
36319-774	広口褐色ガラス瓶(キャップ付)30ml	¥58,800	¥63,900	
36319-775	広口褐色ガラス瓶(キャップ付)60ml	¥74,300	¥80,700	
36319-776	広口褐色ガラス瓶(キャップ付)120ml	¥83,200	¥90,400	
36319-777	広口褐色ガラス瓶(キャップ付)250ml	¥126,200	¥137,100	
36319-778	広口褐色ガラス瓶(キャップ付)500ml	¥75,700	¥82,200	
36319-779	広口褐色ガラス瓶(キャップ付)950ml	¥51,700	¥56,100	
36324-000	保存ビン(テフロン製)	¥68,100	¥76,100	
36324-001	保存ビン(テフロン製)	¥116,500	¥130,300	
36324-002	保存ビン(テフロン製)	¥121,700	¥136,100	
36324-003	保存ビン(テフロン製)	¥222,500	¥239,400	
36324-005	保存ビン(テフロン製)	¥453,600	¥507,500	
36324-006	保存ビン(テフロン製)	¥608,100	¥680,300	
36419-213	メタルタグ(アルミニウム製)79x35mm 48個入	¥34,600	¥38,500	
36433-435	ラベル(危険警告用)	¥32,100	¥44,600	
36433-437	ラベル(危険警告用)	¥32,500	¥41,700	
36434-300	バリケートテープ	¥8,900	¥9,600	
36434-301	バリケートテープ	¥9,100	¥9,800	
36434-302	バリケートテープ	¥19,900	¥21,400	
36434-304	バリケートテープ	¥8,900	¥9,600	
36434-305	バリケートテープ	¥17,300	¥18,500	
36434-307	バリケートテープ	¥30,200	¥32,300	
36434-308	バリケートテープ	¥115,300	¥123,400	
36434-321	バリケートテープ	¥34,500	¥38,600	
36434-327	バリケートテープ	¥30,600	¥32,800	
36434-330	バリケートテープ	¥34,500	¥36,900	
36434-332	バリケートテープ	¥31,300	¥33,500	
36434-334	バリケートテープ	¥35,300	¥37,800	
36436-507	ラベル(ラミネート用/低温耐性)	¥100,400	¥116,400	
36436-530	ラベル(ラミネート用/低温耐性)	¥101,300	¥117,500	
36436-562	ラベル(ラミネート用/低温耐性)	¥105,400	¥122,300	
36440-107	ラベリングテープ (Tape Dispenser)	¥59,000	¥63,100	
36440-446	ラベリングテープ (Labels) 1000入	¥23,500	¥25,400	
36440-640	ラベリングテープ (Labels) 1000入	¥19,800	¥21,100	
36440-843	ラベリングテープ (Labels) 1000入	¥20,700	¥22,100	
36442-248	ラベリングテープ (Tape Roll)50.8mm wide	¥25,300	¥27,700	
36443-003	ラベリングテープ ティスランサー	¥66,300	¥71,000	
36490-003	踏み台(キャスター/手すり付/3段)	¥302,100	¥317,900	
36490-060	踏み台(キャスター付)	¥361,600	¥380,500	
36490-062	踏み台(キャスター付)	¥682,100	¥1,207,000	
36547-009	ラジエーター(赤外線ラジエーター、加熱・乾燥用)	¥188,000	¥353,000	
36548-001	ラジエーター(赤外線ラジエーター、加熱・乾燥用)	¥18,600	¥35,000	
36926-935	Folding Pocket Magnifiers B & L	¥19,700	¥21,100	
36927-000	Measuring Magnifier B & L	¥83,800	¥89,700	
36929-005	Measuring Magnifier B & L	¥45,400	¥48,600	
36929-027	Measuring Magnifier B & L	¥45,400	¥48,600	
36929-049	Measuring Magnifier B & L	¥45,400	¥48,600	
36929-060	Measuring Magnifier B & L	¥45,500	¥48,700	
36933-002	Illuminating Coddington Magnifier	¥33,900	¥61,000	
36934-136	VWR電気伝導率校正用標準液	¥112,100	¥123,600	
36934-160	Platinum Freezer Thermometer	¥58,900	¥110,000	
36937-003	Folding Pocket Magnifiers B & L	¥27,600	¥29,600	
36937-069	ルンバ Coddington Magnifiers B & L	¥38,300	¥41,000	
36955-027	ルンバ Hastings Triplet Magnifiers B&L	¥53,400	¥57,200	
36968-007	ルンバ ハンディ型	¥90,200	¥96,500	
36982-008	KIMAX キャビネット型 100	¥16,400	¥17,800	
36983-000	KIMAX キャビネット型 100	¥12,200	¥13,300	
36984-003	KIMAX キャビネット型 100入	¥14,600	¥15,900	
36993-300	水銀回収容器	¥25,100	¥26,900	
37000-922	Basic Plus Lab Coats 25	¥181,900	¥194,800	
37001-074	靴カバー/フット	¥106,800	¥112,900	
37001-272	PCRプレート96ウェル(PP製)	¥31,500	¥32,100	
37001-286	排液回収タンク 7.6L 1	¥811,700	¥851,900	
37002-000	ワイパー(ワイパー-ルX60/インク追加工済)	¥286,700	¥317,700	
37002-004	ワイパー(ワイパー-ルX60/インク追加工済)	¥205,100	¥215,200	
37002-006	ワイパー(ワイパー-ルX60/インク追加工済)	¥194,200	¥203,900	
37002-008	ワイパー(ワイパー-ルX60/インク追加工済)	¥186,300	¥195,500	
37002-012	ワイパー(ワイパー-ル20キムタ株) 30入	¥49,500	¥66,900	
37002-016	ワイパー(ワイパー-ル/使い捨て用/スタック型)	¥133,400	¥140,000	
37002-020	ワイパー(ワイパー-ルX60/インク追加工済)	¥189,000	¥198,400	
37002-022	ワイパー(キムタ株/再使用可能タイプ)	¥109,500	¥114,900	
37002-024	ワイパー(ワイパー-ルX60/インク追加工済)	¥122,800	¥128,900	
37002-028	ワイパー(ワイパー-ルX70)	¥155,400	¥163,100	
37002-032	ワイパー(ワイパー-ル20キムタ株)	¥143,900	¥151,100	
37002-040	ワイパー(ワイパー-ルX70) 100枚x10	¥187,900	¥212,600	
37002-046	ワイパー(ワイパー-ルX60/インク追加工済)	¥138,000	¥144,900	
37002-048	ワイパー(ワイパー-ル20キムタ株)	¥138,200	¥145,000	
37002-050	ワイパー(キムタ株/再使用可能タイプ)	¥121,900	¥128,000	
37002-052	ワイパー(ワイパー-ルX70)	¥85,400	¥95,200	
37002-058	ワイパー(キムタ株/再使用可能タイプ)	¥72,400	¥76,000	
37650T	10ul マイクロチップ LB 1000/bag	¥5,400	¥5,800	
37660T	10ul マイクロチップ LB S96X10R	¥12,600	¥13,600	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
37699-382	Erem Magic Cutter	¥164,700	¥172,900	
37699-570	ドライバースセット(精密電子機器用)	¥46,200	¥59,000	
37699-716	ドライバースセット	¥52,000	¥58,200	
37777-000	照射容器(マウス) 1	¥194,200	¥220,400	
37777-002	動物用回復容器 1	¥113,500	¥125,400	
37777-004	照射容器キャップ 1	¥228,700	¥247,900	
37820T	10ul クイックラック リーチアップ LB 96X50	¥39,000	¥42,000	
38000T	10ul フィルター リーチアップ LB S96x10R	¥20,300	¥21,900	
39640T	1.7ml LBマイクロチューブ NT 250/bag	¥2,900	¥3,100	
39870-218	防振用グローブ セーフティ&ILJノミクス(革製)	¥24,900	¥26,600	
39870-236	防振用グローブ セーフティ&ILJノミクス(革製)	¥67,700	¥72,500	
39870-286	軽作業用グローブ セーフティ&ILJノミクス	¥16,200	¥17,300	
40000-270	VWRビニール回転式スクラップ (6台収納)	¥81,700	¥88,300	
40000-272	VWRビニール直列式スクラップ (6台収納)	¥46,600	¥50,400	
40000-320	チューブラック(凍結保存用・PP製) 1個	¥10,500	¥11,500	
40000-322	チューブラック(凍結保存用・PP製)	¥7,700	¥8,500	
40000-324	チューブラック(凍結保存用・PP製)	¥18,600	¥20,500	
40000-326	チューブラック(凍結保存用・PP製)	¥7,800	¥8,600	
40000-328	チューブラック(凍結保存用・PP製)	¥7,800	¥8,600	
40000-330	チューブラック(凍結保存用・PP製)	¥7,800	¥8,600	
40000-332	チューブラック(凍結保存用・PP製) 5個	¥51,100	¥70,400	
40000-340	チューブラック(凍結保存用・PP製)	¥153,400	¥169,100	
40000-342	チューブラック(凍結保存用・PP製) 1個	¥126,300	¥139,200	
40000-730	VeryHigh Intensity Microprobe210um	¥279,800	¥302,600	
40101-146	テックスグローブ サテイクリーン(クリーンルーム/ST)	¥255,700	¥273,800	
40101-148	テックスグローブ サテイクリーン(クリーンルーム/ST)	¥256,500	¥255,800	
40101-150	テックスグローブ サテイクリーン(クリーンルーム/ST)	¥255,700	¥273,800	
40101-152	テックスグローブ サテイクリーン(クリーンルーム/ST)	¥255,700	¥273,800	
40101-352	ニトリグローブ サテイクリーン(クリーンルーム用)	¥209,700	¥224,500	
40101-353	ニトリグローブ サテイクリーン(クリーンルーム用)	¥209,700	¥224,500	
40101-354	ニトリグローブ サテイクリーン(クリーンルーム用)	¥209,700	¥224,500	
40101-356	ニトリグローブ サテイクリーン(クリーンルーム用)	¥219,700	¥235,200	
40101-357	ニトリグローブ サテイクリーン(クリーンルーム用)	¥209,700	¥224,500	
40101-358	ニトリグローブ サテイクリーン(クリーンルーム用)	¥209,700	¥224,500	
41013-200	Bitran Liquid-tight Formalin Bags	¥327,600	¥357,800	
414004-002	VWR Disposable Transfer Pipets 500	¥18,700	¥19,700	
414004-004	VWR Disposable Transfer Pipets 500	¥12,400	¥13,700	
414004-005	VWR Disposable Transfer Pipets 500	¥18,000	¥20,900	
414004-014	VWRテイスルサフールトランスファーピペット 500	¥12,800	¥13,500	
414004-015	VWRテイスルサフールトランスファーピペット 500	¥21,200	¥22,400	
414004-021	VWRテイスルサフールトランスファーピペット 500本	¥37,300	¥39,300	
414004-024	VWRテイスルサフールトランスファーピペット 500	¥21,100	¥22,200	
414004-042	VWR Disposable Transfer Pipets 500	¥28,400	¥29,900	
414004-053	VWR分光光度計 0.7mL Flat Lid 1個	¥105,500	¥114,600	
414004-060	分光光度計用皿(石英)	¥305,600	¥331,900	
414004-062	分光光度計用皿(石英)	¥297,700	¥323,400	
414004-069	分光光度計用皿(石英) 2個	¥327,400	¥355,700	
414004-070	VWR分光光度計 0.7mL Flat Lid 2個	¥380,800	¥413,600	
414004-076	VWR分光光度計 0.3mL Flat Lid 1個	¥96,500	¥104,800	
414004-125	VWR Laboratory Bottles 250mL	¥29,100	¥31,300	
414004-144	VWR Graduated Beakers 30mL 100個	¥23,400	¥27,500	
414004-145	VWR Graduated Beakers 50mL 100個	¥25,600	¥27,600	
414004-189	VWR Large Square Bottle PP 3.8L 1	¥25,000	¥26,900	
414004-234	乳鉢(アルミ製) 1	¥93,900	¥102,000	
414004-265	VWR テイスルサフール吸引ピペット 2mL	¥19,400	¥20,400	
414004-266	VWR テイスルサフール吸引ピペット 5mL	¥22,600	¥23,900	
414004-267	VWR テイスルサフール吸引ピペット 10mL	¥23,800	¥25,000	
414004-279	Analytical Funnels 180mm 4	¥20,400	¥21,900	
414004-309	VWR Maximum Protection Coveralls (¥111,200	¥112,000	
414004-414	VWR Boot Covers Universal 200	¥240,700	¥242,500	
414004-498	Anti-Skid ShoeCover Universal 300	¥62,500	¥63,000	
414004-504	VWR Irradiated Bouffant Cap S 1000	¥46,300	¥46,700	
414004-505	VWR Irradiated Bouffant CapsL 1000	¥49,100	¥49,500	
414004-506	VWR Irradiated Beard Cover 1000	¥45,300	¥45,600	
414004-510	VWR Sterile Sleeves Universal	¥195,200	¥196,700	
414004-517	VWRホリステルカバー - 150枚	¥39,500	¥41,000	
414004-666	マスク Blue 50枚x10	¥67,300	¥69,900	
414004-672	Maximum Protection Face Mask 500入	¥97,400	¥101,200	
414004-680	マスク White 50枚x10	¥57,600	¥65,900	
41428-950	Reagentボトル - Vistalab 100mL ST	¥61,200	¥66,100	
41428-952	Reagentボトル - Vistalab 100mL 50入	¥66,300	¥72,100	
41428-958	Reagentボトル - Vistalab 25mL 10x10	¥59,600	¥64,400	
41446-974	Microscopy Illuminator MI-150用	¥653,100	¥1,221,000	
41447-004	Dual Gooseneckホルダー(MI-150)	¥273,000	¥511,000	
43100-064	1860 N95 防塵マスク レディース 20入	¥26,200	¥29,100	
43300-430	96-Well MicroWell Plates 50枚	¥234,800	¥249,200	
43300-942	CSLB MicroSeal Sleeve Seal Tape 1	¥25,300	¥27,400	
45000-672	Cell Culture Grade Water 6本	¥53,400	¥55,500	
45800-016	Anticon100クリーンルーム用ワイパー - 150	¥32,500	¥34,800	
45800-058	クリーンルーム用ワイパー - 150枚/PK	¥22,600	¥24,200	
45800-062	クリーンルーム用ワイパー - 150	¥60,300	¥64,600	
45800-068	Anticon クリーンルーム用ワイパー - 150	¥70,000	¥74,900	
45800-084	ワイパー - PRE MOIST 6% ISO A 50	¥37,100	¥39,700	
45800-134	クリーンルーム用ワイパー - 75枚入	¥60,300	¥64,600	
46000-096	耳栓 QB2HYG	¥49,000	¥55,100	
46100-826	丸底遠心チューブ 30ml用7段階ター 2	¥42,300	¥46,100	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
46100-830	Rubber Adapter Sleeve 15ml 2個	¥39,500	¥43,100	
46300-326	Automation Pipet Tips 960本	¥71,000	¥79,000	
46300-494	マイクロタイテリレーションチューブ 1.2ml 1000	¥53,600	¥63,300	
46300-562	ハック(再密閉可能)	¥89,800	¥96,200	
46300-564	ハック(再密閉可能)	¥93,000	¥99,500	
46300-884	ワイルド(ワイルド/使い捨てハック用)	¥120,900	¥126,300	
46600-032	活栓(コック)	¥57,700	¥64,000	
46600-036	活栓(コック)	¥46,300	¥51,300	
46600-038	活栓(コック)	¥77,300	¥85,600	
46600-040	逆止弁	¥44,700	¥49,600	
46600-050	チューブコネクター	¥21,100	¥23,900	
46600-054	チューブコネクター 5	¥100,900	¥113,900	
46600-056	チューブコネクター(ホリカミ製)	¥34,200	¥15,600	
46600-058	チューブコネクター(ホリカミ製) 20	¥87,900	¥17,800	
46600-060	チューブコネクター(ホリカミ製)	¥51,600	¥23,700	
46600-062	チューブコネクター(ホリカミ製)	¥125,000	¥31,100	
46600-064	チューブコネクター(ホリカミ製)	¥71,200	¥18,200	
46600-066	チューブコネクター(ホリカミ製)	¥78,900	¥20,500	
46600-068	チューブコネクター(ホリカミ製)	¥25,600	¥26,900	
46600-070	チューブコネクター(ホリカミ製)	¥53,800	¥34,100	
46600-072	チューブコネクター(ホリカミ製)	¥43,100	¥39,200	
46600-086	チューブコネクター(ホリカミ製)	¥28,600	¥19,200	
46600-088	チューブコネクター(ホリカミ製)	¥65,600	¥21,300	
46600-112	流量目視計 1個	¥21,900	¥24,700	
46600-114	流量目視計	¥195,500	¥220,600	
46600-122	チューブコネクター(ホリカミ製)	¥33,300	¥11,700	
46600-124	チューブコネクター(ホリカミ製)	¥43,900	¥22,300	
46600-126	チューブコネクター(ホリカミ製)	¥46,000	¥12,100	
46600-128	チューブコネクター(ホリカミ製)	¥69,200	¥13,900	
46600-138	チューブコネクター(ホリカミ製)	¥121,200	¥134,400	
46600-140	チューブコネクター(ホリカミ製)	¥185,400	¥205,300	
46600-146	チューブコネクター(ホリカミ製)	¥90,200	¥99,900	
46600-148	チューブコネクター(ホリカミ製)	¥26,100	¥19,500	
46600-160	TRIlites994Tri-Polymerカートリッジ S 100	¥21,700	¥23,200	
46600-162	TRIlites994Tri-Polymerカートリッジ M 100	¥21,700	¥23,200	
46600-164	TRIlites994Tri-Polymerカートリッジ L 100	¥21,700	¥23,200	
46600-542	トレ(ホリカミ製)	¥359,900	¥377,700	
46600-544	トレ(ホリカミ製)	¥368,400	¥386,600	
46600-546	トレ(ホリカミ製)	¥332,200	¥348,700	
46600-548	トレ(ホリカミ製)	¥285,900	¥304,600	
46600-558	トレ(ホリカミ製)	¥470,000	¥500,700	
46600-678	I-Chem 40ml褐色セロファン 144本	¥171,100	¥198,000	
46610-040	VWR電気伝導率校正用標準液	¥101,600	¥113,100	
46610-290	循環恒温槽(糞便性大腸菌検査用)	¥336,800	¥613,000	
46610-306	防塵マスク 2600 N95	¥21,300	¥22,800	
46610-644	ハック(再密閉可能)	¥63,000	¥67,400	
46610-702	オートサフラー用ハック(1.8mL)	¥55,100	¥59,300	
46620-002	ハック(再密閉可能)	¥37,100	¥39,700	
46620-006	ハック(再密閉可能)	¥49,600	¥53,100	
46620-008	ハック(再密閉可能)	¥37,700	¥40,400	
46620-012	ハック(再密閉可能)	¥118,200	¥126,600	
46620-014	遮光ハック(再密閉可能)	¥26,700	¥28,600	
46620-016	遮光ハック(再密閉可能)	¥33,100	¥35,400	
46620-025	ハック(再密閉可能)10.2x15.2cm 100	¥2,300	¥2,400	
46620-026	ハック(再密閉可能)	¥44,500	¥47,700	
46620-032	ハック(再密閉可能)	¥64,100	¥68,600	
46620-076	ハック(再密閉可能)	¥29,100	¥31,200	
46620-078	ハック(再密閉可能)	¥46,300	¥49,600	
46620-080	ハック(再密閉可能)	¥74,600	¥79,800	
46620-234	実験室用 ムック	¥27,800	¥30,100	
46620-238	実験室用 ムック	¥32,100	¥34,700	
46620-240	実験室用 ムック	¥36,900	¥39,900	
46620-242	実験室用 ムック	¥43,700	¥47,300	
470002-082	VWR Polystyrene Spot Plate 12個入	¥24,500	¥25,700	
470021-368	Zebra Fish Breeding Tank 1	¥30,700	¥30,800	
470104-462	Rubber Policeman 12	¥13,700	¥14,400	
470104-630	Kalliroscope Rheoscope Fluid	¥17,900	¥19,400	
470122-234	SargentWelch Filter Pump/Aspirator	¥71,700	¥76,900	
470123-962	Wash Bottle Carrier 1	¥45,200	¥49,300	
470144-262	Rainbow Laboratory Tape	¥30,500	¥33,300	
470145-118	Disposable Pipets 9inch long, ST	¥116,800	¥128,100	
470145-138	Bel-Bulb Pipettor 1	¥19,600	¥21,200	
470145-298	Pipet Pumps ノール 1&2ML	¥30,100	¥32,600	
470145-876	Glass coverslips 116枚入	¥4,100	¥4,500	
470149-424	ガラスカバー 4x6インチ 50枚入	¥1,700	¥2,000	
470200-930	Glass Depression Slides	¥20,800	¥21,800	
470200-932	Glass Depression Slides	¥24,200	¥26,500	
470225-048	Polyethylene Transfer Pipettes 500	¥29,500	¥30,900	
47727-024	分光光度計(UVビームスケール)	¥77,900	¥142,000	
47727-232	流量モニター 6.4-8.0mm 1	¥19,100	¥20,400	
47727-234	流量モニター 8.0-12.8mm 1	¥19,000	¥20,300	
47727-238	流量モニター 4.8-6.4mm 1	¥19,000	¥20,300	
47727-256	加二重シリンダー - Large 50入	¥88,400	¥94,700	
47727-266	ガラスムック	¥138,300	¥151,700	
47727-268	カートリッジ & 専用容器システム(PS製/滅菌済み)	¥190,000	¥208,400	
47727-270	ガラスムック	¥62,000	¥68,000	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
47727-664	ポットシート 柄付 幅50.8cmx長さ7.3m	¥230,800	¥257,300	
47727-666	ガラススプーン 1mL 18cm 12個	¥19,600	¥20,700	
47727-682	Sterileware Large Sampling Spoon25	¥50,800	¥54,900	
47727-700	ガラススプーン 2.46mL 12個入	¥19,100	¥20,200	
47727-704	ガラススプーン 14.79mL 12個入	¥20,000	¥21,800	
47727-708	万能スプーン (スチール) ポリプロピレン製	¥30,400	¥32,300	
47727-710	万能スプーン (スチール) ポリプロピレン製	¥29,300	¥31,800	
47727-724	葉さし	¥44,500	¥47,200	
47727-822	循環恒温槽(投入型)室温+5℃-99℃	¥159,800	¥300,000	
47727-848	ウォーターバス(プラスチック製)	¥62,000	¥116,000	
47728-756	冷却循環恒温槽 20℃-40℃	¥2,042,400	¥3,824,000	
47729-350	冷却循環恒温槽 20℃-40℃	¥210,200	¥394,000	
47729-566	加液チューブ(ガラス製)	¥141,800	¥161,700	
47729-568	加液チューブ リム無 10x75mm 3ml 1000	¥34,000	¥36,900	
47729-570	加液チューブ リム無 12x75mm 5ml 1000	¥38,800	¥41,300	
47729-572	加液チューブ リム無 13x100mm 9ml 1000	¥46,300	¥50,200	
47729-574	加液チューブ リム無 15x85mm 10ml 1000	¥69,000	¥74,900	
47729-576	加液チューブ リム無 16x100mm 14ml 1000	¥68,000	¥73,800	
47729-578	加液チューブ リム無 16x125mm 18ml 1000	¥75,800	¥82,300	
47729-580	加液チューブ リム無 16x150mm 22ml 1000	¥81,300	¥88,400	
47729-583	加液チューブ リム無 18x150mm 28ml 500	¥60,200	¥65,400	
47729-584	加液チューブ リム無 20x150mm 34ml 500	¥101,700	¥110,400	
47729-586	加液チューブ リム無 25x150mm 60ml 500	¥113,900	¥123,700	
47729-740	オートクレーブ(高圧蒸気滅菌器)	¥8,300	¥15,000	
47729-968	ふるい振とう器	¥142,300	¥160,500	
47729-972	ふるい振とう器	¥140,800	¥158,900	
47729-976	ふるい振とう器	¥140,800	¥158,900	
47729-978	ふるい振とう器	¥140,800	¥158,900	
47729-980	ふるい振とう器	¥140,800	¥158,900	
47729-982	ふるい振とう器	¥140,800	¥158,900	
47729-984	ふるい振とう器	¥149,900	¥169,200	
47729-986	ふるい振とう器	¥149,800	¥169,000	
47729-988	ふるい振とう器	¥166,300	¥187,700	
47729-990	ふるい振とう器	¥181,300	¥204,600	
47729-992	ふるい振とう器	¥198,400	¥223,800	
47729-994	ふるい振とう器	¥230,700	¥260,300	
47729-996	ふるい振とう器	¥284,200	¥320,600	
47729-998	ふるい振とう器	¥297,200	¥335,400	
47730-082	ふるい振とう器	¥138,800	¥156,600	
47730-084	ふるい振とう器	¥140,700	¥158,700	
47730-086	ふるい振とう器	¥142,200	¥160,500	
47730-090	ふるい振とう器	¥142,200	¥160,500	
47730-092	ふるい振とう器	¥138,800	¥156,600	
47730-094	ふるい振とう器	¥139,300	¥157,200	
47730-096	ふるい振とう器	¥138,800	¥156,600	
47730-098	ふるい振とう器	¥147,700	¥166,700	
47730-102	ふるい振とう器	¥165,500	¥186,700	
47730-104	ふるい振とう器	¥176,400	¥199,000	
47730-106	ふるい振とう器	¥198,400	¥223,800	
47730-112	ふるい振とう器	¥292,800	¥330,400	
47730-196	ふるい振とう器	¥142,200	¥160,500	
47730-198	ふるい振とう器	¥142,200	¥160,500	
47730-200	ふるい振とう器	¥138,800	¥156,600	
47730-202	ふるい振とう器	¥141,400	¥159,600	
47730-206	ふるい振とう器	¥138,800	¥156,600	
47730-208	ふるい振とう器	¥138,800	¥156,600	
47730-210	ふるい振とう器	¥142,200	¥160,500	
47730-212	ふるい振とう器	¥149,200	¥168,300	
47730-214	ふるい振とう器	¥151,200	¥170,600	
47730-216	ふるい振とう器	¥165,500	¥186,700	
47730-218	ふるい振とう器	¥178,500	¥201,500	
47730-220	ふるい振とう器	¥200,300	¥226,000	
47730-222	ふるい振とう器	¥229,500	¥259,000	
47730-224	ふるい振とう器	¥280,000	¥315,900	
47730-226	ふるい振とう器	¥294,200	¥331,900	
47730-580	ふるい振とう器	¥107,800	¥121,700	
47730-582	ふるい用黄銅製加液	¥99,000	¥111,700	
47730-730	No-Wire Indexed Test Tube Racks	¥39,500	¥42,400	
47732-446	ポットリザーブ - 15 mL Conical	¥89,300	¥98,400	
47732-448	ポットリザーブ - 35mL(Free-Standing)	¥89,300	¥98,400	
47732-450	ポットリザーブ - 50mL 二カ所 10入	¥88,400	¥97,400	
47732-718	MicroSeal 1200 Wipers 600入	¥122,500	¥131,200	
47733-014	Sterile Mini Tongs 10.8cm 25入	¥35,000	¥38,000	
47733-806	濾過フラスコ(フオートホトム) 75ml	¥112,000	¥121,700	
47733-808	濾過フラスコ(フオートホトム) 150ml	¥117,100	¥127,200	
47733-810	濾過フラスコ(フオートホトム) 350ml	¥129,500	¥140,700	
47733-812	濾過フラスコ(フオートホトム) 75ml	¥88,100	¥95,700	
47733-814	濾過フラスコ(フオートホトム) 150ml	¥94,900	¥103,100	
47734-116	ステンレスチールポット 外径12.7mm	¥79,900	¥86,400	
47734-148	Ultra High Vacuum アルミコップ 1	¥230,600	¥246,900	
47734-150	Ultra High Vacuum アルミコップ 1	¥257,800	¥276,000	
47734-152	Ultra High Vacuum アルミコップ 1	¥253,900	¥271,800	
47734-154	Ultra High Vacuum アルミコップ 1	¥274,600	¥293,900	
47734-692	ピペットチップ 200ul ST 96本x10R	¥71,200	¥76,200	
47734-788	96ウェルプレート 1.0mL ST	¥184,700	¥194,800	
47734-817	シリカフィルム アルミコップ 100枚入	¥96,700	¥102,000	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
47735-396	Funnels FLASH Chromatography 5個入	¥113,500	¥169,300	
47735-406	Silica Gel TLC Plate 25枚	¥172,100	¥170,200	
47735-594	ワスコ リンク	¥17,900	¥19,400	
47735-596	ワスコ リンク	¥19,900	¥21,500	
47735-598	ワスコ リンク	¥22,400	¥24,300	
47735-600	ワスコ リンク	¥24,100	¥26,100	
47735-602	攪拌子 (テフロン樹脂コート)	¥56,300	¥60,100	
47735-678	Amplitude Kappa Sterile LE Wipes	¥383,200	¥410,200	
47742-592	Astral Inoculation Loop 1ul	¥200,200	¥222,800	
47742-594	Astral Inoculation Loop 10ul	¥203,100	¥226,000	
47742-598	Clone-Pik Transfer Device	¥13,600	¥15,200	
47743-834	ウルトラマイクロピペット 8.5mm 100本入	¥82,700	¥89,800	
47743-836	BRAND UVマイクロピペット(ホルリスプレ)1.5-3ml	¥376,200	¥408,700	
47743-838	分光光度計(UVマイクロピペット)	¥72,400	¥132,000	
47743-840	分光光度計(UVマイクロピペット)	¥334,000	¥610,000	
47743-842	BRAND UVマイクロピペット用キャップ(ブルー)100個	¥27,000	¥29,400	
47743-940	Disposable Plastic Cuvettes 500入	¥85,200	¥92,600	
47743-982	96-Well Assay Plates	¥80,000	¥84,400	
47743-996	96ウェルアッセイプレートV底 カラー 50枚	¥79,700	¥84,000	
47744-002	96ウェルアッセイプレートV底 カラー 50枚	¥80,100	¥84,500	
47744-004	96ウェルアッセイプレートV底 白 50枚	¥83,800	¥88,500	
47744-642	BRANDマイクロピペット(ホルリスプレ)1.5-3ml	¥83,200	¥90,400	
47744-644	BRANDマイクロピペット(ホルリスプレ) 100	¥21,100	¥22,900	
47744-652	マイクロピペット(Acrylic)1.5-3ml 100	¥21,100	¥22,900	
47744-998	SpinbarMagneticStirring Bar100	¥235,900	¥255,900	
47745-796	FILTER PAPER #131 SIZE4.7cm 100毎	¥33,400	¥36,800	
47745-880	フィルター-パ- 46cmx50cm 1巻	¥88,500	¥97,500	
47745-882	FILTER PAPER 46x57cm	¥69,200	¥76,300	
47745-914	フィルター-パ-222 サイズ20x20cm 100枚	¥41,900	¥46,100	
47746-930	フィルター-パ- CRCLS BLACK 100枚入	¥89,800	¥98,900	
47746-940	ニトリル製グローブ(N-DEX)Sサイズ 50枚	¥19,000	¥20,300	
47746-942	ニトリル製グローブ(N-DEX)Mサイズ 50枚x20	¥345,100	¥369,500	
47746-944	ニトリル製グローブ(N-DEX)Lサイズ 50枚x20	¥345,100	¥369,500	
47746-948	ニトリル製グローブ(N-DEX)XLサイズ 50枚x20	¥345,100	¥369,500	
47747-396	PVAホビニールアルミコートグローブ 1双	¥28,900	¥31,000	
47749-924	Sealing Films Breathable 100	¥119,000	¥125,500	
47750-040	STOPCOCK PP/TFN ID 2WAY 2MM 1個	¥34,100	¥36,700	
47750-042	STOPCOCK PP/TFN ID 2WAY 4MM 1個	¥31,700	¥33,600	
47750-046	STOPCOCK PP/TFN ID 3WAY 4MM 1個	¥37,400	¥39,700	
47750-094	Scoop Sampling System 60ml 100個入	¥192,400	¥208,000	
47750-146	Magnetic Stirring Bar 1入	¥4,300	¥4,500	
47750-178	Flea Micro Spinbar 2x5mm 1本	¥4,100	¥4,400	
47750-180	Flea Micro Spinbar BLUE PTFE 5x2mm	¥4,000	¥4,300	
47750-194	Flea Micro Spinbar 3x6.5mm 1本	¥4,300	¥4,600	
47750-354	100 Pestles Only 100	¥58,000	¥62,300	
47750-422	Large Weighted Safety Shields 1個	¥578,100	¥618,800	
47750-483	50mL Conical Tube Holders 5個	¥29,100	¥31,200	
47750-528	BEL-BLOTTER 96wellリアクティングツール	¥18,200	¥19,400	
47750-726	SCIENCEWARE Carboys with Spigot 4L	¥90,400	¥96,800	
47750-746	バリエーションドリッパ 12x24個	¥44,500	¥47,600	
47750-814	Table Top Buchner Funnel 45.7cm 1	¥613,300	¥638,600	
47750-926	Leveling Table Acrylic	¥116,100	¥123,600	
47751-004	サングラス(PS製/滅菌済み個別包装)	¥117,200	¥125,400	
47751-006	サングラス(PS製/滅菌済み個別包装)	¥123,500	¥132,200	
47751-008	サングラス(PS製/滅菌済み個別包装)	¥168,300	¥180,100	
47751-090	Disposable Plastic Cuvettes 100入	¥18,300	¥19,900	
47751-730	Cryo-Safe Mini Coolers	¥85,600	¥91,600	
47751-736	Cryo-Safe Maxi Coolers	¥169,400	¥181,400	
47751-756	Cryo-Safe -1℃ Freeze Controller	¥62,100	¥67,300	
47751-770	チューブ用ラック(ホルリスプレ)12mm用	¥50,200	¥53,700	
47751-792	ラットグラス染色用シャー 5.9x11.4 1	¥6,100	¥6,600	
47751-794	Poly-Jaque Lab Support Jack	¥222,900	¥269,000	
47751-888	N-Dex Plus Gloves Mサイズ 50枚	¥17,600	¥18,500	
47751-890	N-Dex Plus Gloves Lサイズ 50枚	¥17,600	¥18,500	
47751-892	N-Dex Plus Gloves XLサイズ 50枚	¥17,600	¥18,500	
48218-804	Prefilled Containers 180ml	¥214,700	¥234,400	
48300-026	VWR顕微鏡用ラットグラス	¥10,900	¥12,000	
48300-364	顕微鏡用ラットグラス(カンテック・チャンバ-)	¥15,600	¥17,200	
48300-400	顕微鏡用ラットグラス(カンテック・チャンバ-)	¥36,400	¥40,200	
48300-411	Vapor Line Integrators 250	¥30,000	¥33,100	
48311-601	顕微鏡用ラットグラス(スリ-70スト)	¥17,700	¥19,100	
48311-703	顕微鏡用ラットグラス(スリ-70スト・70スト)	¥27,900	¥30,000	
48349-058	MicroSlides Frosted Fluorescent 72	¥14,600	¥15,800	
48366-045	顕微鏡用カバーグラス(角型)	¥70,600	¥74,900	
48366-067	顕微鏡用カバーグラス(角型)	¥8,200	¥8,800	
48366-089	顕微鏡用カバーグラス(角型) 1oz	¥9,200	¥9,800	
48366-205	顕微鏡用カバーグラス(角型) 1ml	¥9,800	¥10,400	
48366-227	顕微鏡用カバーグラス(角型) 1ml	¥8,100	¥8,700	
48366-249	顕微鏡用カバーグラス(角型) 1ml	¥8,300	¥8,800	
48366-251	顕微鏡用カバーグラス	¥183,600	¥194,700	
48368-040	顕微鏡用カバーグラス(角型) 18x18mm 1oz	¥7,500	¥7,300	
48368-062	顕微鏡用カバーグラス(角型) 22x22mm	¥58,400	¥62,700	
48376-027	顕微鏡用カバーグラス(ガラス製) 1000	¥33,500	¥36,000	
48376-049	顕微鏡用カバーグラス(ガラス製)	¥30,400	¥32,900	
48380-046	カバーグラス丸 No.1 18mm 1oz	¥35,800	¥38,000	
48380-068	カバーグラス丸 No.1 22mm 1oz	¥32,100	¥34,100	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
48380-080	カバーガラス丸 No.1 25mm 1oz	¥32,600	¥34,600	
48382-063	顕微鏡用カバーガラス(丸型)	¥195,400	¥230,800	
48382-085	顕微鏡用カバーガラス(丸型) 1oz	¥30,800	¥32,600	
48382-117	Microscope Adhesion Slides 72枚入	¥41,200	¥44,600	
48382-136	VWR カバーガラスNo.2 24x50mm 21oz.	¥47,600	¥50,500	
48382-139	VWR カバーガラスNo.2 25x60mm 10oz.	¥117,600	¥124,900	
48382-180	Micro Slides Plain 72枚入	¥57,800	¥62,200	
48393-026	顕微鏡用カバーガラス(角型)	¥67,600	¥71,600	
48393-059	顕微鏡用カバーガラス(角型)	¥61,200	¥65,800	
48393-060	顕微鏡用カバーガラス(角型)	¥63,400	¥68,200	
48393-070	顕微鏡用カバーガラス(角型)	¥60,500	¥64,000	
48393-081	顕微鏡用カバーガラス(角型)	¥56,600	¥60,900	
48393-092	顕微鏡用カバーガラス(角型)	¥67,900	¥71,900	
48393-106	顕微鏡用カバーガラス(角型)	¥58,600	¥62,900	
48393-151	Micro Cover Glasses 1oz	¥10,500	¥11,100	
48393-172	顕微鏡用カバーガラス(角型) 1oz.	¥7,800	¥8,300	
48393-230	顕微鏡用カバーガラス(角型)1oz.	¥7,700	¥8,200	
48393-241	顕微鏡用カバーガラス(角型) 1	¥7,800	¥8,400	
48404-031	顕微鏡用カバーガラス(角型)	¥154,800	¥184,000	
48404-111	顕微鏡用カバーガラス(角型)	¥164,900	¥178,300	
48404-133	顕微鏡用カバーガラス(角型)	¥153,300	¥165,800	
48440-003	Slide Grip Peel-A-Way Plastic 6	¥84,000	¥64,400	
48457-107	スライドガラス用キャビネット	¥19,400	¥21,600	
48711-056	バブブコックトル(0.8%牛乳試験用シールド型)	¥35,200	¥38,400	
48711-100	バブブコックトル(10%アイスクリム試験用)	¥30,100	¥32,900	
48714-055	バブブコックトル(50%乳脂肪・チーズ試験用)	¥26,300	¥28,700	
48717-004	バブブコックトル(0.5%脱脂乳試験用)	¥61,600	¥67,200	
48719-000	バブブコックトル(50%チーズ試験用)	¥99,300	¥108,300	
48754-104	ガラス製アイスホールドガラス(0.8%牛乳試験用) ST 1.1	¥446,700	¥487,100	
48754-110	ガラス製アイスホールドガラス(10%アイスクリム試験用) ST 2.2	¥251,800	¥274,600	
48906-007	ミルシヤー	¥277,700	¥331,600	
48906-040	ミルシヤー	¥473,000	¥953,000	
48906-062	ミルシヤー	¥892,100	¥1,797,000	
48906-084	ミルシヤー	¥1,097,100	¥2,202,000	
48906-109	ミルシヤー	¥1,267,200	¥2,546,000	
48906-110	ミルシヤー	¥1,637,700	¥3,277,000	
48906-131	ミルシヤー	¥2,380,200	¥4,745,000	
48908-068	破碎シリンダー	¥66,500	¥64,200	
48908-080	破碎シリンダー	¥309,000	¥379,800	
490003-230	Low-Adhesion ThermoGraduated Tubes	¥40,300	¥43,100	
490003-692	GeneMate 0.2mL 8連PCRチューブ	¥312,500	¥318,200	
490005-658	Clavies Autoclavable Bags 100枚	¥141,900	¥151,100	
490006-654	Dilux Dilution Reservoir Sterile	¥110,800	¥119,800	
490007-322	AeraSeal AdhesiveSealingFilm ST 50	¥49,900	¥53,700	
4ti-003	BCガラス代(プレート貼付付)750μm(コート付)	¥56,500	¥66,800	
4ti-004	BCガラス代(プレート貼付なし)	¥47,600	¥56,100	
4ti-005	BCガラス代 ダブル 1枚分	¥300	¥400	
4ti-007	バブブコート セットアップ料	¥173,400	¥204,300	
4ti-008	BCガラス代 トリアル 1枚分	¥500	¥600	
4ti-009	BreakAwayプレート専用2DJ-ト印字	¥13,400	¥15,900	
4ti-0110	0.3ml 96well MTP 100plates	¥36,200	¥41,000	
4ti-0116	96well DWP 丸穴/U底 100枚	¥65,000	¥76,600	
4ti-0117	96well DWP 丸穴/V底 100枚	¥65,000	¥76,600	
4ti-0120	1.2ml 96well low profile DWP 50枚	¥74,000	¥87,200	
4ti-0125	96 Deep Well Storage Microplate	¥74,000	¥87,200	
4ti-0126	1.2ml 96well 保存プレートLP(角穴) 100	¥109,300	¥128,700	
4ti-0130	2.2ml 96well 保存プレート(丸穴) 50枚	¥101,700	¥119,900	
4ti-0131	12列リザーブ・タイプウェルプレート 25枚	¥87,400	¥102,900	
4ti-0132	2.0ml 96well 保存プレート(角穴) 50	¥101,700	¥119,900	
4ti-0133	290ml オフ・タイプウェルリザーブ 25枚	¥87,400	¥102,900	
4ti-0135	96 cover mat piercable 100mats	¥132,500	¥156,000	
4ti-0136	2.2ml 96well DeepWellPlate角穴	¥68,600	¥80,900	
4ti-0137	96ウェル用シリコンマット(角穴用) 50枚	¥85,100	¥100,400	
4ti-0138	96ウェル用シリコンマット(丸穴用) 50枚	¥85,100	¥100,400	
4ti-0139	384ウェル用シリコンマット(角穴用) 50枚	¥94,100	¥110,800	
4ti-0147	190ul 384well 保存プレート(角穴)100枚	¥185,200	¥218,200	
4ti-0201	VisionPlate384blackTCtreated7枚付24	¥95,200	¥112,300	
4ti-0203	VisionPlate384 175micron sterile	¥54,500	¥64,200	
4ti-0204	VisionPlate384 black 未コーティング 30	¥51,800	¥61,200	
4ti-0214	UltraVisionPlate 384 UV base 30枚	¥167,100	¥196,700	
4ti-0221	VisionPlate96blackTCtreated7枚付 24	¥95,200	¥112,300	
4ti-0223	VisionPlate 96 190um ST 30	¥52,900	¥62,500	
4ti-0224	VisionPlate 96 190um 30	¥50,000	¥58,900	
4ti-0234	UltraVisionPlate 96 UV base 30枚	¥167,100	¥196,700	
4ti-0241	VisionPlate24blackTCtreated7枚付 24	¥95,200	¥112,300	
4ti-0243	VisionPlate 24 190um ST 30	¥52,900	¥62,500	
4ti-0244	VisionPlate 24 190um 30	¥50,000	¥58,900	
4ti-0254	384ウェル アッセイプレート Solid bottom PS10	¥116,800	¥137,600	
4ti-0262	Assay Plate 24 Black 100	¥162,900	¥191,900	
4ti-0263	96ウェル アッセイプレート Solid bottom PS100	¥122,200	¥143,900	
4ti-0264	384-well Assay Plate, black 100	¥122,200	¥143,900	
4ti-0273	Assay Plate 96 White 100	¥122,200	¥143,900	
4ti-0274	384-well Assay Plate, white 100	¥122,200	¥143,900	
4ti-0280	Vision Plate 384 Lid LP 100	¥60,800	¥71,600	
4ti-0281	Vision Plate 384 Lid LP ST 100	¥62,600	¥73,800	
4ti-0282	Assay Plate 96 Lid LP 80	¥39,800	¥47,000	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
4ti-0283	Assay Plate 96 Lid LP ST 80	¥48,900	¥57,700	
4ti-0284	Assay Plate 24 Lid LP 80	¥39,800	¥47,000	
4ti-0285	Microplate lid ultra-low 100	¥61,400	¥72,300	
4ti-0286	Assay Plate 24 Lid LP ST 80	¥48,900	¥57,700	
4ti-0287	4ti-0960/RIG用 79 50	¥41,100	¥48,500	
4ti-0288	PCR Plate用 79 50	¥41,100	¥48,500	
4ti-0289	4ti-0770用 79 50	¥41,100	¥48,500	
4ti-0290	Universal Microplate用 79 50	¥40,200	¥47,500	
4ti-0291	Auto-Sealing PCR Plate Lid 20	¥54,100	¥63,700	
4ti-0292	FrameStrip Adapter Lid 10	¥16,500	¥19,500	
4ti-0310	PhytoStrip 108strips	¥56,100	¥66,300	
4ti-0371	FrameStrip Adapter 18	¥36,000	¥42,500	
4ti-0373	Non-Skirted PCR Plate Adapter 18	¥36,000	¥42,500	
4ti-0380/C	フレームスタ- 384well PCRプレート480 クリア 50	¥63,700	¥75,100	
4ti-0381	フレームスタ-384wellPCRプレート480用白カバー50	¥72,800	¥86,000	
4ti-0382	4Ti-0381 + 4Ti-0560 50セット	¥113,700	¥134,000	
4ti-0383	4Ti-0380/C + 4Ti-0560 50セット	¥104,600	¥123,100	
4ti-0384	フレームスタ-384PCRプレート白フレーム 50	¥66,400	¥78,100	
4ti-0384/B	フレームスタ-384PCRプレートブルーフレーム 50	¥66,400	¥78,100	
4ti-0384/C	フレームスタ- 384well PCRプレートクリア 50	¥48,700	¥56,700	
4ti-0384/G	フレームスタ-384PCRプレートグリーンフレーム 50	¥66,400	¥78,100	
4ti-0384/R	フレームスタ-384PCRプレートレッドフレーム 50	¥66,400	¥78,100	
4ti-0384/X	フレームスタ-384PCRプレートブラックフレーム 50	¥66,400	¥78,100	
4ti-0385	フレームスタ-384PCRプレートブラックフレーム白カバー50	¥72,800	¥86,000	
4ti-0386	フレームスタ-384PCRプレートブラックフレーム 50	¥86,000	¥101,400	
4ti-0387	フレームスタ-384プレートクリアフレーム/ポストカバー50	¥66,400	¥78,100	
4ti-0394	ポリアミドテープ 1	¥4,700	¥5,600	
4ti-0395	Ergo Freeze ビンク 1	¥21,000	¥24,800	
4ti-0396	Ergo Freeze 100- 1	¥21,000	¥24,800	
4ti-0398	Pierce Plate 1	¥135,200	¥159,300	
4ti-0399	Corbett(Qiagen) CAS robot用アダプター	¥49,600	¥58,400	
4ti-0500	PCR 粘着シール 135x80ミリ 100	¥32,100	¥37,500	
4ti-0500/12	インジ-カット粘着シール 12ストリップタイプ 100	¥39,800	¥47,000	
4ti-0500/8	インジ-カット粘着シール 8ストリップタイプ 100	¥39,800	¥47,000	
4ti-0500FL	Foil Seal Strong(137 x 80mm) 100枚	¥41,300	¥48,800	
4ti-0502	Roller for adhesive seals	¥11,200	¥13,100	
4ti-0503	シールテープリカー 1	¥900	¥1,300	
4ti-0510	保存用粘着シール(非耐熱)130x80ミリ 100	¥13,600	¥16,200	
4ti-0512	Peelable Solvent Resistant Foil 100	¥133,200	¥156,800	
4ti-0513	粘着シールforBrooksXPeel 86mmx200m	¥259,800	¥306,000	
4ti-0516/24	ガス透過性粘着シール 24枚用 ST 50枚	¥78,200	¥92,000	
4ti-0516/384	ガス透過性粘着シール 384枚用ST 50	¥81,800	¥96,300	
4ti-0516/96	ガス透過性粘着シール 96枚用 ST 50枚	¥78,200	¥92,000	
4ti-0517	Air-O-Seal 100枚	¥36,200	¥42,700	
4ti-0517/ST	Air-O-Seal ST 100枚	¥41,800	¥49,300	
4ti-0519	InterSeal 100枚	¥80,200	¥94,600	
4ti-0520	Peel Seal 610mx78mm (4880枚)	¥196,400	¥231,400	
4ti-0520/122	Peel-Seal 1000シール/ロール	¥52,500	¥62,000	
4ti-0520S	Peel Seal サフ-用 5mx78mm	¥8,300	¥9,900	
4ti-0521	Peel-Seal 125x78ミリ 100	¥22,600	¥26,800	
4ti-0522	Peel Seal 500mx115mm (6250枚)	¥206,400	¥243,000	
4ti-0522S	Peel Seal サフ-用 5mx115mm	¥8,300	¥9,900	
4ti-0523	Universal Heat Seal 610mx78mm1roll	¥216,000	¥254,400	
4ti-05231	Universal Heat Seal 125x78mm 100枚	¥24,600	¥29,100	
4ti-0523S	ユニバーサルヒートシール サフ-用5mx78mm	¥8,300	¥9,900	
4ti-0524S	Universal Peel Heat Seal サフ-用1	¥8,300	¥9,900	
4ti-0530	Pierce Seal 610mx78mm (4880枚)	¥196,400	¥231,400	
4ti-0530/122	PierceSeal 1000シール/ロール	¥52,500	¥62,000	
4ti-0530/REMP	Pierce Seal 610mx78mm (4880枚)	¥215,800	¥254,100	
4ti-0530S	Pierce Seal サフ-用 5mx78mm	¥8,300	¥9,900	
4ti-0531	PierceSeal 125x78ミリ 100	¥22,600	¥26,800	
4ti-0531/GR	PierceSeal 125x78ミリ グリット付 100	¥22,600	¥26,800	
4ti-0532	Pierce Seal 500mx115mm (6250枚)	¥206,400	¥243,000	
4ti-0532/REMP	Pierce Seal 500mx115mm (6250枚)	¥226,500	¥266,800	
4ti-0532S	Pierce-Seal サフ-用 5mx115mm	¥8,300	¥9,900	
4ti-0535	Foil Seal 610mx78mm (4880枚)	¥196,400	¥231,400	
4ti-0535/122	Foll Seal 1000シール/ロール	¥52,500	¥62,000	
4ti-0535/REMP	Foil Seal 610mx78mm (4880枚)	¥215,800	¥254,100	
4ti-0535S	Foil-Seal サフ-用 5mx78mm	¥8,300	¥9,900	
4ti-0536	Foil-Seal 125mmx78mm 100枚	¥22,600	¥26,800	
4ti-0537	Foil Seal 500mx115mm (6250枚)	¥231,900	¥273,100	
4ti-0537/REMP	Foil Seal 500mx115mm (6200枚)	¥255,400	¥300,700	
4ti-0537S	Foil Seal サフ-用 5mx115mm	¥8,300	¥9,900	
4ti-0538	ヒ-アスシールストンク 78mm x 610m 1	¥326,200	¥384,100	
4ti-05381	ヒ-アスシールストンク 78mmx125mm 100	¥25,700	¥30,300	
4ti-05381/RA	ラダムアケヒ-アスシールストンク 100	¥78,900	¥93,100	
4ti-0538S	サフ-用ヒ-アスシールストンク 78mmx5m	¥8,300	¥9,900	
4ti-0539	ヒ-アスシールストンク 115mmx500m(Agilent)1	¥385,400	¥453,600	
4ti-0540	Clear Seal 500mx78mm (約4000枚)	¥196,400	¥231,400	
4ti-0540/80	Clear Seal 640シール/ロール	¥46,900	¥55,400	
4ti-0540/REMP	Clear Seal 500mx78mm (4000枚)	¥215,800	¥254,100	
4ti-0540/SLIT	ガス透過クリアヒートシール 500mx78mm(4000枚)	¥288,000	¥339,100	
4ti-0540S	Clear Seal サフ-用 5mx78mm	¥8,300	¥9,900	
4ti-0541	Clear Seal 125x78ミリ 100	¥22,600	¥26,800	
4ti-0541/SLIT	Gas Permeable Clear Seal 100	¥31,300	¥36,900	
4ti-0542	Clear Seal 350mx115ミリ/ロール	¥206,200	¥242,700	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
4ti-0542/REMP	Clear Seal 500mx115mm (6250枚)	¥226,100	¥266,300	
4ti-0542S	Clear Seal サンプル用 5mx115mm	¥8,300	¥9,900	
4ti-0545	PS Foil Seal 610mx78mm (4880枚)	¥215,800	¥254,100	
4ti-0545/122	PS Foil Seal 1000シート用	¥57,600	¥68,000	
4ti-0545S	サンプル PS Foil Seal 78mmx5m	¥8,300	¥9,900	
4ti-0546	PS Foil Seal 500mx115mm (6250枚)	¥231,900	¥273,100	
4ti-0547	Polystyrene Foil Seal 100	¥22,600	¥26,800	
4ti-0548	クリアフィルム プラス 115mmx250m(Agilent)	¥147,500	¥173,700	
4ti-05481	クリアフィルム プラス 125x78mm 100シート	¥22,600	¥26,800	
4ti-0549	クリアフィルム プラス 78mmx250m	¥98,700	¥116,300	
4ti-0550	PCR粘着Foilシート 130x80mm 100	¥29,500	¥34,900	
4ti-0560	qPCR用圧着シート 140x77mm 100	¥80,700	¥95,100	
4ti-0561	qPCR用圧着シート 100mm用	¥423,400	¥498,400	
4ti-0563	Optical Film Compression Pad 5枚	¥23,200	¥27,600	
4ti-0565	Q-Stick粘着シート qPCR用 100	¥47,800	¥56,400	
4ti-0566	Piercable Film 100	¥47,800	¥56,400	
4ti-0573	Clear Weld Mk2 610mx78mm (4880枚)	¥225,200	¥265,300	
4ti-0573/122	Clear Weld Seal Mark II1000シート用	¥60,300	¥71,100	
4ti-0573S	Clear Weld Mk2 サンプル用 5mx78mm	¥8,300	¥9,900	
4ti-0574	Clear Weld Mk2 500mx115mm (6250枚)	¥238,200	¥280,400	
4ti-0574S	Clear Weld Mk2 サンプル用 5mx115mm	¥8,300	¥9,900	
4ti-0575	Clear Weld Seal Mark II 100sheets	¥26,600	¥31,400	
4ti-0580	Clear Seal 3730 610mx78mm (4880枚)	¥196,400	¥231,400	
4ti-0580/122	Clear Seal 3730 1000シート用	¥55,600	¥65,500	
4ti-0580S	サンプル ClearSeal 5mx78mm	¥8,300	¥9,900	
4ti-0581	Clear Seal 3730 1000枚分用	¥159,500	¥187,900	
4ti-0582	Clear Seal 3730 500mx115mm(6250枚)	¥252,000	¥296,900	
4ti-0582S	Clear Seal 3730 サンプル用 5mx115mm	¥8,300	¥9,900	
4ti-0585	DMSO耐性PeelSeal 500mx78mm(4000枚)	¥291,600	¥343,400	
4ti-0585/100	DMSO耐性シート (Peelable) 800シート用	¥69,000	¥81,400	
4ti-0585/REMP	DMSO耐性PeelSeal 500mx78mm	¥319,900	¥376,800	
4ti-0585S	サンプル DMSO ResistantPeelSeal 5m	¥8,300	¥9,900	
4ti-0586	DMSO Resistant Peel Seal V11format	¥423,400	¥498,400	
4ti-0586/REMP	DMSO耐性PeelSeal500mx115mm(6250枚)	¥465,400	¥548,000	
4ti-0586S	サンプルDMSO ResistantPeelSeal 5mx115	¥8,300	¥9,900	
4ti-0587	DMSO Resistant Peel Seal 100sheets	¥30,600	¥36,200	
4ti-0590	Thermal Bond 500mx78mm (4000枚)	¥377,300	¥444,300	
4ti-0590/100	Thermal Bond 800シート用	¥89,600	¥105,700	
4ti-0590S	Thermal Bond サンプル用 5mx78mm	¥8,300	¥9,900	
4ti-0591	Thermal Bond 100sheets	¥30,600	¥36,200	
4ti-0592	Thermal Bond 300mm用	¥377,300	¥444,300	
4ti-0592S	Thermal Bond サンプル用 5mx115mm	¥8,300	¥9,900	
4ti-0597	Gas Permeable Heat Seal Mk 2 100	¥22,600	¥26,800	
4ti-0597/ST	Gas Permeable Heat Seal Mk 2 ST100	¥26,600	¥31,400	
4ti-0598	ガス透過ヒートシールMk2 610mx78mm(4880枚)	¥274,800	¥323,700	
4ti-0598/122	GasPermeable HeatSeal Mk2 1000シート用	¥73,500	¥86,700	
4ti-0598S	ガス透過ヒートシールMk2 サンプル用 5mx78mm	¥8,300	¥9,900	
4ti-0599	ガス透過ヒートシールMk2 500mx115mm(6200枚)	¥330,400	¥389,200	
4ti-0599S	ガス透過ヒートシールMk2 サンプル用 5mx115mm	¥8,300	¥9,900	
4ti-0640	サマルテスフィルム (125mmx80mm) 25	¥149,200	¥175,700	
4ti-0641	サマルテスフィルムスターキット 10枚+シリコンパッド	¥84,200	¥99,400	
4ti-0642	サマルテスフィルム ロールキット 1巻+シリコンパッド	¥80,700	¥95,100	
4ti-0665	オートヒートシーラー a4S	¥4,800,000	オープン価格	
4ti-0665-1	A4S用シールドインキカートリッジ	¥29,900	¥31,000	
4ti-0665-11	A4S 追加1年保証 (計3年保証)	¥1,075,200	¥1,107,000	
4ti-0665-113	A4S 追加1年保証 (計4年保証)	¥1,290,200	¥1,329,000	
4ti-0665-12	A4S 追加2年保証 (計4年保証)	¥2,149,800	¥2,214,000	
4ti-0665-2	A4S用プレートホルダータイプA	¥110,300	¥114,000	
4ti-0665-3	A4S用プレートホルダータイプB	¥136,500	¥141,000	
4ti-0665-4	A4S用ロールホルダーセット	¥143,100	¥148,000	
4ti-0665-41	SILA driver for A4S heat sealer 1	¥315,700	¥325,000	
4ti-0665-5	A4S用バキュームカップ フロント用	¥9,600	¥10,000	
4ti-0665-6	A4S用バキュームカップ リア用	¥9,600	¥10,000	
4ti-0665-7	A4S用ガスシフトワイヤー	¥10,600	¥11,000	
4ti-0665-8	A4S用ガスシフトケーブル	¥67,200	¥69,000	
4ti-0710	プレート96PCRプレートスカート付き 黒 50	¥66,400	¥78,100	
4ti-0710/B	プレート96PCRプレートスカート付き 白 50	¥66,400	¥78,100	
4ti-0710/C	プレート96PCRプレート non skirted 50	¥63,700	¥75,100	
4ti-0710/G	プレート96PCRプレートスカート付き グリーン 50	¥66,400	¥78,100	
4ti-0710/R	プレート96PCRプレートスカート付き レッド 50	¥66,400	¥78,100	
4ti-0711	プレート96PCRプレートスカート付き 黒/白 50	¥72,800	¥86,000	
4ti-0720	プレート96PCRプレートLowProfile 黒 50	¥66,400	¥78,100	
4ti-0720/B	プレート96PCRプレートLowProfile 白 50	¥66,400	¥78,100	
4ti-0720/C	プレート96PCRプレート LowProfile NT50	¥63,700	¥75,100	
4ti-0721	プレート96PCRプレートLowProfile黒/白 50	¥72,800	¥86,000	
4ti-0730	プレート96PCRプレートupstandセミスカート 紫 5	¥66,400	¥78,100	
4ti-0730/B	プレート96PCRプレートupstandセミスカート 青 5	¥66,400	¥78,100	
4ti-0730/C	プレート96PCRプレート upstandセミスカート 50	¥63,700	¥75,100	
4ti-0730/C10P	プレート96PCRプレート upstandセミスカート 10	¥14,700	¥17,400	
4ti-0730/G	プレート96PCRプレートupstandセミスカート 緑 5	¥66,400	¥78,100	
4ti-0730/R	プレート96PCRプレートupstandセミスカート 赤 5	¥66,400	¥78,100	
4ti-0730/X	プレート96PCRプレートupstandセミスカート 黒 5	¥66,400	¥78,100	
4ti-0730C0565	4Ti-0730/C + 4Ti-0565 50セット	¥88,000	¥103,700	
4ti-0735	96 PCRプレートセミスカート クリア 50	¥48,000	¥56,600	
4ti-0736	プレート96PCRプレートセミスカート グラス 50	¥50,700	¥59,900	
4ti-0740	96well skirted PCRプレート 50	¥48,000	¥56,600	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
4ti-0741	96well 8ストリッププレート 初付 50枚入	¥50,700	¥59,900	
4ti-0750	96well PCRプレート 黒印字 50枚	¥48,000	¥56,600	
4ti-0750/16	16 well PCR Plate 初付 300 strips	¥56,300	¥66,500	
4ti-0750/16/B	16 well PCR Plate フル 300 strips	¥56,300	¥66,500	
4ti-0750/16/G	16 well PCR Plate グリーン 300 strips	¥56,300	¥66,500	
4ti-0750/16/P	16 well PCR Plate 白-フル 300strips	¥56,300	¥66,500	
4ti-0750/16/R	16 well PCR Plate レッド 300strips	¥56,300	¥66,500	
4ti-0750/16/Y	16 well PCR Plate イロ 300strips	¥56,300	¥66,500	
4ti-0750/24	24 well PCR Plate 初付 200 strips	¥54,300	¥64,000	
4ti-0750/24/B	24 well PCR Plate フル 200 strips	¥54,300	¥64,000	
4ti-0750/24/G	24 well PCR Plate グリーン 200 strips	¥54,300	¥64,000	
4ti-0750/24/P	24 well PCR Plate 白-フル 200strips	¥54,300	¥64,000	
4ti-0750/24/R	24 well PCR Plate レッド 200strips	¥54,300	¥64,000	
4ti-0750/24/Y	24 well PCR Plate イロ 200strips	¥54,300	¥64,000	
4ti-0750/32	32 well PCR Plate 初付 150 strips	¥54,300	¥64,000	
4ti-0750/32/B	32 well PCR Plate フル 150 strips	¥54,300	¥64,000	
4ti-0750/32/G	32 well PCR Plate グリーン 150 strips	¥54,300	¥64,000	
4ti-0750/32/P	32 well PCR Plate 白-フル 150strips	¥54,300	¥64,000	
4ti-0750/32/R	32 well PCR Plate レッド 150strips	¥54,300	¥64,000	
4ti-0750/32/Y	32 well PCR Plate イロ 150strips	¥54,300	¥64,000	
4ti-0750/48	48 well PCR Plate 初付 100 strips	¥52,100	¥61,400	
4ti-0750/48/B	48 well PCR Plate フル 100 strips	¥52,100	¥61,400	
4ti-0750/48/G	48 well PCR Plate グリーン 100 strips	¥52,100	¥61,400	
4ti-0750/48/P	48 well PCR Plate 白-フル 100strips	¥52,100	¥61,400	
4ti-0750/48/R	48 well PCR Plate レッド 100strips	¥52,100	¥61,400	
4ti-0750/48/Y	48 well PCR Plate イロ 100strips	¥52,100	¥61,400	
4ti-0750/8	8 strip PCR Tubes 初付 600 strips	¥56,300	¥66,500	
4ti-0750/8/B	8 strip PCR tube フル 600 strips	¥56,300	¥66,500	
4ti-0750/8/G	8 strip PCR tube グリーン 600 strips	¥56,300	¥66,500	
4ti-0750/8/P	8 strip PCR tube 白-フル 600strips	¥56,300	¥66,500	
4ti-0750/8/R	8 strip PCR tube レッド 600strips	¥56,300	¥66,500	
4ti-0750/8/Y	8 strip PCR tube イロ 600strips	¥56,300	¥66,500	
4ti-0750/B	96well PCRプレート 黒印字 フル 50枚	¥50,700	¥59,900	
4ti-0750/G	96well PCRプレート 黒印字 グリーン 50枚	¥50,700	¥59,900	
4ti-0750/P	96well PCRプレート 黒印字 白-フル 50枚	¥50,700	¥59,900	
4ti-0750/R	96well PCRプレート 黒印字 レッド 50枚	¥50,700	¥59,900	
4ti-0750/TA	Tear-A-Way 96ウェル PCRプレート ノンスカト50	¥53,800	¥63,500	
4ti-0750/TA12	TearAWay 12x8StripPCRプレートノンスカト50	¥53,800	¥63,500	
4ti-0750/W	96well PCRプレート 黒印字 初付 50枚	¥50,700	¥59,900	
4ti-0750/W/TA	Tear-A-Way 96ウェル PCRプレート ノンスカト50	¥57,200	¥67,500	
4ti-0750/Y	96well PCRプレート 黒印字 イロ 50枚	¥50,700	¥59,900	
4ti-0750-25	96well PCRプレート 黒印字 25	¥32,800	¥38,700	
4ti-0751	8連フラットキャップ 300ストリップ	¥30,200	¥35,700	
4ti-0752	8連ドームキャップ 300ストリップ	¥30,200	¥35,700	
4ti-0753	リム-パール8ウェルチューブストリップ 単体 LP 120	¥26,400	¥31,100	
4ti-0753/757	96ウェルPCRプレート/8連ドームキャップ 50	¥160,600	¥189,100	
4ti-0753/950W	Vari-Strip 初付(プレート/セミカト) 50	¥160,600	¥189,100	
4ti-0753/C	リム-パール8ウェルチューブストリップ 単体 LP 120	¥26,400	¥31,100	
4ti-0753/C757	96ウェルPCRプレート/8連ドームキャップ 50	¥160,600	¥189,100	
4ti-0754	8連LP初付チューブ (4ti-0753)+フラット 120	¥34,000	¥40,200	
4ti-0754/C	8連LP初付チューブ (4ti-0753/C)+フラット 120	¥34,000	¥40,200	
4ti-0755	CrystalStrip 8連フラットキャップ 300	¥31,500	¥37,200	
4ti-0755/120	CrystalStrip 8連フラットキャップ 120	¥14,300	¥16,900	
4ti-0757-F	96ウェルPCRプレート/フレーム ノンスカト 10	¥7,400	¥8,900	
4ti-0760	96well PCRプレート セミカト付 50枚	¥26,000	¥25,300	
4ti-0761	96well PCRプレート セミカト 初付 50枚	¥50,700	¥59,900	
4ti-0762	96 well PCRプレート セミカト ノンスカト 50	¥48,000	¥56,600	
4ti-0770	フレームスタ-96PCRプレート セミカト 白-フル 50	¥66,400	¥78,100	
4ti-0770/B	フレームスタ-96PCRプレート セミカト フル 50	¥66,400	¥78,100	
4ti-0770/C	フレームスタ-0.2mL 96PCRプレート セミカト 50	¥40,200	¥45,500	
4ti-0770/C10P	フレームスタ-96セミカト クリアフレーム 10	¥14,700	¥17,400	
4ti-0770/G	フレームスタ-96PCRプレート セミカト グリーン 50	¥66,400	¥78,100	
4ti-0770/R	フレームスタ-96PCRプレート セミカト レッド 50	¥66,400	¥78,100	
4ti-0770/X	フレームスタ-96PCRプレート セミカト フラット 50	¥66,400	¥78,100	
4ti-0771	フレームスタ-96PCRプレート セミカト付 黒/白50	¥72,800	¥86,000	
4ti-0772	フレームスタ-96PCRプレート セミカト ノンスカト50	¥41,800	¥48,800	
4ti-0772/10P	フレームスタ-96PCRプレート セミカト ノンスカト50	¥15,900	¥18,700	
4ti-0778	Cap Mat for PCR Plates	¥113,900	¥134,300	
4ti-0780	8連0.2ml PCRTubes and ドームcaps 125	¥34,900	¥41,200	
4ti-0781	8連0.2ml PCR tubes 125	¥27,300	¥32,100	
4ti-0782	8 domed caps 125 strips	¥15,900	¥18,700	
4ti-0783	8連フラットキャップ 125ストリップ	¥15,900	¥18,700	
4ti-0784	8連0.2mlPCRチューブ+フラット初付キャップ 12	¥34,900	¥41,200	
4ti-0785/B	FrameStrip 8連初付チューブ+ドームCAP付120	¥34,200	¥40,500	
4ti-0785/G	FrameStrip 8連初付チューブ+ドームCAP付120	¥34,200	¥40,500	
4ti-0785/M	FrameStrip 8連初付チューブ+ドームCAP付120	¥34,200	¥40,500	
4ti-0785/P	FrameStrip 8連初付チューブ+ドームCAP付120	¥34,200	¥40,500	
4ti-0785/R	FrameStrip 8連初付チューブ+ドームCAP付120	¥34,200	¥40,500	
4ti-0785/W/2D	FrameStrip 8連初付チューブ+ドームCap 120	¥40,400	¥47,800	
4ti-0785/X	FrameStrip 8連初付チューブ+ドームCAP付120	¥34,200	¥40,500	
4ti-0785/X/2D	FrameStrip 2DJ+ドーム付8連初付チューブ 120	¥40,400	¥47,800	
4ti-0785/XW	FrameStrip8連初付チューブ+ドームCAP付120	¥37,100	¥43,700	
4ti-0786/B	FrameStrip 8連初付チューブフラットCAP付120	¥34,200	¥40,500	
4ti-0786/G	FrameStrip 8連初付チューブフラットCAP付120	¥34,200	¥40,500	
4ti-0786/M	FrameStrip 8連初付チューブフラットCAP付120	¥34,200	¥40,500	
4ti-0786/P	FrameStrip 8連初付チューブフラットCAP付120	¥34,200	¥40,500	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
4ti-0786/R	FrameStrip 8連カチーフ フラットCAP付120	¥34,200	¥40,500	
4ti-0786/W/2D	FrameStrip 8連カチーフ +フラットCap 120	¥40,400	¥47,800	
4ti-0786/X	FrameStrip 8連カチーフ フラットCAP付120	¥34,200	¥40,500	
4ti-0786/X/2D	FrameStrip 2Dジャケット付8連カチーフ 120	¥40,400	¥47,800	
4ti-0786/XW	FrameStrip8連カチーフ フラットCAP付120	¥37,100	¥43,700	
4ti-0788	12連Opticalチップ 200	¥29,900	¥35,400	
4ti-0790	0.2ml SinglePCR Tube FlatCap 1000	¥17,200	¥20,500	
4ti-0790/2D	0.2mlシングルPCRチューブ フラットCap 960	¥92,900	¥109,500	
4ti-0792	8 strip 200mu PCR Tubes チップ付120	¥30,800	¥36,400	
4ti-0793	8 strip 100mu PCR Tubes チップ付120	¥32,400	¥38,200	
4ti-0794	8stripPCR Tubes w/DomedCap120	¥27,900	¥32,900	
4ti-0795	0.2ml SinglePCR Tube Domed Cap1000	¥17,200	¥20,500	
4ti-0796	4連PCRチューブ cap付 250	¥50,300	¥59,200	
4ti-0900/C	フレームスター96ウェルカート クリアウェル/フレーム 50	¥63,700	¥75,100	
4ti-0901	フレームスター96ウェルカート 白ウェル 50	¥72,800	¥86,000	
4ti-0910/C	フレームスター-Fast 96ウェルカート クリアウェル/フレーム50	¥63,700	¥75,100	
4ti-0910/C10P	フレームスター-Fast 96ウェルカート クリアウェル/フレーム10	¥14,700	¥17,400	
4ti-0910C0565	4Ti-0910/C + 4Ti-0565 50セット	¥88,000	¥103,700	
4ti-0911	フレームスター-Fast96ウェルカート白ウェル/クリアフレーム50	¥72,800	¥86,000	
4ti-0912	フレームスター-Fast 96ウェルカート 白ウェル 50	¥41,800	¥48,800	
4ti-0912/10P	フレームスター-Fast 96ウェルカート 白ウェル 10	¥14,700	¥17,400	
4ti-0950/C	フレームスター480 96PCRプレートウェルカートクリア-50	¥63,700	¥75,100	
4ti-0950W-F	Vari-Strip 480/96 フレーム 白 10	¥7,400	¥8,900	
4ti-0951	フレームスター480 96PCRプレートウェルカート白ウェル50	¥72,800	¥86,000	
4ti-0951/0565	4Ti-0951 + 4Ti-0565 50セット	¥96,300	¥113,500	
4ti-0952	4Ti-0951 + 4Ti-0560 50セット	¥113,700	¥134,000	
4ti-0953	4Ti-0950/C + 4Ti-0560 50セット	¥104,600	¥123,100	
4ti-0954	96PCRプレートウェルカートLC480 エキストラウェル 50	¥78,200	¥92,000	
4ti-0954/0560	96PCRプレートウェルカートEXウェル+qPCRシール各50	¥118,000	¥139,100	
4ti-0955	96PCRプレートウェルカート LC480 白 50	¥50,700	¥59,900	
4ti-0955/0560	96PCRプレートウェルカート 白+qPCRシール 各50	¥93,800	¥110,500	
4ti-0960	フレームスター96 PCRプレートウェルカートフレーム 50	¥66,400	¥78,100	
4ti-0960/0565	4Ti-0960 + 4Ti-0565 50セット	¥90,900	¥107,200	
4ti-0960/B	フレームスター96 PCRプレートウェルフレーム 50	¥66,400	¥78,100	
4ti-0960/C	フレームスター96 PCRプレートウェルカートクリアウェル/フレーム 50	¥63,700	¥75,100	
4ti-0960/C-F	Ergo Plate Stand用アダプター 10個	¥7,400	¥8,900	
4ti-0960/G	フレームスター96 PCRプレートウェルカートフレーム 50	¥66,400	¥78,100	
4ti-0960/R	フレームスター96 PCRプレートウェルカートフレーム 50	¥66,400	¥78,100	
4ti-0960/RA	ランダムアクセス96wellプレート スカート IDVチューブ	¥80,200	¥94,600	
4ti-0960/RA-F	ランダムアクセス96wellプレート付用フレーム 白 10	¥11,200	¥13,100	
4ti-0960/RIG	フレームスター96 PCRプレートウェルカートリジッドカート50	¥66,400	¥78,100	
4ti-0960/W	フレームスター96 PCRプレートウェルカートフレーム 50	¥66,400	¥78,100	
4ti-0960/X	フレームスター96 PCRプレートウェルカートフレーム 50	¥66,400	¥78,100	
4ti-0961	フレームスター96PCRプレートウェルカートフレーム白ウェル50	¥72,800	¥86,000	
4ti-0966	フレームスター96PCRプレートウェルカート黒50	¥86,000	¥101,400	
4ti-0970	フレームスター96wellプレート スカート	¥80,200	¥94,600	
4ti-0970/RA	ランダムアクセス96wellプレート スカート IDVチューブ	¥80,200	¥94,600	
4ti-0970/RA-F	ランダムアクセス96wellプレート付用フレーム 白 10	¥11,200	¥13,100	
4ti-0975/RA	ランダムアクセス96wellプレート スカート IDVチューブ	¥202,000	¥237,900	
4ti-1000/B	FrameStar Break-A-Way 96well Plate	¥66,800	¥78,600	
4ti-1000/C	FrameStar Break-A-Way 96well Plate	¥66,800	¥78,600	
4ti-1000/G	FrameStar Break-A-Way 96well Plate	¥66,800	¥78,600	
4ti-1000/P	FrameStar Break-A-Way 96well Plate	¥66,800	¥78,600	
4ti-1000/R	FrameStar Break-A-Way 96well Plate	¥66,800	¥78,600	
4ti-1000/W	フレームスター-BreakAWay 96well Plate 50	¥66,800	¥78,600	
4ti-1000/X	FrameStar Break-A-Way 96well Plate	¥66,800	¥78,600	
4ti-1001	FrameStar Break-A-Way 96well Plate	¥75,300	¥88,800	
4ti-1200/B	FrameStar Break-A-Way 96 LP 50	¥66,800	¥78,600	
4ti-1200/C	FrameStar Break-A-Way 96 LP 50	¥66,800	¥78,600	
4ti-1200/G	FrameStar Break-A-Way 96 LP 50	¥66,800	¥78,600	
4ti-1200/P	FrameStar Break-A-Way 96 LP 50	¥66,800	¥78,600	
4ti-1200/R	FrameStar Break-A-Way 96 LP 50	¥66,800	¥78,600	
4ti-1200/W	フレームスター-BreakAWay 96well PlateLP 50	¥66,800	¥78,600	
4ti-1200/X	FrameStar Break-A-Way 96 LP 50	¥66,800	¥78,600	
4ti-1201	FrameStar Break-A-Way96 LP 白ウェル50	¥75,300	¥88,800	
4ti-1300/X	FrameStarBreak-2-Ways96wellプレート50	¥66,800	¥78,600	
4ti-1400/X	Break2Ways96wellプレートLowprofile 50	¥66,800	¥78,600	
4ti-4060	ハンディ 2Dハートスキャナー USBケーブル付	¥150,500	¥155,000	
4ti-LB0109	96well Storage Microplate LB	¥52,100	¥61,400	
4ti-LB0125	96well 丸穴DWP for MagStand LB	¥84,200	¥99,400	
4ti-LB0770/C	フレームスター-LB 96PCRプレートウェルカート付 50	¥73,700	¥87,000	
4ti-LB0960RIG	フレームスター-LB96PCRプレートウェルカートリジッドカート50	¥76,400	¥90,000	
4ti-OX730CSBC	フレームスター96 upstandウェルカート BC付 20	¥41,100	¥48,500	
4ti-OX770CSBC	フレームスター96PCRプレートウェルカート BC付 20	¥41,100	¥48,500	
4ti-OX960CSBC	フレームスター96PCRプレートウェルカート BC付 20	¥41,100	¥48,500	
50000-122	Mouse Jacket small	¥38,200	¥40,100	
50000-410	Heavy-Dutyレザーケース Mサイズ	¥356,900	¥374,600	
50000-412	Heavy-Dutyレザーケース Mサイズ	¥356,900	¥374,600	
50000-414	Heavy-Dutyレザーケース Lサイズ	¥365,900	¥384,100	
50405-026	乳鉢(ガラス製)	¥507,200	¥550,900	
50405-048	乳鉢(ガラス製)	¥696,900	¥757,000	
50415-020	乳鉢(ガラス製)	¥13,300	¥14,700	
50415-041	乳鉢(ガラス製)	¥14,100	¥15,500	
50415-063	乳鉢(ガラス製)	¥16,600	¥18,300	
50415-085	乳鉢(ガラス製)	¥22,600	¥24,900	
50415-100	乳鉢(ガラス製)	¥33,000	¥36,400	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
50806-404	白金耳(ル-フ)	¥104,500	¥112,500	
50806-741	白金耳(ル-フ) 7.6cm	¥28,900	¥31,000	
50806-763	白金耳(ル-フ) 7.6cm	¥28,700	¥30,900	
50806-800	白金耳(ル-フ) 7.6cm	¥28,900	¥31,000	
50806-865	白金耳(ル-フ) 7.6cm	¥28,700	¥30,900	
50806-901	白金耳(ル-フ) 7.6cm	¥28,900	¥31,000	
50806-930	加マトクワイ(HPLC用フイタインク)	¥39,600	¥45,400	
50806-932	加マトクワイ(HPLC用フイタインク)	¥44,100	¥48,200	
50807-008	白金耳(ル-フ) プラチ製ワイヤ	¥118,800	¥133,400	
50807-020	白金耳(ル-フ) プラチ製ワイヤ	¥115,700	¥130,000	
50807-041	白金耳(ル-フ) プラチ製ワイヤ	¥128,200	¥144,000	
50807-085	白金耳(ル-フ) プラチ製ワイヤ	¥76,700	¥86,200	
50807-136	ニードル ステンレスチール 22ゲージ(ル-チップ)	¥27,500	¥34,900	
50807-520	加マトクワイ(HPLC用フイタインク)	¥29,400	¥32,100	
50809-020	接種用タンデム(Large)	¥134,800	¥145,000	
50809-022	接種用タンデム(small)	¥107,600	¥116,100	
50814-003	白金耳(ニードル) アルミニウム製ハット付	¥13,300	¥17,600	
50815-028	白金耳(ル-フ) 12本入	¥67,200	¥72,300	
50815-040	白金耳(ル-フ) 12本入	¥44,500	¥59,000	
50815-200	白金耳(ル-フ)	¥12,900	¥13,900	
50817-250	ニードル 刺タ-	¥68,300	¥71,700	
51280-752	クリン-ム用バインダー	¥595,500	¥659,300	
51280-756	クリン-ム用バインダー	¥923,400	¥959,100	
51280-796	クリン-ム用バインダー	¥490,800	¥603,000	
51280-812	クリン-ム用アンダー(フィルムチップ)12本	¥71,500	¥82,500	
51280-826	クリン-ム用バインダー	¥912,900	¥977,300	
51280-874	クリン-ム用アンダー(トラック)10本	¥39,000	¥49,500	
51280-882	クリン-ム用スリット(耐薬品性)	¥200,600	¥336,200	
51280-888	クリン-ム用スリット(耐薬品性)	¥254,800	¥315,200	
51280-916	クリン-ム用スリット(耐薬品性)	¥368,300	¥406,800	
52100-106	ガラスバインダー	¥168,900	¥195,700	
52196-124	CALIBRATION STANDARD 50 PK10	¥98,700	¥155,300	
52196-126	SOLUTION 290 REFERENCE 2ML PK10	¥123,500	¥178,500	
52196-128	STANDARD CALIB 850 MILLIO PK10	¥106,000	¥156,100	
52197-167	Sample Tubes 0.2mL 500本入	¥976,200	¥1,468,900	
52197-416	Sample Tubes and Probe Wipes 250	¥533,100	¥574,000	
52197-452	Osmometer Sampler 20ul	¥672,500	¥1,572,000	
52201-585	真空ポンプ	¥133,400	¥247,000	
52202-564	ウォーターバス(振盪機能付)	¥21,200	¥39,000	
52202-566	ウォーターバス(振盪機能付)	¥26,000	¥48,000	
52202-568	ウォーターバス(振盪機能付)	¥30,500	¥57,000	
52202-570	ウォーターバス(振盪機能付)	¥44,600	¥83,000	
52202-572	ウォーターバス(振盪機能付)	¥40,100	¥74,000	
52202-574	ウォーターバス(振盪機能付)	¥79,100	¥147,000	
52457-010	溶存酸素プローブ(YSI専用)	¥353,900	¥650,000	
52457-064	B.O.D. Bottle Probe Non-stirring	¥334,900	¥615,000	
52457-552	キット(Standard Membranes)	¥39,600	¥73,000	
52457-585	キット(High Sensitivity Membranes)	¥31,200	¥58,000	
52457-610	オリガキット(52457-010/52457-064用)	¥16,000	¥30,000	
52845-009	レンズクリナー(ウェットタイプ) 50シートx12PK	¥15,000	¥17,400	
52846-001	レンズクリナー(ウェットタイプ) 50シートx12PK	¥23,100	¥29,300	
52846-007	レンズクリナー(ウェットタイプ) 100シートx12	¥31,200	¥37,700	
52847-150	Lenx 90 Lens Tissues 1000入	¥75,100	¥80,400	
52848-058	レンズクリナー(プラスチック製BOX入り)	¥142,900	¥153,000	
52848-073	レンズクリナー(金属製BOX入り)	¥61,900	¥66,300	
52857-098	Super Absorbent Blue Soakers 175	¥248,000	¥288,400	
52857-112	Super Absorbent Blue Soakers 1	¥128,500	¥149,600	
52857-116	Super Absorbent Blue Soakers 350	¥248,000	¥288,400	
52857-120	Standard White Soakers	¥233,200	¥270,900	
52859-079	Parafilm M Bemis 50.8cm x 15.2m	¥46,600	¥49,900	
52865-005	けがきペン(Diamond Scriber) 1	¥13,000	¥14,200	
52865-118	けがきペン(Replacement Point)	¥8,700	¥9,500	
52877-490	Securlink Blood Recipient System	¥437,600	¥479,900	
52877-810	マーカー(長期保存用)Assorted	¥21,600	¥23,600	
52877-842	マーカー(長期保存用) Black 6本	¥17,200	¥18,800	
52877-864	マーカー(長期保存用) Red	¥13,100	¥14,400	
52877-886	マーカー(長期保存用) Blue	¥12,500	¥13,600	
52877-900	マーカー(長期保存用) Green	¥12,500	¥13,600	
52934-505	PENETROMETER W/PLUNGER WEIGHTS	¥2,461,100	¥6,105,000	
52934-640	Magnesium Grease Cone	¥819,700	¥2,035,000	
53047-600	ピペット(綿栓付 ガラス/ステンレス)1ml	¥249,800	¥272,400	
53047-602	ピペット(綿栓付 ガラス/ステンレス)5ml	¥301,300	¥328,600	
53047-604	ピペット(綿栓付 ガラス/ステンレス)10ml	¥326,200	¥355,700	
53047-606	ピペット(綿栓付 ガラス/ステンレス)50ml	¥240,400	¥262,200	
53047-608	ピペット(綿栓付 ガラス/ステンレス)25ml	¥511,700	¥558,000	
53047-610	ピペット(綿栓付 ガラス/ステンレス)1ml	¥256,200	¥279,400	
53047-612	ピペット(綿栓付 ガラス/ステンレス)5ml	¥252,600	¥275,500	
53047-614	ピペット(綿栓付 ガラス/ステンレス)10ml	¥298,700	¥325,700	
53047-616	ピペット(綿栓付 ガラス/ステンレス)25ml	¥449,800	¥490,500	
53047-618	ピペット(綿栓付 ガラス/ステンレス)50ml	¥291,900	¥318,300	
53200-100	ドライバイ- 1個	¥41,900	¥44,800	
53200-102	フラスコ(ほうき)	¥44,100	¥47,200	
53281-816	ピペット 5ml/11ml(ガラス製)10ml	¥589,000	¥642,300	
53281-818	ピペット 5ml/10ml(ガラス製) 5ml	¥571,800	¥623,500	
53283-685	加液チューブ 口径 20x125mm 500	¥110,200	¥119,800	
53283-686	加液チューブ 口径 20x150mm 34ml 500	¥120,000	¥130,300	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
53283-800	加チャチューブ 袖口 13x100mm 1000	¥156,600	¥170,100	
53283-802	加チャチューブ 袖口 16x100mm 1000	¥151,300	¥164,400	
53283-806	加チャチューブ 袖口 16x150mm 1000	¥204,700	¥222,400	
53283-916	バースピペット(オートタイプ ガラス製)2ml	¥80,300	¥87,200	
53384-192	ミルクピペット(PS製/滅菌済みタイプ)	¥513,300	¥562,400	
53384-240	ミルクピペット(PS製/滅菌済みタイプ)	¥594,100	¥651,000	
53432-604	Drummond Disposable Micropipets1-5	¥63,800	¥69,000	
53432-706	マイクロピペット(ガラス製) 5ul 250	¥61,100	¥65,100	
53432-728	マイクロピペット(ガラス製) 10ul 250本x4	¥159,400	¥169,800	
53432-740	マイクロピペット(ガラス製)	¥156,400	¥166,500	
53432-761	マイクロピペット(ガラス製)	¥160,900	¥171,400	
53432-783	マイクロピペット(ガラス製)250x4PK	¥160,900	¥171,400	
53432-808	マイクロピペット(ガラス製)	¥159,800	¥170,200	
53432-921	マイクロピペット(ガラス製)100ul 250	¥61,100	¥65,100	
53481-009	ガラスピペット	¥223,800	¥245,300	
53481-064	ガラスピペット	¥225,600	¥247,300	
53481-097	ガラスピペット 10mL	¥232,000	¥256,900	
53481-122	ガラスピペット	¥236,900	¥262,300	
53481-155	ガラスピペット	¥244,600	¥270,900	
53481-188	ガラスピペット	¥248,100	¥271,800	
53481-213	ガラスピペット	¥265,600	¥294,100	
53481-246	ガラスピペット	¥292,800	¥320,900	
53481-280	ガラスピペット	¥157,800	¥172,800	
53481-348	ガラスピペット	¥166,000	¥183,800	
53481-370	ガラスピペット	¥173,600	¥192,200	
53481-406	ガラスピペット	¥173,700	¥192,300	
53481-461	ガラスピペット	¥190,100	¥210,400	
53481-494	ガラスピペット	¥207,800	¥227,800	
53481-520	ガラスピペット	¥223,200	¥247,200	
53497-009	ピペットフイラー	¥16,400	¥17,800	
53497-053	ピペットフイラー(レット)	¥16,200	¥17,500	
53497-055	VWR FLIP ピペットフイラー	¥18,600	¥20,100	
53497-086	ピペットフイラー(ブルー)	¥16,400	¥17,800	
53497-111	ピペットフイラー	¥16,500	¥17,900	
53497-304	Bel-Bulb Pipettor	¥18,500	¥19,900	
53498-103	PORTABLE PIPET-AID Drummond 110V 1	¥223,300	¥400,000	
53498-205	Replacement Filters for Pipet-Aid	¥25,500	¥29,900	
53500-409	連続分注シリンジ	¥28,000	¥30,700	
53500-410	連続分注シリンジ	¥32,700	¥35,800	
53500-420	連続分注シリンジ	¥31,200	¥34,200	
53500-431	連続分注シリンジ	¥48,600	¥53,300	
53500-442	連続分注シリンジ	¥43,400	¥47,500	
53500-497	連続分注シリンジ	¥47,100	¥51,600	
53500-500	連続分注シリンジ	¥50,400	¥55,200	
53500-511	連続分注シリンジ	¥53,100	¥58,200	
53500-522	連続分注シリンジ	¥70,500	¥77,200	
53500-588	連続分注シリンジ	¥16,000	¥17,500	
53500-599	連続分注シリンジ	¥37,200	¥40,700	
53500-602	連続分注シリンジ	¥22,100	¥24,200	
53500-613	連続分注シリンジ	¥20,800	¥22,800	
53500-636	連続分注シリンジ	¥24,400	¥26,800	
53500-637	連続分注シリンジ	¥25,700	¥28,100	
53500-657	連続分注シリンジ	¥17,800	¥19,500	
53501-144	ボトルトップディスペンサー(キャリカレックス)10ml	¥448,600	¥491,600	
53501-650	ボトルトップディスペンサー(モデル506)	¥41,300	¥45,300	
53501-652	ボトルトップディスペンサー(モデル505)	¥25,300	¥27,700	
53501-654	ボトルトップディスペンサー(モデル507)	¥108,000	¥118,300	
53502-211	ピペットホウ 0.2mL	¥15,900	¥17,000	
53502-222	ピペットホウ 2mL	¥13,900	¥14,900	
53502-233	ピペットホウ 10mL	¥15,400	¥17,700	
53502-244	ピペットホウ 25mL	¥20,400	¥21,900	
53502-251	ピペットホウ II 2mL	¥20,200	¥21,600	
53502-253	ピペットホウ II 10mL	¥19,300	¥20,700	
53502-255	ピペットホウ II 25mL	¥26,000	¥27,900	
53506-900	VWRマイクロピペット用交換用チップ	¥13,800	¥14,700	
53506-910	VWRマイクロピペット用交換用チップ	¥14,300	¥15,300	
53506-912	VWRマイクロピペット用交換用チップ	¥14,100	¥15,000	
53506-914	VWRマイクロピペット用交換用チップ	¥14,300	¥15,300	
53506-922	VWRマイクロピペット用交換用チップ	¥14,300	¥15,300	
53506-924	VWRマイクロピペット用交換用チップ	¥14,300	¥15,300	
53506-926	VWRマイクロピペット用交換用チップ	¥14,300	¥15,300	
53506-928	VWRマイクロピペット用交換用チップ	¥14,100	¥15,000	
53507-268	Bulb Assembly for Micropipet	¥11,500	¥12,500	
53507-368	容量固定式マイクロピペット用ディスペンサー 5ul	¥199,200	¥212,300	
53508-375	VWRマイクロピペット用チップリポア	¥11,900	¥12,600	
53508-386	VWRマイクロピペット用チップリポア	¥11,900	¥12,700	
53508-400	マイクロピペット用チップリポア	¥12,100	¥12,900	
53508-422	VWRマイクロピペット用チップリポア 20ul	¥11,400	¥12,100	
53508-535	VWRマイクロピペット用チップリポア 20ul	¥11,900	¥12,700	
53508-590	VWRマイクロピペット用チップリポア 250ul	¥12,000	¥12,700	
53509-300	12連PCRチューブのみ 80本	¥64,200	¥65,300	
53509-302	12連チップのみ 80本	¥29,100	¥29,700	
53509-304	8連/12連PCRチューブ&チップ	¥75,300	¥76,700	
53509-306	12連PCRチューブ+チップ(別袋入) 80	¥74,300	¥75,700	
53547-000	Plastic Syringes 1ML CS7000	¥578,400	¥634,400	
53548-005	ガラスシリンジ 5mL 100本	¥14,100	¥14,800	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
53548-017	HSW Norm-Jectシリンジ (ルア-ロック式)3mL100	¥16,700	¥17,600	
53548-019	HSW Norm-Jectシリンジ ルア-ロック5mL 100入	¥19,700	¥20,600	
53548-023	HSW Norm-Jectシリンジ ルア-ロック 10mL 100	¥22,200	¥23,300	
53548-025	HSW Norm-Jectシリンジ ルア-ロック20mL 100	¥33,300	¥35,000	
53550-116	PCRチップ 8連ストリップ	¥36,500	¥39,300	
53550-118	PCRチップ 8連ストリップ	¥36,000	¥38,700	
53550-122	PCRチップ 8連ストリップ	¥36,200	¥38,900	
53550-126	PCRチップ 8連ストリップ	¥39,300	¥42,300	
53550-128	PCRチップ 8連ストリップ	¥39,700	¥42,700	
53550-130	PCRチップ 8連ストリップ	¥39,700	¥42,700	
53550-132	PCRチップ 8連ストリップ	¥39,700	¥42,700	
53550-134	PCRチップ 8連ストリップ	¥39,700	¥42,700	
53574-100	ピペッティン 1個入	¥30,500	¥32,600	
53576-207	ピペッター	¥162,800	¥176,000	
53577-145	パック(ピペッター) 10本収納	¥78,900	¥84,500	
53577-148	パック(ピペッター) 10本収納	¥69,000	¥73,800	
53600-108	ピペーター(ホリフック)	¥5,100	¥5,400	
53600-141	ピペーター(ホリフック)	¥9,500	¥10,100	
53600-185	ピペーター(ホリフック)	¥14,200	¥15,300	
53600-221	ピペーター(ホリフック) 3000ml	¥19,300	¥20,700	
53600-265	ピペーター(ホリフック)	¥53,500	¥57,200	
53800-008	スプレッダー - BactiCellSpreaders 1	¥13,800	¥14,700	
53800-010	スプレッダー - BactiCellSpreaders 1	¥14,400	¥15,300	
53801-008	攪拌棒	¥3,500	¥3,800	
53873-009	プレス機(小型)	¥916,500	¥1,788,000	
54110-300	実験台保護シート(吸水性)	¥54,400	¥59,600	
54110-560	滅菌物用カバー	¥71,700	¥75,200	
54110-562	滅菌物用カバー	¥91,600	¥96,200	
54110-564	滅菌物用カバー	¥184,900	¥194,100	
54110-568	滅菌物用カバー	¥309,500	¥324,800	
54110-570	滅菌物用カバー	¥576,600	¥605,200	
54110-574	滅菌物用カバー	¥154,100	¥161,700	
54110-576	滅菌物用カバー	¥330,100	¥346,400	
54110-580	滅菌物用カバー	¥88,100	¥94,800	
54110-586	滅菌物用カバー	¥148,100	¥155,500	
54112-100	実験台保護シート(撥水性)の粘着シート	¥457,000	¥521,300	
54112-256	実験台保護シート(撥水性)の粘着シート	¥397,000	¥452,800	
54848-004	液体洗剤(中性/ガラス器具用)19L	¥149,200	¥180,100	
54848-314	トルクレンチ(ター用)	¥321,200	¥387,700	
54848-316	トルクレンチ(ター用)	¥71,300	¥83,000	
54848-318	トルクレンチ(ター用)	¥59,400	¥69,200	
54970-544	真空ポンプ(オイルフリー)	¥1,441,400	¥2,622,000	
54972-709	真空ポンプ(オイル駆動式)	¥1,377,400	¥2,518,000	
54973-075	真空ポンプ(オイル駆動式二段高真空型)	¥1,532,600	¥2,801,000	
54973-084	真空ポンプ(オイル駆動式二段高真空型)	¥1,524,100	¥2,786,000	
54973-133	真空ポンプ(オイル駆動式二段高真空型)	¥897,200	¥1,640,000	
54973-213	真空ポンプ(オイル駆動式二段高真空型)	¥1,424,800	¥2,604,000	
54974-238	真空ポンプ(オイル駆動式二段高真空型)	¥2,141,300	¥3,914,000	
54974-247	真空ポンプ(オイル駆動式二段高真空型)	¥1,555,000	¥2,842,000	
54974-262	真空ポンプ(オイル駆動式二段高真空型)	¥1,919,500	¥3,508,000	
54976-029	真空ポンプ(オイル駆動式二段高真空型)	¥2,196,500	¥4,015,000	
54976-035	真空ポンプ(オイル駆動式二段高真空型)	¥1,654,400	¥3,024,000	
54976-038	真空ポンプ(オイル駆動式二段高真空型)	¥1,697,000	¥3,102,000	
54976-053	真空ポンプ(オイル駆動式二段高真空型)	¥2,075,400	¥3,793,000	
54976-110	真空ポンプ(オイル駆動式二段高真空型)	¥1,424,800	¥2,604,000	
54987-070	真空ポンプ(オイル駆動式二段高真空型)	¥2,201,100	¥4,023,000	
54987-203	真空ポンプ(オイル駆動式)	¥1,663,200	¥3,040,000	
54993-345	真空ポンプ(オイル駆動式二段高真空型)	¥4,352,200	¥7,955,000	
54993-505	真空ポンプ(オイル駆動式二段高真空型)	¥3,550,400	¥6,489,000	
54993-550	真空ポンプ	¥2,594,000	¥4,718,000	
54994-050	真空ポンプ	¥1,927,100	¥3,505,000	
54994-052	真空ポンプ	¥2,717,000	¥4,942,000	
54994-074	交換エレメント	¥90,600	¥177,000	
54994-076	交換エレメント	¥67,600	¥125,000	
54994-084	オイルミトリネーター用ポンプアダプターキット	¥121,500	¥225,000	
54994-086	オイルミトリネーター用ポンプアダプターキット	¥47,800	¥89,000	
54994-092	Acid Neutralization Trap	¥244,900	¥454,000	
54994-094	オイルミトリネーター	¥448,000	¥831,000	
54994-110	真空ポンプ(オイルフリー)	¥511,300	¥930,000	
54996-061	ポンプアダプター(真空ポンプ用)	¥13,200	¥25,000	
54996-082	ポンプアダプター(真空ポンプ用)	¥27,000	¥52,000	
54996-141	ポンプアダプター(真空ポンプ用)	¥69,200	¥133,000	
54996-184	ポンプアダプター(真空ポンプ用) 1L 1入	¥28,000	¥54,000	
55005-007	HyFlex Lite Gloves size 7 12入	¥28,400	¥30,500	
55005-009	HyFlex Lite Gloves size 8 12入	¥28,400	¥30,500	
55007-124	オイルミトリネーター用交換エレメント	¥106,100	¥197,000	
55007-176	ポンプアダプター(トリアイス冷却式)	¥831,000	¥1,541,000	
55008-180	ポンプアダプター(真空計)	¥557,600	¥1,034,000	
55008-186	ポンプアダプター(真空計)	¥94,400	¥175,000	
55009-050	ポンプアダプター(バキュームポンプ用)	¥40,600	¥75,000	
55009-052	ポンプアダプター(バキュームポンプ用)	¥79,500	¥147,000	
55009-054	ポンプアダプター(バキュームポンプ用)	¥82,700	¥153,000	
55009-097	ポンプアダプター(ベンチング・リンク)	¥13,200	¥24,000	
55009-098	ポンプアダプター(バキュームポンプ用)	¥150,500	¥279,000	
55009-101	ポンプアダプター(ベンチング・リンク)	¥20,600	¥38,000	
55009-104	ポンプアダプター(ポンプ)	¥20,200	¥37,000	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
55009-109	ホップアダプター(クランプ)	¥19,500	¥36,000	
55009-115	ホップアダプター(クランプ)	¥13,000	¥24,000	
55009-116	ホップアダプター(セパレート・リンク)	¥19,300	¥36,000	
55009-118	RubberHoseAdapters Welch NW25	¥28,200	¥52,000	
55009-128	ホップアダプター(ドロップ用エレメント)	¥59,600	¥111,000	
55009-131	ホップアダプター(ホース口)	¥25,200	¥47,000	
55009-149	ホップアダプター(ス・アダプター)	¥46,400	¥86,000	
55009-152	ホップアダプター(ホース口)	¥22,700	¥42,000	
55009-153	ホップアダプター(バキュームチューブ&クランプ)	¥25,600	¥47,000	
55009-172	ホップアダプター(ス・アダプター)	¥42,100	¥78,000	
55009-176	ホップアダプター(ス・アダプター)	¥46,400	¥86,000	
55009-191	ホップアダプター(ドロップ用エレメント)	¥64,400	¥119,000	
55009-238	ホップアダプター(リテーナ)	¥55,900	¥104,000	
55009-239	ホップアダプター(ホース口)	¥18,200	¥34,000	
55009-276	ホップアダプター(ドロップ用エレメント)	¥59,600	¥111,000	
55009-280	ホップアダプター(バキュームチューブ&クランプ)	¥27,100	¥50,000	
55009-282	ホップアダプター(バキュームチューブ&クランプ)	¥35,200	¥65,000	
55100-020	排液回収タンク 7.6L 1	¥367,300	¥426,900	
55407-014	分光光度計	¥6,900	¥14,000	
55407-020	分光光度計	¥5,100	¥11,000	
55407-028	分光光度計	¥5,200	¥11,000	
55407-030	分光光度計	¥28,800	¥59,000	
55407-034	分光光度計	¥17,000	¥33,000	
55407-040	分光光度計	¥6,000	¥12,000	
55407-042	分光光度計	¥2,900	¥6,000	
55407-044	分光光度計	¥9,400	¥19,000	
55407-046	分光光度計	¥6,000	¥12,000	
55407-048	分光光度計	¥10,700	¥22,000	
55407-050	分光光度計	¥28,200	¥57,000	
55411-050	VWR Razor Blades 替刃 0.22mm 100枚	¥24,200	¥27,600	
55411-055	VWR Razor Blades 替刃 0.3mm 100枚	¥21,100	¥24,200	
55533-028	耳栓 Laser-Lite1	¥342,000	¥384,700	
55533-030	耳栓 Laser-Lite 200 th ア	¥34,500	¥38,800	
55615-664	耳栓 E-A-R1	¥98,700	¥105,800	
55709-220	凍結保存容器キャニスター型(CCシース)	¥1,029,700	¥1,873,000	
55709-222	凍結保存容器キャニスター型(CCシース)	¥1,041,900	¥1,895,000	
55709-224	凍結保存容器キャニスター型(CCシース)	¥1,267,200	¥2,305,000	
55709-226	液体窒素用テトラ瓶(Vシース)	¥987,800	¥1,797,000	
55709-228	液体窒素用テトラ瓶(Vシース)	¥914,700	¥1,664,000	
55709-230	液体窒素用テトラ瓶(Vシース)	¥2,618,900	¥4,764,000	
55709-234	液体窒素用テトラ瓶(Lシース)	¥514,000	¥935,000	
55709-236	液体窒素用テトラ瓶(Lシース)	¥524,700	¥954,000	
55709-238	液体窒素用テトラ瓶(Lシース)	¥588,000	¥1,069,000	
55709-240	液体窒素用テトラ瓶(Lシース)	¥698,900	¥1,271,000	
55709-242	液体窒素用テトラ瓶(Lシース)	¥871,500	¥1,585,000	
55709-270	液体窒素用テトラ瓶(Lシース)	¥254,400	¥275,000	
55709-272	液体窒素用テトラ瓶(Lシース)	¥252,700	¥273,200	
55709-274	凍結保存容器キャニスター型(CCシース)	¥251,500	¥271,900	
55709-280	液体窒素用テトラ瓶(Lシース)	¥549,300	¥999,000	
55709-282	液体窒素用テトラ瓶(Lシース)	¥549,300	¥999,000	
55709-284	液体窒素用テトラ瓶(Lシース)	¥549,300	¥999,000	
55709-286	液体窒素用テトラ瓶(Lシース)	¥549,300	¥999,000	
55709-288	液体窒素用テトラ瓶(Lシース)	¥176,700	¥191,100	
55709-290	液体窒素用テトラ瓶(Lシース)	¥176,700	¥191,100	
55709-292	液体窒素用テトラ瓶(Lシース)	¥176,700	¥191,100	
55709-294	液体窒素用テトラ瓶(Lシース)	¥178,400	¥192,900	
55709-296	液体窒素用テトラ瓶(Lシース)	¥175,600	¥189,800	
55709-298	液体窒素用テトラ瓶(Vシース)	¥371,500	¥401,600	
55709-300	液体窒素用テトラ瓶(Vシース)	¥358,600	¥652,000	
55710-254	Nunc MegaMAX-100 Box 4個入	¥65,600	¥70,900	
55850-700	VWR Gas Pressure Regulators	¥227,700	¥414,000	
56219-004	防塵マスク 2200 N95 20個	¥24,600	¥26,400	
56219-020	防塵マスク 2300 N95	¥29,900	¥32,100	
56222-100	7800シース FullFacepieceRespiratorsL	¥387,600	¥445,300	
56222-362	防塵マスク 8210 N95 20個	¥27,800	¥29,800	
56222-944	防毒マスク(7700用マスク)	¥323,300	¥349,400	
56222-964	防毒マスク(交換用カートリッジ)	¥629,500	¥680,200	
56223-152	Glove Box Gloves(ハイパー)SZ8.5 15ミル	¥211,900	¥234,600	
56223-210	袖口製グローブ用ハイパー SZ8.5	¥181,700	¥201,200	
56223-250	防毒マスク(交換用カートリッジ)2x18	¥181,600	¥196,100	
56250-024	ワサ台 10to100mm/80mm/32mm	¥5,300	¥5,900	
56250-046	ワサ台 200to500mm/110mm/60mm	¥11,800	¥13,100	
56250-068	ワサ台 1000to3000mm/140mm/90mm	¥15,300	¥17,100	
56250-080	ワサ台 5000mm/170mm/120mm	¥19,300	¥21,600	
56250-104	ワサ台 12000mm/210mm/150mm	¥32,500	¥36,300	
56255-006	ワサ台(袖口製)	¥13,700	¥15,000	
56330-000	ルテイルホップ 気体用	¥14,700	¥15,900	
56335-003	Rubber Bulb Set(ラテックス)200ml	¥35,100	¥36,800	
56340-003	気体採取バッグ	¥14,000	¥15,000	
56341-243	コップ 10mL	¥9,000	¥10,600	
56341-287	コップ 15mL	¥46,200	¥49,900	
56341-367	コップ 35mL 12入	¥51,200	¥55,300	
56341-403	コップ 60mL 12入	¥46,100	¥49,800	
56354-008	マット(作業用)マット	¥61,800	¥88,000	
56355-000	マット(作業用)マット	¥61,900	¥66,300	
56358-009	マット(作業用)マット	¥74,400	¥79,600	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
56359-001	マット(作業用)マット	¥74,400	¥79,600	
56360-120	DURA-STAT SoftWork Surfaces60x90cm	¥48,400	¥51,900	
56360-122	DURA-STAT SoftWork Surfaces	¥61,500	¥65,900	
56360-124	DURA-STAT SoftWork Surfaces	¥1,298,700	¥1,390,400	
56360-126	DURA-STAT SoftWork Surfaces60x90cm	¥46,200	¥49,500	
56360-130	DURA-STAT SoftWork Surfaces	¥567,400	¥607,500	
56510-014	定規(金属製)	¥13,700	¥15,100	
56606-595	ハイルグート 標識板	¥37,400	¥40,100	
56606-613	危険標識板	¥8,700	¥9,300	
56606-619	危険標識板	¥9,100	¥9,700	
56606-622	危険標識板	¥7,800	¥8,400	
56606-627	危険標識板	¥6,600	¥7,100	
56606-629	危険標識板	¥8,700	¥9,300	
56606-630	危険標識板	¥8,500	¥9,100	
56606-660	バル(危険物質輸送用/DOT規格)	¥53,900	¥57,700	
56606-662	バル(危険物質輸送用/DOT規格)	¥51,400	¥55,000	
56606-663	バル(危険物質輸送用/DOT規格)	¥51,400	¥55,000	
56606-669	バル(危険物質輸送用/DOT規格)	¥51,400	¥55,000	
56606-676	危険標識板	¥12,300	¥13,100	
56606-684	危険標識板	¥8,800	¥9,400	
56606-685	危険標識板	¥8,600	¥9,200	
56606-686	危険標識板	¥12,300	¥13,200	
56606-688	危険標識板	¥6,800	¥7,300	
56606-701	バル(危険警告用/EPA/DOT規格)	¥19,800	¥22,200	
56606-703	バル(危険警告用/EPA/DOT規格)	¥72,400	¥81,000	
56606-705	バル(危険警告用/EPA/DOT規格)	¥39,100	¥41,900	
56606-706	バル(危険警告用/EPA/DOT規格)	¥44,600	¥47,700	
56606-714	危険標識板	¥7,600	¥8,100	
56606-715	危険標識板	¥9,100	¥9,700	
56606-716	危険標識板	¥8,900	¥9,500	
56606-718	標識板	¥8,800	¥9,400	
56606-726	標識板	¥5,600	¥5,900	
56606-729	標識板	¥7,800	¥8,400	
56606-730	標識板	¥5,600	¥6,200	
56606-732	標識板	¥7,600	¥8,100	
56606-741	標識板	¥7,600	¥8,100	
56606-750	バル(危険警告用) 5.1x5.1cm	¥14,400	¥16,100	
56606-751	バル(危険警告用)	¥33,700	¥37,700	
56606-756	危険物識別標識板	¥46,700	¥50,000	
56606-757	危険物識別標識板	¥29,400	¥31,400	
56606-763	警告サイン(床置タイプ)	¥28,900	¥30,900	
56606-764	標識板	¥6,800	¥7,200	
56606-765	標識板	¥9,000	¥9,600	
56606-766	標識板	¥8,800	¥9,400	
56606-767	標識板	¥12,100	¥13,000	
56606-780	ハイルグート 標識板	¥7,800	¥8,400	
56606-783	ハイルグート 標識板	¥7,900	¥8,400	
56606-785	標識板	¥5,600	¥5,900	
56606-789	標識板	¥8,800	¥9,400	
56606-822	危険標識板	¥8,700	¥9,300	
56606-823	危険標識板	¥8,500	¥9,100	
56606-847	危険標識板	¥8,800	¥9,400	
56606-850	危険標識板	¥6,600	¥7,100	
56606-852	危険標識板	¥8,800	¥9,400	
56606-853	危険標識板	¥12,500	¥13,400	
56606-857	危険標識板	¥8,800	¥9,400	
56606-858	危険標識板	¥12,100	¥13,000	
56606-860	危険標識板	¥8,800	¥9,400	
56606-863	危険標識板	¥12,500	¥13,300	
56606-866	危険標識板	¥8,700	¥9,300	
56606-867	危険標識板	¥8,500	¥9,100	
56606-868	危険標識板	¥12,100	¥13,000	
56606-870	危険標識板	¥6,500	¥7,000	
56606-872	危険標識板	¥8,500	¥9,100	
56606-873	危険標識板	¥12,500	¥13,300	
56607-054	標識板	¥11,000	¥11,800	
56607-520	危険標識板	¥6,400	¥6,900	
56607-521	危険標識板	¥8,700	¥9,300	
56607-522	危険標識板	¥8,500	¥9,100	
56607-523	危険標識板	¥12,100	¥13,000	
56607-526	危険標識板	¥8,700	¥9,300	
56607-527	危険標識板	¥8,500	¥9,100	
56607-528	危険標識板	¥12,100	¥13,000	
56607-530	危険標識板	¥6,400	¥6,900	
56607-531	危険標識板	¥8,700	¥9,300	
56607-532	危険標識板	¥8,500	¥9,100	
56607-533	危険標識板	¥12,100	¥13,000	
56607-535	危険標識板	¥6,400	¥6,900	
56607-536	危険標識板	¥8,700	¥9,300	
56607-537	危険標識板	¥8,500	¥9,100	
56607-538	危険標識板	¥12,100	¥13,000	
56607-540	危険標識板	¥6,600	¥7,100	
56607-541	危険標識板	¥9,000	¥9,600	
56607-542	危険標識板	¥8,500	¥9,100	
56607-543	危険標識板	¥12,100	¥13,000	
56607-545	危険標識板	¥6,400	¥6,900	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
56607-546	危険標識板	¥8,700	¥9,300	
56607-547	危険標識板	¥8,500	¥9,100	
56607-548	危険標識板	¥12,100	¥13,000	
56607-653	標識板	¥6,400	¥6,900	
56607-654	標識板	¥8,800	¥9,400	
56607-655	標識板	¥8,500	¥9,100	
56607-656	標識板	¥12,100	¥13,000	
56607-657	標識板	¥6,400	¥6,900	
56607-658	標識板	¥9,000	¥9,600	
56607-659	標識板	¥8,500	¥9,100	
56607-661	標識板	¥6,400	¥6,900	
56607-662	標識板	¥8,800	¥9,400	
56607-687	標識板	¥5,400	¥5,700	
56607-696	標識板	¥5,400	¥5,800	
56607-697	標識板	¥7,600	¥8,100	
56607-699	標識板	¥5,400	¥5,800	
56607-702	標識板	¥5,400	¥5,800	
56607-703	標識板	¥7,600	¥8,100	
56607-792	パネル(危険警告用/EPA/DOT規格)	¥67,000	¥73,900	
56607-795	危険標識板	¥8,600	¥9,200	
56607-800	危険標識板	¥6,700	¥7,200	
56607-805	危険標識板	¥6,700	¥7,100	
56607-807	危険標識板	¥8,600	¥9,200	
56607-808	危険標識板	¥12,300	¥13,200	
56607-812	危険標識板	¥8,600	¥9,200	
56607-813	危険標識板	¥12,300	¥13,200	
56607-815	危険標識板	¥6,800	¥7,300	
56607-817	危険標識板	¥8,600	¥9,200	
56607-837	パネル(危険警告用/書き込み式)	¥28,000	¥31,300	
56607-838	パネル(危険警告用/書き込み式)	¥21,500	¥23,000	
56607-839	パネル(危険警告用/書き込み式)	¥28,100	¥30,100	
56607-840	パネル(危険警告用/書き込み式)500	¥14,600	¥15,700	
56607-841	パネル(危険警告用/書き込み式)	¥28,100	¥30,900	
56607-875	危険物識別標識板	¥14,900	¥15,900	
56607-876	危険物識別標識板	¥14,900	¥15,900	
56607-877	危険物識別標識板	¥15,300	¥16,400	
56607-878	危険物識別標識板	¥14,900	¥15,900	
56607-879	危険物識別標識板	¥15,300	¥16,400	
56607-893	危険物識別標識板	¥14,900	¥15,900	
56607-897	危険物識別標識板	¥14,900	¥15,900	
56607-901	危険物識別標識板	¥15,300	¥16,400	
56607-905	危険物識別標識板	¥14,900	¥15,900	
56607-909	危険物識別標識板	¥14,900	¥15,900	
56607-913	危険物識別標識板	¥14,900	¥15,900	
56608-864	ガラス瓶運搬容器(安全キャリ)	¥55,000	¥58,900	
56608-886	ガラス瓶運搬容器(安全キャリ)	¥46,600	¥49,900	
56608-900	ガラス瓶運搬容器(安全キャリ)	¥54,700	¥58,500	
56608-970	ガラス瓶運搬用キャリ	¥34,800	¥53,200	
56608-972	ガラス瓶運搬用キャリ	¥21,400	¥22,500	
56608-975	ガラス瓶運搬用キャリ	¥35,200	¥53,800	
56608-980	ガラス瓶運搬用キャリ	¥32,400	¥54,400	
56608-982	ガラス瓶運搬用キャリ	¥21,400	¥35,900	
56609-187	廃液保管ボトル(安全缶:樹脂製)	¥208,000	¥297,800	
56609-190	液体吸収材	¥233,400	¥249,900	
56609-202	中和処理剤	¥62,700	¥67,100	
56609-204	中和処理剤	¥210,900	¥225,700	
56609-205	Mercury Indicator Powder 250g	¥85,500	¥91,600	
56609-223	液体吸収材 10入	¥58,500	¥62,600	
56609-470	薬災防止缶(安全缶)	¥47,700	¥51,000	
56609-506	薬災防止缶(安全缶)	¥52,900	¥56,600	
56609-517	薬災防止缶(安全缶)	¥28,700	¥29,200	
56609-528	薬災防止缶(安全缶)	¥28,200	¥27,600	
56609-550	薬災防止缶(安全缶)	¥30,700	¥31,300	
56609-561	薬災防止缶(安全缶)	¥61,000	¥65,300	
56609-572	薬災防止缶(安全缶)	¥33,400	¥35,700	
56609-586	薬災防止缶(安全缶)	¥311,000	¥326,400	
56609-600	薬災防止缶(安全缶) 3.8	¥53,300	¥56,000	
56609-749	フェイスシールド 946ml 8.9x23.2cm	¥18,100	¥19,500	
56609-751	Stainless Steel Flame Arrester 1	¥121,800	¥145,900	
56609-765	薬災防止缶(安全缶)	¥116,800	¥125,100	
56609-801	薬災防止缶(安全缶)	¥58,900	¥63,000	
56609-812	薬災防止缶(安全缶)	¥183,200	¥196,000	
56609-834	薬災防止缶(安全缶)	¥193,600	¥207,100	
56609-878	薬災防止缶(安全缶)Coated Metal	¥135,000	¥144,500	
56609-889	薬災防止缶(安全缶) 18.9L	¥380,300	¥685,000	
56609-970	薬災防止缶(安全缶)	¥104,200	¥111,500	
56610-100	ボトルキャット	¥30,200	¥34,400	
56610-102	フェイスシールド 473ml 7.6x17.9cm	¥27,500	¥35,100	
56610-130	廃液保管ボトル(安全缶)	¥570,300	¥610,500	
56610-141	廃液保管ボトル(安全缶)	¥266,700	¥285,300	
56610-143	ドラム缶 保存ケース 30GL	¥35,400	¥26,500	
56610-165	Overpack Dolly 1個	¥501,800	¥510,600	
56610-174	廃液保管ボトル(安全缶:樹脂製)	¥228,500	¥272,600	
56610-175	OILY WASTE CAN POLY RED 10GL	¥89,800	¥96,200	
56610-177	耐火ドラム缶箱(6缶入)	¥75,600	¥80,900	
56610-181	耐火ドラム缶箱(14缶入)	¥99,100	¥106,100	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
56610-190	試薬ボトル(高密度PE製)	¥78,900	¥82,900	
56610-192	試薬ボトル(高密度PE製)	¥148,500	¥166,000	
56610-196	試薬ボトル(高密度PE製)	¥223,100	¥234,200	
56610-212	トランスフュージョン(連結タイプ)	¥257,300	¥275,500	
56610-216	耐火ゴミ箱(6リットル)	¥76,100	¥81,500	
56610-218	耐火ゴミ箱(10リットル)	¥89,400	¥95,700	
56610-220	耐火ゴミ箱(14リットル)	¥99,600	¥106,600	
56610-226	耐火ゴミ箱(10リットル)	¥91,800	¥98,300	
56610-543	可燃物貯蔵容器(自己消化式)	¥136,300	¥180,800	
56610-544	可燃物貯蔵容器(自己消化式)	¥69,800	¥84,600	
56610-545	可燃物貯蔵容器(自己消化式)	¥158,100	¥165,900	
56610-546	可燃物貯蔵容器(自己消化式)	¥162,800	¥142,900	
56610-547	可燃物貯蔵容器(自己消化式)	¥247,400	¥390,800	
56610-548	可燃物貯蔵容器(自己消化式)	¥50,300	¥52,800	
56610-549	可燃物貯蔵容器(自己消化式)	¥77,700	¥81,600	
56610-551	可燃物貯蔵容器(自己消化式)	¥61,700	¥64,700	
56610-552	可燃物貯蔵容器(自己消化式)	¥60,500	¥63,500	
56610-553	可燃物貯蔵容器(自己消化式)	¥162,800	¥170,900	
56610-652	耳栓 Laser-Lite1	¥283,200	¥318,600	
56610-663	耳栓 Pura-Fit 6800 200個	¥31,300	¥33,500	
56610-673	耳栓(付) Pura-Fit 6900 100個	¥33,800	¥36,200	
56611-075	Eyewash bottle 0.9L	¥17,500	¥18,700	
56611-083	Eyewash bottle holder	¥36,600	¥39,700	
56611-199	SCIENCEWARE Fire Blanket 1	¥152,800	¥168,400	
56611-202	消防カバー	¥392,000	¥873,000	
56611-205	消防カバー	¥114,300	¥206,000	
56613-340	Mid-Temperatureクリーンルームグローブ	¥324,900	¥347,800	
56613-472	インターグローブ クリアフィット(クリーンルーム用)M	¥328,900	¥364,300	
56613-474	インターグローブ クリアフィット(クリーンルーム用)XL	¥418,600	¥505,100	
56614-106	ポンプ 16oz ESD Safe Pump Bottle	¥156,000	¥167,000	
56614-284	インターグローブ クリアフィット(クリーンルーム用)S	¥350,500	¥388,200	
56614-288	インターグローブ クリアフィット(クリーンルーム用)L	¥328,900	¥364,300	
56616-018	実験台保護シート(吸水性)	¥63,800	¥70,000	
56616-020	実験台保護シート(吸水性)	¥60,400	¥66,200	
56616-022	実験台保護シート(吸水性)	¥64,200	¥70,400	
56616-030	実験台保護シート(吸水性)	¥48,300	¥53,000	
56616-031	実験台保護シート(吸水性)	¥37,700	¥41,400	
56616-032	実験台保護シート(吸水性)	¥51,700	¥56,700	
56616-079	中和処理剤	¥230,500	¥246,700	
56616-126	中和処理剤	¥187,800	¥201,100	
56616-322	液体吸収材	¥81,200	¥87,000	
56616-324	液体吸収材	¥82,000	¥87,800	
56616-342	液体吸収材	¥105,700	¥113,200	
56616-344	液体吸収材	¥97,500	¥104,400	
56616-346	液体吸収材	¥170,100	¥182,100	
56616-380	液体吸収材	¥50,000	¥53,500	
56616-388	液体吸収材	¥208,800	¥223,600	
56616-390	液体吸収材	¥176,900	¥189,400	
56616-391	液体吸収材	¥150,300	¥160,900	
56616-395	液体吸収材	¥76,100	¥81,500	
56616-397	液体吸収材	¥387,700	¥415,100	
56616-398	液体吸収材	¥263,500	¥282,100	
56616-404	中和処理剤	¥188,000	¥201,300	
56617-000	VWR UNDRPAD DELX23X36 CS150 1	¥63,300	¥75,600	
56617-006	実験台保護シート(吸水性)	¥54,200	¥59,500	
56617-008	実験台保護シート(吸水性)	¥36,400	¥40,700	
56617-014	実験台保護シート(吸水性) 300枚	¥40,400	¥44,300	
56617-016	実験台保護シート(吸水性)	¥40,700	¥45,600	
56617-018	実験台保護シート(吸水性) 10	¥3,000	¥3,200	
56617-329	Double Drum Work Station 59kg	¥1,048,800	¥1,889,000	
56617-330	OILY WASTE CAN 22.7L	¥74,600	¥79,900	
56617-331	廃棄用容器	¥33,300	¥24,900	
56617-335	OILY WASTE CANS 53L	¥121,800	¥130,400	
56617-343	トランスフュージョン	¥102,100	¥109,300	
56617-509	グライファースト用追加ティンラック 20リットル	¥23,400	¥25,100	
56617-513	フューズボックス(丸型 TPX製) 1	¥46,800	¥50,100	
56617-515	フューズボックス(丸型 TPX製) 1	¥50,300	¥53,800	
56617-801	ガラス器具用廃棄容器(床型)	¥52,700	¥60,800	
56617-804	ガラス器具用廃棄容器(卓上型) 6個	¥36,600	¥52,700	
56617-807	ゴミ箱(ハイバート焼却用・タンポール)	¥50,500	¥62,600	
56617-810	ゴミ箱(ハイバート焼却用・タンポール)	¥39,200	¥51,400	
56618-538	トランスフュージョン	¥107,400	¥115,000	
56618-544	トランスフュージョン(連結タイプ)	¥256,900	¥275,000	
56618-546	トランスフュージョン(連結タイプ)	¥460,300	¥492,700	
56618-576	薬災防止缶(安全缶)	¥7,900	¥8,400	
56618-596	薬災防止缶(安全缶)	¥38,100	¥38,700	
56618-598	薬災防止缶(安全缶)	¥38,100	¥38,700	
56621-928	防塵マスク 8233 N100 20個入	¥203,100	¥217,400	
56621-932	防塵マスク 8247 R95 20個	¥78,700	¥90,900	
56621-934	防塵マスク 8271 P95 10入	¥46,800	¥50,100	
56621-936	防塵マスク 8511 N95 10入	¥25,500	¥27,300	
56700-174	スリットフィルター 25x75mm	¥251,500	¥271,800	
56766-364	バッグ(再密閉可能ハイバートバッグ)	¥64,100	¥70,300	
56766-370	再密閉可能ハイバートバッグ 15.2x22.9	¥60,800	¥66,700	
56766-372	バッグ(再密閉可能ハイバートバッグ)	¥73,400	¥80,500	
56766-374	バッグ(再密閉可能ハイバートバッグ)	¥149,300	¥163,700	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
56766-390	再密閉可能な「バグ」-「ト」パック 15.2x22.9	¥55,800	¥61,200	
56775-083	メ「バ」パック (サン) 郵送用封筒	¥146,500	¥166,100	
56920-020	万能スcoop (スクープ) ポリプロピレン製	¥8,300	¥8,800	
56920-025	万能スcoop (スクープ) ポリプロピレン製	¥20,700	¥24,700	
56920-036	万能スcoop (スクープ) ポリプロピレン製	¥32,400	¥34,400	
56920-047	万能スcoop (スクープ) ポリプロピレン製	¥30,700	¥32,500	
56921-006	ステンレススcoop	¥167,000	¥196,100	
56924-478	サンプリンクナイフ	¥219,900	¥233,100	
56924-503	サンプリンクナイフ (PS製/滅菌済み個別包装)	¥90,900	¥101,200	
56924-536	サンプリンクナイフ (PS製/滅菌済み個別包装)	¥116,100	¥129,200	
56924-569	サンプリンクナイフ (PS製/滅菌済み個別包装)	¥174,000	¥184,400	
56924-576	Scoop Sampling System 125ml 100個入	¥202,900	¥219,400	
56924-809	ス「バ」パック (ポリプロピレン製/滅菌済み個別包装)	¥103,600	¥109,800	
56925-109	サンプリンクナイフ (PS製/滅菌済み個別包装)	¥158,400	¥168,000	
56925-131	サンプリンクナイフ (PS製/滅菌済み個別包装)	¥160,400	¥170,000	
56925-164	サンプリンクナイフ (PS製/滅菌済み個別包装)	¥144,700	¥153,400	
56925-197	サンプリンクナイフ (PS製/滅菌済み個別包装)	¥143,200	¥161,800	
57007-000	シーカー	¥4,461,400	¥9,021,000	
57008-002	シーカー	¥426,800	¥788,000	
57008-024	シーカー Crosswise Bar Clamp	¥63,900	¥118,000	
57008-046	シーカー	¥64,400	¥119,000	
57018-775	シーカー	¥19,400	¥37,000	
57018-786	Stainless Steel Clamp	¥18,100	¥32,000	
57018-797	シーカー	¥16,300	¥29,000	
57018-800	シーカー	¥16,700	¥29,000	
57018-811	シーカー	¥16,500	¥29,000	
57018-822	シーカー	¥20,800	¥37,000	
57018-844	シーカー	¥65,300	¥118,000	
57250-020	ハサミ Shears 20.3cm	¥17,500	¥19,000	
57250-043	Shears	¥15,600	¥17,000	
57252-150	Scissors with Footplate	¥7,600	¥8,300	
57320-486	ふるい (黄銅製, 8インチ=直径20.3cm)	¥97,800	¥108,800	
57320-500	ふるい (黄銅製, 8インチ=直径20.3cm)	¥95,900	¥106,700	
57320-522	ふるい (黄銅製, 8インチ=直径20.3cm)	¥78,900	¥87,800	
57320-544	ふるい (黄銅製, 8インチ=直径20.3cm)	¥69,600	¥77,400	
57320-566	ふるい (黄銅製, 8インチ=直径20.3cm)	¥81,800	¥91,100	
57320-588	ふるい (黄銅製, 8インチ=直径20.3cm)	¥82,500	¥91,800	
57320-602	ふるい (黄銅製, 8インチ=直径20.3cm)	¥80,900	¥90,000	
57320-624	ふるい (黄銅製, 8インチ=直径20.3cm)	¥70,300	¥78,200	
57320-646	ふるい (黄銅製, 8インチ=直径20.3cm)	¥83,300	¥92,700	
57320-668	ふるい (黄銅製, 8インチ=直径20.3cm)	¥80,700	¥89,700	
57320-680	ふるい (黄銅製, 8インチ=直径20.3cm)	¥80,200	¥89,200	
57320-726	ふるい (黄銅製, 8インチ=直径20.3cm)	¥81,000	¥90,100	
57320-748	ふるい (黄銅製, 8インチ=直径20.3cm)	¥70,400	¥78,400	
57320-760	ふるい (黄銅製, 8インチ=直径20.3cm)	¥80,100	¥89,200	
57320-781	ふるい (黄銅製, 8インチ=直径20.3cm)	¥70,800	¥78,700	
57320-806	ふるい (黄銅製, 8インチ=直径20.3cm)	¥81,000	¥90,200	
57320-828	ふるい (黄銅製, 8インチ=直径20.3cm)	¥70,800	¥78,800	
57320-840	ふるい (黄銅製, 8インチ=直径20.3cm)	¥78,800	¥87,700	
57320-861	ふるい (黄銅製, 8インチ=直径20.3cm)	¥84,400	¥93,900	
57320-883	ふるい (黄銅製, 8インチ=直径20.3cm)	¥90,900	¥101,100	
57320-908	ふるい (黄銅製, 8インチ=直径20.3cm)	¥107,400	¥119,500	
57320-920	ふるい (黄銅製, 8インチ=直径20.3cm)	¥127,200	¥141,600	
57320-941	ふるい (黄銅製, 8インチ=直径20.3cm)	¥159,500	¥177,400	
57321-503	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥67,000	¥74,600	
57321-525	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥65,800	¥73,200	
57321-547	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥65,800	¥73,200	
57321-569	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥69,200	¥77,000	
57321-580	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥84,800	¥94,300	
57321-605	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥65,800	¥73,200	
57321-649	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥84,800	¥94,300	
57321-660	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥65,800	¥73,200	
57321-682	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥86,600	¥96,400	
57321-729	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥69,200	¥77,000	
57321-740	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥79,600	¥88,600	
57321-762	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥65,800	¥73,200	
57321-784	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥80,700	¥89,800	
57321-809	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥65,800	¥73,200	
57321-842	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥69,200	¥77,000	
57321-864	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥98,100	¥109,100	
57322-007	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥82,700	¥92,000	
57322-029	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥69,200	¥77,000	
57322-040	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥66,500	¥74,000	
57322-062	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥69,200	¥77,000	
57322-084	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥76,700	¥85,300	
57322-109	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥74,100	¥82,500	
57322-120	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥65,800	¥73,200	
57322-142	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥65,800	¥73,200	
57322-164	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥66,300	¥73,700	
57322-186	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥69,200	¥77,000	
57322-200	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥65,800	¥73,200	
57322-222	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥69,200	¥77,000	
57322-244	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥66,900	¥74,500	
57322-266	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥65,800	¥73,200	
57322-288	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥66,300	¥73,700	
57322-302	ふるい (ステンレス鋼製, 8インチ=直径20.3cm)	¥69,200	¥77,000	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
57322-324	ふるい(ステン鋼製、8インチ=直径20.3cm)	¥66,300	¥73,700	
57322-346	ふるい(ステン鋼製、8インチ=直径20.3cm)	¥70,300	¥78,200	
57322-368	ふるい(ステン鋼製、8インチ=直径20.3cm)	¥66,300	¥73,700	
57322-380	ふるい	¥59,500	¥66,200	
57322-404	ふるい(ステン鋼製、8インチ=直径20.3cm)	¥66,900	¥74,500	
57322-426	ふるい(ステン鋼製、8インチ=直径20.3cm)	¥67,600	¥75,200	
57322-448	ふるい(ステン鋼製、8インチ=直径20.3cm)	¥72,900	¥81,100	
57322-460	ふるい(ステン鋼製、8インチ=直径20.3cm)	¥83,300	¥92,700	
57322-481	ふるい(ステン鋼製、8インチ=直径20.3cm)	¥89,300	¥99,300	
57322-506	ふるい(ステン鋼製、8インチ=直径20.3cm)	¥105,800	¥117,700	
57322-528	ふるい(ステン鋼製、8インチ=直径20.3cm)	¥122,200	¥136,000	
57322-540	ふるい(ステン鋼製、8インチ=直径20.3cm)	¥152,100	¥169,200	
57322-561	ふるい(ステン鋼製、8インチ=直径20.3cm)	¥176,000	¥195,800	
57324-002	ふるい(8インチ=直径20.3cm)	¥78,900	¥87,800	
57324-024	ふるい(8インチ=直径20.3cm)	¥79,100	¥88,000	
57324-046	ふるい(8インチ=直径20.3cm)	¥74,600	¥83,000	
57324-068	ふるい(8インチ=直径20.3cm)	¥74,100	¥82,500	
57324-080	ふるい(8インチ=直径20.3cm)	¥75,900	¥84,400	
57324-104	ふるい(8インチ=直径20.3cm)	¥75,700	¥84,200	
57324-126	ふるい(8インチ=直径20.3cm)	¥76,200	¥84,700	
57324-148	ふるい(8インチ=直径20.3cm)	¥79,100	¥88,000	
57324-160	ふるい(8インチ=直径20.3cm)	¥75,700	¥84,200	
57324-181	ふるい(8インチ=直径20.3cm)	¥79,300	¥88,300	
57324-206	ふるい(8インチ=直径20.3cm)	¥76,200	¥84,700	
57324-228	ふるい(8インチ=直径20.3cm)	¥79,300	¥88,300	
57324-240	ふるい(8インチ=直径20.3cm)	¥75,900	¥84,400	
57324-261	ふるい(8インチ=直径20.3cm)	¥79,300	¥88,300	
57324-272	ふるい(8インチ=直径20.3cm)	¥75,700	¥84,200	
57324-283	ふるい(8インチ=直径20.3cm)	¥79,100	¥88,000	
57324-308	ふるい(8インチ=直径20.3cm)	¥79,700	¥88,700	
57324-320	ふるい(8インチ=直径20.3cm)	¥79,500	¥88,400	
57324-363	ふるい(8インチ=直径20.3cm)	¥75,400	¥83,800	
57324-385	ふるい(8インチ=直径20.3cm)	¥75,400	¥83,800	
57324-400	ふるい(8インチ=直径20.3cm)	¥78,200	¥87,000	
57324-421	ふるい(8インチ=直径20.3cm)	¥79,700	¥88,600	
57324-443	ふるい(8インチ=直径20.3cm)	¥87,300	¥97,200	
57324-465	ふるい(8インチ=直径20.3cm)	¥99,500	¥110,700	
57324-487	ふるい(8インチ=直径20.3cm)	¥107,400	¥119,500	
57324-501	ふるい(8インチ=直径20.3cm)	¥122,600	¥136,400	
57324-523	ふるい(8インチ=直径20.3cm)	¥135,100	¥150,300	
57324-545	ふるい(8インチ=直径20.3cm)	¥154,600	¥172,000	
57324-567	ふるい(8インチ=直径20.3cm)	¥193,500	¥215,300	
57331-405	ふるい(黄銅製、8インチ=直径20.3cm)	¥69,900	¥77,700	
57331-427	ふるい(黄銅製、8インチ=直径20.3cm)	¥69,200	¥77,000	
57331-449	ふるい(黄銅製、8インチ=直径20.3cm)	¥69,900	¥77,700	
57331-482	ふるい(黄銅製、8インチ=直径20.3cm)	¥69,900	¥77,700	
57331-507	ふるい(黄銅製、8インチ=直径20.3cm)	¥69,200	¥77,000	
57331-529	ふるい(黄銅製、8インチ=直径20.3cm)	¥70,300	¥78,200	
57331-540	ふるい(黄銅製、8インチ=直径20.3cm)	¥69,200	¥77,000	
57331-562	ふるい(黄銅製、8インチ=直径20.3cm)	¥70,300	¥78,200	
57331-584	ふるい(黄銅製、8インチ=直径20.3cm)	¥69,900	¥77,700	
57331-609	ふるい(黄銅製、8インチ=直径20.3cm)	¥69,600	¥77,400	
57331-620	ふるい(黄銅製、8インチ=直径20.3cm)	¥69,600	¥77,400	
57331-642	ふるい(黄銅製、8インチ=直径20.3cm)	¥69,600	¥77,400	
57331-664	ふるい(黄銅製、8インチ=直径20.3cm)	¥70,300	¥78,200	
57331-700	ふるい(黄銅製、8インチ=直径20.3cm)	¥70,300	¥78,200	
57331-722	ふるい(黄銅製、8インチ=直径20.3cm)	¥82,400	¥91,700	
57331-744	ふるい(黄銅製、8インチ=直径20.3cm)	¥69,200	¥77,000	
57331-766	ふるい(黄銅製、8インチ=直径20.3cm)	¥78,800	¥87,700	
57331-788	ふるい(黄銅製、8インチ=直径20.3cm)	¥83,800	¥93,200	
57331-802	ふるい(黄銅製、8インチ=直径20.3cm)	¥93,400	¥103,900	
57331-846	ふるい(黄銅製、8インチ=直径20.3cm)	¥128,300	¥142,700	
57331-868	ふるい(黄銅製、8インチ=直径20.3cm)	¥158,300	¥176,100	
57332-080	ふるい(ステン鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-084	ふるい(ステン鋼製、12インチ=直径30.5cm)	¥143,300	¥159,500	
57332-086	ふるい(ステン鋼製、12インチ=直径30.5cm)	¥143,300	¥159,500	
57332-088	ふるい(ステン鋼製、12インチ=直径30.5cm)	¥112,600	¥125,300	
57332-094	ふるい(ステン鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-096	ふるい(ステン鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-097	ふるい(ステン鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-098	ふるい(ステン鋼製、12インチ=直径30.5cm)	¥113,100	¥125,800	
57332-099	ふるい(ステン鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-104	ふるい(ステン鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-106	ふるい(ステン鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-108	ふるい(ステン鋼製、12インチ=直径30.5cm)	¥107,300	¥119,300	
57332-110	ふるい(ステン鋼製、12インチ=直径30.5cm)	¥138,500	¥154,100	
57332-112	ふるい(ステン鋼製、12インチ=直径30.5cm)	¥106,500	¥118,600	
57332-114	ふるい(ステン鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-116	ふるい(ステン鋼製、12インチ=直径30.5cm)	¥113,100	¥125,800	
57332-118	ふるい(ステン鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-120	ふるい(ステン鋼製、12インチ=直径30.5cm)	¥113,100	¥125,800	
57332-122	ふるい(ステン鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-124	ふるい(ステン鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-126	ふるい(ステン鋼製、12インチ=直径30.5cm)	¥113,700	¥126,500	
57332-128	ふるい(ステン鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-130	ふるい(ステン鋼製、12インチ=直径30.5cm)	¥112,600	¥125,300	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
57332-132	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-134	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-136	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-138	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-142	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-144	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-146	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-148	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-150	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-152	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-154	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥113,100	¥125,900	
57332-158	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-160	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-162	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-164	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-166	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥106,200	¥118,200	
57332-170	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥120,400	¥134,000	
57332-172	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥120,400	¥134,000	
57332-174	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥131,800	¥146,600	
57332-176	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥131,800	¥146,600	
57332-178	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥139,700	¥155,500	
57332-180	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥143,600	¥159,700	
57332-182	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥158,700	¥176,500	
57332-184	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥157,100	¥174,800	
57332-186	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥183,500	¥204,100	
57332-190	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥226,300	¥251,800	
57332-192	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥225,100	¥250,400	
57332-194	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥237,400	¥264,100	
57332-196	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥235,400	¥261,900	
57332-198	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥276,900	¥308,100	
57332-200	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥279,200	¥310,600	
57332-202	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-204	ふるい(ステンル鋼製、12インチ=直径30.5cm)	¥112,000	¥124,600	
57332-206	ふるい(ステンル鋼製、12インチ=直径30.6cm)	¥150,700	¥167,700	
57332-208	ふるい(ステンル鋼製、12インチ=直径30.7cm)	¥157,700	¥175,400	
57332-210	ふるい(ステンル鋼製、12インチ=直径30.8cm)	¥138,500	¥154,100	
57332-212	ふるい(ステンル鋼製、12インチ=直径30.9cm)	¥142,200	¥158,300	
57332-214	ふるい(ステンル鋼製、12インチ=直径30.10cm)	¥112,000	¥124,600	
57333-965	Mini-Sieve マイクロン・セット PP	¥116,600	¥124,800	
57333-981	ふるい(テフロン製)	¥67,900	¥76,700	
57334-100	ふるい振とう器	¥64,200	¥72,400	
57334-102	ふるい振とう器	¥63,900	¥72,100	
57334-104	ふるい振とう器	¥65,200	¥73,500	
57334-106	ふるい振とう器	¥64,300	¥72,600	
57334-108	ふるい振とう器	¥60,400	¥68,100	
57334-110	ふるい振とう器	¥60,400	¥68,100	
57334-112	ふるい振とう器	¥61,300	¥69,100	
57334-114	ふるい振とう器	¥60,700	¥68,500	
57334-116	ふるい振とう器	¥60,400	¥68,100	
57334-118	ふるい振とう器	¥60,400	¥68,100	
57334-120	ふるい振とう器	¥60,400	¥68,100	
57334-122	ふるい振とう器	¥60,700	¥68,500	
57334-124	ふるい振とう器	¥60,100	¥67,800	
57334-126	ふるい振とう器	¥60,400	¥68,100	
57334-128	ふるい振とう器	¥60,400	¥68,100	
57334-130	ふるい振とう器	¥60,400	¥68,100	
57334-132	ふるい振とう器	¥64,200	¥72,400	
57334-134	ふるい振とう器	¥60,900	¥68,700	
57334-136	ふるい振とう器	¥70,300	¥79,300	
57334-138	ふるい振とう器	¥72,000	¥81,300	
57334-140	ふるい振とう器	¥81,100	¥91,500	
57334-142	ふるい振とう器	¥96,700	¥109,100	
57334-144	ふるい振とう器	¥118,100	¥133,200	
57334-146	ふるい振とう器	¥143,000	¥161,300	
57334-150	ふるい振とう器	¥63,400	¥71,500	
57334-152	ふるい振とう器	¥63,400	¥71,500	
57334-154	ふるい振とう器	¥63,400	¥71,500	
57334-156	ふるい振とう器	¥62,800	¥70,800	
57334-158	ふるい振とう器	¥59,500	¥67,100	
57334-160	ふるい振とう器	¥59,500	¥67,100	
57334-162	ふるい振とう器	¥60,000	¥67,700	
57334-164	ふるい振とう器	¥58,900	¥66,500	
57334-166	ふるい振とう器	¥59,500	¥67,100	
57334-168	ふるい振とう器	¥59,500	¥67,100	
57334-170	ふるい振とう器	¥60,500	¥68,200	
57334-172	ふるい振とう器	¥59,500	¥67,100	
57334-174	ふるい振とう器	¥59,200	¥66,800	
57334-176	ふるい振とう器	¥60,400	¥68,100	
57334-178	ふるい振とう器	¥60,100	¥67,800	
57334-180	ふるい振とう器	¥59,500	¥67,100	
57334-182	ふるい振とう器	¥63,800	¥72,000	
57334-184	ふるい振とう器	¥63,700	¥71,900	
57334-186	ふるい振とう器	¥71,100	¥80,200	
57334-188	ふるい振とう器	¥71,400	¥80,600	
57334-190	ふるい振とう器	¥81,100	¥91,500	
57334-192	ふるい振とう器	¥97,400	¥109,900	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
57334-194	ふるい振とう器	¥117,600	¥132,700	
57334-196	ふるい振とう器	¥144,300	¥162,800	
57334-199	ふるい振とう器	¥80,200	¥90,500	
57334-200	ふるい振とう器	¥46,600	¥52,600	
57334-201	ふるい振とう器	¥70,900	¥80,100	
57334-202	ふるい振とう器	¥46,600	¥52,600	
57334-206	ふるい振とう器	¥46,000	¥51,800	
57334-208	ふるい振とう器	¥45,700	¥51,500	
57334-210	ふるい振とう器	¥46,000	¥51,800	
57334-212	ふるい振とう器	¥46,000	¥51,800	
57334-214	ふるい振とう器	¥46,000	¥51,800	
57334-216	ふるい振とう器	¥46,400	¥52,400	
57334-220	ふるい振とう器	¥46,800	¥52,800	
57334-222	ふるい振とう器	¥45,900	¥51,800	
57334-224	ふるい振とう器	¥47,100	¥53,200	
57334-226	ふるい振とう器	¥46,400	¥52,300	
57334-228	ふるい振とう器	¥47,100	¥53,100	
57334-230	ふるい振とう器	¥47,100	¥53,200	
57334-232	ふるい振とう器	¥46,400	¥52,300	
57334-234	ふるい振とう器	¥47,100	¥53,200	
57334-236	ふるい振とう器	¥47,100	¥53,200	
57334-238	ふるい振とう器	¥46,400	¥52,300	
57334-240	ふるい振とう器	¥47,100	¥53,200	
57334-242	ふるい振とう器	¥46,700	¥52,700	
57334-244	ふるい振とう器	¥47,100	¥53,100	
57334-246	ふるい振とう器	¥47,800	¥54,000	
57334-248	ふるい振とう器	¥47,100	¥53,200	
57334-250	ふるい振とう器	¥47,100	¥53,200	
57334-252	ふるい振とう器	¥47,100	¥53,200	
57334-254	ふるい振とう器	¥47,100	¥53,200	
57334-256	ふるい振とう器	¥47,100	¥53,200	
57334-258	ふるい振とう器	¥47,100	¥53,200	
57334-260	ふるい振とう器	¥47,100	¥53,200	
57334-262	ふるい振とう器	¥46,800	¥52,800	
57334-264	ふるい振とう器	¥47,100	¥53,200	
57334-266	ふるい振とう器	¥46,900	¥53,000	
57334-268	ふるい振とう器	¥47,100	¥53,200	
57334-270	ふるい振とう器	¥47,100	¥53,200	
57334-272	ふるい振とう器	¥47,100	¥53,200	
57334-274	ふるい振とう器	¥47,100	¥53,200	
57334-276	ふるい振とう器	¥47,100	¥53,200	
57334-278	ふるい振とう器	¥46,900	¥53,000	
57334-280	ふるい振とう器	¥46,900	¥53,000	
57334-282	ふるい振とう器	¥47,900	¥54,000	
57334-284	ふるい振とう器	¥47,100	¥53,200	
57334-286	ふるい振とう器	¥47,000	¥58,300	
57334-288	ふるい振とう器	¥47,500	¥53,700	
57334-290	ふるい振とう器	¥51,900	¥58,600	
57334-292	ふるい振とう器	¥56,300	¥63,500	
57334-294	ふるい振とう器	¥63,300	¥78,600	
57334-296	ふるい振とう器	¥71,400	¥80,600	
57334-298	ふるい振とう器	¥86,700	¥97,900	
57334-300	ふるい振とう器	¥105,500	¥125,000	
57334-302	ふるい振とう器	¥131,000	¥147,800	
57334-304	ふるい振とう器	¥155,500	¥175,500	
57334-306	ふるい振とう器	¥174,300	¥196,600	
57334-308	ふるい振とう器	¥429,100	¥484,200	
57334-310	ふるい振とう器	¥46,000	¥51,800	
57334-314	ふるい振とう器	¥45,800	¥51,700	
57334-324	ふるい振とう器	¥45,500	¥51,400	
57334-326	ふるい振とう器	¥46,000	¥51,800	
57334-328	ふるい振とう器	¥45,100	¥50,900	
57334-330	ふるい振とう器	¥46,400	¥52,400	
57334-334	ふるい振とう器	¥46,400	¥52,300	
57334-336	ふるい振とう器	¥45,600	¥51,400	
57334-339	ふるい振とう器	¥47,100	¥53,100	
57334-340	ふるい振とう器	¥46,400	¥52,300	
57334-342	ふるい振とう器	¥46,700	¥52,700	
57334-344	ふるい振とう器	¥46,400	¥52,300	
57334-346	ふるい振とう器	¥47,100	¥53,100	
57334-348	ふるい振とう器	¥47,900	¥54,000	
57334-350	ふるい振とう器	¥47,100	¥53,200	
57334-352	ふるい振とう器	¥45,900	¥51,800	
57334-354	ふるい振とう器	¥46,800	¥52,800	
57334-356	ふるい振とう器	¥47,100	¥53,200	
57334-358	ふるい振とう器	¥46,800	¥52,800	
57334-360	ふるい振とう器	¥47,900	¥54,000	
57334-364	ふるい振とう器	¥47,900	¥54,000	
57334-366	ふるい振とう器	¥51,600	¥58,200	
57334-368	ふるい振とう器	¥55,300	¥62,400	
57334-370	ふるい振とう器	¥63,300	¥71,500	
57334-372	ふるい振とう器	¥70,900	¥80,000	
57334-374	ふるい振とう器	¥86,100	¥97,200	
57334-376	ふるい振とう器	¥106,500	¥120,200	
57334-378	ふるい振とう器	¥128,600	¥145,000	
57334-380	ふるい振とう器	¥155,800	¥175,800	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
57334-382	ふるい振とう器	¥171,400	¥193,400	
57334-384	ふるい振とう器	¥427,100	¥481,900	
57334-385	ふるい振とう器	¥98,300	¥110,900	
57334-386	ふるい振とう器	¥67,700	¥76,300	
57334-387	ふるい振とう器	¥92,000	¥103,800	
57334-388	ふるい振とう器	¥66,900	¥75,500	
57334-392	ふるい振とう器	¥66,700	¥75,200	
57334-394	ふるい振とう器	¥67,900	¥76,700	
57334-396	ふるい振とう器	¥67,700	¥76,300	
57334-398	ふるい振とう器	¥66,700	¥75,200	
57334-400	ふるい振とう器	¥68,300	¥77,100	
57334-402	ふるい振とう器	¥66,900	¥75,500	
57334-404	ふるい振とう器	¥66,100	¥74,500	
57334-406	ふるい振とう器	¥68,000	¥76,700	
57334-408	ふるい振とう器	¥68,300	¥77,100	
57334-410	ふるい振とう器	¥67,700	¥76,400	
57334-412	ふるい振とう器	¥66,700	¥75,200	
57334-414	ふるい振とう器	¥67,700	¥76,300	
57334-416	ふるい振とう器	¥68,000	¥76,700	
57334-418	ふるい振とう器	¥68,300	¥77,100	
57334-420	ふるい振とう器	¥67,700	¥76,400	
57334-422	ふるい振とう器	¥67,600	¥76,200	
57334-424	ふるい振とう器	¥67,300	¥76,000	
57334-426	ふるい振とう器	¥67,700	¥76,400	
57334-428	ふるい振とう器	¥66,900	¥75,500	
57334-430	ふるい振とう器	¥67,300	¥76,000	
57334-432	ふるい振とう器	¥68,000	¥76,700	
57334-434	ふるい振とう器	¥67,700	¥76,400	
57334-436	ふるい振とう器	¥68,000	¥76,700	
57334-438	ふるい振とう器	¥69,000	¥77,900	
57334-440	ふるい振とう器	¥68,000	¥76,700	
57334-442	ふるい振とう器	¥68,000	¥76,700	
57334-444	ふるい振とう器	¥68,000	¥76,700	
57334-446	ふるい振とう器	¥68,000	¥76,700	
57334-448	ふるい振とう器	¥67,700	¥76,400	
57334-450	ふるい振とう器	¥68,000	¥76,700	
57334-452	ふるい振とう器	¥67,000	¥75,700	
57334-454	ふるい振とう器	¥68,000	¥76,800	
57334-456	ふるい振とう器	¥67,700	¥76,400	
57334-458	ふるい振とう器	¥67,700	¥76,400	
57334-460	ふるい振とう器	¥67,700	¥76,400	
57334-462	ふるい振とう器	¥68,000	¥76,700	
57334-464	ふるい振とう器	¥67,700	¥76,400	
57334-466	ふるい振とう器	¥68,000	¥76,700	
57334-468	ふるい振とう器	¥68,000	¥76,700	
57334-470	ふるい振とう器	¥68,000	¥76,700	
57334-472	ふるい振とう器	¥70,100	¥79,100	
57334-474	ふるい振とう器	¥71,300	¥80,500	
57334-476	ふるい振とう器	¥77,800	¥87,800	
57334-478	ふるい振とう器	¥85,500	¥96,500	
57334-480	ふるい振とう器 1個	¥92,400	¥104,200	
57334-482	ふるい振とう器	¥103,600	¥116,900	
57334-484	ふるい振とう器	¥121,400	¥137,000	
57334-486	ふるい振とう器	¥135,700	¥160,700	
57334-488	ふるい振とう器	¥174,700	¥207,000	
57334-490	ふるい振とう器	¥212,100	¥239,300	
57334-492	ふるい振とう器	¥226,200	¥255,200	
57334-494	ふるい振とう器	¥518,700	¥585,200	
57334-504	ふるい振とう器	¥66,700	¥75,200	
57334-516	ふるい振とう器	¥67,700	¥76,400	
57334-518	ふるい振とう器	¥66,700	¥75,200	
57334-520	ふるい振とう器	¥68,000	¥76,700	
57334-522	ふるい振とう器	¥68,300	¥77,100	
57334-526	ふるい振とう器	¥67,400	¥76,000	
57334-528	ふるい振とう器	¥67,700	¥76,300	
57334-530	ふるい振とう器	¥66,700	¥75,200	
57334-532	ふるい振とう器	¥67,700	¥76,400	
57334-534	ふるい振とう器	¥68,300	¥77,100	
57334-536	ふるい振とう器	¥66,400	¥74,900	
57334-538	ふるい振とう器	¥67,300	¥76,000	
57334-540	ふるい振とう器	¥67,300	¥76,000	
57334-542	ふるい振とう器	¥67,600	¥76,200	
57334-544	ふるい振とう器	¥67,900	¥76,600	
57334-546	ふるい振とう器	¥67,600	¥76,200	
57334-548	ふるい振とう器	¥67,900	¥76,600	
57334-550	ふるい振とう器	¥68,100	¥76,900	
57334-552	ふるい振とう器	¥66,900	¥75,500	
57334-554	ふるい振とう器	¥68,000	¥76,700	
57334-556	ふるい振とう器	¥68,000	¥76,700	
57334-558	ふるい振とう器	¥67,700	¥76,400	
57334-560	ふるい振とう器	¥69,000	¥77,900	
57334-562	ふるい振とう器	¥67,700	¥76,400	
57334-564	ふるい振とう器	¥68,000	¥76,700	
57334-566	ふるい振とう器	¥68,000	¥76,700	
57334-568	ふるい振とう器	¥68,000	¥76,700	
57334-570	ふるい振とう器	¥67,700	¥84,000	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
57334-572	ふるい振とう器	¥68,000	¥76,700	
57334-574	ふるい振とう器	¥67,700	¥76,400	
57334-576	ふるい振とう器	¥68,000	¥76,700	
57334-578	ふるい振とう器	¥67,700	¥76,400	
57334-580	ふるい振とう器	¥68,000	¥84,400	
57334-582	ふるい振とう器	¥69,800	¥86,600	
57334-584	ふるい振とう器	¥72,100	¥81,300	
57334-586	ふるい振とう器	¥77,800	¥87,800	
57334-588	ふるい振とう器	¥86,800	¥98,000	
57334-590	ふるい振とう器	¥92,400	¥104,200	
57334-592	ふるい振とう器	¥104,600	¥118,000	
57334-594	ふるい振とう器	¥120,600	¥136,100	
57334-596	ふるい振とう器	¥137,600	¥155,200	
57334-598	ふるい振とう器	¥173,500	¥195,800	
57334-600	ふるい振とう器	¥207,200	¥233,800	
57334-602	ふるい振とう器	¥220,300	¥248,500	
57334-604	ふるい振とう器	¥522,100	¥589,100	
57334-606	ふるい振とう器	¥46,400	¥52,300	
57334-610	ふるい振とう器	¥46,400	¥52,300	
57334-612	ふるい振とう器	¥46,400	¥52,300	
57334-614	ふるい振とう器	¥46,400	¥52,300	
57334-616	ふるい振とう器	¥47,100	¥53,100	
57334-618	ふるい振とう器	¥45,900	¥51,800	
57334-620	ふるい振とう器	¥46,400	¥52,300	
57334-622	ふるい振とう器	¥46,400	¥52,300	
57334-624	ふるい振とう器	¥46,400	¥52,300	
57334-626	ふるい振とう器	¥46,400	¥52,300	
57334-628	ふるい振とう器	¥46,400	¥52,300	
57334-630	ふるい振とう器	¥47,100	¥53,200	
57334-632	ふるい振とう器	¥47,100	¥53,200	
57334-634	ふるい振とう器	¥46,400	¥52,300	
57334-636	ふるい振とう器	¥45,900	¥51,800	
57334-640	ふるい用アセリ-	¥19,300	¥21,800	
57334-642	ふるい用アセリ-	¥19,300	¥21,800	
57334-644	ふるい用アセリ-	¥31,700	¥35,700	
57334-646	ふるい用アセリ-	¥32,600	¥36,800	
57334-648	ふるい用アセリ-	¥12,900	¥14,500	
57334-650	ふるい用アセリ-	¥16,200	¥18,300	
57334-652	ふるい用アセリ-	¥35,000	¥39,500	
57334-654	ふるい用アセリ-	¥35,600	¥40,100	
57334-656	ふるい用アセリ-	¥43,100	¥48,600	
57334-658	ふるい用アセリ-	¥43,100	¥48,600	
57334-660	ふるい用アセリ-	¥15,900	¥17,900	
57334-662	ふるい用アセリ-	¥23,900	¥27,000	
57355-006	ふるい用受け皿	¥24,500	¥27,300	
57355-028	ふるい用受け皿	¥26,000	¥28,900	
57355-040	ふるい用受け皿	¥44,500	¥49,500	
57355-061	ふるい用受け皿	¥73,600	¥81,900	
57356-020	ふるい用蓋	¥36,100	¥40,200	
57356-042	ふるい用蓋	¥23,900	¥26,600	
57357-001	ふるい用蓋	¥19,200	¥21,400	
57377-009	ふるい(染料・塗料分別用)	¥113,700	¥123,500	
57378-001	ふるい(染料・塗料分別用)	¥46,200	¥50,200	
57379-004	ふるい(染料・塗料分別用)	¥39,400	¥42,800	
57420-061	Sieve Sound Enclosure for Ro-Tap	¥699,600	¥1,370,000	
57462-006	実験台保護シート(撥水性珪素樹脂シート)	¥493,700	¥581,700	
57497-306	ピペットライナー 12入	¥72,100	¥77,100	
57901-106	Quaver Spatula	¥36,500	¥38,700	
57901-140	Quaver Teflon Spatula	¥49,000	¥51,900	
57930-025	ルードラバースパチラ	¥24,300	¥28,200	
57930-047	ルードラバースパチラ	¥26,900	¥31,200	
57930-069	ルードラバースパチラ	¥45,100	¥52,300	
57949-033	マイクロスパチラ	¥26,000	¥28,200	
57949-037	ステンレスマイクロスパチラ	¥5,800	¥6,300	
57949-039	ステンレスマイクロスパチラ	¥5,700	¥6,200	
57949-041	ステンレスマイクロスパチラ	¥5,700	¥6,200	
57952-311	ステンレスマイクロスポンジ	¥46,200	¥50,200	
58017-890	フィルタ(紫外線除菌用)1mm	¥210,900	¥235,100	
58017-895	フィルタ(紫外線除菌用)2mm 50個	¥215,800	¥240,500	
58017-900	フィルタ(紫外線除菌用)	¥215,800	¥240,500	
58017-929	Flatpack Chambers 50	¥354,900	¥395,300	
58019-508	ワイヤ(フォワード/電話機用ワイヤタイプ)	¥248,600	¥266,100	
58019-512	ワットケル(光学機器用ワイヤ) 100	¥46,400	¥49,700	
58020-014	Vakuwash Cuvette Washer 1個	¥139,500	¥150,800	
58020-107	RACK SPECTRO CELL POPYPROPYLENE	¥27,700	¥29,700	
58400-022	固定棚	¥1,073,300	¥1,580,000	
58540-047	ワイヤーストック(セロ製)	¥34,300	¥38,000	
58540-069	ワイヤーストック(セロ製)	¥27,200	¥33,100	
58542-246	ワイヤーストック(セロ製)	¥65,900	¥73,000	
58542-360	ワイヤーストック(セロ製)	¥101,200	¥112,100	
58544-106	ワッシュ(ステンレス製洗浄用)2個入	¥21,400	¥22,800	
58578-039	Replacement Cartridge for 1C200	¥85,700	¥95,800	
58615-000	オートレーブ(ホールド型)	¥316,400	¥575,000	
58618-009	オートレーブ(ホールド型)	¥273,800	¥498,000	
58752-861	滅菌用ラップ 30.5x30.5	¥142,100	¥156,600	
58752-920	滅菌用ラップ	¥232,400	¥256,100	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
58753-106	滅菌用パ〇子(蒸気/カ〇滅菌兼用)	¥86,000	¥94,800	
58753-139	滅菌用パ〇子(蒸気/カ〇滅菌兼用)	¥114,900	¥131,400	
58753-161	滅菌用パ〇子(蒸気/カ〇滅菌兼用)	¥165,800	¥182,700	
58753-194	滅菌用パ〇子(蒸気/カ〇滅菌兼用)	¥211,800	¥233,400	
58753-202	滅菌用パ〇子(蒸気/カ〇滅菌兼用)	¥150,700	¥166,100	
58753-220	滅菌用パ〇子(蒸気/カ〇滅菌兼用)	¥189,300	¥208,700	
58815-176	シール	¥815,900	¥1,513,000	
58815-206	Microtube Foam Inserts pack of 2	¥28,400	¥53,000	
58815-212	Platform Head Only, 15.2cm	¥23,900	¥44,000	
58815-252	20-Place Vortexer Attachment 1個	¥66,300	¥132,000	
58939-990	オートクレーブ (高圧蒸気滅菌器)	¥73,100	¥137,000	
58939-994	オートクレーブ (高圧蒸気滅菌器)	¥114,900	¥215,000	
58940-004	オートクレーブ (高圧蒸気滅菌器)	¥1,304,200	¥2,447,000	
58947-106	攪拌子(テロン樹脂コト)22.2x7.9mm	¥2,100	¥2,200	
58947-114	攪拌子(テロン樹脂コト)41.3x7.9mm	¥2,900	¥3,100	
58947-120	攪拌子(テロン樹脂コト)12.7x9.5mm	¥4,500	¥4,700	
58947-122	攪拌子(テロン樹脂コト)	¥2,800	¥2,900	
58947-124	攪拌子(テロン樹脂コト)22.2x9.5mm	¥2,800	¥3,000	
58947-128	攪拌子(テロン樹脂コト)34.9x9.5mm	¥3,400	¥4,000	
58947-132	攪拌子(テロン樹脂コト)47.6x9.5mm	¥5,600	¥6,700	
58947-136	攪拌子(テロン樹脂コト)	¥5,500	¥7,000	
58947-138	攪拌子(テロン樹脂コト) 38.1x12.7mm	¥14,000	¥19,200	
58947-140	攪拌子(テロン樹脂コト/テ〇イスク)	¥55,500	¥58,800	
58947-142	攪拌子(テロン樹脂コト/テ〇イスク)	¥79,000	¥84,500	
58947-144	攪拌子(テロン樹脂コト/テ〇イスク)	¥95,800	¥103,000	
58947-146	攪拌子(テロン樹脂コト/テ〇イスク)	¥131,700	¥142,400	
58947-150	攪拌子(テロン樹脂)	¥4,600	¥5,000	
58947-412	攪拌機用羽根(テロン)	¥5,500	¥10,000	
58947-820	攪拌子(テロン樹脂コト/Spinplus)	¥6,100	¥6,500	
58947-822	攪拌子(テロン樹脂コト/Spinplus)	¥7,200	¥7,700	
58947-824	攪拌子(テロン樹脂コト/Spinplus)	¥7,100	¥7,600	
58947-826	攪拌子(テロン樹脂コト/Spinplus)	¥8,000	¥8,600	
58947-828	攪拌子(テロン樹脂コト/Spinplus)	¥8,100	¥9,000	
58947-830	攪拌子(テロン樹脂コト/Spinplus)	¥38,400	¥41,000	
58947-849	攪拌子(テロン樹脂コト/Circulus)初付	¥8,700	¥10,500	
58947-860	攪拌子(テロン樹脂コト/Circulus)初付	¥6,800	¥8,200	
58947-907	攪拌子(テロン樹脂コト/Circulus)セット	¥6,800	¥8,200	
58947-929	攪拌子(テロン樹脂コト/Circulus)10-	¥6,200	¥6,600	
58947-962	攪拌子(テロン樹脂コト/Circulus)7-	¥8,000	¥8,500	
58948-003	攪拌子セット(テロン樹脂コト)	¥47,900	¥60,400	
58948-025	攪拌子セット(Octagon)	¥15,300	¥16,400	
58948-069	攪拌子(テロン樹脂コト/小型)	¥16,700	¥17,700	
58948-070	攪拌子(テロン樹脂コト)	¥52,900	¥64,900	
58948-072	攪拌子(テロン樹脂コト)	¥72,200	¥83,300	
58948-074	攪拌子(テロン樹脂コト/カ〇ル〇)	¥127,800	¥169,100	
58948-078	攪拌子セット(テロン樹脂コト/タ〇型)	¥84,200	¥102,000	
58948-080	攪拌子(テロン樹脂コト)22.2x4.8mm	¥3,100	¥3,600	
58948-091	攪拌子(テロン樹脂コト)12.7x3mm	¥4,200	¥4,400	
58948-116	攪拌子(テロン樹脂コト)12.7x7.9mm	¥4,100	¥4,400	
58948-138	攪拌子(テロン樹脂コト)25.4x7.9mm	¥3,800	¥4,100	
58948-150	攪拌子(テロン樹脂コト)38.1x7.9mm	¥4,800	¥5,400	
58948-171	攪拌子(テロン樹脂コト)50.8x7.9mm	¥5,700	¥6,100	
58948-193	攪拌子(テロン樹脂コト)76.2x12.7mm	¥15,900	¥16,800	
58948-218	攪拌子(テロン樹脂コト)15.9x7.9mm	¥1,900	¥2,000	
58948-230	攪拌子(テロン樹脂コト)25.4x9.5mm	¥4,300	¥4,600	
58948-251	攪拌子(テロン樹脂コト)38.1x9.5mm	¥3,100	¥3,600	
58948-273	攪拌子(テロン樹脂コト)50.8x9.5mm	¥5,300	¥6,200	
58948-280	攪拌子回収セット	¥25,300	¥27,000	
58948-309	攪拌子(テロン樹脂コト)63.5x7.9mm	¥6,000	¥6,400	
58948-353	攪拌子(テロン樹脂コト/小型) 1	¥3,100	¥3,300	
58948-375	攪拌子(テロン樹脂コト/小型)	¥110,300	¥116,800	
58948-377	VWR SPINBAR MAGNETIC 5x2mm	¥4,600	¥4,900	
58948-397	攪拌子(テロン樹脂コト/小型)	¥109,800	¥116,300	
58948-400	攪拌子(テロン樹脂コト/小型)	¥123,100	¥131,800	
58948-411	攪拌子(テロン樹脂コト/小型)	¥100,000	¥105,900	
58948-422	攪拌子 (カ〇製)22.2x6.4mm	¥19,100	¥20,400	
58948-444	攪拌子 (カ〇製)	¥19,100	¥20,300	
58948-466	攪拌子 (カ〇製)	¥19,200	¥20,400	
58948-488	攪拌子 (カ〇製)	¥20,100	¥21,400	
58948-950	攪拌子(テロン樹脂コト)28.6x7.9mm	¥2,200	¥2,400	
58948-952	攪拌子(テロン樹脂)	¥27,600	¥33,400	
58948-954	攪拌子(テロン樹脂)	¥4,800	¥5,200	
58948-956	攪拌子(テロン樹脂)	¥2,800	¥2,900	
58948-960	攪拌子(テロン樹脂)	¥3,700	¥3,900	
58948-962	攪拌子(テロン樹脂)	¥31,500	¥39,600	
58948-966	攪拌子(テロン樹脂)	¥5,100	¥5,500	
58948-968	攪拌子(テロン樹脂)	¥6,200	¥6,500	
58948-972	攪拌子(テロン樹脂コト)	¥62,000	¥66,000	
58948-974	攪拌子(テロン樹脂コト)	¥75,900	¥80,900	
58948-976	攪拌子(テロン樹脂コト/小型)	¥98,100	¥105,000	
58948-978	攪拌子(テロン樹脂コト)	¥76,900	¥81,900	
58948-980	攪拌子(テロン樹脂コト)50.8x9.5mm	¥4,500	¥4,800	
58948-981	攪拌子(テロン樹脂コト)76.2x12.7mm	¥10,100	¥13,300	
58948-982	攪拌子(テロン樹脂コト)12.7x7.9mm	¥2,300	¥2,400	
58948-983	攪拌子(テロン樹脂コト)25.4x9.5mm	¥2,600	¥3,000	
58948-984	攪拌子(テロン樹脂コト)15.9x7.9mm	¥2,100	¥2,300	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
58948-985	攪拌子(テロン樹脂コート)22.2x9.5mm	¥2,700	¥2,900	
58948-986	攪拌子(テロン樹脂コート)63.5x12.7mm	¥12,000	¥15,800	
58948-988	攪拌子(テロン樹脂コート)38.1x9.5mm	¥3,100	¥3,600	
58948-998	攪拌子(テロン樹脂コート/Double Spinfan)	¥7,900	¥8,400	
58949-002	攪拌子(テロン樹脂コート/タコ型)	¥13,700	¥19,400	
58949-004	攪拌子(テロン樹脂コート/Double Spinfan)	¥9,300	¥13,100	
58949-006	攪拌子(テロン樹脂コート/タコ型)	¥4,700	¥5,000	
58949-008	攪拌子(テロン樹脂コート/タコ型)	¥15,000	¥23,400	
58949-010	攪拌子(テロン樹脂コート/タコ型)	¥6,800	¥7,200	
58949-012	攪拌子(テロン樹脂コート/タコ型)	¥5,100	¥5,400	
58949-016	攪拌子(テロン樹脂コート/タコ型)	¥6,900	¥7,400	
58949-034	攪拌子(テロン樹脂コート)	¥14,300	¥15,300	
58949-036	攪拌子(テロン樹脂コート)	¥12,600	¥13,500	
58949-038	攪拌子(テロン樹脂コート)	¥17,200	¥18,600	
58949-040	攪拌子(テロン樹脂コート)	¥12,700	¥13,600	
58949-042	攪拌子(テロン樹脂コート/カーハック)	¥18,500	¥19,600	
58949-044	攪拌子(テロン樹脂コート/カー)	¥2,500	¥3,200	
58949-046	攪拌子(テロン樹脂コート/カー)	¥2,500	¥3,200	
58949-048	攪拌子(テロン樹脂コート/カー)	¥5,200	¥6,200	
58949-050	攪拌子(テロン樹脂コート/カー)	¥2,500	¥2,700	
58949-052	攪拌子(テロン樹脂コート/カー)	¥4,000	¥4,800	
58949-054	攪拌子(テロン樹脂コート/カー)	¥4,000	¥4,800	
58949-056	攪拌子(テロン樹脂コート/カー)	¥3,800	¥5,000	
58949-058	攪拌子(テロン樹脂コート/カー)	¥3,800	¥5,000	
58949-060	攪拌子(テロン樹脂コート/カー)	¥3,800	¥5,000	
58949-061	攪拌子(テロン樹脂)	¥2,400	¥2,600	
58949-062	攪拌子(テロン樹脂コート/カー)	¥3,100	¥3,800	
58949-064	攪拌子(テロン樹脂コート/カー)	¥3,100	¥3,800	
58949-066	攪拌子(テロン樹脂コート/カー)	¥3,000	¥3,800	
58949-083	攪拌子(テロン樹脂)	¥2,600	¥2,800	
58949-094	攪拌子(テロン樹脂)	¥4,500	¥4,700	
58949-107	Magnetic Stirring Bar 1個	¥6,900	¥7,400	
58949-108	攪拌子(テロン樹脂)	¥4,400	¥4,700	
58949-130	攪拌子(テロン樹脂)	¥5,000	¥5,300	
58949-152	攪拌子(テロン樹脂)	¥12,800	¥14,800	
58949-196	攪拌子(テロン樹脂コート/タコ型)	¥6,300	¥6,600	
58949-200	攪拌子(テロン樹脂コート/タコ型)	¥5,800	¥6,200	
58949-210	攪拌子(テロン樹脂コート/タコ型)	¥7,900	¥8,500	
58949-232	攪拌子(テロン樹脂コート/くさび型)	¥3,500	¥3,800	
58949-243	攪拌子(テロン樹脂コート/タコ型)	¥11,900	¥12,700	
58949-254	攪拌子(テロン樹脂コート/くさび型)	¥6,100	¥8,600	
58949-265	攪拌子(テロン樹脂コート/くさび型)	¥5,800	¥8,200	
58949-270	攪拌子(テロン樹脂コート/Spinvane)	¥6,300	¥6,700	
58949-272	攪拌子(テロン樹脂コート/Spinvane)	¥4,100	¥4,400	
58949-276	攪拌子(テロン樹脂コート/Spinvane)	¥6,100	¥6,500	
58949-280	攪拌子回収口 305mm	¥14,000	¥14,900	
58949-287	攪拌子回収口 457mm	¥11,600	¥12,400	
58949-700	プラスチック廃棄用ゴミ箱タンク・ホルダー	¥53,500	¥60,300	
58976-921	Sample Bags(Standard)	¥88,500	¥98,500	
58976-929	ポリバケツ 33x18.5cm 500枚入	¥238,800	¥261,000	
58981-009	ミキサー用500mlバケツ容器	¥270,200	¥499,000	
58981-020	ミキサー交換用ステンレスプレート	¥80,500	¥163,000	
58983-139	ステンレス製ミニコンテナ 37-110mL 1個	¥320,400	¥346,400	
59060-025	攪拌棒(ガラス製)	¥102,600	¥111,900	
59060-047	攪拌棒(ガラス製) 100本入	¥49,300	¥53,700	
59060-069	攪拌棒(ガラス製) 200x5	¥82,000	¥89,400	
59060-080	攪拌棒(ガラス製)	¥188,100	¥205,100	
59060-105	攪拌棒(ガラス製)	¥119,000	¥129,800	
59060-127	攪拌棒(ガラス製)	¥179,500	¥195,800	
59063-001	攪拌棒(PP製)	¥8,900	¥9,600	
59075-002	踏み台(キャスター付)	¥188,800	¥216,000	
59075-005	踏み台(キャスター付)交換用ゴム	¥8,700	¥9,400	
59075-022	踏み台(ステップ スワivelキャスター付) 1個	¥125,400	¥135,800	
59075-077	踏み台(キャスター付)	¥201,400	¥211,400	
59075-412	踏み台(キャスター付)	¥176,500	¥201,900	
59075-428	踏み台(キャスター付) 1台	¥194,900	¥223,000	
59075-438	踏み台(キャスター付)	¥192,400	¥208,000	
59075-440	踏み台(キャスター付) 1台	¥177,400	¥201,100	
59097-102	Replacement Stopcocks with Plug	¥48,500	¥53,400	
59-2000	ntellixseal SA セミオートヒートシーラー	¥950,000	オープン価格	
59-2001	Standard 96&384 Plate用アダプター	¥22,700	¥23,000	
59-2002	Roche 1536Plate アダプター 1	¥57,500	¥59,000	
59-2003	96well Plate用アダプター 1	¥83,100	¥86,000	
59-2004	384well Plate用アダプター 1	¥43,800	¥45,000	
59-2005	4s3 ラックアダプター用プレートホルダーアダプター 1	¥98,500	¥102,000	
59-2007	プレートホルダーアダプター 1	¥134,500	¥139,000	
59-2008	Weighted Sealing Platen 1	¥20,500	¥21,000	
59-2009	シーリングゴム 1	¥20,500	¥21,000	
59344-055	Haigh Vacuum Grease 150g 1本	¥28,400	¥30,400	
59580-047	ゴム栓	¥48,200	¥52,600	
59580-069	ゴム栓	¥33,200	¥36,200	
59580-080	ゴム栓	¥29,800	¥32,500	
59580-105	ゴム栓	¥31,700	¥34,600	
59580-127	ゴム栓	¥32,300	¥35,200	
59580-149	ゴム栓	¥38,600	¥42,200	
59580-160	ゴム栓	¥40,100	¥43,800	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
59580-182	コム栓	¥42,200	¥46,100	
59580-207	コム栓	¥41,000	¥44,700	
59580-240	コム栓	¥45,100	¥49,200	
59580-262	コム栓	¥39,400	¥43,000	
59580-284	コム栓	¥44,600	¥48,700	
59580-309	コム栓	¥39,600	¥43,200	
59580-342	コム栓	¥43,000	¥47,000	
59580-364	コム栓	¥36,600	¥40,000	
59580-386	コム栓	¥42,800	¥46,700	
59580-400	コム栓	¥38,700	¥42,200	
59580-422	コム栓	¥42,500	¥46,400	
59580-444	コム栓	¥48,300	¥52,800	
59580-466	コム栓	¥45,400	¥49,600	
59580-488	コム栓	¥43,100	¥47,100	
59580-502	コム栓	¥21,300	¥23,300	
59580-524	コム栓	¥22,600	¥24,600	
59580-546	コム栓	¥34,900	¥38,100	
59581-221	コム栓(フック/サイズ4) 36入	¥45,100	¥49,300	
59581-287	コム栓(フック/サイズ6) 23入	¥46,800	¥51,100	
59581-301	コム栓(フック/サイズ6.5) 21入	¥67,000	¥73,200	
59585-164	コム栓(穴あき)	¥44,800	¥48,900	
59585-200	コム栓(穴あき)	¥39,100	¥42,600	
59585-288	コム栓(穴あき) サイズ9	¥36,200	¥39,500	
59585-302	コム栓(穴あき)	¥38,500	¥42,000	
59585-346	コム栓(穴あき)	¥55,100	¥60,200	
59586-157	栓(穴あき)	¥35,300	¥38,600	
59586-158	栓(穴あき)	¥35,000	¥38,200	
59586-159	栓(穴あき)	¥35,300	¥38,600	
59586-161	栓(穴あき)	¥34,300	¥37,400	
59586-162	栓(穴あき)	¥34,300	¥37,400	
59586-163	栓(穴あき)	¥34,100	¥37,200	
59586-164	栓(穴あき)	¥35,300	¥38,600	
59589-074	コム栓	¥53,500	¥58,500	
59589-096	コム栓	¥36,200	¥39,500	
59589-110	コム栓	¥48,200	¥52,600	
59589-132	コム栓	¥51,000	¥55,700	
59589-154	コム栓	¥50,800	¥55,400	
59589-176	コム栓	¥51,800	¥56,600	
59589-198	コム栓	¥48,100	¥52,500	
59589-212	コム栓	¥51,200	¥55,900	
59589-234	コム栓	¥47,400	¥51,800	
59589-278	コム栓	¥46,400	¥50,700	
59589-290	コム栓	¥47,600	¥52,000	
59589-314	コム栓	¥49,700	¥54,200	
59589-336	コム栓	¥48,200	¥52,600	
59589-358	コム栓	¥56,900	¥62,100	
59589-370	コム栓	¥46,700	¥51,000	
59589-374	コム栓	¥44,200	¥48,300	
59590-109	刃コム栓	¥14,400	¥13,900	
59590-120	刃コム栓	¥17,800	¥17,300	
59590-142	刃コム栓	¥18,900	¥18,300	
59590-164	刃コム栓	¥29,700	¥28,800	
59590-186	刃コム栓	¥38,000	¥36,800	
59590-200	刃コム栓	¥45,500	¥44,100	
60002-246	Synthware平底ワシ 1L 24/40 1個	¥29,000	¥31,900	
60002-248	Synthware平底ワシ 1L 29/42 1個	¥30,700	¥33,700	
60003-288	ReactionTube Flsk w/High Vac Valve	¥48,700	¥53,500	
60004-930	Ribbed Sleeves 24/40 1個	¥14,800	¥16,300	
60004-970	Glass Cutter 1本	¥33,400	¥36,700	
60005-092	Synthware Cannula 18G	¥34,000	¥37,400	
60079-000	栓-トプレート	¥60,000	¥65,300	
60079-010	VWR Talon Connector	¥8,900	¥10,000	
60079-053	マルチポートコネクタ	¥17,200	¥18,900	
60079-100	VWR Talon End-to-End Connector	¥12,000	¥13,100	
60079-106	VWR Talon Horizontal Mounting Bars	¥17,200	¥18,900	
60079-107	VWR Talon Horizontal Mounting Bars	¥19,100	¥20,900	
60079-166	フックコネクタ	¥6,200	¥6,800	
60079-257	フレーム用安定台	¥12,600	¥13,800	
60079-461	マルチポート 1524mm	¥26,900	¥29,500	
60079-483	マルチポート 1829mm	¥28,700	¥31,300	
60079-530	ステンレスポート	¥5,200	¥5,600	
60079-531	ステンレスポート	¥5,200	¥5,700	
60079-532	ステンレスポート	¥8,300	¥9,000	
60079-533	ステンレスポート	¥14,000	¥15,300	
60079-534	ステンレスポート	¥18,400	¥20,000	
60079-535	ステンレスポート	¥23,300	¥25,500	
60079-536	ステンレスポート	¥32,300	¥35,400	
60079-537	ステンレスポート	¥41,400	¥45,200	
60079-538	ステンレスポート	¥49,400	¥54,000	
60079-539	ステンレスポート	¥58,000	¥63,600	
60079-540	ステンレスポート	¥75,300	¥81,400	
60097-055	フレームコネクタ	¥126,200	¥138,400	
60097-056	VWR Talon S-Connector	¥12,600	¥13,800	
60100-105	支持ポート	¥197,900	¥216,900	
60119-001	フレームリフト	¥154,000	¥161,700	
60141-000	運搬用具(ホッチャ用リフト)	¥365,400	¥454,000	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
60141-124	クランプ (ナイロンベルト) ストラップ 付	¥53,500	¥57,900	
60141-126	クランプ (ナイロンベルト) ストラップ・チェーン付	¥117,500	¥128,600	
60142-003	クランプ (ナイロンベルト)	¥42,300	¥47,000	
60142-004	交換用ストラップ	¥55,600	¥59,700	
60142-006	クランプ (ナイロンベルト)	¥73,700	¥80,800	
60142-007	交換用チェーン	¥55,600	¥60,300	
60315-000	表面張力測定装置	¥2,914,800	¥5,534,000	
60315-002	表面張力測定装置	¥16,000	¥31,000	
60315-004	表面張力測定装置	¥554,900	¥1,065,000	
60315-008	表面張力測定装置	¥3,578,000	¥6,793,000	
60315-010	表面張力測定装置	¥16,300	¥31,000	
60315-012	表面張力測定装置	¥552,200	¥1,060,000	
60371-000	マイクロシリンダ 700(Luer Tip) 10ul	¥45,500	¥49,200	
60371-011	マイクロシリンダ 700(Luer Tip) 25ul	¥57,800	¥62,400	
60371-022	マイクロシリンダ 700(Luer Tip) 50ul	¥57,100	¥61,700	
60371-044	マイクロシリンダ 700(Luer Tip) 100ul	¥61,700	¥66,700	
60371-066	マイクロシリンダ 700(Luer Tip) 250ul	¥66,400	¥71,800	
60371-088	マイクロシリンダ 700(Luer Tip) 500ul	¥75,400	¥81,500	
60372-025	マイクロシリンダ 700 26S 10ul	¥130,300	¥140,800	
60372-091	マイクロシリンダ 700 26S 10ul	¥28,900	¥31,200	
60372-105	マイクロシリンダ 700 22S 25ul	¥38,800	¥42,000	
60373-028	マイクロシリンダ 700 22S 50ul	¥47,800	¥51,700	
60373-040	マイクロシリンダ 700 22S 100ul	¥63,100	¥68,300	
60373-061	マイクロシリンダ 700 22S 250ul	¥56,100	¥60,600	
60373-083	マイクロシリンダ 700 22S 500ul	¥65,100	¥70,400	
60373-323	マイクロシリンダ 700シリーズ 500ul	¥85,500	¥92,400	
60373-345	マイクロシリンダ 700シリーズ 250ul	¥74,600	¥80,600	
60373-361	マイクロシリンダ 600シリーズ 2.5ul	¥96,200	¥104,000	
60373-370	マイクロシリンダ 700 26S 5ul	¥42,400	¥45,800	
60373-375	マイクロシリンダ 700 26S 5ul	¥51,800	¥56,000	
60373-402	マイクロシリンダ 7673A 5ul	¥56,000	¥60,600	
60373-404	マイクロシリンダ 7673A 5ul	¥300,200	¥324,500	
60373-405	マイクロシリンダ 7673A 5ul	¥302,600	¥327,100	
60373-408	マイクロシリンダ 7673A 5ul	¥53,400	¥57,700	
60373-409	マイクロシリンダ 7673A 5ul	¥301,300	¥325,700	
60373-425	Guide Assembly(5-10ulシリンダ用)	¥36,800	¥39,800	
60373-470	マイクロシリンダ 700 26S 10ul	¥55,400	¥59,900	
60373-471	マイクロシリンダ 700 26S 10ul	¥248,900	¥269,100	
60373-475	マイクロシリンダ 7673A 10ul	¥65,800	¥71,100	
60373-476	マイクロシリンダ 7673A 10ul	¥65,500	¥70,800	
60373-477	マイクロシリンダ 7673A 10ul	¥49,700	¥53,800	
60373-505	マイクロシリンダ 700 22S 25ul	¥51,400	¥55,500	
60373-538	マイクロシリンダ 700 22S 50ul	¥59,300	¥64,200	
60373-560	マイクロシリンダ 700 22S 100ul	¥66,400	¥71,800	
60373-593	マイクロシリンダ 700 22S 250ul	¥75,800	¥81,900	
60373-629	マイクロシリンダ 700 22シリーズ 500ul	¥80,800	¥87,300	
60373-946	マイクロシリンダ 700シリーズ 10ul	¥33,200	¥35,800	
60373-950	マイクロシリンダ 7673A 10ul	¥36,100	¥39,100	
60373-952	マイクロシリンダ 7673A 10ul	¥37,900	¥40,900	
60373-953	マイクロシリンダ 7673A 10ul	¥201,400	¥217,800	
60373-954	マイクロシリンダ 7673A 10ul	¥211,500	¥228,600	
60373-956	マイクロシリンダ 7673A 10ul	¥208,100	¥224,900	
60373-957	マイクロシリンダ 700シリーズ 25ul	¥35,000	¥37,800	
60373-958	マイクロシリンダ 7673A 10ul	¥37,300	¥40,400	
60373-960	マイクロシリンダ 600シリーズ 5.0ul	¥101,400	¥109,600	
60373-968	マイクロシリンダ 700シリーズ 50ul	¥42,400	¥45,800	
60373-979	マイクロシリンダ 700シリーズ (22S/100uL)	¥49,900	¥53,900	
60373-985	マイクロシリンダ 1700シリーズ 10ul	¥49,300	¥53,300	
60373-987	マイクロシリンダ 1700シリーズ 25ul	¥59,100	¥63,900	
60373-988	マイクロシリンダ 1700シリーズ 50ul	¥51,500	¥55,600	
60373-990	マイクロシリンダ 1700シリーズ 100ul	¥60,000	¥64,900	
60373-992	マイクロシリンダ 1700シリーズ 250ul	¥53,500	¥57,800	
60373-994	マイクロシリンダ 1700シリーズ 500ul	¥65,400	¥70,700	
60374-004	マイクロシリンダ 7000 0.5ul	¥95,500	¥103,200	
60374-006	マイクロシリンダ 7000 22S/1ul	¥84,000	¥90,800	
60374-008	マイクロシリンダ 7000 2ul	¥90,400	¥97,700	
60374-010	マイクロシリンダ 7000 2ul	¥93,600	¥101,200	
60374-012	マイクロシリンダ 7000 5ul	¥82,800	¥89,500	
60374-016	Chaney Adapter(マイクロシリンダ用)	¥84,600	¥91,500	
60374-018	Guide Assembly(マイクロシリンダ用)	¥35,000	¥37,800	
60374-024	Chaney Adapter(5-10ulシリンダ用)	¥100,600	¥108,700	
60374-030	マイクロシリンダ 1700 26S 10ul	¥60,900	¥65,900	
60374-036	マイクロシリンダ 700 26S 10ul	¥93,000	¥100,600	
60375-045	マイクロシリンダ 7000 22S/1ul	¥70,100	¥75,700	
60375-158	マイクロシリンダ 7000 25S/1ul	¥106,700	¥115,300	
60375-373	マイクロシリンダ 1700 50ul	¥49,500	¥53,600	
60375-409	マイクロシリンダ 1700 100ul	¥54,500	¥59,000	
60375-431	マイクロシリンダ 1700 250ul	¥53,000	¥57,300	
60375-442	マイクロシリンダ 1700 22S 250ul	¥63,800	¥69,000	
60375-464	マイクロシリンダ 1700 500ul	¥50,500	¥54,600	
60375-497	マイクロシリンダ 1000/S2 1ml	¥61,600	¥66,600	
60375-522	マイクロシリンダ 1000/S2 2.5ml	¥63,300	¥68,400	
60375-555	マイクロシリンダ 1000/S2 5ml	¥86,100	¥93,100	
60375-588	マイクロシリンダ 1000/S2 10ml	¥75,700	¥81,900	
60375-624	マイクロシリンダ 1000/S2 25ml	¥266,300	¥287,900	
60375-668	マイクロシリンダ 1000/S2 50ml	¥375,600	¥406,100	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
60375-670	マイクロシリンジ 1000/S2 100ml	¥522,700	¥565,100	
60376-026	マイクロシリンジ 1700 50ul	¥47,200	¥51,100	
60376-048	マイクロシリンジ 1700 100ul	¥49,300	¥53,300	
60376-060	マイクロシリンジ 1700 250ul	¥48,700	¥52,700	
60376-081	マイクロシリンジ 1700 500ul	¥46,400	¥50,200	
60376-106	マイクロシリンジ 1000/S2(Luer Tip) 1ml	¥53,100	¥57,400	
60376-128	マイクロシリンジ 1000/S2(Luer Tip) 2.5ml	¥46,300	¥55,100	
60376-140	マイクロシリンジ 1000/S2(Luer Tip) 5ml	¥56,400	¥60,900	
60376-161	マイクロシリンジ 1000/S2(Luer Tip) 10ml	¥67,500	¥72,900	
60376-172	マイクロシリンジ 1700 26S 10ul	¥44,100	¥47,700	
60376-177	マイクロシリンジ 7673A 10ul	¥56,800	¥61,400	
60376-178	マイクロシリンジ 7673A 10ul	¥61,300	¥66,200	
60376-179	マイクロシリンジ 7673A 10ul	¥60,400	¥65,300	
60376-180	マイクロシリンジ 7673A 10ul	¥51,000	¥55,100	
60376-182	マイクロシリンジ 7673A 10ul	¥320,200	¥346,200	
60376-183	マイクロシリンジ 1700 22S 25ul	¥56,000	¥60,500	
60376-194	マイクロシリンジ 1700 22S 50ul	¥51,100	¥55,300	
60376-208	マイクロシリンジ 1700 100ul	¥57,700	¥62,300	
60376-219	マイクロシリンジ 1700 500ul	¥60,400	¥71,800	
60376-220	マイクロシリンジ 1700 26S 10ul	¥52,000	¥56,200	
60376-230	マイクロシリンジ 1700 22S 25ul	¥48,300	¥52,200	
60376-241	マイクロシリンジ 1700 22S 50ul	¥51,900	¥56,100	
60376-252	マイクロシリンジ 1700 22S 100ul	¥47,300	¥51,100	
60376-263	マイクロシリンジ 1700 22S 250ul	¥50,300	¥54,400	
60376-274	マイクロシリンジ 1700 500ul	¥47,900	¥51,700	
60376-285	マイクロシリンジ 1000/S2 1ml	¥57,400	¥62,000	
60376-300	マイクロシリンジ 1000/S2 2.5ml	¥48,900	¥52,800	
60376-321	マイクロシリンジ 1000/S2 5ml	¥61,000	¥65,900	
60376-354	マイクロシリンジ 1000/S2 10ml	¥72,400	¥78,300	
60376-412	マイクロシリンジ 1000/S2 2.5ml	¥78,400	¥84,700	
60376-423	マイクロシリンジ 1000/S2 5ml	¥76,000	¥82,200	
60376-434	マイクロシリンジ 1000/S2 10ml	¥93,400	¥100,900	
60376-438	マイクロシリンジ 1000/S5 5ml	¥71,400	¥77,200	
60376-551	マイクロシリンジ 1000/S5 0.5ml	¥66,800	¥72,200	
60376-558	マイクロシリンジ 1000/S2 1ml	¥92,800	¥100,300	
60376-561	マイクロシリンジ 1000/S5 1ml	¥54,900	¥59,400	
60376-565	マイクロシリンジ 1000/S5 2.5ml	¥57,000	¥61,600	
60463-107	Rectangular Tanks 22.7L	¥186,700	¥209,700	
60490-100	ジョイント テープ	¥19,500	¥20,900	
60600-310	加針チューブ (カーラール付 ガラス製)	¥22,200	¥24,200	
60600-316	加針チューブ (カーラール付 ガラス製)	¥22,200	¥24,100	
60600-320	加針チューブ (カーラール付 ガラス製)	¥27,700	¥30,200	
60600-326	加針チューブ (カーラール付 ガラス製)	¥27,900	¥30,500	
60600-334	加針チューブ (カーラール付 ガラス製)	¥27,300	¥29,800	
60782-425	スリッパチップ 2mL ST 自立型 50x10	¥81,900	¥88,100	
60799-008	Test Paper 100stripx25vial	¥68,000	¥79,200	
60818-270	加針チューブ (ガラス製) 5mL 12x75mm	¥60,700	¥65,900	
60819-854	加針チューブ 6x50mm 100本	¥8,900	¥9,700	
60820-114	Capillary Tube Sealant	¥43,600	¥47,900	
60820-116	Funnel PS 55mm 100個入	¥35,000	¥37,700	
60820-140	Tiered Plug Stopper Closure 1000	¥15,500	¥17,000	
60825-402	加針チューブ (ガラス製)	¥34,800	¥37,800	
60825-406	加針チューブ (ガラス製)	¥37,200	¥40,400	
60825-414	加針チューブ (ガラス製)	¥46,300	¥50,300	
60825-425	加針チューブ (ガラス製) 16x100 1000本	¥86,500	¥93,900	
60825-430	加針チューブ (ガラス製)	¥71,600	¥77,800	
60825-435	加針チューブ (ガラス製)	¥75,800	¥82,300	
60825-447	加針チューブ (ガラス製)	¥83,300	¥90,500	
60825-500	KIMAXワンド 加針チューブ 4ml	¥49,900	¥54,400	
60825-502	加針チューブ (ガラス製) 6ml 12x75mm	¥66,300	¥72,300	
60825-504	KIMAX ティスト 加針チューブ 20ml 1000本	¥109,300	¥119,200	
60825-510	KIMAXワンド 加針チューブ 16x100mm	¥91,800	¥100,100	
60825-837	加針チューブ ガラス 16x100mm 15ml 1000	¥233,900	¥255,000	
60825-979	KIMAX加針チューブ 24mL 16x150cm	¥109,500	¥119,400	
60825-980	KIMAX加針チューブ 28mL 500本入	¥100,700	¥109,800	
60826-110	加針チューブ (カーラール付 ガラス製)	¥84,700	¥92,300	
60826-144	加針チューブ (カーラール付 ガラス製)	¥164,300	¥179,200	
60826-304	遠心チューブ (ガラス製+ガラス)	¥103,800	¥113,100	
60826-315	遠心チューブ (ガラス製+ガラス)	¥111,800	¥121,900	
60826-825	ティスト 加針チューブ 4ml	¥62,800	¥67,800	
60827-453	加針チューブ 9ml 13x100mm 144本	¥307,900	¥332,800	
60828-680	スプレッダー - BactiCell Spreaders	¥43,400	¥48,300	
60828-682	スプレッダー - BactiCell Spreaders	¥138,500	¥154,100	
60828-684	スプレッダー - BactiCell Spreaders	¥42,800	¥47,600	
60828-688	スプレッダー - BactiCell Spreaders	¥138,500	¥154,100	
60828-720	VWR Safe-T-Flex Caps 1000個入	¥33,400	¥37,000	
60828-722	VWR Safe-T-Flex Caps 1000個入	¥75,700	¥83,900	
60828-724	VWR Safe-T-Flex Caps 1000個入	¥33,400	¥37,000	
60828-726	VWR Safe-T-Flex Caps 1000個入	¥75,700	¥92,300	
60828-728	VWR Safe-T-Flex Caps 1000個入	¥33,500	¥37,100	
60828-730	VWR Safe-T-Flex Caps 1000個入	¥33,200	¥36,800	
60828-732	VWR Safe-T-Flex Caps 1000個入	¥33,500	¥37,100	
60828-760	VWR Safe-T-Flex Caps 1000個入	¥33,100	¥36,600	
60828-766	VWR Safe-T-Flex Caps 1000個入	¥33,400	¥37,000	
60869-056	チップのみ (スリッパチップ サイズ 1) 500個	¥28,100	¥30,600	
60869-089	チップのみ (スリッパチップ サイズ 2) 500個	¥20,800	¥22,700	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
60872-310	Square Petri Dish with Grid ST 500	¥217,000	¥233,400	
60872-488	Disposable Base Molds 1000	¥121,200	¥127,200	
60872-492	Biopsy Foam Pads 1000枚入	¥26,300	¥27,600	
60872-628	PCRチューブ(ト-ムキャップ付)0.2ml 1000本	¥23,700	¥25,500	
60872-688	Rack with 96 Plain Individual Tube	¥48,500	¥52,100	
60872-716	1.2ml Tube Polypropylene 4800本入	¥110,700	¥119,100	
60872-736	8連PCRチューブ(キャップ付)0.2ml 125本	¥82,500	¥88,800	
60872-738	8連PCRチューブ(キャップ付)0.2ml 125本	¥82,500	¥88,800	
60872-742	8連PCRチューブ(キャップ付)0.2ml 125本	¥82,500	¥88,800	
60872-744	8連PCRチューブ(キャップ付)0.2ml 125本	¥82,500	¥88,800	
60882-178	ホ-リスレンパ-イアル	¥68,600	¥75,300	
60882-189	ホ-リスレンパ-イアル	¥58,500	¥64,200	
60882-205	Plastic Foam Plugs 50	¥39,100	¥42,900	
60882-207	Plastic Foam Plugs 80	¥55,400	¥60,800	
60914-210	トル-イキハ-ター用)	¥38,500	¥42,000	
60915-404	20mlチューブ用ラック 24.6x10.4x6.4cm 1	¥13,500	¥14,500	
60915-410	20mlチューブ用ラック 24.6x10.4x6.4cm 1	¥19,500	¥21,300	
60915-422	No-Wire Racks 21本立 1個入	¥20,300	¥22,400	
60916-260	チューブキャップ	¥27,300	¥29,500	
60917-202	ミチューブ96穴ラック	¥236,000	¥257,600	
60917-207	チューブキャップ (PS製/滅菌済み個別包装)	¥70,400	¥76,100	
60917-250	電気泳動用プレートラック	¥74,800	¥80,100	
60917-462	チューブ用ラック 24穴	¥33,900	¥36,200	
60941-062	VWR シーリングフィルム(Plate Film)100枚	¥35,000	¥37,700	
60941-064	VWR シーリングフィルム 14.6x7.9cm 滅菌済	¥47,700	¥51,300	
60941-066	VWR シーリングフィルム(Strip Film)200枚	¥27,500	¥29,600	
60941-068	VWR シーリングフィルム(StripFilm)ST200枚	¥36,800	¥39,600	
60941-084	VWR Rayon Films 100枚	¥56,300	¥60,500	
60941-088	VWR Pierceable Polyethylene Films	¥37,500	¥40,300	
60941-104	VWR シーリングフィルム(ア-ト/マ-キカ-リア付)	¥26,400	¥28,400	
60941-118	VWR Plate Roller 1個	¥23,000	¥27,200	
60985-454	チューブ用ラック 10-13mm 72本立	¥16,500	¥17,600	
60985-476	チューブ用ラック 13-16mm 40本立	¥15,400	¥16,400	
60985-498	チューブ用ラック 16-20mm 40本立	¥15,500	¥16,600	
60986-340	チューブ用ラック(丸型 TPX製) 1	¥12,500	¥13,300	
60986-342	チューブ用ラック(丸型 TPX製) 加-プレート無	¥13,900	¥14,800	
60986-358	TPXラック BBL RND 20穴 1個	¥10,700	¥11,400	
60986-402	チューブ用ラック(ホ-リスレン製) PS製 1個	¥80,000	¥85,700	
60986-427	チューブ用ラック(ホ-リスレン製)	¥80,800	¥86,500	
60995-955	乾燥ラック(卓上・壁掛け兼用型)	¥76,000	¥81,300	
60996-009	乾燥ラック(ワ-型)	¥201,000	¥217,100	
60996-012	乾燥ラック(軽量型)	¥292,600	¥322,400	
61000-110	Pikovskayas Agar 500g	¥74,200	¥84,100	
61161-036	デジタル温度計(テ-イアル)	¥32,500	¥34,900	
61161-319	乾燥剤	¥45,800	¥51,100	
61161-321	乾燥剤	¥46,200	¥50,200	
61161-364	Refrigerator/Freezer Thermometer	¥45,400	¥50,000	
611A46	Blood DNA Auto Plate 96回用	¥69,900	¥71,900	
611CDS1	ア-ゲル-ア-リカム611用 リニテ-ションキット 1ヶ入	¥6,400	¥6,600	
611CEL1	ア-ゲル-ア-リカム611用 UVラップ	¥39,500	¥40,300	
611S46	Blood DNA Auto Tube 96回用	¥69,900	¥71,900	
61220-601	Digital Thermometer with Probe	¥432,500	¥805,000	
61222-686	温度計保護ケース	¥11,400	¥12,200	
61222-702	温度計保護ケース	¥23,900	¥25,600	
61279-934	温度記録計(温度幅5種 選択可能)	¥560,400	¥1,063,000	
61280-105	温度記録計	¥36,700	¥69,000	
61280-132	Red Pens(温度記録計用) 6	¥47,400	¥58,100	
61280-321	温度記録計	¥41,600	¥78,000	
61280-323	温度記録計	¥41,600	¥78,000	
61280-338	テ-加カ-(温度計)	¥31,100	¥60,000	
61280-340	テ-加カ-(温度計)	¥35,600	¥41,900	
612A46	Tissue DNA Auto Plate 96回用	¥49,800	¥51,300	
612S46	Tissue DNA Auto Tube 96回用	¥49,800	¥51,300	
613A46	Plant DNA Auto Plate 96回用	¥46,600	¥47,900	
613A46-SE	Plant DNA Auto Plate 96回用	¥50,500	¥51,900	
613S46	Plant DNA Auto Tube 96回用	¥46,600	¥47,900	
613S46-SE	Plant DNA Auto Tube 96回用	¥50,500	¥51,900	
619A45	Rice DNA Auto Plate 48回用	¥31,100	¥31,900	
619S45	Rice DNA Auto Tube 48回用	¥31,100	¥31,900	
61CA46	OptiPure cfDNA Auto Plate 96回用	¥136,000	¥139,800	
61CS46	OptiPure cfDNA Auto Tube 96回用	¥136,000	¥139,800	
61EA46	OptiPure BloodDNA AutoPlate 96回用	¥50,500	¥51,900	
61EA46-BS	Dried Blood Spot DNA Plate 96回用	¥73,800	¥75,900	
61ES46	OptiPure BloodDNA AutoTube 96回用	¥50,500	¥51,900	
61ES46-BS	Dried Blood Spot DNA Tube 96回用	¥73,800	¥75,900	
61FA46	TANBead Fungi DNA AutoPlate 96回用	¥46,600	¥47,900	
61FS46	TANBead Fungi DNA Auto Tube 96回用	¥46,600	¥47,900	
61GA46	Gram Bacteria DNA AutoPlate 96回用	¥50,500	¥51,900	
61GS46	Gram Bacteria DNA Auto Tube 96回用	¥50,500	¥51,900	
61HA46	HPV DNA Auto Plate 96回用	¥62,200	¥63,900	
61HS46	HPV DNA Auto Tube 96回用	¥62,200	¥63,900	
61PA46	FFPE DNA Auto Plate 96回用	¥73,800	¥75,900	
61PS46	FFPE DNA Auto Tube 96回用	¥73,800	¥75,900	
622A46	OptiPure Leukocyte RNA Auto Plate	¥111,900	¥115,000	
622S46	OptiPure Leukocyte RNA Auto Tube	¥111,900	¥115,000	
62362-000	ホ-リスレン(ア-ゲル)	¥114,800	¥177,000	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
62363-003	ラベタイマー(アナログ)60 Minute	¥112,600	¥206,000	
62364-185	タイマー Interval Clock 5 Channel	¥43,900	¥47,400	
62371-040	ラベタイマー	¥342,900	¥864,000	
62373-007	タイマー(暗室用)	¥437,900	¥1,002,000	
62373-109	ラベタイマー	¥361,600	¥859,000	
62378-204	Wall Clocks	¥28,400	¥52,000	
62400-595	Dounce Tissue Glinders 1mL 2個入	¥167,900	¥183,100	
62401-966	Jacketed Celstir(ガラス製)75ml	¥479,700	¥523,200	
62404-012	Head Friction Drive TT RD65	¥162,800	¥180,200	
62404-014	Head Twin-Disk Test-Tube RD10	¥495,100	¥548,300	
62546-048	Enclosed Variable Transformer	¥472,800	¥950,000	
62660-348	取納ホックス	¥34,600	¥37,100	
62660-508	取納ホックス	¥67,000	¥71,700	
62660-585	取納ホックス	¥119,100	¥127,500	
62660-767	取納ホックス	¥129,000	¥138,200	
62660-790	取納ホックス	¥88,700	¥94,900	
62660-880	取納ホックス	¥61,000	¥65,300	
62660-916	取納ホックス	¥74,100	¥79,300	
62660-949	取納ホックス	¥26,200	¥28,100	
62661-455	取納トレ 24個入	¥62,100	¥74,600	
62661-458	取納トレ	¥61,900	¥74,400	
62661-464	取納トレ	¥64,500	¥77,500	
62661-466	取納トレ	¥99,200	¥119,200	
62661-476	取納トレ	¥11,400	¥13,700	
62661-478	取納トレ	¥12,300	¥14,700	
62661-480	取納トレ	¥21,100	¥25,300	
62661-490	取納トレ	¥117,900	¥76,900	
62661-498	取納トレ	¥156,400	¥169,000	
62661-504	取納トレ	¥5,600	¥1,900	
62661-506	取納トレ	¥10,000	¥3,600	
62661-508	取納トレ	¥27,700	¥9,000	
62661-522	取納トレ	¥58,100	¥69,800	
62661-524	取納トレ	¥58,200	¥69,900	
62661-532	取納トレ	¥10,700	¥3,600	
62661-534	取納トレ	¥10,600	¥3,600	
62662-401	トレ(ホリフビ製)	¥25,400	¥29,000	
62662-423	トレ(ホリフビ製)	¥43,200	¥48,500	
62662-445	トレ(ホリフビ製)	¥58,400	¥63,000	
62662-489	トレ(ホリフビ製)	¥98,500	¥105,500	
62662-503	トレ(ホリフビ製)	¥115,300	¥123,400	
62662-966	トレ(ホリフビ製)	¥86,200	¥111,400	
62663-186	トレ(ホリフビ製)加付	¥20,500	¥22,000	
62663-222	トレ(ホリフビ製)加付	¥37,800	¥40,500	
62663-266	トレ(ホリフビ製)加付	¥85,300	¥91,300	
62687-050	トレ(ステンレス製)	¥36,600	¥45,100	
62687-069	トレ(ステンレス製)	¥36,100	¥42,400	
62687-070	トレ(ステンレス製)	¥58,000	¥69,400	
62860-020	Nalgene Tubing Connectors 12個入	¥49,900	¥56,100	
62866-107	マイクロレスタクト(角度調節可)	¥123,400	¥131,500	
62866-140	加ニクシリタ Assorted 3サイズx各40	¥76,200	¥81,600	
62866-176	Self Zeroing Burette 25ml	¥184,400	¥197,400	
62866-178	Self Zeroing Burette 10mL 1	¥161,600	¥173,000	
62866-192	ハトリテ イソ用ラック	¥67,300	¥72,100	
62866-206	ハトリテ イソ用ラック	¥26,900	¥29,400	
62866-550	Centrifuge Tubes Rimless 90mL PK12	¥40,100	¥45,600	
62991-905	コムチーブ	¥282,600	¥296,600	
62993-014	ネオレチーブ	¥33,400	¥31,400	
62993-018	ネオレチーブ	¥35,400	¥33,300	
62993-019	ネオレチーブ	¥55,900	¥60,400	
62993-023	ネオレチーブ	¥38,200	¥35,100	
62993-024	ネオレチーブ	¥59,000	¥54,300	
62993-057	ネオレチーブ	¥72,200	¥78,100	
62993-058	ネオレチーブ	¥106,500	¥115,200	
62993-059	ネオレチーブ	¥183,700	¥198,600	
62993-061	ネオレチーブ	¥44,800	¥41,100	
62993-093	ネオレチーブ	¥50,600	¥54,700	
62993-114	ネオレチーブ	¥101,200	¥109,400	
62993-115	ネオレチーブ	¥139,400	¥128,200	
62993-118	ネオレチーブ	¥94,900	¥102,600	
62993-148	ネオレチーブ	¥171,600	¥185,500	
62993-163	ネオレチーブ	¥195,900	¥211,700	
62993-172	ネオレチーブ	¥253,800	¥274,400	
62993-690	コムチーブ(耐圧/耐真空)	¥279,600	¥293,500	
62993-726	コムチーブ(耐圧/耐真空)	¥167,000	¥175,200	
62993-759	コムチーブ(耐圧/耐真空)	¥227,100	¥238,300	
62993-781	コムチーブ(耐圧/耐真空)	¥265,200	¥278,400	
62993-817	コムチーブ(耐圧/耐真空)	¥209,100	¥219,400	
62993-840	コムチーブ(耐圧/耐真空)	¥523,700	¥549,600	
62993-908	コムチーブ(耐圧/耐真空)	¥541,500	¥568,300	
62994-219	コムチーブ(耐圧/耐真空)	¥60,200	¥63,200	
62994-241	コムチーブ(耐圧/耐真空)	¥85,700	¥90,000	
62995-059	コムチーブ(耐圧/耐真空)	¥116,900	¥122,700	
62995-062	コムチーブ(耐圧/耐真空)	¥171,600	¥180,100	
62995-095	コムチーブ(耐圧/耐真空)	¥670,900	¥704,100	
62995-120	コムチーブ(耐圧/耐真空)	¥630,300	¥661,500	
62995-153	コムチーブ(耐圧/耐真空)	¥705,100	¥740,100	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
62995-335	コトフチ-7 (耐圧/耐真空)	¥40,900	¥43,000	
62995-368	コトフチ-7 (耐圧/耐真空)	¥122,000	¥128,100	
62995-390	コトフチ-7 (耐圧/耐真空)	¥174,300	¥182,900	
62995-426	コトフチ-7 (耐圧/耐真空)	¥189,000	¥198,300	
62995-608	コトフチ-7	¥23,900	¥25,100	
62995-618	コトフチ-7	¥26,700	¥28,000	
62995-630	コトフチ-7	¥37,800	¥39,600	
62995-663	コトフチ-7	¥41,300	¥43,400	
62995-670	コトフチ-7	¥55,800	¥58,600	
62995-680	コトフチ-7	¥73,200	¥76,800	
62995-696	コトフチ-7	¥45,200	¥47,400	
62995-721	コトフチ-7	¥92,100	¥96,600	
62996-145	コトフチ-7	¥106,900	¥112,200	
62996-178	コトフチ-7	¥73,000	¥76,600	
62996-203	コトフチ-7	¥141,600	¥148,600	
62996-350	コトフチ-7	¥23,700	¥24,900	
62996-382	コトフチ-7	¥35,900	¥37,700	
62996-418	コトフチ-7	¥35,500	¥37,300	
62996-440	コトフチ-7	¥49,900	¥52,400	
62996-473	コトフチ-7	¥40,400	¥42,400	
62996-509	コトフチ-7	¥56,600	¥59,400	
62996-531	コトフチ-7	¥70,400	¥73,900	
62996-564	コトフチ-7	¥45,200	¥47,400	
62996-597	コトフチ-7	¥64,100	¥67,300	
62996-622	コトフチ-7	¥64,400	¥67,600	
62996-655	コトフチ-7	¥57,800	¥60,600	
62996-666	コトフチ-7	¥92,900	¥97,500	
62996-688	コトフチ-7	¥110,700	¥116,100	
62996-699	コトフチ-7	¥158,900	¥166,700	
62999-371	シリコンチューブ	¥69,900	¥71,300	
62999-406	シリコンチューブ	¥84,300	¥88,500	
62999-439	シリコンチューブ 3/16x3/8	¥115,700	¥118,000	
62999-530	シリコンチューブ	¥154,700	¥162,400	
62999-850	シリコンチューブ	¥65,100	¥68,300	
62999-852	シリコンチューブ 0.020x0.037	¥165,100	¥173,300	
62999-854	シリコンチューブ 0.025x0.047	¥138,100	¥145,000	
62999-856	シリコンチューブ	¥156,600	¥164,400	
62999-860	シリコンチューブ	¥67,400	¥70,700	
62999-862	シリコンチューブ	¥66,200	¥69,500	
62999-864	シリコンチューブ	¥117,000	¥122,800	
62999-866	シリコンチューブ	¥136,000	¥142,800	
62999-868	シリコンチューブ	¥104,200	¥109,400	
62999-870	シリコンチューブ	¥87,500	¥91,800	
62999-874	シリコンチューブ	¥129,000	¥135,400	
62999-878	シリコンチューブ	¥194,000	¥203,600	
62999-880	シリコンチューブ	¥231,200	¥242,600	
62999-884	シリコンチューブ	¥338,500	¥355,200	
62999-886	シリコンチューブ	¥383,500	¥402,500	
62999-888	シリコンチューブ	¥342,300	¥359,300	
62999-890	シリコンチューブ	¥417,900	¥438,600	
62999-892	シリコンチューブ	¥424,300	¥445,300	
62999-896	シリコンチューブ 0.625x0.875/0.125in 1	¥527,700	¥553,900	
62999-900	シリコンチューブ	¥141,700	¥148,800	
62999-902	シリコンチューブ	¥128,900	¥135,300	
62999-904	シリコンチューブ	¥202,300	¥212,400	
62999-906	シリコンチューブ	¥241,600	¥253,600	
62999-908	シリコンチューブ	¥381,100	¥400,000	
62999-910	シリコンチューブ	¥436,600	¥458,200	
62999-912	シリコンチューブ	¥519,000	¥544,700	
62999-914	シリコンチューブ	¥553,600	¥581,100	
62999-926	シリコンチューブ 0.125x0.250	¥78,400	¥82,300	
62999-934	HelixMark シリコンチューブ 50フィート/巻	¥236,600	¥248,400	
62999-952	シリコンチューブ	¥615,400	¥645,900	
62999-954	シリコンチューブ 50フィート	¥130,000	¥136,400	
62999-956	シリコンチューブ	¥189,800	¥199,200	
62999-958	シリコンチューブ	¥314,700	¥330,300	
62999-960	シリコンチューブ	¥528,900	¥555,100	
62999-962	シリコンチューブ	¥840,900	¥882,600	
62999-964	シリコンチューブ 0.030x0.065	¥56,300	¥59,100	
62999-990	シリコンチューブ 0.012x0.025	¥53,200	¥55,900	
63000-122	シリコンチューブ	¥140,600	¥147,600	
63000-124	シリコンチューブ	¥104,100	¥109,200	
63000-126	シリコンチューブ	¥109,000	¥114,400	
63000-128	シリコンチューブ	¥122,200	¥128,300	
63000-130	シリコンチューブ	¥186,300	¥195,600	
63000-134	シリコンチューブ	¥237,200	¥248,900	
63000-136	シリコンチューブ	¥271,700	¥285,200	
63000-138	シリコンチューブ	¥380,300	¥399,100	
63007-319	コトフチ-7 (高品質)	¥88,900	¥96,100	
63007-329	コトフチ-7 (高品質)	¥107,300	¥116,000	
63007-900	滅菌用バフ(蒸気/ガス滅菌兼用 0.7kg)	¥9,100	¥9,500	
63007-902	滅菌用バフ(蒸気/ガス滅菌兼用 0.7kg)	¥14,600	¥15,300	
63007-904	滅菌用バフ(蒸気/ガス滅菌兼用 0.7kg)	¥16,100	¥16,900	
63007-906	滅菌用バフ(蒸気/ガス滅菌兼用 0.7kg)	¥14,100	¥14,800	
63007-908	滅菌用バフ(蒸気/ガス滅菌兼用) 1	¥17,900	¥18,800	
63007-910	滅菌用バフ(蒸気/ガス滅菌兼用 0.7kg)	¥124,000	¥130,200	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
63008-942	ホリカタチューブ	¥266,100	¥287,700	
63008-944	ホリカタチューブ	¥441,000	¥476,800	
63008-946	ホリカタチューブ	¥461,000	¥498,400	
63008-950	ホリカタチューブ	¥691,800	¥747,900	
63009-252	タロンチューブ 3350	¥75,300	¥81,400	
63009-254	タロンチューブ 3350	¥85,300	¥92,300	
63009-256	タロンチューブ 3350	¥116,700	¥126,100	
63009-258	タロンチューブ 3350	¥94,800	¥102,500	
63009-260	タロンチューブ 3350	¥134,700	¥145,600	
63009-262	タロンチューブ 3350	¥98,800	¥106,800	
63009-265	タロンチューブ 3350	¥138,000	¥149,200	
63009-267	タロンチューブ 3350	¥103,900	¥112,300	
63009-269	タロンチューブ 3350	¥103,600	¥112,000	
63009-271	タロンチューブ 3350	¥198,500	¥214,600	
63009-273	タロンチューブ 3350	¥284,800	¥307,900	
63009-275	タロンチューブ 3350	¥372,100	¥402,300	
63009-277	タロンチューブ 3350	¥163,000	¥176,200	
63009-279	タロンチューブ 3350	¥208,900	¥225,900	
63009-281	タロンチューブ 3350	¥321,300	¥347,300	
63009-283	タロンチューブ 3350	¥452,300	¥489,000	
63009-285	タロンチューブ 3350	¥252,000	¥272,400	
63009-287	タロンチューブ 3350	¥370,000	¥400,000	
63009-289	タロンチューブ 3350	¥287,700	¥311,000	
63009-291	タロンチューブ 3350	¥427,600	¥462,300	
63009-293	タロンチューブ 3350	¥492,200	¥532,100	
63009-297	タロンチューブ 3350	¥463,600	¥501,200	
63009-299	タロンチューブ 3350	¥379,800	¥410,600	
63009-301	タロンチューブ 3350	¥525,900	¥568,500	
63009-303	タロンチューブ 3350	¥698,100	¥754,700	
63009-307	タロンチューブ 3350	¥651,900	¥704,700	
63009-309	タロンチューブ 3350	¥915,200	¥989,400	
63009-311	タロンチューブ 3350	¥923,600	¥998,500	
63009-313	タロンチューブ 3350	¥1,097,500	¥1,186,500	
63010-224	タロンチューブ SE-200	¥205,500	¥222,100	
63010-406	タロンチューブ F-4040-A	¥54,700	¥59,200	
63010-408	Glove Box Glove 30ミル 非製	¥293,800	¥325,400	
63010-417	タロンチューブ F-4040-A	¥72,100	¥77,900	
63010-428	タロンチューブ F-4040-A	¥101,400	¥113,900	
63010-436	タロンチューブ F-4040-A	¥105,900	¥114,500	
63010-440	タロンチューブ F-4040-A	¥125,100	¥135,300	
63010-450	タロンチューブ F-4040-A	¥154,000	¥166,500	
63010-461	タロンチューブ F-4040-A	¥310,900	¥336,100	
63010-468	タロンチューブ F-4040-A	¥360,200	¥389,400	
63014-305	プラスチックチューブ (耐腐食性)	¥123,800	¥133,800	
63014-382	プラスチックチューブ (耐腐食性)	¥720,100	¥778,500	
63014-407	プラスチックチューブ (耐腐食性)	¥875,200	¥946,100	
63014-422	タロンチューブ 3370 I.B.	¥870,700	¥941,300	
63014-424	タロンチューブ 3370 I.B.	¥1,001,500	¥1,082,700	
63014-426	タロンチューブ 3370 I.B.	¥1,482,900	¥1,603,000	
63014-428	タロンチューブ 3370 I.B.	¥1,625,400	¥1,757,200	
63014-432	タロンチューブ 3370 I.B.	¥2,298,800	¥2,485,100	
63014-690	珪素樹脂チューブ (ケミコ)	¥45,800	¥50,200	
63014-695	珪素樹脂チューブ (ケミコ)	¥115,300	¥124,600	
63014-696	珪素樹脂チューブ (ケミコ)	¥282,400	¥305,300	
63014-699	珪素樹脂チューブ (ケミコ)	¥438,800	¥474,300	
63014-800	珪素樹脂チューブ (ケミコ)	¥101,500	¥109,800	
63014-805	珪素樹脂チューブ (ケミコ)	¥328,600	¥355,300	
63014-809	珪素樹脂チューブ (ケミコ)	¥220,000	¥237,800	
63014-815	珪素樹脂チューブ (ケミコ)	¥362,600	¥392,000	
63014-850	珪素樹脂チューブ (ケミコ)	¥296,000	¥320,000	
63014-860	珪素樹脂チューブ (ケミコ)	¥850,600	¥919,500	
63014-861	珪素樹脂チューブ (ケミコ)	¥957,200	¥1,034,800	
63014-900	珪素樹脂チューブ (ケミコ)	¥54,100	¥58,500	
63014-905	珪素樹脂チューブ (ケミコ)	¥222,400	¥240,400	
63014-909	珪素樹脂チューブ (ケミコ)	¥104,400	¥112,900	
63014-914	珪素樹脂チューブ (ケミコ)	¥306,000	¥330,800	
63014-915	珪素樹脂チューブ (ケミコ)	¥349,500	¥377,800	
63014-925	珪素樹脂チューブ (ケミコ)	¥548,000	¥592,400	
63014-934	珪素樹脂チューブ (ケミコ)	¥707,300	¥764,700	
63014-950	珪素樹脂チューブ (ケミコ)	¥165,000	¥178,300	
63014-960	珪素樹脂チューブ (ケミコ)	¥173,100	¥187,100	
63014-961	珪素樹脂チューブ (ケミコ)	¥230,300	¥248,900	
63022-100	チューブキャップ	¥23,400	¥26,500	
63022-102	チューブキャップ用替刃 24	¥265,400	¥284,200	
63022-150	VACUUM PRESSURE INDICATOR	¥37,800	¥44,900	
63022-403	チューブクランプ (オートクレーブ可/シングル)	¥30,800	¥32,700	
63022-405	チューブクランプ (オートクレーブ可/15ホ°シヨソ)	¥12,900	¥13,600	
63022-407	チューブクランプ (オートクレーブ可/12ホ°シヨソ)	¥21,200	¥24,900	
63300-102	TruCLEAN Polyfoam Sponges 25個	¥118,800	¥127,200	
63300-105	TruCLEAN ホリカタカラム用 25入	¥58,100	¥62,200	
63395-005	Lab-Grade Dewar Flasks 350ml	¥135,500	¥152,400	
63397-000	Lab-Grade Dewar Flasks	¥165,200	¥186,800	
63400-114	1000mL Powder Flask (IIRホ°用)	¥177,000	¥322,000	
63400-116	2000mL Powder Flask (IIRホ°用)	¥216,600	¥394,000	
63400-118	500mL Powder Flask (IIRホ°用)	¥162,600	¥296,000	
63405-109	液体窒素用テフロン瓶 (ロケツの球状)	¥1,450,500	¥2,754,000	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
63415-000	サモワスコ(ステン製)	¥427,200	¥490,200	
63415-033	サモワスコ(ステン製)	¥454,500	¥498,900	
63422-028	液体窒素用テラー瓶 350ml	¥140,900	¥156,700	
63422-040	液体窒素用テラー瓶 665ml	¥141,600	¥157,500	
63422-061	液体窒素用テラー瓶 1000ml	¥157,000	¥174,700	
63422-108	液体窒素用テラー瓶 1900ml	¥373,100	¥415,100	
63422-141	液体窒素用テラー瓶 4300ml	¥446,900	¥497,200	
63550-055	GLOVES BUTYL 14IN 7MIL 1双	¥36,600	¥40,500	
63550-057	GLOVES BUTYL 14IN 7MIL 1双	¥36,600	¥40,500	
635A46	Viral Auto Plate 96回用	¥50,500	¥51,900	
635S46	Viral Auto Tube 96回用	¥50,500	¥51,900	
63884-961	液体窒素用テラー瓶(浅型真空)	¥202,700	¥225,600	
63884-963	液体窒素用テラー瓶(浅型真空) 275mL	¥257,300	¥285,800	
63884-965	液体窒素用テラー瓶(浅型真空) 600mL	¥296,100	¥329,400	
63884-969	液体窒素用テラー瓶(浅型真空)	¥303,100	¥337,600	
63884-971	液体窒素用テラー瓶(浅型真空)	¥278,900	¥310,300	
63884-973	液体窒素用テラー瓶(浅型真空)	¥395,300	¥439,300	
63884-975	液体窒素用テラー瓶(浅型真空)	¥496,300	¥928,000	
65000-000	スリッパ - (ホリア®製)	¥137,700	¥149,000	
65000-002	スリッパ - (ホリア®製)	¥104,600	¥112,300	
65000-004	スリッパ - (ホリア®製)	¥112,200	¥120,100	
65000-006	スリッパ - (ホリア®製)	¥73,900	¥79,200	
65000-008	スリッパ - (ホリア®製)	¥124,100	¥133,200	
65000-010	スリッパ - (ホリア®製)	¥85,300	¥91,400	
65000-012	スリッパ - (ホリア®製)	¥79,200	¥84,800	
65-54001	簡易リムーバー(8連) 1	¥7,100	¥5,300	
66001-104	PESTLE MIXER W/MATCH 1.5MLTU	¥97,800	¥104,900	
66001-148	洗浄瓶(広口) Methanol 500ml 3本入	¥17,100	¥18,300	
66001-152	洗浄瓶(広口)	¥17,100	¥18,300	
66001-200	Cryo-Safe Cold Box(加付)12本立	¥58,300	¥62,400	
66001-236	オートサンワラー用ハイル	¥32,500	¥34,800	
66001-239	オートサンワラー用ハイル	¥32,200	¥34,400	
66001-240	オートサンワラー用ハイル	¥30,900	¥33,100	
66001-242	オートサンワラー用ハイル	¥29,800	¥31,900	
66004-004	STOPPER BUTYL/PTFE 20MM 100個入	¥101,400	¥110,600	
66004-006	STOPPER BUTYL STRAIGHT 20MM 1000個	¥250,600	¥273,300	
66004-008	ガラスバイアル	¥236,500	¥257,900	
66008-783	Cryogenic Vial Holder 50-Position	¥43,900	¥47,400	
66009-453	スロアプレイク・アンプル	¥41,200	¥44,900	
66009-454	スロアプレイク・アンプル	¥76,100	¥83,000	
66009-458	スロアプレイク・アンプル(クリアガラス) 5mL 144本	¥94,800	¥103,300	
66009-460	スロアプレイク・アンプル	¥39,300	¥42,900	
66009-462	スロアプレイク・アンプル	¥75,700	¥82,500	
66009-546	ガラスバイアル	¥142,600	¥157,900	
66009-548	ガラスバイアル	¥140,700	¥155,800	
66009-550	ガラスバイアル	¥55,000	¥60,900	
66009-551	ガラスバイアル	¥48,700	¥53,900	
66009-553	ガラスバイアル	¥65,900	¥73,000	
66009-554	ガラスバイアル	¥70,800	¥78,300	
66009-555	ガラスバイアル	¥107,000	¥118,500	
66009-556	ガラスバイアル	¥112,700	¥124,800	
66009-557	ガラスバイアル	¥104,500	¥115,800	
66009-558	ガラスバイアル	¥108,200	¥119,800	
66009-559	ガラスバイアル	¥129,600	¥143,500	
66009-560	ガラスバイアル	¥130,000	¥144,000	
66009-561	ガラスバイアル	¥169,300	¥187,500	
66009-562	ガラスバイアル	¥88,200	¥97,700	
66009-563	ガラスバイアル	¥107,300	¥118,800	
66009-565	ガラスバイアル	¥101,500	¥112,400	
66009-566	ガラスバイアル	¥49,400	¥54,600	
66009-567	ガラスバイアル	¥59,900	¥66,300	
66009-575	ガラスバイアル	¥123,700	¥137,000	
66009-870	オートサンワラー用ハイルキット(広口)	¥238,900	¥194,300	
66009-882	オートサンワラー用ハイルキット	¥31,000	¥33,300	
66009-902	オートサンワラー用ハイル 200本	¥63,800	¥68,600	
66009-924	オートサンワラー用ハイル 200本	¥47,400	¥50,900	
66009-968	オートサンワラー用ハイル	¥348,800	¥374,900	
66009-981	ガラスバイアル 4mL	¥130,300	¥142,100	
66009-984	ガラスバイアル 8mL	¥194,200	¥211,800	
66009-985	ガラスバイアル 12ml 19x67mm 144本	¥217,200	¥236,900	
66009-986	ガラスバイアル 16mL	¥246,900	¥269,200	
66009-987	ガラスバイアル 24mL	¥288,200	¥314,200	
66009-988	ガラスバイアル 2mL	¥210,600	¥229,600	
66009-990	オートサンワラー用ハイル	¥417,600	¥455,400	
66010-243	ガラスバイアル	¥905,200	¥987,100	
66010-247	ガラスバイアル	¥744,300	¥811,600	
66010-249	ガラスバイアル	¥626,700	¥683,400	
66010-251	ガラスバイアル	¥558,700	¥609,300	
66010-253	ガラスバイアル	¥550,400	¥600,200	
66010-255	ガラスバイアル	¥696,200	¥759,300	
66010-260	ガラスバイアル	¥875,100	¥954,300	
66010-265	ガラスバイアル	¥794,400	¥866,400	
66010-267	ガラスバイアル	¥823,300	¥897,800	
66010-268	ガラスバイアル	¥579,800	¥632,300	
66010-270	ガラスバイアル	¥486,500	¥530,600	
66010-272	ガラスバイアル	¥561,300	¥612,100	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
66010-274	ガラスバール	¥451,400	¥492,300	
66010-276	ガラスバール	¥344,900	¥376,100	
66010-286	ガラスバール	¥217,300	¥236,900	
66010-300	ガラスバール	¥216,900	¥236,500	
66010-322	ガラスバール	¥249,000	¥271,600	
66010-408	ガラスバール	¥669,600	¥730,100	
66010-410	ガラスバール	¥420,400	¥458,400	
66010-423	ガラスバール	¥142,500	¥155,400	
66010-425	ガラスバール	¥99,600	¥108,600	
66010-426	ガラスバール	¥129,800	¥141,600	
66010-427	ガラスバール	¥128,400	¥140,000	
66010-428	ガラスバール	¥161,200	¥175,800	
66010-429	ガラスバール	¥203,100	¥221,500	
66010-437	サンプルバール(キャップ付)11drams 432本	¥591,500	¥645,000	
66010-442	ガラスバール	¥808,400	¥881,600	
66010-445	ガラスバール	¥656,100	¥715,500	
66010-446	ガラスバール	¥261,700	¥285,400	
66010-448	ガラスバール	¥273,100	¥297,900	
66010-449	ガラスバール	¥502,300	¥547,800	
66010-451	ガラスバール	¥586,300	¥639,300	
66010-452	ガラスバール	¥287,300	¥313,300	
66010-482	オートサンワラー用バール	¥489,900	¥534,200	
66010-484	オートサンワラー用バール	¥575,200	¥627,200	
66010-490	オートサンワラー用バール	¥123,400	¥134,500	
66010-491	オートサンワラー用バール	¥127,000	¥138,500	
66010-492	オートサンワラー用バール	¥129,200	¥140,900	
66010-493	オートサンワラー用バール	¥141,500	¥154,300	
66010-494	オートサンワラー用バール	¥161,800	¥176,400	
66010-495	オートサンワラー用バール	¥209,200	¥228,100	
66010-496	オートサンワラー用バール	¥226,800	¥247,300	
66010-497	オートサンワラー用バール	¥228,600	¥249,300	
66010-498	オートサンワラー用バール	¥250,700	¥273,400	
66010-522	ガラスバール	¥987,800	¥1,077,300	
66010-560	オートサンワラー用バール	¥304,000	¥331,500	
66010-564	オートサンワラー用バール	¥364,600	¥397,700	
66010-570	オートサンワラー用バール	¥424,300	¥462,700	
66010-581	オートサンワラー用バール	¥496,500	¥541,500	
66010-642	オートサンワラー用バール	¥100,700	¥108,300	
66010-672	バールキャップ	¥18,800	¥20,600	
66010-683	バールキャップ	¥23,000	¥25,100	
66010-694	バールキャップ	¥68,100	¥74,300	
66010-720	バールシール	¥51,200	¥55,900	
66010-730	バールシール	¥59,100	¥64,400	
66010-741	バールシール	¥61,200	¥66,800	
66010-751	バールシール	¥112,900	¥123,100	
66010-774	バールシール	¥51,900	¥56,600	
66010-785	バールシール	¥68,400	¥74,600	
66010-795	ガラスバール	¥74,500	¥81,200	
66010-796	バールシール	¥92,000	¥100,300	
66010-797	ガラスバール	¥76,100	¥82,900	
66010-798	ガラスバール	¥65,400	¥71,300	
66010-799	ガラスバール	¥68,000	¥74,200	
66010-865	Gray Butyl Rubber Stoppers	¥119,000	¥129,800	
66010-887	バールキャップ	¥27,400	¥29,900	
66010-890	オートサンワラー用バール	¥125,800	¥137,200	
66010-894	オートサンワラー用バール	¥149,700	¥162,600	
66010-912	バールシール	¥69,400	¥75,700	
66010-946	バールシール	¥231,500	¥252,500	
66010-947	バールシール	¥193,600	¥211,100	
66011-020	ガラスバール	¥31,600	¥34,300	
66011-041	ガラスバール	¥49,400	¥53,700	
66011-063	ガラスバール	¥27,600	¥30,000	
66011-100	ガラスバール	¥28,600	¥31,100	
66011-121	ガラスバール	¥59,100	¥64,200	
66011-143	ガラスバール	¥52,300	¥56,800	
66011-165	ガラスバール	¥41,200	¥44,800	
66011-188	バールシール 100	¥18,900	¥20,200	
66011-202	オートサンワラー用バールインサート	¥63,000	¥67,500	
66011-204	オートサンワラー用バールインサート	¥52,900	¥56,600	
66011-206	オートサンワラー用バール(クリア)4mL	¥19,700	¥21,100	
66011-256	サンプルバール(キャップ付)2ml 288本入	¥144,300	¥157,400	
66011-512	EC サンプルバール 24ml 200個入	¥96,300	¥105,000	
66012-022	ガラスバール	¥69,300	¥75,300	
66012-044	ガラスバール	¥43,800	¥47,600	
66012-066	ガラスバール	¥48,800	¥58,400	
66012-088	ガラスバール	¥49,100	¥53,300	
66012-140	ガラスバール	¥240,600	¥262,400	
66012-142	ガラスバール	¥268,000	¥292,300	
66012-146	ガラスバール	¥295,600	¥322,300	
66013-025	ガラスバール	¥40,300	¥43,900	
66013-069	ガラスバール	¥64,900	¥70,800	
66013-080	ガラスバール	¥78,800	¥85,900	
66013-105	ガラスバール	¥133,800	¥145,900	
66013-120	ガラスバール	¥502,500	¥547,900	
66013-122	ガラスバール	¥616,900	¥672,700	
66013-127	ガラスバール	¥345,700	¥362,800	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
66013-129	ガラスバイアル	¥637,900	¥695,700	
66013-131	ガラスバイアル	¥605,600	¥660,400	
66013-143	ガラスバイアル	¥714,300	¥778,900	
66013-150	ガラスバイアル	¥533,100	¥581,300	
66013-152	ガラスバイアル	¥600,700	¥655,000	
66013-156	ガラスバイアル	¥697,500	¥760,600	
66014-004	KIMAXメソイボトル 125mL 48本入	¥259,000	¥282,400	
66014-006	KIMAXメソイボトル 250mL 48本入	¥322,500	¥351,700	
66014-301	Vial EPA WAT ANYL 20ml 432本	¥298,500	¥325,500	
66014-303	ガラスバイアル	¥228,600	¥249,300	
66014-304	ガラスバイアル	¥292,000	¥318,500	
66014-346	ガラスバイアル 45ml without cap 432本	¥368,400	¥401,800	
66014-358	ガラスバイアル	¥59,500	¥64,900	
66014-359	ガラスバイアル	¥57,100	¥62,200	
66014-361	ガラスバイアル	¥100,600	¥120,700	
66014-364	ガラスバイアル	¥68,000	¥74,100	
66014-380	バイアルキャップ	¥20,300	¥22,200	
66015-450	スチロール棒瓶	¥47,000	¥50,500	
66015-483	スチロール棒瓶	¥47,200	¥50,800	
66015-519	スチロール棒瓶	¥44,800	¥55,300	
66015-600	スチロール棒瓶	¥79,100	¥85,100	
66015-632	スチロール棒瓶	¥189,100	¥203,500	
66017-366	ガラスバイアル	¥289,600	¥315,800	
66018-574	ガラスバイアル	¥126,200	¥139,700	
66018-665	ガラスバイアル	¥174,600	¥193,300	
66020-190	ガラスバイアル	¥40,600	¥44,300	
66020-191	ガラスバイアル	¥36,500	¥39,800	
66020-192	ガラスバイアル	¥44,800	¥48,800	
66020-194	ガラスバイアル	¥61,700	¥67,200	
66020-196	ガラスバイアル	¥67,900	¥74,000	
66020-197	ガラスバイアル	¥92,700	¥101,100	
66020-198	ガラスバイアル	¥145,500	¥158,700	
66020-257	液体シリンジョン用バイアル 4mL 1000本	¥218,400	¥239,300	
66020-800	オートサンワラー用バイアル	¥114,600	¥123,200	
66020-801	オートサンワラー用バイアル	¥100,500	¥108,100	
66020-802	オートサンワラー用バイアル	¥122,200	¥131,300	
66020-804	オートサンワラー用バイアル	¥135,300	¥145,400	
66020-806	オートサンワラー用バイアル	¥154,500	¥166,100	
66020-902	オートサンワラー用バイアル	¥158,800	¥170,600	
66020-903	オートサンワラー用バイアル	¥170,500	¥183,300	
66020-904	オートサンワラー用バイアル	¥162,700	¥156,100	
66020-905	オートサンワラー用バイアル	¥180,200	¥172,900	
66021-307	サンワールバイアル(キャップ付)6ml 1000本入	¥299,000	¥326,100	
66021-497	シリンジョンバイアル(スクリーキャップ付)20ml	¥224,700	¥245,100	
66021-978	クワイバールバイアル(アウターキャップ付/丸底/自立型)	¥285,900	¥298,600	
66021-983	Corning Brand 96well UV Plates 50	¥433,900	¥453,100	
66022-004	液体シリンジョン用バイアル 500	¥133,300	¥139,900	
66022-060	液体シリンジョン用バイアル	¥124,300	¥130,500	
66022-081	液体シリンジョン用バイアル	¥125,500	¥131,700	
66022-106	液体シリンジョン用バイアル	¥126,300	¥132,600	
66022-128	液体シリンジョン用バイアル	¥150,100	¥157,500	
66022-140	液体シリンジョン用バイアル	¥184,900	¥201,600	
66022-150	液体シリンジョン用バイアル	¥372,800	¥406,500	
66022-503	液体シリンジョン用バイアル用ラック	¥101,200	¥108,300	
66023-540	ラック (ワイヤ無し) 16mm	¥22,900	¥26,800	
66024-200	液体シリンジョン用バイアル	¥298,600	¥325,600	
66030-040	バイアルキャップ	¥17,700	¥19,000	
66030-042	バイアルキャップ	¥29,300	¥31,300	
66030-160	オートサンワラー用バイアルキット(広口)	¥174,400	¥186,700	
66030-392	バイアルキャップ	¥12,300	¥13,200	
66030-412	バイアルキャップ 100	¥9,900	¥10,600	
66030-420	Caps with Septa OpenTop 8-425 100	¥54,800	¥58,700	
66030-450	オートサンワラー用バイアルキット	¥74,200	¥79,400	
66030-478	オートサンワラー用バイアルキット	¥75,800	¥81,200	
66030-480	オートサンワラー用バイアルキット	¥78,000	¥83,500	
66030-698	オートサンワラー用バイアルキット	¥141,800	¥151,800	
66030-837	フィルター(水系シリンジョン用) 100	¥193,600	¥219,300	
66030-841	Target シリンジョン用フィルター 0.45um 30mm	¥193,100	¥208,700	
66030-847	Target Syringe Filter PVDF 0.45um 100	¥193,500	¥209,200	
66030-859	Target シリンジョン用フィルター 0.45um 17mm	¥134,400	¥145,300	
66030-863	Target シリンジョン用フィルター 0.45um 17mm	¥137,700	¥148,800	
66030-867	Target シリンジョン用フィルター 0.45um 17mm	¥151,500	¥163,800	
66030-870	フィルター(水系シリンジョン用)	¥136,300	¥147,400	
66030-872	フィルター(水系シリンジョン用)	¥135,100	¥146,000	
66030-874	フィルター(非水系シリンジョン用)	¥121,400	¥131,200	
66030-878	フィルター(非水系シリンジョン用)	¥122,000	¥131,900	
66030-928	バイアルシール	¥26,100	¥27,900	
66030-970	セパレーター(PreSlit PTFE/シリコン) 100	¥39,500	¥42,300	
66030-982	8-425 Screw-Thread Vial PP 100本	¥34,900	¥37,400	
66033-022	オートサンワラー用バイアルキット	¥46,300	¥49,600	
66064-096	Target DP9-425 スクリューバイアル用セパレーター	¥28,600	¥30,600	
66064-264	Conical Glass Insert	¥43,100	¥46,200	
66064-395	フィルター(水系シリンジョン用)5.0um/25mm	¥176,300	¥190,500	
66064-410	フィルター(水系/非水系兼用シリンジョン用)	¥198,700	¥216,900	
66064-412	フィルター(水系シリンジョン用)	¥135,100	¥146,000	
66064-414	フィルター(水系シリンジョン用)	¥136,500	¥147,600	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
66064-435	Target シリカフィルター 0.45um 17mm	¥152,300	¥164,600	
66064-450	Micro-Centrifuge Filters PVDF0.2um	¥168,700	¥182,400	
66064-468	25mL Centrifugal Filters PVDF 50	¥189,600	¥205,000	
66064-766	ルアーロックリッパ(オール スリット)10ml 100	¥77,600	¥83,900	
66064-770	フェイスシールド 3ML ルアーロック 100本	¥61,200	¥66,100	
66064-772	フェイスシールド 5ML ルアーロック 100本	¥69,600	¥75,200	
66064-823	Target シリカフィルター 0.45um 17mm	¥171,800	¥185,700	
66064-834	フィルター(水系シリカフィルター)	¥134,200	¥145,100	
66064-836	フィルター(水系シリカフィルター)	¥122,600	¥132,500	
66064-838	フィルター(水系/非水系兼用シリカフィルター)	¥118,600	¥128,200	
66064-840	フィルター(溶解試験用シリカフィルター)	¥207,600	¥224,400	
66064-842	フィルター(溶解試験用シリカフィルター)	¥209,600	¥226,600	
66065-142	VIAL 4ML CLEAR AP2000 PK100	¥263,700	¥290,600	
66100-036	貯蔵缶(ハダール式・ステンレス鋼製)	¥1,717,700	¥3,094,000	
66103-241	ゴミ箱(取り外し式キャスター付き)	¥80,500	¥86,100	
66103-300	ゴミ箱(取り外し式キャスター付き)	¥106,000	¥109,000	
66103-423	ゴミ箱(取り外し式キャスター付き)	¥178,200	¥170,200	
66104-110	ちりとり(ダストボックス用)	¥13,100	¥13,900	
66112-027	時計皿(浅型)	¥26,300	¥28,600	
66112-049	時計皿(浅型)	¥24,400	¥26,500	
66112-060	時計皿(浅型)	¥22,500	¥24,500	
66112-082	時計皿(浅型)	¥21,300	¥23,100	
66112-107	時計皿(浅型)	¥23,700	¥25,700	
66112-129	時計皿(浅型)	¥24,900	¥27,100	
66112-140	時計皿(浅型)	¥26,200	¥28,400	
66112-162	時計皿(浅型)	¥27,100	¥29,400	
66112-184	時計皿(浅型)	¥30,100	¥32,700	
66112-209	時計皿(浅型)	¥41,900	¥45,500	
66112-220	時計皿(浅型)	¥46,200	¥50,200	
66112-242	時計皿(浅型)	¥62,200	¥67,600	
66121-968	導電率計用電極 YSI 3400シリーズ	¥203,800	¥373,000	
66121-980	Self-Stirring B.O.D. Probe	¥542,300	¥996,000	
66175-100	サンプル採取器	¥96,500	¥103,700	
66175-101	サンプル採取器	¥35,000	¥36,800	
66175-612	サンプルバッグ 207mL 500枚入	¥54,900	¥57,600	
66176-507	KIMAX Nessler Color Comparison Tub	¥360,600	¥393,300	
66190-080	Vaccu-Pette96 ピペティングデバイス	¥29,500	¥31,600	
66197-008	顕微鏡用スライドガラス(カウティング・チャンバー)	¥254,200	¥281,500	
66198-000	顕微鏡用スライドガラス(カウティング・チャンバー)	¥28,000	¥31,000	
66238-332	Rocker II Platform Mixer Boekel	¥52,000	¥95,000	
66238-336	Rocker II Platform Mixer Boekel	¥512,200	¥932,000	
66238-422	Gable Lid(ウォーターバ 16L用)	¥169,900	¥319,000	
66238-424	Gable Lid(ウォーターバ 14L用)	¥177,200	¥332,000	
66238-426	Gable Lid(ウォーターバ 28L用)	¥246,900	¥462,000	
66238-440	Flat Lid(ウォーターバ 14L用/穴開き)	¥363,500	¥681,000	
66238-596	冷却循環恒温槽 20℃-40℃	¥168,100	¥315,000	
665A46	OptiPure Viral Auto Plate 96回用	¥58,300	¥59,900	
665S46	OptiPure Viral Auto Tube 96回用	¥58,300	¥59,900	
67100-734	有機カセット	¥1,398,200	¥2,519,000	
67100-736	有機カセット	¥873,700	¥1,574,000	
68310-252	CONSTIX Seald Foam Swab 500本入	¥49,500	¥53,000	
68310-254	綿棒 PP 500本入	¥85,100	¥91,100	
68310-355	PROSAT滅菌ワイパー(ホリパ 30x48PK)	¥677,700	¥725,500	
68400-018	NovaClean Lab and Glass Cleaner 1	¥35,200	¥41,800	
68400-020	NovaClean Lab and Glass Cleaner	¥185,000	¥200,000	
68400-422	CSLB MicroSeal Sleeve Seal Tape 1	¥26,800	¥28,900	
68400-682	NovaHol Cleanroom Detergent 3.78L	¥289,900	¥313,400	
68400-684	NovaHol Cleanroom Detergent 0.95L	¥337,400	¥364,800	
68400-776	CuffStrips ST 200	¥642,600	¥1,169,000	
685A46	Virapid Virus Auto Plate 96回用	¥46,600	¥47,900	
685S46	Virapid Virus Auto Tube 96回用	¥46,600	¥47,900	
68888-000	Silicone Adapter for macro 1入	¥7,000	¥7,600	
69100-210	VWR 滅菌ゴキブリ(Indirect) 50個入	¥431,000	¥441,000	
6960BP01	アセリ-ダストカバー	¥9,000	¥10,000	2024/3/14 追加
6960BP02	アセリ-ダストカバー	¥9,000	¥10,000	2024/3/14 追加
6960CP01	アセリ-ダストカバー LE台はかり表示部用	¥8,500	¥9,500	2024/3/14 追加
6960CP02	アセリ-ダストカバー LE丸皿ひょう量用	¥10,000	¥11,100	2024/3/14 追加
6960CP03	アセリ-ダストカバー LE角皿ひょう量用	¥9,500	¥10,600	2024/3/14 追加
6960CP04	アセリ-ダストカバー 12kg未満表示部用	¥10,000	¥11,100	2024/3/14 追加
6960CP05	アセリ-ダストカバー LE分析用	¥8,500	¥9,500	2024/3/14 追加
6960FB02	アセリ-ダストカバー	¥10,000	¥11,100	2024/3/14 追加
6960LC01	アセリ-プラスチックダストカバー	¥12,000	¥13,400	2024/3/14 追加
6960LP03	アセリ-ダストカバー	¥10,000	¥11,100	2024/3/14 追加
6960MS03	アセリ-ダストカバー	¥8,500	¥9,500	2024/3/14 追加
6960SE03	アセリ-ダストカバー	¥6,500	¥7,200	2024/3/14 追加
6960TE01	アセリ-ダストカバー	¥8,500	¥9,500	2024/3/14 追加
6960TE02	アセリ-ダストカバー	¥11,000	¥12,300	2024/3/14 追加
6960TE03	アセリ-ダストカバー	¥11,000	¥12,300	2024/3/14 追加
69898530A	接続ケーブル(3m) 1	¥89,000	¥97,900	
6EMA46	Microbiome DNA Auto Plate 96回用	¥77,700	¥79,900	
6EMS46	Microbiome DNA Auto Tube 96回用	¥77,700	¥79,900	
6GMA46	Food and Feed DNA AutoPlate 96回用	¥62,200	¥63,900	
6GMS46	Food and Feed DNA AutoTube 96回用	¥62,200	¥63,900	
6K2A46	Tissue RNA Auto Plate 96回用	¥62,200	¥63,900	
6K2S46	Tissue RNA Auto Tube 96回用	¥62,200	¥63,900	
6K3A46	Plant RNA Auto Plate 96回用	¥62,200	¥63,900	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
6K3S46	Plant RNA Auto Tube 96回用	¥62,200	¥63,900	
6PEA46	プラスミド抽出キット/プレート/一時販売停止	¥58,300	¥59,900	
6PES46	プラスミド抽出キット/チューブ/一時販売停止	¥58,300	¥59,900	
6T2A46	Tissue Total DNA Auto Plate 96回用	¥50,500	¥51,900	
6T2S46	Tissue Total DNA Auto Tube 96回用	¥50,500	¥51,900	
70230-010	フラスコ(ワイド用) 1	¥5,400	¥6,100	
70230-012	フラスコ(ワイド用) 1	¥4,500	¥4,700	
70230-014	フラスコ(ワイド用) 1	¥3,900	¥4,400	
70230-018	フラスコ(ボトル用)	¥11,000	¥12,300	
70380-018	リポリター 48ウェルプレート用	¥308,500	¥333,500	
70680-046	Fillpour Dialysis Tube Funnels 10	¥33,100	¥36,100	
70680-093	SCIENCEWARE Flaskup Flask Holder 3	¥17,300	¥18,400	
70680-110	Fritware Porous Polyethylen Sheet	¥87,700	¥95,500	
70680-112	Fritware Porous Polyethylene Sheet	¥133,700	¥163,200	
70680-168	Round Lead Rings 1個	¥17,200	¥18,700	
70680-172	栓(テフロン製)	¥101,500	¥109,700	
70680-180	栓(テフロン製)	¥79,000	¥85,400	
70680-186	栓(テフロン製)	¥69,200	¥74,700	
70680-192	栓(テフロン製)	¥66,200	¥70,600	
70680-270	Covamat Clear Bench Liner 1枚	¥80,800	¥86,500	
70680-298	薬さじ	¥34,400	¥37,300	
70680-314	マイクロチューブカッター 1個	¥20,700	¥22,400	
70680-316	Replacement Blade 1個	¥8,000	¥8,400	
70680-322	Screw Style Tubing Clamp 3個入	¥14,200	¥15,100	
70680-342	バリエーションチューブ 24x36インチ	¥119,600	¥128,000	
70680-344	バリエーションチューブ 24x30インチ	¥113,500	¥121,500	
70680-346	バリエーションチューブ 12x24インチ	¥44,300	¥47,400	
70680-376	VIAL W/SNAP CAP 237ml 12本入	¥11,600	¥12,700	
70680-392	栓(テフロン製)	¥65,900	¥70,200	
70680-394	栓(テフロン製)	¥60,400	¥64,300	
70680-396	栓(テフロン製)	¥52,700	¥56,100	
70680-398	栓(テフロン製)	¥46,700	¥50,400	
70680-600	SCIENCEWARE Flaskup Flask Holder 3	¥16,100	¥17,100	
70690-018	TRAY MICROSCOPE SLIDE PK6	¥64,700	¥69,300	
70690-104	バリエーションチューブ (HDPE製) 1000入	¥171,200	¥183,300	
70690-120	Swing Jaw Tubing Clamp 3入	¥41,400	¥43,900	
70690-124	Swing Jaw Tubing Clamp 3入	¥44,200	¥46,800	
70690-140	Conway Diffusion Cell 3	¥33,900	¥36,300	
70690-152	セフェリウムカラム W/CORD 1	¥25,700	¥27,500	
70690-178	薬さじ	¥15,600	¥17,000	
70690-182	水銀回収容器用交換用スポンジパッド 12	¥23,300	¥24,900	
70690-202	Industrial Drum Funnel PE 25L	¥56,100	¥60,100	
70690-276	Fritware Immersion Filters 125mm	¥61,700	¥67,400	
71000-280	Aspirator Pump	¥17,200	¥23,100	
71000-312	Aseptic Filling Bell 1入	¥82,600	¥97,200	
71000-318	Aseptic Filling Bell 1入	¥44,400	¥52,200	
71001-642	PolarSafe Cardboard Cryo Box 81wel	¥2,500	¥2,700	
72310-388	マイクロシリンジ 1700 25ul	¥53,000	¥57,300	
72310-390	マイクロシリンジ 1700 10ul	¥67,900	¥73,400	
72310-392	マイクロシリンジ 1000/S5 10ml	¥88,900	¥96,100	
72350-004	血球計算盤(ケース付)	¥500,600	¥554,300	
72460-044	凍結乾燥容器 ガラス製 150mL 1入	¥107,200	¥78,600	
72460-046	凍結乾燥容器 ガラス製 300mL	¥119,500	¥100,000	
72460-048	Flask Freeze Dry 900mL bottom 1個	¥165,400	¥137,800	
72460-050	凍結乾燥容器 ガラス製 1.2L	¥192,500	¥160,600	
72690-016	加液チューブ ガラス 20x150mm 30ml 250	¥292,500	¥319,000	
72690-018	加液チューブ ガラス 20x125mm 25ml 250	¥131,800	¥143,700	
72690-020	加液チューブ ガラス 16x150mm 20ml 1000	¥372,500	¥406,200	
72690-022	加液チューブ ガラス 16x125mm 16ml 250	¥155,900	¥170,000	
72690-024	加液チューブ ガラス 13x100mm 8ml 500	¥293,600	¥320,200	
72760-048	Center Well 10mm Diameter 144本入	¥99,200	¥108,200	
72760-050	STOPPER TOP CS144	¥136,800	¥149,200	
73000-202	小動物用ケージ (エコアミカケージ)	¥394,400	¥414,000	
73000-214	小動物用ケージ (エコケージ/フェイス) 1	¥11,000	¥11,600	
73320-084	Foam Mailer for Five 3/5ml Tubes	¥96,900	¥100,700	
73320-122	メイキング(郵送用梱包器材セット)	¥440,600	¥472,100	
73520-486	Pinchbar Design 96well Plate 180枚	¥340,500	¥361,200	
73960-370	Cryovial 5ml SS O-Ring	¥645,800	¥1,258,000	
73990-016	トランスファースティック(スプレッド)	¥291,800	¥323,200	
73990-026	トランスファースティック(スプレッド) 500入	¥44,500	¥49,200	
73990-078	トランスファースティック S Extライナー 400本	¥65,800	¥72,900	
73990-088	トランスファースティック 500本	¥64,400	¥71,400	
73990-092	トランスファースティック(スプレッド) 400	¥50,100	¥55,500	
73990-098	トランスファースティック(スプレッド) 400	¥47,100	¥52,200	
74280-016	Coplin Staining Jar with Cover	¥9,500	¥10,500	
74900-012	洗浄瓶(広口) 500 mL 6本入	¥28,400	¥30,400	
74950-090	マルチ用 PP製	¥27,200	¥29,100	
74950-172	Mini Sampler Spoons 0.1mL 1個	¥14,700	¥15,600	
74950-181	Economy Sterile Sample Spoons	¥64,500	¥71,100	
74950-286	Magnetic Stir Bar 1入	¥4,200	¥4,400	
74950-288	VWR Spinbar Magnetic Stir Bars	¥4,300	¥4,500	
74950-290	VWR Spinbar Magnetic Stir Bars	¥4,400	¥4,700	
74950-350	Secador Mini Desiccator Cabinets	¥138,700	¥149,500	
74950-366	PLATE POROUS MEDIUM REPLACEMENT	¥194,000	¥211,700	
75799-934	Cell Scrapers 滅菌済 100入	¥44,900	¥48,300	
75813-140	VWR Drosophila Stock Bottle 6oz	¥100,700	¥108,300	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
75832-404	タリオンチューブ 3355-L	¥447,900	¥484,200	
75838-670	Pre-Gasketed Sanitary Fittings 10	¥91,500	¥101,700	
75838-672	Pre-Gasketed Sanitary Fittings 10	¥80,000	¥88,900	
75838-678	SteriEnz-II/AT Std Sanitary Fittin	¥90,000	¥97,100	
75845-208	VWR FlipCap CentrifugeTubes ST15ml	¥178,300	¥196,500	
75845-210	VWR FlipCap CentrifugeTubes ST50ml	¥239,600	¥251,500	
75845-212	VWR FlipCap CentrifugeTubes ST15ml	¥200,800	¥210,800	
75856-764	KIMAX Volumetric Pipets 1mL 6本	¥92,000	¥100,300	
75856-804	KIMAX Volumetric Pipets 5mL 6本	¥124,700	¥136,000	
75856-852	KIMAX Volumetric Pipets 10mL 6本	¥151,600	¥181,800	
75856-856	KIMAX Volumetric Pipets 15mL 3本	¥87,200	¥95,100	
75856-862	KIMAX Volumetric Pipets 30mL 3本	¥109,600	¥119,500	
75856-866	KIMAX Volumetric Pipets 50mL 3本	¥124,300	¥149,100	
76008-422	Basic Protection Lab Coat 30	¥12,700	¥14,600	
76008-426	Basic Protection Lab Coat 30	¥12,800	¥13,500	
76036-006	Cell Lifter with Narrow Blade 100	¥52,500	¥55,100	
76051-210	Nitrile Exam Gloves 50	¥10,500	¥11,300	
76203-706	Protected Ambient Shipper 20	¥94,300	¥107,900	
76203-710	Protected Ambient Shipper 12	¥98,600	¥110,700	
76203-714	Ambient Gel Wraps 32	¥31,500	¥36,100	
76203-884	Protected Ambient Shipper 12	¥84,500	¥94,900	
76210-524	96-Well Deep Well Plates 50	¥124,100	¥133,500	
76210-676	Dippers, HDPE 1	¥52,900	¥57,000	
76293-642	鋳鉄 サポートリング 4φ	¥5,900	¥6,500	
76320-276	Disposable Serological Pipets 400	¥136,800	¥143,600	
76320-278	Disposable Serological Pipets 400	¥136,900	¥143,700	
76320-280	Disposable Serological Pipets 150	¥127,500	¥133,800	
76320-282	Disposable Serological Pipets 100	¥159,600	¥167,500	
76355-906	カバーガラス丸型 12MM 100	¥27,400	¥30,100	
76357-172	STOPPER STOPPER RUBBER 13mm 250	¥196,100	¥213,000	
76457-428	VWR Single-EdgedRazorBlades 100	¥6,100	¥11,000	
76457-540	Large Animal Dissecting Kit 1	¥24,400	¥28,100	
77100-278	カマトク プライ (HPLC用フイタインク)	¥365,100	¥379,500	
77100-286	カマトク プライ (HPLC用フイタインク)	¥330,000	¥344,000	
77100-290	カマトク プライ (HPLC用フイタインク)	¥327,400	¥340,800	
77100-294	カマトク プライ (HPLC用フイタインク)	¥381,800	¥395,800	
77100-628	カマトク プライ (HPLC用フイタインク)	¥315,400	¥344,900	
77100-630	カマトク プライ (HPLC用フイタインク)	¥315,400	¥344,900	
77100-644	カマトク プライ (HPLC用フイタインク)	¥332,200	¥360,700	
77100-646	カマトク プライ (HPLC用フイタインク)	¥332,200	¥360,700	
77100-700	カマトク プライ (HPLC用フイタインク)	¥305,600	¥332,000	
77100-722	カマトク プライ (HPLC用フイタインク)	¥318,200	¥344,900	
77100-724	カマトク プライ (HPLC用フイタインク)	¥318,200	¥344,900	
77100-788	カマトク プライ (HPLC用フイタインク)	¥314,800	¥343,100	
77100-790	カマトク プライ (HPLC用フイタインク)	¥314,800	¥343,100	
77100-814	カマトク プライ (HPLC用フイタインク)	¥338,500	¥367,200	
77100-816	カマトク プライ (HPLC用フイタインク)	¥338,500	¥367,200	
77100-874	カマトク プライ (HPLC用フイタインク)	¥315,300	¥343,700	
77100-876	カマトク プライ (HPLC用フイタインク)	¥315,300	¥343,700	
77100-946	カマトク プライ (HPLC用フイタインク)	¥227,400	¥237,900	
77100-948	カマトク プライ (HPLC用フイタインク)	¥227,400	¥237,900	
77100-970	カマトク プライ (HPLC用フイタインク)	¥188,900	¥196,700	
77100-972	カマトク プライ (HPLC用フイタインク)	¥191,100	¥199,600	
77101-026	カマトク プライ (HPLC用フイタインク)	¥435,200	¥473,800	
77101-028	カマトク プライ (HPLC用フイタインク)	¥435,200	¥473,800	
77101-030	カマトク プライ (HPLC用フイタインク)	¥436,000	¥474,600	
77101-114	カマトク プライ (HPLC用フイタインク)	¥331,400	¥359,400	
77101-116	カマトク プライ (HPLC用フイタインク)	¥331,400	¥359,400	
77101-118	カマトク プライ (HPLC用フイタインク)	¥331,900	¥360,100	
77101-120	カマトク プライ (HPLC用フイタインク)	¥331,400	¥359,400	
77101-122	カマトク プライ (HPLC用フイタインク)	¥331,400	¥359,400	
77101-294	カマトク プライ (HPLC用フイタインク)	¥332,100	¥328,000	
77101-310	カマトク プライ (HPLC用フイタインク)	¥282,800	¥305,000	
77101-312	カマトク プライ (HPLC用フイタインク)	¥282,800	¥310,400	
77101-338	カマトク プライ (HPLC用フイタインク)	¥362,300	¥387,800	
77101-340	カマトク プライ (HPLC用フイタインク)	¥366,700	¥393,300	
77101-366	カマトク プライ (HPLC用フイタインク)	¥365,300	¥388,700	
77101-368	カマトク プライ (HPLC用フイタインク)	¥365,300	¥388,700	
77776-338	シリコンチューブ	¥353,100	¥370,600	
77776-750	ビーカー用キャップ (ホ-ル紙製)	¥36,900	¥40,100	
77776-752	ビーカー用キャップ (ホ-ル紙製)	¥52,100	¥56,600	
77776-754	ビーカー用キャップ (ホ-ル紙製)	¥82,300	¥88,600	
77776-756	ビーカー用キャップ (ホ-ル紙製)	¥98,700	¥106,200	
77776-758	ビーカー用キャップ (ホ-ル紙製)	¥70,000	¥75,300	
77776-774	加-ニング シリカ - small 50入	¥64,800	¥69,400	
77776-776	加-ニング シリカ - Medium 50入	¥52,200	¥55,900	
77776-794	電気泳動用プレート	¥74,100	¥80,400	
77788-016	BCR クリ-ル-ム用ホ-ト-ハ-ル- 250枚x5	¥269,000	¥288,000	
80060-606	TLCスプレ-ティング用キャビ-リ-チューブ	¥187,800	¥219,900	
80060-608	TLCスプレ-ティング用キャビ-リ-チューブ	¥170,300	¥203,500	
80060-610	TLC SPOTTING CAPILLARY TUBES	¥84,000	¥98,200	
80060-664	Replacement stopcock	¥41,300	¥48,400	
80060-914	ワ-スプレ-ティング用 (For 13-425 Vials)	¥17,700	¥21,200	
80060-916	ワ-スプレ-ティング用 (For 15-425 Vials)	¥18,900	¥20,500	
80060-918	ワ-スプレ-ティング用 (For 18-425 Vials)	¥20,500	¥22,300	
80060-920	ワ-スプレ-ティング用 (For 20-400 Vials)	¥21,500	¥25,600	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
80060-922	フラスコ用フタ(For 22-400 Vials)	¥21,500	¥24,500	
80060-924	フラスコ用フタ(For 24-400 Vials)	¥27,900	¥33,400	
80060-926	コネクタ	¥144,600	¥157,100	
80061-030	SLEEVE PTFE ECONOMY 3入	¥47,300	¥51,400	
80061-180	オリフク PTFE Covered 100mm	¥33,200	¥36,100	
80061-850	Solid Cap with PTFE Seal #25 1	¥2,700	¥2,900	
80061-984	Hose Connector GL14 Bent	¥2,500	¥2,700	
80061-986	Hose Conn GL14 Straight	¥3,200	¥3,300	
80062-304	Stirrer Shaft 480mm OAL	¥76,700	¥82,900	
80062-412	Stirrer Paddle PTFE 50mm D	¥177,700	¥192,200	
80063-072	ホップアダプター(バキュームトラップ)	¥720,700	¥1,317,000	
80063-100	O-Ring 15.2/30 12個	¥12,900	¥14,000	
80064-216	MAINFOLD DOUBLE O-RING 3-PORTS	¥859,100	¥1,563,000	
80064-242	ADAPTER AF STYLE B AF 10MM	¥22,200	¥24,100	
80064-270	ADAPTER O-8MM 1入	¥159,200	¥173,000	
80064-292	BURBLER CAJON CONN. COMPLETE 1	¥124,700	¥135,500	
80064-380	Vacuum Trap O-Ring Body	¥79,000	¥85,800	
80064-394	STOPPER O-RING JOINT #15	¥19,400	¥21,000	
80064-438	VACUUM TRAP ASSEMBLY ONLY	¥474,800	¥515,800	
80064-440	VACUUM TRAP STOPPER ONLY	¥362,600	¥393,900	
80064-530	ガラスジョイント	¥99,000	¥107,600	
80064-532	ガラスジョイント	¥99,200	¥107,800	
80064-534	ガラスジョイント	¥114,200	¥124,000	
80064-536	ガラスジョイント	¥113,100	¥122,900	
80064-558	ガラスジョイント バッファ	¥50,700	¥60,500	
80064-586	ガラスケルダール形	¥69,500	¥75,500	
80064-588	ガラスケルダール形	¥76,700	¥83,300	
80064-590	ガラスケルダール形	¥79,400	¥86,300	
80064-592	ガラスケルダール形	¥83,000	¥90,200	
80064-594	ガラスケルダール形	¥90,200	¥98,000	
80064-596	ガラスケルダール形	¥88,800	¥96,500	
80064-598	ガラスケルダール形	¥94,400	¥102,600	
80064-600	ガラスケルダール形	¥101,700	¥110,500	
80064-626	Flask Storage 500ml 24/40	¥145,700	¥158,300	
80064-646	Flask Storage 100ml 24/40	¥152,600	¥165,800	
80064-656	Flask Reaction Glass 500ml	¥122,500	¥133,100	
80064-682	ガラス(Reaction/AirFree)250ml	¥132,100	¥143,500	
80064-684	FLASK REACTION 500ML 24/40 1	¥132,100	¥143,500	
80064-686	ガラス(Reaction/AirFree)100ml	¥119,700	¥130,000	
80064-738	TUBE REACTION 100ML 14/20 1	¥106,900	¥116,100	
80064-740	TUBE REACTION 100ML 24/40 1	¥108,400	¥117,800	
80064-854	FILTER TUBE 20MM C 14/20	¥118,400	¥128,700	
80064-856	FILTER TUBE 20MM M 14/20	¥118,500	¥128,700	
80064-862	FILTER TUBE 30MM M 14/20	¥126,700	¥137,700	
80064-866	FILTER TUBE 40MM C 24/40	¥130,400	¥141,700	
80064-868	FILTER TUBE 40MM M 24/40	¥130,400	¥141,600	
80065-488	ガラスジョイントアダプター(Reducing)	¥135,500	¥147,200	
80065-490	ガラスジョイントアダプター(Reducing)	¥137,000	¥148,800	
80065-500	ガラスジョイントアダプター(Reducing)	¥39,700	¥43,100	
80065-522	ガラスジョイントアダプター(Reducing)	¥47,100	¥51,100	
80065-524	ガラスジョイントアダプター(Reducing)	¥34,000	¥36,900	
80065-544	ガラスジョイントアダプター(Enlarging)	¥38,100	¥41,400	
80065-574	ガラスジョイントアダプター(Enlarging)	¥38,400	¥41,700	
80065-710	ガラスジョイント	¥32,700	¥37,400	
80065-712	ガラスジョイント	¥24,500	¥26,600	
80065-720	ガラスジョイント	¥27,800	¥30,300	
80065-752	ガラスジョイント	¥37,100	¥40,300	
80065-762	ガラスジョイント	¥42,500	¥46,100	
80065-772	ガラスジョイント	¥128,200	¥139,200	
80065-798	ガラスジョイント	¥88,400	¥100,800	
80065-802	ガラスジョイント	¥60,200	¥65,400	
80065-900	ガラスジョイント	¥93,800	¥101,900	
80065-902	ガラスジョイント	¥79,500	¥86,300	
80065-908	ガラスジョイント	¥85,700	¥93,100	
80065-972	ガラスジョイント #24/40	¥32,100	¥34,800	
80065-980	ガラスジョイント	¥31,500	¥34,300	
80065-994	Adapter Universal Inlet 29/42	¥34,000	¥36,900	
80066-052	ADAPTER TAKE-OFF LONG STEM 1	¥55,900	¥60,700	
80066-054	ガラスジョイント	¥68,200	¥74,100	
80066-058	ガラスジョイント	¥47,200	¥56,400	
80066-086	ガラスジョイント	¥56,800	¥61,700	
80066-090	ガラスジョイント	¥70,800	¥80,700	
80066-100	ガラスジョイント	¥82,400	¥89,500	
80066-462	ガラスジョイント(TLC展開用)	¥101,100	¥115,300	
80066-464	ガラスジョイント(TLC展開用)	¥142,000	¥169,700	
80066-466	ガラスジョイント(TLC展開用)	¥196,900	¥213,900	
80066-732	Reservoir 250ml 24/40 JT	¥68,200	¥74,100	
80066-818	Column Chromatography 24/40	¥141,800	¥154,000	
80066-820	Column Chromatography	¥163,900	¥178,100	
80066-826	Column Chromatography 1	¥140,900	¥153,000	
80066-828	Column Chromatography	¥163,900	¥178,000	
80066-834	ガラスジョイント 24/40 1個	¥175,500	¥190,600	
80066-836	Column Chromatography	¥199,200	¥216,400	
80066-842	ガラスジョイント 24/40 1個	¥182,300	¥198,100	
80066-844	Column Chromatography	¥211,200	¥229,400	
80066-850	Column Flask Chromatography	¥227,400	¥247,100	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
80066-852	Column Chromatography	¥248,300	¥269,700	
80066-854	Column Chromatography	¥290,400	¥315,500	
80066-878	Reservoir 250mL 24/40JT 1個	¥83,700	¥90,900	
80066-880	Reservoir 500ml	¥85,900	¥93,300	
80066-882	リザーブ 1000mL 1個	¥92,500	¥100,500	
80066-896	Adapter Flow Control 24/40	¥115,600	¥125,600	
80066-902	Column with Reservoir 100ml	¥122,800	¥133,400	
80066-908	Column with Reservoir 250ml	¥166,000	¥180,300	
80066-920	Column with Reservoir 1L	¥226,000	¥245,500	
80067-112	Flow Control Adapter	¥113,400	¥123,200	
80067-124	加圧ろ過フィルター 250mL 1個	¥149,600	¥162,500	
80067-220	Condenser Allihn 300mm	¥114,900	¥124,800	
80067-264	蛇管冷却器(コネクター)	¥215,900	¥234,500	
80067-278	蛇管冷却器(コネクター)	¥221,900	¥241,100	
80067-290	蛇管冷却器(コネクター)	¥192,700	¥209,300	
80067-292	蛇管冷却器(コネクター)	¥184,100	¥200,000	
80067-364	蛇管冷却器(コネクター)	¥195,700	¥212,600	
80067-368	蛇管冷却器(コネクター)	¥138,700	¥150,700	
80067-398	蛇管冷却器(コネクター)	¥78,200	¥85,000	
80067-572	多孔質沸騰石	¥175,700	¥190,900	
80067-574	多孔質沸騰石	¥219,500	¥238,400	
80067-582	多孔質沸騰石	¥223,200	¥242,500	
80067-584	多孔質沸騰石	¥238,800	¥259,400	
80067-614	多孔質沸騰石	¥375,400	¥407,800	
80067-874	Distilling Receiver 12ml 14/	¥260,800	¥283,300	
80067-904	Distilling Receiver w/Hooks 1	¥54,100	¥58,800	
80067-960	Drying Apparatus Complete Receiver	¥646,200	¥701,900	
80067-962	Drying Chamber Only	¥219,600	¥238,600	
80067-964	Desiccant Flask Only 34/45	¥270,100	¥293,500	
80068-042	ロトリバポレータフタ(ガラス瓶)	¥46,400	¥53,000	
80068-046	ロトリバポレータフタ(ガラス瓶)	¥58,100	¥63,100	
80068-048	ロトリバポレータフタトラップ、自動洗浄型	¥141,400	¥153,600	
80068-050	ロトリバポレータフタトラップ、自動洗浄型	¥141,000	¥153,100	
80068-052	ロトリバポレータフタトラップ、自動洗浄型	¥145,400	¥157,900	
80068-054	ロトリバポレータフタトラップ、自動洗浄型	¥121,500	¥132,000	
80068-060	ロトリバポレータフタトラップ	¥134,000	¥160,200	
80068-062	ロトリバポレータフタトラップ	¥127,000	¥138,000	
80068-064	ロトリバポレータフタトラップ	¥98,300	¥117,500	
80068-068	ロトリバポレータフタトラップ	¥116,100	¥126,100	
80068-080	ロトリバポレータフタトラップ	¥138,100	¥150,000	
80068-082	ロトリバポレータフタトラップ	¥127,000	¥144,800	
80068-084	ロトリバポレータフタトラップ	¥126,200	¥137,100	
80068-086	ロトリバポレータフタトラップ	¥121,800	¥145,500	
80068-090	ロトリバポレータフタトラップ	¥120,600	¥131,000	
80068-102	ロトリバポレータフタトラップ	¥128,800	¥139,900	
80068-104	ロトリバポレータフタトラップ	¥137,700	¥149,600	
80068-106	ロトリバポレータフタトラップ	¥143,300	¥155,700	
80068-108	ロトリバポレータフタトラップ	¥149,800	¥162,700	
80068-112	ロトリバポレータフタトラップ	¥122,800	¥133,400	
80068-130	ロトリバポレータフタトラップ 29/	¥89,200	¥96,900	
80068-388	ロトリバポレータフタ(ガラス瓶)	¥126,100	¥137,000	
80068-392	ロトリバポレータフタ(ガラス瓶)	¥130,500	¥141,800	
80068-394	ロトリバポレータフタ(ガラス瓶)	¥112,300	¥122,000	
80068-398	ロトリバポレータフタ(ガラス瓶)	¥144,400	¥156,800	
80068-400	ロトリバポレータフタ(ガラス瓶)	¥163,800	¥178,000	
80068-402	ロトリバポレータフタ(ガラス瓶)	¥141,900	¥154,200	
80068-404	ロトリバポレータフタ(ガラス瓶)	¥185,400	¥201,400	
80068-406	ロトリバポレータフタ(ガラス瓶)	¥223,600	¥242,900	
80068-408	ロトリバポレータフタ(ガラス瓶)	¥181,000	¥196,600	
80068-416	ロトリバポレータフタ(ガラス瓶)	¥214,300	¥232,800	
80068-418	ロトリバポレータフタ(ガラス瓶)	¥213,300	¥231,700	
80068-428	ロトリバポレータフタ(ガラス瓶)	¥260,000	¥282,500	
80068-430	ロトリバポレータフタ(ガラス瓶)	¥263,100	¥285,800	
80068-440	ロトリバポレータフタ(ガラス瓶)	¥87,200	¥94,700	
80068-442	ロトリバポレータフタ(ガラス瓶)	¥88,300	¥96,000	
80068-444	ロトリバポレータフタ(ガラス瓶)	¥86,300	¥93,800	
80068-452	ロトリバポレータフタ(ガラス瓶)	¥94,200	¥102,300	
80068-454	ロトリバポレータフタ(ガラス瓶)	¥115,800	¥125,800	
80068-456	ロトリバポレータフタ(ガラス瓶)	¥90,700	¥98,600	
80068-472	ロトリバポレータフタ(ガラス瓶)	¥92,700	¥100,700	
80068-474	ロトリバポレータフタ(ガラス瓶)	¥95,800	¥104,000	
80068-756	丸底フラスコ(共栓)	¥31,000	¥33,600	
80068-762	丸底フラスコ(共栓) 100ml	¥33,400	¥36,300	
80068-770	丸底フラスコ(共栓)	¥34,100	¥37,000	
80068-774	丸底フラスコ(共栓)	¥34,000	¥36,900	
80068-778	丸底フラスコ(共栓)	¥43,800	¥52,100	
80068-782	丸底フラスコ(共栓)	¥38,000	¥41,300	
80068-786	丸底フラスコ(共栓)	¥48,800	¥53,000	
80068-790	丸底フラスコ(共栓)	¥150,900	¥164,000	
80068-794	丸底フラスコ(共栓)	¥154,700	¥184,900	
80068-828	丸底フラスコ(共栓)	¥37,400	¥40,700	
80068-852	丸底フラスコ(共栓)	¥24,200	¥26,300	
80068-856	丸底フラスコ(共栓) 10ml	¥26,500	¥28,800	
80068-860	丸底フラスコ(共栓)	¥24,500	¥26,600	
80068-864	丸底フラスコ(共栓)	¥28,700	¥31,200	
80068-868	丸底フラスコ(共栓) 50ml	¥25,300	¥27,500	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
80068-870	丸底フラスコ(共栓)	¥28,300	¥30,800	
80068-872	丸底フラスコ(共栓)	¥28,600	¥33,900	
80068-874	丸底フラスコ(共栓)	¥23,400	¥25,400	
80068-876	丸底フラスコ(共栓)	¥28,500	¥31,000	
80068-878	丸底フラスコ(共栓)	¥35,300	¥38,300	
80068-880	丸底フラスコ(共栓)	¥36,600	¥39,800	
80068-964	共通摺合なす型フラスコ	¥36,100	¥39,200	
80068-968	共通摺合なす型フラスコ	¥40,300	¥48,200	
80068-972	共通摺合なす型フラスコ	¥48,800	¥58,300	
80068-976	共通摺合なす型フラスコ	¥71,000	¥77,100	
80068-978	共通摺合なす型フラスコ	¥64,000	¥69,500	
80068-980	共通摺合なす型フラスコ	¥78,400	¥85,200	
80068-998	共通摺合なす型フラスコ	¥136,500	¥148,300	
80069-006	共通摺合なす型フラスコ	¥30,200	¥32,800	
80069-010	共通摺合なす型フラスコ	¥30,200	¥32,800	
80069-014	共通摺合なす型フラスコ	¥29,900	¥32,500	
80069-018	共通摺合なす型フラスコ	¥28,800	¥31,300	
80069-022	共通摺合なす型フラスコ	¥34,100	¥37,000	
80069-134	共通摺合二口フラスコ	¥88,000	¥95,600	
80069-138	共通摺合二口フラスコ	¥105,000	¥114,100	
80069-146	共通摺合二口フラスコ	¥85,500	¥92,900	
80069-154	共通摺合二口フラスコ	¥105,000	¥114,100	
80069-174	共通摺合二口フラスコ	¥56,300	¥61,200	
80069-178	共通摺合二口フラスコ 25ml	¥71,600	¥77,800	
80069-184	共通摺合二口フラスコ	¥73,800	¥80,100	
80069-190	共通摺合二口フラスコ	¥38,200	¥41,500	
80069-196	共通摺合二口フラスコ	¥70,900	¥77,000	
80069-448	共通摺合三口フラスコ	¥119,000	¥129,200	
80069-452	共通摺合三口フラスコ	¥141,000	¥153,200	
80069-454	共通摺合三口フラスコ	¥153,200	¥166,400	
80069-456	共通摺合三口フラスコ	¥135,400	¥147,100	
80069-610	共通摺合三口フラスコ	¥67,000	¥72,800	
80069-612	共通摺合三口フラスコ	¥67,800	¥73,600	
80069-618	共通摺合三口フラスコ	¥65,500	¥71,100	
80069-624	共通摺合三口フラスコ	¥67,400	¥73,200	
80069-630	Round-Bottom Boiling Flasks 250ml	¥82,400	¥89,500	
80069-666	フラスコ 2L	¥302,200	¥328,300	
80069-806	フラスコ 25ML 1	¥99,400	¥108,000	
80070-186	梨型フラスコ(共栓) 500ml 1	¥57,000	¥68,100	
80070-198	梨型フラスコ(共栓) 5mL 14/20	¥32,100	¥34,800	
80070-200	梨型フラスコ(共栓) 5mL 19/22	¥25,200	¥27,300	
80070-636	付加口付(目盛付き)ガラス製	¥206,400	¥246,700	
80070-638	付加口付(目盛付き)ガラス製	¥214,800	¥233,400	
80070-640	付加口付(目盛付き)ガラス製	¥246,700	¥268,000	
80070-642	付加口付(目盛付き)ガラス製	¥270,600	¥293,900	
80070-646	付加口付(目盛付き)ガラス製	¥322,800	¥350,700	
80070-654	付加口付(目盛付き)ガラス製	¥179,600	¥195,000	
80070-660	付加口付(目盛付き)ガラス製	¥177,000	¥192,300	
80070-666	付加口付(目盛付き)ガラス製	¥179,600	¥195,000	
80070-672	付加口付(目盛付き)ガラス製	¥177,800	¥193,200	
80070-844	粉末口付(ガラス製)	¥40,600	¥44,100	
80070-846	粉末口付(ガラス製)	¥45,800	¥49,700	
80070-858	粉末口付(ガラス製)	¥30,400	¥33,000	
80070-886	粉末口付(ガラス製 変型)	¥82,600	¥89,700	
80070-890	粉末口付(ガラス製 変型)	¥96,100	¥104,400	
80070-900	粉末口付(ガラス製 変型)	¥113,200	¥123,000	
80071-368	Peptide synthesis Vessel 10mL	¥169,400	¥184,000	
80071-370	ガラス製合成容器 25ml	¥186,700	¥202,800	
80071-376	ガラス製合成容器 250ml	¥236,700	¥257,100	
80071-378	Peptide Synthesis Vessel 10ml	¥187,500	¥203,700	
80071-380	ガラス製合成容器 25ml	¥205,600	¥223,300	
80071-382	ガラス製合成容器 50ml	¥221,000	¥240,000	
80071-410	濾過フラスコ	¥60,400	¥65,600	
80071-412	濾過フラスコ	¥61,700	¥67,000	
80071-414	濾過フラスコ	¥66,900	¥72,600	
80071-418	濾過フラスコ	¥80,300	¥87,200	
80071-420	濾過フラスコ 120ml	¥84,900	¥92,300	
80071-422	濾過フラスコ 75ml	¥84,700	¥92,000	
80071-424	濾過フラスコ 150ml	¥91,400	¥109,300	
80071-426	濾過フラスコ 350ml	¥105,600	¥114,700	
80071-440	濾過フラスコ(フット付) 75ml	¥54,300	¥59,000	
80071-442	濾過フラスコ(フット付) 150ml	¥59,200	¥64,400	
80071-444	濾過フラスコ(フット付) 350ml	¥76,200	¥82,800	
80071-798	Stirrer Bearing Only 29/42	¥267,100	¥288,700	
80072-148	Sublimator Apparatus Complete 150ml	¥597,100	¥648,600	
80072-468	バルブ(真空マニフォルド)	¥646,900	¥1,182,000	
80072-470	バルブ(真空マニフォルド)	¥799,000	¥1,460,000	
80072-472	バルブ(真空マニフォルド)	¥903,300	¥1,651,000	
80072-474	マニフォルド(活栓付き2列型)ガラス製	¥617,000	¥1,122,000	
80072-476	マニフォルド(活栓付き2列型)ガラス製	¥711,000	¥1,300,000	
80072-478	マニフォルド(活栓付き2列型)ガラス製	¥780,500	¥1,427,000	
80072-548	ポンプアセサリー(バキュームトラップ)	¥106,600	¥195,000	
80072-550	ポンプアセサリー(バキュームトラップ)	¥120,700	¥221,000	
80072-552	ポンプアセサリー(バキュームトラップ)	¥137,100	¥251,000	
80072-554	ポンプアセサリー(バキュームトラップ)	¥179,300	¥328,000	
80072-578	ポンプアセサリー(バキュームトラップ)	¥195,500	¥357,000	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
80073-174	Joint O-Ring Seal 15mm 1個	¥10,200	¥11,000	
80073-440	Fritted Filter Discs Borosilicate	¥8,600	¥9,400	
80073-448	Fritted Filter Discs Borosilicate	¥8,500	¥9,200	
80073-456	Fritted Filter Discs Borosilicate	¥8,500	¥9,300	
80073-514	Filter Disc Fritted 60mm XC 1	¥33,100	¥35,900	
80074-032	濾過ワスリ	¥34,100	¥37,000	
80074-034	濾過ワスリ	¥48,200	¥52,400	
80074-058	Replacement 4mm Glass Plug	¥20,800	¥22,600	
80074-160	Replacement 4mm Teflon Plug Only	¥23,100	¥25,100	
80074-164	PLUG REPLACEMENT 2mm 11/25 1個	¥14,100	¥15,300	
80076-786	DuPont タハツカハール S 25枚	¥263,400	¥281,900	
80076-788	DuPont タハツカハール M 25枚	¥270,100	¥289,200	
80076-790	DuPont タハツカハール L 25枚	¥263,400	¥281,900	
80076-792	DuPont タハツカハール XL 25枚	¥274,000	¥293,400	
80076-794	DuPont タハツカハール XXL 25枚	¥291,400	¥311,900	
80076-874	DuPont タハツカハール M 100枚	¥291,600	¥312,100	
80076-876	DuPont タハツカハール L 100枚	¥297,400	¥318,400	
80078-782	Heating Tapes 115V 2.5cmx3.1m	¥215,500	¥238,800	
80079-876	VWR Serum Stoppers 20mm 1000個入	¥1,342,400	¥1,458,200	
80079-904	VWR Flip Cap Aluminium Seals 13mm	¥167,400	¥181,900	
80079-938	Sterile Dippas 滅菌済 100入	¥300,700	¥326,600	
80081-188	VWR ティンクサール スリット Standard 21cm	¥24,200	¥26,000	
80081-190	ティンクサール スリット 300	¥24,200	¥26,000	
80081-192	VWR ティンクサール スリット Macro 31cm	¥31,500	¥33,800	
80081-194	VWR ティンクサール スリット Micro ST 300	¥15,600	¥16,800	
80086-976	VWR Cabinet Slide Storage	¥82,800	¥86,900	
80087-062	BRAND 96well PP Plate 330uL	¥177,200	¥190,700	
80087-110	BRAND 0.2mL PCR Tubes	¥53,600	¥57,700	
80087-272	Nasal Cytology Curette 100入	¥253,000	¥270,800	
80089-644	Cell Culture Aspiration Station 115V	¥1,002,300	¥1,823,000	
80089-648	Cell Culture Aspiration Station 110V	¥743,200	¥1,352,000	
80089-736	Vacuum Aspirator Bottle 1個	¥157,800	¥185,800	
80090-036	NSK-24 Nitri-Solve Knit Gloves XL 12	¥66,100	¥70,800	
80092-830	Ansell Therm-A-Knit グローブ 12双	¥66,400	¥71,100	
80094-624	Cytology Brush 20.3cm 500個	¥178,000	¥195,300	
82001-070	排液回収タック用 Carbon filter KT/POLY	¥234,300	¥250,800	
82001-072	排液回収タック用 Carbon filter KT/SS	¥332,900	¥356,400	
82002-304	シリコン ルー 20ml 個包装 160入	¥113,900	¥118,300	
82003-414	VWR Microscope スライド ホック 1個	¥14,200	¥15,400	
82003-416	VWR マイナス 20 スライド ホック 1個	¥12,000	¥13,000	
82003-632	排液回収タック 7.6L 1	¥780,000	¥882,200	
82003-820	VWR Light-Duty Tissue Wipers 280	¥2,200	¥2,400	
82004-218	Ball Valve Lockouts 1.3-3.2cm	¥30,100	¥32,200	
82004-222	Ball Valve Lockouts 5.1-20.3cm	¥58,800	¥62,900	
82004-224	Ball Valve Lockouts 5.1-20.3cm	¥61,500	¥65,900	
82004-226	Ball Valve Lockouts 5.1-20.3cm	¥61,500	¥65,900	
82005-332	Economy NMR Sample Tube 5	¥20,900	¥23,000	
82005-426	Petroff-Hausser Cover Glass 1	¥36,400	¥40,300	
82006-730	CleanMark Thermal Label 2.2x1.3cm 1	¥74,800	¥80,100	
82006-731	CleanMark Thermal Label 2.5x7.6cm 1	¥160,500	¥171,800	
82007-142	Mechanical Cryogenic freezer box	¥4,100	¥4,900	
82007-144	Mechanical Cryogenic Freezer Box 1	¥4,400	¥4,800	
82007-146	Fiberboard Dividers 16穴 1入	¥1,600	¥1,700	
82007-148	Fiberboard Dividers 25穴 1入	¥1,600	¥1,700	
82007-158	VWR LN2 Cryogenic freezer box 1	¥4,600	¥5,000	
82007-160	VWR LN2 Cryogenic freezer box 1	¥4,500	¥4,800	
82007-162	VWR CryoPro 紙製 ホック 81本立	¥4,400	¥4,700	
82007-178	ホック リスタル 950ml 1000入	¥1,424,000	¥1,600,700	
82007-208	Removable Needles (22S) 6本入	¥47,300	¥51,100	
82007-220	VWR Magnetized Pipet Rack 1	¥53,700	¥59,000	
82007-226	VWR Pipetter Holder 1	¥14,100	¥16,900	
82010-460	Standard Metal Hub Luer Lock Needle 6	¥57,200	¥61,800	
82013-160	VWR セーフティチップ サイズ 38-430 96入	¥112,500	¥132,400	
82013-208	排液回収タック 3.8L 1	¥330,300	¥346,700	
82017-632	VWR Floating Foam Tube Racks 5	¥8,600	¥9,400	
82017-634	VWR Floating Foam Tube Racks	¥8,800	¥9,500	
82017-696	Kit HPLC manifold 1	¥352,100	¥369,500	
82018-048	VWR スリット-リア遠心管 50ml 50本x10袋	¥86,300	¥90,600	
82018-050	遠心チューブ 50ml (自立型・PP製) 500本	¥83,900	¥96,900	
82018-052	VWR スリット-リア遠心管 50ml 500本/袋	¥64,400	¥67,600	
82020-064	デジタル温度計 (ホック 取り付け型)	¥25,200	¥45,000	
82020-808	Sample Vials 3.8mL 144本	¥109,700	¥121,500	
82020-845	VWR UNDERPAD TISSUE 50入	¥11,100	¥12,100	
82020-852	Pens Series 39 Disp	¥127,200	¥144,800	
82020-854	ティンクサール スリット-リア	¥143,100	¥163,000	
82020-900	クリーニング用ワスリ (MF製/洗浄済) 150枚	¥41,400	¥43,700	
82020-913	VWR マイナス スリット 20枚用	¥19,200	¥20,800	
82021-118	Fiberboard Dividers 16穴 12枚入	¥8,900	¥9,600	
82021-120	Fiberboard Dividers 25穴 12枚入	¥8,900	¥9,600	
82021-133	VWR Laboratory Notebooks 1冊	¥18,600	¥20,400	
82021-135	Laboratory Notebooks 200pages 1入	¥21,200	¥23,200	
82021-172	VWR デジタル温度計	¥21,200	¥22,700	
82024-452	VWR 8ケル チューブ ラック 1個	¥12,600	¥13,700	
82024-454	VWR リバース チューブ ラック アウト 5個	¥43,400	¥47,200	
82024-471	マイナス チューブ ラック 80本立 1個- 5個	¥19,700	¥21,400	
82024-494	マイナス チューブ ラック 96本立 アウト ホック 5個	¥41,400	¥45,100	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
82024-510	VWR Floating Microtube Racks 4	¥15,300	¥16,700	
82024-548	Large Pipet Drawer Organizer 1	¥24,900	¥27,000	
82024-554	VWR Small Latex Bulbs 2ml 72	¥27,300	¥29,700	
82024-584	VWR Microscope SlideBoxes 50-Place	¥91,900	¥99,500	
82024-598	VWR 5-Place Slide Mailers 25	¥17,700	¥19,200	
82024-606	VWR Economy Microscope Slide Box 1	¥8,700	¥9,400	
82026-407	Laboratory Notebooks 200pages	¥208,000	¥228,300	
82027-378	解剖用ピンセット ラウトタイプ 152mm	¥4,900	¥5,400	
82027-382	VWRピンセット(ステン製) 1本	¥12,200	¥13,400	
82027-384	VWR Dissecting Forceps Curved	¥7,300	¥8,000	
82027-386	VWR Dissecting Forceps Fine Tip	¥7,500	¥8,200	
82027-390	VWR Cover Glass Forceps Bent 114mm	¥6,400	¥7,000	
82027-398	ピンセット(解剖用) 114mm	¥6,100	¥6,700	
82027-436	Specimen Forceps 1	¥11,500	¥13,900	
82027-438	Specimen Forceps 1	¥7,200	¥8,700	
82027-490	Round Tapered Spatulas 1	¥6,700	¥7,200	
82027-504	VWRスチールハンドル(w/PVC Handles) 1	¥7,600	¥8,200	
82027-510	VWRスチールハンドル(w/PVC Handles) 1	¥11,100	¥12,000	
82027-520	VWR Round/Tapered Micro Spatulas	¥4,900	¥5,300	
82027-522	VWR Hayman Style Spatulas	¥4,800	¥5,200	
82027-524	VWR ステンレススチール スチール 16.5cm	¥4,800	¥5,200	
82027-526	VWR ステンレススチール スチール 16.5cm 12本	¥32,800	¥39,200	
82027-534	テフロンコートマイクロナード	¥8,100	¥8,800	
82027-578	VWR Dissecting Scissors 1本	¥6,000	¥6,600	
82027-584	VWRハサミ(先端カーブ)ステン製 1本	¥5,800	¥6,400	
82027-588	VWR Dissecting Scissors 140mm	¥9,200	¥10,100	
82027-594	VWR Dissecting Scissors 165mm	¥11,200	¥12,300	
82027-604	VWRステンレススチール製ピンセット滅菌缶(角型)	¥14,900	¥16,100	
82027-606	VWRステンレススチール製ピンセット滅菌缶(角型)	¥15,200	¥16,500	
82027-612	VWR Zahn Buo-Type Viscometers 1本	¥63,800	¥68,900	
82027-614	VWR Zahn Buo-Type Viscometers 1本	¥98,000	¥105,900	
82027-616	VWR Zahn Buo-Type Viscometers 1本	¥97,500	¥105,400	
82027-618	VWR Zahn Buo-Type Viscometers 1本	¥63,800	¥68,900	
82027-620	VWR Zahn Buo-Type Viscometers 1本	¥97,100	¥104,900	
82027-688	Clear Bags with Hang Hole 4mil 5x8	¥53,700	¥57,500	
82028-028	General-Purpose Pads 100入	¥70,000	¥74,900	
82028-474	Carboy 10L PC w/handle CS2	¥121,600	¥136,400	
82029-064	Puritan DoubleCotton Applicator 25	¥1,800	¥2,000	
82029-734	PolarPack保冷剤 12個入	¥8,800	¥9,400	
82029-900	Operating Scissors 140mm	¥31,400	¥33,500	
82030-468	Disposable Hemacytometer 50入	¥76,200	¥83,600	
82030-472	血球計算盤(Fuchs Rosenthal)	¥76,900	¥84,300	
82030-474	血球計算盤(Fuchs Rosenthal)	¥30,900	¥33,900	
82030-482	血球計算盤(Burker Turk)	¥27,400	¥30,100	
82030-864	Heavy Absorbency Pad 200入	¥13,100	¥14,500	
82031-192	FUL-BAF Flask(キャップなし) 2.5L 1本	¥34,900	¥37,900	
82031-356	VWR Centrifugal Filters 0.2um 100	¥138,300	¥149,900	
82031-360	VWR Centrifugal Filters 0.45um 100	¥142,400	¥154,400	
82031-368	96穴平底プレート HighFlange 0.4mL	¥332,800	¥353,800	
82031-432	ShipSafeサードラール(Class9トライイ)	¥123,500	¥132,300	
82031-436	CuffStrips ST 400	¥583,400	¥1,061,000	
82050-716	96マイクロプレート(平底)MICROLON200 40	¥60,200	¥49,800	
82050-788	UV-Star 96マイクロプレート Chimney Style	¥240,500	¥167,200	
82051-212	ティンク-サークルピンセット(紙/ガラス)50ml	¥163,600	¥170,000	
82051-244	96ウェルプレートマイクロプレート ST	¥192,300	¥151,300	
82051-348	UV-Star 384マイクロプレート クリア 40入	¥540,900	¥386,800	
83007-374	VWR 96ウェルPCRプレート NT Half-Skirted	¥23,200	¥23,600	
83007-630	ボトルトップディスプレイ(容量可変)	¥592,500	¥637,500	
83007-646	ボトルトップディスプレイ 2-10ml	¥472,400	¥508,500	
83007-648	ボトルトップディスプレイ 10ml	¥475,900	¥512,100	
83007-654	ボトルトップディスプレイ(容量可変)	¥637,300	¥685,600	
83007-662	ボトルトップディスプレイ 50ml	¥946,600	¥1,018,200	
83007-708	TUBE DISPENSE (CERAMUS)	¥23,600	¥25,400	
83007-744	TUBE SUCTION	¥12,600	¥13,500	
83007-794	Adapter Air Hose	¥9,400	¥10,200	
83007-804	pipetus用 17マイクロ 5個	¥33,600	¥61,000	
83007-810	HOUSING NOSEPIECE 1個	¥28,000	¥30,100	
83008-710	Measuring Scoops HDPE 500mL 6個	¥48,500	¥52,600	
83009-302	W4 Critical Task Wipers 500枚	¥132,200	¥151,000	
83009-672	VWR 96well RealTime PCR Plate 10	¥23,200	¥23,600	
83009-678	Individual Tubes 1.2mL 1000本入	¥28,400	¥30,400	
87002-476	HISTOSETTE I Tissue Cassette 1500	¥140,400	¥147,300	
87002-574	生検飼料包埋用チップ	¥332,800	¥365,000	
87003-290	0.65ml 遠心チューブ 1000本	¥42,200	¥45,400	
87003-294	1.7ml マイクロプレート 500本	¥21,900	¥23,600	
87003-923	Prequalified Nomadic Shipper LRG	¥153,500	¥172,800	
87004-016	VWR PureStep Adhesive Mats	¥97,100	¥104,000	
87004-068	VWR PureStep粘着マット 30枚重ね	¥176,600	¥189,100	
87004-136	PharMed コレクション 25マイクロ	¥72,000	¥77,900	
87004-144	PharMed コレクション 25マイクロ	¥86,100	¥93,000	
87004-421	ThermoSafe Insulated Shippers	¥54,000	¥57,900	
89000-012	マイクロ遠心チューブ 0.5ml 1,000本	¥26,900	¥28,200	
89000-028	マイクロ遠心チューブ 1.5ml NT 500本	¥16,300	¥17,200	
89000-030	マイクロ遠心チューブ 1.5ml アンバー 500本	¥16,000	¥16,800	
89000-032	マイクロ遠心チューブ 1.5ml YL 500本	¥15,900	¥16,700	
89000-034	マイクロ遠心チューブ 1.5ml OR 500本	¥16,000	¥16,800	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
89000-036	マイクロ遠心チューブ 1.5ml RD 500本	¥16,000	¥18,400	
89000-038	マイクロ遠心チューブ 1.5ml PL 500本	¥15,900	¥16,700	
89000-040	マイクロ遠心チューブ 1.5ml BL 500本	¥17,300	¥18,100	
89000-042	マイクロ遠心チューブ 1.5ml GR 500本	¥15,900	¥16,700	
89000-044	マイクロ遠心チューブ 1.5ml アウト 500本	¥15,400	¥16,100	
89000-634	WypAll X80タワシ 初付 Pop-Up Box	¥89,900	¥96,300	
89000-692	ワイルド L40 41.7x24.9cm 900枚	¥103,000	¥115,100	
89000-750	KIMAX BerzeliusBeakers 50-500ml	¥65,700	¥71,700	
89000-860	加チャーチューブ CAP付 16x150mm 20ml 72	¥177,400	¥193,400	
89000-862	加チャーチューブ CAP付 25x200mm 70ml 36	¥179,700	¥195,900	
89000-926	KIMAXR GL45 Media/Storage Bottle	¥84,600	¥92,300	
89000-928	KIMAXR GL45 Media/Storage Bottle	¥98,200	¥107,100	
89000-936	KIMAXR GL45 Media/Storage Bottle	¥134,800	¥147,000	
89000-964	KIMAX RAY-SORB ガラスボトル 500mL	¥36,400	¥39,700	
89000-966	KIMAX RAY-SORB ガラスボトル 1L	¥42,100	¥45,900	
89001-066	KIMAX Griffin Beakers Glass 600ml	¥203,600	¥244,200	
89001-070	KIMAXR Griffin Beaker 200ml 4個	¥83,600	¥91,200	
89001-420	KIMAX Powder Funnels 100x18	¥196,300	¥214,100	
89001-454	KIMAX Culture Tube Durham 0.7ml 72	¥27,000	¥29,400	
89001-458	KIMAX Culture Tube Plain PK72	¥115,300	¥125,700	
89001-462	加チャーチューブ CAP付 13x100mm 8ml 72	¥145,300	¥158,500	
89001-464	加チャーチューブ CAP付 16x100mm 12ml 72	¥210,200	¥229,200	
89001-466	加チャーチューブ CAP付 16x125mm 16ml 72	¥208,200	¥227,000	
89001-468	加チャーチューブ CAP付 20x125mm 25ml 48	¥166,600	¥181,600	
89001-470	加チャーチューブ CAP付 20x150mm 30ml 48	¥190,100	¥207,200	
89001-472	加チャーチューブ CAP付 25x150mm 50ml 36	¥153,900	¥167,900	
89001-474	加チャーチューブ CAP付 25x200mm 70ml 36	¥178,500	¥194,700	
89001-476	加チャーチューブ CAP付 37x200mm 160ml 24	¥309,600	¥337,600	
89001-478	加チャーチューブ CAP付 13x100mm 8ml 72	¥204,900	¥223,400	
89001-480	加チャーチューブ CAP付 16x100mm 12ml 72	¥213,100	¥232,400	
89001-482	加チャーチューブ CAP付 16x125mm 16ml 72	¥232,800	¥253,900	
89001-484	加チャーチューブ CAP付 16x150mm 20ml 72	¥247,800	¥270,300	
89001-486	加チャーチューブ CAP付 20x125mm 25ml 48	¥211,300	¥230,400	
89001-488	加チャーチューブ CAP付 20x150mm 30ml 48	¥220,200	¥240,200	
89001-490	加チャーチューブ CAP付 25x150mm 50ml 36	¥199,600	¥217,700	
89001-520	KIMAX SolidGlassBeads 6mm 0.45kg 1	¥177,800	¥193,900	
89001-522	Sample/Serum ボトル 200ml 48	¥168,600	¥183,900	
89001-544	Milk Dilution ボトル 目盛付 48	¥238,300	¥259,900	
89001-546	Milk Dilution ボトル 目盛付 48	¥338,300	¥368,900	
89001-692	KIMAX ReverseScaleCylinder 250ml 1	¥53,200	¥58,000	
89002-066	KIMAX Volumetric Flasks 900mL 1	¥178,500	¥194,700	
89002-194	KIMAX Centrifuge Tubes 15mL 12本	¥168,800	¥184,100	
89002-198	KIMAX Centrifuge Tubes 15mL 12本入	¥268,400	¥292,700	
89002-200	KIMAX遠心管 50mL(目盛付/V底)12本入	¥383,500	¥418,200	
89002-214	KIMAX Centrifuge Tubes 10mL 12本	¥317,800	¥346,600	
89002-216	KIMAX Centrifuge Tubes 40mL 12本	¥414,000	¥451,500	
89002-222	遠心チューブ 丸底 スクリューキャップ 50ml 12	¥211,700	¥230,900	
89002-238	KIMAX Graduates 2.5oz 16oz/500ml	¥188,000	¥205,000	
89003-002	WypAll X80タワシ レット Pop-Up Box	¥54,800	¥63,700	
89003-004	WypAll X80タワシ ボール BRAG Box	¥22,900	¥24,000	
89003-006	WypAll X80タワシ ボール ジェンボロール	¥52,700	¥64,000	
89003-010	WypAll X80タワシ レット ジェンボロール	¥46,200	¥49,500	
89003-064	W5 Critical Task Wipers 100枚入	¥26,000	¥27,300	
89003-384	Reagent Reservoirs 50mL ST 50個	¥67,100	¥72,400	
89003-516	スリット 5ml 12本入	¥142,500	¥155,400	
89003-530	25ml スリット 5ml (ガラス製)	¥122,400	¥133,500	
89003-552	スクリューキャップ ミニ 13mm用 13-415 1000	¥91,200	¥99,400	
89003-554	スクリューキャップ ミニ 15mm用 15-415 1000	¥107,800	¥117,600	
89003-556	キャップ 18-415 20mm用 フェノール 1000	¥124,100	¥135,400	
89003-558	スクリューキャップ ライン付 PP 1000	¥70,100	¥76,400	
89003-842	WypAll X80タワシ レット クォーターフォルト	¥45,500	¥54,400	
89004-282	VWR スクリューキャップ 遠心チューブ 0.5mL ST	¥83,100	¥87,200	
89004-284	凍結保存チューブ 0.5ml/自立/未滅菌	¥56,300	¥59,100	
89004-286	凍結保存チューブ 0.5ml/自立/滅菌済	¥78,200	¥82,100	
89004-290	Screw-Cap tube ST 500入	¥82,400	¥86,400	
89004-298	2.0ml Screw Cap tube 500本	¥82,300	¥86,400	
89004-300	VWR スクリューキャップ 遠心チューブ 2mL 500	¥56,300	¥65,000	
89004-320	凍結保存チューブ 自立 1.5ml 500入	¥58,600	¥61,500	
89004-384	WypAll X80タワシ 初付 クォーターフォルト 50	¥15,200	¥18,300	
89004-418	VWR Sterile sample bags 150ml	¥70,900	¥78,100	
89004-554	Kimble Serum Bottle Clear 200ml	¥110,100	¥120,100	
89004-556	WypAll X80タワシ 初付 ジェンボロール	¥68,300	¥76,300	
89004-570	VWR ハンフトップ PCRクワイック 2個入	¥64,500	¥65,600	
89005-130	VWR Laboratory Notebooks	¥25,000	¥28,400	
89005-238	Write-On ラベリングテープ	¥13,500	¥15,200	
89005-240	Write-On ラベリングテープ	¥63,600	¥70,000	
89005-324	UV protectant アルバニウム 1000枚	¥50,600	¥53,900	
89005-330	VWR Reclosable UV7 紫外線遮光フィルム	¥73,300	¥78,100	
89005-356	ドキュメントポーチ (Document Pouch)	¥56,700	¥62,200	
89005-360	SAFEGUARD Specimen Bags	¥72,200	¥79,200	
89005-362	SAFEGUARD Specimen Bags	¥72,500	¥79,500	
89005-544	Prepacked Lens Station w/SILI	¥49,700	¥53,200	
89005-580	Mini Tube System 12-Tube Strip 960	¥81,200	¥84,800	
89005-596	Graduated Transport チューブ 1000本	¥235,000	¥254,000	
89005-758	Ultra-MicroBIO-CERT Cuvettes3.5mm	¥97,900	¥106,300	
89009-200	Poxygrid 50mL Tube Slant Rack	¥52,200	¥55,900	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
89009-656	Precut Pierceable Vinyl Films 100	¥68,500	¥73,700	
89009-666	VWR HDPEコップ 500mL 100個	¥117,900	¥124,300	
89009-668	VWR HDPEコップ 1L 100個	¥204,900	¥216,100	
89011-724	VWR Glove Box Holders	¥23,500	¥25,400	
89012-138	KIMAX Media Bottle Graduated 5L	¥98,100	¥115,600	
89012-366	VWR Glass Burets Stopcock付 50ml	¥122,100	¥132,600	
89012-608	Activated Charcoal Filter VetEquip	¥12,800	¥14,600	
89012-688	Induction Chambers 2L 1入	¥345,500	¥364,500	
89012-816	Hooded Coveralls with Elastic Ankles	¥217,800	¥233,100	
89013-156	VWRグローブ (DuPont製) Large	¥151,800	¥162,500	
89013-158	VWRグローブ (DuPont製) X-Large	¥180,500	¥193,300	
89020-076	保護メガネ Clear 1	¥4,200	¥4,600	
89022-356	Harris Uni-Core Punch 3mm	¥63,000	¥72,400	
89023-804	ハサミ(個別包装) 25本	¥105,500	¥117,100	
89024-859	オートクレーブライナー 1000枚	¥96,100	¥105,600	
89024-876	Chamber Cooling Each 0.2mL	¥375,300	¥410,000	
89024-892	Chambers Specialty EA 0.2 PCR	¥675,800	¥737,800	
89025-820	Volumetric Flasks(スリット付)10ml	¥42,600	¥46,600	
89026-142	Funnels, PTFE	¥166,000	¥169,300	
89026-378	フラスコ(ラット) 85ml 40x100mm 1	¥37,800	¥38,700	
89026-386	Tissue Grinders Glass Vessel 2mL	¥49,800	¥54,100	
89026-398	Tissue Grinders Plain Plunger 2mL	¥24,600	¥25,100	
89026-490	VWR PTFE Thin Sleeves 10	¥58,700	¥60,700	
89026-602	広口透明ガラスボトル 2L 6個	¥244,500	¥277,700	
89026-604	広口褐色ガラスボトル 1L 12個	¥180,700	¥205,200	
89026-606	広口褐色ガラスボトル 2.5L 4個	¥164,100	¥185,700	
89026-608	広口透明ガラスボトル 500mL 12個	¥178,800	¥203,000	
89027-106	Absorber Light-Duty Roll	¥198,500	¥212,500	
89027-108	Spillyter Soaker(吸収材) 250枚入	¥91,500	¥98,000	
89027-110	Absorber Heavy-Duty Roll	¥158,100	¥169,200	
89027-112	Absorber Heavy-Duty Pad	¥118,700	¥127,100	
89027-280	O-Ring Perfluoro Size 110	¥21,400	¥22,900	
89027-322	O-RING PERFLUORO Size210 1	¥28,900	¥31,300	
89027-334	O-Ring Perfluoro Size 220	¥46,500	¥50,200	
89027-702	Shaker Platform 33x27.9cm	¥104,700	¥273,000	
89027-858	キャップ付広口ビュ PP 30mL 800本	¥352,100	¥378,800	
89028-166	ADAPTER#15 O-RING JOINT TO	¥93,400	¥101,400	
89028-170	ADAPTER THERMOCOUPLE GAUGE	¥237,800	¥257,000	
89028-172	GAUGE ADAPTER COMPLETE #25	¥168,100	¥181,700	
89028-174	GAUGE ADAPTER COMPLETE #30	¥246,900	¥266,900	
89028-176	BUSHING #303 STAINLESS STEEL	¥63,700	¥68,900	
89028-424	ロータリーポンプ用ガスケット(アルミ) 1個	¥135,400	¥145,800	
89028-692	FUNNEL SEP SQUIBB 30ML 1個	¥75,300	¥88,100	
89028-694	FUNNEL SEPARATOR SQUIBB 60ML 1個	¥77,900	¥84,600	
89028-696	FUNNEL SEP SQUIBB 125ML 1個	¥88,600	¥96,300	
89028-698	FUNNEL SEP SQUIBB 250ML 1個	¥104,600	¥113,700	
89028-700	FUNNEL SEPARATOR SQUIBB 500ML 1個	¥122,100	¥145,900	
89028-802	フラスコ(ラット) 28ml	¥103,200	¥112,100	
89028-804	フラスコ(ラット) 65ml	¥44,900	¥47,100	
89028-806	フラスコ(ラット) 200ml	¥120,900	¥131,300	
89028-808	フラスコ(ラット) 420ml	¥133,500	¥145,000	
89028-810	Pressure Vessel only 28mL 1個	¥86,100	¥93,600	
89028-812	Pressure Vessel only 65mL 1個	¥92,400	¥100,400	
89028-814	Pressure Vessel only 200mL 1個	¥104,900	¥114,000	
89028-816	Pressure Vessel only 420mL 1個	¥116,900	¥127,000	
89029-233	Viton O-Rings サイズ 210 6個	¥15,800	¥18,100	
89029-235	Viton O-Rings サイズ 212	¥14,300	¥15,500	
89029-247	Viton O-Rings サイズ 220	¥21,200	¥24,200	
89029-406	VACUUM MANIFOLD 4-PORT 1個	¥718,000	¥1,312,000	
89029-408	VACUUM MANIFOLD 5-PORT 1個	¥877,000	¥1,603,000	
89029-590	Pressure Relief Valve Complete 1	¥85,700	¥93,100	
89029-592	Spring Loaded Teflon Plug Only 1	¥21,800	¥23,700	
89029-726	Corrosion Disc	¥8,900	¥9,800	
89029-728	Repture Disc Inconel	¥8,900	¥9,800	
89029-766	VESSEL ACID DIGESTION 45ML 1	¥826,000	¥923,100	
89029-919	1-200ml用 1.0L用 960本	¥43,200	¥45,800	
89030-552	Magnetic Rare Earth Stir Bar 10入	¥68,700	¥77,100	
89030-781	Komet Magnetic Stirring Bars 21mm	¥39,900	¥46,600	
89031-270	VWR Critical Swab Cotton Head 1000	¥17,800	¥18,900	
89031-272	Critical Swab Large Cotton Head 500	¥11,800	¥12,500	
89031-458	VWR Rota-Rack Tube Racks	¥28,000	¥33,600	
89031-950	Silver Shield 4H Gloves サイズ 10	¥530,500	¥587,500	
89032-092	VWRヒーター(Standard) 500W	¥1,142,000	¥2,077,000	
89032-104	ヒーター 120V Model5000	¥2,500,700	¥4,547,000	
89032-164	Replacement Power Supply 1個	¥28,600	¥52,000	
89032-170	Mop Head Refill w/フリーストール 1個	¥562,600	¥619,900	
89032-174	VWR Mop Head Refills	¥78,000	¥85,900	
89032-180	VWR Mop Head Refill 35.6cm 12個	¥228,800	¥251,900	
89032-182	ホリスティックフリーストール 35.6cm 12個	¥208,600	¥197,100	
89033-116	1/2" 製作用(綿ビロード) 36枚入	¥60,700	¥65,000	
89037-480	乳鉢(メー製) 25x20	¥98,000	¥106,400	
89037-484	乳鉢(メー製) 51x38	¥153,800	¥167,100	
89038-038	VWR Gooch Crucibles 13mL	¥37,800	¥41,100	
89038-040	VWR Gooch Crucibles 25mL	¥40,000	¥43,400	
89038-048	ろつぼ 1	¥29,300	¥31,800	
89038-050	ろつぼ 1入	¥29,300	¥31,800	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
89038-134	0-ト(フ)型 ろ過用 陶器製)	¥195,400	¥212,300	
89038-144	VWR乳鉢 50mL 1個入	¥10,700	¥11,700	
89038-146	VWR乳鉢 65mL 1個入	¥11,700	¥12,700	
89038-148	VWR 乳鉢 145mL 1個入	¥15,800	¥17,200	
89038-154	VWR乳鉢 750mL 1個入	¥43,800	¥47,600	
89038-160	VWR乳棒 114mm 1本入	¥12,200	¥13,200	
89038-162	VWR乳棒 130mm 1本入	¥12,500	¥13,600	
89038-164	VWR 乳棒(89038-148用) 1本入	¥15,100	¥16,400	
89038-170	VWR乳棒 222mm 1本入	¥32,300	¥35,100	
89038-196	ろつぽ 1	¥30,600	¥33,200	
89038-756	VWR Visual Alarm タイマー 1	¥10,800	¥11,900	
89039-656	VWR遠心チューブ 50mL ST 500本	¥57,900	¥60,800	
89039-658	50ml チューブ フラットキャップ 500本入	¥67,700	¥93,200	
89039-660	VWR遠心チューブ 50mL ST 500本	¥58,200	¥61,100	
89039-664	VWR遠心チューブ 15mL 500本	¥48,500	¥50,900	
89039-668	VWR遠心チューブ 15mL ST 500本	¥51,900	¥54,400	
89039-690	SlideFolder Microscope Slide Holde	¥101,700	¥106,700	
89040-062	Pipet Tips 1000ul 96/177x167mm	¥109,100	¥116,800	
89041-324	CAPS ALUMINUM SEAL 20mm 200	¥164,900	¥179,100	
89041-330	Crimper, Hand Operated 20mm 1	¥187,700	¥203,900	
89041-332	Pliers Decapper 20mm 1	¥56,100	¥67,000	
89041-338	ADAPTER MICROWAVE VIAL 24/40	¥159,000	¥172,700	
89041-776	粘着マット 60枚入	¥172,100	¥199,600	
89042-240	Stir Bar Restrainer 1	¥9,600	¥10,300	
89042-314	ディスプレイホルダー(ガラス) 1個	¥131,800	¥150,100	
89044-508	VWR Graduated Beakers 25mL 12個	¥75,200	¥82,900	
89044-510	VWR Graduated Beakers 50mL 12個	¥76,100	¥83,800	
89045-760	Straight-Sided Jars 12入	¥107,600	¥117,300	
89045-848	Straight-Sided Jars 24入	¥123,700	¥134,900	
89045-850	Straight-Sided Jars 24入	¥136,100	¥148,500	
89045-854	Straight-Sided Jars 12入	¥141,000	¥153,800	
89045-856	Straight-Sided Jars 12入	¥156,600	¥170,800	
89045-866	Straight-Sided Jars 24入	¥61,500	¥67,100	
89046-978	Straight-Sided Vials 4mL 250本	¥299,200	¥326,200	
89047-206	VWR 注射器用シリンジ 1mm 90ul	¥112,700	¥118,500	
89047-208	VWR 注射器用シリンジ 2mm 400ul	¥112,300	¥118,100	
89047-296	VWR Traceable Spatula balance	¥36,800	¥41,400	
89049-012	Magnetic Stirrer 1	¥48,200	¥52,100	
89049-034	ColdStorageAdhesiveSealingFoil 50	¥41,900	¥45,100	
89049-142	Vacuum Filtration Adapter	¥70,600	¥76,700	
89049-148	Vacuum Filtration Adapter	¥57,700	¥62,600	
89049-150	Vacuum Filtration Adapter	¥66,800	¥72,500	
89049-170	遠心チューブ 15ml PP 滅菌済 500本	¥72,100	¥75,700	
89049-172	VWRマルチ遠心チューブ 15ml ST 500本	¥84,200	¥88,400	
89049-176	VWRマルチ遠心チューブ 50ml ST	¥103,700	¥108,900	
89049-178	VWR 96-Well Thermal Cycling Plate	¥16,300	¥16,600	
89049-670	UNIFROST Slides 1440枚	¥176,700	¥193,800	
89049-982	PolarPack Foam Brick Gel Packs 158	¥34,400	¥36,800	
89049-986	PolarPack Foam Brick 72	¥58,400	¥62,600	
89049-990	PolarPack Gel Packs(保冷剤)24個	¥88,800	¥95,200	
89051-098	24/40 Adapter Inner	¥26,500	¥31,200	
89051-194	24/25 Adapter Twin	¥75,400	¥88,700	
89051-882	5555チューブホルダー 250ml	¥125,000	¥147,100	
89051-884	5555チューブホルダー 250ml	¥106,800	¥125,700	
89052-284	Jacketed Recycling Column 1セット	¥662,000	¥1,314,000	
89052-844	7mm PTFE Plug	¥22,700	¥26,700	
89055-514	Chimney only ST 24/40 #5 Joint	¥219,300	¥258,000	
89055-516	Adapter Tube only #5 Joint	¥60,500	¥71,200	
89055-518	Distillation Apparatus	¥328,900	¥388,100	
89057-716	Gas Washing Bottle 250mL 1	¥226,300	¥266,400	
89059-306	Midget Impinger(Complete) 1個	¥78,500	¥92,400	
89059-376	Impinger Air Sampling (Complete) 1	¥86,100	¥101,300	
89060-256	Pinch Clamp w/Screw-Locking Device	¥14,600	¥16,700	
89060-606	O-Rings Viton 12入	¥11,700	¥15,100	
89060-626	O-Rings Viton 12入	¥14,900	¥19,100	
89060-728	O-Rings Viton 12入	¥20,200	¥26,100	
89060-892	O-Rings Viton 3入	¥19,500	¥25,100	
89060-898	O-Rings Kalrez 1入	¥37,400	¥48,100	
89060-926	O-Rings Kalrez 1入	¥100,100	¥128,500	
89061-626	Inner Top Only 50mm	¥175,500	¥206,400	
89063-334	Pressure Tube 120mL 1	¥47,100	¥55,400	
89064-926	Septa, Sleeve Type, Rubber 72入	¥89,400	¥105,200	
89064-976	Septum for 24/40 joints	¥20,800	¥24,500	
89064-978	Septum for 24/40 joints	¥20,700	¥24,300	
89065-280	UltraSenseECI [®] カリブコントロール 100	¥21,800	¥22,400	
89065-290	Steraffirm Bowie-Dick Test Pack 20	¥300,300	¥341,600	
89066-380	Surgical Platform	¥513,300	¥628,500	
89067-546	Rodent Grain-Based Diet	¥66,500	¥69,800	
89067-562	Rodent Purified Diet	¥72,800	¥76,400	
89067-850	Mouse Igloo 1個	¥3,200	¥3,400	
89067-852	Mouse Igloo 1個	¥3,300	¥3,400	
89067-854	Fast-Trac 1個	¥4,800	¥5,000	
89067-856	Fast-Trac 1個	¥4,800	¥5,000	
89067-866	Rodent Retreats Rat EnrichmentDevi	¥4,100	¥4,300	
89067-868	Mouse Tunnels 1個	¥3,900	¥4,100	
89067-874	Rat Tunnels 1個	¥4,100	¥4,300	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
89067-876	Rat Tunnels 1個	¥4,100	¥4,300	
89068-044	Gloves Fluonit 468 サイズ 8	¥96,000	¥102,800	
89068-468	VWR Silicone Tubing 1	¥68,500	¥71,900	
89068-474	VWR Silicone Tubing 1	¥74,000	¥77,700	
89073-624	Scrapers&Spatulas Rubbermaid	¥5,500	¥5,900	
89073-844	液体吸収材(オイル専用/0-10リットル)	¥87,300	¥93,400	
89074-286	Vacuum Take-Off Long Stem Adapter	¥49,300	¥58,000	
89074-898	ガラス 25ml 14/20 1	¥22,200	¥26,100	
89075-592	Stainless Steel Beaker 1	¥86,000	¥105,800	
89079-154	Volumetricガラス 10ml 6	¥148,000	¥161,400	
89079-158	Volumetricガラス 25ml 6	¥154,700	¥168,700	
89079-160	Volumetricガラス 50ml 6	¥166,900	¥182,000	
89079-162	KIMAX Volumetric Flasks 栓付100ml	¥203,200	¥221,600	
89079-176	KIMAX Volumetric Flasks 10ml	¥157,000	¥171,200	
89079-182	KIMAX Volumetric Flasks 50ml	¥173,100	¥188,800	
89079-184	KIMAX Volumetric Flasks 100ml	¥209,800	¥228,700	
89079-192	KIMAX Volumetric Flasks 1L	¥94,700	¥103,200	
89079-220	Volumetricガラス 10ml 6	¥136,300	¥148,700	
89079-228	KIMAX Volumetric Flasks 100ml	¥203,800	¥222,200	
89079-404	Microwave Reaction Vials 2-5mL	¥382,500	¥421,500	
89079-535	Light Sensitive遠心ボトル 50ml ST 50	¥15,600	¥16,300	
89079-537	Light Sensitive遠心ボトル 50ml ST 25	¥10,100	¥10,600	
89079-974	VWR Single Channel Pipettors 1	¥216,000	¥191,900	
89080-300	Dissolution Cannula Filters 1000	¥491,500	¥526,200	
89080-532	24-Well Deep Well Plates 5	¥32,900	¥34,700	
89080-534	24-Well Deep Well Plate 10mL	¥51,400	¥54,200	
89082-458	マグネティックスターバール-アソートメント 11個入	¥46,400	¥49,400	
89083-242	Clamp, Dual Adjustment 1	¥31,400	¥37,500	
89083-356	加圧シリンダー - 125個	¥70,300	¥82,700	
89083-358	加圧シリンダー - 125個	¥84,600	¥99,500	
89083-360	加圧シリンダー - 125個	¥59,600	¥70,200	
89083-362	加圧シリンダー - 125個	¥61,700	¥72,600	
89083-364	加圧シリンダー - 125個	¥68,500	¥80,600	
89083-974	Tapered Five-Neck Flasks 125mL 1	¥444,700	¥523,200	
89083-984	Tapered Five-Neck Flasks 50mL 1	¥444,700	¥523,200	
89084-914	Clear Vial 4mL 200本	¥146,000	¥176,300	
89085-160	Volumetric Pipets Class A 6入	¥55,400	¥60,500	
89085-162	Volumetric Pipets Class A 6入	¥61,000	¥66,600	
89085-224	Volumetric Pipets Class A 6入	¥128,200	¥139,900	
89085-226	Volumetric Pipets Class A 6入	¥137,500	¥149,900	
89085-300	ディスク-サンプルプレート(紙/75)2ml	¥248,900	¥261,300	
89085-399	VWR Frosted Microscope Slide72x20P	¥130,600	¥143,200	
89085-518	PolarSafe 81-Place Freezer Box 5個	¥25,700	¥31,300	
89085-524	PolarSafe 81-Place Freezer Box 5個	¥25,500	¥33,800	
89087-398	カートリッジ	¥13,300	¥16,100	
89088-038	Vacuum Aspirator Bottle 1個	¥205,400	¥219,900	
89088-054	Carton Cover (56617-801用)	¥24,300	¥26,600	
89089-120	Transparent Poly Bags 100	¥26,000	¥28,600	
89089-922	HTS Transwell Support Systems PC	¥240,300	¥256,000	
89090-484	Crystallizing Dish 90x50mm 6入	¥100,700	¥108,900	
89091-138	PYREX Separatory Funnels125mL 6個	¥427,800	¥449,000	
89091-492	PYREX Glass Beads 5mm径 0.45kgx4P	¥572,100	¥600,500	
89091-610	PYREX Stopper 24/40 40mm 6本入	¥127,500	¥137,800	
89092-718	Packer Bottles 15mL 48	¥71,600	¥77,800	
89092-870	Plastic Foam Plugs 40	¥39,900	¥43,800	
89093-134	バケツ(樹脂製安全バケツ) 0.6ガロン 30入	¥89,200	¥95,100	
89093-138	バケツ(樹脂製安全バケツ) 3.5ガロン 6入	¥52,200	¥57,100	
89093-140	バケツ(樹脂製安全バケツ) 18.9L 4個	¥49,500	¥53,800	
89093-142	バケツ(樹脂製安全バケツ) 6.5ガロン 2個	¥26,700	¥29,100	
89093-450	Poxygrid 500mL Centrifuge TubeRack	¥25,300	¥27,400	
89093-590	96ウェルプレート V底 5枚x20PK	¥113,900	¥122,600	
89093-592	96ウェルプレート(ヒートシールド/フタ付)100枚	¥115,500	¥124,300	
89093-610	96ウェルプレート(ヒートシールド/フタ付)ST50枚	¥126,300	¥135,800	
89093-612	96ウェルプレート(ヒートシールド/フタ付)ST50枚	¥127,400	¥137,100	
89093-680	96ウェルプレート フタ付 ST 50枚	¥140,300	¥150,900	
89093-690	96well plate ST	¥385,700	¥414,900	
89093-830	VWR TraceClean 20mL NPC Vials 72本	¥95,700	¥103,900	
89093-838	VWR TraceClean 20mL QA Vials 72本	¥190,900	¥207,300	
89093-848	VWR TraceClean 20mL NPC Vials 72本	¥66,200	¥71,900	
89093-850	ガラスバケツ 40ml Clear 72個	¥87,600	¥95,200	
89093-872	TraceClean 60 mL Vials, 72	¥115,000	¥125,000	
89093-924	ガラス瓶 ホスト丸瓶(オートクレーブ付)	¥21,600	¥23,500	
89093-938	VWR TraceClean広口ガラス瓶 60mL 24	¥65,800	¥71,500	
89093-944	VWR TraceClean 広口ガラス瓶 250mL 12	¥36,600	¥39,800	
89093-952	VWR TraceClean 広口ガラス瓶 250mL 12	¥40,400	¥43,800	
89093-960	VWR TraceClean 広口ガラス瓶 250mL 12	¥41,300	¥44,900	
89093-962	VWR TraceClean褐色広口ボトル 946ml	¥45,900	¥49,900	
89093-964	VWR TraceClean褐色広口ボトル 946ml	¥53,300	¥57,900	
89093-966	VWR TraceClean褐色広口ボトル 946ml	¥56,100	¥61,000	
89093-970	VWR TraceClean広口ガラス瓶125mL 24本	¥59,100	¥64,200	
89093-994	VWR TraceClean 広口ガラス瓶 250mL 12	¥31,100	¥33,800	
89094-002	VWR TraceClean 広口ガラス瓶 250mL 12	¥35,800	¥38,900	
89094-010	VWR TraceClean 広口ガラス瓶 250mL 12	¥41,200	¥44,800	
89094-018	VWR TraceClean広口ガラス瓶 120mL24個	¥53,800	¥58,500	
89094-024	広口ガラス瓶 120ml 24個	¥65,400	¥74,500	
89094-034	VWR TraceClean褐色広口ボトル 120ml	¥49,000	¥53,200	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
89094-036	VWR TraceClean褐色広口ボトル 250ml	¥63,500	¥69,000	
89094-038	VWR TraceClean褐色広口ボトル 500ml	¥33,700	¥36,600	
89094-040	VWR TraceClean褐色広口ボトル 946ml	¥41,400	¥45,000	
89094-042	VWR TraceClean褐色広口ボトル 120ml	¥61,300	¥70,400	
89094-044	VWR TraceClean褐色広口ボトル 250ml	¥71,500	¥77,700	
89094-046	VWR TraceClean褐色広口ボトル 500ml	¥45,400	¥49,300	
89094-048	VWR TraceClean褐色広口ボトル 946ml	¥49,600	¥53,900	
89094-050	VWR TraceClean褐色広口ボトル 120ml	¥68,200	¥77,700	
89094-052	VWR TraceClean褐色広口ボトル 250ml	¥88,700	¥96,400	
89094-054	VWR TraceClean褐色広口ボトル 500ml	¥49,000	¥53,300	
89094-056	VWR TraceClean褐色広口ボトル 946ml	¥58,700	¥63,800	
89094-108	VWR TraceClean広口ガラス瓶500mL 24個	¥113,600	¥122,200	
89094-656	VWRテイスホ-サ-アルサ-ル- 100mL ST	¥64,900	¥68,200	
89094-658	VWRテイスホ-リサ-ル- 100mL ST 100個	¥80,700	¥84,700	
89094-660	VWRテイスホ-リサ-ル- 100mL 50個x8	¥209,800	¥143,100	
89094-662	VWRテイスホ-サ-アルサ-ル- 25mL ST 200	¥64,300	¥67,500	
89094-664	VWRテイスホ-サ-アルサ-ル- 25mL ST 100	¥84,800	¥87,200	
89094-668	VWRテイスホ-サ-アルサ-ル- 25mL 100個	¥30,000	¥30,800	
89094-670	VWRテイスホ-リサ-ル- 25mL 5個x40	¥81,200	¥81,000	
89094-672	VWRテイスホ-リサ-ル- 25mL ST 100個	¥72,700	¥72,500	
89094-674	VWRテイスホ-リサ-ル- 50mL 50個x8	¥129,400	¥135,800	
89094-676	VWRテイスホ-リサ-ル- 50mL 5個x40	¥61,100	¥64,200	
89094-680	VWRテイスホ-リサ-ル- 50mL ST 5個x40袋	¥58,600	¥61,500	
89094-682	VWRテイスホ-リサ-ル- 50mL ST 100個	¥72,100	¥71,900	
89094-688	VWRテイスホ-リサ-ル- 50mL ST 200個	¥63,200	¥66,400	
89094-738	VWR Traceable Type-K Thermometer	¥103,600	¥193,000	
89094-966	Tips Vistarak FTR 250ul 960本	¥82,200	¥89,500	
89094-974	Tips Vistarak FTR 1250ul 960本	¥92,700	¥100,300	
89095-258	VWR Erlenmeyer Flasks 125mL	¥260,400	¥271,900	
89095-260	VWR Erlenmeyer Flasks 125mL standa	¥214,300	¥223,800	
89095-266	VWR Erlenmeyer Flasks 250mL vente	¥247,100	¥258,100	
89095-268	VWR Erlenmeyer Flasks 250mL standa	¥240,300	¥250,900	
89095-274	VWR Erlenmeyer Flasks 500mL vented	¥192,700	¥201,300	
89095-276	VWR ErlenmeyerFlasks 500mLstandard	¥189,200	¥197,500	
89095-282	VWR Erlenmeyer Flasks 1L Vented	¥273,100	¥285,200	
89095-284	VWR Erlenmeyer Flasks 1L Standard	¥265,700	¥277,500	
89095-594	VWR Thermometers 1	¥19,200	¥20,900	
89095-630	VWR Certified Thermometers	¥112,300	¥122,500	
89095-670	VWR Freezer Thermometer 1本	¥10,200	¥11,600	
89095-672	VWR Freezer Thermometer 1本	¥24,200	¥27,400	
89095-686	VWR Bottle Thermometers 15-50°C 1	¥10,300	¥11,700	
89095-702	VWR Pocket Thermometers 1	¥6,600	¥7,500	
89095-704	VWR Pocket Thermometers 1	¥4,600	¥5,200	
89095-996	Handle Truclean ADJ 72-144巾	¥73,700	¥78,900	
89095-998	Handle Truclean ADJ 100-216巾	¥93,500	¥100,100	
89096-004	TruCLEAN Cover Perforated PK25	¥106,100	¥124,900	
89096-006	Wringer Truclean Stainless Stl	¥291,300	¥311,900	
89096-030	TruCLEAN Mop Fused Cover 25入	¥329,800	¥353,100	
89096-036	Mop Antimicrobial Microfiber 12入	¥138,300	¥148,000	
89096-038	Mop Cleanroom 12入	¥202,200	¥216,500	
89096-040	Mop Microfiber Blue 12入	¥138,000	¥147,700	
89096-722	リアルタイムPCR用8連1-7° +7ラットキャップ° 120	¥59,900	¥60,900	
89097-388	Cylinder TD Kimcote CS1	¥116,500	¥127,100	
89097-534	VWR Sleeve Stoppers 100個	¥20,700	¥22,600	
89097-540	VWR Sleeve Stoppers 100個	¥31,400	¥34,300	
89097-556	VWR Sleeve Stoppers 19/22 100個	¥48,900	¥53,400	
89097-558	VWR Sleeve White Stoppers 100個	¥70,800	¥77,300	
89097-816	VWRテイスホ-サ-アルサ-ル- 117 ST 100	¥61,500	¥66,200	
89097-818	VWRテイスホ-リサ-ル- 1000枚/巻	¥6,400	¥7,800	
89097-822	VWRテイスホ-リサ-ル- 1000枚/巻	¥6,400	¥7,800	
89097-824	VWRテイスホ-リサ-ル- 1000枚/巻	¥6,400	¥7,800	
89097-826	VWRテイスホ-リサ-ル- 1000枚/巻	¥6,400	¥7,800	
89097-828	VWRテイスホ-リサ-ル- 1000枚/巻	¥6,400	¥7,800	
89097-830	VWRテイスホ-リサ-ル- 1000枚/巻	¥6,400	¥7,800	
89097-882	VWR Tamper-Evident Labels 1巻(500	¥16,600	¥18,100	
89097-892	VWR Protective Labels 500枚	¥30,100	¥38,300	
89097-902	VWRテイスホ-リサ-ル- 長さ12.7m 6本	¥17,200	¥18,700	
89097-912	VWRテイスホ-リサ-ル- 長さ12.7m 6本	¥17,200	¥18,700	
89097-916	VWRテイスホ-リサ-ル- 初付 長さ12.7m 6本	¥17,200	¥18,700	
89097-920	テイスホ-リサ-ル- (5本用)レイホ-ル° 24本	¥55,300	¥60,100	
89097-926	VWRテイスホ-リサ-ル- 1000枚/巻	¥7,700	¥8,400	
89097-928	VWRテイスホ-リサ-ル- 117- 長さ12.7m 6本	¥17,200	¥20,300	
89097-930	VWRテイスホ-リサ-ル- グリーン 6巻	¥17,200	¥20,300	
89097-934	VWRテイスホ-リサ-ル- 初付 6巻	¥17,200	¥20,300	
89097-936	VWRテイスホ-リサ-ル- ブルー 6巻	¥17,000	¥18,500	
89097-938	VWRテイスホ-リサ-ル- ピンク 6巻	¥17,200	¥20,300	
89097-954	VWRテイスホ-リサ-ル- 初付 6巻入	¥38,400	¥41,800	
89097-958	VWRテイスホ-リサ-ル- グリーン 6巻入	¥38,400	¥41,800	
89097-974	VWRテイスホ-リサ-ル- 4巻	¥16,200	¥17,700	
89097-976	VWRテイスホ-リサ-ル- 4巻	¥16,200	¥17,700	
89097-982	VWRテイスホ-リサ-ル- 4巻	¥16,200	¥17,700	
89097-984	VWRテイスホ-リサ-ル- 4巻	¥16,200	¥17,700	
89097-986	VWRテイスホ-リサ-ル- 4巻	¥16,100	¥17,500	
89097-990	VWRテイスホ-リサ-ル- 幅19mm×長12.7m	¥55,700	¥66,000	
89097-996	VWRテイスホ-リサ-ル- 117- 4巻	¥16,200	¥17,600	
89097-998	VWRテイスホ-リサ-ル- 4巻	¥16,200	¥17,600	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
89098-000	VWRラベリングテープ 4巻	¥16,200	¥17,600	
89098-002	VWRラベリングテープ 4巻	¥16,200	¥17,600	
89098-010	VWRラベリングテープ 4巻	¥16,200	¥17,700	
89098-016	VWRラベリングテープ 幅19mm×全長54.9m 4	¥31,800	¥37,100	
89098-018	VWRラベリングテープ 幅19mm×全長54.9m 4	¥31,800	¥37,100	
89098-024	VWRラベリングテープ 初付 4巻	¥31,700	¥36,900	
89098-028	VWRラベリングテープ グリーン 4巻	¥31,800	¥37,100	
89098-030	VWRラベリングテープ 幅19mm×全長54.9m 4	¥31,700	¥36,900	
89098-032	VWRラベリングテープ 緑 4巻	¥31,800	¥37,100	
89098-034	VWRラベリングテープ ブルー 4巻	¥31,800	¥37,100	
89098-036	VWRラベリングテープ 幅19mm×全長54.9m 4個	¥31,700	¥36,900	
89098-054	VWRラベリングテープ アクア 3巻	¥15,700	¥17,000	
89098-058	VWRラベリングテープ 初付 幅25.4mm	¥15,500	¥18,400	
89098-062	General-Purpose 赤 用ラベリングテープ	¥55,100	¥65,200	
89098-068	VWRラベリングテープ 白 3巻	¥15,600	¥16,900	
89098-070	VWRラベリングテープ グリーン 3巻	¥15,700	¥17,000	
89098-072	General-Purpose Labeling Tape 3	¥15,700	¥17,000	
89098-074	VWRラベリングテープ 緑 3巻	¥15,600	¥18,500	
89098-076	VWRラベリングテープ ブルー 3巻	¥15,600	¥16,900	
89098-078	VWRラベリングテープ ピンク 3巻	¥15,700	¥17,000	
89098-082	VWRラベリングテープ サモン 3巻	¥15,700	¥17,000	
89098-092	VWR ラベリングテープ 54.9m/幅25.4mm	¥32,900	¥35,900	
89098-094	VWRラベリングテープ 幅25.4mm 長さ54.9m	¥32,400	¥38,600	
89098-096	VWR ラベリングテープ 54.9m/幅25.4mm	¥32,600	¥35,600	
89098-098	VWR ラベリングテープ 54.9m/幅25.4mm	¥32,900	¥35,900	
89098-104	VWR ラベリングテープ 54.9m/幅25.4mm	¥32,900	¥35,900	
89098-106	VWR ラベリングテープ 54.9m/幅25.4mm	¥32,800	¥35,700	
89107-628	Tissue Grinder Sterile CS10	¥169,100	¥185,300	
89107-632	VWR DisposablePetriDishes 100x25mm	¥53,600	¥56,300	
89107-766	Cleanroom Marker Black	¥88,800	¥96,000	
89107-768	VWR Protection SPP Bouffant Caps	¥60,700	¥61,200	
89107-770	VWR Basic Protection Cap 1000枚	¥71,400	¥72,000	
89107-778	VWR Basic Protection Cover 1000枚	¥43,900	¥44,200	
89108-148	Cryotube 5ml Selfstand	¥320,100	¥366,300	
89125-616	DuPont Clean-Processedカバーオール L	¥353,000	¥377,900	
89125-646	クリーニング用袖カバー ST	¥200,000	¥214,100	
89125-846	VWR SpecimenContainer 237ml ST 100	¥66,100	¥73,200	
89126-578	VWR TraceClean ストレートカウストル 60mL	¥63,000	¥68,500	
89126-580	VWR TraceClean ストレートカウストル 120mL	¥79,800	¥86,700	
89126-584	VWR TraceClean ストレートカウストル 500mL	¥49,000	¥53,200	
89126-694	96-Well PCR Plates 100	¥125,700	¥134,500	
89126-722	Microcentrifuge Tubes w/caps 250	¥12,100	¥12,900	
89126-754	Whirl-Pak Filter Bags 2.72L 100	¥128,400	¥134,700	
89127-172	DuPont タイック IsoClean カバーオール 25	¥319,500	¥342,000	
89127-174	DuPont タイック IsoClean カバーオール 25	¥338,700	¥362,600	
89127-176	DuPont タイック IsoClean カバーオール 25	¥273,900	¥293,200	
89127-178	DuPont タイック IsoClean カバーオール 25	¥273,900	¥293,200	
89127-212	DuPont Clean-Processedカバーオール XXL	¥372,800	¥399,100	
89127-218	DuPont Clean-Processedカバーオール 3XL	¥433,400	¥463,900	
89127-224	DuPont Clean-Processedカバーオール 4XL	¥431,400	¥461,800	
89127-226	DuPont Clean-Processedカバーオール 5XL	¥465,900	¥498,800	
89127-230	タイック IsoClaen カバーオール ST L 25枚	¥286,000	¥306,100	
89127-234	タイック IsoClaen カバーオール ST M 25枚	¥303,700	¥325,200	
89127-236	DuPont Clean-Processedカバーオール M	¥353,000	¥377,900	
89127-240	DuPont Clean-Processedカバーオール S	¥338,200	¥362,100	
89127-244	タイック IsoClaen カバーオール ST XL 25枚	¥293,100	¥313,700	
89127-246	DuPont Clean-Processedカバーオール XL	¥353,000	¥377,900	
89127-324	タイック IsoClaen フックカバー ST 100枚	¥202,800	¥217,100	
89127-326	タイック IsoClaen フックカバー 100枚	¥124,200	¥132,900	
89127-328	DuPontタイックIsoCleanフックカバー ST100	¥201,300	¥215,500	
89127-336	DuPontタイックフックカバー L 45.7cm 100	¥229,500	¥245,700	
89127-340	DuPontタイックフックカバー M 45.7cm 100	¥211,300	¥226,200	
89127-350	DuPontタイック腕カバー 45.7cm 100枚入	¥106,100	¥113,600	
89127-354	DuPontタイック腕カバー 45.7cm 100枚入	¥243,200	¥260,400	
89129-470	バリエーショントピック (PP製厚さ2ミル)	¥107,600	¥129,800	
89129-598	DuPont タイック IsoClean カバーオール 25	¥184,500	¥197,500	
89129-600	DuPont タイック IsoClean カバーオール 25	¥254,100	¥272,000	
89129-614	DuPont タイック IsoClean カバーオール 25	¥301,400	¥322,700	
89129-686	DuPontタイックフックカバー XL 45.7cm 100	¥186,800	¥200,000	
89130-774	Whirl-Pak Filter Bags 2.04L 250	¥177,800	¥205,300	
89130-776	サブスクリプション 5 441mL 100枚	¥100,900	¥105,900	
89130-828	滅菌済サブスクリプション 2.04L 19x38.1 250	¥116,700	¥122,500	
89132-028	滅菌済サブスクリプション 2.72L 25.4x38.1cm	¥112,100	¥117,700	
89132-320	Kelly Forceps ストレート 14.0cm	¥28,200	¥29,600	
89133-752	Pur-Wraps Long Cotton Swab 200本	¥14,000	¥14,800	
89133-910	Tube PCR Reaction Strip w/cap 125	¥74,300	¥79,900	
89134-164	BluntNeedles w/LuerStubAdapters 25	¥20,800	¥21,800	
89134-166	BluntNeedles w/LuerStubAdapter 100	¥50,900	¥53,400	
89134-190	30G BLUNT NEEDLE 0.5IN 100BK 100	¥50,800	¥53,300	
89134-428	VWR Polyolefin Films w/Silicone 10	¥126,300	¥135,800	
89135-370	拭き取り用拭き布 3 スティック 20リットル用 SZ6 3	¥393,000	¥420,700	
89135-430	PCRプレート付 2ml 1000	¥94,300	¥101,400	
89135-450	PCR Plate box(rack) PK5	¥74,200	¥79,900	
89136-072	ACCU-JET PRO用フィルター 1個	¥22,400	¥41,000	
89136-568	Biotix Filter Pipet Tips300ul 960	¥79,700	¥85,300	
89136-576	Biotix Filter Pipet Tips 10ul 960	¥79,300	¥84,900	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
89136-580	Biotix Filter Pipet Tips100ul 960	¥80,900	¥86,600	
89136-652	Mouse Restrainer Scapular Acc 1個	¥45,900	¥50,600	
89136-654	Kit with Regular Swab and Cap	¥486,200	¥510,200	
89137-864	シリコンチューブ	¥71,300	¥74,900	
89137-880	Pharma50 シリコンチューブ 50フィート×1巻	¥171,300	¥179,800	
89140-162	Neptune100ulシリコンチューブ 96本×10R	¥66,900	¥70,900	
89140-164	Neptune100ulシリコンチューブ 96本×10R	¥63,700	¥67,500	
89140-196	VWR Tracable Flashing LED Timer	¥17,100	¥32,000	
89140-598	Scoop Sampling System 60ml 10個入	¥36,200	¥39,100	
89140-800	Self-Seal SterilizationPouches 200	¥47,000	¥51,600	
89140-818	STERLIN Nitrile Gloves S 200枚	¥31,400	¥30,100	
89140-820	STERLIN Nitrile Gloves M 200枚	¥31,400	¥30,100	
89141-020	BOTTLE WASHVENT ETOH 500ML 3本	¥23,500	¥25,100	
89141-046	VWR広口ボトル LDPE製(Ethanol)500mL	¥20,400	¥21,900	
89141-048	VWR Safety-Vented 広口ボトル Ethano	¥14,000	¥14,900	
89167-106	Round Cover Slips 1oz	¥84,500	¥102,000	
89167-108	Round Cover Slips 1oz	¥74,300	¥81,500	
89167-170	TUBE W/CAP 16x125mm PK100	¥61,800	¥73,600	
89167-176	STOPPER BUTYL RUBBER PK100	¥27,800	¥33,200	
89167-178	Anaerobic Culture Tube 100本入	¥328,700	¥357,000	
89167-180	18x150mm Glass Tube Only 100本入	¥317,700	¥345,100	
89167-184	Aluminium Seal(封筒) 100個入	¥6,900	¥7,500	
89168-688	RH/TEMP SD CARD DATALOGGER	¥230,600	¥250,500	
89168-870	冷凍/冷蔵庫用 温度計 1	¥15,400	¥16,400	
89170-926	吸収シート 50mL	¥43,000	¥48,300	
89170-934	4-Bay Aquic-Pak Absorbent Pak 900	¥244,400	¥274,500	
89170-936	6-Bay Aquic-Pak Absorbent Pak 900	¥359,500	¥403,700	
89170-948	95kPa 2層キャットパッド(吸収パッド付)100	¥109,800	¥123,300	
89170-950	95kPa 2層キャットパッド 100	¥230,700	¥259,100	
89170-952	95kPa 2層キャットパッド 100	¥170,600	¥191,500	
89171-106	Ambient Gel Wraps 16	¥29,200	¥29,600	
89171-108	Insulated Foil Envelope Mailer 32	¥120,000	¥134,800	
89173-086	BlackDragonZeroパッド-フリーテックXS	¥12,800	¥13,700	
89173-088	BlackDragonZeroパッド-フリーテックS	¥12,800	¥13,700	
89173-090	BlackDragonZeroパッド-フリーテックM	¥12,800	¥13,700	
89173-092	BlackDragonZeroパッド-フリーテックL	¥12,800	¥13,700	
89173-094	BlackDragonZeroパッド-フリーテックXL	¥12,800	¥13,700	
89173-616	パッド-フリーテックロー XS 100入	¥21,800	¥22,400	
89174-468	Centrifuge Tubes Sterile 15mL 500	¥187,400	¥201,000	
89174-474	Centrifuge Tubes Sterile 50mL 300	¥161,300	¥172,800	
89174-491	Luer Slip Bulk Syringe 1mL 100	¥19,200	¥21,000	
89174-502	Syringe 5mL Luer Lock 100本/袋	¥30,500	¥36,700	
89176-472	Prep-Safe MCT Tube Floating Rack	¥9,500	¥10,200	
89182-906	96ウェル ティーグウェルプレート 560ul ST 25入	¥234,000	¥246,800	
89194-872	Flocked Swab Standard Tip ST 50入	¥13,600	¥14,900	
89194-874	Flocked Swab, Mini Tip ST 10入	¥12,300	¥13,500	
89194-893	Puritan Flocked Swabs 50本	¥11,300	¥12,400	
89198-888	Tubing Pharma65 50FT	¥194,400	¥204,000	
89198-890	Tubing Pharma65 50FT	¥224,600	¥235,700	
89199-064	Wide Mouth Clear Bottle 1L 10本入	¥235,800	¥284,700	
89199-068	Wide Mouth Clear Bottle 10L 1本入	¥423,800	¥511,600	
89199-070	Wide Mouth Clear Bottle 20L 1本入	¥917,200	¥1,107,200	
89200-976	Gosselin Dippers w/handle	¥198,900	¥205,600	
89202-966	VWR MX ControllerCirculatingBath 7	¥1,551,800	¥2,740,000	
89203-942	拭き取り紙 A5 パッド - S/M 200枚	¥534,900	¥572,700	
89203-944	拭き取り紙 A5 パッド - Univ 200枚	¥569,400	¥609,600	
89203-946	拭き取り紙 A5 パッド - XL/XXL 200枚	¥575,300	¥615,900	
89205-946	LabMarker Gamma Irradiated Black	¥55,400	¥59,900	
89209-458	Self-Wringing Cleanroom Sponge Mop	¥189,700	¥259,300	
89213-042	Greenbox Shipping box 1	¥316,800	¥339,500	
89215-572	V30 Nemesis Clear Lens	¥81,200	¥87,000	
89215-758	Test Tube Racks Nylon 13mm 1	¥18,200	¥19,600	
89215-762	Test Tube Racks Nylon 13mm 1	¥18,300	¥19,700	
89215-766	Test Tube Racks Nylon 16mm 1	¥18,800	¥20,200	
89215-768	Test Tube Racks Nylon 16mm 1	¥18,800	¥20,300	
89215-770	Test Tube Racks Nylon 16mm 1	¥18,800	¥20,300	
89215-772	Test Tube Racks Nylon 16mm 1	¥18,800	¥20,300	
89215-774	Test Tube Racks Nylon 16mm 1	¥18,800	¥20,300	
89215-776	Test Tube Racks Nylon 16mm 1	¥18,900	¥20,400	
89215-820	Test Tube Racks 13mm 1	¥15,100	¥16,200	
89217-468	BT Barrier Pipet Tips 1000ul 96本x	¥72,300	¥76,700	
89218-057	KIMWIPES Delicate Task Wipers 8400	¥188,600	¥204,200	
89218-066	SteriPettes Non-Plugged	¥54,400	¥59,400	
89218-068	SteriPettes Cotton-Plugged	¥93,100	¥101,500	
89218-294	384-Well PCR Plates 30ul 10	¥27,000	¥27,500	
89218-750	Double Ended Pestle 12	¥28,200	¥30,300	
89219-270	Tri-bevel Needles 102mm,滅菌 100本	¥41,300	¥43,400	
89219-272	Vet Premium Hypodermic Needles 25	¥36,100	¥37,900	
89219-276	Tri-bevel Needles 89mm, 滅菌 100本	¥23,600	¥24,700	
89219-589	PYREX Fiber Glass Roving	¥865,400	¥935,600	
89220-070	Fused-Insert Vials Screw Caps 100	¥176,000	¥193,500	
89224-424	Citrate Azidetween Carbon 2.5kg	¥285,700	¥320,800	
89225-738	Pikovskayas Agar 2.5kg	¥228,800	¥256,700	
89225-740	Pikovskayas Agar 5kg	¥451,900	¥497,600	
89230-600	BlackDragonZeroパッド-フリーテックロー	¥10,600	¥8,600	
89230-602	BlackDragonZeroパッド-フリーテックロー	¥10,300	¥11,000	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
89230-604	BlackDragonZeroパ [®] ウ [®] -フリ-ニトリルカ [®] ロー [®]	¥10,300	¥11,000	
89230-606	BlackDragonZeroパ [®] ウ [®] -フリ-ニトリルカ [®] ロー [®]	¥10,300	¥11,000	
89230-608	BlackDragonZeroパ [®] ウ [®] -フリ-ニトリルカ [®] ロー [®]	¥10,300	¥11,000	
89230-610	BlackDragonZeroパ [®] ウ [®] -フリ-ニトリルカ [®] ロー [®]	¥10,300	¥11,000	
89230-648	80Well Reversible Tube Rack 5	¥43,600	¥46,600	
89233-128	VWR Disposable PP Spatulas 300	¥15,000	¥17,800	
89233-356	Mini Ice Pan with Lid 4L	¥85,100	¥90,700	
89233-358	Mini Ice Pan with Lid 4L	¥84,700	¥90,200	
89233-360	Mini Ice Pan with Lid 4L	¥84,700	¥90,200	
89233-488	給餌針(滅菌済/テ [®] イス [®] -サ [®] ブル) 20G	¥436,200	¥457,800	
89233-508	給餌針(滅菌済/テ [®] イス [®] -サ [®] ブル) 14G	¥360,200	¥378,100	
89233-510	給餌針(滅菌済/テ [®] イス [®] -サ [®] ブル) 14G	¥418,000	¥438,800	
89233-524	給餌針(滅菌済/テ [®] イス [®] -サ [®] ブル) 15G	¥316,700	¥332,400	
89233-802	Protexis テ [®] ツク [®] 滅菌カ [®] ロー [®] 100枚	¥227,400	¥243,400	
89233-818	テ [®] ツク [®] 滅菌カ [®] ロー [®] 200双	¥916,100	¥865,400	
89233-820	テ [®] ツク [®] 滅菌カ [®] ロー [®] 200双	¥916,100	¥865,400	
89233-986	Sense-able Detectable Scoops	¥355,200	¥378,400	
89233-994	Liquid Nitrogen Cooled Mini Mortar	¥322,500	¥345,900	
89234-142	Laboratory Hypodermic Needles 33G	¥196,600	¥219,900	
89234-146	Laboratory Hypodermic Needles 33G	¥215,400	¥241,000	
89234-192	Deflected Point Needles 12入	¥208,300	¥233,100	
89234-224	Deflected Point Needles 12入	¥187,500	¥196,800	
89234-824	液体洗剤 Contrex AL 950mL	¥18,900	¥20,300	
89236-002	VWRフリ-サ [®] -ホ [®] ツク [®] スチ [®] ル 5個入	¥28,100	¥30,200	
89236-110	ス [®] チ [®] ル(ホ [®] ツク [®] 製/滅菌済み個別包装)	¥55,300	¥58,700	
89236-678	テ [®] イス [®] -サ [®] ブル [®] ハ [®] ット 100mL ST 50本	¥181,600	¥204,000	
89238-034	Trigger Sprayer with Adapter	¥22,400	¥24,000	
89259-898	Aluminum Test Tube Rack 24本用 1	¥12,400	¥13,300	
89259-944	DissectingForceps w/StraightTips 1	¥3,900	¥4,200	
89259-954	VWR Filter Forceps 114mm 1	¥2,500	¥2,700	
89260-602	Switch-Grid Test Tube Racks 1個	¥16,100	¥17,200	
89260-788	Qlicksmart Blade Removal System 1	¥39,300	¥39,600	
89305-304	VWR Disposable Paper Lab Funnel 25	¥12,600	¥14,300	
89370-260	カ [®] ラス管カ [®] ッター (Griffin Type)	¥21,100	¥24,100	
89399-924	Pipet Pasteur Boro 9inch	¥65,500	¥71,400	
89400-380	Acid Gas Cartridge Pack of 2	¥17,200	¥18,500	
89401-562	50ml チ [®] ュー [®] フラットキャ [®] ップ [®] ST 500本入	¥68,500	¥71,900	
89402-798	73-701 グ [®] -マ [®] シ [®] ルト [®] カ [®] ロー [®] 200双	¥351,100	¥375,800	
89402-800	73-701 グ [®] -マ [®] シ [®] ルト [®] カ [®] ロー [®] 200双	¥349,900	¥374,600	
89402-802	73-701 グ [®] -マ [®] シ [®] ルト [®] カ [®] ロー [®] 200双	¥343,900	¥368,100	
89402-804	73-701 グ [®] -マ [®] シ [®] ルト [®] カ [®] ロー [®] 200双	¥343,900	¥368,200	
89402-806	73-701 グ [®] -マ [®] シ [®] ルト [®] カ [®] ロー [®] 200双	¥343,900	¥368,200	
89402-808	73-701 グ [®] -マ [®] シ [®] ルト [®] カ [®] ロー [®] 200双	¥349,900	¥374,600	
89402-810	73-701 グ [®] -マ [®] シ [®] ルト [®] カ [®] ロー [®] 200双	¥349,900	¥374,600	
89403-312	Reclosable Biohazard Bags 1000	¥66,200	¥70,500	
89403-846	タ [®] イ [®] ン [®] チ [®] ュー [®] 1巻	¥21,300	¥23,000	
89404-162	タ [®] イ [®] ン [®] チ [®] ュー [®] 1巻	¥14,000	¥15,100	
89404-300	タ [®] イ [®] ン [®] チ [®] ュー [®] 1巻	¥236,500	¥248,200	
89408-048	Trypic Soy Broth 100mL	¥98,700	¥109,300	
89412-072	Mouse Igloo 1個	¥3,200	¥3,400	
89415-020	ニトリルカ [®] ロー [®] S 1000枚	¥224,200	¥228,200	
89415-022	ニトリルカ [®] ロー [®] M 1000枚	¥224,200	¥228,200	
89415-024	ニトリルカ [®] ロー [®] L 1000枚	¥224,200	¥228,200	
89415-026	ニトリルカ [®] ロー [®] XL 1000枚	¥224,200	¥228,200	
89415-028	ニトリルカ [®] ロー [®] 2XL 1000枚	¥224,200	¥228,200	
89427-828	O-Rings Kalrez 1入	¥42,400	¥59,000	
89428-072	Funnel Addition 24/40 500ml 1個	¥378,700	¥445,600	
89428-610	KIMTECHPURE SleeveProtector ST 100	¥246,900	¥264,300	
89428-758	Tubing Sizer, Metric Gauge 1	¥34,500	¥37,000	
89429-048	48ツル Thin-Walled PCR Plate 50枚入	¥53,000	¥57,100	
89493-012	Integrated Hood&Mask ST 75	¥480,900	¥514,800	
89495-444	VWR Cuff Tape 1巻入	¥27,000	¥29,000	
89495-448	VWR Extreme Temperature Tape	¥19,400	¥20,800	
89496-702	Insulated Shipper case of 2	¥134,600	¥144,200	
89496-718	Insulated Shipper case of 2	¥128,900	¥138,200	
89496-746	Insulated Shipper case of 4	¥150,900	¥161,700	
89497-632	Fine Tip Transfer Pipet 400本x8入	¥147,700	¥174,800	
89497-676	ト [®] ランス [®] ア [®] -ピ [®] ハ [®] ット 400本	¥82,300	¥88,600	
89497-720	Exact Volume Transfer Pipets 75ul	¥157,200	¥169,100	
89497-728	Self Standing Sample Tubes 5mL 100	¥124,200	¥136,300	
89497-738	Self Standing SampleTubes 5mL 1000	¥171,700	¥188,300	
89497-740	Self Standing Sample Tubes 5mL 500	¥100,400	¥110,100	
89498-054	Low Profile Grip Rack Blue 1	¥18,100	¥19,500	
89498-056	Low Profile Grip Rack Orange 1	¥18,200	¥19,500	
89498-062	Centrifuge Tube Rack	¥15,900	¥17,100	
89498-258	Spectainer Urine Containers 300個	¥389,000	¥426,600	
89498-366	Specimen Container 60ML 500個	¥173,200	¥186,800	
89499-096	MicrewtubeMicrocentrifugeTubes1000	¥80,800	¥86,900	
89499-176	MicrewtubeMicrocentrifugeTubes1000	¥116,200	¥125,000	
89511-912	Cryo-Industrial Apron 1入	¥98,000	¥109,500	
89522-664	Eyewash Station	¥43,600	¥46,700	
89522-672	BasicProtectionBouffantCaps 1000	¥73,100	¥73,600	
90000-004	ハ [®] トリ [®] イ [®] ツ [®] PS ST カ [®] -付 20枚x25袋	¥314,700	¥345,200	
93000-032	ス [®] -ル [®] -カ [®] リア [®] 遠 [®] 心 [®] チ [®] ュー [®] 50ml ST 50本x10袋	¥79,600	¥83,500	
93000-036	ス [®] -ル [®] -カ [®] リア [®] 遠 [®] 心 [®] チ [®] ュー [®] 50ml ST 25本x20ツ	¥100,200	¥105,100	
93000-680	VWRテ [®] イス [®] -サ [®] ブル [®] ハ [®] ット 10ml	¥183,300	¥199,800	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
93000-690	VWR [®] イソホ-サ [®] フ [®] 肋 [®] ラ [®] ス [®] ビ [®] ハ [®] ヲ [®] ヲ [®] 1mL	¥265,000	¥289,000	
93000-692	VWR [®] イソホ-サ [®] フ [®] 肋 [®] ラ [®] ス [®] ビ [®] ハ [®] ヲ [®] ヲ [®] 1mL	¥171,300	¥184,100	
93000-696	VWR [®] イソホ-サ [®] フ [®] 肋 [®] ラ [®] ス [®] ビ [®] ハ [®] ヲ [®] ヲ [®] 5mL	¥114,400	¥124,700	
93000-698	VWR [®] イソホ-サ [®] フ [®] 肋 [®] ラ [®] ス [®] ビ [®] ハ [®] ヲ [®] ヲ [®] 10mL	¥137,900	¥150,300	
93000-704	VWR [®] イソホ-サ [®] フ [®] 肋 [®] ラ [®] ス [®] ビ [®] ハ [®] ヲ [®] ヲ [®] 5ml	¥106,800	¥116,500	
93000-706	VWR [®] イソホ-サ [®] フ [®] 肋 [®] ラ [®] ス [®] ビ [®] ハ [®] ヲ [®] ヲ [®] 10ml	¥119,700	¥130,600	
93000-928	Sample Vial Screw Thread 15mL 144	¥130,300	¥144,300	
93000-936	Packer Bottles 250mL 84	¥133,200	¥144,700	
93001-118	VWR リアルタイムPCR用8連1-7 +7フットキャップ	¥49,200	¥50,100	
93001-762	Septa 24mm 100個入	¥33,900	¥36,800	
94002-802	Roll-O-Matic ステンレススチールチップ [®] (ストレート)	¥210,000	¥237,500	
94023-908	STERLING Nitrile-XTRA Sサイズ [®]	¥28,100	¥25,500	
94023-910	STERLING Nitrile-XTRA Mサイズ [®]	¥28,100	¥25,500	
94023-912	STERLING Nitrile-XTRA Lサイズ [®]	¥28,100	¥25,500	
95029-584	Monoject Tuberculin Syringes 240	¥247,200	¥256,900	
95029-702	Self-Contained Biological Indicator	¥97,900	¥107,300	
95029-704	Self-Contained Biological Indicator	¥254,800	¥279,500	
95035-042	Polysorbate 20 500ml	¥75,500	¥83,200	
95039-092	Disposable Biopsy Punches 25個	¥94,700	¥99,400	
95039-094	Disposable Biopsy Punches 25個	¥81,900	¥85,900	
95039-098	Disposable Biopsy Punches 25個	¥73,000	¥76,600	
95039-100	Disposable Biopsy Punches 25個	¥71,600	¥75,200	
95039-102	Disposable Biopsy Punches 25個	¥81,000	¥78,200	
95039-810	Formalin Absorbent Pads Case of 6	¥181,200	¥196,800	
95041-302	Mini Nasopharyngeal Swab ST 500	¥332,000	¥351,900	
95041-714	Alcohol Prep Pads 4000	¥86,000	¥94,300	
95042-206	VWR Double-Slide Holders/Mailers25	¥48,200	¥52,900	
95042-548	オートクレーブ [®] バ [®] ック [®] 20.3x30.5cm 100枚	¥35,400	¥38,800	
95042-552	VWRオートクレーブ [®] バ [®] ック [®] PP 100枚入	¥36,200	¥39,600	
95042-554	オートクレーブ [®] バ [®] ック [®] 30.5x61cm 100枚	¥48,300	¥55,600	
95042-556	オートクレーブ [®] バ [®] ック [®] 100枚	¥106,900	¥117,200	
95042-563	VWRオートクレーブ [®] バ [®] ック [®] PP 25.4x38.1cm	¥46,400	¥50,900	
95042-564	オートクレーブ [®] バ [®] ック [®] 20.3x25.4cm 100枚	¥37,400	¥41,100	
95042-716	Surgical Instrument Stain Remover 12	¥84,500	¥88,700	
95054-974	Normadic Shipping Box	¥128,300	¥137,500	
95054-990	Encore Economy Table Paper	¥49,200	¥53,200	
95057-225	VACUETTE EDTA [®] 1-7 [®] (K3EDTA)4mL 50	¥5,900	¥6,400	
95057-838	Premium Frozen Section Compound	¥19,100	¥20,900	
95059-670	バ [®] イ [®] ル [®] 3-ト [®] バ [®] ック [®] (3-Wall)	¥60,500	¥66,400	
95062-775	Flexible Tip Digital Thermometer	¥9,900	¥10,800	
97000-030	電極 PH ORION REFILLABLE 120mm	¥156,900	¥168,300	
97000-076	Pressure&Vacuum Station 115V	¥246,500	¥448,000	
97000-584	VWR [®] リス [®] フ [®] レ [®] キ [®] ャ [®] ハ [®] ヲ [®] ヲ [®] (マイ [®] ク [®])2.5-4.5ml	¥67,200	¥78,000	
97000-586	VWR [®] リス [®] フ [®] レ [®] キ [®] ャ [®] ハ [®] ヲ [®] ヲ [®] (セミ [®] マイ [®] ク [®])1.5-3.0ml	¥71,600	¥77,800	
97000-588	VWR [®] ア [®] リ [®] ル [®] キ [®] ャ [®] ハ [®] ヲ [®] ヲ [®] (マイ [®] ク [®])2.5-4.5ml	¥78,700	¥85,500	
97000-590	VWR [®] ア [®] リ [®] ル [®] キ [®] ャ [®] ハ [®] ヲ [®] ヲ [®] (セミ [®] マイ [®] ク [®])1.5-3.0ml	¥77,200	¥83,900	
97001-248	Sure-Grip EX Safety Cabinets 15L	¥625,500	¥656,500	
97001-314	キャ [®] ビ [®] ネット用スチール棚 34x33cm	¥26,100	¥24,200	
97001-362	防火キャ [®] ビ [®] ネット(シングルドア)57L	¥836,100	¥1,477,000	
97001-384	EXセ [®] フ [®] ィ [®] ヤ [®] キ [®] ビ [®] ネット 170L(手動ト [®] ア [®] /2シ [®] ル [®])	¥630,800	¥1,054,000	
97004-376	Jitterbug-2(7 [®] レ [®] ト [®] 2枚用)	¥1,235,600	¥2,248,000	
97004-378	Jitterbug-4(7 [®] レ [®] ト [®] 4枚用)	¥2,099,400	¥3,819,000	
97005-270	ステンレスカート 39.4x61x81.3cm 1	¥317,700	¥343,400	
97006-114	防火キャ [®] ビ [®] ネット(シングルドア)45L	¥768,200	¥1,410,000	
97006-122	Sure-GripEXセ [®] フ [®] ィ [®] ヤ [®] キ [®] ビ [®] ネット(加 [®] タ [®] ト [®] ヲ [®])	¥625,500	¥1,105,000	
97006-126	Compac キャ [®] ビ [®] ネット (one door) 1	¥815,100	¥855,500	
97006-128	防火キャ [®] ビ [®] ネット(シングルドア)57L	¥907,500	¥1,603,000	
97006-920	Centrifugal Filt 0.45UM Ptfе	¥191,700	¥207,200	
97008-196	VWR耐水性クイ [®] カ [®] ト [®] ク [®] 0-7 [®] S	¥47,400	¥50,700	
97008-198	VWR耐水性クイ [®] カ [®] ト [®] ク [®] 0-7 [®] M	¥47,400	¥50,700	
97008-202	VWR耐水性クイ [®] カ [®] ト [®] ク [®] 0-7 [®] XL	¥47,100	¥50,400	
97008-204	VWR耐水性クイ [®] カ [®] ト [®] ク [®] 0-7 [®] S	¥51,600	¥55,200	
97008-206	VWR耐水性クイ [®] カ [®] ト [®] ク [®] 0-7 [®] M	¥51,600	¥55,200	
97008-208	VWR耐水性クイ [®] カ [®] ト [®] ク [®] 0-7 [®] L	¥51,600	¥55,200	
97008-210	VWR耐水性クイ [®] カ [®] ト [®] ク [®] 0-7 [®] XL	¥51,600	¥55,200	
97008-214	VWR耐水性クイ [®] カ [®] ト [®] ク [®] 0-7 [®] M	¥58,700	¥62,900	
97008-216	VWR耐水性クイ [®] カ [®] ト [®] ク [®] 0-7 [®] L	¥58,700	¥62,900	
97008-218	VWR耐水性クイ [®] カ [®] ト [®] ク [®] 0-7 [®] XL	¥58,700	¥62,900	
97008-220	VWR耐水性クイ [®] カ [®] ト [®] ク [®] 0-7 [®] S	¥69,300	¥74,200	
97008-222	VWR耐水性クイ [®] カ [®] ト [®] ク [®] 0-7 [®] M	¥69,300	¥74,200	
97008-224	VWR耐水性クイ [®] カ [®] ト [®] ク [®] 0-7 [®] L	¥69,800	¥74,700	
97008-226	VWR耐水性クイ [®] カ [®] ト [®] ク [®] 0-7 [®] XL	¥69,300	¥74,200	
97008-228	VWRクイ [®] カ [®] ト [®] ク [®] 0-7 [®] (ウォ [®] タ [®] ー [®] フ [®] ル [®]) S	¥50,500	¥54,100	
97008-230	VWRクイ [®] カ [®] ト [®] ク [®] 0-7 [®] (ウォ [®] タ [®] ー [®] フ [®] ル [®]) M	¥50,500	¥54,100	
97008-232	VWRクイ [®] カ [®] ト [®] ク [®] 0-7 [®] (ウォ [®] タ [®] ー [®] フ [®] ル [®]) L	¥50,500	¥54,100	
97008-234	VWRクイ [®] カ [®] ト [®] ク [®] 0-7 [®] (ウォ [®] タ [®] ー [®] フ [®] ル [®]) XL	¥50,500	¥54,100	
97008-236	VWRクイ [®] カ [®] ト [®] ク [®] 0-7 [®] (ウォ [®] タ [®] ー [®] フ [®] ル [®]) S	¥54,700	¥58,500	
97008-238	VWRクイ [®] カ [®] ト [®] ク [®] 0-7 [®] (ウォ [®] タ [®] ー [®] フ [®] ル [®]) M	¥53,300	¥57,600	
97008-240	VWRクイ [®] カ [®] ト [®] ク [®] 0-7 [®] (ウォ [®] タ [®] ー [®] フ [®] ル [®]) L	¥54,700	¥58,500	
97008-242	VWRクイ [®] カ [®] ト [®] ク [®] 0-7 [®] (ウォ [®] タ [®] ー [®] フ [®] ル [®]) XL	¥54,700	¥58,500	
97008-244	VWRクイ [®] カ [®] ト [®] ク [®] 0-7 [®] (ウォ [®] タ [®] ー [®] フ [®] ル [®]) S	¥62,200	¥66,600	
97008-246	VWRクイ [®] カ [®] ト [®] ク [®] 0-7 [®] (ウォ [®] タ [®] ー [®] フ [®] ル [®]) M	¥63,500	¥68,000	
97008-248	VWRクイ [®] カ [®] ト [®] ク [®] 0-7 [®] (ウォ [®] タ [®] ー [®] フ [®] ル [®]) L	¥62,600	¥67,000	
97008-250	VWRクイ [®] カ [®] ト [®] ク [®] 0-7 [®] (ウォ [®] タ [®] ー [®] フ [®] ル [®]) XL	¥62,600	¥67,000	
97008-252	VWRクイ [®] カ [®] ト [®] ク [®] 0-7 [®] (ウォ [®] タ [®] ー [®] フ [®] ル [®]) S	¥76,400	¥81,800	
97008-254	VWRクイ [®] カ [®] ト [®] ク [®] 0-7 [®] (ウォ [®] タ [®] ー [®] フ [®] ル [®]) M	¥76,400	¥81,800	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
97008-256	VWRクワイカートクローブ (ウォータープルーフ) L	¥76,900	¥82,300	
97008-258	VWRクワイカートクローブ (ウォータープルーフ) XL	¥76,900	¥82,300	
97012-550	O-Ring Thermal Desorber Trap	¥45,500	¥83,000	
97014-322	VWR LN2 Cryogenic freezer box 12	¥39,000	¥42,200	
97017-726	VWR Medium-Duty Plastic Carts	¥473,800	¥233,100	
97017-728	VWR Medium-Duty Plastic Carts	¥437,100	¥233,100	
97021-800	Buffer Powder Pillows pH6.00 50個	¥26,700	¥34,900	
97023-356	Venusil XBP C18 HPLC Columns 5um	¥1,299,900	¥2,364,000	
97034-978	HPLC Columns 2.1x150mm	¥581,200	¥1,093,000	
97035-040	HPLC Columns 100x2.1mm	¥513,600	¥966,000	
97035-046	HPLC Columns 2.1x150mm	¥581,800	¥1,094,000	
97035-534	クランプ (セーフティメッシュ表記ストラップ付)	¥69,500	¥75,900	
97035-536	クランプ (セーフティメッシュストラップ/チェーン付) 1	¥78,900	¥82,800	
97041-236	Anti-Skid ShoeCovers Universal300	¥25,900	¥26,900	
97043-938	ULTRASONIC CLEANER 5.7L 220V	¥374,600	¥681,000	
97044-508	96-Well Microplates Inserts	¥22,800	¥24,500	
97044-562	EVA Plate Covers, Pierceable	¥7,900	¥8,500	
97044-568	EVA Plate Covers, Pierceable	¥5,500	¥5,900	
97044-944	96well Microplate 350ul	¥8,100	¥8,700	
97045-086	ASM CL 350ul INS NAT FEP/SIL 100	¥85,500	¥91,900	
97045-322	96well Microplate 350ul	¥6,900	¥7,400	
97045-378	Plate with Glass Insert 2.2mL 1枚	¥34,800	¥37,400	
97045-788	Kit 1.0ml Blue Mat	¥81,400	¥88,000	
97046-020	Seal 8mm T/S/T	¥30,600	¥32,900	
97046-054	Vial 8mm Clear Round 6x32mm 100本	¥36,600	¥39,300	
97046-200	9mm Yellow Screw Thread Cap&Septa	¥23,700	¥25,400	
97046-948	Septa 10mm T/S/T 100入	¥34,200	¥36,700	
97047-012	Limited Volume Vial Inserts 100	¥123,700	¥133,000	
97047-104	MicroLiter 11mm Crimp-Top Seals100	¥35,800	¥38,500	
97047-330	MicroLiter 11mm Crimp-Top Vials100	¥25,200	¥27,100	
97047-988	Clear Vials with Beveled Top 20ml	¥40,300	¥43,300	
97054-980	Vials 10mL	¥34,300	¥38,400	
97054-990	20mm Crimp Top Vials 10mL	¥36,900	¥41,200	
97058-256	Rheodyne Injector valve 7725 1個	¥1,414,400	¥1,573,700	
97058-340	2/6 MAN SYR W/POS SENSOR	¥1,281,900	¥2,400,000	
97060-174	Horizontal 50ml Tube Holder 1	¥77,600	¥144,000	
97061-022	EDTA Disodium Salt Dihydrate 500g	¥78,000	¥89,200	
97061-250	Dimethyl Sulfoxide 500mL	¥59,900	¥66,700	
97062-968	Destaining Bags 25	¥82,400	¥91,600	
97063-144	Carbenicillin Disodium Salt 1g	¥37,400	¥45,500	
97063-146	Carbenicillin Disodium Salt 250mg	¥25,500	¥28,200	
97066-206	VWR Vacuum Filtration System 12	¥87,500	¥94,600	
98000-200	平面温度計	¥25,100	¥45,000	
A311-00	CCK-8 Cell Counting Kit 50回用④	¥2,400	¥3,400	
AA36645K7	Water reagent 4L	¥42,700	¥48,200	
AA4549406	Silver(I) behenate 5g	¥39,000	¥41,000	
AA4549414	Silver(I) behenate 25g	¥159,900	¥180,100	
AS-08091-01	NE-NG040用 アダプター 2枚入	¥3,200	¥3,500	
AS-08091-02	NE-NG040用 アダプター 2枚入	¥3,200	¥3,500	
AS-08091-03	NE-NG040用 リターナル	¥14,300	¥17,500	
AS6123020002	NE-NG040用 Lid	¥11,300	¥12,200	
B-002	LittleGenius フラワー	¥2,300	¥2,400	
B-003	LittleGene ヲケット	¥1,900	¥2,000	
B-005	Life/GeneTouch ヒューズ	¥1,900	¥2,000	
B-006	LifePro/ECO/LittleG/GeneQ温度ヒューズ	¥1,900	¥2,000	
B-007	LittleGene 液晶	¥9,700	¥9,800	
B-008	GenePro/GeneQ 液晶カバー	¥87,000	¥88,400	
B-010	LifePro 液晶カバー	¥110,900	¥112,700	
B-011	Life/GeneTouch タッチパネル	¥25,100	¥25,500	
B-012	LifePro/Touch/GenePro 電源	¥185,900	¥188,900	
B-013	ファン LifePro/Touch GenePro/Touch	¥7,700	¥7,900	
B-014	LittleGene TE	¥42,500	¥43,200	
B-015	LittleGenius フラワーセンサー	¥23,200	¥23,600	
B-016	LifeECO/GeneTouch ネットワークセンサー	¥11,600	¥11,800	
B-017	LifeTouch/GeneTouch フラワーセンサー	¥22,000	¥22,400	
B-018	タッチパネル	¥1,900	¥2,000	
B-019	LittleGene2 フラワーモジュール24	¥119,800	¥121,700	
B-020	LittleGene2 ネットワークモジュール	¥41,000	¥41,600	
B-021	LittleGene電源基板	¥260,500	¥264,700	
B-022	LittleGene カバーユニット	¥83,100	¥84,400	
B-024	MB-101 フラワー基板	¥1,900	¥2,000	
B-025	LifePro TE (新:薄)	¥40,200	¥40,800	
B-029	GeneQ TE	¥59,500	¥60,500	
B-030	GeneQ フラワーカバー	¥13,500	¥13,700	
B-031	GeneQ ヒューズ	¥8,900	¥9,000	
B-032	GeneQ カバー	¥16,200	¥16,500	
B-033	GeneQ 薄膜スイッチ	¥9,300	¥9,400	
B-034	LittleGenius サポートボード	¥5,800	¥5,900	
B-035	LittleGenius ネットワーク	¥1,500	¥1,600	
B-036	LittleGene フラワースイッチ	¥12,000	¥12,200	
B-037	LittleGene ヒューズ	¥8,900	¥9,000	
B-039	LittleGene フラワーカバー	¥2,100	¥2,200	
B-040	LittleGene 薄膜スイッチ	¥12,000	¥12,200	
B-041	GenePro 薄膜スイッチ	¥22,400	¥22,800	
B-043	LifePro ヒューズ	¥9,700	¥9,800	
B-045	LifePro フラワーカバー	¥4,600	¥4,700	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
B-047	LifePro TE (旧:厚)	¥40,200	¥40,800	
B-048	LifeECO フロントパネル 8ピン	¥231,300	¥235,000	
B-049	LifeECO 接続ケーブル ワイヤレスなし	¥106,300	¥108,000	
B-050	LifeTouch/ECO/GenoTouch CPUボード	¥368,700	¥374,700	
B-051	LifeECO 電源(ハーフサイズ)	¥152,700	¥155,100	
B-052	LifeECO CPUトランスファボード	¥4,600	¥4,700	
B-053	LifeECO TE	¥38,600	¥39,300	
B-054	LifeECO ヒーター	¥9,700	¥9,800	
B-055	ファン LifeECO	¥7,700	¥7,900	
B-057	ThermocellMixingBlock デバイスレ-	¥30,900	¥31,400	
B-058	GeneQ フロントパネル24	¥135,300	¥137,500	
B-059	GeneQ 接続ケーブル	¥70,000	¥71,100	
B-060	GeneQ CPUボード	¥132,600	¥134,700	
B-061	GeneQ 電源基板	¥299,900	¥304,800	
B-062	GeneQ フロントパネル18	¥133,300	¥135,500	
B-063	LittleGene フロントパネル16	¥115,900	¥117,800	
B-064	LittleGene/Genius接続ケーブル	¥41,000	¥41,600	
B-065	LittleGenius フロントパネル	¥116,300	¥118,200	
B-066	LittleGene CPUボード	¥100,500	¥102,100	
B-067	LifePro/GenoPro 駆動基板	¥282,100	¥286,700	
B-068	GenePro CPUボード	¥253,100	¥257,200	
B-069	GenePro フロントパネル96	¥328,500	¥333,800	
B-070	GenePro 接続ケーブル96	¥121,700	¥123,700	
B-071	GenePro フロントパネル48	¥344,000	¥349,500	
B-074	LifePro CPUボード	¥253,100	¥257,200	
B-077	LifePro/Touch 接続ケーブル	¥106,300	¥108,000	
B-078	Life/GenoTouch ドライバボード	¥282,100	¥286,700	
B-079	GenoTouch CPUボード	¥268,200	¥272,600	
B-080	LifeTouch CPUボード	¥268,200	¥272,600	
B-082	GeneQ ホットジョイモーター	¥1,500	¥1,600	
B-085	B-384GA 接続ケーブル	¥7,700	¥7,900	
B-086	LittleGenius リッドシール	¥5,800	¥5,900	
B-088	GenoTouch フロントパネル96	¥328,500	¥333,800	
B-089	GenoTouch 接続ケーブル96	¥121,700	¥123,700	
B-090	GenoTouch フロントパネル48	¥344,000	¥349,500	
B-091	GenoTouch 接続ケーブル48(左)	¥94,700	¥96,200	
B-092	GenoTouch 接続ケーブル48(右)	¥94,700	¥96,200	
B-093	LifeTouch フロントパネル96	¥251,200	¥255,300	
B-098	LifeECO ドライバボードコネクタ	¥282,100	¥286,700	
B-100	LittleGene2 TE	¥42,500	¥43,200	
B-105	LifePro/LifeTouch 接続ケーブル	¥9,700	¥9,800	
B-106	LifeECO フロントパネル	¥30,900	¥31,400	
B-107	LifeECO タッチディスプレイ	¥37,900	¥38,500	
B-109	接続ケーブル カバープレート	¥13,500	¥13,700	
B-110	LineGene96 Plus フロントパネル	¥2,521,700	¥2,562,500	
B-111	LineGene96 Plus 光電部品	¥3,250,000	¥3,302,700	
B-112	LineGene96 Plus光電部品 (PMT不含)	¥2,099,900	¥2,133,900	
B-113	LineGene96 Plus PMT	¥1,150,200	¥1,168,800	
B-114	LineGene96 Plus 光ファイバー	¥145,400	¥147,700	
B-115	LineGene96 Plus ステップモーター	¥238,200	¥242,000	
B-116	LineGene96 Plus 接続ケーブル	¥168,400	¥171,100	
B-117	LineGene96 Plus CPUボード	¥462,000	¥469,400	
B-118	LineGene96 Plus ドライバボード	¥432,000	¥439,000	
B-119	LineGene96 Plus 光電制御基板	¥276,400	¥280,800	
B-120	LineGene96 Plus TE	¥209,400	¥212,800	
B-121	LineGene96 Plus フロントパネルミスタ	¥22,000	¥22,400	
B-122	LineGene96 Plus 接続ケーブルミスタ	¥11,600	¥11,800	
B-123	LifeECO LCDディスプレイ	¥306,100	¥311,000	
B-126	GeneQ サーマットボード	¥11,600	¥11,800	
B-127	GeneQ 耐熱スポンジ	¥1,900	¥2,000	
B-128	LifeECO用ヘルシイ	¥37,700	¥38,300	
B-129	LifeECO用接続ケーブル PCB	¥1,500	¥1,600	
B-130	LifePro用ヒンジパーツ	¥25,600	¥26,000	
B-134	LifeECO カバーフィルム	¥1,400	¥1,500	
B-135	LittleGene2 フラッシュカバー	¥14,500	¥14,700	
B-136	LittleGene2 装飾カバー	¥2,400	¥2,500	
B-137	LifeECO接続ケーブル 8ピン	¥112,600	¥114,400	
B-138	LifeECO 上カバーシール	¥39,600	¥40,300	
B-139	GenePro 接続ケーブル384	¥165,200	¥167,900	
B-140	GenoTouch 384well用接続ケーブル	¥165,200	¥167,900	
B-145	LifeECO シリコンパッド 2.4.5mm厚	¥4,800	¥4,900	
B-146	LifeECO シリコンパッド 1	¥4,800	¥4,900	
B-147	GenoTouch 384well用フロントパネル	¥355,100	¥360,800	
B-149	MB-101用Cooling&Heatingユニット	¥115,900	¥117,800	
B-150	LifeECO用温度ヒューズ (12A)	¥1,900	¥2,000	
B-151	LifeECO/LifeTouch ヒンジパーツ(New)	¥25,600	¥26,000	
B-155	GenoTouch ヘルシイ	¥7,200	¥7,400	
B-157	LifeECO接続ケーブル 10ピン	¥112,600	¥114,400	
B-158	LifeECO フロントパネル 10ピン	¥231,300	¥235,000	
B-159	GenoTouch Adapter plate 1	¥32,400	¥32,900	
B-160	GenoTouchPlus フロントパネル 1	¥328,500	¥333,800	
B-161	Life/GenoTouch 液晶パネル	¥222,200	¥225,800	
B-162	シフトタッチ+用接続ケーブル A(対面左)	¥118,400	¥120,300	
B-163	シフトタッチ+用接続ケーブル B(対面右)	¥118,400	¥120,300	
B-164	LifeECO フラッシュカバー 1	¥36,200	¥36,800	
B-165	GeneExplorer ドライバボード	¥198,100	¥201,300	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
B-166	GeneExplorer 電源	¥253,600	¥257,700	
B-167	GeneExplorer デバイス	¥208,700	¥212,100	
B-168	GeneExplorer デバイスレイト	¥257,000	¥261,200	
B-169	GeneExplorer CPU コントロールボード	¥172,900	¥175,700	
B-170	GeneExplorer デバイスポート(96)	¥384,100	¥390,300	
B-171	GeneExplorer ネットワークポート(96)	¥236,700	¥240,500	
B-172	GeneExplorer ネットワーク用ヒューズ	¥2,400	¥2,500	
B-173	GeneExplorer ネットワーク用サーミスタ	¥12,100	¥12,300	
B-174	GeneExplorer ヒーティングファン	¥37,200	¥37,800	
B-175	GeneExplorer デバイスポート用サーミスタ	¥27,500	¥28,000	
B-176	GeneExplorer ハードディスク	¥50,200	¥51,100	
B-177	GeneExplorer USB ケーブル	¥15,500	¥15,700	
B-178	LifeECO デバイスモジュール Ver3.0用	¥299,500	¥304,400	
B-179	LifeECO ネットワークモジュール Ver3.0用	¥140,100	¥142,400	
B-183	GeneExplorer ネットワークアダプタ	¥58,000	¥58,900	
B-184	GeneExplorer デバイスポート	¥62,800	¥63,800	
B-185	LifeECO 電源スイッチ 1	¥2,400	¥2,500	
B-187	LifeECO用発泡スチロール緩衝材(左のみ)1	¥3,600	¥3,700	
B-188	LifeECO用発泡スチロール緩衝材(右のみ)1	¥3,600	¥3,700	
B-189	GE用ハードウェア一式(底面&背面)	¥19,300	¥19,600	
B-190	LifeECO パワーケーブル	¥4,800	¥4,900	
B-191	LifeECO パワーフィルターキット(10A)	¥12,800	¥13,000	
B-48DA	■ GeneTouch48x0.2ml デバイス	¥230,000	¥260,000	
B-48DA+	■ GeneTouch Plus48x0.2ml デバイス	¥230,000	¥260,000	
B-96GA	■ GeneTouch96x0.2ml デバイス	¥230,000	¥260,000	
BCE124I-1SJP	電子天秤 ひょう量皿付 1	¥210,000	¥232,000	2024/3/14 追加
BCF0001	Optical Calibration fixture for BP	¥12,900	¥13,700	
BDF1503	BP Gelカセット 1.5% DF 250bp-1.5kb 3	¥61,700	¥65,800	
BDF1510	BP Gelカセット 1.5%Dye-Free0.3-1.5kb 10	¥166,100	¥177,100	
BDF1550	BP Gelカセット 1.5%Dye-Free0.3-1.5kb 50	¥707,500	¥754,300	
BDF2003	BP Gelカセット 2% DF 100-600bp 3	¥61,700	¥65,800	
BDF2010	BP Gelカセット 2% DF 100-600bp 10	¥166,100	¥177,100	
BDF2050	BP Gelカセット 2% Dye-Free100-600bp 50	¥707,500	¥754,300	
BDH11685GL	Water, ASTM Type II 5gal.	¥47,100	¥52,400	
BDH5092500ML	BDH Reference Standard Buffer-12.45PH	¥13,100	¥12,900	
BDQ3003	BP Gelカセット 3% Dye-Free 100-250bp 3	¥61,700	¥65,800	
BDQ3010	BP Gelカセット 3%Dye-Free 100-250bp 10	¥166,100	¥177,100	
BEF2003	BP Gelカセット 2% Dye-Free 100-600bp 3	¥61,700	¥65,800	
BEF2010	BP Gelカセット 2% Dye-Free100-600bp 10	¥166,100	¥177,100	
BEF2050	BP Gelカセット 2% Dye-Free100-600bp 50	¥707,500	¥754,300	
BH710931	BioHit MidiPlus デバイス	¥48,000	¥50,000	
BH711002.01	交換用バッテリー ニッケル(NiMH)6V	¥9,000	¥11,000	
BH720010	Proline 0.1-2.5ul シングル	¥23,000	¥25,000	
BH720015	Proline 0.5-10ul シングル	¥23,000	¥25,000	
BH720025	Proline 5-50ul シングル	¥23,000	¥25,000	
BH720050	Proline 10-100ul シングル	¥23,000	¥25,000	
BH720060	Proline 100-1000ul シングル	¥23,000	¥25,000	
BH720070	Proline 20-200ul シングル	¥23,000	¥25,000	
BH720080	Proline 2-20ul シングル	¥23,000	¥25,000	
BH720110	Proline 1-5ml シングル	¥23,000	¥25,000	
BH720210	Proline 0.5-10ul 8ch	¥80,000	¥85,000	
BH720220	Proline 5-50ul 8ch	¥80,000	¥85,000	
BH720240	Proline 50-300ul 8ch	¥80,000	¥85,000	
BH720310	Proline 0.5-10ul 12ch	¥102,000	¥108,000	
BH720320	Proline 5-50ul 12ch	¥102,000	¥108,000	
BH720340	Proline 50-300ul 12ch	¥102,000	¥108,000	
BH721005	交換用セーフティフィルターセット2000.5000	¥3,800	¥3,900	
BH721006	交換用セーフティフィルターセット500.1000.1200	¥3,800	¥3,900	
BH721007	交換用セーフティフィルターセット100.200	¥3,800	¥3,900	
BH721008	親水性フィルターセット (50個入り)	¥3,800	¥3,900	
BH721014	親水性フィルターセット (50個入り)	¥3,800	¥3,900	
BH721015	疎水性フィルターセット 50個入 φ6.73mm	¥6,000	¥6,200	
BH721016	交換用セーフティフィルターセット500.1000.1200	¥6,000	¥6,200	
BH721017	疎水性フィルターセット 50個入 φ3.15mm	¥6,000	¥6,200	
BH721018	疎水性フィルターセット 50個入 φ2.51mm	¥6,000	¥6,200	
BH721670	Silicon Grease 1.5ml	¥11,000	¥11,500	
BH722001	デバイスチップ 固定容量 5ul	¥19,000	¥20,000	
BH722004	デバイスチップ 固定容量 10ul	¥19,000	¥20,000	
BH722010	デバイスチップ 固定容量 20ul	¥19,000	¥20,000	
BH722015	デバイスチップ 固定容量 25ul	¥19,000	¥20,000	
BH722020	デバイスチップ 固定容量 50ul	¥19,000	¥20,000	
BH722025	デバイスチップ 固定容量 100ul	¥19,000	¥20,000	
BH722030	デバイスチップ 固定容量 200ul	¥19,000	¥20,000	
BH722035	デバイスチップ 固定容量 250ul	¥19,000	¥20,000	
BH722040	デバイスチップ 固定容量 500ul	¥19,000	¥20,000	
BH722045	デバイスチップ 固定容量 1000ul	¥19,000	¥20,000	
BH722050	デバイスチップ 固定容量 2000ul	¥19,000	¥20,000	
BH722055	デバイスチップ 固定容量 5ml	¥19,000	¥20,000	
BH725010	BIOHIT m-Line(シングル) 0.1-3ul	¥38,000	¥39,000	
BH725020	BIOHIT m-Line(シングル) 0.5-10ul	¥38,000	¥39,000	
BH725030	BIOHIT m-Line(シングル) 2-20ul	¥38,000	¥39,000	
BH725050	BIOHIT m-Line(シングル) 10-100ul	¥38,000	¥39,000	
BH725060	BIOHIT m-Line(シングル) 20-200ul	¥38,000	¥39,000	
BH725070	BIOHIT m-Line(シングル) 100-1000ul	¥38,000	¥39,000	
BH725080	BIOHIT m-Line(シングル) 500-5000ul	¥38,000	¥39,000	
BH725090	BIOHIT m-Line(シングル) 1-10ml	¥38,000	¥39,000	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
BH725120	BIOHIT m-Line 0.5-10ul 8ch	¥112,000	¥120,000	
BH725130	BIOHIT m-Line 5-100ul 8ch	¥112,000	¥120,000	
BH725140	BIOHIT m-Line 30-300ul 8ch	¥112,000	¥120,000	
BH725220	BIOHIT m-Line 0.5-10ul 12ch	¥129,000	¥135,000	
BH725230	BIOHIT m-Line 5-100ul 12ch	¥129,000	¥135,000	
BH725240	BIOHIT m-Line 30-300ul 12ch	¥129,000	¥135,000	
BH725620	ピエット用リアスタッド 1	¥15,200	¥15,700	
BH728010	Prolineアース(マニュアル) 1ch 0.1-3ul	¥30,000	¥31,000	
BH728020	Prolineアース(マニュアル) 1ch 0.5-10ul	¥30,000	¥31,000	
BH728030	Prolineアース(マニュアル) 1ch 2-20ul	¥30,000	¥31,000	
BH728050	Prolineアース(マニュアル) 1ch 10-100ul	¥30,000	¥31,000	
BH728060	Prolineアース(マニュアル) 1ch 20-200ul	¥30,000	¥31,000	
BH728070	Prolineアース(マニュアル)1ch 100-1000ul	¥30,000	¥31,000	
BH728080	Prolineアース(マニュアル)1ch 500-5000ul	¥30,000	¥31,000	
BH728090	Prolineアース(マニュアル) 1ch 1-10ml	¥30,000	¥31,000	
BH728120	Prolineアース(マニュアル) 8ch 0.5-10ul	¥100,000	¥108,000	
BH728130	Prolineアース(マニュアル) 8ch 10-100ul	¥100,000	¥108,000	
BH728140	Prolineアース(マニュアル) 8ch 30-300ul	¥100,000	¥108,000	
BH728220	Prolineアース(マニュアル) 12ch 0.5-10ul	¥112,000	¥123,000	
BH728230	Prolineアース(マニュアル) 12ch 10-100ul	¥112,000	¥123,000	
BH728240	Prolineアース(マニュアル) 12ch 30-300ul	¥112,000	¥123,000	
BH728520	Prolineアース(マニュアル)1ch 10ul	¥25,600	¥26,000	
BH728590	Prolineアース(マニュアル)1ch 10mL	¥25,600	¥26,000	
BH730981	eLINE シングル充電スタンド	¥30,000	¥32,000	
BH730991	eLINE 4台掛回転充電スタンド	¥63,000	¥66,000	
BH731009	e-line電動ピエットバッテリー 4.8V	¥15,000	¥16,000	
BH780082	BIOHIT ロングチップ 1000ul 1000本	¥11,500	¥11,600	
BH780300	BIOHIT Tip 5000ul バルク 100本	¥2,900	¥3,000	
BH780300CTN	スタンダードチップ バルク 5ml 100x16	¥39,040	¥40,640	
BH780304	BIOHIT Tip 5000ul オートクループ ラック 50	¥2,600	¥2,700	
BH780304CTN	スタンダードチップ シングルラック 5ml 50x10	¥21,200	¥23,000	
BH780305	BIOHIT Tip 5000ul ST 50本	¥3,200	¥3,300	
BH780305CTN	スタンダードチップ シングルラックST 5ml 50x10	¥28,600	¥28,800	
BH780308	BIOHIT Tip 5000ul バルク 1000本	¥24,000	¥25,000	
BH780310	BIOHIT Tip 10ml 100本	¥8,800	¥9,000	
BH780310CTN	スタンダードチップ バルク 10ml 100x8	¥65,600	¥62,400	
BH783201	フィルターチップ 10ul/46mm 96本x10R	¥19,000	¥20,000	
BH783210	スタンダードチップ 96本x10R	¥11,400	¥11,600	
BH790010	BIOHIT リフィルシステムTip 10ul シングルパック	¥11,000	¥11,200	
BH790010CTN	スタンダードチップ 10ul 96x10Rx8	¥73,600	¥78,400	
BH790011	スタンダードチップ ST 10ul 96本x10R	¥15,000	¥15,600	
BH790011CTN	スタンダードチップ ST 10ul 96x10Rx8	¥104,800	¥108,800	
BH790011F	フィルターチップ 10ul 96本x10R	¥19,000	¥20,000	
BH790012	BIOHIT Tip 10ul リフィル 96本x10	¥8,300	¥8,600	
BH790012CTN	スタンダードチップ リフィル 10ul 96x10x19	¥136,800	¥143,640	
BH790013	スタンダードチップ 10ul リフィルパック S96x20	¥19,500	¥20,000	
BH790014	スタンダードチップ 0.2-10ul 1000本	¥5,500	¥5,700	
BH790021F	フィルターチップ 20ul 96本x10R	¥19,000	¥20,000	
BH790101F	フィルターチップ 120ul 96本x10R	¥19,000	¥20,000	
BH790200	BIOHIT Tip 200ul シングルラック 96x10	¥10,000	¥10,200	
BH790201	BIOHIT Tip 200ul 96x10 ST	¥12,700	¥13,200	
BH790201F	フィルターチップ 200ul 96本x10R	¥19,000	¥20,000	
BH790202	BIOHIT Tip 200ul リフィル 96本x10	¥7,000	¥7,200	
BH790301F	フィルターチップ 300ul 96本x10R	¥19,000	¥20,000	
BH790350	BIOHIT Tip 350ul シングルラック	¥10,000	¥10,200	
BH790350CTN	スタンダードチップ 350ul 96x10Rx8	¥71,040	¥73,600	
BH790351	BIOHIT Tip 350ul シングルラック ST	¥12,700	¥13,200	
BH790351CTN	スタンダードチップ ST 350ul 96x10Rx8	¥87,840	¥90,400	
BH790352	BIOHIT Tip 350ul ツインパック 96x10	¥7,000	¥7,200	
BH790352CTN	スタンダードチップ リフィル 350ul 96x10x19	¥110,580	¥115,900	
BH790353	BIOHIT Tip 350ul リフィルパック 96x15R	¥13,800	¥14,200	
BH791000	BIOHIT Tip 1000ul シングルラック	¥10,000	¥10,200	
BH791001	BIOHIT Tip 1000ul シングルラック ST	¥12,700	¥13,200	
BH791001CTN	スタンダードチップ ST 1000ul 96x10Rx4	¥43,920	¥47,200	
BH791001F	フィルターチップ 1000ul 96本x10R	¥19,000	¥20,000	
BH791002	BIOHIT Tip 1000ul リフィルパック 96x10	¥7,500	¥7,700	
BH791002CTN	スタンダードチップ リフィル 1000ul 96x10x4	¥25,520	¥27,200	
BH791003	BIOHIT Tip 1000ul リフィルパック滅菌96x10	¥9,800	¥10,100	
BH791020	ワイドネット 1000ul (96本x10ラック)	¥11,000	¥11,200	
BH791021	ワイドネット ST 1000ul 96本x10R	¥14,300	¥14,900	
BH791200	BIOHIT Tip 1200ul シングルラック	¥12,000	¥12,600	
BH791200CTN	スタンダードチップ 1200ul 96x10Rx4	¥43,120	¥44,800	
BH791201	BIOHIT Tip 1200ul シングルラック ST	¥15,300	¥16,000	
BH791201CTN	スタンダードチップ ST 1200ul 96x10Rx4	¥53,600	¥56,800	
BH791202	BIOHIT Tip 1200ul ツインパック 96x10	¥9,100	¥9,500	
BH791202CTN	スタンダードチップ リフィル 1200ul 96x10x4	¥30,960	¥34,000	
BH791210	ロングチップ 1200ul (96本x10ラック)	¥13,000	¥13,600	
BH791211	ロングチップ ST 1200ul 96本x10R	¥16,800	¥17,500	
BH791211F	フィルターチップ 1200ul 96本x10R	¥20,000	¥22,000	
BH791212	ロングチップ 1200ul リフィル 96本x10R	¥9,800	¥10,300	
BH791213	ロングチップ ST 1200ul XT 96本x10R	¥13,800	¥14,200	
BHLH-711018	XL・MidiPlus10mlチップ用アダプターセット	¥12,000	¥12,500	
BHLH-711019	アダプターセット(Midiアース用) 標準	¥12,000	¥12,500	
BHLH-712149	交換用バッテリー(Midiアース用) 3	¥8,500	¥9,000	
BHLH-712153	シリコンアダプター	¥8,500	¥9,000	
BHLH-7129120	交換用リフィル(XL Midiアース用)未滅菌 5	¥7,100	¥7,300	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
BHLH-7129131	交換用フィルター(XL Midi7 用)滅菌済 5	¥9,500	¥9,800	
BHLH-728671	Prolineアース 3本パック 1	¥78,000	¥80,000	
BHLH-729080	Tacta 手動ピペット 1ch 500-5000ul	¥44,000	¥48,000	
BHLH-729120	Tacta 手動ピペット 8ch 0.5-10ul	¥123,000	¥130,000	
BHLH-729672	Tacta4 pack20 4本パック	¥149,000	¥157,000	
BHLH-729673	Tacta4 pack100 4本パック	¥149,000	¥157,000	
BHLH-729674	Tacta5 pack10 5本パック	¥184,000	¥192,000	
BHLH-729675	Tacta5 pack20 5本パック	¥184,000	¥192,000	
BHLH-736139	ACFアプター-Picusリ-ス 本体用 1	¥17,000	¥18,000	
BHLH-780314	スタンダードチップ ラック 10ml 35本入 1	¥4,300	¥4,400	
BHLH-780316	スタンダードチップ ハルク 10ml 250本	¥12,400	¥12,600	
BHLH780316CTN	スタンダードチップ ハルク 10ml 250x8	¥81,760	¥84,000	
BHLH-795001F	フィルターチップ 100-5000µL (50本)	¥5,900	¥6,000	
BHZ7503	BP Gelカセット 0.75% DF 50kb 3	¥76,400	¥81,500	
BHZ7510	BP Gelカセット 0.75% Dye-Free 50kb 10	¥183,800	¥195,900	
BHZ7550	BP Gelカセット 0.75% Dye-Free 50kb 50	¥798,300	¥851,100	
BIO-21040	BIOTAQ 500u dNTP付 ㊟	¥18,000	¥18,800	
BIO-21040/5	BIOTAQ 500ux5 dNTP付 ㊟	¥85,200	¥89,300	
BIO-21047	BIOTAQ-HS 500u dNTP付 ㊟	¥31,000	¥32,800	
BIO-21047/5	BIOTAQ-HS 500ux5 dNTP付㊟	¥149,500	¥155,800	
BIO-21052	Accuzyme 500u㊟	¥89,700	¥98,300	
BIO-21083	MangoTaq 1000u㊟	¥51,500	¥56,800	
BIO-21105	MyTaq DNA 赤リミターゼ 500units ㊟	¥28,800	¥31,600	
BIO-21106	MyTaq DNA 赤リミターゼ 2500units ㊟	¥129,100	¥141,500	
BIO-21107	MyTaq DNA 赤リミターゼ 5000units ㊟	¥248,700	¥272,700	
BIO-21108	MyTaq Red DNA 赤リミターゼ 500unit㊟	¥34,800	¥38,400	
BIO-21109	MyTaq Red DNA赤リミターゼ 2500unit㊟	¥129,100	¥141,500	
BIO-21110	MyTaq Red DNA赤リミターゼ 5000unit㊟	¥248,700	¥272,700	
BIO-21111	MyTaq HS DNA 赤リミターゼ 250unit㊟	¥36,700	¥40,400	
BIO-21112	MyTaq HS DNA赤リミターゼ 1000unit㊟	¥124,600	¥130,000	
BIO-21113	MyTaq HS DNA赤リミターゼ 2500unit㊟	¥300,000	¥320,000	
BIO-21115	MyTaqHS RedDNA赤リミターゼ 1000unit㊟	¥136,700	¥149,900	
BIO-21116	MyTaqHS RedDNA赤リミターゼ 2500unit㊟	¥325,200	¥356,600	
BIO-21117	MyFi DNA 赤リミターゼ 250units ㊟	¥46,200	¥50,800	
BIO-21118	MyFi DNA 赤リミターゼ 500units ㊟	¥86,300	¥94,800	
BIO-21119	MyFi DNA 赤リミターゼ 2500units ㊟	¥403,200	¥441,800	
BIO-21126	MyTaq Extract-PCR Kit 100 ㊟	¥27,600	¥30,400	
BIO-21127	MyTaq Extract-PCR Kit 500 ㊟	¥108,300	¥118,700	
BIO-25006	2xBIOMix RED 500㊟	¥91,200	¥100,300	
BIO-25012	2xBIOMix 500㊟	¥86,700	¥95,200	
BIO-25020	2xImmoMix(BIOTAQ HS) 500㊟	¥132,100	¥145,100	
BIO-25022	2xImmoMix RED 500㊟	¥153,700	¥168,700	
BIO-25028	2xAccuzyme Mix 500㊟	¥176,000	¥193,100	
BIO-25033	MangoMix 250(1.25mlx5)㊟	¥37,500	¥41,200	
BIO-25034	MangoMix 1000(1.25mlx20)㊟	¥137,000	¥150,300	
BIO-25041	MyTaq Mix 200㊟	¥36,700	¥40,400	
BIO-25042	MyTaq Mix 1000㊟	¥154,100	¥169,100	
BIO-25043	MyTaq Red Mix 200㊟	¥36,700	¥40,400	
BIO-25044	MyTaq RED Mix 1000㊟	¥154,100	¥169,100	
BIO-25045	MyTaq HS Mix 200㊟	¥55,700	¥61,200	
BIO-25046	MyTaq HS Mix 1000㊟	¥226,000	¥247,900	
BIO-25047	MyTaq HS RED Mix 200rxn㊟	¥52,200	¥57,600	
BIO-25048	MyTaq HS RED Mix 1000㊟	¥226,000	¥247,900	
BIO-25049	MyFi Mix 2x 100Reactions ㊟	¥37,500	¥41,200	
BIO-25050	MyFi Mix 2x 500Reactions ㊟	¥172,300	¥189,100	
BIO-33025	HyperLadder 1kb 200Lane㊟	¥20,400	¥22,800	
BIO-33026	HyperLadder1kb 500Lane㊟	¥45,800	¥50,400	
BIO-33029	HyperLadder100bp 200Lane㊟	¥34,100	¥37,600	
BIO-33030	HyperLadder100bp 500Lanes ㊟	¥74,200	¥81,600	
BIO-33031	HyperLadder25bp 200Lanes ㊟	¥31,800	¥35,200	
BIO-33032	HyperLadder25bp 500Lanes ㊟	¥72,700	¥80,000	
BIO-33039	HyperLadder50bp 200 Lanes ㊟	¥32,900	¥36,400	
BIO-33040	HyperLadder50bp 500 Lanes ㊟	¥73,100	¥80,400	
BIO-33045	Easy LadderI 200Lanes ㊟	¥25,700	¥28,400	
BIO-33046	Easy LadderI 500Lanes ㊟	¥56,400	¥62,000	
BIO-33053	HyperLadder 1kb 100Lanes ㊟	¥16,700	¥18,400	
BIO-33054	HyperLadder 50bp 100Lanes ㊟	¥22,300	¥24,800	
BIO-33056	HyperLadder 100bp 100Lanes ㊟	¥20,100	¥22,000	
BIO-33057	HyperLadder 25bp 100Lanes ㊟	¥22,300	¥24,800	
BIO-37025	NH4 Buffer 10x 3x1.2ml ㊟	¥11,400	¥12,800	
BIO-37026	MgCl2 Solution 50mM 3x1.2m㊟	¥11,400	¥12,800	
BIO-37035	X-GAL 1g ㊟	¥29,200	¥32,000	
BIO-37036	IPTG 5g ㊟	¥27,300	¥30,000	
BIO-37037	Proteinase K 100mg ㊟	¥31,000	¥34,400	
BIO-37045	DNA Loading Buffer Blue 5x ㊟	¥12,900	¥14,400	
BIO-37084	Proteinase K Solution 5ml ㊟	¥34,800	¥38,400	
BIO-39025	dNTP Set 100mM 4X250ul ㊟	¥29,200	¥32,000	
BIO-39026	dNTP Set 100mM 4X4X250ul ㊟	¥104,100	¥114,300	
BIO-39028	100mM dNTPMix 50umol 500ul㊟	¥20,400	¥22,800	
BIO-39029	100mM dNTPMix 200umol 4X500ul㊟	¥67,400	¥74,000	
BIO-39035	100mMdUTPindividual25umol250ul㊟	¥21,600	¥24,000	
BIO-39036	100mMdATPindividual25umol250ul㊟	¥13,300	¥14,800	
BIO-39037	100mMdGTPindividual25umol250ul㊟	¥13,300	¥14,800	
BIO-39038	100mMdCTPindividual25umol250ul㊟	¥13,300	¥14,800	
BIO-39039	100mMdTTPindividual25umol250ul㊟	¥13,300	¥14,800	
BIO-39041	dUTP Mix 25umol (1 x 500ul)㊟	¥19,700	¥21,600	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
BIO-39043	40mM dNTPMix 20umol 500ul ^②	¥13,600	¥15,200	
BIO-39044	10mM dNTPMix 10umol 1.0ml ^②	¥12,900	¥14,400	
BIO-39049	dNTP Set 100mM 4X1ml ^②	¥102,200	¥112,300	
BIO-39053	10mM dNTPMix 100umol 10x1ml ^②	¥69,700	¥76,400	
BIO-41026	Agarose Molecular Grade 100g	¥22,300	¥24,800	
BIO-41027	Agarose Tablet 300g	¥76,100	¥83,600	
BIO-65048	MyTaq One-Step RT-PCR Kit 25 ^②	¥51,100	¥56,400	
BIO-65049	MyTaq One-Step RT-PCR Kit 100 ^②	¥188,200	¥206,300	
BJ9006	加マトク フライ (固相抽出法)	¥312,700	¥350,800	
BJ9008	加マトク フライ (固相抽出法)	¥320,600	¥250,400	
BJ9054	加マトク フライ (固相抽出法)	¥341,700	¥273,300	
BJ9056	加マトク フライ (固相抽出法)	¥338,600	¥263,900	
BLF7503	BP Gelカセット 0.75% DF 1-10kb 3	¥76,400	¥81,500	
BLF7510	BP Gelカセット 0.75% Dye-Free2-10kb 10	¥183,800	¥195,900	
BLF7550	BP Gelカセット 0.75% Dye-Free2-10kb 50	¥798,300	¥851,100	
BLU0001	■ Blue Pippin 本体 (モニター付)	¥3,900,000	¥4,200,000	
BMF7503	BPカセット0.75% DF18-27kb 3	¥76,400	¥81,500	
BPLUS03	BP Gelカセット High Pass Plus 3	¥79,400	¥84,600	
BPLUS10	BP Gelカセット High Pass Plus 10	¥213,200	¥227,300	
BPLUS50	BP Gelカセット High Pass Plus 50	¥873,300	¥931,100	
BR0501402	4X プローブ マスターミックスキット(No Rox) ^②	¥64,000	¥73,000	
BR0501502	4X プローブ マスターミックスキット(LRox) ^②	¥64,000	¥73,000	
BR0501602	4X プローブ マスターミックスキット(HRox) ^②	¥64,000	¥73,000	
BR0501701	4X SYBR Green マスターミックスキット(No Rox) ^②	¥12,900	¥14,000	
BR0501801	4X SYBR Green マスターミックスキット(LRox) ^②	¥12,900	¥14,000	
BR0501901	4X SYBR Green マスターミックスキット(HRox) ^②	¥12,900	¥14,000	
BSP0161-A	フルオート入Wメンブリ 0.2μ 26cmX3.3m 〇-	¥43,400	¥46,400	2024/4/1追加
BUF7503	BP Gelカセット0.75% Dye-Free10-18kb 3	¥76,400	¥81,500	
BUF7510	BP Gelカセット0.75% Dye-Free10-18kb 10	¥183,800	¥195,900	
BUF7550	BP Gelカセット0.75% Dye-Free10-18kb 50	¥798,300	¥851,100	
C1	G100 プローブ 0.2mlx96ウェル	¥19,600	¥20,400	
C115-01	ClonExpressカセットワンステップ クロニング キット ^②	¥30,000	¥46,600	
C115-02	ClonExpressカセットワンステップ クロニング キット ^②	¥56,000	¥75,900	
C2	G100 プローブ 0.5mlx58ウェル	¥19,600	¥20,400	
C3	G100 プローブ 1.5mlx39ウェル	¥19,600	¥20,400	
C4	G100 プローブ 2.0mlx39ウェル	¥19,600	¥20,400	
C5	G100 プローブ 5.0mlx18ウェル	¥19,600	¥20,400	
CDF1503	Gelカセット 1.5% DF 250bp-1.5kb 3	¥61,700	¥65,800	
CDF1510	Gelカセット 1.5% DF 250bp-1.5kb 10	¥166,100	¥177,100	
CDF1550	Gelカセット 1.5% DF 250bp-1.5kb 50	¥707,500	¥754,300	
CDF2003	Gelカセット 2% DF 100-600bp 3	¥61,700	¥65,800	
CDF2010	Gelカセット 2% DF 100-600bp 10	¥166,100	¥177,100	
CDF2050	Gelカセット 2% DF 100-600bp 50	¥707,500	¥754,300	
CDH1504	HT コントロールDNA Gelカセット1.5%用 ④	¥27,600	¥29,400	
CDH2004	HT コントロールDNA Gelカセット2%用 ④	¥27,600	¥29,400	
CDH3004	HT コントロールDNA Gelカセット3%用 ④	¥27,600	¥29,400	
CDH7504	HT コントロールDNA Gelカセット0.75%用 ④	¥55,100	¥58,800	
CDP3003	Gelカセット 3% DF 100-250bp 3	¥61,700	¥65,800	
CDP3010	Gelカセット 3% DF 100-250bp 10	¥166,100	¥177,100	
CEF2003	Gelカセット 2% EF 100-600bp 3	¥61,700	¥65,800	
CEF2010	Gelカセット2%EtBr-FreeFor100-600bp 10	¥166,100	¥177,100	
CEF2050	Gelカセット2%EtBr-FreeFor100-600bp 50	¥707,500	¥754,300	
CELECO1	HLS E.coliCellSuspension kit	¥14,700	¥15,700	
CELMWB1	HLS Cell suspension kit	¥14,700	¥15,700	
CHP7504	0.75%PacBio ValidationDNA ④	¥30,100	¥32,100	
CIS1504	コントロールDNA1.5%用internal standard④	¥30,100	¥32,100	
CIS2004	コントロールDNA2%用internal standard④	¥30,100	¥32,100	
CIS3004	コントロールDNA3%用internal standard④	¥30,100	¥32,100	
CON1504	コントロールDNA Gelカセット1.5%用 ④	¥30,100	¥32,100	
CON2004	コントロールDNA Gelカセット2%用 ④	¥30,100	¥32,100	
CON3004	コントロールDNA Gelカセット3%用 ④	¥30,100	¥32,100	
CON7504	コントロールDNA Gelカセット0.75%用 ④	¥30,100	¥32,100	
CON7604	コントロールDNA Gelカセット0.75%用	¥30,100	¥32,100	
CONTROL01DNAE	Control DNA Template exo ^②	¥19,300	¥20,600	
CONTROL01EXO	Control Primer & Probe Mix exo ^②	¥19,300	¥20,600	
CS-04-001	パソパソカ 20mlx5本 ④	¥13,000	¥15,000	
CS-06-001	パソパソカ タルク 20ml ④	¥5,000	¥9,000	
CS-07-001	パソパソカ-hRM 20ml ④	¥8,000	¥9,000	
CS-09-001	パソパソカ-DMSO7リ 20ml ④	¥7,000	¥9,000	
CSM0026	フリカ用紙(YJ-360T-S用) 20	¥32,200	¥34,900	
CV5550	不織布タレ	¥146,400	¥159,800	
D101-01	MycoBlueマイコブルー マスターミックス 20回用 ^②	¥14,000	¥20,200	
D101-02	MycoBlueマイコブルー マスターミックス 50回用 ^②	¥31,300	¥45,200	
D15	GA150 プローブ 1.5mlx24ウェル	¥19,600	¥20,400	
D16	GA150 プローブ 2.0mlx24ウェル	¥19,600	¥20,400	
D17	GA150 プローブ 0.2mlx96ウェル	¥22,500	¥23,400	
D6	GA150 プローブ ウェル径15mmx12ウェル	¥19,600	¥20,400	
DC-8CH	マイクロ8/1chチップホルダー - プラチナチップ用	¥16,900	¥19,200	
DE111-01	RNase A (100 mg/ml) 500ul ^②	¥14,900	¥21,600	
EK1725852	Electron Microscope	¥419,800	¥460,500	
ELD2004	ELF DNA Gelカセット 2% w/リカ 4	¥42,300	¥45,100	
ELD2010	ELF DNA Gelカセット 2% w/リカ 10	¥84,900	¥90,500	
ELD4004	ELF DNA Gelカセット 0.75% w/リカ 4	¥42,300	¥45,100	
ELD4010	ELF DNA Gelカセット 0.75% w/リカ 10	¥84,900	¥90,500	
ELD7504	ELF DNA Gelカセット 0.75% 4	¥42,300	¥45,100	
ELD7510	ELF DNA Gelカセット 0.75% 10	¥84,900	¥90,500	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
ELF0001	■ ELF 本体	¥5,200,000	¥5,500,000	
ELP3004	ELF SDS Protein Gelカセット 3% 4	¥64,700	¥69,000	
ELP3010	ELF SDS Protein Gelカセット 3% 10	¥161,700	¥172,400	
ELP5004	ELF SDS Protein Gelカセット 5% 4	¥64,700	¥69,000	
ELP5010	ELF SDS Protein Gelカセット 5% 10	¥161,700	¥172,400	
EN402-01	DNase I RNase-free(50U/ul)1000U ²⁰	¥18,600	¥26,800	
EP13920C	ガラスバイアル 20ml	¥157,700	¥179,700	
EP13920CCT	ガラスバイアル	¥191,700	¥218,800	
EP14040C	ガラスバイアル	¥169,900	¥193,700	
EP14040CCT	ガラスバイアル	¥212,300	¥241,900	
EP14040CDB	ガラスバイアル	¥240,400	¥274,500	
EP14040CEP	ガラスバイアル	¥176,600	¥201,300	
EP14040CTS	ガラスバイアル	¥215,800	¥246,100	
EP14060C	ガラスバイアル	¥410,400	¥475,400	
EP14140A	ガラスバイアル	¥170,600	¥194,400	
EP14140ACT	ガラスバイアル	¥206,700	¥235,500	
EP14140ADB	ガラスバイアル	¥218,300	¥249,200	
EP14140AEP	ガラスバイアル	¥179,200	¥204,300	
EP14140ATS	ガラスバイアル	¥219,600	¥250,300	
EP33920A	ガラスバイアル	¥73,500	¥83,700	
EP33920C	ガラスバイアル	¥73,500	¥83,800	
EP34040C	ガラスバイアル	¥90,000	¥102,500	
EP34040CDB	ガラスバイアル	¥177,200	¥201,900	
EP34060C	ガラスバイアル	¥282,400	¥326,500	
EP34140A	ガラスバイアル	¥90,000	¥102,500	
EP34140ADB	ガラスバイアル	¥177,100	¥201,900	
EP34160A	ガラスバイアル	¥252,800	¥292,400	
EPC2002APF	広口褐色ガラスボトル 60mL 24個	¥162,600	¥184,600	
EPC2004APF	広口褐色ガラスボトル 125mL 24個	¥159,600	¥181,400	
EPSCASV22	バイアル瓶(テフロンガラス製)	¥165,600	¥186,500	
EPST10020	Depyrogenated Vial 100ml ST 50本	¥362,000	¥411,100	
EPST1020	Depyrogenated Vial 10ml ST 50本	¥140,000	¥159,000	
EPST2020	Depyrogenated Vial 20ml ST 50本	¥333,400	¥378,500	
EPST213	Depyrogenated Vial 2ml ST 100本	¥214,500	¥243,900	
EPST3020	Depyrogenated Vial 30ml ST 50本	¥357,800	¥406,300	
EPST5020	Vial 50ml ST50本	¥533,600	¥605,700	
EPST520	Depyrogenated Vial 5ml ST 50本	¥136,800	¥155,400	
EPSTT13100S	バイアル瓶(テフロンガラス製)	¥421,800	¥478,300	
EPSTT16100S	バイアル瓶(テフロンガラス製)	¥499,800	¥567,600	
EVA13305000	ポンプアセサリー(オイル蒸気逆流防止トラップ)	¥389,200	¥708,000	
EVA1638	ポンプアセサリー(真空計アダプター)	¥274,500	¥499,000	
EVA1641	ポンプアセサリー(真空計アダプター)	¥246,600	¥449,000	
EVA46220000	ポンプアセサリー(オイルミスト除去フィルター)	¥279,600	¥508,000	
EVA46226000	ポンプアセサリー(オイルミスト除去フィルター)	¥383,700	¥740,000	
EVA46229000	ポンプアセサリー(オイルミスト除去フィルター)	¥557,500	¥1,014,000	
EVC10511311	ポンプアセサリー(バキューム・コネクション)	¥24,600	¥45,000	
EVC10511366	ポンプアセサリー(バキューム・コネクション)	¥9,600	¥17,000	
EVC10511396	ポンプアセサリー(O-リング類)	¥6,500	¥12,000	
EVC10511410	ポンプアセサリー(エルボ)	¥37,200	¥68,000	
EVC10511645	ポンプアセサリー(ホース口)	¥22,700	¥41,000	
EVC10512401	ポンプアセサリー(O-リング類)	¥6,700	¥12,000	
EVC10514311	ポンプアセサリー(バキューム・コネクション)	¥22,700	¥41,000	
EVC10514366	ポンプアセサリー(バキューム・コネクション)	¥11,000	¥20,000	
EVC10514401	ポンプアセサリー(O-リング類)	¥15,600	¥28,000	
EVC10514410	ポンプアセサリー(エルボ)	¥46,100	¥84,000	
EVC10514645	ポンプアセサリー(ホース口)	¥36,800	¥67,000	
EVC10516366	ポンプアセサリー(バキューム・コネクション)	¥17,800	¥32,000	
EVC10516401	ポンプアセサリー(O-リング類)	¥11,200	¥20,000	
EVC10516410	ポンプアセサリー(エルボ)	¥73,900	¥134,000	
EVC10516611	ポンプアセサリー(バキューム・コネクション)	¥29,400	¥53,000	
EVC10520050	ポンプアセサリー(バキューム・コネクション)	¥55,400	¥101,000	
EVC31105000	ポンプアセサリー(真空弁)	¥263,300	¥479,000	
EVD35610000	ポンプアセサリー(真空計)	¥672,200	¥1,223,000	
EVD40001010	ポンプアセサリー(ケーブル)	¥29,500	¥54,000	
EVD40001050	ポンプアセサリー(ケーブル)	¥73,200	¥133,000	
EVD40001100	ポンプアセサリー(ケーブル)	¥104,800	¥191,000	
EVH02124013	ポンプアセサリー(O-リング類)	¥7,700	¥14,000	
EVH02124032	ポンプアセサリー(O-リング類)	¥19,500	¥36,000	
EVH02124035	ポンプアセサリー(O-リング類)	¥22,700	¥41,000	
EVH11025015	Oil Vacuum Ultragrade19 1L	¥29,000	¥53,000	
FG-0008DC	PCR8連トームチップ [®] 小分梱包12 x 10	¥6,800	¥7,500	
FG-0008FC	PCR8連フラットチップ [®] 小分梱包12 x 10	¥6,800	¥7,500	
FG-0018WF/SE	0.1ml PCR8連トームチップ [®] + フラットチップ [®] 一体型 整列 12x10	¥14,400	¥14,900	
FG-0088WD/SE	0.2ml PCR8連トームチップ [®] + トームチップ [®] 一体型 整列 12x10	¥14,400	¥14,900	
FG-0088WF/SE	0.2ml PCR8連トームチップ [®] + フラットチップ [®] 一体型 整列 12x10	¥14,400	¥14,900	
FG-008DC	PCR8連トームチップ [®] 120	¥4,500	¥6,000	
FG-008FC	PCR8連フラットチップ [®] 120	¥5,500	¥6,000	
FG-008FCP	PCR8連フラットチップ [®] PLUS 120	¥5,500	¥6,000	
FG-011F	0.1ml PCRシグナルチップ [®] Flat 1000	¥4,500	¥4,700	
FG-012FC	PCR 12連フラットチップ [®] 8ストリップ x 25袋	¥12,600	¥13,000	
FG-018WF	0.1ml PCR8連トームチップ [®] + フラットチップ [®] 一体型 120	¥11,300	¥11,700	
FG-021D	0.2ml PCRシグナルチップ [®] Dome 1000	¥4,500	¥4,700	
FG-021D/SP	0.2ml PCRシグナルチップ [®] Dome小分 100x10	¥6,700	¥6,900	
FG-021F	0.2ml PCRシグナルチップ [®] Flat 1000	¥4,500	¥4,700	
FG-021F/SP	0.2ml PCRシグナルチップ [®] Flat小分 100x10	¥6,700	¥6,900	
FG-028DC	0.2ml PCR8連トームチップ [®] + トームチップ [®] 120	¥9,000	¥10,000	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
FG-028FC	0.2ml PCR8連1-7°+フラットチップ° 120	¥10,200	¥10,000	
FG-0384MH	384well PCR Plate用マッドハート° 10	¥7,800	¥8,000	
FG-0384P	384well PCR/QPCR°プレート 50	¥29,800	¥30,700	
FG-07	白色プレート(FAS-DigiPRO/FAS5専用)	¥33,000	¥41,500	
FG-088WD	0.2ml PCR8連1-7°+トームチップ°一体型 120	¥11,300	¥11,700	
FG-088WF	0.2ml PCR8連1-7°+フラットチップ°一体型 120	¥11,300	¥11,700	
FG-102	10ulチップ° NT 1000本	¥2,500	¥2,600	
FG-102R	10ulチップ° P.Pチップ入 NT 96本x10パック	¥8,100	¥8,300	
FG-102RS	10ulチップ° P.Pチップ入 NT ST 96本x10パック	¥10,200	¥10,500	
FG-103	10ulチップ° NT 1000本	¥2,500	¥2,600	
FG-103R	10ulチップ° PPチップ入 NT 96本x10R	¥8,100	¥8,300	
FG-103RS	10ulチップ° PPチップ入 ST 96本x10R	¥10,200	¥10,500	
FG-1700	96wellPCR°プレートUltraEasyCut ST5枚x5	¥13,600	¥14,000	
FG-1702	96well PCR°プレート UltraEasyCut 5枚x5	¥10,700	¥11,000	
FG-1740	96wellPCR°プレートノズカト ST 5枚x5	¥18,900	¥19,500	
FG-1742	96well PCR°プレートノズカト 5枚x5	¥15,900	¥16,400	
FG-2040	10ulチップ° 滅菌済 96本 x 50パック	¥46,200	¥47,600	
FG-2040-1	10ulチップ° 滅菌済 96本 x 10パック	¥12,600	¥13,000	
FG-2074	10ulチップ° 1000本	¥3,500	¥3,600	
FG-301	200ulチップ° 1000本	¥2,500	¥2,600	
FG-301R	200ulチップ° P.Pチップ入 96本x10パック	¥8,100	¥8,300	
FG-301RS	200ulチップ° P.Pチップ入 ST 96本x10パック	¥10,200	¥10,500	
FG-302	200ulチップ° 1000本	¥2,500	¥2,600	
FG-302R	200ulチップ° P.Pチップ入 96本x10パック	¥9,200	¥9,400	
FG-302RS	200ulチップ° P.Pチップ入ST 96本x10パック	¥11,300	¥11,600	
FG-402	1000ul FG°チップ° 1000/bag	¥4,400	¥4,700	
FG-402R	1000ul FG°チップ° 100X10R	¥9,000	¥9,700	
FG-402RS	1000ul FG°チップ° S100X10R	¥10,400	¥11,100	
FG-AatII	FG Restriction Enzyme Aat II [®]	¥8,200	¥10,000	
FG-AccI	FG Restriction Enzyme Acc I [®]	¥8,800	¥10,000	
FG-AccIII	FG Restriction Enzyme Acc III [®]	¥8,800	¥10,000	
FG-AcuI	FG Restriction Enzyme Acu I [®]	¥8,200	¥10,000	
FG-AflII	FG Restriction Enzyme Afl II [®]	¥8,800	¥10,000	
FG-AgeI	FG Restriction Enzyme Age I [®]	¥8,200	¥15,000	
FG-Alw26I	FG Restriction Enzyme Alw26 I [®]	¥11,000	¥13,000	
FG-AlwI	FG Restriction Enzyme AlwI [®]	¥9,900	¥11,000	
FG-ApaI	FG Restriction Enzyme Apa [®]	¥5,800	¥7,000	
FG-ApaLI	FG Restriction Enzyme ApaL [®]	¥8,200	¥10,000	
FG-ApoI	FG Restriction Enzyme ApoI [®]	¥9,900	¥11,000	
FG-AscI	FG Restriction Enzyme Asc [®]	¥8,200	¥10,000	
FG-AvaI	FG Restriction Enzyme Ava I [®]	¥8,200	¥10,000	
FG-AvaII	FG Restriction Enzyme Ava II [®]	¥8,200	¥10,000	
FG-AvrII	FG Restriction Enzyme Avr II [®]	¥8,200	¥10,000	
FG-BalI	FG Restriction Enzyme Bal [®]	¥19,400	¥23,000	
FG-BamHI	FG Restriction Enzyme BamH I [®]	¥4,100	¥5,000	
FG-BclI	FG Restriction Enzyme Bcl I [®]	¥6,800	¥8,000	
FG-BglI	FG Restriction Enzyme Bgl I [®]	¥6,800	¥8,000	
FG-BglII	FG Restriction Enzyme Bgl II [®]	¥4,800	¥6,000	
FG-BsaI	FG Restriction Enzyme Bsa I [®]	¥8,200	¥10,000	
FG-BsaWI	FG Restriction Enzyme BsaW [®]	¥8,200	¥10,000	
FG-BsiWI	FG Restriction Enzyme BsiW I [®]	¥8,200	¥10,000	
FG-BsmBI	FG Restriction Enzyme BsmB I [®]	¥8,800	¥10,000	
FG-BsoBI	FG Restriction Enzyme BsoB I [®]	¥8,200	¥10,000	
FG-BspEI	FG Restriction Enzyme BspE I [®]	¥8,200	¥10,000	
FG-BsrFI	FG Restriction Enzyme BsrF I [®]	¥8,200	¥10,000	
FG-BstYI	FG Restriction Enzyme BstY I [®]	¥8,200	¥10,000	
FG-BtsCI	FG Restriction Enzyme BtsC I [®]	¥6,100	¥7,000	
FG-CDC01	FastGene Cap easy	¥17,800	¥18,400	
FG-CDC02	CapEasy v2 (0.1mL/0.2mL兼用)	¥21,000	¥21,600	
FG-Cfr10I	FG Restriction Enzyme Cfr10 I [®]	¥8,700	¥10,000	
FG-Cfr42I	FG Restriction Enzyme Cfr42 I [®]	¥6,400	¥7,000	
FG-Cfr9I	FG Restriction Enzyme Cfr9 I [®]	¥8,200	¥10,000	
FG-CH01	FastGene Western ECL Kit 50mL ^④	¥12,100	¥14,000	
FG-CH05	FastGene Block & Go 500mL	¥21,100	¥24,000	
FG-ClaI	FG Restriction Enzyme Cla I [®]	¥6,100	¥7,000	
FG-CRY-05RC	SBSチップ FG-CRY0.5mL用 10	¥11,500	¥11,800	
FG-CRY-10RC	SBSチップ FG-CRY1.0mL用 10	¥11,500	¥11,800	
FG-CRY-10S	1.0mL クライオチップ° アクター ST 500	¥32,500	¥33,500	
FG-CRY-20RC	SBSチップ FG-CRY2.0mL用 10	¥11,500	¥11,900	
FG-CRY-20S	2.0mL クライオチップ° アクター ST 500	¥32,500	¥33,500	
FG-CRY-2D	2Dインサート FG-CRY用 500	¥23,100	¥23,800	
FG-CRY-In-10S	1.0mL クライオチップ° インナー ST 500	¥32,500	¥33,500	
FG-CRY-In-20S	2.0mL クライオチップ° インナー ST 500	¥32,500	¥33,500	
FG-CviAI	FG Restriction Enzyme CviAI [®]	¥10,000	¥11,000	
FG-DdeI	FG Restriction Enzyme Dde I [®]	¥8,200	¥10,000	
FG-DpnI	FG Restriction Enzyme Dpn I [®]	¥6,800	¥8,000	
FG-DpnII	FG Restriction Enzyme Dpn II [®]	¥7,400	¥9,000	
FG-DraI	FG Restriction Enzyme Dra I [®]	¥6,800	¥8,000	
FG-EagI	FG Restriction Enzyme Eag [®]	¥7,600	¥9,000	
FG-Eco47I	FG Restriction Enzyme Eco47 I [®]	¥5,100	¥6,000	
FG-EcoNI	FG Restriction Enzyme EcoNI [®]	¥9,200	¥11,000	
FG-EcoO109I	FG Restriction Enzyme EcoO109 I [®]	¥8,200	¥10,000	
FG-EcoRI	FG Restriction Enzyme EcoR I [®]	¥8,200	¥10,000	
FG-EcoRV	FG Restriction Enzyme EcoR V [®]	¥6,800	¥8,000	
FG-EcoT38I	FG Restriction Enzyme EcoT38 I [®]	¥8,200	¥10,000	
FG-EGL010	FG Taq Mix with Rox 5x1ml [®]	¥35,000	¥41,000	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
FG-Esp3I	FG Restriction Enzyme Esp3 I ②	¥8,800	¥10,000	
FGF-10002	1000ulロングフィルターチップシステムキットST96x10	¥19,900	¥20,500	
FGF-10003	1000ul リフィル フィルターチップ ST 96x10R	¥17,300	¥17,800	
FGF-2042	10ulフィルターチップ(ラックPP)滅菌済 96x10R	¥16,100	¥16,600	
FGF-2043	10ul ショートリフィル フィルターチップ ST 96x10R	¥14,700	¥15,200	
FGF-2072	10ulロングフィルターチップシステムキット(ST)96x10	¥16,100	¥16,600	
FGF-2073	10ul ショートリフィル フィルターチップ ST 96x10R	¥14,900	¥15,200	
FGF-20S2	20ul フィルターチップシステムキット(PP) ST96x10	¥14,900	¥15,400	
FGF-20S3	20ul リフィル フィルターチップ ST 96x10R	¥13,600	¥14,000	
FGF-7032	200ulフィルターチップシステムキット(PP) ST96x10	¥14,900	¥15,400	
FGF-7033	200ul リフィル フィルターチップ ST 96x10R	¥13,600	¥14,000	
FG-FokI	FG Restriction Enzyme Fok I ②	¥8,200	¥10,000	
FG-FP-01	FastGene フェルトバック(アルコール耐性) 1	¥1,500	¥1,600	
FG-FspI	FG Restriction Enzyme Fsp ②	¥8,200	¥10,000	
FG-GDS4	DNA Saveプレート-2.0 4サクル100入	¥124,000	¥163,000	
FG-GTD02-100	FastGene DNA Save チューブ 100	¥149,000	¥128,000	
FG-GTD02-3	FastGene DNA Save チューブ トライアル 3	¥4,500	¥5,000	
FG-GTR05-100	FastGene RNA Save チューブ 100	¥165,000	¥146,000	
FG-GTR05-3	FastGene RNA Save チューブ トライアル 3	¥5,500	¥5,000	
FG-HaeII	FG Restriction Enzyme Hae II ②	¥6,800	¥8,000	
FG-HaeIII	FG Restriction Enzyme Hae III ②	¥8,200	¥10,000	
FG-HgaI	FG Restriction Enzyme Hga I ②	¥8,200	¥10,000	
FG-HincII	FG Restriction Enzyme Hinc II ②	¥6,800	¥8,000	
FG-HindII	FG Restriction Enzyme Hind II ②	¥6,800	¥8,000	
FG-HindIII	FG Restriction Enzyme Hind III ②	¥8,200	¥10,000	
FG-HinfI	FG Restriction Enzyme Hinf ②	¥6,800	¥8,000	
FG-HinPII	FG Restriction Enzyme HinPI ②	¥8,200	¥10,000	
FG-HpaI	FG Restriction Enzyme Hpa ②	¥8,200	¥10,000	
FG-HpaII	FG Restriction Enzyme Hpa II ②	¥7,600	¥9,000	
FG-HphI	FG Restriction Enzyme Hph I ②	¥7,600	¥9,000	
FG-Hpy188I	FG Restriction Enzyme Hpy188 I ②	¥8,200	¥10,000	
FG-Hpy99I	FG Restriction Enzyme Hpy99 I ②	¥8,200	¥10,000	
FG-HpyCH4V	FG Restriction Enzyme HpyCH4 V ②	¥8,200	¥10,000	
FG-Kpn2I	FG Restriction Enzyme Kpn2 ②	¥7,600	¥9,000	
FG-KpnI	FG Restriction Enzyme Kpn ②	¥4,600	¥5,000	
FG-LK30	FG Kickspeed DNA Ligation kit ②	¥8,100	¥9,000	
FG-LK60	FG Kickspeed DNA Ligation kit ②	¥14,500	¥17,000	
FG-LM50	FG Kickspeed 2X DNA Ligation Mix ②	¥13,700	¥16,000	
FG-Lsp1109I	FG Restriction Enzyme Lsp1109 I ②	¥5,500	¥6,000	
FG-MboI	FG Restriction Enzyme Mbo ②	¥9,200	¥11,000	
FG-MboII	FG Restriction Enzyme Mbo II ②	¥8,200	¥10,000	
FG-MCT0155C	FG 1.5mL マイクロチューブ(リフトアップ) 1000	¥4,700	¥4,900	
FG-MCT015C	FastGene 1.5mL マイクロチューブ 1000	¥5,200	¥5,400	
FG-MCT015CLB	FG 1.5mL 核酸低吸着マイクロチューブ 100	¥3,000	¥3,200	
FG-MCT020C	FastGene 2.0mL マイクロチューブ 500	¥4,500	¥4,700	
FG-MCT020CLB	FG 2.0mL 核酸低吸着マイクロチューブ 100	¥3,000	¥3,200	
FG-MluI	FG Restriction Enzyme Mlu I ②	¥6,800	¥8,000	
FG-MM003PG	FG 2X PG ProbeMix w/enhancer ②	¥55,700	¥66,000	
FG-MnlI	FG Restriction Enzyme Mnl ②	¥8,800	¥10,000	
FG-MseI	FG Restriction Enzyme Mse ②	¥8,200	¥10,000	
FG-MspAII	FG Restriction Enzyme MspA I ②	¥8,800	¥10,000	
FG-MspI	FG Restriction Enzyme Msp ②	¥6,800	¥8,000	
FG-MunI	FG Restriction Enzyme Mun I ②	¥4,600	¥5,000	
FG-NaeI	FG Restriction Enzyme Nae ②	¥6,800	¥8,000	
FG-NcoI	FG Restriction Enzyme Nco I ②	¥6,100	¥7,000	
FG-NdeI	FG Restriction Enzyme Nde I ②	¥6,800	¥8,000	
FG-NgoMIV	FG Restriction Enzyme NgoM IV ②	¥8,200	¥10,000	
FG-NheI	FG Restriction Enzyme Nhe I ②	¥6,800	¥8,000	
FG-NlaIV	FG Restriction Enzyme Nla IV ②	¥8,200	¥10,000	
FG-NotI	FG Restriction Enzyme Not ②	¥6,100	¥7,000	
FG-NP01	■ FG NanoSpec Photometer	¥1,750,000	¥1,890,000	
FG-NPSC01	NanoSpec用 Calibration solution	¥20,100	¥21,800	
FG-NPSC02	NanoView用 Calibration solution	¥20,100	¥21,800	
FG-NruI	FG Restriction Enzyme Nru I ②	¥6,800	¥8,000	
FG-NtBstNBI	FG Restriction Enzyme Nt.BstNB I ②	¥6,800	¥8,000	
FG-OP008FC	FG 0.2mL 8連PCRチップ(フロッグ)	¥5,500	¥6,000	
FG-OP028FC	0.2ml PCR8連チューブ(frosted)+フロッグチップ 120	¥10,200	¥10,000	
FG-PaeR7I	FG Restriction Enzyme PaeR7 I ②	¥6,800	¥8,000	
FG-PflMI	FG Restriction Enzyme PflM I ②	¥6,800	¥8,000	
FG-PleI	FG Restriction Enzyme Ple I ②	¥8,200	¥10,000	
FG-PluTI	FG Restriction Enzyme PluT ②	¥8,200	¥10,000	
FG-PspGI	FG Restriction Enzyme PspG I ②	¥6,300	¥7,000	
FG-PstI	FG Restriction Enzyme Pst I ②	¥4,700	¥6,000	
FG-PvuI	FG Restriction Enzyme Pvu ②	¥8,200	¥10,000	
FG-PvuII	FG Restriction Enzyme Pvu II ②	¥6,800	¥8,000	
FG-REB1	10X FastGene Buffer I 1ml ②	¥1,400	¥2,000	
FG-REB2	10X FastGene Buffer II 1ml ②	¥1,400	¥2,000	
FG-REB3	10X FastGene Buffer III 1ml ②	¥1,400	¥2,000	
FG-REB4	10X FastGene Buffer IV 1ml ②	¥1,400	¥2,000	
FG-REBAccIII	10X FastGene Buffer Acc III 1ml ②	¥1,400	¥2,000	
FG-REBBalI	10X FastGene Buffer Bal I 1ml ②	¥1,400	¥2,000	
FG-REBBamHI	10X FastGene Buffer BamH I 1ml ②	¥1,400	¥2,000	
FG-REBDpnII	10X FastGene Buffer Dpn II 1ml ②	¥1,400	¥2,000	
FG-REBEcoRI	10X FastGene Buffer EcoR I 1ml ②	¥1,400	¥2,000	
FG-REBHF	10X FastCut Buffer 1 ml ②	¥1,800	¥2,000	
FG-REDiA	FastGene Diluent A 1ml ②	¥1,400	¥2,000	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
FG-REDiB	FastGene Diluent B 1ml ^②	¥1,400	¥2,000	
FG-REDiC	FastGene Diluent C 1ml ^②	¥1,400	¥2,000	
FG-RsaI	FG Restriction Enzyme Rsa I ^②	¥6,800	¥8,000	
FG-SA5600	凍結用 ^① 用 ^② ハル 1.5-2ml用 40x12mm	¥27,100	¥28,100	
FG-SacI	FG Restriction Enzyme Sac I ^②	¥6,800	¥8,000	
FG-SacII	FG Restriction Enzyme Sac II ^②	¥6,800	¥8,000	
FG-SalI	FG Restriction Enzyme Sal ^②	¥10,200	¥12,000	
FG-Sau96I	FG Restriction Enzyme Sau96 I ^②	¥6,100	¥7,000	
FG-SbfI	FG Restriction Enzyme Sbf ^②	¥8,200	¥10,000	
FG-ScaI	FG Restriction Enzyme Sca I ^②	¥8,200	¥10,000	
FG-SCR-15LB	FG 1.5mL スクリューキャップ ^① DLB ^② 100	¥7,800	¥8,200	
FG-SCR-20LB	FG 2.0mL スクリューキャップ ^① DLB ^② 100	¥7,800	¥8,200	
FG-SCRB-OL05S	0.5mL オリジナル無スクリューキャップ ^① チューブ ^② S 500	¥26,700	¥27,500	
FG-SCRG-OL20S	1.8mL オリジナル無スクリューキャップ ^① チューブ ^② S 500	¥26,700	¥27,500	
FG-SCRLP-05S	FG 0.5mL スクリューキャップ ^① チューブ ^② 500	¥29,700	¥31,200	
FG-SCRLP-20S	FG 2.0mL スクリューキャップ ^① チューブ ^② 500	¥29,700	¥31,200	
FG-SCT05-BS	FG 0.5mL凍結保存チューブ ^① 青(キャップ付)	¥23,000	¥23,800	
FG-SCT05-GS	FG 0.5mL凍結保存チューブ ^① 緑(キャップ付)	¥23,000	¥23,800	
FG-SCT05-RS	FG 0.5mL凍結保存チューブ ^① 赤(キャップ付)	¥23,000	¥23,800	
FG-SCT05-S	FG 0.5mL凍結保存チューブ ^① NT(キャップ付)	¥23,000	¥23,800	
FG-SCT05-SHS	FG 0.5mL凍結保存チューブ ^① 遮光(キャップ付)	¥25,000	¥25,800	
FG-SCT05-YS	FG 0.5mL凍結保存チューブ ^① 黄(キャップ付)	¥23,000	¥23,800	
FG-SCT20-BS	FG 2.0mL凍結保存チューブ ^① 青(キャップ付)	¥23,000	¥23,800	
FG-SCT20-GS	FG 2.0mL凍結保存チューブ ^① 緑(キャップ付)	¥23,000	¥23,800	
FG-SCT20-RS	FG 2.0mL凍結保存チューブ ^① 赤(キャップ付)	¥23,000	¥23,800	
FG-SCT20-S	FG 2.0mL凍結保存チューブ ^① NT(キャップ付)	¥23,000	¥23,800	
FG-SCT20-SHS	FG 2.0mL凍結保存チューブ ^① 遮光(キャップ付)	¥25,000	¥25,800	
FG-SCT20-YS	FG 2.0mL凍結保存チューブ ^① 黄(キャップ付)	¥23,000	¥23,800	
FG-SdaI	FG Restriction Enzyme Sda I ^②	¥6,100	¥7,000	
FG-SgrAI	FG Restriction Enzyme SgrA I ^②	¥7,800	¥9,000	
FG-SmaI	FG Restriction Enzyme Sma I ^②	¥6,800	¥8,000	
FG-SnaBI	FG Restriction Enzyme SnaB ^②	¥6,100	¥7,000	
FG-SpeI	FG Restriction Enzyme Spe ^②	¥6,800	¥8,000	
FG-SphI	FG Restriction Enzyme Sph I ^②	¥5,800	¥7,000	
FG-Sse9I	FG Restriction Enzyme Sse9 ^②	¥8,200	¥10,000	
FG-SspI	FG Restriction Enzyme Ssp I ^②	¥6,800	¥8,000	
FG-StuI	FG Restriction Enzyme Stu I ^②	¥6,100	¥7,000	
FG-StyD4I	FG Restriction Enzyme StyD4 I ^②	¥8,200	¥10,000	
FG-SwaI	FG Restriction Enzyme Swa I ^②	¥6,100	¥7,000	
FG-SWv2-AB01	FG SeeWell 96穴マルチ ^① 0.1mL用	¥26,000	¥29,000	
FG-SWv2-AB02	FG SeeWell 96穴マルチ ^① 0.2mL用	¥26,000	¥29,000	
FG-T4	FastGene T4 DNA Ligase(20,000 U) ^②	¥11,300	¥13,000	
FG-T4BP	FastGene T4 DNA Ligase(100,000U) ^②	¥45,100	¥53,000	
FG-T4HC	FastGene T4 DNA Ligase(100,000U) ^②	¥45,100	¥53,000	
FG-TaqI	FG Restriction Enzyme Taq I ^②	¥6,800	¥8,000	
FG-TspMI	FG Restriction Enzyme TspM I ^②	¥8,200	¥10,000	
FG-Tth111I	FG Restriction Enzyme Tth111 I ^②	¥6,100	¥7,000	
FG-XbaI	FG Restriction Enzyme Xba I ^②	¥6,100	¥7,000	
FG-XhoI	FG Restriction Enzyme Xho ^②	¥4,700	¥6,000	
FG-XmaI	FG Restriction Enzyme Xma ^②	¥6,800	¥8,000	
G100	G100 トライバ ^① スインキ ^② ター本体	¥48,700	¥50,600	
GA150-1	■ GA150 トライバ ^① スインキ ^② ター本体	¥69,100	¥71,900	
GA150-1A	■ GA150 トライバ ^① スインキ ^② ター本体(浅型)	¥71,600	¥74,400	
GE-96G	■ GeneExplorer サーマライザー 96x0.2mL	¥560,000	¥620,000	
GM60207	バ ^① イ ^② ク ^③ イ ^④ B ^⑤ ク ^⑥ ア ^⑦ ル ^⑧ BNBZF3RT 0.45u 3m	¥50,900	¥53,200	2024/4/1追加
GP-07LED	■ FAS-Digi PRO	¥875,000	¥1,000,000	
GP-07LED/PC	■ FAS-Digi PRO (PC付)	¥1,100,000	¥1,150,000	
GP-07LED/TAB	■ FAS-Digi PRO (タブレット付)	¥1,000,000	¥1,150,000	
GP-07LEDAF	FAS-Digi PRO用 ^① ア ^② ル ^③ フィルタ ^④ 1	¥16,100	¥17,500	
GP-08LED	■ FAS-Digi Compact	¥600,000	¥660,000	
GP-08LEDUPJP	■ FAS-Digi Compact ^① ア ^② ル ^③ プレート ^④ セット	¥290,000	¥330,000	
GP-FAS-X	■ FAS-X	¥1,675,000	¥1,800,000	
GVGT2001	FastDryer 1to24tubes or SiglePlate	¥251,000	¥228,000	
H2O-A-PACK	ア ^① ク ^② ム ^③ H20システム型 UV・VF用 ^④ カ ^⑤ トリ ^⑥ ジ ^⑦ キット	¥102,700	¥106,000	
H2O-B-PACK	ア ^① ク ^② ム ^③ H20システム型 UF用 ^④ カ ^⑤ トリ ^⑥ ジ ^⑦ キット	¥102,700	¥106,000	
H2O-CPR	ア ^① ク ^② ム ^③ Mini Plus用 ^④ フ ^⑤ トリ ^⑥ ト ^⑦ メント ^⑧ カ ^⑨ トリ ^⑩ ジ ^⑪ キット	¥46,100	¥47,500	
H2O-MA-T	ア ^① ク ^② ム ^③ Mini Plus 水道水供給 ^④ タイ ^⑤ プ	¥602,000	¥617,000	
H2O-MA-UV-T	ア ^① ク ^② ム ^③ Mini Plus UV 水道水供給 ^④ タイ ^⑤ プ	¥672,000	¥689,000	
H2O-MM-T	ア ^① ク ^② ム ^③ Mini 純水 ^④ マ ^⑤ ニ ^⑥ リ ^⑦ 供給 ^⑧ タイ ^⑨ プ	¥410,000	¥420,000	
H2O-MM-UV-T	ア ^① ク ^② ム ^③ Mini UV ^④ マ ^⑤ ニ ^⑥ リ ^⑦ 供給 ^⑧ タイ ^⑨ プ	¥494,000	¥506,000	
H2O-MU-T	ア ^① ク ^② ム ^③ Mini Essential 純水 ^④ マ ^⑤ ニ ^⑥ リ ^⑦ 供給 ^⑧ タイ ^⑨ プ	¥480,000	¥492,000	
H2O-MU-UV-T	ア ^① ク ^② ム ^③ Mini Essential UV 純水 ^④ マ ^⑤ ニ ^⑥ リ ^⑦ 供給 ^⑧ タイ ^⑨ プ	¥580,000	¥595,000	
H2O-S-PACK	ア ^① ク ^② ム ^③ Mini用 Scientific カ ^④ トリ ^⑤ ジ ^⑥ キット	¥37,100	¥38,200	
H2O-U-PACK	ア ^① ク ^② ム ^③ UV・VF・DI・UF用 ^④ カ ^⑤ トリ ^⑥ ジ ^⑦ キット	¥77,900	¥99,300	
HEX0004	HLS HMW DNAextraction kit 4	¥117,600	¥125,400	
HEX0012	HLS HMW DNAextraction kit 12	¥338,100	¥360,500	
HIT0004	HLS HMW DNA isolation kit 4	¥141,100	¥150,500	
HIT0012	HLS HMW DNAisolation kit 12	¥388,100	¥413,800	
HLT0001	■ HLS2 本体	¥5,200,000	¥5,550,000	
HPE7504	HT Gel ^① セット ^② 0.75% High-Pass 4	¥126,800	¥135,200	
HPE7510	HT Gel ^① セット ^② 0.75% High-Pass 10	¥255,400	¥272,300	
HPF7504	HT Gel ^① セット ^② 0.75% High-Pass 4	¥126,800	¥135,200	
HPF7510	HT Gel ^① セット ^② 0.75% High-Pass 10	¥255,400	¥272,300	
HSS0004	HLS DNA Size Selction Kit 4	¥102,900	¥110,000	
HSS0012	HLS DNA Size Selction Kit 12	¥294,000	¥310,000	
HTC1504	HT Gel ^① セット ^② 1.5% 4	¥113,900	¥121,500	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
HTC1510	HT Gelカセット 1.5% 10	¥213,200	¥227,300	
HTC2004	HT Gelカセット 2% 100-600bp 4	¥113,900	¥121,500	
HTC2010	HT Gelカセット 2% 100-600bp 10	¥213,200	¥227,300	
HTG3004	HT Gelカセット 3% 100-250bp 4	¥113,900	¥121,500	
HTG3010	HT Gelカセット 3% 100-250bp 10	¥213,200	¥227,300	
HTG3050	HT Gelカセット 3% 100-250bp 50	¥818,500	¥872,700	
HTP0001	■PippinHT 本体	¥6,200,000	¥6,950,000	
ib10310	Nikon用K7用マウント/取付サビレス 1	¥129,000	¥200,000	
ib10310-T1	Nikon用K7用マウント+T1取付サビレス 1	¥129,000	¥200,000	
ib10310-Ti2	Nikon用K7用マウント+Ti2取付サビレス 1	¥129,000	¥200,000	
ib10720	■BL Stage Top Incubation System	¥3,134,000	¥4,870,000	
ib10722	■BL Stage Top Incubation System	¥3,506,000	¥5,448,000	
ib10724	■BL Stage Top Incubation System	¥3,401,000	¥5,284,000	
ib10726	■BL Stage Top Incubation System	¥3,773,000	¥5,862,000	
ib10802	Elbow Luerコネクタ PP ST25	¥18,400	¥23,700	
ib10811	ガラスカラム-スリッパ 24x60mm No.1.5 100	¥18,400	¥23,700	
ib10812	ガラスカラム-スリッパ 25x75mm No.1.5 100	¥40,900	¥52,900	
ib10813	ガラス-スリッパ 未コーティング No.1.5 100	¥40,900	¥52,900	
ib10814	ガラス-スリッパ ibiTreat No.1.5 100	¥40,900	¥52,900	
ib10815	ガラスカラム-スリッパ ヒートアップツール付 100	¥33,100	¥42,700	
ib10816	ガラスカラム-スリッパ Grid-500 No.1.5 100	¥106,600	¥137,600	
ib10817	ガラスカラム-スリッパ Grid-50 No.1.5 100	¥106,600	¥137,600	
ib10820	インラインアークインジェクションホースST	¥18,400	¥23,700	
ib10821	ホースクリップ	¥17,300	¥22,400	
ib10822	ルアープラグ (ス) ST	¥18,400	¥23,700	
ib10823	ルアーコネクタ (ス) ST	¥18,400	¥23,700	
ib10824	ルアーコネクタ Male PP ST25	¥18,400	¥23,700	
ib10825	ルアーコネクタ (ス) ST	¥18,400	¥23,700	
ib10826	ルアーコネクタ (ス) ST	¥18,400	¥23,700	
ib10827	0.8mmID Yチューブフィッティング ST25	¥18,400	¥23,700	
ib10828	1.6mmID Yチューブフィッティング ST25	¥18,400	¥23,700	
ib10829	フィッティングレギュレータ-0.5-1.6ST	¥18,400	¥23,700	
ib10830	uSlide用シリアルコネクタ ST2x6	¥28,900	¥37,300	
ib10831	チューブアダプターセット ST6x2	¥28,900	¥37,300	
ib10840	シリコンチューブ 0.5mm IDx5mST	¥21,000	¥27,100	
ib10841	シリコンチューブ 0.8mm IDx5mST	¥21,000	¥27,100	
ib10842	シリコンチューブ 1.5mm IDx5mST	¥20,900	¥27,000	
ib10871	ibiSeal 20x56 ST60	¥24,300	¥31,400	
ib10872	ibiSeal 22x48 ST60	¥24,300	¥31,400	
ib10873	ibiSeal 25x52 ST60	¥24,300	¥31,400	
ib10874	ibiSeal 76x114 ST30	¥24,300	¥31,400	
ib10875	ibiSeal 71x107 ST30	¥24,300	¥31,400	
ib10876	ibiSeal 7x56 ST60	¥24,300	¥31,400	
ib10877	ibiSeal 8x8 ST120	¥24,300	¥31,400	
ib10878	ibiSeal 35 ST60	¥24,300	¥31,400	
ib10902	■ ibidi pump system	¥3,860,000	¥4,990,000	
ib10903	■ fluidic unit(ibidi pump system用)	¥527,000	¥593,000	
ib10905	■ ibidi air pressure pump	オープン価格	¥4,280,000	
ib10906	■ ibidi pump system(Fluidic Quad)	¥5,090,000	¥6,750,000	
ib10910	■ ibidi pump system:4Fluidic Units	¥5,090,000	¥6,750,000	
ib10915	■ BL Heating system 96 ウェルプレート用	¥769,000	¥1,194,000	
ib10918	■ BL ヒーティングシステム Universal Fit	¥1,466,000	¥2,278,000	
ib10927	■ BL Heating System, Universal Fit	¥1,618,000	¥2,514,000	
ib10928	BL HeatedPlate in MultiWell Format	¥481,000	¥747,000	
ib10929	■ BL ヒーティングシステム 8ウェルプレートHKFrame	¥1,732,000	¥2,692,000	
ib10929-PL	BL Heated Plate & Heated Lid 1	¥1,198,000	¥1,861,000	
ib10932	BL ヒーティング インサート(uDish 35mm Low)	¥91,000	¥141,000	
ib10933	BL ヒーティング インサート(all ibidi uslide)	¥91,000	¥141,000	
ib10934	BL ヒーティング インサート(uDish 35mmHigh)	¥91,000	¥141,000	
ib10936	BL ヒーティング インサート アダプター	¥46,000	¥72,000	
ib10937	BL ヒーティング インサート(LabTekカラムケース)	¥91,000	¥141,000	
ib10961	Perfusion Setブルー 15cm ID0.8mm 4	¥95,000	¥111,700	
ib10962	Perfusion Setレッド 15cm ID1.6mm 4	¥95,000	¥111,700	
ib10963	Perfusion Set緑 50cm ID0.8mm 4	¥95,000	¥111,700	
ib10964	Perfusion Set黄緑 50cm ID1.6mm 4	¥95,000	¥111,700	
ib10965	Perfusion Set白 15cm ID 0.5mm 4	¥95,000	¥111,700	
ib10966	Perfusion Setブラック 50cm ID 0.5mm 4	¥95,000	¥111,700	
ib10967	Perfusion Setブルー 100cm ID0.5mm 4	¥95,000	¥111,700	
ib10968	Perfusion Setグレー 100cm ID0.8mm 4	¥95,000	¥111,700	
ib10969	Perfusion Setブラック 100cm ID1.6mm 4	¥95,000	¥111,700	
ib10971	Filter/Reservoir set 10ml ST10	¥66,700	¥75,000	
ib10972	Filter/Reservoir set 2ml ST10	¥66,700	¥75,000	
ib10974	Filter/Reservoir set 50 ml ST10	¥66,700	¥75,000	
ib10976	Reservoir Holder 10ml	¥25,000	¥34,200	
ib10977	Reservoir Holder 2 ml	¥25,000	¥34,200	
ib10978	Reservoir Holder 50 ml	¥25,000	¥34,200	
ib10982	Pump Demo Set	¥40,000	¥43,000	
ib11920	■BL Gas Incubation System for CO2	¥1,669,000	¥2,592,000	
ib11921	■SL Gas Incubation System CO2	¥1,669,000	¥2,503,000	
ib11922	■BL Gas Incubation System CO2&O2	¥2,040,000	¥3,170,000	
ib11923	■SL Gas Incubation System CO2/O2	¥2,040,000	¥3,030,000	
ib11929-115	■BL I7° レジスタ-Generator 1bar ready	¥211,000	¥328,000	
ib11929-230	■BL I7° レジスタ-Generator 1bar ready	¥211,000	¥328,000	
ib11930	BLModificationHeatedLidTypeA,w2129	¥204,000	¥317,000	
ib11933	SL Modification of Heated Lid	¥206,000	¥320,000	
ib11934	SL Modification of Heated Lid	¥412,000	¥640,000	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
ib12110	■SL Heating System Slide/Dish	¥2,014,000	¥2,989,000	
ib12111	■SL Incubation Chamber Slide/Dish	¥1,393,000	¥2,165,000	
ib12116	SilverLine Slide Holder	¥109,000	¥169,000	
ib12117	SilverLine Dish Holder	¥144,000	¥224,000	
ib12130	■SL Heating System 4 Slides	¥1,786,000	¥2,775,000	
ib12131	■SL Incubation Chamber 4 Slides	¥1,165,000	¥1,811,000	
ib12150	■SL Heating System マルチウェルプレート	¥1,887,000	¥2,793,000	
ib12151	■SL Incubation Chamber マルチウェルプレート	¥1,267,000	¥1,968,000	
ib12170	SL Objective Heater Universal	¥132,000	¥205,000	
ib12720	■SL Stage Top Incubator	¥3,535,000	¥5,492,000	
ib12722	■SL Stage Top Incubator	¥3,874,000	¥6,019,000	
ib12724	■SL Stage Top Incubator	¥3,408,000	¥5,296,000	
ib12726	■SL Stage Top Incubator	¥3,747,000	¥5,823,000	
IB1651512	オートラジオグラフィ用フィルム 35x43cm 50	¥1,365,600	¥1,497,700	
IB1660760	オートラジオグラフィ用フィルム 13x18cm 50	¥590,600	¥647,800	
ib50001	Mounting Medium 15ml	¥11,200	¥14,400	
ib50001-4	Mounting Medium 4x15ml	¥39,500	¥50,700	
ib50011	Mounting Medium DAPI 15ml付 ④	¥12,600	¥16,200	
ib50011-4	Mounting Medium DAPI 4x15ml付 ④	¥44,800	¥57,500	
ib50051	Anti-Evaporation Oil 125ml	¥39,400	¥50,500	
ib50101	Immersion Oil 15ml	¥14,700	¥18,900	
ib50201	コラーゲンタイプ1 rat 5mg/ml 1x5ml ②	¥40,400	¥48,400	
ib50202	コラーゲンタイプ1 rat 5mg/ml 4x5ml ②	¥113,100	¥135,500	
ib50203	コラーゲンタイプ1 rat 5mg/ml 100ml ②	¥484,500	¥580,900	
ib50204	コラーゲンタイプ1 rat 10mg/ml 1x5ml ②	¥76,000	¥91,100	
ib50205	コラーゲンタイプ1 rat 10mg/ml 4x5ml ②	¥212,800	¥255,100	
ib50206	コラーゲンタイプ1 rat 10mg/ml 100ml ②	¥912,000	¥1,093,400	
ib50301	コラーゲンタイプ1 bovine 5mg/ml 1x5ml ④	¥40,400	¥48,400	
ib50302	コラーゲンタイプ1 bovine 5mg/ml 4x5ml ④	¥113,100	¥135,500	
ib50303	コラーゲンタイプ1 bovine 5mg/ml 100ml ④	¥484,500	¥580,900	
ib80003	uSlide専用ツク MetalRack	¥29,400	¥37,900	
ib80022	freezing medium 5x20ml ④	¥46,900	¥60,200	
ib80023	freezing medium 120ml ④	¥47,400	¥60,800	
ib80030	ibidi uSlide microscopy rack	¥55,100	¥71,200	
ib80031	Magnetic Lid for uSlides 4	¥80,800	¥104,400	
ib80035	ibidi uDish microscopy rack base	¥50,700	¥78,800	
ib80036	magnetic lid for uDish 35 mm low	¥17,700	¥27,400	
ib80037	magnetic lid for uDish 35 mm high	¥17,700	¥27,400	
ib80040	Clamp for stickySlides 1	¥18,400	¥23,700	
ib80042	Adapter for stickySlide I Luer 1	¥91,900	¥118,600	
ib80043	Adapter for stickySlide VI 0.4 1	¥91,900	¥118,600	
ib80044	Adapter for stickySlide Chemotaxis1	¥91,900	¥118,600	
ib80045	Adapter stickySlide 8/8h/18h用 1	¥91,900	¥118,600	
ib80050	Optional DIC lid ST15	¥27,900	¥36,000	
ib80055	DIC Lid for uSlides ST15	¥27,900	¥36,000	
ib80106	uSlide I ibiTreat ST15	¥47,200	¥61,000	
ib80126	uSlide y-shaped ibiTreat ST15	¥51,900	¥67,000	
ib80131	uDish 35mm low 未コティング ST60	¥53,500	¥69,100	
ib80136	uDish 35mm low ibiTreat ST60	¥53,500	¥69,100	
ib80137	uDish 35mm low カラ底 ST60	¥59,300	¥76,600	
ib80156	uDish35mm low grid500 ibiTreat ST60	¥58,100	¥75,000	
ib80161	uSlide I 0.2Luer 未コティング ST15	¥47,200	¥61,000	
ib80166	uSlide I 0.2 Luer ibiTreat ST15	¥47,200	¥61,000	
ib80167	uSlide I 0.2Luer カラ底 ST15	¥63,400	¥81,800	
ib80168	Sticky-Slide I 0.2Luer ST 15	¥55,600	¥71,800	
ib80171	uSlide I 0.4Luer 未コティング ST15	¥47,200	¥61,000	
ib80172	uSlide I 0.4Luer コラーゲンIV ST15	¥76,200	¥98,400	
ib80176	uSlide I 0.4Luer ibiTreat ST15	¥47,200	¥61,000	
ib80177	uSlide I 0.4Luer カラ底 ST15	¥63,400	¥81,800	
ib80178	Sticky-Slide I 0.4Luer ST15	¥55,600	¥71,800	
ib80181	uSlide I 0.6Luer 未コティング S15	¥47,200	¥61,000	
ib80186	uSlide I 0.6 Luer ibiTreat S15	¥47,200	¥61,000	
ib80187	uSlide I 0.6Luer カラ底 ST15	¥63,400	¥81,800	
ib80188	Sticky-Slide I 0.6Luer ST 15	¥55,600	¥71,800	
ib80191	uSlide I 0.8Luer 未コティング S15	¥47,200	¥61,000	
ib80196	uSlide I 0.8 Luer ibiTreat S15	¥47,200	¥61,000	
ib80197	uSlide I 0.8Luer カラ底 ST15	¥63,400	¥81,800	
ib80198	Sticky-Slide I 0.8Luer ST 15	¥55,600	¥71,800	
ib80206	Culture-Insert 2ウェルuDish low ST 30	¥66,100	¥85,400	
ib80209	25 Culture-Inserts 2ウェル ST 25	¥35,200	¥45,400	
ib80209-150	25 Culture-Inserts 2ウェル ST 6x25	¥182,100	¥235,100	
ib80242	Culture Insert 24ウェル ST 3	¥123,300	¥148,800	
ib80281	uSlide 2ウェル ST ノコト 15	¥36,300	¥46,900	
ib80286	uSlide 2ウェル ST ibiTreat 15	¥36,300	¥46,900	
ib80287	uSlide 2ウェル カラ底 ST 15	¥54,600	¥70,500	
ib80296	uSlide2ウェル Ph+ ST ibiTreat15	¥39,900	¥51,500	
ib80297	uSlide2ウェル Ph+ カラ底 ST 15	¥54,600	¥70,500	
ib80311	uSlideⅢ3in1未コティング ST15	¥58,800	¥75,900	
ib80316	uSlideⅢ3in1ibiTreat ST15	¥58,800	¥75,900	
ib80322	uSlide Chemotaxis コラーゲンIV 10	¥103,600	¥125,000	
ib80326	uSlide Chemotaxis ibiTreat 10	¥103,600	¥125,000	
ib80328	sticky-Slide Chemotaxis ST 10	¥113,800	¥137,300	
ib80350	uSlideSpheroidPerfusionBioinert15	¥110,000	¥132,700	
ib80351	uSlideSpheroidPerfusion未コティング 15	¥83,600	¥107,900	
ib80356	uSlideSpheroidPerfusionibiTreat15	¥83,600	¥107,900	
ib80366	Culture-Insert3ウェル in uDish ST 30	¥68,200	¥88,100	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
ib80369	25 Culture-Inserts 3well ST 25	¥35,700	¥46,100	
ib80371	uSlide III 3D Perfusion 未コーティング	¥69,800	¥84,200	
ib80376	uSlide III 3D Perfusion ibiTreat	¥69,800	¥84,200	
ib80381	3 well Chamber ST 15	¥52,500	¥63,300	
ib80406	Dish 35mm High ibiTreat S30	¥73,500	¥94,900	
ib80409	25 Micro-Insert 4well ST	¥42,500	¥54,900	
ib80411	uDish Quad 未コーティング ST 60	¥103,900	¥134,200	
ib80416	uDish 35mm Quad ibiTreat ST 60	¥103,900	¥134,200	
ib80420	uSlide 4 well Bioinert ST 15	¥83,600	¥107,900	
ib80421	uSlide 4well ST ノット 15	¥36,300	¥46,900	
ib80424	uSlide 4well ST ホールタイプ 15	¥98,700	¥127,400	
ib80426	uSlide 4well ST ibiTreat 15	¥36,300	¥46,900	
ib80427	uSlide 4well ガラス底 ST15	¥54,600	¥70,500	
ib80446	uSlide 4well Ph+ ST ibiTreat15	¥39,900	¥51,500	
ib80447	uSlide 4well Ph+ ガラス底 ST15	¥54,600	¥70,500	
ib80466	Culture-Insert 4well in uDish ST 30	¥75,600	¥97,600	
ib80469	25 Culture-Inserts 4well ST 25	¥43,600	¥56,200	
ib80486	micro-Insert 4well FulTrac in uDish	¥95,500	¥123,300	
ib80489	25 micro-Inserts 4well FulTrac ST	¥58,300	¥75,200	
ib80600	uSlide VI 0.4 Bioinert ST 15	¥92,300	¥119,100	
ib80601	uSlide VI flow 未コーティング ST15	¥53,000	¥68,400	
ib80602	uSlide VI flow コーティング ST15	¥98,700	¥127,400	
ib80604	uSlide VI 0.4 Poly-L-Lysine ST15	¥98,700	¥127,400	
ib80606	uSlide VI 0.4 flow ibiTreat ST15	¥53,000	¥68,400	
ib80606-90	uSlide VI 0.4 flow ibiTreat ST 90	¥276,600	¥349,000	
ib80607	uSlide VI 0.5 ガラス底 ST 15	¥73,000	¥94,200	
ib80608	Sticky-Slide VI 0.4 ST 15	¥61,900	¥80,000	
ib80609	uSlide VI 0.4 コーティング I ST 15	¥98,700	¥127,400	
ib80621	uSlide VI flat 未コーティング ST15	¥53,000	¥68,400	
ib80626	uSlide VI flat ibiTreat ST15	¥53,000	¥68,400	
ib80661	uSlide VI 0.1 未コーティング ST15	¥58,300	¥75,200	
ib80666	uSlide VI 0.1 ibiTreat ST15	¥58,300	¥75,200	
ib80800	uSlide 8well high Bioinert ST15	¥83,600	¥107,900	
ib80801	uSlide 8well high 未コーティング ST 15	¥36,300	¥46,900	
ib80802	uSlide 8well high コーティング S 15	¥98,700	¥127,400	
ib80804	uSlide 8well high Poly-L-Lysine S15	¥98,700	¥127,400	
ib80806	uSlide 8well high ibiTreat 15	¥36,300	¥46,900	
ib80806-90	uSlide 8well high ibiTreat 90	¥181,100	¥239,300	
ib80806-G500	uSlide 8well high ibiTreat Grid500	¥60,400	¥77,900	
ib80807	uSlide 8well high ガラス底 15	¥54,600	¥70,500	
ib80807-90	uSlide 8well high ガラス底 90	¥287,100	¥359,400	
ib80808	sticky-Slide 8well high 15	¥50,900	¥65,700	
ib80809	uSlide 8well high コーティング I ST 15	¥98,700	¥127,400	
ib80816	Cedex XS Analyzer用 Smart Slide S15	¥49,900	¥60,200	
ib80821	uSlide 8well 未コーティング ST 15	¥36,300	¥46,900	
ib80822	uSlide 8well コーティング IV ST15	¥98,700	¥127,400	
ib80824	uSlide 8well Poly-L-Lysine ST15	¥98,700	¥127,400	
ib80826	uSlide 8well ibiTreat ST 15	¥36,300	¥46,900	
ib80826-90	uSlide 8well ibiTreat ST 90	¥181,100	¥239,300	
ib80826-G500	uSlide 8well ibiTreat ST Grid500 15	¥60,400	¥77,900	
ib80827	uSlide 8well ガラス底 ST15	¥54,600	¥70,500	
ib80827-90	uSlide 8well ガラス底 ST 90	¥287,100	¥359,400	
ib80828	Sticky-Slide 8well ST 15	¥50,900	¥65,700	
ib80829	uSlide 8well コーティング I ST 15	¥98,700	¥127,400	
ib80841	8 well Chamber ST 15	¥52,500	¥63,300	
ib81101	uSlide I Luerロック 未コーティング S16	¥53,000	¥68,400	
ib81106	uSlide I Luerロック ibiTreat ST16	¥53,000	¥68,400	
ib81128	Sticky-Slide I 0.1 Luer ST 15	¥55,600	¥71,800	
ib81148	uDish 35mm High ガラス底 Grid-50 30	¥75,600	¥97,600	
ib81150	uDish 35mm High Bioinert ST 30	¥73,000	¥88,000	
ib81151	uDish 35mm High 未コーティング ST60	¥53,500	¥69,100	
ib81156	uDish 35mm High ibiTreat ST 60	¥53,500	¥69,100	
ib81156-400	uDish 35mm High ibiTreat ST 400	¥311,800	¥391,700	
ib81158	uDish 35mm High ガラス底 60	¥59,300	¥76,600	
ib81158-400	uDish 35mm High ガラス底 400	¥343,800	¥433,900	
ib81161	uDish 35mm High grid500 未コーティング	¥58,100	¥75,000	
ib81166	uDish 35mm High grid500 ibiTreat	¥58,100	¥75,000	
ib81168	uDish 35mm High ガラス底 Grid-500 30	¥75,600	¥97,600	
ib81176	Culture-Insert 2well uDish high ST30	¥66,100	¥85,400	
ib81191	uDish 35mm High ESS 28kPa 未コーティング ST30	¥63,000	¥81,300	
ib81199	uDish 35mm High ESS 未コーティング ST 30	¥70,300	¥90,800	
ib81201	uChamber 12well ST 15	¥52,500	¥63,300	
ib81218-200	ガラス底 Dish 35mm ST 200	¥113,400	¥136,800	
ib81218-800	ガラス底 Dish 35mm ST 800	¥435,700	¥465,100	
ib81291	uDish 35mm High ESS 1.5kPa 未コーティング ST30	¥63,000	¥81,300	
ib81391	uDish 35mm High ESS 15kPa 未コーティング ST30	¥63,000	¥81,300	
IB8142358	オートクレーブ用フィルム 1入	¥120,100	¥131,700	
ib81501	uSlide 15well 3D 未コーティング ST15	¥66,100	¥85,400	
ib81506	uSlide 15well 3D ibiTreat ST15	¥66,100	¥85,400	
ib81507	uSlide 15well 3D ガラス底 ST15	¥66,100	¥85,400	
ib81806	uSlide 2x9well ibiTreat ST 15	¥39,700	¥51,200	
ib81811	uSlide 18well 未コーティング ST 15	¥36,300	¥46,900	
ib81816	uSlide 18well ibiTreat ST 15	¥36,300	¥46,900	
ib81817	uSlide 18well ガラス底 ST 15	¥54,600	¥70,500	
ib81818	Sticky-Slide 18well ST 15	¥50,900	¥65,700	
ib81821	uSlide 18well flat 未コーティング ST15	¥39,900	¥51,500	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
ib81826	uSlide187μm flat ibiTreatST15	¥39,900	¥51,500	
ib82101	uSlide 17μm 未コーティング ST 15	¥36,300	¥46,900	
ib82106	uSlide 17μm ibiTreat ST 15	¥36,300	¥46,900	
ib82107	uSlide 17μm カラ底 ST 15	¥53,900	¥69,600	
ib82421	uPlate 247μm 未コーティング ST 15	¥85,000	¥109,800	
ib82426	uPlate 247μm ibiTreat ST 15	¥85,000	¥109,800	
ib83601	uSlide VI 0.4 uPattern sqr20 10	¥151,100	¥182,300	
ib83601S	uSlide VI 0.4 uPattern sqr20 2	¥34,000	¥41,000	
ib83602	uSlide VI 0.4 uPattern 10	¥151,100	¥182,300	
ib83602S	uSlide VI 0.4 uPattern 2	¥34,000	¥41,000	
ib83603	uSlide VI 0.4 uPattern sqr30,10	¥151,100	¥182,300	
ib83603S	uSlide VI 0.4 uPattern sqr30 2	¥34,000	¥41,000	
ib83651	uSlide VI0.4 uPattern テストパターン1 2口	¥44,800	¥54,100	
ib83652	uSlide VI0.4 uPattern テストパターン2, 2	¥44,800	¥54,100	
ib83653	uSlide VI0.4 uPattern テストパターン3, 2	¥44,800	¥54,100	
ib83801	uSlide 87μm high uPattern sqr20 10	¥151,100	¥182,300	
ib83801S	uSlide 87μm high uPattern sqr20 2	¥34,000	¥41,000	
ib83802	uSlide 87μm high uPattern 10	¥151,100	¥182,300	
ib83802S	uSlide 87μm high uPattern 2	¥34,000	¥41,000	
ib83803	uSlide 87μm high uPattern sqr30,10	¥151,100	¥182,300	
ib83803S	uSlide 87μm high uPattern sqr30,2	¥34,000	¥41,000	
ib83851	uSlide 87μm high uPattern テストパターン1 2	¥44,800	¥54,100	
ib83852	uSlide 87μm high uPattern テストパターン2,2	¥44,800	¥54,100	
ib83853	uSlide 87μm high uPattern テストパターン3,2	¥44,800	¥54,100	
ib85216	uSlide ibiPore SiN 0.5um/20% 10	¥217,300	¥262,200	
ib85216S	uSlide ibiPoreSiN 0.5um/20% Trial 2	¥51,700	¥62,300	
ib85226	uSlide ibiPore SiN 3.0um/5% 10	¥217,300	¥262,200	
ib85226S	uSlide ibiPoreSiN 3.0um/5% Trial 2	¥51,700	¥62,300	
ib85236	uSlide ibiPore SiN 5.0um/5% 10	¥217,300	¥262,200	
ib85236S	uSlide ibiPore SiN 5.0um/5% Trial 2	¥51,700	¥62,300	
ib85246	uSlide ibiPore SiN 8.0um/5% 10	¥217,300	¥262,200	
ib85246S	uSlide ibiPore SiN 8.0um/5% Trial 2	¥51,700	¥62,300	
IB8646770	オートクレーブ用フィルム	¥203,400	¥223,100	
ib87171	uSlide I Luer 3D 未コーティング ST 15	¥73,000	¥94,200	
ib87176	uSlide I Luer 3D ibiTreat ST 15	¥73,000	¥94,200	
ib88407	uPlate 3847μm Blackカラ底 S 15	¥165,200	¥199,300	
ib89601	uPlate 967μm Round 未コーティング ST 15	¥85,000	¥109,800	
ib89606	uPlate 967μm Round ibiTreat ST 15	¥85,000	¥109,800	
ib89621	uPlate 967μm Square 未コーティング ST 15	¥85,000	¥109,800	
ib89626	uPlate 967μm Square ibiTreat ST 15	¥85,000	¥109,800	
ib89627	uPlate 967μm Blackカラ底 S 15	¥118,400	¥142,900	
ib89646	uPlate 967μm 3D ibiTreat ST15	¥139,600	¥168,500	
IC7667194	Cleaning Solution ES-7X 4本入	¥185,100	¥128,200	
IR1200060	広口ガラス瓶(キャップ付き)60ml Short	¥53,600	¥62,100	
IR1200125	広口ガラス瓶(キャップ付き)125ml Short	¥58,500	¥68,000	
IR1200250	広口ガラス瓶(キャップ付き)250ml Short	¥35,900	¥42,800	
IR1200500	広口ガラス瓶(キャップ付き)500ml Short	¥38,500	¥44,900	
IR1202000	広口ガラス瓶(キャップ付き)2 000ml Short	¥42,000	¥48,600	
IR1210125	広口ガラス瓶(キャップ付き)125ml Tall	¥29,900	¥34,700	
IR1210250	広口ガラス瓶(キャップ付き)250ml Tall	¥34,000	¥40,300	
IR1210500	広口ガラス瓶(キャップ付き)500ml Tall	¥37,000	¥43,900	
IR1211000	広口ガラス瓶(キャップ付き)1000ml Tall	¥44,900	¥52,200	
IR1212000	広口ガラス瓶(キャップ付き)2 000ml Tall	¥40,200	¥45,400	
IR1214000	広口ガラス瓶(キャップ付き)4 000ml Tall	¥44,100	¥51,100	
IR1290125	ガラス瓶 ホスト丸瓶(キャップ付)透明125ml	¥25,200	¥29,200	
IR1290250	ガラス瓶 ホスト丸瓶(キャップ付)透明250ml	¥22,100	¥25,600	
IR1290500	ガラス瓶 ホスト丸瓶(キャップ付)透明500ml	¥39,000	¥45,100	
IR1291000	ガラス瓶 ホスト丸瓶(キャップ付)透明1000ml	¥37,300	¥44,400	
IR1400060	広口ガラス瓶(キャップ付き)60ml Short	¥53,300	¥61,800	
IR1400120	広口ガラス瓶(キャップ付き)120ml Short	¥109,700	¥127,300	
IR1410060	広口ガラス瓶(キャップ付き)60ml 褐色	¥53,300	¥61,800	
IR1410120	広口ガラス瓶(キャップ付き)120ml 褐色	¥29,600	¥34,300	
IR1410250	広口ガラス瓶(キャップ付き)250ml 褐色	¥31,800	¥36,700	
IR1410500	広口ガラス瓶(キャップ付き)500ml 褐色	¥36,400	¥42,800	
IR1410950	広口ガラス瓶(キャップ付き)950ml 褐色	¥43,100	¥49,900	
IR1490125	ガラス瓶 ホスト丸瓶(キャップ付)褐色125ml	¥20,300	¥23,500	
IR1490250	ガラス瓶 ホスト丸瓶(キャップ付)褐色250ml	¥24,800	¥27,900	
IR1490500	ガラス瓶 ホスト丸瓶(キャップ付)褐色500ml	¥34,900	¥40,300	
IR1491000	ガラス瓶 ホスト丸瓶(キャップ付)褐色1000ml	¥42,800	¥49,500	
IR2200060	広口ガラス瓶(キャップ付き)60ml Short	¥60,800	¥70,300	
IR2200125	広口ガラス瓶(キャップ付き)125ml Short	¥68,600	¥79,700	
IR2200250	広口ガラス瓶(キャップ付き)250ml Short	¥38,300	¥44,400	
IR2200500	広口ガラス瓶(キャップ付き)500ml Short	¥40,800	¥47,300	
IR2202000	広口ガラス瓶(キャップ付き)2 000ml Short	¥49,600	¥55,500	
IR2210125	広口ガラス瓶(キャップ付き)125ml Tall	¥34,700	¥40,200	
IR2210250	広口ガラス瓶(キャップ付き)250ml Tall	¥38,600	¥44,700	
IR2210500	広口ガラス瓶(キャップ付き)500ml Tall	¥41,600	¥48,200	
IR2211000	広口ガラス瓶(キャップ付き)1 000ml Tall	¥50,000	¥57,900	
IR2212000	広口ガラス瓶(キャップ付き)2 000ml Tall	¥39,100	¥45,300	
IR2214000	広口ガラス瓶(キャップ付き)4 000ml Tall	¥52,100	¥60,400	
IR2290125	ガラス瓶 ホスト丸瓶(キャップ付)透明125ml	¥26,700	¥31,000	
IR2290250	ガラス瓶 ホスト丸瓶(キャップ付)透明250ml	¥25,900	¥29,100	
IR2290500	ガラス瓶 ホスト丸瓶(キャップ付)透明500ml	¥38,900	¥45,000	
IR2291000	ガラス瓶 ホスト丸瓶(キャップ付)透明1000ml	¥41,000	¥47,500	
IR2400060	広口ガラス瓶(キャップ付き)60ml Short	¥98,000	¥113,100	
IR2400120	広口ガラス瓶(キャップ付き)120ml Short	¥87,700	¥102,000	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
IR2410060	広口ガラス瓶(キャップ付き褐色)60ml	¥62,600	¥72,200	
IR2410120	広口ガラス瓶(キャップ付き褐色)120ml	¥34,600	¥40,000	
IR2410250	広口ガラス瓶(キャップ付き褐色)250ml	¥38,000	¥44,000	
IR2410500	広口ガラス瓶(キャップ付き褐色)500ml	¥41,500	¥48,100	
IR2410950	広口ガラス瓶(キャップ付き褐色)950ml	¥49,600	¥57,500	
IR2490125	ガラス瓶 ホストン丸瓶(キャップ付)褐色125ml	¥25,200	¥29,100	
IR2490250	ガラス瓶 ホストン丸瓶(キャップ付)褐色250ml	¥25,200	¥29,100	
IR2490500	ガラス瓶 ホストン丸瓶(キャップ付)褐色500ml	¥35,200	¥40,700	
IR2491000	ガラス瓶 ホストン丸瓶(キャップ付)褐色1000ml	¥44,700	¥51,700	
IR2880022	Septa Tef Silicone 24個入	¥25,500	¥29,500	
IR2880022BP	Septa バル 2000個入	¥813,900	¥942,000	
IR2972400UB	Closure Open Septa 24個入	¥29,100	¥34,600	
IR2972414	オープントップキャップのみ 24入	¥25,800	¥29,400	
IR3110500	ポリ瓶(キャップ付き高密度ポリエチレン製)	¥82,400	¥95,600	
IR3111000	ポリ瓶(キャップ付き高密度ポリエチレン製)	¥51,300	¥59,800	
IR3120250	ポリ瓶(キャップ付き高密度ポリエチレン製)	¥46,300	¥52,300	
IR3120500	ポリ瓶(キャップ付き高密度ポリエチレン製)	¥68,100	¥76,800	
IR3120950	ポリ瓶(キャップ付き高密度ポリエチレン製)	¥39,900	¥45,100	
IR3130250	ポリ瓶(キャップ付き高密度ポリエチレン製)	¥45,100	¥50,900	
IR3130500	ポリ瓶(キャップ付き高密度ポリエチレン製)	¥55,600	¥62,800	
IR3131000	ポリ瓶(キャップ付き高密度ポリエチレン製)	¥32,300	¥37,400	
IR3140001	液体コンテナー	¥80,700	¥91,600	
IR3140005	液体コンテナー	¥45,200	¥51,000	
IR3140025	液体コンテナー	¥100,100	¥113,200	
IR3141000	液体コンテナー	¥55,400	¥62,500	
IR3190125	ポリ瓶(キャップ付き高密度ポリエチレン製)	¥103,000	¥116,400	
IR3191000	ポリ瓶(キャップ付き高密度ポリエチレン製)	¥47,600	¥53,800	
IR3200060	広口ガラス瓶(キャップ付き) 60ml Short	¥69,600	¥80,800	
IR3200125	広口ガラス瓶(キャップ付き)125ml Short	¥77,700	¥89,400	
IR3200250	広口ガラス瓶(キャップ付き)250ml Short	¥43,800	¥50,600	
IR3200500	広口ガラス瓶(キャップ付き)500ml Short	¥46,900	¥54,300	
IR3202000	広口ガラス瓶(キャップ付き)2 000ml Short	¥53,800	¥62,300	
IR3210125	広口ガラス瓶(キャップ付き)125ml Tall	¥39,300	¥45,500	
IR3210250	広口ガラス瓶(キャップ付き)250ml Tall	¥43,300	¥50,100	
IR3210500	広口ガラス瓶(キャップ付き)500ml Tall	¥47,400	¥54,800	
IR3211000	広口ガラス瓶(キャップ付き)1 000ml Tall	¥57,400	¥66,200	
IR3212000	広口ガラス瓶(キャップ付き)2 000ml Tall	¥44,400	¥51,300	
IR3214000	広口ガラス瓶(キャップ付き)4 000ml Tall	¥58,100	¥67,500	
IR3290125	ガラス瓶 ホストン丸瓶(キャップ付)透明125ml	¥30,900	¥35,800	
IR3290250	ガラス瓶 ホストン丸瓶(キャップ付)透明250ml	¥29,400	¥34,000	
IR3291000	ガラス瓶 ホストン丸瓶(キャップ付)透明1000ml	¥47,300	¥54,800	
IR3400060	広口ガラス瓶(キャップ付き)60ml Short	¥69,200	¥80,500	
IR3400120	広口ガラス瓶(キャップ付き)120ml Short	¥161,800	¥186,200	
IR3410060	広口ガラス瓶(キャップ付き褐色)60ml	¥71,100	¥82,500	
IR3410120	広口ガラス瓶(キャップ付き褐色)120ml	¥39,500	¥45,700	
IR3410250	広口ガラス瓶(キャップ付き褐色)250ml	¥43,800	¥50,700	
IR3410500	広口ガラス瓶(キャップ付き褐色)500ml	¥46,300	¥53,500	
IR3410950	広口ガラス瓶(キャップ付き褐色)950ml	¥58,000	¥66,900	
IR3490125	ガラス瓶 ホストン丸瓶(キャップ付)褐色125ml	¥25,500	¥29,500	
IR3490250	ガラス瓶 ホストン丸瓶(キャップ付)褐色250ml	¥28,600	¥33,100	
IR3490500	ガラス瓶 ホストン丸瓶(キャップ付)褐色500ml	¥36,200	¥41,800	
IR3491000	ガラス瓶 ホストン丸瓶(キャップ付)褐色1000ml	¥48,300	¥58,200	
IR5030004	バル(環境検査用/カート)	¥70,500	¥81,900	
IR5030006	バル(環境検査用/カート)	¥70,800	¥82,100	
IR5030007	バル(環境検査用/カート)	¥71,200	¥82,600	
IR5030008	バル(環境検査用/カート)	¥70,100	¥81,300	
IR5030009	バル(環境検査用/カート)	¥70,500	¥81,700	
IR5030010	バル(環境検査用/カート)	¥70,200	¥81,400	
IRC1260020	ガラスバケツ	¥74,200	¥85,900	
IRC1360040	ガラスバケツ	¥82,800	¥95,800	
IRC1460040	ガラスバケツ	¥106,400	¥123,100	
IRC2260020	ガラスバケツ	¥106,800	¥123,300	
IRC2360040	ガラスバケツ	¥120,300	¥138,800	
IRC2460040	ガラスバケツ	¥131,700	¥152,000	
IRC3260020	ガラスバケツ	¥127,300	¥147,100	
IRC3360040	ガラスバケツ	¥151,900	¥175,600	
IRC3460040	ガラスバケツ	¥153,100	¥177,000	
IRLB2360040	ガラスバケツ	¥119,500	¥137,900	
IRS1260020	ガラスバケツ	¥70,200	¥81,400	
IRS1290250	ガラス瓶 ホストン丸瓶(オープントップキャップ付)	¥25,900	¥29,600	
IRS1360040	ガラスバケツ	¥77,700	¥89,400	
IRS1360060	ガラスバケツ	¥101,900	¥117,800	
IRS1460040	ガラスバケツ	¥92,100	¥106,000	
IRS1490125	ガラス瓶 ホストン丸瓶(オープントップキャップ付)	¥24,200	¥27,100	
IRS2260020	ガラスバケツ	¥99,200	¥114,700	
IRS2360040	ガラスバケツ	¥111,400	¥129,300	
IRS2460040	ガラスバケツ	¥114,500	¥132,800	
IRS3260020	ガラスバケツ	¥122,500	¥141,300	
IRS3290250	ガラス瓶 ホストン丸瓶(オープントップキャップ)	¥31,700	¥36,600	
IRS3360040	ガラスバケツ	¥134,200	¥155,300	
IRS3460040	ガラスバケツ	¥134,200	¥155,300	
IRS3490125	ガラス瓶 ホストン丸瓶(オープントップキャップ)	¥28,400	¥32,900	
IRS3490250	ガラス瓶 ホストン丸瓶(オープントップキャップ)	¥32,100	¥37,200	
IRT1360040	ガラスバケツ	¥101,300	¥117,300	
IRT2360040	ガラスバケツ	¥145,700	¥168,700	
IRT2460040	ガラスバケツ	¥144,200	¥166,900	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
IRT3360040	ガラスバイアル	¥174,600	¥202,200	
IRT3460040	ガラスバイアル	¥175,600	¥203,300	
IRV2200125	広口ガラス瓶(キャップ付き)125ml Short	¥80,900	¥94,000	
IRV2200250	広口ガラス瓶(キャップ付き)250ml Short	¥44,900	¥51,900	
IRV2200500	広口ガラス瓶(キャップ付き)500ml Short	¥47,300	¥54,700	
IRV2210125	広口ガラス瓶(キャップ付き)125ml Tall	¥80,600	¥93,600	
IRV2211000	広口ガラス瓶(キャップ付き)1 000ml Tall	¥59,400	¥68,900	
IRV3200125	広口ガラス瓶(キャップ付き)125ml Short	¥93,700	¥108,300	
IRV3200250	広口ガラス瓶(キャップ付き)250ml Short	¥50,300	¥58,100	
IRV3200500	広口ガラス瓶(キャップ付き)500ml Short	¥55,400	¥64,100	
IRV3210250	広口ガラス瓶(キャップ付き)250ml Tall	¥50,400	¥58,200	
IRV3211000	広口ガラス瓶(キャップ付き)1 000ml Tall	¥67,500	¥78,100	
J1-L56003	アセリ- プリンタ用ロール紙x2/リールx2	¥6,300	¥6,600	2024/3/14 追加
J1-L56006	アセリ- プリンタ用ロール紙x2/リールx2	¥6,300	¥6,600	2024/3/14 追加
J1-L56009N	プリンタ用粘着テープ付ロール紙x6/リールx1	¥12,500	¥12,900	2024/3/14 追加
J1-L56010	アセリ- YDP03用インクカートリッジ 2本	¥2,700	¥3,000	2024/3/14 追加
J1-L56011	MMA30用ロール紙セット 5巻	¥6,300	¥6,600	2024/3/14 追加
J1-L90002CUB	動物ひょう量皿(ステンレス製) 円錐型 1	¥49,000	¥51,000	2024/3/14 追加
J1-LMMA1000	MMA30用秤量皿(ガラス製) 1000	オープン価格	¥32,000	2024/3/14 追加
J1-MD8BASE	SS製サニタリーフラスコ・フィルターアダプター-1	¥81,000	¥89,100	
JT05401	加圧カートリッジ(低圧カラム用充填剤)	¥130,200	¥144,200	
JT30031	Ammonium Bicarbonate 500g	¥34,300	¥38,200	
JT32487	Potassium Phosphate 12kg	¥678,000	¥754,300	
JT36056	Sodium carbonate 3kg	¥337,900	¥375,900	
JT40061	Starch Soluble BAKER ANALYZED 1	¥106,900	¥122,300	
JT40094	Sodium Acetate 1kg	¥169,100	¥189,900	
JT42182	Water for HPLC 1L	¥33,900	¥33,500	
JT46311	加圧カートリッジ(低圧カラム用充填剤)	¥256,400	¥280,500	
KBB1001	Pippin Pulse用10X KBB Buffer	¥27,600	¥29,400	
KT1615002942	ガラスジョイント	¥97,000	¥105,800	
KT1695002440	KONTES Distilling Adapters	¥95,700	¥110,800	
KT1850002440	ガラスジョイントジョック付	¥96,500	¥105,200	
KT1990002440	ガラスジョイントジョック付	¥42,600	¥46,500	
KT2112501420	ガラスジョイント 3方	¥123,000	¥133,700	
KT2112752440	Septum Inlet Adapter PTFE Plug	¥126,300	¥137,800	
KT2113001420	ガラスジョイントゴム管付	¥87,800	¥95,800	
KT2113002440	ガラスジョイントゴム管付	¥90,100	¥98,300	
KT2131000114	KONTES Airless-ware Flasks 10mL	¥93,600	¥101,600	
KT2131000214	KONTES Airless-ware Flasks 25mL	¥104,600	¥114,100	
KT2131000514	KONTES Airless-ware Flasks 50mL	¥107,800	¥117,600	
KT2131001014	Kontes Airless-Ware Flask	¥110,500	¥120,500	
KT2132101000	フラスコ(Airless/Straus) 1000ML	¥160,300	¥174,800	
KT2160500000	マニホールド(活栓付き)	¥500,000	¥545,200	
KT2160500005	マニホールド(活栓付き)	¥689,300	¥1,265,000	
KT2160600000	真空ジョック(ガラス製)	¥944,600	¥1,030,100	
KT2160600003	真空ジョック(ガラス製)	¥680,000	¥741,500	
KT2161000000	圧力解放弁 1	¥211,800	¥230,900	
KT2411600000	KONTES Arsine Generator 125ml	¥143,400	¥155,700	
KT2737500000	ガラスジョイント	¥83,700	¥90,900	
KT2737600000	ガラスジョイント	¥77,900	¥85,000	
KT2834000000	Condenser Dewar	¥172,200	¥187,800	
KT2871000000	Distilling head	¥243,900	¥266,000	
KT2930500323	Funnel Buchner 30ml C	¥92,400	¥100,800	
KT2940000010	丸底フラスコ(共栓)	¥27,700	¥30,300	
KT2940000015	丸底フラスコ(共栓)	¥31,500	¥34,400	
KT2940000025	丸底フラスコ(共栓)	¥26,100	¥28,400	
KT2940000035	丸底フラスコ(共栓)	¥26,200	¥28,600	
KT2940000050	丸底フラスコ(共栓)	¥28,800	¥31,500	
KT2940000100	丸底フラスコ(共栓)	¥22,300	¥24,300	
KT2940000250	丸底フラスコ(共栓)	¥34,300	¥37,400	
KT2940000500	丸底フラスコ(共栓)	¥37,300	¥40,700	
KT2940100010	Flask Round Bottom 10mL 1本	¥26,900	¥29,300	
KT2940100025	Flask Round Bottom 25mL 1本	¥26,600	¥29,000	
KT2940100050	Flask Round Bottom 50mL 1本	¥27,000	¥29,400	
KT2940100100	Flask Round Bottom 100mL 1本	¥28,600	¥31,000	
KT2942500010	梨型フラスコ	¥33,100	¥36,100	
KT2942500025	梨型フラスコ	¥33,500	¥36,500	
KT2942500050	梨型フラスコ	¥35,200	¥38,400	
KT2942500075	梨型フラスコ	¥34,800	¥38,000	
KT2942500100	梨型フラスコ	¥36,900	¥40,300	
KT2955000050	フラスコ 50ml	¥75,400	¥81,900	
KT2955500025	フラスコ 25ml	¥70,100	¥76,400	
KT3170000050	ビーカー(シヤケット型) 50ml	¥109,400	¥119,300	
KT3170000100	ビーカー(シヤケット型) 100ml	¥124,800	¥136,100	
KT3170000250	ビーカー(シヤケット型) 250ml	¥161,300	¥175,900	
KT3170000500	ビーカー(シヤケット型) 500ml	¥220,600	¥240,600	
KT4100900015	遠心チューブ 15ml w/screw cap 12本	¥163,600	¥178,400	
KT4161701020	KONTES Tank with Glass Lid	¥183,100	¥196,800	
KT4161750000	加圧カートリッジ(TLC Aluminum Rack)	¥168,600	¥181,200	
KT4161920000	加圧カートリッジ(TLC平衡化装置)	¥155,300	¥167,000	
KT4201600000	FUNNEL NMR POLYPRO 200MM PK100	¥142,000	¥154,800	
KT4201634503	活栓(HDPE製3方向ルーロック付き)	¥68,400	¥74,600	
KT4204000720	KONTES Flex-Column 8mL 10本入	¥186,900	¥203,800	
KT4204000750	KONTES Flex-Column 19mL 10本入	¥256,800	¥280,000	
KT4204010730	Flex-Column ID 0.7cm x Length 30cm	¥37,300	¥40,600	
KT4204011020	KONTES FLEX COLUMN 16mL	¥31,000	¥33,800	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
KT4204011030	Flex-Column Chromatography Column 1	¥31,000	¥33,800	
KT4204012520	KONTES FLEX COLUMN	¥55,700	¥60,800	
KT4204072003	ADAPTER F-Luer 1個	¥18,700	¥20,400	
KT4204080000	Adapter Luer	¥30,300	¥33,000	
KT4206100430	Chromatography Column w/Reservoir	¥153,700	¥165,200	
KT4206100630	Chromatography Column w/Reservoir	¥209,900	¥225,600	
KT4208230018	Kontes PTFE Tubing	¥41,200	¥44,900	
KT4208230116	Kontes PTFE Tubing	¥26,500	¥28,900	
KT4208361620	Flow Adapter 2.5cm 1個	¥360,700	¥393,300	
KT4208370000	Packing Reservoirs 150ML	¥201,500	¥219,700	
KT4208370020	Packing Reservoirs 500ML	¥396,900	¥432,800	
KT4225300025	TLC Reagent Sprayer 25mL	¥188,900	¥203,000	
KT4225300125	TLC Reagent Sprayer 125mL	¥187,900	¥201,900	
KT4577502425	Condensers Dewar Joints 24/40	¥361,300	¥394,000	
KT4599600000	CUVETTE WASHER 1個	¥252,000	¥274,800	
KT5030000121	Column Three-Chamber24/40	¥139,000	¥151,600	
KT5700010250	冷却濃縮容器(目盛付)	¥80,500	¥87,800	
KT5700010500	冷却濃縮容器(目盛付)	¥95,500	¥104,100	
KT5700110250	冷却濃縮容器(目盛付)	¥90,400	¥98,600	
KT5700250250	冷却濃縮容器(目盛付)	¥269,700	¥294,100	
KT5700250500	冷却濃縮容器(目盛付)	¥253,600	¥276,600	
KT5700350250	冷却濃縮容器(目盛付)	¥279,600	¥305,000	
KT5700350500	冷却濃縮容器(目盛付)	¥279,100	¥304,300	
KT5700370500	冷却濃縮容器(目盛付)	¥88,000	¥95,900	
KT5700501025	冷却濃縮容器(目盛付)	¥53,800	¥58,700	
KT5700511025	冷却濃縮容器(目盛付)	¥54,300	¥59,200	
KT5700512526	冷却濃縮容器(目盛付)	¥64,500	¥70,300	
KT5702000114	トラップ O-リングホーラー	¥126,500	¥137,900	
KT5702000119	トラップ O-リングホーラー	¥126,500	¥137,900	
KT5702002514	トラップ O-リングホーラー	¥126,800	¥138,300	
KT5702002519	トラップ O-リングホーラー	¥126,900	¥138,400	
KT5702002524	トラップ O-リングホーラー	¥130,800	¥142,700	
KT5702005014	トラップ O-リングホーラー	¥126,100	¥137,500	
KT5702005024	トラップ O-リングホーラー	¥131,700	¥143,600	
KT5702602524	Trap for Rotary Evap 250ml	¥129,900	¥141,600	
KT6015000524	フラスコ(フラットボトム/ヒール)500ML	¥42,600	¥46,500	
KT6015000624	フラスコ(フラットボトム/ヒール)1000ML	¥67,200	¥73,300	
KT6060201129	フラスコ 500ml Round Bottom	¥107,100	¥116,800	
KT6060201529	フラスコ 1000ml Round Bottom	¥129,400	¥141,100	
KT6230200300	Wide-Mouth Volumetric Flasks,300mL	¥87,200	¥95,100	
KT6230200400	Wide-Mouth Volumetric Flasks,400mL	¥94,800	¥103,400	
KT6572502523	Gas Washing Bottle 250ml 1	¥321,800	¥350,900	
KT7365003515	リクワリタ 化繊・シリコン・チューブ (35x150mm)	¥121,200	¥132,100	
KT7365003530	リクワリタ 化繊・シリコン・チューブ	¥138,800	¥151,300	
KT7375520000	気体採取器 イピソジヤ	¥27,700	¥30,200	
KT7495101500	KONTES Microtube only 1.5mL NT 100	¥22,600	¥24,800	
KT7495101590	ホシノナゲル 用マイクロチューブ 1.5ml	¥41,800	¥45,100	
KT7880000021	攪拌機アセンブリ	¥32,000	¥34,600	
KT7880200021	攪拌機(Arrow Junior)	¥114,900	¥124,100	
KT8264000004	バルブ(O-ringが追加) HI-VAC付)	¥66,900	¥72,900	
KT8265000004	バルブ(真空マニフォルド)	¥58,000	¥63,200	
KT8265000008	バルブ(真空マニフォルド)	¥101,600	¥110,800	
KT8265000012	バルブ(真空マニフォルド)	¥119,700	¥130,500	
KT8265010004	VALVE PLUG SIZE4	¥20,100	¥21,900	
KT8265100004	Right Angle Valves 1個	¥58,800	¥64,100	
KT8266000004	VALVE HI VAC PTFE SZ4 1個	¥31,700	¥34,500	
KT8266010008	Valve Plug Size8 1個	¥27,500	¥30,000	
KT8266100004	VALVE HI VAC PTFE Size4 1個	¥58,400	¥63,700	
KT8266200004	HI-VAC Low Holdup Valves 1個	¥39,400	¥43,000	
KT8501000008	KONTES Solid Pennyheadストッパー - サイズ8	¥10,200	¥11,200	
KT8518002440	Stopper Glass 24/40	¥28,600	¥31,200	
KT8817000000	Flask Moisture 2.2ML	¥79,600	¥86,800	
KT8823000010	FRENMAYER FLASK 10ML PK12	¥108,600	¥118,500	
KT8823600010	REACTION FLASK 10ML PK12	¥188,100	¥204,300	
KT8833500009	KONTES ティッシュ - サンプル スリッパ 250	¥25,800	¥28,100	
KT8850000002	ホシノナゲル	¥82,900	¥90,400	
KT8850000007	ホシノナゲル - 7mL	¥93,300	¥101,800	
KT8850000015	ホシノナゲル	¥102,700	¥111,900	
KT8850000040	ホシノナゲル	¥140,800	¥153,500	
KT8853000002	ホシノナゲル - 2mL 1	¥140,500	¥153,200	
KT8853000007	ホシノナゲル - 7ml	¥124,500	¥135,800	
KT8853000015	ホシノナゲル - 15ml	¥139,700	¥152,300	
KT8853000040	ホシノナゲル -(Dounce Tissue Grinders)	¥160,700	¥175,300	
KT8853000100	ホシノナゲル	¥347,500	¥378,900	
KT8853010002	ホシノナゲル - 2ml	¥42,800	¥46,700	
KT8853010007	ホシノナゲル - 7ml	¥51,900	¥56,600	
KT8853010015	ホシノナゲル - 15ml	¥54,000	¥58,900	
KT8853010040	Tissue Grind Pestle Large 40mL 1	¥83,200	¥90,700	
KT8853010100	ホシノナゲル - 100ml	¥114,400	¥124,800	
KT8853020002	Tissue Grind Pestle	¥45,100	¥49,200	
KT8853020007	ホシノナゲル - 7ml	¥48,800	¥53,200	
KT8853020015	ホシノナゲル - 15ml	¥56,500	¥61,700	
KT8853020040	Tissue Grind Pestle Small 40mL 1	¥81,600	¥88,900	
KT8853020100	ホシノナゲル - 100ml	¥114,500	¥124,800	
KT8853030002	Tissue Grind Tube 2mL	¥62,300	¥67,900	
KT8853030007	ホシノナゲル - 7ml	¥56,400	¥61,500	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
KT8853030015	Tissue Grind Tube 15mL	¥60,700	¥66,200	
KT8853030040	Tissue Grind Tube SZ 40mL	¥69,600	¥75,900	
KT8853030100	Tissue Grind Tube 100mL 1	¥129,900	¥141,700	
KT8854500021	KONTES DUALL 滅菌ナイフ - サイズ 21 3mL	¥125,400	¥136,800	
KT8854500024	KONTES DUALL 滅菌ナイフ - サイズ 24 30mL	¥156,000	¥170,100	
KT8854500025	滅菌ナイフ - Ä ユスル 50ml	¥198,000	¥215,900	
KT8854900021	Tissue Grinder 5mL 1個	¥134,400	¥146,600	
KT8855000019	滅菌ナイフ -	¥107,500	¥117,200	
KT8855000021	滅菌ナイフ -	¥85,800	¥93,500	
KT8855000022	滅菌ナイフ -	¥95,300	¥104,000	
KT8855000023	滅菌ナイフ -	¥101,800	¥111,000	
KT8855000024	滅菌ナイフ -	¥105,500	¥115,100	
KT8855100020	TISSUE GRIND COMP SZ20 1個	¥135,700	¥147,900	
KT8855100023	TISSUE GRIND COMP SZ23 1個	¥161,100	¥175,600	
KT8855110020	Tissue Grind Pestle Size20 1本	¥53,700	¥58,600	
KT8855110023	Tissue Grind Pestle Size23 1本	¥88,500	¥96,500	
KT8855120020	Tissue Grind Tube SZ20	¥83,400	¥90,900	
KT8855200100	TISSUE GRIND COMP 100ML 1個	¥174,400	¥190,200	
KT8855210100	Tissue Grind Pestle 100mL	¥107,200	¥116,900	
KT8857520018	TISSUE GRIND TUBE Size18 1本	¥55,600	¥60,700	
KT8860000018	KONTES Tissue Grinders サイズ 18 0.5mL	¥121,400	¥132,400	
KT8860000019	KONTES Tissue Grinders サイズ 19 1mL	¥100,400	¥109,500	
KT8860000020	KONTES Tissue Grinders サイズ 20 3mL	¥106,900	¥116,600	
KT8860000021	KONTES Tissue Grinders サイズ 21 5mL	¥104,000	¥113,400	
KT8860000022	KONTES Tissue Grinders サイズ 22 8mL	¥100,200	¥109,300	
KT8860000023	KONTES Tissue Grinders サイズ 23 17mL	¥116,800	¥127,400	
KT8860000024	KONTES Tissue Grinders サイズ 24 45mL	¥130,800	¥142,600	
KT8968200000	Vial Hydrolysis	¥250,800	¥273,500	
KT8968608910	Tubes Vacuum Hydrolysis 8mm	¥241,000	¥262,800	
KT9537530000	加圧クワイ (HPLC用ろ過器)	¥91,800	¥100,100	
KT9537530090	ALUMINUM CLAMP 90MM	¥137,600	¥150,100	
KT9537550000	フィルター-Glass Fannel/support	¥190,600	¥207,800	
KT9537600000	フィルター(フィルターホルダ)	¥65,100	¥71,000	
KT9537630801	STOPPER	¥16,600	¥18,100	
KT9537800000	フィルター-PTFE Faced Funnnel/Support	¥269,300	¥293,600	
KT9538010001	GLASS BASE ONLY SS 47mm 1	¥83,100	¥90,600	
KT9538050000	フィルター(フィルターホルダ)	¥226,300	¥246,800	
KT9538100000	SS SUPPORT SCREEN 47mm 1	¥190,400	¥207,600	
KT9538110000	TFE GASKET 47mm 25入	¥74,400	¥81,100	
KT9538201047	加圧クワイ (HPLC用ろ過器)	¥474,000	¥516,900	
KT9538204047	加圧クワイ (HPLC用ろ過器)	¥555,000	¥605,200	
KT9538250000	フィルター(フィルターホルダ)	¥444,100	¥484,300	
KT9538250090	90mm Filtration Apparatus	¥944,500	¥1,733,000	
KT9538260000	加圧クワイ (HPLC用ろ過器)	¥245,700	¥267,900	
KT9538280000	加圧クワイ (HPLC用ろ過器)	¥239,600	¥261,300	
KT9538290000	加圧クワイ (HPLC用ろ過器)	¥277,300	¥302,400	
KT9538300000	フィルター(フィルターホルダ)	¥40,800	¥44,400	
KT9539012002	Kontes Ultra-Ware 加工底皿 -B- -2L	¥263,400	¥287,300	
KT9539015002	Kontes Ultra-Ware 加工底皿 -B- -5L	¥575,800	¥628,000	
KT9539060047	加圧クワイ (HPLC用ろ過器)	¥145,000	¥158,200	
KT9539070000	ADAPTER SOLVENT BOTTLE 1	¥78,400	¥85,500	
KT9539090000	HighTemperature ScrewThread Cap 1	¥13,200	¥14,400	
KT9539135000	CAP HPLC ECONOMY 3 HOLE 1	¥46,700	¥50,900	
KT9541000623	Funnel Buchner 60ml C	¥144,800	¥157,900	
KT9711000047	SPE Glassware Set 47mm 1	¥362,700	¥395,500	
KT9711000090	SPE Glassware Set 90 mm 1	¥690,600	¥746,600	
KT9711010347	SPE FUNNEL 47mm 300ml 1個	¥95,100	¥103,800	
KT9711020090	SPE SUPPORT BASE 90MM	¥318,500	¥344,300	
KT9711034290	SPE SCREEN KEL-F 90MM	¥525,700	¥521,400	
LH-712146	Midi Plus用ノズル(黒)	¥8,000	¥8,500	
LH-725630	6台掛け回転式クワイ	¥24,500	¥25,500	
LH-729010	Tacta 手動ピペット 1ch 0.1-3ul	¥44,000	¥48,000	
LH-729030	Tacta 手動ピペット 1ch 2-20ul	¥44,000	¥48,000	
LH-729060	Tacta 手動ピペット 1ch 20-200ul	¥44,000	¥48,000	
LH-729070	Tacta 手動ピペット 1ch 100-1000ul	¥44,000	¥48,000	
LH-729140	Tacta 手動ピペット 8ch 30-300ul	¥112,000	¥130,000	
LH-729220	Tacta 手動ピペット 12ch 0.5-10ul	¥140,000	¥145,000	
LH-729230	Tacta 手動ピペット 12ch 5-100ul	¥140,000	¥145,000	
LH-729240	Tacta 手動ピペット 12ch 30-300ul	¥140,000	¥145,000	
LH-731610	eLINE直結充電用ICカードACアダプター	¥17,000	¥18,000	
LH-B790204	FlexiBulk 0.5-200ul 960本	¥4,800	¥5,000	
LH-B790354	FlexiBulk 350ul 960本	¥4,750	¥4,900	
LH-B791004	BIOHIT Tip 1000ul バルブ 480本	¥2,500	¥2,600	
LH-B791024	FlexiBulk 10-1000ul 480本	¥3,000	¥3,100	
LH-B791204	BIOHIT Tip 1200ul バルブ 480本	¥3,200	¥3,300	
LI1000-103	2Dフューン 1000ul キヤップなし 96x10ウェル	¥142,200	¥194,100	
M1	Miniシリス アンプル 0.2mlx40ウェル	¥13,100	¥13,600	
M2	Miniシリス アンプル 0.5mlx24ウェル	¥13,100	¥13,600	
M3	Miniシリス アンプル 1.5mlx15ウェル	¥13,100	¥13,600	
M4	Miniシリス アンプル 2.0mlx15ウェル	¥13,100	¥13,600	
M5	Miniシリス アンプル 1.8mlx8ウェル	¥13,100	¥13,600	
M6	Miniシリス アンプル 15mlx6ウェル	¥13,100	¥13,600	
M611A46	Blood DNA Auto Plate 96回用	¥69,900	¥71,900	
M611S46	Blood DNA Auto Tube 96回用	¥69,900	¥71,900	
M612A46	Tissue DNA Auto Plate 96回用	¥49,800	¥51,300	
M612S46	Tissue DNA Auto Tube 96回用	¥49,800	¥51,300	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
M613A46	Plant DNA Auto Plate 96回用	¥46,600	¥47,900	
M613A46-SE	Plant DNA AutoPlate 96回用	¥50,500	¥51,900	
M613S46	Plant DNA Auto Tube 96回用	¥46,600	¥47,900	
M613S46-SE	Plant DNA AutoTube 96回用	¥50,500	¥51,900	
M61CA46	OptiPure cfDNA Auto Plate 96回用	¥155,400	¥159,700	
M61CS46	OptiPure cfDNA Auto Tube 96回用	¥155,400	¥159,700	
M61EA46	OptiPure BloodDNA AutoPlate 96回用	¥50,500	¥51,900	
M61ES46	OptiPure Blood DNA AutoTube 96回用	¥50,500	¥51,900	
M61FA46	TANBead Fungi DNA AutoPlate 96回用	¥46,600	¥47,900	
M61FS46	TANBead Fungi DNA AutoTube 96回用	¥46,600	¥47,900	
M61GA46	Gram Bacteria DNA AutoPlate 96回用	¥50,500	¥51,900	
M61GS46	Gram Bacteria DNA Auto Tube 96回用	¥50,500	¥51,900	
M61HA46	TANBead HPV DNA Auto Plate 96回用	¥62,200	¥63,900	
M61HS46	TANBead HPV DNA Auto Tube 96回用	¥62,200	¥63,900	
M61PA46	OptiPure FFPE DNA AutoPlate 96回用	¥73,800	¥75,900	
M61PS46	OptiPure FFPE DNA AutoTube 96回用	¥73,800	¥75,900	
M622A46	OptiPure Leukocyte RNA Auto Plate	¥111,900	¥115,000	
M622S46	OptiPure Leukocyte RNA Auto Tube	¥111,900	¥115,000	
M635A46	Viral Auto Plate 96回用	¥50,500	¥51,900	
M635S46	Viral Auto Tube 96回用	¥50,500	¥51,900	
M665A46	OptiPure Viral Auto Plate 96回用	¥58,300	¥59,900	
M665S46	OptiPure Viral Auto Tube 96回用	¥58,300	¥59,900	
M685A46	Virapid Virus Auto Plate 96回用	¥46,600	¥47,900	
M685S46	Virapid Virus Auto Tube 96回用	¥46,600	¥47,900	
M6EMA46	Microbiome DNA Auto Plate 96回用	¥77,700	¥79,900	
M6EMS46	Microbiome DNA Auto Tube 96回用	¥77,700	¥79,900	
M6GMA46	Food and Feed DNA AutoPlate 96回用	¥62,200	¥63,900	
M6GMS46	Food and Feed DNA AutoTube 96回用	¥62,200	¥63,900	
M6K2A46	Tissue RNA Auto Plate 96回用	¥62,200	¥63,900	
M6K2S46	Tissue RNA Auto Tube 96回用	¥62,200	¥63,900	
M6K3A46	Plant RNA Auto Plate 96回用	¥62,200	¥63,900	
M6K3S46	Plant RNA Auto Tube 96回用	¥62,200	¥63,900	
M6PEA46	アズミト抽出キット/プレート/一時販売停止	¥58,300	¥59,900	
M6PES46	アズミト抽出キット/チューブ/一時販売停止	¥58,300	¥59,900	
M6T2A46	Tissue Total DNA Auto Plate 96回用	¥50,500	¥51,900	
M6T2S46	Tissue Total DNA Auto Tube 96回用	¥50,500	¥51,900	
M7	Miniシリウス アンプル 50 mLx2ケル	¥13,100	¥13,600	
M8	Miniシリウス アンプル 5.0mLx12ケル	¥13,100	¥13,600	
M1500-101	2Dチューブ 500ul キャップなし 96x10ラック	¥142,200	¥194,100	
M1500-103	2Dチューブ 500ul キャップなし 96x10ラック	¥88,000	¥99,700	
M1500-104	2Dチューブ 500ul キャップ白バルク ST96x10	¥196,400	¥236,300	
M1500-NCNSP	2Dtube 450ul アンプル96x10R cap無	¥108,000	¥151,900	
M1500-RC-PL	2D空ラック 500ul アンプル 10	¥15,800	¥30,700	
MILENIA01	Milenia HybriDetect 1 ④	¥269,600	¥293,200	
Mini-7KS	Mini-7KS ミニ遠心機本体	¥47,000	¥48,900	
MiniC-100	■ MiniC-100 Coolingミニ冷蔵庫本体	¥77,300	¥80,400	
MP-2500	■ MP-2500 マイクロプレート遠心機本体	¥77,700	¥80,800	
NE-AG09-0.5	FG MGA7カラム-スタックレット(TBE) 100	¥39,000	¥40,000	
NE-AG10	FG MGA7カラム-スタックレット(TAE) 75	¥40,800	¥42,000	
NE-AG11	FG MGA7カラム-スタックレット(バルク)100	¥31,600	¥33,000	
NE-AG12	FG MGX7カラム-スタックレット(バルク)100	¥36,000	¥37,000	
NE-AG13	FG MGX7カラム-スタックレット(TAE) 100	¥36,000	¥37,000	
NE-AS01	NG002用レンズ	¥1,600	¥1,700	
NE-AS02	NG002用チューブ 1.5/2.0ml	¥3,200	¥3,500	
NE-AS03	NG002用チューブ 8連用	¥3,200	¥3,500	
NE-AS04	NG002用チューブ ストックケース	¥3,200	¥3,500	
NE-AS05	アタプター0.5-0.6ml tubes 1個	¥1,600	¥1,700	
NE-AS06	アタプター0.2ml tubes 1個	¥1,600	¥1,700	
NE-FG-12	■ FastGene Blue/White LED Tab 1	¥141,500	¥153,600	
NE-FG-QS1	FastGene Q-Stain 1L ④	¥32,500	¥34,000	
NE-GPG	ポンジフィルター-ジョーグ 1	¥11,800	¥12,000	
NE-LS05	DNAreleasy Advance 10preps②	¥5,000	¥5,500	
NE-LS06	DNAreleasy Advance 50preps②	¥15,400	¥16,500	
NE-LS33	FastGene BAC free HS TAQ 500U ②	¥39,000	¥41,000	
NE-LS33S	FastGene BAC free HS TAQ 25U ②	¥2,700	¥3,000	
NE-LS36	FastGene Hifi 2x HS Master Mix②	¥35,100	¥36,000	
NE-LS52	FG Scriptase Basic 100rxn ②	¥25,900	¥27,000	
NE-LS52S	FG ScriptaseBasic サンプルキット 10rxn②	¥2,600	¥3,000	
NE-LS53	FG Scriptase II 100rxn ②	¥52,200	¥54,500	
NE-LS53S	FG Scriptase II サンプルキット 10rxn ②	¥5,000	¥6,000	
NE-LS62	FG ScriptaseBasic cDNA 100rxn ②	¥52,200	¥54,500	
NE-LS62S	FG ScriptaseBasic cDNAサンプルキット10rxn②	¥5,000	¥6,000	
NE-LS63	FG Scriptase II cDNA Kit 100rxn②	¥76,800	¥81,000	
NE-LS63S	FG ScriptaseII cDNAサンプルキット10rxn②	¥7,300	¥8,500	
NE-LS64	FG ScriptaseII R-Mix(5X) 100rxn ②	¥57,000	¥60,000	
NE-LS64S	FG ScriptaseII R-Mix(5X) 10rxn ②	¥5,800	¥7,000	
NE-LS65	FG ScriptaseII R-Mix(5X) 100rxn ②	¥57,000	¥60,000	
NE-LS65S	FG ScriptaseII R-Mix(5X) 10rxn ②	¥5,800	¥7,000	
NE-MG03	ミドリグリーン アンプル 50ul ④	¥2,200	¥3,000	
NE-MG04	ミドリグリーン アンプル 1mL ④	¥26,500	¥27,000	
NE-MG05	ミドリグリーン グレイト 50ul ④	¥2,200	¥3,000	
NE-MG06	ミドリグリーン グレイト 1mL ④	¥26,500	¥27,000	
NE-MG09	ミドリグリーン イクストラ 50ul ④	¥2,200	¥3,000	
NE-MG10	ミドリグリーン イクストラ 1ml ④	¥30,000	¥31,000	
NE-MWD100P	100bp DNA Ladder PLUS, 500ul④	¥8,300	¥9,000	
NE-MWD1P	FG 1kb DNA ladder PLUS ④	¥8,300	¥9,000	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
NE-MWD50	DNA ladder 50bp ④	¥10,100	¥11,000	
NE-MWP03	BLUEStar [®] ブルートリート [®] マーカー 500ul ^②	¥26,300	¥28,400	
NE-MWP03-8	BLUEStar [®] ブルートリート [®] マーカー 4000ul ^②	¥180,000	¥195,000	
NE-MWP04	BlueStar PLUS 500ul ②	¥28,500	¥31,000	
NE-MWP06	BlueEasy [®] ブルートリート [®] マーカー 500ul ②	¥25,000	¥27,000	
NE-MWP07	FG ブルートリート [®] マーカー 500ul ②	¥17,500	¥19,000	
NE-NG002B	FastGene Mini Centrifuge Blue	¥40,200	¥43,600	
NE-NG040	■ FastGene Plate Centrifuge	¥146,800	¥165,000	
NE-PG-MOPS10	FG MOPS [®] バッファ [®] 10	¥14,500	¥24,000	
NE-PGS012	FG PAGE [®] ゲル 10x8 12% 12well 10 ④	¥31,100	¥40,000	
NE-PGS012S	サッポ [®] PAGE [®] ゲル 10x8 12% 12well 1 ④	¥3,700	¥6,000	
NE-PGS412	FG PAGE [®] ゲル 10x8 4-12% 12well 10 ④	¥31,100	¥40,000	
NE-PGS412S	サッポ [®] PAGE [®] ゲル 10x8 4-12% 12well 1 ④	¥3,700	¥6,000	
NE-PGS420	FG PAGE [®] ゲル 10x8 4-20% 12well 10 ④	¥31,100	¥40,000	
NE-PGS420S	サッポ [®] PAGE [®] ゲル 10x8 4-20% 12well 1 ④	¥3,700	¥6,000	
NE-PGS816	FG PAGE [®] ゲル 10x8 8-16% 12well 10 ④	¥31,100	¥40,000	
NE-PGS816S	サッポ [®] PAGE [®] ゲル 10x8 8-16% 12well 1 ④	¥3,700	¥6,000	
NE-VX2	Mini Vortexer ONIGIRI	¥41,000	¥45,800	
PAC20KB	BP Gel [®] 0.75% DF PacBio-20kb 10	¥183,800	¥195,900	
PAC30KB	BP Gel [®] 0.75% DF PacBio-40kb 10	¥183,800	¥195,900	
PC96-BL	ブッシュキップ [®] マット 青 96x10枚	¥29,800	¥33,800	
PC96-COL-NS	ブッシュキップ [®] マット アウト 96x10枚	¥16,400	¥19,200	
PC96-XT-BL	ブッシュキップ [®] マット 青 イクスター [®] 96x10枚	¥29,800	¥33,800	
PC96-XT-YE	ブッシュキップ [®] マット 黄 イクスター [®] 96x10枚	¥29,800	¥33,800	
PC96-YE	ブッシュキップ [®] マット 黄 96x10枚	¥29,800	¥33,800	
PCF0001	Optical Calibration fixture Pippin	¥12,900	¥13,700	
PGB1000	Sage Pulse Gel Box	¥216,800	¥230,000	
PH2Y1600	テックス製 マルチプレックススクリーン(クリンルーム) 5.5	¥453,600	¥403,600	
PH2Y1601	テックス製 マルチプレックススクリーン(クリンルーム) 6	¥421,000	¥374,700	
PH2Y1602	テックス製 マルチプレックススクリーン(クリンルーム) 6.5	¥442,200	¥393,500	
PH2Y1603	テックス製 マルチプレックススクリーン(クリンルーム) 7	¥442,200	¥393,500	
PH2Y1604	テックス製 マルチプレックススクリーン(クリンルーム) 7.5	¥442,200	¥393,500	
PH2Y1605	テックス製 マルチプレックススクリーン(クリンルーム) 8	¥442,200	¥393,500	
PH2Y1606	テックス製 マルチプレックススクリーン(クリンルーム) 8.5	¥442,200	¥393,500	
PH2Y1607	テックス製 マルチプレックススクリーン(クリンルーム) 9	¥494,500	¥440,000	
PH2Y1624	テックス製 CR100(クリンルーム用/ST) サイズ 6	¥425,400	¥349,000	
PH2Y1625	テックス製 CR100(クリンルーム用/ST) 6.5	¥493,400	¥404,800	
PH2Y1626	テックス製 CR100(クリンルーム用/ST) 7	¥490,600	¥402,500	
PH2Y1627	テックス製 CR100(クリンルーム用/ST) 7.5	¥493,400	¥404,800	
PH2Y1628	テックス製 CR100(クリンルーム用/ST) 8	¥482,400	¥395,800	
PH2Y1629	テックス製 CR100(クリンルーム用/ST) 8.5	¥493,400	¥404,800	
PH2Y1630	テックス製 CR100(クリンルーム用/ST) 9	¥493,400	¥404,800	
PH2Y1640	テックス製 CR100(クリンルーム用/ST) 5.5	¥732,400	¥600,900	
PH2Y1730	テックス製 CP100BT(クリンルーム用) S	¥516,800	¥488,100	
PH2Y1731	テックス製 CP100BT(クリンルーム用) M	¥516,800	¥553,200	
PH2Y1732	テックス製 CP100BT(クリンルーム用) L	¥516,800	¥553,200	
PH2Y1733	テックス製 CP100BT(クリンルーム用) XL	¥516,800	¥553,200	
PH2Y1850	ニトリ製 CP100(クリンルーム用) S	¥601,100	¥575,900	
PH2Y1851	ニトリ製 CP100(クリンルーム用) M	¥601,100	¥575,900	
PH2Y1852	ニトリ製 CP100(クリンルーム用) L	¥601,100	¥575,900	
PH2Y1853	ニトリ製 CP100(クリンルーム用) XL	¥601,100	¥575,900	
PIPO001	■ Pippin Prep 本体	¥2,500,000	¥2,650,000	
PPI0200	■ Pippin Pulse 本体	¥970,000	¥1,020,000	
PRACTUM124-1S	ニ-Practum 天びん 120g	¥249,000	¥299,000	2024/3/14 追加
PRACTUM412-1S	ニ-Practum 天びん 410g	¥124,000	¥149,000	2024/3/14 追加
PRACTUM51011S	ニ-Practum 天びん 5100g	¥149,000	¥179,000	2024/3/14 追加
PVM020C-099	フルホトランスメンブレン0.2u 26cmx3.3m ロール	¥50,000	¥52,200	2024/4/1追加
Q712-02	TaqProUniv SYBRqPCRマスタミックス500回 ^②	¥26,100	¥34,800	
QUINTIX1102-1	電子天秤 Quintix 0.01g 1100g	¥190,000	¥240,000	2024/3/14 追加
QUINTIX124-1S	電子天秤 Quintix 0.1mg 120g	¥314,000	¥397,000	2024/3/14 追加
QUINTIX224-1S	電子天秤 Quintix 0.1mg/220g	¥346,000	¥437,000	2024/3/14 追加
QUINTIX313-1S	電子天秤 Quintix 0.001g 310g	¥234,000	¥296,000	2024/3/14 追加
QUINTIX412-1S	電子天秤 Quintix 0.01g 410g	¥160,000	¥202,000	2024/3/14 追加
QUINTIX51021S	電子天秤 Quintix 0.01g 100g	¥294,000	¥372,000	2024/3/14 追加
QUINTIX61011S	電子天秤 Quintix 0.1g 6100g	¥224,000	¥283,000	2024/3/14 追加
QUINTIX61021S	電子天秤 Quintix 0.01g 6100g	¥313,000	¥396,000	2024/3/14 追加
R333-00	HiScript III RT SuperMix 10回用 ^②	¥5,800	¥8,400	
R333-01	HiScript III RT SuperMix 100回用 ^②	¥58,300	¥83,500	
R501-01	RNA KeeperTissueStabilizer 100ml	¥15,700	¥22,700	
RBE2010	試薬キットM1 (BEF2010用) ④	¥30,100	¥32,100	
RBF1510	試薬キットR2 (BDF1510用) ④	¥30,100	¥32,100	
RBF2010	試薬キットV2 (BDF2010用) ④	¥40,400	¥43,100	
RBF7510	試薬キットS1 (BLF7510用) ④	¥40,400	¥43,100	
RBQ3010	試薬キットQ3 (BDQ3010/3003用) ④	¥30,100	¥32,100	
RC161816	Standard Buffer pH12.45 500ml	¥25,400	¥28,200	
RC91501	Water Reagent Grade ACS Type1	¥19,000	¥20,900	
RC915025	Water Reagent Grade ACS Type1	¥37,900	¥41,700	
RC91505	Water Reagent Grade ACS Type1	¥49,800	¥54,700	
RCK-SX96	2D空ラック スライド 10	¥40,600	¥42,200	
RED2010	試薬キットO2 (ELD2010用) ④	¥30,100	¥32,100	
RED4010	試薬キット40 (ELD4010用) ④	¥30,100	¥32,100	
RED7510	試薬キット75 (ELD7510用) ④	¥30,100	¥32,100	
REK2010	試薬キットE (CEF2010用) ④	¥30,100	¥32,100	
REX0004	試薬キット(HEX0004用) ④	¥17,600	¥18,800	
REX0012	試薬キット(HEX00012用) ④	¥44,100	¥47,000	
RHF7510	試薬キット75F (HPF7510用) ④	¥40,400	¥43,100	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
RHT1510	試薬キット15C (HTC1510用) ④	¥30,100	¥32,100	
RHT2010	試薬キット20B (HTC2010用) ④	¥30,100	¥32,100	
RIK1510	試薬キットK (CDF1510用) ④	¥30,100	¥32,100	
RIK2010	試薬キットL (CDF2010用) ④	¥30,100	¥32,100	
RIT0004	試薬キット(HIT0004用) ④	¥17,600	¥18,800	
RIT0012	試薬キット(HIT0012用) ④	¥44,100	¥47,000	
RUK7510	試薬キットU1 (BUF7510用) ④	¥40,400	¥43,100	
S6534FMGUK	ミニプレート 0.2u 28mm LL型 γ線ST 50	¥18,200	¥20,000	
S6534FMOSK	ミニプレート 0.2u 28mm ルア-ロック EOG滅菌 50	¥16,200	¥18,000	
S6534FMQ	ミニプレート 0.2u 28mm ルア-ロック 500入	¥142,000	¥152,000	
S6555FMGUK	ミニプレート 0.45u 28mm LL型 γ線ST 50	¥18,200	¥20,000	
S6555FMOSK	ミニプレート 0.45u 28mm ルア-ロック EOG滅菌 50	¥16,200	¥18,000	
S6555FMQ	ミニプレート 0.45u 28mm ルア-ロック 500入	¥142,000	¥152,000	
S6596FMOSK	ミニプレートHY 0.2u 26mm ルア-ロック S 50	¥21,800	¥23,800	
S6596FMQ	ミニプレートHY 0.2u 26mm ルア-ロック 500	¥163,000	¥178,000	
S7575FXOSK	ミニプレートSRP 0.2u 25mm ルア-スリッパ T 50	¥37,000	¥38,000	
S7594FMOSK	ミニプレート 5.0u 28mm ルア-ロック S 50入	¥16,200	¥18,000	
S7597FXOSK	ミニプレート 0.2u 28mm 酢酸セルロース ST 50	¥16,200	¥18,000	
S7597FXQ	ミニプレート 0.2u 28mm ルア-スリッパ 500	¥142,000	¥152,000	
S7598FXOSK	ミニプレート 0.45um 28mm 酢酸セルロース ST 50	¥16,200	¥18,000	
S7598FXQ	ミニプレート 0.45u 28mm ルア-スリッパ 500	¥142,000	¥152,000	
SCIT-BK-NS	スクリーンキャップ (イナ) 黒 100枚 96x10	¥74,500	¥84,400	
SCIT-WH-NS	スクリーンキャップ (イナ) 白 100枚 96x10	¥74,500	¥84,400	
SDLC01	プレートラック マルチステーション	¥360,000	¥400,000	
SDLV100040E0E	プレートラック DY 5本 V1000 12	¥76,800	¥84,480	
SECURA225D-1S	天秤 Secura 0.01/0.1mg//120/220g	¥579,000	¥695,000	2024/3/14 追加
SLA-E13200	■ SLA-E13200自動核酸抽出装置	¥2,680,000	¥2,880,000	
SP-0078	0.3ml 96well MTP 100plates	¥37,700	¥43,900	
SpinTips	Spin Tips 96本	¥3,900	¥4,000	
SR16511	ポリカーボネート製フィルムカバー 16511 1	¥12,500	¥12,800	
SR17801	ナイロン製フィルムカバー 17801 1	¥27,000	¥27,500	
SX300-003	2Dチューブ 300ul キャップなし 96x10ラック	¥108,300	¥151,900	
TABAS03KIT	TwistAmp Basic [®]	¥157,100	¥170,800	
TAEXO02KIT	TwistAmp exo [®]	¥157,100	¥170,800	
TALQBAS01	TwistAmp Liquid Basic [®]	¥110,700	¥120,400	
TALQEXO01	TwistAmp Liquid exo [®]	¥114,900	¥125,000	
TCPP10RS	FG-102用空ラック 10μl PP	¥1,300	¥1,400	
TCPP200RS	FG-302用空ラック 200μl PP	¥1,600	¥1,700	
TCPPRS	FG-301用空ラック 200μl PP	¥1,600	¥1,700	
TWTX1004	クリーニング用ワイパー (MF製/洗浄済)	¥71,400	¥75,300	
TWTX1008	AlphaLite Wipers 150枚	¥46,100	¥48,600	
TWTX1009	クリーニング用ワイパー (MF製/洗浄済)	¥52,500	¥55,400	
TWTX1010	クリーニング用ワイパー (MF製/洗浄済) 100枚	¥67,300	¥71,000	
TWTX1012	クリーニング用ワイパー MF製/洗浄済 100枚入	¥103,800	¥109,400	
TWTX1013	クリーニング用ワイパー (MF製/洗浄済)	¥54,700	¥57,700	
TWTX1025	WIPER ALPHHAWIPE POLY 60x112cm 25P	¥62,400	¥65,800	
TWTX1050	クリーニング用ワイパー (MF製/洗浄済) 100枚	¥87,300	¥92,000	
TWTX1060	クリーニング用ワイパー (MF製/洗浄済) 100枚	¥655,200	¥691,000	
TWTX1109	クリーニング用ワイパー (MF&セルロース製)	¥42,600	¥45,300	
TWTX1112	クリーニング用ワイパー (MF&セルロース製)	¥40,700	¥43,200	
TWTX1118	クリーニング用ワイパー (MF&セルロース製)	¥41,900	¥45,100	
TWTX2064	クリーニング用ワイパー (ホリステル製) 300枚入	¥31,400	¥33,100	
TWTX2069	クリーニング用ワイパー (ホリステル製) 150枚入	¥38,300	¥42,500	
TWTX2412	クリーニング用ワイパー (100%ホリステル) 30.5x30.5	¥95,800	¥101,000	
TWTX2418	クリーニング用ワイパー (100%ホリステル) 45.7x45.7	¥198,800	¥230,700	
TWTX3042	Sterile Vertex Dry HS 30cm 100/PK	¥84,700	¥89,300	
TWTX3212	クリーニング用ワイパー (MF製/ST/洗浄済)	¥525,800	¥554,400	
TWTX3215	SterileWipe AS10 Wiper 100枚入	¥149,300	¥157,600	
TWTX3340	TOC Cleaning Validation Kit 1	¥121,000	¥127,400	
TWTX3342	TOC Cleaning Validation Bulk Kit	¥434,900	¥458,700	
TWTX4004	クリーニング用ワイパー (ナイド製/洗浄済) 600入	¥134,800	¥142,900	
TWTX4009	クリーニング用ワイパー (ナイド製/洗浄済) 150入	¥115,800	¥123,400	
TWTX4012	クリーニング用ワイパー (ナイド製/洗浄済) 200入	¥204,400	¥217,600	
TWTX404	クリーニング用ワイパー (MF製) 1200入	¥68,700	¥79,800	
TWTX409	クリーニング用ワイパー (MF製) 300枚入	¥58,900	¥62,100	
TWTX49	Vertex Dry Wiper HS 23cm 150/PK	¥61,400	¥64,800	
TWTX5708	TexWrite22 Cleanroom Notebooks	¥42,400	¥47,100	
TWTX5800	クリーニング用無塵紙 (クラス10対応)	¥369,900	¥429,200	
TWTX5812	クリーニング用無塵紙 (耐熱)	¥211,100	¥239,900	
TWTX5814	クリーニング用無塵紙 (耐熱)	¥241,800	¥267,800	
TWTX5815	クリーニング用無塵紙 (耐熱)	¥247,100	¥267,200	
TWTX5816	クリーニング用無塵紙 (耐熱)	¥266,500	¥286,400	
TWTX5820	接着テープ	¥71,500	¥77,900	
TWTX5831	クリーニング用無塵紙 (耐熱)	¥280,200	¥357,400	
TWTX5832	クリーニング用無塵紙 (裏書対応)	¥156,300	¥187,700	
TWTX604	クリーニング用ワイパー (MF&セルロース製) 1200枚	¥41,200	¥43,500	
TWTX606	クリーニング用ワイパー (MF&セルロース製) 600枚	¥38,200	¥40,300	
TWTX609	クリーニング用ワイパー (MF&セルロース製) 300	¥228,100	¥185,200	
TWTX612	クリーニング用ワイパー (MF&セルロース製) 150入	¥29,100	¥30,700	
TWTX705	CleanTips Mid-Size Swab 1000本	¥28,400	¥37,100	
TWTX706A	CleanTips Swab 100本入	¥39,600	¥46,000	
TWTX707A	CleanTips Rect Swab 100本入	¥33,300	¥35,100	
TWTX709A	CleanTips Swab 100本入	¥36,600	¥38,600	
TWTX7102	ClipperMop 9x9インチ w/EXT Handle	¥152,800	¥186,500	
TWTX7103	ClipperMop 12x12インチ w/EXT Handle	¥172,000	¥194,100	
TWTX7108	AlphaMop (ヘッド、ハンドル別付) 1個	¥104,300	¥121,000	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
TWTX710A	CleanTips Swab LargeFlexible 100本	¥36,500	¥38,500	
TWTX714A	CleanTips 糸リステリアル綿棒 100本	¥45,700	¥48,200	
TWTX714K	LowTOCAAlpha PolyesterSwab 20本入	¥37,700	¥39,800	
TWTX716	CleanTips HPLC SamplingSwabLarge20	¥31,000	¥32,700	
TWTX720B	CleanTips Cotton Foam Swab 500本	¥87,200	¥92,300	
TWTX730	SWAB CLEANTIP MICROSUB 7IN PK50	¥28,300	¥37,400	
TWTX740B	CleanTips Swab Medium 500本入	¥96,100	¥101,400	
TWTX742B	CleanTips Swab 500本入	¥63,000	¥69,800	
TWTX743B	CleanTips PolyesterAlphaSwab 500入	¥73,800	¥77,800	
TWTX754B	CleanTipsPolyesterAlphaMiniSwab500	¥72,400	¥76,400	
TWTX758B	CleanTips Swab POLY/MICRO 500本入	¥71,700	¥80,100	
TWTX758E	Transplex MicroAlpha ESDswab 500入	¥76,600	¥81,200	
TWTX759B	CleanTipPolyesterAbsorbondSwab 500	¥71,500	¥83,000	
TWTX761	CleanTips 糸リステリアル綿棒 100入	¥30,900	¥32,700	
TWTX761K	CleanTipsLowTOC Small Flex Swab 20	¥36,400	¥38,400	
TWTX762	CleanTipPolyesterAbsorbondSwab 100	¥29,300	¥31,000	
USTAR01	Lateral Flow Detection Units	¥141,900	¥154,300	
VC1001	ヒパルシオン100 10 000MWCO PES 2入	¥19,600	¥21,560	
VC1002	ヒパルシオン100 10 000MWCO PES 10入	¥69,800	¥76,780	
VC1011	ヒパルシオン100 5 000MWCO PES 2入	¥19,600	¥21,560	
VC1012	ヒパルシオン100 5 000MWCO PES 10入	¥69,800	¥76,780	
VC1021	ヒパルシオン100 30 000MWCO PES 2入	¥19,600	¥21,560	
VC1022	ヒパルシオン100 30 000MWCO PES 10入	¥69,800	¥76,780	
VC1031	ヒパルシオン100 50 000MWCO PES 2入	¥19,600	¥21,560	
VC1032	ヒパルシオン100 50 000MWCO PES 10入	¥69,800	¥76,780	
VC1041	ヒパルシオン100 100 000MWCO PES 2入	¥19,600	¥21,560	
VC1042	ヒパルシオン100 100 000MWCO PES 10入	¥69,800	¥76,780	
VC1051	ヒパルシオン100 300 000MWCO PES 2入	¥19,600	¥21,560	
VC1052	ヒパルシオン100 300 000MWCO PES 10入	¥69,800	¥76,780	
VC1061	ヒパルシオン100 1 000 000MWCO PES 2入	¥19,600	¥21,560	
VC1062	ヒパルシオン100 1 000 000MWCO PES 10入	¥69,800	¥76,780	
VC1072	ヒパルシオン100 0.2u PES 10入	¥69,800	¥76,780	
VCA002	加圧システム 空気圧コントローラ 1	¥76,600	¥84,260	
VCA005	加圧遠心ハット用チューブホルダ 1	¥15,600	¥17,160	
VCA010	加圧システム マニピュレータ(加圧用) 1	¥11,000	¥12,100	
VCA011	加圧システム マニピュレータ(加圧遠心用) 1	¥8,800	¥9,680	
VCA012	加圧システム 外径4mm耐圧チューブ(3m) 1	¥2,800	¥3,080	
VCA014	ヒパルセル100用サトノリ製交換フィル 10	¥4,600	¥5,060	
VCA200	ヒパルシオン20 加圧ハット 1	¥24,000	¥26,400	
VCA800	ヒパルセル100 加圧ハット 交換フィル5枚付	¥41,000	¥45,100	
VF05H0	ヒパルシオン50R 10000MWCO RC 1入	¥42,800	¥45,000	
VF05H1	ヒパルシオン50R 5000MWCO RC 1入	¥42,800	¥45,000	
VF05H2	ヒパルシオン50R 30000MWCO RC 1入	¥42,800	¥45,000	
VF05H4	ヒパルシオン50R 100 000MWCO RC 1入	¥42,800	¥45,000	
VF05P0	ヒパルシオン50 10 000MWCO PES 2入	¥25,800	¥30,000	
VF05P1	ヒパルシオン50 5 000MWCO PES 2入	¥25,800	¥30,000	
VF05P2	ヒパルシオン50 30 000MWCO PES 2入	¥25,800	¥30,000	
VF05P3	ヒパルシオン50 50 000MWCO PES 2入	¥25,800	¥30,000	
VF05P4	ヒパルシオン50 100 000MWCO PES 2入	¥25,800	¥30,000	
VF05P6	ヒパルシオン50 1,000,000 MWCO PES 2入	¥25,800	¥30,000	
VF05P7	ヒパルシオン50 0.2u PES 2入	¥25,800	¥30,000	
VF05P9	ヒパルシオン50 3 000MWCO PES 2入	¥25,800	¥30,000	
VF20H0	ヒパルシオン200 10000MWCO Hydrosart 1	¥68,800	¥78,000	
VF20H2	ヒパルシオン200 30000MWCO Hydrosart 1	¥68,800	¥78,000	
VF20H4	ヒパルシオン200 100000MWCO Hydrosart 1	¥68,800	¥78,000	
VF20H9	ヒパルシオン200 2000MWCO Hydrosart 1	¥68,800	¥78,000	
VF20P0	ヒパルシオン200 10000MWCO PES 1入	¥68,800	¥78,000	
VF20P1	ヒパルシオン200 5000MWCO PES 1入	¥68,800	¥78,000	
VF20P2	ヒパルシオン200 30000MWCO PES 1入	¥68,800	¥78,000	
VF20P3	ヒパルシオン200 50000MWCO PES 1入	¥68,800	¥78,000	
VF20P4	ヒパルシオン200 100000MWCO PES 1入	¥68,800	¥78,000	
VF20P7	ヒパルシオン200 0.2um PES 1入	¥68,800	¥78,000	
VF20P9	ヒパルシオン200 3000MWCO PES 1入	¥68,800	¥78,000	
VFA004	ヒパルシオン- タンク用チューブ LG16 3m 1	¥6,800	¥8,000	
VFA005	ヒパルシオン-50 Y字リコネクター 1	¥2,800	¥3,400	
VFA006	ヒパルシオン-50 タンク用フィルトレーションカートリッジ-バ 1	¥28,800	¥32,000	
VFA009	ヒパルシオン-流量リコネクター-セット 1	¥2,200	¥2,400	
VFA016	ヒパルシオン-50 タンク 1	¥35,000	¥36,000	
VFA020	ヒパルシオン-50 圧力指示計 1	¥13,200	¥14,520	
VFA030	ヒパルシオン-50 T字リコネクター(並列用) 2	¥2,000	¥2,200	
VFA031	ヒパルシオン-50 連結用リコネクター 6	¥4,200	¥4,620	
VFA032	ヒパルシオン-50 スクリュー フィッティング 10	¥2,700	¥2,970	
VK01P042	ヒパルシオン 0.8μ/0.4ml 100	¥31,800	¥38,000	
VN01H01	ヒパルシオン500 10 000MWCO 1.5ml 25入	¥18,000	¥18,600	
VN01H02	ヒパルシオン500 10 000MWCO 1.5ml 100入	¥58,000	¥59,600	
VN01H21	ヒパルシオン500 30 000MWCO 1.5ml 25入	¥18,000	¥18,600	
VN01H22	ヒパルシオン500 30 000MWCO 1.5ml 100入	¥58,000	¥59,600	
VN01H31	ヒパルシオン500 50 000MWCO 1.5ml 25入	¥18,000	¥18,600	
VN01H32	ヒパルシオン500 50000 MWCO 1.5ml 100入	¥58,000	¥59,600	
VN01H41	ヒパルシオン500 100 000MWCO 1.5ml 25入	¥18,000	¥18,600	
VN01H42	ヒパルシオン500 100 000MWCO 1.5ml 100入	¥58,000	¥59,600	
VN01H81	ヒパルシオン500 125000MWCO 1.5ml 25入	¥18,000	¥18,600	
VN01H91	ヒパルシオン500 2 000MWCO 1.5ml 25入	¥18,000	¥18,600	
VN01H92	ヒパルシオン500 2 000MWCO 1.5ml 100入	¥58,000	¥59,600	
VN02H01	ヒパルシオン2 10 000MWCO 15ml 25入	¥28,000	¥30,800	
VN02H02	ヒパルシオン2 10 000MWCO 15ml 100入	¥100,000	¥110,000	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
VN02H21	ヒパソ2 30 000MWCO 15ml 25入	¥28,000	¥30,800	
VN02H22	ヒパソ2 30 000MWCO 15ml 100入	¥100,000	¥110,000	
VN02H31	ヒパソ2 50 000MWCO 15ml 25入	¥28,000	¥30,800	
VN02H32	ヒパソ2 50 000MWCO 15ml 100入	¥100,000	¥110,000	
VN02H41	ヒパソ2 30 000MWCO 15ml 25入	¥28,000	¥30,800	
VN02H42	ヒパソ2 30 000MWCO 15ml 100入	¥100,000	¥110,000	
VN02H81	ヒパソ2 25入	¥28,000	¥30,800	
VN02H91	ヒパソ2 2 000MWCO 15ml 25入	¥28,000	¥30,800	
VN02H92	ヒパソ2 2 000MWCO 15ml 100入	¥100,000	¥110,000	
VNCT01	ヒパソ1-7 100	¥10,600	¥11,660	
VS0101	ヒパソ500 10 000MWCO 25入	¥16,200	¥18,000	
VS0102	ヒパソ500 10 000MWCO 100入	¥51,000	¥56,800	
VS0111	ヒパソ500 5 000MWCO 25入	¥16,200	¥18,000	
VS0112	ヒパソ500 5 000MWCO 100入	¥51,000	¥56,800	
VS0121	ヒパソ500 30 000MWCO 25入	¥16,200	¥18,000	
VS0122	ヒパソ500 30 000MWCO 100入	¥51,000	¥56,800	
VS0131	ヒパソ500 50 000MWCO 25入	¥16,200	¥18,000	
VS0132	ヒパソ500 50 000MWCO 100入	¥51,000	¥56,800	
VS0141	ヒパソ500 100 000MWCO 25入	¥16,200	¥18,000	
VS0142	ヒパソ500 100 000MWCO 100入	¥51,000	¥56,800	
VS0151	ヒパソ500 300 000MWCO 25入	¥16,200	¥18,000	
VS0152	ヒパソ500 300 000MWCO 100入	¥51,000	¥56,800	
VS0161	ヒパソ500 1 000 000MWCO 25入	¥16,200	¥18,000	
VS0162	ヒパソ500 1 000 000MWCO 100入	¥51,000	¥56,800	
VS0172	ヒパソ500 0.2u 100入	¥51,000	¥56,800	
VS0191	ヒパソ500 3 000MWCO 25入	¥16,200	¥18,000	
VS0192	ヒパソ500 3 000MWCO 100入	¥51,000	¥56,800	
VS0201	ヒパソ2 10 000MWCO PES 25入	¥28,000	¥30,800	
VS0202	ヒパソ2 10 000MWCO PES 100入	¥86,000	¥94,600	
VS0211	ヒパソ2 5 000MWCO PES 25入	¥28,000	¥30,800	
VS0212	ヒパソ2 5 000MWCO PES 100入	¥86,000	¥94,600	
VS0221	ヒパソ2 30 000MWCO PES 25入	¥28,000	¥30,800	
VS0222	ヒパソ2 30 000MWCO PES 100入	¥86,000	¥94,600	
VS0231	ヒパソ2 50 000MWCO PES 25入	¥28,000	¥30,800	
VS0232	ヒパソ2 50 000MWCO PES 100入	¥86,000	¥94,600	
VS0241	ヒパソ2 100 000MWCO PES 25入	¥28,000	¥30,800	
VS0242	ヒパソ2 100 000MWCO PES 100入	¥86,000	¥94,600	
VS0252	ヒパソ2 300 000MWCO PES 100入	¥86,000	¥94,600	
VS0291	ヒパソ2 3 000MWCO PES 25入	¥28,000	¥30,800	
VS0292	ヒパソ2 3 000MWCO PES 100入	¥86,000	¥94,600	
VS02H01	ヒパソ2 10 000MWCO HD 25入	¥28,000	¥30,800	
VS02H02	ヒパソ2 10 000MWCO HD 100入	¥86,000	¥94,600	
VS02H11	ヒパソ2 5 000MWCO HD 25入	¥28,000	¥30,800	
VS02H12	ヒパソ2 5 000MWCO HD 100入	¥86,000	¥94,600	
VS02H21	ヒパソ2 30 000MWCO HD 25入	¥28,000	¥30,800	
VS02H22	ヒパソ2 30 000MWCO HD 100入	¥86,000	¥94,600	
VS02H91	ヒパソ2 2 000MWCO HD 25入	¥28,000	¥30,800	
VS02H92	ヒパソ2 2 000MWCO HD 100入	¥86,000	¥94,600	
VS02V2	ヒパソ2 10 000MWCO CTA 100入	¥86,000	¥94,600	
VS02X2	ヒパソ2 10 000MWCO CTA 100入	¥86,000	¥94,600	
VS04T01	ヒパソ4 MWCO 10k PES	¥32,800	¥34,800	
VS04T02	ヒパソ4 MWCO 10k PES	¥99,600	¥106,000	
VS04T11	ヒパソ4 MWCO 5k PES	¥32,800	¥34,800	
VS04T12	ヒパソ4 MWCO 5k PES	¥99,600	¥106,000	
VS04T21	ヒパソ4 MWCO 30k PES	¥32,800	¥34,800	
VS04T22	ヒパソ4 MWCO 30k PES	¥99,600	¥106,000	
VS04T31	ヒパソ4 MWCO 50k PES	¥32,800	¥34,800	
VS04T32	ヒパソ4 MWCO 50k PES	¥99,600	¥106,000	
VS04T41	ヒパソ4 MWCO 100k PES	¥32,800	¥34,800	
VS04T42	ヒパソ4 MWCO 100k PES	¥99,600	¥106,000	
VS04T91	ヒパソ4 MWCO 3k PES	¥32,800	¥34,800	
VS04T92	ヒパソ4 MWCO 3k PES	¥99,600	¥106,000	
VS0601	ヒパソ6 10 000MWCO PES 25入	¥31,800	¥32,800	
VS0602	ヒパソ6 10 000MWCO PES 100入	¥96,000	¥100,000	
VS0611	ヒパソ6 5 000MWCO PES 25入	¥31,800	¥32,800	
VS0612	ヒパソ6 5 000MWCO PES 100入	¥96,000	¥100,000	
VS0621	ヒパソ6 30 000MWCO PES 25入	¥31,800	¥32,800	
VS0622	ヒパソ6 30 000MWCO PES 100入	¥96,000	¥100,000	
VS0631	ヒパソ6 50 000MWCO PES 25入	¥31,800	¥32,800	
VS0632	ヒパソ6 50 000MWCO PES 100入	¥96,000	¥100,000	
VS0641	ヒパソ6 100 000MWCO PES 25入	¥31,800	¥32,800	
VS0642	ヒパソ6 100 000MWCO PES 100入	¥96,000	¥100,000	
VS0651	ヒパソ6 300 000MWCO PES 25入	¥31,800	¥32,800	
VS0652	ヒパソ6 300 000MWCO PES 100入	¥96,000	¥100,000	
VS0661	ヒパソ6 1 000 000MWCO PES 25入	¥31,800	¥32,800	
VS0662	ヒパソ6 1 000 000MWCO PES 100入	¥96,000	¥100,000	
VS0672	ヒパソ6 0.2u PES 100入	¥96,000	¥100,000	
VS0691	ヒパソ6 3 000MWCO PES 25入	¥31,800	¥32,800	
VS0692	ヒパソ6 3 000MWCO PES 100入	¥96,000	¥100,000	
VS15RH01	ヒパソ15R 10 000MWCO 12個入	¥21,400	¥22,800	
VS15RH02	ヒパソ15R 10 000MWCO 48個入	¥65,400	¥72,000	
VS15RH11	ヒパソ15R 5 000MWCO 12個入	¥21,400	¥22,800	
VS15RH12	ヒパソ15R 5 000MWCO 48個入	¥65,400	¥72,000	
VS15RH21	ヒパソ15R 30 000MWCO 12個入	¥21,400	¥22,800	
VS15RH22	ヒパソ15R 30 000MWCO 48個入	¥65,400	¥72,000	
VS15RH91	ヒパソ15R 2 000MWCO 12個入	¥21,400	¥22,800	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
VS15RH92	ヒパルゲン15R 2 000MWCO 48個入	¥65,400	¥72,000	
VS15RXETO	ヒパルゲンイントレスト MWCO20K 12個	¥25,000	¥27,000	
VS15T01	ヒパルゲン15turbo 10 000MWCO PES 12	¥21,400	¥22,800	
VS15T02	ヒパルゲン15turbo 10 000MWCO PES 48	¥65,400	¥72,000	
VS15T11	ヒパルゲン15turbo 5 000MWCO PES 12入	¥21,400	¥22,800	
VS15T12	ヒパルゲン15turbo 5 000MWCO PES 48入	¥65,400	¥72,000	
VS15T21	ヒパルゲン15turbo 30 000MWCO PES 12	¥21,400	¥22,800	
VS15T22	ヒパルゲン15turbo 30 000MWCO PES 48	¥65,400	¥72,000	
VS15T31	ヒパルゲン15turbo 50 000MWCO PES 12	¥21,400	¥22,800	
VS15T32	ヒパルゲン15turbo 50 000MWCO PES 48	¥65,400	¥72,000	
VS15T41	ヒパルゲン15turbo 100 000MWCO PES 12	¥21,400	¥22,800	
VS15T42	ヒパルゲン15turbo 100 000MWCO PES 48	¥65,400	¥72,000	
VS15T91	ヒパルゲン15turbo 3 000MWCO PES 12入	¥21,400	¥22,800	
VS15T92	ヒパルゲン15turbo 3 000MWCO PES 48入	¥65,400	¥72,000	
VS15TR01	ヒパルゲン15turbo 10,000 MWCO 12入	¥21,400	¥22,800	
VS15TR21	ヒパルゲンターホ15 RC 30,000 MWCO 12入	¥21,400	¥22,800	
VS15TR22	ヒパルゲンターホ15 RC 30,000 MWCO 48入	¥65,400	¥72,000	
VS15TR32	ヒパルゲンターホ15 RC 50,000 MWCO 48入	¥65,400	¥72,000	
VS2001	ヒパルゲン20 10 000MWCO PES 12入	¥22,400	¥23,800	
VS2002	ヒパルゲン20 10 000MWCO PES 48入	¥71,600	¥74,000	
VS2011	ヒパルゲン20 5 000MWCO PES 12入	¥22,400	¥23,800	
VS2012	ヒパルゲン20 5 000MWCO PES 48入	¥71,600	¥74,000	
VS2021	ヒパルゲン20 30 000MWCO PES 12入	¥22,400	¥23,800	
VS2022	ヒパルゲン20 30 000MWCO PES 48入	¥71,600	¥74,000	
VS2031	ヒパルゲン20 50 000MWCO PES 12入	¥22,400	¥23,800	
VS2032	ヒパルゲン20 50 000MWCO PES 48入	¥71,600	¥74,000	
VS2041	ヒパルゲン20 100 000MWCO PES 12入	¥22,400	¥23,800	
VS2042	ヒパルゲン20 100 000MWCO PES 48入	¥71,600	¥74,000	
VS2051	ヒパルゲン20 300 000MWCO PES 12入	¥22,400	¥23,800	
VS2052	ヒパルゲン20 300 000MWCO PES 48入	¥71,600	¥74,000	
VS2061	ヒパルゲン20 1 000 000MWCO PES 12	¥22,400	¥23,800	
VS2062	ヒパルゲン20 1 000 000MWCO PES 48	¥71,600	¥74,000	
VS2072	ヒパルゲン20 0.2u PES 48入	¥71,600	¥74,000	
VS2091	ヒパルゲン20 3 000MWCO PES 12入	¥22,400	¥23,800	
VS2092	ヒパルゲン20 3 000MWCO PES 48入	¥71,600	¥74,000	
VSA005	ヒパルゲン20 タイフイルトレーションキャップ 12	¥5,000	¥5,600	
W651446	2.0mL 凍乾アンプル 567/box	¥194,300	¥207,400	
W651450	11mm グレー プラストール - 144/box	¥115,100	¥122,900	
W651469	5mL Cryoアンプル 293/box	¥146,400	¥156,300	
W651506	2.0mL 凍乾アンプル(Pre-scored)567/box	¥218,700	¥233,500	
W651600	クイオアイル リーザーホックス 2ml用 15	¥51,700	¥54,100	
W651600-XL	クイオアイルXL リーザーホックス5ml用 15	¥56,900	¥59,500	
W651601	クイオアイル リーザーホックス 2ml用 15	¥51,700	¥54,100	
W651601-XL	クイオアイルXL リーザーホックス 3-5ml用 15	¥56,900	¥59,500	
W651602	クイオアイル リーザーホックス 2ml用 15	¥51,700	¥54,100	
W651602-XL	クイオアイルXL リーザーホックス 3-5ml用 15	¥56,900	¥59,500	
W651603	クイオアイル リーザーホックス 2ml用 15	¥51,700	¥54,100	
W651603-XL	クイオアイルXL リーザーホックス 3-5ml用 15	¥56,900	¥59,500	
W651604	クイオアイル リーザーホックス 2ml用 15	¥51,700	¥54,100	
W651604-XL	クイオアイルXL リーザーホックス 3-5ml用 15	¥56,900	¥59,500	
W651605	クイオアイル リーザーホックス 2ml用 15	¥51,700	¥54,100	
W651605-XL	クイオアイルXL リーザーホックス 3-5ml用 15	¥56,900	¥59,500	
W651610-B	Tissueバール用ホックス 青 15/box	¥50,000	¥52,400	
W651610-G	Tissueバール用ホックス 緑 15/box	¥50,000	¥52,400	
W651610-P	Tissueバール用ホックス ピンク 15/box	¥50,000	¥52,400	
W651610-R	Tissueバール用ホックス 赤 15/box	¥50,000	¥52,400	
W651610-W	Tissueバール用ホックス 白 15/box	¥50,000	¥52,400	
W651610-Y	Tissueバール用ホックス 黄 15/box	¥50,000	¥52,400	
W651611	FTAカード用ホックス 2/box	¥11,800	¥12,400	
W651612	CryoFile Storage Box 4	¥23,000	¥24,100	
W651702-B	KeepIT-25 リーザーホックス青 10/box	¥49,200	¥51,500	
W651702-G	KeepIT-25 リーザーホックス緑 10/box	¥49,200	¥51,500	
W651702-P	KeepIT-25 リーザーホックスピンク 10/box	¥49,200	¥51,500	
W651702-R	KeepIT-25 リーザーホックス赤 10/box	¥49,200	¥51,500	
W651702-W	KeepIT-25 リーザーホックス白 10/box	¥49,200	¥51,500	
W651702-Y	KeepIT-25 リーザーホックス黄 10/box	¥49,200	¥51,500	
W651703-B	KeepIT-81 リーザーホックス青 10/box	¥67,300	¥70,400	
W651703-G	KeepIT-81 リーザーホックス緑 10/box	¥67,300	¥70,400	
W651703-P	KeepIT-81 リーザーホックスピンク 10/box	¥67,300	¥70,400	
W651703-R	KeepIT-81 リーザーホックス赤 10/box	¥67,300	¥70,400	
W651703-W	KeepIT-81 リーザーホックス白 10/box	¥67,300	¥70,400	
W651703-Y	KeepIT-81 リーザーホックス黄 10/box	¥67,300	¥70,400	
W651704-B	KeepIT-100 リーザーホックス青 10/box	¥67,300	¥70,400	
W651704-G	KeepIT-100 リーザーホックス緑 10/box	¥67,300	¥70,400	
W651704-P	KeepIT-100 リーザーホックスピンク 10/box	¥67,300	¥70,400	
W651704-R	KeepIT-100 リーザーホックス赤 10/box	¥67,300	¥70,400	
W651704-W	KeepIT-100 リーザーホックス白 10/box	¥67,300	¥70,400	
W651704-Y	KeepIT-100 リーザーホックス黄 10/box	¥67,300	¥70,400	
W985100	5mL Cryo Tissue Vial 白 ST250	¥162,900	¥170,400	
W985810	Cryoアンプル 50本立 5/box	¥82,400	¥87,100	
W985850	2ml Cryoアンプル 丸底(W) NT 1000	¥190,000	¥198,900	
W985851	2ml Cryoアンプル 丸底 NT 1000	¥155,200	¥162,400	
W985852	2ml Cryoアンプル 自立(W) NT 1000	¥206,100	¥215,700	
W985853	2ml Cryoアンプル 自立 NT 1000	¥168,400	¥176,200	
W985860	1.2ml Cryoアンプル 丸底(W) NT ST500	¥85,800	¥89,800	
W985861	2ml Cryoアンプル 丸底(W) NT ST500	¥98,800	¥103,400	

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
W985862	1.2ml Cryoアタ- 自立(WS) NT ST500	¥91,200	¥95,400	
W985863	2ml Cryoアタ- 自立(WS) 白 ST500	¥107,500	¥112,600	
W985863-100	2ml Cryoアタ- 自立(WS) 白 ST100	¥29,600	¥31,000	
W985863-BC	2ml Cryoアタ-ブリカト自立(WS)白 500	¥290,500	¥304,000	
W985864	2ml Cryoアタ- 自立(WS) 赤 ST500	¥107,500	¥112,600	
W985864-100	2ml Cryoアタ- 自立(WS) 赤 ST100	¥29,600	¥31,000	
W985864-BC	2ml Cryoアタ-ブリカト自立(WS)赤 500	¥290,500	¥304,000	
W985865	2ml Cryoアタ- 自立(WS) ピンク ST500	¥107,500	¥112,600	
W985865-100	2ml Cryoアタ- 自立(WS) ピンク ST100	¥29,600	¥31,000	
W985865-BC	2mlCryoアタ-ブリカト自立(WS)ピンク500	¥290,500	¥304,000	
W985866	2ml Cryoアタ- 自立(WS) 黄 ST500	¥107,500	¥112,600	
W985866-100	2ml Cryoアタ- 自立(WS) 黄 ST100	¥29,600	¥31,000	
W985866-BC	2ml Cryoアタ-ブリカト自立(WS)黄 500	¥290,500	¥304,000	
W985867	2ml Cryoアタ- 自立(WS) 緑 ST500	¥107,500	¥112,600	
W985867-100	2ml Cryoアタ- 自立(WS) 緑 ST100	¥29,600	¥31,000	
W985867-BC	2ml Cryoアタ-ブリカト自立(WS)緑 500	¥290,500	¥304,000	
W985868	2ml Cryoアタ- 自立(WS) 青 ST500	¥107,500	¥112,600	
W985868-100	2ml Cryoアタ- 自立(WS) 青 ST100	¥29,600	¥31,000	
W985868-BC	2ml Cryoアタ-ブリカト自立(WS)青 500	¥290,500	¥304,000	
W985869	3ml Cryoアタ- 自立(WS) NT ST500	¥128,800	¥134,800	
W985870	4ml Cryoアタ- 自立(WS) NT ST500	¥138,500	¥145,000	
W985871	5ml Cryoアタ- 自立(WS) NT ST500	¥145,000	¥151,800	
W985872	2ml Cryoアタ- 自立(WS) NT ST500	¥107,500	¥112,600	
W985872-100	2ml Cryoアタ- 自立(WS) NT ST100	¥29,600	¥31,000	
W985874	0.5ml Cryoアタ- 自立(WS) NT ST500	¥108,300	¥113,400	
W985881	2DJ-ド イサト 500	¥141,500	¥148,000	
W985883	1.2ml Cryoアタ- 自立(WS) 白 ST500	¥89,200	¥93,400	
W985884	1.2ml Cryoアタ- 自立(WS) 赤 ST500	¥91,200	¥95,400	
W985885	1.2ml Cryoアタ- 自立(WS) ピンク ST500	¥91,200	¥95,400	
W985886	1.2ml Cryoアタ- 自立(WS) 黄 ST500	¥89,200	¥93,400	
W985887	1.2ml Cryoアタ- 自立(WS) 緑 ST500	¥89,400	¥93,500	
W985888	1.2ml Cryoアタ- 自立(WS) 青 ST500	¥89,400	¥93,500	
W985900	2ml Cryoイナ- 丸底(WS) NT 1000	¥190,000	¥198,900	
W985901	2ml Cryoイナ- 丸底 NT 1000	¥155,200	¥162,400	
W985902	2ml Cryoイナ- 自立(WS) NT 1000	¥206,100	¥215,700	
W985903	2ml Cryoイナ- 自立 NT 1000	¥168,400	¥176,200	
W985910	1.2ml Cryoイナ- 丸底(WS) NT ST500	¥85,900	¥89,900	
W985911	2ml Cryoイナ- 丸底(WS) NT ST500	¥98,200	¥102,700	
W985912	3ml Cryoイナ- 丸底(WS) NT ST500	¥126,300	¥132,200	
W985913	4ml Cryoイナ- 丸底(WS) NT ST500	¥135,900	¥142,300	
W985914	5ml Cryoイナ- 丸底(WS) NT ST500	¥142,400	¥149,000	
W985915	1.2ml Cryoイナ- 自立(WS) NT ST500	¥90,600	¥94,900	
W985916	2ml Cryoイナ- 自立(WS) 白 ST500	¥107,600	¥112,600	
W985917	2ml Cryoイナ- 自立(WS) 赤 ST500	¥107,600	¥112,600	
W985918	2ml Cryoイナ- 自立(WS) ピンク ST500	¥107,600	¥112,600	
W985919	2ml Cryoイナ- 自立(WS) 黄 ST500	¥107,600	¥112,600	
W985920	2ml Cryoイナ- 自立(WS) 緑 ST500	¥107,600	¥112,600	
W985921	2ml Cryoイナ- 自立(WS) 青 ST500	¥107,600	¥112,600	
W985922	2ml Cryoイナ- 自立(WS) NT ST500	¥107,600	¥112,600	
W985923	3ml Cryoイナ- 自立(WS) NT ST500	¥126,300	¥132,200	
W985924	4ml Cryoイナ- 自立(WS) NT ST500	¥135,900	¥142,300	
W985925	5ml Cryoイナ- 自立(WS) NT ST500	¥142,400	¥149,000	
WFT10011	FilterTips 10uL 透明 ST96本x10R	¥13,800	¥16,700	
WFT10020	FilterTips 20uL ST 96本x10パック	¥13,800	¥17,600	
WFT10100	FilterTips 100uL ST 96本x10パック	¥13,800	¥17,600	
WFT10200	FilterTips 200uL ST 96本x10パック	¥13,800	¥17,600	
WFT11250	FilterTips 1250uL ST 96本x10パック	¥17,000	¥19,500	
WH-10	サ-モスタットウォーターバ-ス 10L	¥146,800	¥152,700	
WLT20010	フィルタ-チップ° 10uL LB ST 96本x10R	¥15,200	¥19,500	
WLT20011	フィルタ-チップ° 10uLEXONカ° LB ST96本x10R	¥15,200	¥19,500	
WLT20012	フィルタ-チップ° 10uLカトカ° LB ST96本x10R	¥15,900	¥20,400	
WLT20020	フィルタ-チップ° 20uL LB ST 96本x10R	¥15,900	¥20,400	
WLT20100	フィルタ-チップ° 100uL LB ST 96本x10R	¥15,900	¥20,400	
WLT20200	フィルタ-チップ° 200uL EXONカ° LB ST96本x10R	¥15,900	¥20,400	
WLT21000	フィルタ-チップ° 1000uL LB ST96本x10R	¥17,300	¥22,300	
WLT21250	フィルタ-チップ° 1000uLEXONカ° LB ST96本x10R	¥17,300	¥22,300	
YAW31	アタ-リ- ピンセット (1mg-5g用) 115mm	¥6,400	¥7,800	2024/3/14 追加
YAW32	アタ-リ- ピンセット (1g-200g用) 160mm	¥10,600	¥12,800	2024/3/14 追加
YCC01-0041M3	プリンタ-接続ケーブル 1	¥34,000	¥37,400	
YCC01-USBM2	アタ-リ- USB変換出力ケーブル 1.8m	¥29,500	¥32,000	2024/3/14 追加
YCW0121	アタ-リ- 校正用分銅 1mg	¥11,300	¥11,600	2024/3/14 追加
YCW1121	アタ-リ- 校正用分銅 10mg	¥11,300	¥11,600	2024/3/14 追加
YCW113	アタ-リ- 校正用分銅 10mg F1	¥6,300	¥6,900	2024/3/14 追加
YCW1221	アタ-リ- 校正用分銅 20mg	¥11,300	¥11,600	2024/3/14 追加
YCW2121	アタ-リ- 校正用分銅 100mg E2	¥11,300	¥11,600	2024/3/14 追加
YCW312	アタ-リ- 校正用分銅 1g E2	¥19,000	¥19,800	2024/3/14 追加
YCW313	アタ-リ- 校正用分銅 1g F1	¥11,300	¥12,400	2024/3/14 追加
YCW413	アタ-リ- 校正用分銅 10g F1	¥12,600	¥15,000	2024/3/14 追加
YCW414	アタ-リ- 校正用分銅 10g F2	¥11,900	¥12,600	2024/3/14 追加
YCW453	アタ-リ- 校正用分銅 50g F1	¥20,000	¥20,600	2024/3/14 追加
YCW512	アタ-リ- 校正用分銅 100g E2	¥32,500	¥36,800	2024/3/14 追加
YCW513	アタ-リ- 校正用分銅 100g F1	¥23,100	¥28,200	2024/3/14 追加
YCW5228	アタ-リ- 校正用分銅 200g	¥28,100	¥31,600	2024/3/14 追加
YCW5238	アタ-リ- 校正用分銅 200g	¥18,200	¥23,600	2024/3/14 追加
YCW524	アタ-リ- 校正用分銅 F2級 200g	¥26,300	¥28,200	2024/3/14 追加
YCW5538	アタ-リ- 校正用分銅 500g	¥25,000	¥40,000	2024/3/14 追加

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。

2024年4月 日本ジェネティクス株式会社 取扱商品 価格改訂対象品一覧

型番	商品名	2024年3月31日まで 旧定価 (税抜)	2024年4月1日より 新定価 (税抜)	備考
YCW6138	アセリ- 校正用分銅 1kg	¥35,000	¥44,000	2024/3/14 追加
YCW6238	アセリ- 校正用分銅 2kg	¥60,000	¥70,400	2024/3/14 追加
YCW624	アセリ- 校正用分銅 F2級 2kg	¥77,000	¥84,200	2024/3/14 追加
YCW6528	アセリ- 校正用分銅 5kg E2	¥138,800	¥161,000	2024/3/14 追加
YCW6538	アセリ- 校正用分銅 5kg	¥96,200	¥101,400	2024/3/14 追加
YDK01	アセリ- 比重測定キット 0.1mm用	¥212,000	¥236,000	2024/3/14 追加
YDP20-OCE	アセリ- 天秤用アリンタ	¥134,000	¥157,000	2024/3/14 追加
YEPS01-PS4	ACアダプター 1	¥10,000	¥11,000	2024/3/14 追加
YJ-360T-S	Photometer用アリンタ(サ-マルアリンタミニサイズ)	¥74,400	¥87,300	
YPE01RC	アセリ- フットスイッチ (Tコネクタ付)	¥108,000	¥118,000	2024/3/14 追加
YWT09	アセリ- 天びん台	¥410,000	¥448,000	2024/3/14 追加

2024年4月1日現在の情報です。取扱状況、価格は予告なく変更する場合がございます。